



Iris
Biotech



LIGATION
TECHNOLOGIES



Version: IB7_3

Empowering Peptide Innovation

With this guiding theme in mind, Iris Biotech's mission is to support researchers by supplying

- innovative technologies,
- rare compounds,
- as well as a broad portfolio on standard consumables,

available in flexible quantities from small scale to bulk quantities. To fulfill our dedication "Empowering Peptide Innovation", we are attending various conferences, symposia, and exhibitions each year. This allows us to remain in direct contact with scientists all over the world, both from academia and industry, to exchange knowledge, and to gather new ideas to tackle your current challenges.

Guided by our dedication to provide

- competent service,
- as well as novel substances and
- latest technologies,

Iris Biotech is your trusted partner for the world of peptides, while having strong expertise in associated disciplines. Thus, our portfolio comprises reagents and tools for the synthesis and modification of peptides, e.g. amino acids, resins and solvents but also for related technologies such as Drug Delivery, Linkerology® and Life Sciences.



Amino Acids



Building Blocks



Life Sciences



Drug Delivery



Reagents



Resins



Linkerology®



Click Chemistry

Owed to the growing demand for tailor-made compounds, our portfolio is fine-tuned by our Custom Synthesis Service at Iris Biotech Laboratories. Our skilled scientists offer profound expertise in

- *de novo* route development,
- upscaling towards larger scale production,
- as well as synthesis optimization for increased efficiency.

Examples are the synthesis of rare chiral building blocks, unnatural amino acid derivatives, sophisticated orthogonal protecting groups, heterocycles, building blocks for nucleotides, PEGs and PEG-analogues as well as specific linkers for controlled drug delivery and release.


Portfolio Overview

Peptide Synthesis and Modification	Linkerology® and Drug Delivery	Life Sciences
(Protected) Amino Acids Standards such as Fmoc-D/L-AAA and Boc-D/L-AAA, Smoc amino acids for peptide synthesis in water, variety of protecting groups (e.g. Pbf, Trt, ^t Bu, Bzl, Acm, Mob, SIT, Phacm, Allocam, Mmt), unusual amino acids, fluorinated derivatives, substituted prolines, arginine analogues	Linkers for Solid Phase Peptide Synthesis Cleavable Linkers Val-Ala based, Val-Cit based, disulfide-based, Dde-helping hands	Biotinylation Reagents Carbohydrates Galactose, Glucose, Maltose, Mannose, Xylose and others
Building Blocks Amino alcohols, amino aldehydes, diamines and hydrazines, (pseudoproline) dipeptides, polyamines and spermines, fatty acid derivatives	Photo-Activatable Linkers Functionalized Linkers Clickable linkers, trifunctional linkers, linkers with maleimide function, cross-linkers, selective N-term acylation and biotinylation	Peptides Substrates & Inhibitors E.g. protein kinase inhibitors, substrates for fusion (Halo/Snap/Clip)-tagged proteins
Reagents Coupling reagents, solvents and scavengers, protecting groups	PROTACs Ligands, linkers & modules	Natural Products Dyes and Fluorescent Labels E.g. ICG, AMC, DAPI
Resins Preloaded resins (e.g. based on Trityl, TCP, TentaGel, Methoxybenzhydryl, Merrifield, PAM, Rink, Wang), scavenger resins, hydrazone resins	Fullerenes, Poly(2-oxazolines) & Dextrans Poly-Amino Acids Poly-Arg, Poly-Glu, Poly-Lys, Poly-Orn, Poly-Sar	Maillard & Amadori Reaction Products Large portfolio of derivatives useful as standards for food, pharma and cosmetics industry
	PEGylation Branched PEGylating reagents, (amino-)PEG-acids, PEG-amines & hydrazides & guanidines, reagents for Click-conjugation, Biotin-PEG-reagents, PEG-thiols, PEG-maleimides, other PEGylating reagents	Vitamins

Custom Synthesis

Your project requires a compound not listed in our portfolio?
Get in contact and inquire about our custom synthesis capabilities.

Our experienced scientists are excited to accept your synthetic challenge!
In such cases, your request undergoes the following stages:



Step-by-Step Analysis

- Customer's demands

Process Evaluation

- Detailed literature review
- Synthetic possibilities

Strategy Development

- Protocol development
- Method development and validation
- Customized synthesis

Quality Consistency

- Identity confirmation
- Purity verification

Our Service Promise

All our services are based on high standards, transparency & documentation, trust, honesty & confidentiality, as well as the required know-how.

High Standards

- Values: sustainability & responsibility
- State-of-the-art equipment & latest technologies
- High quality standards
- Qualified suppliers & regular audits

Transparency & Documentation

- Talk to our specialists – customer care
- Certificates of analysis & impurity profiling
- Analytical and process reports

Trust, Honesty & Confidentiality

- Intergenerational business valuing partnerships
- Meeting the customer's expectations
- Integrity towards our customers

Our Know-How

- One-step reactions & complex multi-step synthesis
- Scalability from mg to kg quantities
- Route scouting



Table of contents

Ligation Technologies

1. Description of the Technology	1
<hr/>	
2. Tools for the Formation of Thioesters	6
2.1. Dawson Linkers	6
2.2. Hydrazone Resins	8
<hr/>	
3. Building Blocks for Ligation	14
3.1. Cysteine Building Blocks for Ligation	14
3.2. Ligation at the Position of Glycine	29
3.3. Building Blocks and Reagents for Ligation with Seleno Amino Acids	31
3.4. Oxime Ligation using Aminooxy-Amino Acids	35
3.5. Activated Cysteine-Based Protein Ligation (ACPL)	36
<hr/>	
4. Enzyme-mediated Ligation	38
<hr/>	
5. Examples	40
5.1. Kinetically Controlled Four-Fragment Ligation for the Convergent Chemical Synthesis of Proteins	40
5.2. Three-Fragment Synthesis of a 116mer with Desulfurization	41
<hr/>	
Code of Conduct	43
Terms and Conditions of Sales	45
Index	49


1. Description of the Technology

Background

Peptides are typically synthesized via solid phase (peptide) synthesis (SP(PS)). The main advantage in comparison to classical solution phase synthesis is the fast and easy separation of the desired product from excess reagents by filtration. Typically, in SP(PS), the molecule being synthesized (e.g., a growing peptide chain) is attached to an insoluble solid support, while reagents are added to the suspension. This setup enables the removal of excess reagents and dissolved byproducts via simple disposal of the reaction solution. Filtration followed by washing of the resin with various solvents leaves the purified protected peptide on the carrier. Consequently, an excess of reagents can be employed in order to shorten reaction times and allow for a quantitative turnover of the substrate, which in turn leads to higher yields and less side products.

The Main Advantages of SPPS

1. Fast work-up through easy separation of solid support from dissolved reactants and by products by filtration, and multiple, rapid washing steps.
2. Improved reaction times, turnover and yield by use of excess amounts of reagents.
3. The syntheses can easily be automatized.
4. Toxic or hazardous materials can be handled safely while being attached to the resin.
5. Minimal physical product loss.
6. Pseudo-dilution phenomena on an individual bead can enable cyclization and avoid formation of dimers.



Protocol of solid phase synthesis of peptides:

Attachment of the first amino acid to a resin
Repetitive cycles of:

- Deprotection
- Washing
- Coupling
- Washing
- ...
- Deprotection

Final cleavage from the resin and purification.

# of Steps	Yield per Step/Final Yield			
	98,0%	99,0%	99,5%	99,9%
10	82%	90%	95%	99%
20	67%	82%	90%	98%
30	55%	74%	86%	97%
40	45%	67%	82%	96%
50	36%	61%	78%	95%
60	30%	55%	74%	94%
70	24%	49%	70%	93%
80	20%	45%	67%	92%
90	16%	40%	64%	91%
100	13%	37%	61%	90%
110	11%	33%	58%	90%
120	9%	30%	55%	89%
130	7%	27%	52%	88%
140	6%	24%	50%	87%
150	5%	22%	47%	86%
160	4%	20%	45%	85%
170	3%	18%	43%	84%
180	3%	16%	41%	84%
190	2%	15%	39%	83%
200	2%	13%	37%	82%

Fig. 1: Concept of solid phase peptide synthesis (SPPS).

Two reaction steps have to be carried out per amino acid of a peptide chain. I.e., the synthesis of a 50mer requires 100 chemical steps. Assuming, that the yield per chemical step is 98%, which in general means is rather high, there is only 13% yield of the target sequence left (*Fig. 1*). Nevertheless, within a long sequence, the probability that difficult fragments are appearing is rather high, which will lower the overall yield of the synthesis. In contrast, an extremely high average yield of 99.9% per coupling step will result in 90% global yield. As this is rather unrealistic, new chemical methodologies such as fragment condensation and native chemical ligation must be developed to facilitate access to longer peptides of 100 amino acids and more.

Native Chemical Ligation

Native chemical ligation (NCL) of unprotected peptide segments involves the reaction between a first peptide fragment α -thioester and a second peptide fragment, which carries a cysteine on the N-terminus, to yield a product with a native amide bond at the ligation site. Peptide- α -thioalkyl esters are commonly used because of their ease of preparation. These thioalkyl esters are rather unreactive, so the ligation reaction is catalyzed by *in situ* transthioesterification with thiol additives, which are either a mixture of thiophenol/benzyl mercaptan or other alkanethiols. Despite the use of this thiol additive, ligation reactions typically take 24–48 h.




Fig. 2: Mechanism of native chemical ligation.

back to content ↑

Ligation Technologies

Requirements

1. Unprotected peptide fragment as C-terminal thioester.
2. Cysteine on N-terminus of a second unprotected peptide fragment.

Limitations

1. Cysteine appears very rarely in natural sequences.
2. Multiple, consecutive ligations are hardly possible.
3. Cysteine is prone to racemization.
4. Slow kinetics, as ligations normally take 24-48 h; hence a good thiol additive, like 4-mercaptophenylacetic acid (MPAA), is necessary.

The main advantage of this ligation technology vs. fragment condensation is the usually good solubility of unprotected peptide fragments compared to fully protected peptides, which have very limited solubility. One main challenge, however, is that the presence of cysteine within a peptide sequence is a prerequisite for NCL, but at the same time cysteine is the least abundant of all proteinogenic amino acids since it occurs with an average frequency of only 1.4% within natural sequences. This becomes even more of an obstacle if multiple consecutive ligations are necessary to construct the target sequence. Hence, sophisticated variations must be implemented in the absence of cysteine.




Fig. 3: The role of MPAA as an additive in native chemical ligation.

Element	Electronegativity	Atomic Radius	Electron Configuration
O	3.4	0.60 ⁰ Å, 1.38 ⁻² ₄ Å	[He] 2s ² 2p ⁴
S	2.6	1.04 ⁰ Å, 1.84 ⁻² ₆ Å	[Ne] 3s ² 3p ⁴
Se	2.6	1.16 ⁰ Å, 1.98 ⁻² ₆ Å	[Ar] 3d ¹⁰ 4s ² 4p ⁴

Structure	Leaving Group	pK _a
	Alkyl Alcohol	16-18
	Alkyl Thiol	10-11
	Alkyl Selenol	7-8
	Phenol	9.89
	Thiophenol	6.6
	Cys S-H	8.3
	Sec Se-H	5.2
	Alkyl Carboxylic Acids	~4

The mechanistic details of NCL have been studied by several groups in order to find appropriate additives. Substituted thiophenols with pK_a > 6 were found to best combine the ability to exchange rapidly and completely with thioalkyl esters, and to then act as efficient leaving groups in the reaction of the peptide thioester with the thiol side chain of an N-terminal cysteine. 4-mercaptophenylacetic acid (MPAA), a non-malodorous, water-soluble thiol, proved to be a highly effective and practical additive. MPAA reacts orders of magnitude faster than other thiols in model studies of NCL and in the synthesis of small proteins, as demonstrated in the synthesis of turkey ovomucoid third domain (OMTKY3). MPAA should find broad use in native chemical ligation and in the total synthesis of proteins.

The need for a cysteine residue on the N-terminus of one of the reacting peptide fragments has motivated the development of β-, γ-, and δ-thiolated variants of other amino acids, as well as

thiol-containing auxiliaries that can be used as Cys surrogates in ligation chemistry. After ligation reactions at these residues, the thiol auxiliary is desulfurized (usually by means of radical-based protocols) to yield native polypeptide products. However, this transformation is not chemoselective in the presence of other unprotected cysteine residues that might be found elsewhere in the sequence. This limitation of desulfurization chemistry has led chemists to expand the native chemical ligation toolbox to include the 21st amino acid, selenocysteine (Sec), as well as various non-natural selenoamino acids (to date specifically based on Pro and Phe). The key advantage of carrying out ligation chemistry with selenoamino acids rather than thioamino acids is that chemoselective deselenization can be performed under mild conditions (typically with a phosphine reductant and a hydrogen atom source) that do not affect unprotected Cys residues. Payne *et al.* and Dery *et al.* have demonstrated that ligation products can be subjected to oxidative deselenization to afford Ser in place of Sec at the ligation junction.




Fig. 4: Differences in properties and capabilities using either sulfur or selenium in native chemical ligation.

Ligation Technologies





References:

- [13] Protein synthesis by chemical ligation of unprotected peptides in aqueous solution; T. W. Muir, P. E. Dawson, S. B. H. Kent; **Solid-Phase Peptide Synthesis** 1997; **289**: Academic Press, 266-298.
- Accelerated Protein Synthesis via One-Pot Ligation-Deselenization Chemistry; N. J. Mitchell, J. Sayers, S. S. Kulkarni, D. Clayton, A. M. Goldys, J. Ripoll-Rozada, P. J. Barbosa Pereira, B. Chan, L. Radom, R. J. Payne; **Chem** 2017; **2**: 703-715. <https://doi.org/10.1016/j.chempr.2017.04.003>
- Chemical synthesis of proteins; B. L. Nilsson, M. B. Soellner, R. T. Raines; **Annu Rev Biophys Biomol Struct** 2005; **34**: 91-118. <https://doi.org/10.1146/annurev.biophys.34.040204.144700>
- Construction of Challenging Proline-Proline Junctions via Diselenide-Selenoester Ligation Chemistry; J. Sayers, P. M. T. Karpati, N. J. Mitchell, A. M. Goldys, S. M. Kwong, N. Firth, B. Chan, R. J. Payne; **J Am Chem Soc** 2018; **140**: 13327-13334. <https://doi.org/10.1021/jacs.8b07877>
- Expressed protein ligation: a general method for protein engineering; T. W. Muir, D. Sondhi, P. A. Cole; **Proc Natl Acad Sci U S A** 1998; **95**: 6705-6710. <https://doi.org/10.1073/pnas.95.12.6705>
- Free-radical-based, specific desulfurization of cysteine: a powerful advance in the synthesis of polypeptides and glycopolypeptides; Q. Wan, S. J. Danishefsky; **Angew Chem Int Ed Engl** 2007; **46**: 9248-9252. <https://doi.org/10.1002/anie.200704195>
- Insights into the mechanism and catalysis of the native chemical ligation reaction; E. C. Johnson, S. B. Kent; **J Am Chem Soc** 2006; **128**: 6640-6646. <https://doi.org/10.1021/ja058344i>
- Native Chemical Ligation of Peptides and Proteins; P. A. Cistrome, M. J. Bird, D. T. Flood, A. P. Silvestri, J. C. J. Hintzen, D. A. Thompson, P. E. Dawson; **Curr Protoc Chem Biol** 2019; **11**: e61. <https://doi.org/10.1002/cpch.61>
- Oxidative Deselenization of Selenocysteine: Applications for Programmed Ligation at Serine; L. R. Malins, N. J. Mitchell, S. McGowan, R. J. Payne; **Angew Chem Int Ed Engl** 2015; **54**: 12716-12721. <https://doi.org/10.1002/anie.201504639>
- Protein synthesis by native chemical ligation: expanded scope by using straightforward methodology; T. M. Hackeng, J. H. Griffin, P. E. Dawson; **Proc Natl Acad Sci U S A** 1999; **96**: 10068-10073. <https://doi.org/10.1073/pnas.96.18.10068>
- Rapid additive-free selenocystine-selenoester peptide ligation; N. J. Mitchell, L. R. Malins, X. Liu, R. E. Thompson, B. Chan, L. Radom, R. J. Payne; **J Am Chem Soc** 2015; **137**: 14011-14014. <https://doi.org/10.1021/jacs.5b07237>
- Selenopeptide chemistry; M. Muttenthaler, P. F. Alewood; **J Pept Sci** 2008; **14**: 1223-1239. <https://doi.org/10.1002/psc.1075>
- Syntheses and biological activities of selenium analogs of alpha-rat atrial natriuretic peptide; T. Koide, H. Itoh, A. Otaka, M. Furuya, Y. Kitajima, N. Fujii; **Chem Pharm Bull (Tokyo)** 1993; **41**: 1596-1600. <https://doi.org/10.1248/cpb.41.1596>
- Synthesis of proteins by native chemical ligation; P. E. Dawson, T. W. Muir, I. Clark-Lewis, S. B. Kent; **Science** 1994; **266**: 776-779. <https://doi.org/10.1126/science.7973629>
- Synthesis of selenocysteine peptides and their oxidation to diselenide-bridged compounds; D. Besse, L. Moroder; **Journal of Peptide Science** 1997; **3**: 442-453. [https://doi.org/10.1002/\(sici\)1099-1387\(199711\)3:6<442::Aid-psc122>3.0.co;2-2](https://doi.org/10.1002/(sici)1099-1387(199711)3:6<442::Aid-psc122>3.0.co;2-2)
- Synthetic study on selenocystine-containing peptides; T. Koide, H. Itoh, A. Otaka, H. Yasui, M. Kuroda, N. Esaki, K. Soda, N. Fujii; **Chem Pharm Bull (Tokyo)** 1993; **41**: 502-506. <https://doi.org/10.1248/cpb.41.502>
- Total chemical synthesis of enzymes; S. Kent; **J Pept Sci** 2003; **9**: 574-593. <https://doi.org/10.1002/psc.475>
- Über Peptidsynthesen. 8. Mitteilung Bildung von S-haltigen Peptiden durch intramolekulare Wanderung von Aminoacylresten; T. Wieland, E. Bokelmann, L. Bauer, H. U. Lang, H. Lau; **Justus Liebigs Annalen der Chemie** 1953; **583**: 129-149. <https://doi.org/10.1002/jlac.19535830110>

2. Tools for the Formation of Thioesters


2.1. Dawson Linkers

The straightforward C-terminal modification of peptides assembled on a solid support used to be a significant challenge in peptide and protein chemistry. This was also true for C-terminal thioester peptides, which are important intermediates in the generation of active esters, amides, and hydrazides. Moreover, they are also an essential component of many synthetic strategies for protein synthesis. The most efficient approach for the synthesis of peptidyl thioesters used to be the *in situ* neutralization protocol for Boc solid-phase peptide synthesis using thioester linkage.

	Product details
FAA3165 Fmoc-Dbz-OH 3-[(9-Fluorenylmethyloxycarbonyl)amino]-4-amino-benzoic acid CAS-No. 1071446-05-3 Formula C ₂₂ H ₁₈ N ₂ O ₄ Mol. weight 374,39 g/mol	 
FAA3166 Fmoc-MeDbz-OH 3-[(9-Fluorenylmethyloxycarbonyl)amino]-4-(methylamino)benzoic acid CAS-No. 1788861-35-7 Formula C ₂₃ H ₂₀ N ₂ O ₄ Mol. weight 388,44 g/mol	 
FAA3167 Fmoc-Dbz(o-Boc)-OH 3-(Fmoc-amino)-4-(Boc-amino)-benzoic acid CAS-No. 1823479-63-5 Formula C ₂₇ H ₂₆ N ₂ O ₆ Mol. weight 474,51 g/mol	 
FAA3168 Fmoc-Dbz(o-Alloc)-OH 3-(Fmoc-amino)-4-(Alloc-amino)-benzoic acid CAS-No. 2143465-53-4 Formula C ₂₆ H ₂₂ N ₂ O ₆ Mol. weight 458,46 g/mol	 

[back to content ↑](#)

Ligation Technologies

		Product details
FAA3169	Fmoc-Dbz(<i>o</i> -NO ₂)-OH 4-(Fmoc-amino)-3-Nitro-benzoic acid CAS-No. 1342864-48-5 Formula C ₂₂ H ₁₆ N ₂ O ₆ Mol. weight 404,37 g/mol	 

Importantly, the linker is a stable amide during chain assembly and the key activation step utilizes the most robust reaction in solid phase peptide synthesis: the acylation of an amine. As a result, the method is compatible with amino acid side chains and protecting groups commonly used in peptide synthesis.




Fig. 5: Formation of thioesters via Dawson linker.

References:

- Alternative chemistries for the synthesis of thrombospondin-1 type 1 repeats; T. K. Tiefenbrunn, J. Blanco-Canosa, P. E. Dawson; **Biopolymers** 2010; **94**: 405-413. <https://doi.org/10.1002/bip.21486>
- An Efficient Fmoc-SPPS Approach for the Generation of Thioester Peptide Precursors for Use in Native Chemical Ligation; J. B. Blanco-Canosa, P. E. Dawson; **Angew Chem** 2008; **120**: 6957-6961. <https://doi.org/10.1002/ange.200705471>
- Protein synthesis assisted by native chemical ligation at leucine; Z. Harpaz, P. Siman, K. S. Kumar, A. Brik; **Chembiochem : a European journal of chemical biology** 2010; **11**: 1232-1235. <https://doi.org/10.1002/cbic.201000168>
- A reversible protection strategy to improve Fmoc-SPPS of peptide thioesters by the N-Acylurea approach; S. K. Mahto, C. J. Howard, J. C. Shimko, J. J. Ottesen; **Chembiochem : a European journal of chemical biology** 2011; **12**: 2488-2494. <https://doi.org/10.1002/cbic.201100472>
- A semisynthesis platform for investigating structure-function relationships in the N-terminal domain of the anthrax Lethal Factor; B. L. Pentelute, A. P. Barker, B. E. Janowiak, S. B. Kent, R. J. Collier; **ACS Chem Biol** 2010; **5**: 359-364. <https://doi.org/10.1021/cb100003r>


2.2. Hydrazone Resins

Product details

PYV1000 Fmoc-NHN=Pyv Resin

Fmoc-hydrazoneo-pyruvyl-aminomethylpolystyrene resin


Mesh Size 100-200 mesh
Loading > 0.3 mmol/g
DVB 1% DVB



BR-5279 Fmoc-NHNH-2CT Resin

Fmoc-hydrazine-2-chlorotriptyl resin


Mesh Size 100-200 mesh
Loading 0,4-1,4 mmol/g
DVB 1%



BR-5280 Fmoc-NHNH-Trt-PEG Resin

Fmoc-hydrazine-trityl polyethyleneglycol resin


Loading 0,10-0,30 mmol/g
DVB 1% DVB



PYV1100 Fmoc-L-Ala-NHN=Pyv Resin

Fmoc-L-alanyl-hydrazoneo-pyruvyl-aminomethylpolystyrene resin


Mesh Size 100-200 mesh
Loading > 0.3 mmol/g
DVB 1% DVB



PYV1110 Fmoc-L-Arg(Pbf)-NHN=Pyv Resin







Fmoc-*N'*-2,2,4,6,7-pentamethylidihydrobenzofuran-5-sulfonyl-L-arginyl-hydrazoneo-pyruvyl-aminomethylpolystyrene resin

Mesh Size 100-200 mesh
Loading > 0.3 mmol/g
DVB 1% DVB



[back to content ↑](#)

Ligation Technologies


			Product details
PYV1120	Fmoc-L-Asn(Trt)-NHN=Pyv Resin	<p>Fmoc-N-beta-trityl-L-asparaginyl-hydrazone-pyruvyl-aminomethylpolystyrene resin</p> <p>Mesh Size 100-200 mesh Loading > 0.3 mmol/g DVB 1% DVB</p>	 
PYV1130	Fmoc-L-Asp(OtBu)-NHN=Pyv Resin	<p>Fmoc-L-aspartyl-beta-t-butyl ester-alpha-hydrazono-pyruvyl-aminomethylpolystyrene resin</p> <p>Mesh Size 100-200 mesh Loading > 0.3 mmol/g DVB 1% DVB</p>	 
PYV1140	Fmoc-L-Cys(Trt)-NHN=Pyv Resin	<p>Fmoc-S-trityl-L-cysteinyl-hydrazone-pyruvyl-amino-methylpolystyrene resin</p> <p>Mesh Size 100-200 mesh Loading > 0.3 mmol/g DVB 1% DVB</p>	 
PYV1150	Fmoc-L-Glu(tBu)-NHN=Pyv Resin	<p>Fmoc-L-glutamyl-gamma-t-butyl ester-alpha-hydrazono-pyruvyl-aminomethylpolystyrene resin</p> <p>Mesh Size 100-200 mesh Loading > 0.3 mmol/g DVB 1% DVB</p>	 
PYV1160	Fmoc-L-Gln(Trt)-NHN=Pyv Resin	<p>Fmoc-N-gamma-trityl-L-glutaminyl-hydrazone-pyruvyl-aminomethylpolystyrene resin</p> <p>Mesh Size 100-200 mesh Loading > 0.3 mmol/g DVB 1% DVB</p>	 
PYV1170	Fmoc-Gly-NHN=Pyv Resin	<p>Fmoc-glycyl-hydrazone-pyruvyl-aminomethylpolystyrene resin</p> <p>Mesh Size 100-200 mesh Loading > 0.3 mmol/g DVB 1% DVB</p>	 

Product details

PYV1180 Fmoc-L-His(Trt)-NHN=Pyv Resin


Fmoc-N-trityl-L-histidyl-hydrazone-pyruvyl-amino-methylpolystyrene resin

Mesh Size 100-200 mesh
 Loading > 0.3 mmol/g
 DVB 1% DVB


PYV1190 Fmoc-L-Ile-NHN=Pyv Resin


Fmoc-L-isoleucyl-hydrazone-pyruvyl-aminomethylpolystyrene resin

Mesh Size 100-200 mesh
 Loading > 0.3 mmol/g
 DVB 1% DVB


PYV1200 Fmoc-L-Leu-NHN=Pyv Resin


Fmoc-L-leucyl-hydrazone-pyruvyl-aminomethylpolystyrene resin

Mesh Size 100-200 mesh
 Loading > 0.3 mmol/g
 DVB 1% DVB


PYV1210 Fmoc-L-Lys(Boc)-NHN=Pyv Resin


Fmoc-N-epsilon-t-butyloxycarbonyl-L-lysyl-hydrazone-pyruvyl-aminomethylpolystyrene resin

Mesh Size 100-200 mesh
 Loading > 0.3 mmol/g
 DVB 1% DVB


PYV1220 Fmoc-L-Met-NHN=Pyv Resin


Fmoc-L-methionyl-hydrazone-pyruvyl-aminomethylpolystyrene resin

Mesh Size 100-200 mesh
 Loading > 0.3 mmol/g
 DVB 1% DVB








PYV1230 Fmoc-L-Phe-NHN=Pyv Resin

Fmoc-L-phenylalanyl-hydrazone-pyruvyl-aminomethylpolystyrene resin

Mesh Size 100-200 mesh
 Loading > 0.3 mmol/g
 DVB 1% DVB


[back to content ↑](#)

Ligation Technologies

		Product details
PYV1240	Fmoc-L-Pro-NHN=Pyv Resin	 Fmoc-L-prolinyl-hydrazone-pyruvyl-aminomethylpolystyrene resin Mesh Size 100-200 mesh Loading > 0.3 mmol/g DVB 1% DVB
PYV1250	Fmoc-L-Ser(tBu)-NHN=Pyv Resin	 Fmoc-O-t-butyl-L-seryl-hydrazone-pyruvyl-aminomethylpolystyrene resin Mesh Size 100-200 mesh Loading > 0.3 mmol/g DVB 1% DVB
PYV1260	Fmoc-L-Thr(tBu)-NHN=Pyv Resin	 Fmoc-O-t-butyl-L-threonyl-hydrazone-pyruvyl-aminomethylpolystyrene resin Mesh Size 100-200 mesh Loading > 0.3 mmol/g DVB 1% DVB
PYV1270	Fmoc-L-Trp(Boc)-NHN=Pyv Resin	 Fmoc-N-t-butyloxycarbonyl-L-tryptophyl-hydrazono-pyruvyl-aminomethylpolystyrene resin Mesh Size 100-200 mesh Loading > 0.3 mmol/g DVB 1% DVB
PYV1280	Fmoc-L-Tyr(tBu)-NHN=Pyv Resin	 Fmoc-O-t-butyl-L-tyrosyl-hydrazone-pyruvyl-aminomethylpolystyrene resin Mesh Size 100-200 mesh Loading > 0.3 mmol/g DVB 1% DVB
PYV1290	Fmoc-L-Val-NHN=Pyv Resin	 Fmoc-L-valyl-hydrazone-pyruvyl-aminomethylpolystyrene resin Mesh Size 100-200 mesh Loading > 0.3 mmol/g DVB 1% DVB

Peptide hydrazides can be easily synthesized using hydrazone resins, which are polystyrene resins functionalized with the hydrazone linker (Fig. 6). The hydrazone linker is completely stable in the course of standard Fmoc-SPPS. Moreover, it tolerates treatment with 5% TFA/DCM, thus permitting selective removal of Mtt or similar acid-labile protecting groups. Subsequent application of cleavage cocktails containing neat TFA permits to obtain the desired peptides in high yields and purity.




Fig. 6: Synthesizing peptide hydrazides utilizing hydrazone resin.

Synthesized peptide hydrazides can be applied as building blocks for the chemoselective conjugation to an aldehyde or ketone group on a second peptide (or other biomacromolecule) using the hydrazone ligation technique (Fig. 7). The resulting conjugate shows a non-natural hydrazone linkage replacing the natural amide bond. This, however, is in many cases tolerated, and the conjugate displays the same properties as a natural protein fragment.



Fig. 7: Hydrazone ligation.

A peptide hydrazide can be converted into a peptide thioester. Two different protocols (Fig. 8 and Fig. 9) have been published recently. These methods include, but are not limited to, solution phase oxidation of C-terminal hydrazides to the corresponding acylazides. The required oxidants – typically used in large excess – cannot be employed in concert with N-terminal thiazolidines, a common peptide protecting strategy utilized in the synthesis of proteins requiring multiple ligations, aryl amines, or other redox-sensitive residues which can be incorporated via SPPS.




Fig. 8: Synthesizing peptide thioesters utilizing hydrazone resin via oxidation.

A more robust method for the generation of Fmoc-SPPS peptide thioesters compatible with multiple ligations utilizes mild and selective conditions at stoichiometric quantities of reactants. This process is suitable for multiple ligations to be performed sequentially in one pot, hence, has the potential for a general route to synthesize larger and complex proteins via multiple chemical ligation steps.




Fig. 9: General route for synthesizing peptide thioesters utilizing hydrazone resin.

back to content ↑

Ligation Technologies

Peptide hydrazides exposed to acetylacetone (acac, 2.5 eq.) in a solution of 6 M guanidinium chloride and 2% thiophenol at acidic conditions will initially form a hydrazine. This cyclizes to pyrazole at low pH and can be replaced with thiols like MPAA resulting in the corresponding peptide thioester.

References:

- *Chemical synthesis of proteins using peptide hydrazides as thioester surrogates; J. S. Zheng, S. Tang, Y. K. Qi, Z. P. Wang, L. Liu; Nat Protoc 2013; 8: 2483-2495. <https://doi.org/10.1038/nprot.2013.152>*
- *Convenient method of peptide hydrazide synthesis using a new hydrazone resin; P. S. Chelushkin, K. V. Polyanichko, M. V. Leko, M. Y. Dorosh, T. Bruckdorfer, S. V. Burov; Tetrahedron Letters 2015; 56: 619-622. <https://doi.org/10.1016/j.tetlet.2014.12.056>*
- *Leveraging the Knorr Pyrazole Synthesis for the Facile Generation of Thioester Surrogates for use in Native Chemical Ligation; D. T. Flood, J. C. J. Hintzen, M. J. Bird, P. A. Cistrome, J. S. Chen, P. E. Dawson; Angew Chem Int Ed Engl 2018; 57: 11634-11639. <https://doi.org/10.1002/anie.201805191>*
- *Protein chemical synthesis by ligation of peptide hydrazides; G. M. Fang, Y. M. Li, F. Shen, Y. C. Huang, J. B. Li, Y. Lin, H. K. Cui, L. Liu; Angew Chem Int Ed Engl 2011; 50: 7645-7649. <https://doi.org/10.1002/anie.201100996>*
- *Selective Coupling at the α -Amino Group of Cysteine Using Transfer Active-ester-condensation Technology to Synthesize a Linear Octadecapeptide; Y. Liao, Y. Kong, N. Hu, Z. Jin, P. Wang; Chem Lett 2010; 39: 196-197. <https://doi.org/10.1246/cl.2010.196>*



Iris
Biotech



Any Questions or Suggestions?

We are there for you – simply choose one of the numerous possibilities to get in touch!

📞 +49 (0) 9231 97121-0
📠 +49 (0) 9231 97121-99
✉️ info@iris-biotech.de
🌐 www.iris-biotech.de

3. Building Blocks for Ligation


3.1. Cysteine Building Blocks for Ligation

[Product details](#)

BAA1135 Boc-L-Thz-OH

(R)-N-t-Butyloxycarbonyl-thiazolidine-4-carboxylic acid


CAS-No. 51077-16-8
 Formula C₉H₁₅NO₄S
 Mol. weight 233,29 g/mol



BAA1186 Boc-D-Thz-OH

(S)-N-(t-Butyloxycarbonyl)-thiazolidine-4-carboxylic acid


CAS-No. 63091-82-7
 Formula C₉H₁₅NO₄S
 Mol. weight 233,29 g/mol



FAA1427 Fmoc-L-Thz-OH

(R)-N-(9-Fluorenylmethyloxycarbonyl)-thiazolidine-L-4-carboxylic acid


CAS-No. 133054-21-4
 Formula C₁₉H₁₇NO₄S
 Mol. weight 355,42 g/mol



FAA1495 Fmoc-D-Thz-OH

(S)-N-alpha-(9-Fluorenylmethyloxycarbonyl)-thiazolidine-4-carboxylic acid


CAS-No. 198545-89-0
 Formula C₁₉H₁₇NO₄S
 Mol. weight 355,42 g/mol



FAA1437 Fmoc-L-Thz(Me₂)-OH

(R)-N-(9-Fluorenylmethyloxycarbonyl)-2,2-dimethyl-thiazolidine-4-carboxylic acid

CAS-No. 873842-06-9
 Formula C₂₁H₂₁NO₄S
 Mol. weight 383,46 g/mol


[back to content ↑](#)

Ligation Technologies

Product details


FAA3160 Fmoc-D-Thz(Me2)-OH

(S)-N-(9-Fluorenylmethyloxycarbonyl)-2,2-dimethyl-thiazolidine-4-carboxylic acid

CAS-No. 1932198-36-1

Formula C₂₁H₂₁NO₄S

Mol. weight 383,46 g/mol



1,3-Thiazolidine-4-carboxylic acid (Thz) is a protected cysteine. It can be converted into free cysteine by reaction with aqueous 0.8 M methoxylamine (MeONH₂).

Product details


HAA1017 H-D-Cys-OH*HCl*H₂O

D-Cysteine Hydrochloride

CAS-No. 32443-99-5

Formula C₃H₇NO₂S*HCl*H₂O

Mol. weight 121,2*36,45*18,01 g/mol




HAA1574 H-L-Cys(Bzl)-OH

S-Benzyl-L-cysteine

CAS-No. 3054-01-1

Formula C₁₀H₁₃NO₂S

Mol. weight 211,29 g/mol




HAA6110 H-D-Cys(Bzl)-OH

S-Benzyl-D-cysteine

CAS-No. 23032-53-3

Formula C₁₀H₁₃NO₂S

Mol. weight 211,29 g/mol




HAA6070 H-L-Cys(Acm)-OH*HCl

S-(Acetyl-aminomethyl)-L-cysteine hydrochloride

CAS-No. 28798-28-9


Formula C₆H₁₂N₂O₃S*HCl

Mol. weight 192,24*36,45 g/mol




[Product details](#)
HAA6080 H-L-Cys(Bzl)-OMe*HCl


S-Benzyl-L-cysteine methyl ester hydrochloride

CAS-No. 16741-80-3
Formula C₁₁H₁₅NO₂S*HCl
Mol. weight 225,31*36,45 g/mol

HAA6090 H-L-Cys(MBzl)-OH


S-(4-Methylbenzyl)-L-cysteine

CAS-No. 42294-52-0
Formula C₁₁H₁₅NO₂S
Mol. weight 225,3 g/mol

HAA9270 H-L-Cys(MDNPE)-OH


1-[4',5'-(methylenedioxy)-2'-nitrophenyl]ethyl-L-cysteine

CAS-No. 1551078-43-3
Formula C₁₂H₁₄N₂O₆S
Mol. weight 314,31 g/mol

HAA3500 H-D-Cys(Mmt)-OH


S-p-methoxytrityl-D-cysteine

CAS-No. 926935-33-3
Formula C₂₃H₂₃NO₃S
Mol. weight 393,5 g/mol

HAA6100 H-L-Cys(Mob)-OH







S-(4-Methoxybenzyl)-L-cysteine

CAS-No. 2544-31-2
Formula C₁₁H₁₅NO₃S
Mol. weight 241,3 g/mol

HAA3510 H-L-Cys(Npys)-OH*HCl

S-(3-nitro-2-pyridylthio)-L-cysteine hydrochloride

CAS-No. 108807-66-5
Formula C₈H₉N₃O₄S₂*HCl
Mol. weight 275,30*36,45 g/mol

[back to content ↑](#)

Ligation Technologies


		Product details
HAA6160	H-L-Cys(Trt)-OH S-Trityl-L-cysteine	<p>CAS-No. 2799-07-7 Formula C₂₂H₂₁NO₂S Mol. weight 363,48 g/mol</p>  
HAA1995	H-L-Cys(Trt)-OtBu*HCl S-Trityl-L-cysteine t-butyl ester hydrochloride	<p>CAS-No. 158009-03-1 Formula C₂₆H₂₉NO₂S*HCl Mol. weight 419,58*36,45 g/mol</p>  
HAA2100	H-D-Cys(Trt)-OtBu*HCl S-Trityl-D-cysteine t-butyl ester hydrochloride	<p>CAS-No. 439089-10-8 Formula C₂₆H₂₉NO₂S*HCl Mol. weight 419,58*36,45 g/mol</p>  
HAA6140	H-L-Cys(StBu)-OH S-Thio-t-butyl-L-cysteine	<p>CAS-No. 30044-51-0 Formula C₇H₁₅NO₂S₂ Mol. weight 209,32 g/mol</p>  
HAA6150	H-L-Cys(tBu)-OH*HCl S-t-Butyl-L-cysteine hydrochloride	<p>CAS-No. 2481-09-6 Formula C₇H₁₅NO₂S*HCl Mol. weight 177,26*36,45 g/mol</p>  
HAA9320	H-L-Cys(DMNB)-OH S-(4,5-dimethoxy-2-nitrobenzyl)-L-cysteine	<p>CAS-No. 214633-68-8 Formula C₁₂H₁₆N₂O₆S Mol. weight 316,33 g/mol</p>  

[Product details](#)
AAA2015 Aloc-L-Cys(Trt)-OH
N-alpha-Allyloxycarbonyl-S-trityl-L-cysteine

CAS-No. 96865-72-4

Formula C₂₆H₂₅NO₄S


Mol. weight 447,55 g/mol


HAA2810 N₃-L-Cys(Trt)-OH*CHA
(R)-2-azido-3-(tritylthio)propanoic acid cyclohexylamine

CAS-No. 1286670-90-3

Formula C₂₂H₁₉N₃O₂S*C₆H₁₃N


Mol. weight 389,47*99,17 g/mol


BAA1083 Boc-L-Cys-OH
N-alpha-t-Butyloxycarbonyl-L-cysteine

CAS-No. 20887-95-0

Formula C₈H₁₅NO₄S


Mol. weight 221,27 g/mol


BAA1170 Boc-D-Cys-OH
N-alpha-t-Butyloxycarbonyl-D-cysteine

CAS-No. 149270-12-2

Formula C₈H₁₅NO₄S


Mol. weight 221,27 g/mol


BAA1078 Boc-L-Cys(Acm)-OH
N-alpha-t-Butyloxycarbonyl-S-(acetyl-amino-methyl)-L-cysteine

CAS-No. 19746-37-3

Formula C₁₁H₂₀N₂O₅S


Mol. weight 292,36 g/mol


BAA1510 Boc-L-Cys(Acm,O)-OH
N-alpha-t-Butyloxycarbonyl-S-(acetyl-amino-methyl)-S-oxo-L-cysteine

CAS-No. 75893-04-8

Formula C₁₁H₂₀N₂O₆S

Mol. weight 308,35 g/mol


[back to content ↑](#)

Ligation Technologies

Product details


BAA1079 Boc-L-Cys(Bzl)-OH

N-alpha-t-Butyloxycarbonyl-S-benzyl-L-cysteine

CAS-No. 5068-28-0

Formula C₁₅H₂₁NO₄S

Mol. weight 311,38 g/mol




BAA5410 Boc-D-Cys(Bzl)-OH

N-alpha-t-Butyloxycarbonyl-S-benzyl-D-cysteine

CAS-No. 102830-49-9

Formula C₁₅H₂₁NO₄S

Mol. weight 311,38 g/mol




BAA5510 Boc-L-Cys(Fm)-OH

N-alpha-t-Butyloxycarbonyl-S-(9-fluorenyl-methyl)-L-cysteine

CAS-No. 84888-35-7

Formula C₂₂H₂₅NO₄S

Mol. weight 399,51 g/mol




BAA1080 Boc-L-Cys(MBzl)-OH

N-alpha-t-Butyloxycarbonyl-S-(4-methyl-benzyl)-L-cysteine

CAS-No. 61925-77-7

Formula C₁₆H₂₃NO₄S

Mol. weight 325,43 g/mol




BAA5420 Boc-D-Cys(MBzl)-OH

N-alpha-t-Butyloxycarbonyl-S-(4-methyl-benzyl)-D-cysteine

CAS-No. 61925-78-8

Formula C₁₆H₂₃NO₄S

Mol. weight 325,43 g/mol




BAA1081 Boc-L-Cys(Mob)-OH

N-alpha-t-Butyloxycarbonyl-S-(4-methoxy-benzy-l)-L-cysteine

CAS-No. 18942-46-6

Formula C₁₆H₂₃NO₅S


Mol. weight 341,43 g/mol



[Product details](#)
BAA5430 Boc-D-Cys(Mob)-OH


N-alpha-t-Butyloxycarbonyl-S-(4-methoxy-benzylo)-D-cysteine

CAS-No. 58290-35-0
 Formula C₁₆H₂₃NO₅S
 Mol. weight 341,43 g/mol


BAA1860 Boc-L-Cys(Npys)-OH


N-alpha-t-Butyloxycarbonyl-S-(3-nitro-2-pyridylthio)-L-cysteine

CAS-No. 76880-29-0
 Formula C₁₃H₁₇N₃O₆S₂
 Mol. weight 375,42 g/mol


BAA6390 Boc-L-Cys(Phacm)-OH


N-alpha-t-Butyloxycarbonyl-S-(Phenylacetylaminomethyl)-L-cysteine

CAS-No. 57084-73-8
 Formula C₁₇H₂₄N₂O₅S
 Mol. weight 368,45 g/mol


BAA3140 Boc-L-Cys(Ph)-OH


N-alpha-t-Butyloxycarbonyl-S-phenyl-L-cysteine

CAS-No. 163705-28-0
 Formula C₁₄H₁₉NO₄S
 Mol. weight 297,37 g/mol


BAA1082 Boc-L-Cys(tBu)-OH


N-alpha-t-Butyloxycarbonyl-S-t-butyl-L-cysteine

CAS-No. 56976-06-8
 Formula C₁₂H₂₃NO₄S
 Mol. weight 277,37 g/mol








BAA1084 Boc-L-Cys(Trt)-OH

N-alpha-t-Butyloxycarbonyl-S-trityl-L-cysteine

CAS-No. 21947-98-8
 Formula C₂₇H₂₉NO₄S
 Mol. weight 463,59 g/mol


[back to content ↑](#)

Ligation Technologies


		Product details
BAA5000	Boc-D-Cys(Trt)-OH N-alpha-t-Butyloxycarbonyl-S-trityl-D-cysteine CAS-No. 87494-13-1 Formula C ₂₇ H ₂₉ NO ₄ S Mol. weight 463,59 g/mol	 
IAD1040	Boc-L-Ser[Fmoc-L-Cys(Trt)]-OH O-(N-(((9H-fluoren-9-yl)methoxy)carbonyl)-S-trityl-L-cysteinyl)-N-(tert-butoxycarbonyl)-L-serine Formula C ₄₅ H ₄₄ N ₂ O ₈ S Mol. weight 772,9	 
IAD2040	Boc-L-Thr[Fmoc-L-Cys(Trt)]-OH O-(N-(((9H-fluoren-9-yl)methoxy)carbonyl)-S-trityl-L-cysteinyl)-N-(tert-butoxycarbonyl)-L-threonine CAS-No. 944283-30-1 Formula C ₄₆ H ₄₆ N ₂ O ₈ S Mol. weight 786,93 g/mol	 
FAA1362	Fmoc-L-Cys-OH*H ₂ O N-alpha-(9-Fluorenylmethyloxycarbonyl)-L-cysteine monohydrat CAS-No. 135248-89-4 net Formula C ₁₈ H ₁₇ NO ₄ S*H ₂ O Mol. weight 343,40*18,01 g/mol	 
FAA1470	Fmoc-D-Cys-OH*H ₂ O N-alpha-(9-Fluorenylmethyloxycarbonyl)-D-cysteine monohydrat CAS-No. 157355-80-1 net Formula C ₁₈ H ₁₇ NO ₄ S*H ₂ O Mol. weight 343,4*18,01 g/mol	 
FAA5150	Fmoc-L-Cys(Aapam)-OH N-alpha-(9-Fluorenylmethyloxycarbonyl)-S-((4-(allyloxycarbonylamino)phenylacetylaminomethyl)-L-cysteine CAS-No. 1946783-89-6 Formula C ₃₁ H ₃₁ N ₃ O ₇ S Mol. weight 589,66 g/mol	 

Product details

FAA1506 Fmoc-L-Cys(Acm)-OH


N-alpha-(9-Fluorenylmethyloxycarbonyl)-S-(acetylaminomethyl)-L-cysteine

CAS-No. 86060-81-3
 Formula C₂₁H₂₂N₂O₅S
 Mol. weight 414,48 g/mol


FAA6230 Fmoc-D-Cys(Acm)-OH


N-alpha-(9-Fluorenylmethyloxycarbonyl)-S-(acetylaminomethyl)-D-cysteine

CAS-No. 168300-88-7
 Formula C₂₁H₂₂N₂O₅S
 Mol. weight 414,48 g/mol


FAA4751 Fmoc-L-Cys(Ac-OtBu)-OH*DCHA


N-alpha-(9-Fluorenylmethyloxycarbonyl)-S-(t-butoxycarbonylmethyl)-L-cysteine dicyclohexylamine

CAS-No. 269730-62-3 net
 Formula C₂₄H₂₇NO₆S*C₁₂H₂₃N
 Mol. weight 457,54*181,32 g/mol


FAA7610 Fmoc-L-Cys(Allocam)-OH


N-alpha-(9-Fluorenylmethyloxycarbonyl)-S-((allyloxy carbonylamino)methyl)-L-cysteine

CAS-No. 232953-09-2
 Formula C₂₃H₂₄N₂O₆S
 Mol. weight 456,51 g/mol


FAA6270 Fmoc-L-Cys(Bzl)-OH


N-alpha-(9-Fluorenylmethyloxycarbonyl)-S-benzyl-L-cysteine

CAS-No. 53298-33-2
 Formula C₂₅H₂₃NO₄S
 Mol. weight 433,52 g/mol


FAA6940 Fmoc-L-Cys(Ddm)-OH

N-alpha-(9-Fluorenylmethyloxycarbonyl)-S-di(4-methoxyphenyl)methyl-L-cysteine

CAS-No. 1403825-56-8
 Formula C₃₃H₃₁NO₆S
 Mol. weight 569,67 g/mol


[back to content ↑](#)


Ligation Technologies

Product details

FAA6950 Fmoc-D-Cys(Ddm)-OH

N-alpha-(9-Fluorenylmethyloxycarbonyl)-S-di(4-methoxyphenyl)methyl-D-cysteine


Formula $C_{33}H_{31}NO_6S$
Mol. weight 569,67 g/mol



FAA3190 Fmoc-L-Cys(Dpm)-OH

N-alpha-(9-Fluorenylmethyloxycarbonyl)-S-diphenylmethyl-L-cysteine


CAS-No. 247595-29-5
Formula $C_{31}H_{27}NO_4S$
Mol. weight 509,62 g/mol



FAA5650 Fmoc-D-Cys(Dpm)-OH

N-alpha-(9-Fluorenylmethyloxycarbonyl)-S-diphenylmethyl-D-cysteine


Formula $C_{31}H_{27}NO_4S$
Mol. weight 509,62 g/mol



FAA1714 Fmoc-L-Cys(MBzl)-OH

N-alpha-(9-Fluorenylmethyloxycarbonyl)-S-(4-methylbenzyl)-L-cysteine


CAS-No. 136050-67-4
Formula $C_{36}H_{25}NO_4S$
Mol. weight 447,53 g/mol



FAA1030 Fmoc-L-Cys(Mmt)-OH

N-alpha-(9-Fluorenylmethyloxycarbonyl)-S-p-methoxytrityl-L-cysteine


CAS-No. 177582-21-7
Formula $C_{38}H_{33}NO_5S$
Mol. weight 615,74 g/mol



FAA1614 Fmoc-D-Cys(Mmt)-OH

N-alpha-(9-Fluorenylmethyloxycarbonyl)-S-p-methoxytrityl-D-cysteine

CAS-No. 1198791-73-9
Formula $C_{38}H_{33}NO_5S$
Mol. weight 615,74 g/mol




Product details

FAA1715 Fmoc-L-Cys(Mob)-OH


N-alpha-(9-Fluorenylmethyloxycarbonyl)-S-(4-methoxybenzyl)-L-cysteine

CAS-No. 141892-41-3
 Formula C₂₆H₂₅NO₅S
 Mol. weight 463,55 g/mol


FAA4155 Fmoc-L-Cys(Msbh)-OH


N-alpha-(9-Fluorenylmethyloxycarbonyl)-S-(4,4'-dimethylsulfonylbenzhydryl)-L-cysteine

CAS-No. 1584646-97-8
 Formula C₃₃H₃₁NO₆S₃
 Mol. weight 633,80 g/mol


FAA8420 Fmoc-L-Cys(NDBF)-OH


N-alpha-(9-Fluorenylmethyloxycarbonyl)-S-(1-(3-nitro-dibenzofuran-2-yl)-ethyl)-L-cysteine

CAS-No. 1895883-28-9
 Formula C₃₂H₂₆N₂O₇S
 Mol. weight 582,62 g/mol


FAA8150 Fmoc-D-Cys(Msbh)-OH


N-alpha-(9-Fluorenylmethyloxycarbonyl)-S-(4,4'-dimethylsulfonylbenzhydryl)-D-cysteine

Formula C₃₃H₃₁NO₆S₃
 Mol. weight 633,80 g/mol


FAA1975 Fmoc-L-Cys(Npys)-OH


N-alpha-(9-Fluorenylmethyloxycarbonyl)-S-(3-nitro-2-pyridylthio)-L-cysteine

CAS-No. 159700-51-3
 Formula C₂₃H₁₉N₃O₆S₂
 Mol. weight 497,54 g/mol


FAA3970 Fmoc-L-Cys(oNv)-OH

N-alpha-(9-Fluorenylmethyloxycarbonyl)-S-(2-nitroveratryl)-L-cysteine

CAS-No. 214633-71-3
 Formula C₂₇H₂₆N₂O₈S
 Mol. weight 538,57 g/mol


[back to content ↑](#)

Ligation Technologies

Product details


FAA6910 Fmoc-L-Cys(Phacm)-OH

N-alpha-(9-Fluorenylmethyloxycarbonyl)-S-Phenylacetylaminomethyl-L-cysteine

CAS-No. 159680-21-4

Formula C₂₇H₂₆N₂O₅S

Mol. weight 490,57 g/mol




FAA3710 Fmoc-D-Cys(Phacm)-OH

N-alpha-(9-Fluorenylmethyloxycarbonyl)-S-(phenylacetylaminomethyl)-D-cysteine

CAS-No. 1565818-55-4

Formula C₂₇H₂₆N₂O₅S

Mol. weight 490,57 g/mol




FAA8495 Fmoc-L-Cys(SIT)-OH

N-alpha-(9-Fluorenylmethyloxycarbonyl)-S-(sec-isoamyl mercaptan)-L-cysteine

CAS-No. 2545642-31-5

Formula C₂₃H₂₇NO₄S₂

Mol. weight 445,59 g/mol




FAA1575 Fmoc-L-Cys(StBu)-OH

N-alpha-(9-Fluorenylmethyloxycarbonyl)-S-(t-butylthio)-L-cysteine

CAS-No. 73724-43-3

Formula C₂₂H₂₅NO₄S₂

Mol. weight 431,57 g/mol




FAA1965 Fmoc-D-Cys(StBu)-OH

N-alpha-(9-Fluorenylmethyloxycarbonyl)-S-(t-butylthio)-D-cysteine

CAS-No. 501326-55-2

Formula C₂₂H₂₅NO₄S₂

Mol. weight 431,57 g/mol




FAA1716 Fmoc-L-Cys(tBu)-OH

N-alpha-(9-Fluorenylmethyloxycarbonyl)-S-t-butyldimethyl-L-cysteine

CAS-No. 67436-13-9

Formula C₂₂H₂₅NO₄S


Mol. weight 399,51 g/mol



[Product details](#)
FAA4160 Fmoc-L-Cys(Thp)-OH


N-alpha-(9-Fluorenylmethyloxycarbonyl)-S-tetrahydro-pyranyl-L-cysteine

CAS-No. 1673576-83-4
 Formula C₂₃H₂₅NO₅S
 Mol. weight 427,15 g/mol


FAA1040 Fmoc-L-Cys(Trt)-OH


N-alpha-(9-Fluorenylmethyloxycarbonyl)-S-trityl-L-cysteine

CAS-No. 103213-32-7
 Formula C₃₇H₃₁NO₄S
 Mol. weight 585,71 g/mol


FAA1035 Fmoc-D-Cys(Trt)-OH


N-alpha-(9-Fluorenylmethyloxycarbonyl)-S-trityl-D-cysteine

CAS-No. 167015-11-4
 Formula C₃₇H₃₁NO₄S
 Mol. weight 585,71 g/mol


FAA5670 Fmoc-L-Cys(Trt)-OMe


N-alpha-(9-Fluorenylmethyloxycarbonyl)-S-trityl-L-cysteine methyl ester

CAS-No. 245088-56-6
 Formula C₃₈H₃₃NO₄S
 Mol. weight 599,74 g/mol


TAA1508 Trt-L-Cys(Trt)-OH*DEA


N-alpha-S-Bistrityl-L-cysteine diethylamine

CAS-No. 27486-88-0
 Formula C₄₁H₃₅NO₂S*C₄H₁₁N
 Mol. weight 678,9 g/mol


ZAA1161 Z-L-Cys(Bzl)-OH

N-alpha-Benzylloxycarbonyl-S-benzyl-L-cysteine

CAS-No. 3257-18-9
 Formula C₁₈H₁₉NO₄S
 Mol. weight 345,42 g/mol


[back to content ↑](#)

Ligation Technologies

Product details


ZAA1310 Z-L-Cys(Trt)-OH

N-alpha-Benzyloxycarbonyl-S-trityl-L-cysteine

CAS-No. 26311-04-6

Formula C₃₀H₂₇NO₄S

Mol. weight 497,60 g/mol




PSI1450 Fmoc-L-Ala-L-Cys[PSI(Dmp,H)pro]-OH

(S)-3-(N-(9-Fluorenylmethyloxycarbonyl)-L-alanyl)-2-(2,4-dimethoxyphenyl)thiazolidine-4-carboxylic acid

CAS-No. 2022956-37-0

Formula C₃₀H₃₀N₂O₇S

Mol. weight 562,63 g/mol




PSI1460 Fmoc-L-Asn(Trt)-L-Cys[PSI(Dmp,H)pro]-OH

(S)-3-(N-(9-Fluorenylmethyloxycarbonyl)-beta-trityl-L-asparaginyl)-2-(2,4-dimethoxyphenyl)thiazolidine-4-carboxylic acid

Formula C₅₀H₄₅N₃O₈S

Mol. weight 847,97 g/mol




PSI1470 Fmoc-L-Asp(tBu)-L-Cys[PSI(Dmp,H)pro]-OH

(S)-3-(N-(9-Fluorenylmethyloxycarbonyl)-beta-t-butyl-L-aspartyl)-2-(2,4-dimethoxyphenyl)thiazolidine-4-carboxylic acid

CAS-No. 1359754-16-7

Formula C₃₅H₃₉N₂O₉S

Mol. weight 662,75 g/mol




PSI1580 Fmoc-L-Cys(Trt)-L-Cys[PSI(Dmp,H)pro]-OH

(R)-3-(N-(9-Fluorenylmethyloxycarbonyl)-S-trityl-L-cysteinyl)-2-(2,4-dimethoxyphenyl)thiazolidine-4-carboxylic acid

CAS-No. 2022956-75-6

Formula C₄₉H₄₄N₂O₇S₂

Mol. weight 837,01 g/mol




PSI1480 Fmoc-L-Gln(Trt)-L-Cys[PSI(Dmp,H)pro]-OH

(S)-3-(N-(9-Fluorenylmethyloxycarbonyl)-gamma-trityl-L-glutaminyl)-2-(2,4-dimethoxyphenyl)thiazolidine-4-carboxylic acid

Formula C₅₁H₄₇N₃O₈S


Mol. weight 862,00 g/mol



PSI1490 Fmoc-L-Glu(tBu)-L-Cys[PSI(Dmp,H)pro]-OH


(S)-3-(N-(9-Fluorenylmethyloxycarbonyl)-gamma-t-butyl-L-glutamyl)-2-(2,4-dimethoxyphenyl)thiazolidine-4-carboxylic acid

Formula C₃₆H₄₀N₂O₉S
Mol. weight 676,78 g/mol

**PSI1440 Fmoc-Gly-L-Cys[PSI(Dmp,H)pro]-OH**


(S)-3-(N-(9-Fluorenylmethyloxycarbonyl)-glycyl)-2-(2,4-dimethoxyphenyl)thiazolidine-4-carboxylic acid

CAS-No. 1926163-05-4
Formula C₂₉H₂₈N₂O₇S
Mol. weight 548,61 g/mol

**PSI1500 Fmoc-L-Ile-L-Cys[PSI(Dmp,H)pro]-OH**


(S)-3-(N-(9-Fluorenylmethyloxycarbonyl)-L-isoleucyl)-2-(2,4-dimethoxyphenyl)thiazolidine-4-carboxylic acid

Formula C₃₃H₃₆N₂O₈S
Mol. weight 604,71 g/mol

**PSI1510 Fmoc-L-Leu-L-Cys[PSI(Dmp,H)pro]-OH**


(S)-3-(N-(9-Fluorenylmethyloxycarbonyl)-L-leucyl)-2-(2,4-dimethoxyphenyl)thiazolidine-4-carboxylic acid

CAS-No. 1926163-06-5
Formula C₃₃H₃₆N₂O₈S
Mol. weight 604,71 g/mol

**PSI1520 Fmoc-L-Lys(Boc)-L-Cys[PSI(Dmp,H)pro]-OH**


(S)-3-(N-(9-Fluorenylmethyloxycarbonyl)-N-epsi-lon-t-butyloxycarbonyl-L-lysyl)-2-(2,4-dimethoxyphenyl)thiazolidine-4-carboxylic acid

CAS-No. 1926163-07-6
Formula C₃₈H₄₅N₃O₉S
Mol. weight 719,84 g/mol





**PSI1530 Fmoc-L-Phe-L-Cys[PSI(Dmp,H)pro]-OH**

(S)-3-(N-(9-Fluorenylmethyloxycarbonyl)-L-phenylalanyl)-2-(2,4-dimethoxyphenyl)thiazolidine-4-carboxylic acid

Formula C₃₆H₃₄N₂O₈S
Mol. weight 638,73 g/mol




Ligation Technologies

		Product details
PSI1550	Fmoc-L-Thr(<i>t</i> Bu)-L-Cys[PSI(Dmp,H)pro]-OH (S)-3-(N-(9-Fluorenylmethyloxycarbonyl)-O- <i>t</i> -butyl-L-threonyl)-2-(2,4-dimethoxyphenyl)thiazolidine-4-carboxylic acid Formula C ₃₅ H ₄₀ N ₂ O ₈ S Mol. weight 648,77 g/mol	 
PSI1560	Fmoc-L-Tyr(<i>t</i> Bu)-L-Cys[PSI(Dmp,H)pro]-OH (S)-3-(N-(9-Fluorenylmethyloxycarbonyl)-O- <i>t</i> -butyl-L-thyrosyl)-2-(2,4-dimethoxyphenyl)thiazolidine-4-carboxylic acid Formula C ₃₆ H ₄₂ N ₂ O ₈ S Mol. weight 710,84 g/mol	 
PSI1570	Fmoc-L-Val-L-Cys[PSI(Dmp,H)pro]-OH (S)-3-(N-(9-Fluorenylmethyloxycarbonyl)-L-valyl)-2-(2,4-dimethoxyphenyl)thiazolidine-4-carboxylic acid CAS-No. 1926163-08-7 Formula C ₃₂ H ₃₄ N ₂ O ₈ S Mol. weight 590,69 g/mol	 
SAA1110	Smoc-L-Cys(Trt)-OH N-(((2,7-disulfo-9H-fluoren-9-yl)methoxy)carbonyl)-S-trityl-L-cysteine potassium salt CAS-No. 2442552-68-1 (net) Formula C ₃₇ H ₂₉ K ₂ NO ₁₀ S ₃ Mol. weight 822,01 g/mol	 

3.2. Ligation at the Position of Glycine

On average, cysteine occurs at a very low frequency in natural peptide sequences, a fact that frequently causes challenges when planning NCL strategies. On the other hand, glycine is a much more common amino acid, and being the simplest amino acid in nature, lends itself ideally for the construction of an NCL auxiliary.

		Product details
PAA2000	tBu-SS-Photo(Fmoc)-Gly-OH Photocleavable-NCL-auxiliary-Gly-OH CAS-No. 1994388-93-0 Formula C ₃₆ H ₄₄ N ₄ O ₉ S ₂ Mol. weight 740,89 g/mol	 

The building block tBu-SS-Photo(Fmoc)-Gly-OH is based on a glycine scaffold that is N-alkylated with a photolabile thiol-bearing auxiliary. This auxiliary mimics the action of an N-terminal cysteine's sulphydryl group in NCL, and can be tracelessly removed following ligation, leaving only a glycine behind. Being able to perform NCL at the position of glycine residues enables the synthetic chemist with a much higher flexibility for the design of peptides and proteins. Following SPPS, this building block is attached to the N-terminus of a peptide sequence in lieu of a glycine residue. The auxiliary's Fmoc-protected amino functionality can subsequently be deprotected and functionalized, e.g., with a PEG, which is useful for increasing the solubility of peptide fragments, and for facilitating their purification by precipitation with EtOH/Et₂O. This is particularly valuable if the peptide's amino acid side chains are supposed to be further derivatized post-SPPS.




Fig. 10: Photo-cleavable glycine auxiliary building block enabling native chemical ligation at the position of glycine. It has additional capabilities to be conjugated with a solubilizing tag, e.g. PEG, in order to improve the solubility of hydrophobic peptide sequences.

back to content ↑

Ligation Technologies

Reference:

- A PEGylated photocleavable auxiliary mediates the sequential enzymatic glycosylation and native chemical ligation of peptides; C. Bello, S. Wang, L. Meng, K. W. Moremen, C. F. Becker; *Angew Chem Int Ed Engl* 2015; **54**: 7711-7715. <https://doi.org/10.1002/anie.201501517>


3.3. Building Blocks and Reagents for Ligation with Seleno Amino Acids

Our portfolio contains a variety of selenocysteine (Sec) derivatives as well as selenazolidine carboxylic acids (Sez derivatives). Sez can be deprotected and converted to Sec by treatment with O-methylhydroxylamine (MeONH_2) at pH 4 or by using Cu(II) salts.

Product details


HAA9255 H-L-Sec(DMNB)-OH*TFA

Dimethoxynitrobenzyl selenocysteine TFA salt
CAS-No. 1644398-13-9
Formula $\text{C}_{12}\text{H}_{16}\text{N}_2\text{O}_6\text{Se}^*\text{CF}_3\text{COOH}$
Mol. weight 363,24*114,02 g/mol




HAA9360 H-L-Sec(MDNPE)-OH

Se-(Methyl-*o*-nitropiperonyl)-selenocysteine
CAS-No. 2235373-47-2
Formula $\text{C}_{12}\text{H}_{14}\text{N}_2\text{O}_6\text{Se}$
Mol. weight 361,21 g/mol



HAA9230 H-L-Sec(MDNPE)*TFA

(2*R*)-2-amino-3-((1-(6-nitrobenzo[d][1,3]dioxol-5-yl)ethyl)selanyl)propanoic acid trifluoroacetate
CAS-No. 2235373-48-3
Formula $\text{C}_{12}\text{H}_{14}\text{N}_2\text{O}_6\text{Se}^*\text{CF}_3\text{CO}_2\text{H}$
Mol. weight 361,22*114,02 g/mol



Product details


HAA9350 (H-L-Sec-OH)2

L-Selenocystine, (H-Sec)2, (H-L-Sec)2

CAS-No. 29621-88-3

 Formula C₆H₁₂N₂O₄Se₂

Mol. weight 334,11 g/mol



BAA3680 (Boc-L-Sec)₂

N-alpha-t-Butyloxycarbonyl-L-selenocystine

CAS-No. 877754-71-7

 Formula C₁₆H₂₈N₂O₈Se₂

Mol. weight 534,35 g/mol



BAA3760 Boc-L-Sec(Mob)-OH

N-alpha-t-Butyloxycarbonyl-Se-(4-methoxybenzyl)-L-selenocysteine

CAS-No. 959415-39-5

 Formula C₁₆H₂₃NO₅Se

Mol. weight 388,32 g/mol



BAA4830 Boc-L-Sec(Mob)-OPfp

N-alpha-tert-Butoxycarbonyl-4-methoxybenzyl-L-selenocysteine pentafluorophenyl ester

CAS-No. 1257525-48-6

 Formula C₂₂H₂₂F₅NO₅Se

Mol. weight 554,38 g/mol



FAA8705 Fmoc-L-Sec(Mob)-OH

N-alpha-(9-Fluorenylmethyloxycarbonyl)-Se-(4-methoxybenzyl)-L-selenocysteine

CAS-No. 150308-80-8

 Formula C₂₆H₂₅NO₅Se

Mol. weight 510,46 g/mol



FAA8760 Fmoc-L-Sec(Mob)-OPfp

N-alpha-(9-Fluorenylmethyloxycarbonyl)-L-4-methoxybenzyl selenocysteine pentafluorophenyl ester





CAS-No. 939431-43-3

 Formula C₃₂H₂₄F₅NO₅Se


Mol. weight 676,51 g/mol


[back to content ↑](#)

Ligation Technologies

		Product details
FAA8710 Fmoc-D-Sec(Mob)-OH	<p>N-alpha-(9-Fluorenylmethyloxycarbonyl)-Se-(4-methoxybenzyl)-D-selenocysteine</p> <p>Formula $C_{26}H_{25}NO_5Se$ Mol. weight 510,46 g/mol</p>	 
FAA8465 Fmoc-L-Sec(Xan)-OH	<p>Fmoc-Se-xanthyl-L-selenocysteine</p> <p>CAS-No. 1639843-35-8 Formula $C_{31}H_{25}NO_5Se$ Mol. weight 570,49 g/mol</p>	 
FAA8600 Fmoc-D-Sec(Xan)-OH	<p>Fmoc-Se-xanthyl-D-selenocysteine</p> <p>Formula $C_{31}H_{25}NO_5Se$ Mol. weight 570,49 g/mol</p>	 
FAA4205 Fmoc-L-Selenomethionine	<p>N-alpha-(9-Fluorenylmethyloxycarbonyl)-L-selenomethionine</p> <p>CAS-No. 1217852-49-7 Formula $C_{20}H_{21}NO_4Se$ Mol. weight 418,35 g/mol</p>	 

Besides Sec, our portfolio contains a variety of selenazolidine carboxylic acids (Sez derivatives). Sez can be deprotected and converted to Sec by treatment with O-methylhydroxylamine ($MeONH_2$) at pH 4 or by using Cu(II) salts.

		Product details
BAA4880 Boc-L-Sez-OH	<p>Boc selenazolidine carboxylic acid</p> <p>CAS-No. 1841180-44-6 Formula $C_9H_{15}NO_4Se$ Mol. weight 280,19 g/mol</p>	 


FAA8860 Fmoc-L-Sez-OH

Fmoc selenazolidine carboxylic acid

CAS-No. 1985651-74-8

Formula C₁₉H₁₇NO₄Se

Mol. weight 402,31 g/mol

**General Procedure for the Synthesis of Selenoesters**

The peptide residue after Fmoc-SPPS is dissolved in anhydrous DMF and cooled to 0 °C. Diphenyl diselenide (DPDS) (30 eq. in DMF) and Bu₃P (30 eq.) are subsequently added. After 3 h at 0 °C, the solvent is removed *in vacuo*.




Fig. 11: Formation of selenoesters.

One-Pot Additive-Free Diselenide-Selenoester Ligation-Deselenization Reactions**Conditions of Additive-free Ligation:**

2.5 mM final concentration of diselenide dimer in 6 M Gdn*HCl, and 0.1 M Na₂HPO₄ (pH 7.2; reduced to 6.2–6.5 upon addition to peptide fragments).



Conditions of One-Pot Deselenization:

Hexane extraction (×5) followed by the addition of 0.25 M TCEP, 25 mM DTT in 6 M Gdn*HCl, and 0.1 M Na₂HPO₄ (pH 5–6).

Reference:

- Accelerated Protein Synthesis via One-Pot Ligation-Deselenization Chemistry; N. J. Mitchell, J. Sayers, S. S. Kulkarni, D. Clayton, A. M. Goldys, J. Ripoll-Rozada, P. J. Barbosa Pereira, B. Chan, L. Radom, R. J. Payne; *Chem* 2017; **2**: 703–715. <https://doi.org/10.1016/j.chempr.2017.04.003>

Ligation Technologies

		Product details
LS-3405	TCEP*HCl	<p>3-[bis(2-carboxyethyl)phosphanyl]propanoic acid, hydrochloride</p> <p>CAS-No. 51805-45-9 Formula C₉H₁₅O₆P*HCl Mol. weight 250,19*36,45 g/mol</p> 
RL-1020	DTT (racemic)	<p>DL-Dithiothreitol</p> <p>CAS-No. 3483-12-3 Formula C₆H₁₀O₂S₂ Mol. weight 154,25 g/mol</p> 

3.4. Oxime Ligation using Aminooxy-Amino Acids

Replacing the primary amino group of an amino acid by an aminooxy moiety leads to an increase in nucleophilicity. Thus, after completion of the peptide synthesis and deprotection of the aminooxy-function, chemoselective reactions with carbonyl compounds under formation of a kinetically stable oxime bond can be performed. Compared to imines, oximes display much higher stability toward hydrolysis. This increased stability is explained by the α -effect provided by the heteroatom adjacent to the sp^2 nitrogen.

This simple reaction can be utilized for the synthesis of large proteins by fragment assembly = oxime ligation. Successive deprotection of aldehyde and/or aminooxy moieties allows for the controlled and sequential ligation of peptide fragments.







Fig. 12: Illustration of oxime ligation.

The oxime linkage is reported as good peptidomimetic instead of the amide bond. The low usage of this ligation type is most probably due to the possible hydrolysis of the oxime linkage as well as the formation of a mixture of products due to *E/Z*-isomerization.

		Product details
FAA8450	Fmoc-AAHA(Boc)-OH (S)	 
(S)-2-(Fmoc-amino)-6-(Boc-aminoxy)hexanoic acid		
CAS-No.	357278-11-6	
Formula	C ₂₆ H ₃₂ N ₂ O ₇	
Mol. weight	484,54 g/mol	
FAA8455	Fmoc-L-trans-Hyp(NHBoc)-OH	 
Fmoc-4-(Boc-amino)oxy-proline (2S,4R)		
CAS-No.	1015426-45-5	
Formula	C ₂₅ H ₂₈ N ₂ O ₇	
Mol. weight	468,50 g/mol	
FAA8460	Fmoc-L-cis-Hyp(NHBoc)-OH	 
Fmoc-4-(Boc-amino)oxy-proline (2S,4S)		
CAS-No.	1015426-31-9	
Formula	C ₂₅ H ₂₈ N ₂ O ₇	
Mol. weight	468,50 g/mol	

References:

- Site-specific cross-linking of proteins to DNA via a new biorthogonal approach employing oxime ligation; S. S. Pujari, Y. Zhang, S. Ji, M. D. Distefano, N. Y. Tretyakova; **Chem. Commun.** 2018; **54**: 6296-6299. <https://doi.org/10.1039/C8CC01300D>
- 44. Amino-oxy-derivatives. Part I. Some a-amino-oxy-acids and a-amino-oxy-hydrazides. D. McHale, J. Green, P. Mamalis; **J. Chem. Soc.** 1960; 225-229. <https://doi.org/10.1039/JR9600000225>
- SAR by Oxime-Containing Peptide Libraries: Application to Tsg101 Ligand Optimization; F. Liu, A. G. Stephen, A. A. Waheed, M. J. Aman, E. O. Freed, R. J. Fisher, T. R. Burke; **ChemBioChem** 2008; **9(12)**: 2000-2004. <https://doi.org/10.1002/cbic.200800281>
- A Versatile Set of Aminoxy Amino Acids for the Synthesis of Neoglycopeptides; M. R. Carrasco, R. T. Brown; **J. Org. Chem.** 2003; **68**: 8853-8858. <https://doi.org/10.1021/jo034984x>

3.5. Activated Cysteine-Based Protein Ligation (ACPL)

This method is based on the direct activation of Cysteine by the small molecule cyanating reagent 2-nitro-5-thiocyanatobenzoic acid (NTCB). The transferred cyanide then forms a thiocyanate which undergoes a reversible intramolecular addition with the Cys N-amide to generate a five-membered 1-acyl-2-iminothiazolidine intermediate which can be reacted with nucleophiles. When the nucleophile is hydrazine, the afforded peptide hydrazide can then undergo ligation, either *via* reaction with an aldehyde/ketone (hydrazone ligation) or *via* transformation to a peptide thioester and subsequent chemical ligation, both as described above. Without addition of a nucleophilic amine, the formed five-membered ring undergoes hydrolysis.

back to content ↑

Ligation Technologies





Fig. 13: Schematic illustration of the activated Cysteine-based protein ligation.

		Product details
RL-4080	NTCB	 

Reference:

- Expressed Protein Ligation without Intein; Y. Qiao, G. Yu, K. C. Kratch, X. Aria Wang, W. Wei Wang, S. Z. Leeuwon, S. Xu, J. S. Morse, W. Ray Liu; *J. Am. Chem. Soc.* 2002; **142**: 7057-7054.
<https://doi.org/10.1021/jacs.0c00252>

4. Enzyme-mediated Ligation

The transpeptidase activity of sortase can be used to produce fusion proteins *in vitro*. The enzyme recognition motif (LPXTG) is added to the C-terminus of a protein of interest, while an oligo (glycine) motif is added to the N-terminus of the second protein. Upon addition of sortase A, the two peptides are covalently linked through a native peptide bond while losing one glycine residue.






Fig. 14: Enzymatically supported conjugation utilizing sortase A.

Here we present several building blocks, which can be used for convenient N-terminal oligo-glycine design of ligation fragments.

		Product details
HAA2860	N₃-Gly-Gly-Gly-Gly-OH	
CAS-No.	2250433-77-1	
Formula	C ₁₀ H ₁₅ N ₇ O ₆	
Mol. weight	329,27 g/mol	
HAA2840	N₃-Gly-Gly-Gly-OH	
Azido-glycylglycylglycine		
CAS-No.	1993176-75-2	
Formula	C ₆ H ₉ N ₅ O ₄	
Mol. weight	215,17 g/mol	
HAA2850	N₃-Gly-Gly-OH*DCHA	
Azido-glycylglycine dicyclohexylamine		
CAS-No.	855750-87-7 net	
Formula	C ₄ H ₆ N ₄ O ₃ *C ₁₂ H ₂₃ N	
Mol. weight	158,12*181,32 g/mol	

back to content ↑

Ligation Technologies

		Product details
FDP1030	Fmoc-Gly-Gly-OH	 
N-alpha-(9-Fluorenylmethyloxycarbonyl)-glycinyl-glycine	CAS-No. 35665-38-4 Formula C ₁₉ H ₁₈ N ₂ O ₅ Mol. weight 354,36 g/mol	
HAA2870	H-(Gly)₃-Lys(N₃)-OH*HCl	 
Triglycyl-epsilon-azido-L-lysine hydrochloride	CAS-No. 2250437-45-5 net Formula C ₁₂ H ₂₁ N ₃ O ₅ *HCl Mol. weight 343,34*36,45 g/mol	

References:

- Assembly of Oligoglycine Layers on Mica Surface; S. V. Tsygankova, A. A. Chinarev, A. B. Tuzikov, I. S. Zaitsev, N. Severin, A. A. Kalachev, J. P. Rabe, N. V. Bovin; *Journal of Biomaterials and Nanobiotechnology* 2011; **02**: 91-97. <https://doi.org/10.4236/jbnb.2011.21012>
- Biantennary oligoglycines and glyco-oligoglycines self-associating in aqueous medium; S. V. Tsygankova, A. A. Chinarev, A. B. Tuzikov, N. Severin, A. A. Kalachev, J. P. Rabe, A. S. Gambaryan, N. V. Bovin; *Beilstein J Org Chem* 2014; **10**: 1372-1382. <https://doi.org/10.3762/bjoc.10.140>
- DNA-bending finger: artificial design of 6-zinc finger peptides with polyglycine linker and induction of DNA bending; M. Imanishi, Y. Hori, M. Nagaoka, Y. Sugiura; *Biochemistry* 2000; **39**: 4383-4390. <https://doi.org/10.1021/bi992989b>
- Effects of the length of a glycine linker connecting the N-and C-termini of a circularly permuted dihydrofolate reductase; M. Iwakura, T. Nakamura; *Protein Eng* 1998; **11**: 707-713. <https://doi.org/10.1093/protein/11.8.707>
- Linkers in the structural biology of protein-protein interactions; V. P. Reddy Chichili, V. Kumar, J. Sivaraman; *Protein Sci* 2013; **22**: 153-167. <https://doi.org/10.1002/pro.2206>
- Optimizing the stability of single-chain proteins by linker length and composition mutagenesis; C. R. Robinson, R. T. Sauer; *Proc Natl Acad Sci U S A* 1998; **95**: 5929-5934. <https://doi.org/10.1073/pnas.95.11.5929>
- Sortase-mediated protein ligation: a new method for protein engineering; H. Mao, S. A. Hart, A. Schink, B. A. Pollok; *J Am Chem Soc* 2004; **126**: 2670-2671. <https://doi.org/10.1021/ja039915e>
- Sortase-tag expressed protein ligation: combining protein purification and site-specific bioconjugation into a single step; R. Warden-Rothman, I. Caturegli, V. Popik, A. Tsourkas; *Anal Chem* 2013; **85**: 11090-11097. <https://doi.org/10.1021/ac402871k>

5. Examples

5.1. Kinetically Controlled Four-Fragment Ligation for the Convergent Chemical Synthesis of Proteins

The convergent synthesis of proteins by the chemical ligation of multiple peptide segments is achieved by a kinetically controlled four-fragment ligation. This was made possible by the controlled extension of the key Cys-peptide intermediate by ligation at either the N- or C-termini. For more information, please refer to the following publication:

Reference:

- Kinetically controlled ligation for the convergent chemical synthesis of proteins; D. Bang, B. L. Pentelute, S. B. Kent; *Angew Chem Int Ed Engl* 2006; **45**: 3985–3988.




Fig. 15: Kinetically controlled four-fragment ligation.

back to content ↑

5.2. Three-Fragment Synthesis of a 116mer with Desulfurization

The discovery of native chemical ligation, the regio- and chemoselective coupling of two unprotected peptide segments, facilitated the synthesis of polypeptides with more than 200 amino acids. However, this approach initially relied on the presence of at least one cysteine residue in the sequence at a convenient position. Hence, postligation-desulfurization protocols were developed that allowed for ligation at amino acid residues other than cysteine. Please refer to the following review for a comprehensive overview of “the development and recent progress on the chemical synthesis of peptides and proteins encompassing postligation-desulfurization at alanine, valine, lysine, threonine, leucine, proline, arginine, aspartic acid, glutamate, phenylalanine, glutamine, and tryptophan”.

Reference:

→ Postligation-Desulfurization: A General Approach for Chemical Protein Synthesis; J. Ma, J. Zeng, Q. Wan; **Protein Ligation and Total Synthesis II** 2015; (Ed.: L. Liu): Springer International Publishing, Cham, 57-101.




Fig. 16: Example of a synthesis of a peptide via ligation of three segments and thiazolidine as temporarily protected cysteine.

Notes

back to content ↑

Code of Conduct

As business activity of Iris Biotech GmbH impacts people's lives and health, it must be operated in ethical and correct manner and act with integrity and responsibility. To ensure high ethical standards and fair business practices, Iris Biotech GmbH applies an integrated policy known as its Code of Conduct.

In 2001 Iris Biotech GmbH was founded just at the beginning of the Biotech movement and the first remarkable breakthrough of biotech pharma products. Although the biotech field is rather young compared to other industries we believe on long-term business, a good partnership between our business partners and Iris Biotech GmbH and a good reputation. It is our duty as well as our responsibility to maintain and to extend this over the next generations – based on the principles of an honourable and prudent tradesman which based upon the concept of honourable entrepreneurship.

This Code of Conduct has been developed following the "Voluntary Guidelines for Manufacturers of Fine Chemical Intermediates and Active Ingredients" issued by AIME (Agrochemical & Intermediates Manufacturers in Europe) and the requirements of some of our business associates.

Iris Biotech GmbH commits to hold this Code of Conduct and to include and apply its principles in the management system and the company policies.

Ethics

Iris Biotech GmbH undertakes business in an ethical manner and acts with integrity. All corruption, extortion and embezzlement are prohibited. We do not pay or accept bribes or participate in other illegal inducements in business or government relationships. We conduct our business in compliance with all applicable anti-trust laws. Employees are encouraged to report concerns or illegal activities in the workplace, without threat of reprisal, intimidation or harassment.

Labour

Iris Biotech GmbH is committed to uphold the human rights of workers and to treat them with dignity and respect. Child labour, workplace harassment, discrimination, and harsh and inhumane treatment are prohibited. Iris Biotech GmbH respects the rights of the employees to associate freely, join or not join labour unions, seek representation and join workers' councils. Employees are paid and their working timetable is established according to applicable wage and labour laws. Employees are able to communicate openly with management regarding working conditions without threat of reprisal, intimidation or harassment.

General Policies

Contracts and Secrecy Agreements are binding and the confidential information received is only used for intended purposes. Clear management and organizational structures exist to provide efficient normal working and to address problems quickly. Know-how is protected and intellectual property is respected.

Health and Safety

Iris Biotech GmbH provides a safe and healthy working environment to the employees and protects them from overexposure to chemical and physical hazards. Products are produced, stored and shipped under the guidelines of the relevant chemical and safety legislation. Risks and emergency scenarios are identified and evaluated, and their possible impact is minimized by implementing emergency plans and written procedures. Safety information regarding hazardous materials is available to educate, train and protect workers from hazards. Preventive equipment and facilities maintenance is performed at suitable periods to reduce potential hazards. Employees are regularly trained in health and safety matters and are informed about product properties and risk classification when it is required.

Environment

Iris Biotech GmbH operates in an environmentally responsible and efficient manner, minimizing adverse impacts on the environment. Waste streams are managed to ensure a safe handling, movement, storage, recycling and reuse, before and after being generated. Systems to prevent and mitigate accidental spills and releases to the environment are in place. All required environmental permits and licenses are obtained and their operational and reporting requirements are complied with.

Production and Quality Management

A quality management system following the Good Distribution Practices (GDP rules) of Active Pharmaceutical Ingredients is established covering all the aspects of the worldwide distribution of products. Regular audits are performed to evaluate the efficiency and fulfilling of the quality system. Process controls to provide reproducible product quality are established. There are preventive maintenance procedures to ensure plant reliability and the lowest risk of failure. Staff is trained periodically about GMP and GDP rules. Procedures are established and installations are designed to avoid cross contamination. Batch and analytical records are kept for inspection and audit purposes for suitable periods according guidelines.

Research and Development

Research and development staff education is appropriate to their functional activity and they are trained to develop, optimize and scale-up the processes. Intellectual property is respected and know-how protected. Development of manufacturing processes reflects the principles of the Green Chemistry according to the American Chemical Society Green Chemistry Institute. Animal testing is not used unless alternatives are not scientifically valid or accepted by regulators. If animal testing is carried out, animals are treated so that pain and stress are minimized.

back to content ↑

Terms and Conditions of Sales

All orders placed by a buyer are accepted and all contracts are made subject to the terms which shall prevail and be effective notwithstanding any variations or additions contained in any order or other document submitted by the buyer. No modification of these terms shall be binding upon Iris Biotech GmbH unless made in writing by an authorised representative of Iris Biotech GmbH.

Placing of Orders

Every order made by the buyer shall be deemed an offer by the buyer to purchase products from Iris Biotech GmbH and will not be binding on Iris Biotech GmbH until a duly authorised representative of Iris Biotech GmbH has accepted the offer made by the buyer. Iris Biotech GmbH may accept orders from commercial, educational or government organisations, but not from private individuals and Iris Biotech GmbH reserves the right to insist on a written order and/or references from the buyer before proceeding.

There is no minimum order value. At the time of acceptance of an order Iris Biotech GmbH will either arrange prompt despatch from stock or the manufacture/acquisition of material to satisfy the order. In the event of the latter Iris Biotech GmbH will indicate an estimated delivery date. In addition to all its other rights Iris Biotech GmbH reserves the right to refuse the subsequent cancellation of the order if Iris Biotech GmbH expects to deliver the product on or prior to the estimated delivery date. Time shall not be of the essence in respect of delivery of the products. If Iris Biotech GmbH is unable to deliver any products by reason of any circumstances beyond its reasonable control („Force Majeure“) then the period for delivery shall be extended by the time lost due to such Force Majeure. Details of Force Majeure will be forwarded by Iris Biotech GmbH to the buyer as soon as reasonably practicable.

Prices, Quotations and Payments

Prices are subject to change. For the avoidance of doubt, the price advised by Iris Biotech GmbH at the time of the buyer placing the order shall supersede any previous price indications. The buyer must contact the local office of Iris Biotech GmbH before ordering if further information is required. Unless otherwise agreed by the buyer and Iris Biotech GmbH, the price shall be for delivery ex-works. In the event that the buyer requires delivery of the products otherwise than ex-works the buyer should contact the local office of Iris Biotech GmbH in order to detail its requirements. Iris Biotech GmbH shall, at its discretion, arrange the buyer's delivery requirements including, without limitation, transit insurance, the mode of transit (Iris Biotech GmbH reserves the right to vary the mode of transit if any regulations or other relevant considerations so require) and any special packaging requirements (including cylinders). For the avoidance of doubt all costs of delivery and packaging in accordance with the buyer's requests over and above that of delivery in standard packaging ex-works shall be for the buyer's account unless otherwise agreed by both parties. Incoterms 2020 shall apply. Any tax, duty or charge imposed by governmental authority or otherwise and any other applicable taxes, duties or charges shall be for the buyer's account. Iris Biotech GmbH may, on request and where possible, provide quotations for multiple packs or bulk quantities, and non-listed items. Irrespective of the type of request or means of response all quotations must be accepted by the buyer without condition and in writing before an order will be accepted by Iris Biotech GmbH. Unless agreed in writing on different terms, quotations are valid for 30 days from the date thereof. Payment terms are net 30 days from invoice date unless otherwise agreed in writing. Iris Biotech GmbH reserves the right to request advance payment at its discretion. For overseas transactions the buyer shall pay all the banking charges of Iris Biotech GmbH. The buyer shall not be entitled to withhold or set-off payment for the products for any reason whatsoever. Government/

Corporate Visa and MasterCard (and other such credit cards) may be accepted on approved accounts for payment of the products. Personal credit cards are not acceptable. Failure to comply with the terms of payment of Iris Biotech GmbH shall constitute default without reminder. In these circumstances Iris Biotech GmbH may (without prejudice to any other of its rights under these terms) charge interest to accrue on a daily basis at the rate of 2% per month from the date upon which payment falls due to the actual date of payment (such interest shall be paid monthly). If the buyer shall fail to fulfil the payment terms in respect of any invoice of Iris Biotech GmbH Iris Biotech GmbH may demand payment of all outstanding balances from the buyer whether due or not and/or cancel all outstanding orders and/or decline to make further deliveries or provision of services except upon receipt of cash or satisfactory securities. Until payment by the buyer in full of the price and any other monies due to Iris Biotech GmbH in respect of all other products or services supplied or agreed to be supplied by Iris Biotech GmbH to the buyer (including but without limitation any costs of delivery) the property in the products shall remain vested in Iris Biotech GmbH.

Shipping, Packaging and Returns

The buyer shall inspect goods immediately on receipt and inform Iris Biotech GmbH of any shortage or damage within five days. Quality problems must be notified within ten days of receipt. Goods must not be returned without prior written authorisation of Iris Biotech GmbH. Iris Biotech GmbH shall at its sole discretion replace the defective products (or parts thereof) free of charge or refund the price (or proportionate price) to buyer. Opened or damaged containers cannot be returned by the buyer without the written prior agreement of Iris Biotech GmbH. In the case of agreed damaged containers which cannot be so returned, the buyer assumes responsibility for the safe disposal of such containers in accordance with all applicable laws.

Product Quality, Specifications and Technical Information

Products are analysed in the Quality Control laboratories of Iris Biotech GmbH's production partners by methods and procedures which Iris Biotech GmbH considers appropriate. In the event of any dispute concerning reported discrepancies arising from the buyer's analytical results, determined by the buyer's own analytical procedures, Iris Biotech GmbH reserves the right to rely on the results of own analytical methods of Iris Biotech GmbH. Certificates of Analysis or Certificates of Conformity are available at the discretion of Iris Biotech GmbH for bulk orders but not normally for prepack orders. Iris Biotech GmbH reserves the right to make a charge for such certification. Specifications may change and reasonable variation from any value listed should not form the basis of a dispute. Any supply by Iris Biotech GmbH of bespoke or custom product for a buyer shall be to a specification agreed by both parties in writing. Technical information, provided orally, in writing, or by electronic means by or on behalf of Iris Biotech GmbH, including any descriptions, references, illustrations or diagrams in any catalogue or brochure, is provided for guidance purposes only and is subject to change.

Safety

All chemicals should be handled only by competent, suitably trained persons, familiar with laboratory procedures and potential chemical hazards. The burden of safe use of the products of Iris Biotech GmbH vests in the buyer. The buyer assumes all responsibility for warning his employees, and any persons who might reasonably be expected to come into contact with the products, of all risks to person and property in any way connected with the products and for instructing them in their safe handling and use. The buyer also assumes the responsibility for the safe disposal of all products in accordance with all applicable laws.

back to content ↑

Uses, Warranties and Liabilities

All products of Iris Biotech GmbH are intended for laboratory research purposes and unless otherwise stated on product labels, in the catalogue and product information sheet of Iris Biotech GmbH or in other literature furnished to the buyer, are not to be used for any other purposes, including but not limited to use as or as components in drugs for human or animal use, medical devices, cosmetics, food additives, household chemicals, agricultural or horticultural products or pesticides. Iris Biotech GmbH offers no warranty regarding the fitness of any product for a particular purpose and shall not be responsible for any loss or damage whatsoever arising there from. No warranty or representation is given by Iris Biotech GmbH that the products do not infringe any letters patent, trademarks, registered designs or other industrial rights. The buyer further warrants to Iris Biotech GmbH that any use of the products in the United States of America shall not result in the products becoming adulterated or misbranded within the meaning of the Federal Food, Drug and Cosmetic Act (or such equivalent legislation in force in the buyer's jurisdiction) and shall not be materials which may not, under sections 404, 505 or 512 of the Act, be introduced into interstate commerce. The buyer acknowledges that, since the products of Iris Biotech GmbH are intended for research purposes, they may not be on the Toxic Substances Control Act 1976 („TSCA“) inventory. The buyer warrants that it shall ensure that the products are approved for use under the TSCA (or such other equivalent legislation in force in the buyer's jurisdiction), if applicable. The buyer shall be responsible for complying with any legislation or regulations governing the use of the products and their importation into the country of destination (for the avoidance of doubt to include, without limitation, the TSCA and all its amendments, all EINECS, ELINCS and NONS regulations). If any licence or consent of any government or other authority shall be required for the acquisition, carriage or use of the products by the buyer the buyer shall obtain the same at its own expense and if necessary produce evidence of the same to Iris Biotech GmbH on demand. Failure to do so shall not entitle the buyer to withhold or delay payment. Any additional expenses or charges incurred by Iris Biotech GmbH resulting from such failure shall be for the buyer's account. Save for death or personal injury caused by negligence of Iris Biotech GmbH, sole obligation of Iris Biotech GmbH and buyer's exclusive remedy with respect to the products proved to the satisfaction of Iris Biotech GmbH to be defective or products incorrectly supplied shall be to accept the return of said products to Iris Biotech GmbH for refund of the actual purchase price paid by the buyer (or proportionate part thereof), or replacement of the defective product (or part thereof) with alternative product. Iris Biotech GmbH shall have no liability to the buyer under or arising directly or indirectly out of or otherwise in connection with the supply of products by Iris Biotech GmbH to the buyer and/or their re-sale or use by the buyer or for any product, process or services of the buyer which in any way comprises the product in contract tort (including negligence or breach of statutory duty) or otherwise for pure economic loss, loss of profit, business, reputation, depletion of brand, contracts, revenues or anticipated savings or for any special indirect or consequential damage or loss of any nature except as may otherwise be expressly provided for in these terms. All implied warranties, terms and representations in respect of the products (whether implied by statute or otherwise) are excluded to the fullest extent permitted by law. The buyer shall indemnify Iris Biotech GmbH for and against any and all losses, damages and expenses, including legal fees and other costs of defending any action, that Iris Biotech GmbH may sustain or incur as a result of any act or omission by the buyer, its officers, agents or employees, its successors or assignees, its customers or all other third parties, whether direct or indirect, in connection with the use of any product. For the avoidance of doubt and in the event that Iris Biotech GmbH supplies bespoke or custom product to the buyer's design or specification, this indemnity shall extend to include any claim by a third party that the manufacture of the product for the buyer or the use of the product by the buyer infringes the intellectual property rights of any third party.

General

Iris Biotech GmbH shall be entitled to assign or sub-contract all or any of its rights and obligations hereunder. The buyer shall not be entitled to assign, transfer, sub-contract or otherwise delegate any of its rights or obligations hereunder. Any delay or forbearance by Iris Biotech GmbH in exercising any right or remedy under these terms shall not constitute a waiver of such right or remedy. If any provision of these terms is held by any competent authority to be invalid or unenforceable in whole or in part the validity of the other provisions of these terms and the remainder of the provision in question shall not be affected. These terms shall be governed by German Law and the German Courts shall have exclusive jurisdiction for the hearing of any dispute between the parties save in relation to enforcement where the jurisdiction of the German Courts shall be non-exclusive.

Notes

Index

Product code	Product name	Page	Product code	Product name	Page
BAA3680	(Boc-L-Sec)2	32	FAA8600	Fmoc-D-Sec(Xan)-OH	33
HAA9350	(H-L-Sec-OH)2	32	FAA3160	Fmoc-D-Thz(Me2)-OH	15
AAA2015	Aloc-L-Cys(Trt)-OH	18	FAA1495	Fmoc-D-Thz-OH	14
BAA5410	Boc-D-Cys(Bzl)-OH	19	FAA3168	Fmoc-Dbz(<i>o</i> -Alloc)-OH	6
BAA5420	Boc-D-Cys(MBzl)-OH	19	FAA3167	Fmoc-Dbz(<i>o</i> -Boc)-OH	6
BAA5430	Boc-D-Cys(Mob)-OH	20	FAA3169	Fmoc-Dbz(<i>o</i> -NO ₂)-OH	7
BAA5000	Boc-D-Cys(Trt)-OH	21	FAA3165	Fmoc-Dbz-OH	6
BAA1170	Boc-D-Cys-OH	18	FDP1030	Fmoc-Gly-Gly-OH	39
BAA1186	Boc-D-Thz-OH	14	PSI1440	Fmoc-Gly-L-Cys[PSI(Dmp,H)pro]-OH	28
BAA1078	Boc-L-Cys(Acm)-OH	18	PYV1170	Fmoc-Gly-NHN=Pyv Resin	9
BAA1510	Boc-L-Cys(Acm,O)-OH	18	PSI1450	Fmoc-L-Ala-L-Cys[PSI(Dmp,H)pro]-OH	27
BAA1079	Boc-L-Cys(Bzl)-OH	19	PYV1100	Fmoc-L-Ala-NHN=Pyv Resin	8
BAA5510	Boc-L-Cys(Fm)-OH	19	PYV1110	Fmoc-L-Arg(Pbf)-NHN=Pyv Resin	8
BAA1080	Boc-L-Cys(MBzl)-OH	19	PSI1460	Fmoc-L-Asn(Trt)-L-Cys[PSI(Dmp,H)pro]-OH	27
BAA1081	Boc-L-Cys(Mob)-OH	19	PYV1120	Fmoc-L-Asn(Trt)-NHN=Pyv Resin	9
BAA1860	Boc-L-Cys(Npys)-OH	20	PYV1130	Fmoc-L-Asp(tBu)-NHN=Pyv Resin	9
BAA3140	Boc-L-Cys(Ph)-OH	20	PSI1470	Fmoc-L-Asp(tBu)-L-Cys[PSI(Dmp,H)pro]-OH	27
BAA6390	Boc-L-Cys(Phacm)-OH	20	FAA8460	Fmoc-L- <i>cis</i> -Hyp(NHBoc)-OH	36
BAA1082	Boc-L-Cys(tBu)-OH	20	FAA5150	Fmoc-L-Cys(Aapam)-OH	21
BAA1084	Boc-L-Cys(Trt)-OH	20	FAA4751	Fmoc-L-Cys(Ac-tBu)-OH*DCHA	22
BAA1083	Boc-L-Cys-OH	18	FAA1506	Fmoc-L-Cys(Acm)-OH	22
BAA3760	Boc-L-Sec(Mob)-OH	32	FAA7610	Fmoc-L-Cys(Allocam)-OH	22
BAA4830	Boc-L-Sec(Mob)-OPfp	32	FAA6270	Fmoc-L-Cys(Bzl)-OH	22
IAD1040	Boc-L-Ser[Fmoc-L-Cys(Trt)]-OH	21	FAA6940	Fmoc-L-Cys(Ddm)-OH	22
BAA4880	Boc-L-Sez-OH	33	FAA3190	Fmoc-L-Cys(Dpm)-OH	23
IAD2040	Boc-L-Thr[Fmoc-L-Cys(Trt)]-OH	21	FAA1714	Fmoc-L-Cys(MBzl)-OH	23
BAA1135	Boc-L-Thz-OH	14	FAA1030	Fmoc-L-Cys(Mmt)-OH	23
RL-1020	DTT (racemic)	35	FAA1715	Fmoc-L-Cys(Mob)-OH	24
FAA8450	Fmoc-AAHA(Boc)-OH (S)	36	FAA4155	Fmoc-L-Cys(Msbh)-OH	24
FAA6230	Fmoc-D-Cys(Acm)-OH	22	FAA8420	Fmoc-L-Cys(NDBF)-OH	24
FAA6950	Fmoc-D-Cys(Ddm)-OH	23	FAA1975	Fmoc-L-Cys(Npys)-OH	24
FAA5650	Fmoc-D-Cys(Dpm)-OH	23	FAA3970	Fmoc-L-Cys(<i>o</i> Nv)-OH	24
FAA1614	Fmoc-D-Cys(Mmt)-OH	23	FAA6910	Fmoc-L-Cys(Phacm)-OH	25
FAA8150	Fmoc-D-Cys(Msbh)-OH	24	FAA8495	Fmoc-L-Cys(SIT)-OH	25
FAA3710	Fmoc-D-Cys(Phacm)-OH	25	FAA1575	Fmoc-L-Cys(StBu)-OH	25
FAA1965	Fmoc-D-Cys(StBu)-OH	25	FAA1716	Fmoc-L-Cys(tBu)-OH	25
FAA1035	Fmoc-D-Cys(Trt)-OH	26	FAA4160	Fmoc-L-Cys(Thp)-OH	26
FAA1470	Fmoc-D-Cys-OH*H ₂ O	21	PSI1580	Fmoc-L-Cys(Trt)-L-Cys[PSI(Dmp,H)pro]-OH	27
FAA8710	Fmoc-D-Sec(Mob)-OH	33	PYV1140	Fmoc-L-Cys(Trt)-NHN=Pyv Resin	9

Product code	Product name	Page
FAA1040	Fmoc-L-Cys(Trt)-OH	26
FAA5670	Fmoc-L-Cys(Trt)-OMe	26
FAA1362	Fmoc-L-Cys-OH*H ₂ O	21
PSI1480	Fmoc-L-Gln(Trt)-L-Cys[PSI(Dmp,H)pro]-OH	27
PYV1160	Fmoc-L-Gln(Trt)-NHN=Pyv Resin	9
PSI1490	Fmoc-L-Glu(tBu)-L-Cys[PSI(Dmp,H)pro]-OH	28
PYV1150	Fmoc-L-Glu(tBu)-NHN=Pyv Resin	9
PYV1180	Fmoc-L-His(Trt)-NHN=Pyv Resin	10
PSI1500	Fmoc-L-Ile-L-Cys[PSI(Dmp,H)pro]-OH	28
PYV1190	Fmoc-L-Ile-NHN=Pyv Resin	10
PSI1510	Fmoc-L-Leu-L-Cys[PSI(Dmp,H)pro]-OH	28
PYV1200	Fmoc-L-Leu-NHN=Pyv Resin	10
PSI1520	Fmoc-L-Lys(Boc)-L-Cys[PSI(Dmp,H)pro]-OH	28
PYV1210	Fmoc-L-Lys(Boc)-NHN=Pyv Resin	10
PYV1220	Fmoc-L-Met-NHN=Pyv Resin	10
PSI1530	Fmoc-L-Phe-L-Cys[PSI(Dmp,H)pro]-OH	28
PYV1230	Fmoc-L-Phe-NHN=Pyv Resin	10
PYV1240	Fmoc-L-Pro-NHN=Pyv Resin	11
FAA8705	Fmoc-L-Sec(Mob)-OH	32
FAA8760	Fmoc-L-Sec(Mob)-OPfp	32
FAA8465	Fmoc-L-Sec(Xan)-OH	33
FAA4205	Fmoc-L-Selenomethionine	33
PYV1250	Fmoc-L-Ser(tBu)-NHN=Pyv Resin	11
FAA8860	Fmoc-L-Sez-OH	34
PSI1550	Fmoc-L-Thr(tBu)-L-Cys[PSI(Dmp,H)pro]-OH	29
PYV1260	Fmoc-L-Thr(tBu)-NHN=Pyv Resin	11
FAA1437	Fmoc-L-Thz(Me2)-OH	14
FAA1427	Fmoc-L-Thz-OH	14
FAA8455	Fmoc-L-trans-Hyp(NHBoc)-OH	36
PYV1270	Fmoc-L-Trp(Boc)-NHN=Pyv Resin	11
PSI1560	Fmoc-L-Tyr(tBu)-L-Cys[PSI(Dmp,H)pro]-OH	29
PYV1280	Fmoc-L-Tyr(tBu)-NHN=Pyv Resin	11
PSI1570	Fmoc-L-Val-L-Cys[PSI(Dmp,H)pro]-OH	29
PYV1290	Fmoc-L-Val-NHN=Pyv Resin	11
FAA3166	Fmoc-MeDbz-OH	6
PYV1000	Fmoc-NHN=Pyv Resin	8
BR-5279	Fmoc-NHNH-2CT Resin	8
BR-5280	Fmoc-NHNH-Trt-PEG Resin	8
HAA2870	H-(Gly)3-Lys(N ₂)-OH*HCl	39

Product code	Product name	Page
HAA6110	H-D-Cys(Bzl)-OH	15
HAA3500	H-D-Cys(Mmt)-OH	16
HAA2100	H-D-Cys(Trt)-OtBu*HCl	17
HAA1017	H-D-Cys-OH*HCl*H ₂ O	15
HAA6070	H-L-Cys(Acm)-OH*HCl	15
HAA1574	H-L-Cys(Bzl)-OH	15
HAA6080	H-L-Cys(Bzl)-OMe*HCl	16
HAA9320	H-L-Cys(DMNB)-OH	17
HAA6090	H-L-Cys(MBzl)-OH	16
HAA9270	H-L-Cys(MDNPE)-OH	16
HAA6100	H-L-Cys(Mob)-OH	16
HAA3510	H-L-Cys(Npys)-OH*HCl	16
HAA6140	H-L-Cys(StBu)-OH	17
HAA6150	H-L-Cys(tBu)-OH*HCl	17
HAA6160	H-L-Cys(Trt)-OH	17
HAA1995	H-L-Cys(Trt)-OtBu*HCl	17
HAA9255	H-L-Sec(DMNB)-OH*TFA	31
HAA9230	H-L-Sec(MDNPE)*TFA	31
HAA9360	H-L-Sec(MDNPE)-OH	31
HAA2860	N ₃ -Gly-Gly-Gly-Gly-Gly-OH	38
HAA2840	N ₃ -Gly-Gly-Gly-OH	38
HAA2850	N ₃ -Gly-Gly-OH*DCHA	38
HAA2810	N ₃ -L-Cys(Trt)-OH*CHA	18
RL-4080	NTCB	37
SAA1110	Smoc-L-Cys(Trt)-OH	29
PAA2000	tBu-SS-Photo(Fmoc)-Gly-OH	29
LS-3405	TCEP*HCl	35
TAA1508	Trt-L-Cys(Trt)-OH*DEA	26
ZAA1161	Z-L-Cys(Bzl)-OH	26
ZAA1310	Z-L-Cys(Trt)-OH	27

back to content ↑



Get in Contact

**Iris**
Biotech

Iris Biotech GmbH
Adalbert-Zoellner-Str. 1
95615 Marktredwitz
Germany

+49 (0) 9231 97121-0
+49 (0) 9231 97121-99
info@iris-biotech.de
www.iris-biotech.de

Distribution Partners

The list contains the current distributors of Iris Biotech in different regions of the world. The latest list of distribution partners and contact details is available at: www.iris-biotech.de/distribution-partner

China:

Chengdu Yoo Technology Co., Ltd.

Japan:

BizCom Japan, Inc.

Shigematsu & Co., Ltd

Cosmo Bio Co., Ltd.

USA & Canada:

Peptide Solutions LLC

India, Bangladesh, Oman, Sri Lanka, United Arab Emirates:

Sumit Biosciences Pvt Ltd.

Singapore, Thailand, Malaysia, Indonesia, Vietnam, Philippines, Brunei, Burma (Myanmar), Laos, Cambodia, Timor-Leste:

SciClix Pte. Ltd.

Empowering Peptide Innovation