




Iris
Biotech



CLICK CHEMISTRY



Version: IB1_4

Empowering Peptide Innovation

With this guiding theme in mind, Iris Biotech's mission is to support researchers by supplying

- innovative technologies,
- rare compounds,
- as well as a broad portfolio on standard consumables,

available in flexible quantities from small scale to bulk quantities. To fulfill our dedication "Empowering Peptide Innovation", we are attending various conferences, symposia, and exhibitions each year. This allows us to remain in direct contact with scientists all over the world, both from academia and industry, to exchange knowledge, and to gather new ideas to tackle your current challenges.

Guided by our dedication to provide

- competent service,
- as well as novel substances and
- latest technologies,

Iris Biotech is your trusted partner for the world of peptides, while having strong expertise in associated disciplines. Thus, our portfolio comprises reagents and tools for the synthesis and modification of peptides, e.g., amino acids, resins and solvents but also for related technologies such as drug delivery, linkerology® and life sciences.

Owed to the growing demand for tailor-made compounds, our portfolio is fine-tuned by our custom synthesis service at Iris Biotech Laboratories. Our skilled scientists offer profound expertise in

- *de novo* route development,
- upscaling towards larger scale production,
- as well as synthesis optimization for increased efficiency.

Examples are the synthesis of rare chiral building blocks, unnatural amino acid derivatives, sophisticated orthogonal protecting groups, heterocycles, building blocks for nucleotides, PEGs and PEG-analogs as well as specific linkers for controlled drug delivery and release.



Amino Acids



Building Blocks



Life Sciences



Drug Delivery



Reagents



Resins



Linkerology®



Click Chemistry

Portfolio Overview

Peptide Synthesis and Modification

(Protected) Amino Acids

Standards such as Fmoc-D/L-AAA and Boc-D/L-AAA, Smoc amino acids for peptide synthesis in water, variety of protecting groups (e.g., Pbf, Trt, ^tBu, Bzl, Acn, Mob, SIT, Phacn, Allocam, Mmt), unusual amino acids, fluorinated derivatives, substituted prolines, arginine analogs

Building Blocks

Amino alcohols, amino aldehydes, diamines and hydrazines, (pseudoproline) dipeptides, polyamines and spermines, fatty acid derivatives, peptide nucleic acids (PNAs)

Reagents

Coupling reagents, solvents and scavengers, protecting groups

Resins

Preloaded resins (e.g., based on Trityl, TCP, TentaGel, Methoxybenzhydryl, Merrifield, PAM, Rink, Wang), scavenger resins, hydrazone resins, poly(acrylamide) resins, Cyclover

Linkerology® and Drug Delivery

Linkers for Solid Phase Peptide Synthesis

Cleavable Linkers

Val-Ala-based, Val-Cit-based, disulfide-based, Dde-helping hands, pH-sensitive linkers

Photo-Activatable Linkers

Functionalized Linkers

Clickable linkers, trifunctional linkers, linkers with maleimide function, cross-linkers, selective N-term acylation and biotinylation, 5HP2O

PROTACs

Ligands, linkers & modules

Fullerenes, Poly(2-oxazolines), Dextrans & Plant-Derived Cholesterol

Superparamagnetic Iron Oxide Nanoparticles

Poly-Amino Acids

Poly-Arg, Poly-Glu, Poly-Lys, Poly-Orn, Poly-Sar

PEGylation

Branched PEGylating reagents, (amino-)PEG-acids, PEG-amines & hydrazides & guanidines, reagents for Click-conjugation, Biotin-PEG-reagents, PEG-thiols, PEG-maleimides, other PEGylating reagents

Life Sciences

Biotinylation Reagents

Carbohydrates

Galactose, Glucose, Mannose, Xylose and others

Drug Metabolites

Peptides

Substrates & Inhibitors

E.g., protein kinase inhibitors, substrates for fusion (Halo/Snap/Clip)-tagged proteins

Natural Products

Dyes and Fluorescent Labels

E.g., ICG, AMC, DAPI

Maillard & Amadori Reaction Products

Large portfolio of derivatives useful as standards for food, pharma and cosmetics industry

Vitamins

Custom Synthesis

Your project requires a compound not listed in our portfolio?
Get in contact and inquire about our custom synthesis capabilities.

Our experienced scientists are excited to accept your synthetic challenge!

In such cases, your request undergoes the following stages:



Step-by-Step Analysis

- Customer's demands

Process Evaluation

- Detailed literature review
- Synthetic possibilities

Strategy Development

- Protocol development
- Method development and validation
- Customized synthesis

Quality Consistency

- Identity confirmation
- Purity verification

Our Service Promise

All our services are based on high standards, transparency & documentation, trust, honesty & confidentiality, as well as the required know-how.

High Standards

- Values: sustainability & responsibility
- State-of-the-art equipment & latest technologies
- High quality standards
- Qualified suppliers & regular audits

Transparency & Documentation

- Talk to our specialists – customer care
- Certificates of analysis & origin
- Impurity profiling
- Safety data sheets
- Analytical and process reports

Trust, Honesty & Confidentiality

- Intergenerational business valuing partnerships
- Meeting the customer's expectations
- Integrity towards our customers

Our Know-How

- One-step reactions & complex multi-step synthesis
- Scalability from mg to kg quantities
- Route scouting





Table of Contents


1. The Click Reaction	1
1.1. The 1 st Generation Click Reaction: CuAAC	1
1.2. Catalyst-free Click Reactions: 2 nd and 3 rd Generation Click Chemistry	3
1.2.1. 2 nd Generation: Strain-Promoted Azide-Alkyne Cycloadditions (SPAAC)	3
1.2.2. 3 rd Generation: Inverse Electron-Demand Diels-Alder (IEDDA) Reactions	17
1.2.3. Custom Synthesis of DBCO, Tetrazine and TCO Derivatives	26
2. Amino Acid Derivatives and Related Building Blocks for Click Chemistry	28
2.1. Recombinant Incorporation of Amino Acids into Proteins	28
2.2. Peptide Synthesis with Azido and Alkyne Amino Acids	31
2.3. Azido Amino Acids and Related Derivatives	32
2.4. Alkyne Amino Acids and Related Derivatives	66
3. Spermines and Amines for Click Chemistry	83
4. Click Reagents for Drug Delivery	85
4.1. Principles of Polymer Therapeutics	85
4.2. PEGylation Improves Drug Delivery and Pharmacokinetics	87
4.3. Azido-PEG Derivatives for Click Chemistry	90
4.4. Alkyne-PEG Derivatives for Click Chemistry	106
4.5. Poly(Amino Acids) and Poly(Oxazolines) for Click Chemistry	114

5. Click Chemistry Tools for Proteomics	117
5.1. Indocyanine Green Dyes for Click Chemistry	117
5.2. Clickable Linkers for Selective Protein Labeling	119
6. Carbohydrates for Click Chemistry	126
6.1. Galactose Derivatives	127
6.2. Glucose Derivatives	128
6.3. Mannose Derivatives	129
6.4. Lactose Derivatives	130
7. Proteolysis Targeting Chimeras (PROTACs®)	132
Index	139
Code of Conduct	147
Terms and Conditions of Sales	149

1. The Click Reaction

1.1. The 1st Generation Click Reaction: CuAAC

Alkynes and azides can undergo a Cu(I)-catalyzed azide-alkyne 1,3-dipolar cycloaddition (CuAAC) to afford 1,4-disubstituted 1,2,3-triazoles. Developed by K. Barry Sharpless and Morton Meldal, this type of chemical transformation was quickly dubbed “Click chemistry”. It has since become a widely used reaction that is orthogonal to many other types of chemical transformations and is used in various kinds of applications. Due to its high thermodynamic driving force, which is usually greater than 20 kcal/mol, the Click reaction rapidly proceeds to completion in almost all cases. Moreover, while the thermal Huisgen 1,3-dipolar cycloaddition affords a mixture of both the 1,4-disubstituted and the 1,5-disubstituted regioisomers, the CuAAC is highly selective for the 1,4-disubstituted isomer only (Fig. 1). Worth noting is the fact that ruthenium is also able to catalyze a 1,3-dipolar cycloaddition between an azide and an alkyne affording the 1,5-disubstituted product instead.



Cycloaddition reactions such as the [3+2] azide-alkyne and the [4+2] Diels-Alder reaction, have become common conjugation techniques. Applications range from imaging and drug design to the development of sensors, thereby covering such diverse fields as chemical biology, material science, surface and polymer chemistry.

Tris(benzyltriazolylmethyl)amine (TBTA; [RL-2010 on page 82](#)) is stabilizing copper(I) towards oxidation in solution by forming a complex and effectively catalyzes quantitative and regioselective Click cycloaddition reactions in a variety of aqueous and organic solvents. Among scientists, CuAAC has found widespread use as a biochemical tool for the site-specific labeling of peptides, proteins, and other biomolecules.

THPTA ([RL-2210 on page 82](#)) is a water-soluble alternative to TBTA ([RL-2010 on page 82](#)) and a highly efficient ligand for Click chemistry in partially organic and particularly in completely aqueous reactions. The benefits of a completely aqueous reaction include the biological labelling of live cells or the labeling of proteins without the concern of denaturing secondary structures. THPTA complexes Cu(I) and thus blocks its bioavailability. This mitigates potentially toxic effects while maintaining the catalytic effectiveness in Click conjugations. Successful Click reactions with oligonucleotides can be found in many publications.

A variety of azido and alkyne building blocks are available from Iris Biotech. Some of those compounds can be incorporated into peptides and proteins by recombinant syntheses, particularly by non-neutral protein translation using the amber-suppression-based orthogonal system, while others are suitable for solid phase peptide synthesis. The presence of an azido or alkyne function at a particular position of a peptide sequence opens up the possibility for the site-selective conjugation of other biomolecules (e.g., carbohydrates), labels or APIs.

References:

- A Stepwise Huisgen Cycloaddition Process: Copper(I)-Catalyzed Regioselective “Ligation” of Azides, Terminal Alkynes; V. V. Rostovtsev, L. G. Green, V. V. Fokin, K. B. Sharpless; **Angew. Chem. Int. Ed.** 2002; **41**: 2596-2599. [https://doi.org/10.1002/1521-3773\(20020715\)41:14<2596::Aid-anie2596>3.0.Co;2-4](https://doi.org/10.1002/1521-3773(20020715)41:14<2596::Aid-anie2596>3.0.Co;2-4)
- Peptidotriazoles on solid phase: [1,2,3]-triazoles by regiospecific copper(i)-catalyzed 1,3-dipolar cycloadditions of terminal alkynes to azides; C. W. Tornøe, C. Christensen, M. Meldal; **J Org Chem** 2002; **67**: 3057-64. <https://doi.org/10.1021/jo011148j>
- Click Chemistry: Diverse Chemical Function from a Few Good Reactions; H. C. Kolb, M. G. Finn, K. B. Sharpless; **Angew. Chem. Int. Ed.** 2001; **40**: 2004-2021. [https://doi.org/10.1002/1521-3773\(20010601\)40:11<2004::aid-anie2004>3.0.co;2-5](https://doi.org/10.1002/1521-3773(20010601)40:11<2004::aid-anie2004>3.0.co;2-5)
- The growing impact of click chemistry on drug discovery; H. C. Kolb, K. B. Sharpless; **Drug Discov Today** 2003; **8**: 1128-37. [https://doi.org/10.1016/S1359-6446\(03\)02933-7](https://doi.org/10.1016/S1359-6446(03)02933-7)
- Polytriazoles as copper(I)-stabilizing ligands in catalysis; T. R. Chan, R. Hilgraf, K. B. Sharpless, V. V. Fokin; **Org Lett** 2004; **6**: 2853-5. <https://doi.org/10.1021/ol0493094>
- CuI-Catalyzed Alkyne-Azide “Click” Cycloadditions from a Mechanistic and Synthetic Perspective; V. D. Bock, H. Hiemstra, J. H. van Maarseveen; **Eur. J. Org. Chem.** 2006; **1**: 51-68. <https://doi.org/10.1002/ejoc.200500483>
- A3-type star polymers via click chemistry; O. Altintas, B. Yankul, G. Hizal, U. Tunca; **J. Polym. Sci. A. Polym. Chem.** 2006; **44**: 6458-6465. <https://doi.org/10.1002/pola.21728>
- Preparation of alumina supported copper nanoparticles and their application in the synthesis of 1,2,3-triazoles; M. L. Kantam, V. S. Jaya, B. Sreedhar, M. M. Rao, B. M. Choudary; **J. Mol. Catal. A Chem.** 2006; **256**: 273-277. <https://doi.org/10.1016/j.molcata.2006.04.054>
- A rapid and versatile method to label receptor ligands using „click“ chemistry: Validation with the muscarinic M1 antagonist pirenzepine; D. Bonnet, B. Ilien, J. L. Galzi, S. Riche, C. Antheaune, M. Hibert; **Bioconjug Chem** 2006; **17**: 1618-23. <https://doi.org/10.1021/bc060140j>
- Alkyne-azide click reaction catalyzed by metallic copper under ultrasound; P. Cintas, A. Barge, S. Tagliapietra, L. Boffa, G. Cravotto; **Nat Protoc** 2010; **5**: 607-16. <https://doi.org/10.1038/nprot.2010.1>
- Synthesis of a DOTA--biotin conjugate for radionuclide chelation via Cu-free click chemistry; M. K. Schultz, S. G. Parameswarappa, F. C. Pigge; **Org Lett** 2010; **12**: 2398-401. <https://doi.org/10.1021/ol100774p>
- „Click“ cycloaddition catalysts: copper(I) and copper(II) tris(triazolylmethyl)amine complexes; P. S. Donnelly, S. D. Zanatta, S. C. Zammit, J. M. White, S. J. Williams; **Chem Commun** 2008: 2459-61. <https://doi.org/10.1039/b719724a>
- Click Chemistry and Radiochemistry: The First 10 Years; J. P. Meyer, P. Adumeau, J. S. Lewis, B. M. Zeglis; **Bioconjug Chem** 2016; **27**: 2791-2807. <https://doi.org/10.1021/acs.bioconjchem.6b00561>

1.2. Catalyst-free Click Reactions: 2nd and 3rd Generation Click Chemistry

Introduced in 2002, the copper-catalyzed variant of the azide-alkyne cycloaddition (CuAAC) reaction has found broad applicability in various fields and is as such currently the most widely used conjugation technique. The presence of copper, however, limits *in vivo* applications of this reaction for several reasons:


- High cell toxicity
- Undesired oxidation of proteins and
- Inhibition of luminescence properties of nanocrystals

To allow for fast and efficient *in vivo* conjugations, new methodologies were developed that do not require the use of a metal catalyst while still making use of bioorthogonal functional groups. The most commonly used approaches can be classified into two categories.

1.2.1. 2nd Generation: Strain-Promoted Azide-Alkyne Cycloadditions (SPAAC)

As early as 1961, Wittig and Krebs noted the propensity of cyclooctyne to strongly react with phenyl azide via a 1,3-dipolar cycloaddition, forming a triazole product. This finding stood in stark contrast to previous research that found slow kinetics for Huisgen 1,3-dipolar cycloadditions of azides with unstrained, linear alkynes. The latter reaction can be drastically accelerated by copper catalysis. The use of this metal, however, is linked with several drawbacks as noted above.

This property of cyclooctynes was exploited by Bertozzi *et al.* in the design of SPAAC reagents for bio-orthogonal couplings to azide-bearing biomolecules in live cells or organisms such as *C. elegans*, zebrafish or mice. By modifying the cyclooctyne core structure of SPAAC reagents with heteroatoms, fluorine substituents and fused rings, key properties such as cycloaddition kinetics, stability, solubility, and pharmacokinetics could be optimized.



In Fig. 4, various strained cyclooctynes and cyclononynes are depicted with their corresponding reactivities, as determined by their reaction with benzyl azide as a model compound. In general, the presence of atoms with high electronegativity next to the alkyne function, i.e. good σ -acceptors, leads to increased reactivity. A higher reactivity also correlates with increased ring strain, as exemplified by dibenzo-fused cyclooctynes (DiBO, DBCO) and bicyclo[6.1.0]non-4-yne (BCN).

A relatively new addition to this ensemble is 4,8-diazacyclononyne (DACN). While exhibiting a reactivity twice as high as OCT, DACN is also more hydrophilic than most cyclooctynes, highly stable (both thermal and chemical stability), and highly selective towards ynophiles. Additionally, the two endocyclic nitrogens in DACN may serve as additional attachment points for further conjugation, rendering the compound functionally versatile.

CliCr® - an innovative Click reagent

CliCr® is based on the small-molecule TMT-H-Sulfoximine (TMT-HSI). It constitutes a superior class of reagents for metal-free click strain promoted cycloaddition-conjugation with azides. The imine in the 7-membered CliCr® ring can be conveniently functionalized with a variety of linkers, e.g., *via* acylation, sulfonylation, N-alkylation, or carbamoylation. CliCr® reagents can be used in diverse applications, including the construction of antibody-drug conjugates (ADCs), small molecule-drug conjugates, oligonucleotide conjugates as well as for diagnostic labelling of a variety of nanoparticles and other agents. CliCr® can also be used in the conjugation of larger proteins (including optionality for native peptide release) or for *ex vivo* cell modification (e.g., glycocalyx).




Fig. 3: Chemical structure of the CliCr® base compound and its derivatization possibilities.




Fig. 4: Reaction rate constants for different strained cycloheptyne (CliCr), cyclooctyne and cyclononyne (DACN) derivatives.

[↑ back to content](#)

Key benefits:

More attractive CoG via faster click reactions:

Shorter reaction times than all other marketed copper-free click reagents, providing greater chemical yield.

Generation of (biodegradable) bioconjugates:

Highly stable reagents yielding (biodegradable) linked drug products.


Greater variety of click reactions:

Due to a larger variation of linkers that can easily be attached to the seven-membered ring.

Broad applicability:

Next to straightforward bioconjugation, a plethora of additional important applications in biochemical, aqueous environments is envisioned, such as surface plasmon resonance (SPR) applications and conjugation of chelator moieties for radio-active isotope incorporation in theragnostic applications.

CliCr® is shown to be > 5 times more reactive than BCN. The reaction progress of 5 mM CliCr® (blue line) or BCN-OH (red line), respectively, with 1.3 eq. of benzylazide in CDCl₃ at room temperature was monitored by MS. The reaction conversion was measured based on the increase of the triazole signals (see <https://doi.org/10.1039/d0sc03477k>).







The increased hydrophilicity as observed for antibody CliCr®-derived conjugates as compared to DBCO-linked constructs, translates into superior *in vivo* biodistribution, namely less liver and spleen uptake. Alike, also significantly faster and more cell labeling was detected upon using CliCr® as compared to DBCO.

Within our portfolio, we offer a selection of CliCr® derivatives. The CliCr® derivatives can be clicked to azide compounds with high efficiency and excellent stability. For further derivatives, large scale production or GMP grade, please get in contact!

CliCr® is provided under an intellectual property license from Cristal Therapeutics. The trademark CliCr® is the property of Cristal Therapeutics. For information on purchasing a license of CliCr® reagents, contact Cristal Therapeutics via Oxfordlaan 55, 6229 EV Maastricht (The Netherlands) or via info@crystaltherapeutics.com.

CliCr® Products

	Product details
RL-4180 CliCr® base compound TMTH-Sulfoximine CAS-No. 2408481-82-1 Formula C ₁₀ H ₁₇ NOS Mol. weight 199,31 g/mol	 
RL-4190 CliCr®-beta-Ala-NH ₂ *TFA TMTH-sulfoximine beta-alanine amide TFA salt CAS-No. 3038495-92-7 net Formula C ₁₃ H ₂₂ N ₂ O ₂ S*CF ₃ COOH Mol. weight 270,39*114,02 g/mol	 
RL-4200 CliCr®-Suc TMTH-sulfoximine succinic acid CAS-No. 2479971-29-2 Formula C ₁₄ H ₂₁ NO ₄ S Mol. weight 299,39 g/mol	 
RL-4330 CliCr®-OSu TMTH-sulfoximine succinimidyl ester CAS-No. 2408481-89-8 Formula C ₁₅ H ₂₀ N ₂ O ₅ S Mol. weight 340,39 g/mol	 
RL-4205 CliCr® Feasibility Service This paid service package consists of 20 hours of consulting and advice by the CliCr® experts of Cristal Therapeutics for the period of 1 year, paid in advance, no refunds. This Service Package will help to maximize your results when testing CliCr® - the innovative click reagent. For more information, please get in contact!	 

[↑ back to content](#)

The Click Reaction

Amino Acid Derivatives and Related Building Blocks for Click Chemistry

Spermines and Amines for Click Chemistry

Click Reagents for Drug Delivery

Click Chemistry Tools for Proteomics

Carbohydrates for Click Chemistry

Proteolysis Targeting Chimeras (PROTACs®)

Index

References:

- *TMTHSI, a superior 7-membered ring alkyne containing reagent for strain-promoted azide-alkyne cycloaddition reactions*; J. Weterings, C. J. F. Rijcken, H. Veldhuis, T. Meulemans, D. Hadavi, M. Timmers, M. Honing, H. Ippel, R. M. J. Liskamp; **Chem. Sci.** 2020; **11**: 9011-9016. <https://doi.org/10.1039/d0sc03477k>
- *Exploring the Chemical Properties and Medicinal Applications of Tetramethylthiocycloheptyne Sulfoximine Used in Strain-Promoted Azide-Alkyne Cycloaddition Reactions*; M. Timmers, A. Kipper, R. Frey, S. Notermans, M. Voievudskyi, C. Wilson, N. Hentzen, M. Ringle, C. Bovino, B. Stump, C. J. F. Rijcken, T. Vermonden, I. Dijkgraaf, R. Liskamp; **Pharmaceuticals** 2023; **16**: 1155. <https://doi.org/10.3390/ph16081155>
- *Specific N-terminal attachment of TMTHSI linkers to native peptides and proteins for strain-promoted azide alkyne cycloaddition*; M. Timmers, W. Peeters, N. J. Hauwert, C. J. F. Rijcken, T. Vermonden, I. Dijkgraaf, R. M. J. Liskamp; **Chem. Commun.** 2023; <https://doi.org/10.1039/d3cc03397j>
- *Versatile click linker enabling native peptide release from nanocarriers upon redox trigger*; E. R. Hebels, S. Dietl, M. Timmers, J. Hak, A. van den Dikkenberg, C. J. F. Rijcken, W. E. Hennink, R. M. J. Liskamp, T. Vermonden; **Bioconj. Chem.** 2023, **34**: 2375-2386. <https://doi.org/10.1021/acs.bioconjchem.3c00484>











You want more details about CliCr®?

Watch the recording of our online workshop!



Products with DBCO

		Product details
RL-2870	ICG-DBCO	
Indocyanine green dibenzoazacyclooctyne		
CAS-No.	3024705-52-7	
Formula	C ₆₃ H ₆₄ N ₄ O ₅ S	
Mol. weight	989,27 g/mol	
		
PEG7465	Me-PEG(12)-DBCO	
Methyl-12(ethylene glycol)-amido-dibenzoazacyclooctyne		
Formula	C ₄₅ H ₆₈ N ₂ O ₁₄	
Mol. weight	861,04 g/mol	
		

		Product details
RL-2120	DBCO-NH₂	
Dibenzocyclooctyne-amine		
CAS-No.	1255942-06-3	
Formula	C ₁₈ H ₁₆ N ₂ O	
Mol. weight	276,33 g/mol	
		
RL-2430	DBCO-COOH	
Dibenzoazacyclooctyne-carboxylic acid		
CAS-No.	1425485-72-8	
Formula	C ₂₁ H ₁₉ NO ₃	
Mol. weight	333,38 g/mol	
		
RL-2440	DBCO-NHS	
Dibenzoazacyclooctyne-carboxylic acid succinimidyl ester		
CAS-No.	1384870-47-6	
Formula	C ₂₅ H ₂₂ N ₃ O ₅	
Mol. weight	430,45 g/mol	
		
RL-2490	DBCO-mal	
Dibenzoazacyclooctyne-maleimide		
CAS-No.	1395786-30-7	
Formula	C ₂₅ H ₂₁ N ₃ O ₄	
Mol. weight	427,45 g/mol	
		
RL-4020	DBCO-C6-Alkyne	
N-(propargylamidoadipoyl)-dibenzoazacyclooctyne		
Formula	C ₂₄ H ₂₂ N ₂ O ₂	
Mol. weight	370,45 g/mol	
		
ADC1620	DBCO-cyclobutane-1,1-dicarboxamide-Ala-PAB	
dibenzoazacyclooctyne-cyclobutane-1,1-dicarboxamide-alanyl-(4-aminobenzyl alcohol)		
CAS-No.	2576471-46-8	
Formula	C ₃₄ H ₃₄ N ₄ O ₅	
Mol. weight	578,66 g/mol	
		


ADC1630 DBCO-cyclobutane-1,1-dicarboxamide-Ala-PAB-PNP

dibenzoazacyclooctyne-cyclobutane-1,1-dicarboxamide-alanyl-(4-aminobenzyl)-(4-nitrophenyl)-carbonate

CAS-No. 2576471-43-5

Formula $C_{41}H_{37}N_5O_9$

Mol. weight 743,76 g/mol


**ADC1520 DBCO-cyclobutane-1,1-dicarboxamide-Cit-PAB**

dibenzoazacyclooctyne-cyclobutane-1,1-dicarboxamide-citrullyl-(4-aminobenzyl alcohol)

CAS-No. 2576471-51-5

Formula $C_{37}H_{40}N_6O_6$

Mol. weight 664,75 g/mol


**ADC1530 DBCO-cyclobutane-1,1-dicarboxamide-Cit-PAB-PNP**

dibenzoazacyclooctyne-cyclobutane-1,1-dicarboxamide-citrullyl-(4-aminobenzyl)-(4-nitrophenyl)-carbonate

CAS-No. 2576471-34-4

Formula $C_{44}H_{43}N_{10}O_{10}$

Mol. weight 829,85 g/mol

**RL-2480 DBCO-PEG(3)-BisSulfonThiol-Linker**

Dibenzoazacyclooctyne-PEG(3)-BisSulfon-Thiol-Linker

Formula $C_{59}H_{69}N_3O_{14}S_3$

Mol. weight 1140,39 g/mol


**RL-3670 Halo-DBCO**

N-[2-[2-[(6-chlorohexyl)oxy]ethoxy]ethyl]-gamma-oxo-dibenz[b,f]azocine-5(6H)-butanamide

CAS-No. 1808119-16-5

Formula $C_{29}H_{35}ClN_2O_4$

Mol. weight 511,06 g/mol


**LS-4270 Biotin-DBCO**







(3aS,4S,6aR)-N-[3-(11,12-Didehydrodibenz[b,f]azocin-5(6H)-yl)-3-oxopropyl]hexahydro-2-oxo-1H-thieno[3,4-d]imidazole-4-pentanamide

CAS-No. 1418217-95-4

Formula $C_{28}H_{30}N_4O_3S$

Mol. weight 502,63 g/mol



			Product details
RL-2520	Biotin-PEG(4)-DBCO		
Dibenzoazacyclooctyne-tetra(ethylene glycol)-biotin			
CAS-No.	1255942-07-4		
Formula	C ₃₉ H ₅₁ N ₅ O ₈ S		
Mol. weight	749,92 g/mol		
PEG8140	Biotin-PEG(4)-Dde-DBCO		
N-(15-(4,4-dimethyl-2,6-dioxocyclohexylidene)-19-oxo-19-(azadibenzocyclooctyn-1-yl)-3,6,9,12-tetraoxa-16-azanonadecyl)-5-(((3aS,4S,6aR)-2-oxohexahydro-1H-thieno[3,4-d]imidazol-4-yl)pentanamide			
CAS-No.	1807512-43-1		
Formula	C ₄₇ H ₆₁ N ₅ O ₉ S		
Mol. weight	872,08 g/mol		
PEG8120	Biotin-PEG(4)-SS-DBCO		
N-(2-((3-(3-(azadibenzocyclooctyn-1-yl)-3-oxopropylamino)-3-oxopropyl)disulfanyl)ethyl)-1-(5-(((3aS,4S,6aR)-2-oxohexahydro-1H-thieno[3,4-d]imidazol-4-yl)pentanamido)-3,6,9,12-tetraoxapentadecan-15-amide			
CAS-No.	1807512-43-1		
Formula	C ₄₄ H ₆₀ N ₆ O ₉ S ₃		
Mol. weight	913,18 g/mol		
RL-2420	DBCO-PEG(4)-NH₂*TFA		
Dibenzoazacyclooctyne-tetra(ethylene glycol)-amine trifluoro acetic acid salt			
CAS-No.	1255942-08-5		
Formula	C ₂₉ H ₃₇ N ₃ O ₆ *C ₂ F ₃ HO ₂		
Mol. weight	523,62*114,02 g/mol		
RL-2421	DBCO-Sulfo-PEG(4)-NH₂		
Dibenzoazacyclooctyne-tetra(ethylene glycol)amine			
CAS-No.	2055198-05-3		
Formula	C ₃₂ H ₄₂ N ₄ O ₁₀ S		
Mol. weight	674,76 g/mol		
RL-2510	DBCO-PEG(4)-OH		
Dibenzoazacyclooctyne-tetra(ethylene glycol)			
CAS-No.	1416711-60-8		
Formula	C ₂₉ H ₃₆ N ₂ O ₆		
Mol. weight	508,61 g/mol		

The Click Reaction

Amino Acid Derivatives and Related Building Blocks for Click Chemistry

Spermines and Amines for Click Chemistry

Click Reagents for Drug Delivery

Click Chemistry Tools for Proteomics

Carbohydrates for Click Chemistry

Proteolysis Targeting Chimeras (PROTACs*)

Index

[↑ back to content](#)


PEG6790 Bromoacetamido-PEG(4)-DBCO

Bromoacetamido-tetra(ethylene glycol)-amido-dibenzoazacyclooctyne

CAS-No. 2735663-70-2

Formula $C_{31}H_{38}BrN_3O_7$

Mol. weight 644,55 g/mol


**PEG6740 DBCO-PEG(4)-TFP**

Dibenzoazacyclooctyne-tetra(ethylene glycol)-propionyl 2,3,5,6-tetrafluorophenol ester

CAS-No. 2247993-79-7

Formula $C_{37}H_{38}F_4N_2O_8$

Mol. weight 714,7 g/mol


**RL-2500 DBCO-PEG(4)-mal**

Dibenzoazacyclooctyne-tetra(ethylene glycol)-maleimide

CAS-No. 1480516-75-3

Formula $C_{36}H_{42}N_4O_9$

Mol. weight 674,74 g/mol


**RL-2450 DBCO-PEG(5)-COOH**

Dibenzoazacyclooctyne-penta(ethylene glycol)-propionic acid

CAS-No. 1870899-46-9

Formula $C_{32}H_{40}N_2O_9$

Mol. weight 596,67 g/mol


**RL-2460 DBCO-PEG(5)-NHS**

Dibenzoazacyclooctyne-penta(ethylene glycol)-propionic acid succinimidyl ester

CAS-No. 1378531-80-6

Formula $C_{36}H_{43}N_3O_{11}$

Mol. weight 693,74 g/mol


**PEG6800 Bromoacetamido-PEG(12)-DBCO**







Bromoacetamido-dodeca(ethylene glycol)-amido-dibenzoazacyclooctyne







CAS-No. 2852742-40-4






Formula $C_{47}H_{70}BrN_3O_{17}$

Mol. weight 996,97 g/mol





		Product details
PEG6810	Bromoacetamido-PEG(24)-DBCO	
Bromoacetamido-24(ethylene glycol)-amido-dibenzoazacyclooctyne		
CAS-No.	2852742-74-4	
Formula	C ₇₁ H ₁₁₈ BrN ₃ O ₂₇	
Mol. weight	1525,6 g/mol	
		
PEG6750	DBCO-PEG(12)-TFP	
Dibenzoazacyclooctyne-dodeca(ethylene glycol)-propionyl 2,3,5,6-tetrafluorophenol ester		
CAS-No.	3038549-18-4	
Formula	C ₅₃ H ₇₀ F ₄ N ₂ O ₁₆	
Mol. weight	1067,12 g/mol	
		
PEG6770	DBCO-PEG(12)-MAL	
Dibenzoazacyclooctyne-dodeca(ethylene glycol)-maleimide		
CAS-No.	2011777-01-6	
Formula	C ₅₂ H ₇₄ N ₄ O ₁₇	
Mol. weight	1027,16 g/mol	
		
PEG6830	DBCO-PEG(12)-(5)6-carboxyfluorescein	
Dibenzoazacyclooctyne-dodeca(ethylene glycol)-(5)6-carboxyfluorescein		
Formula	C ₄₆ H ₃₉ N ₃ O ₁₀	
Mol. weight	1234,34 g/mol	
		
PEG7460	DBCO-PEG(24)-OMe	
alpha-Methoxy-24(ethylene glycol)-amido-dibenzoazacyclooctyne		
CAS-No.	3023098-71-4	
Formula	C ₆₈ H ₁₁₄ N ₂ O ₂₆	
Mol. weight	1375,63 g/mol	
		
PEG6760	DBCO-PEG(24)-TFP	
Dibenzoazacyclooctyne-24(ethylene glycol)-propionyl 2,3,5,6-tetrafluorophenol ester		
CAS-No.	2754372-40-0	
Formula	C ₇₇ H ₁₁₈ F ₄ N ₂ O ₂₈	
Mol. weight	1595,75 g/mol	
		

			Product details
PEG6780	DBCO-PEG(24)-MAL		
Dibenzoazacyclooctyne-24(ethylene glycol)-maleimide			
CAS-No.	2924872-84-2		
Formula	C ₇₆ H ₁₂₂ N ₄ O ₂₉		
Mol. weight	1555,79 g/mol		
PEG6765	DBCO-PEG(36)-TFP		
Dibenzoazacyclooctyne-36(ethylene glycol)-propionyl 2,3,5,6-tetrafluorophenol ester			
CAS-No.	2924873-16-3		
Formula	C ₁₀₁ H ₁₆₆ F ₄ N ₂ O ₄₀		
Mol. weight	2124,41 g/mol		
RL-2530	DBCO-mPEG (5kDa)		
alpha-Dibenzoazacyclooctyne-omega-methoxy-poly(ethylene glycol)			
CAS-No.	2262541-53-5		
Mol. weight	5000 Da		
RL-2540	DBCO-mPEG (10kDa)		
alpha-Dibenzoazacyclooctyne-omega-methoxy-poly(ethylene glycol)			
CAS-No.	2262541-53-5		
Mol. weight	10000 Da		
RL-2550	DBCO-mPEG (20kDa)		
alpha-Dibenzoazacyclooctyne-omega-methoxy-poly(ethylene glycol)			
CAS-No.	2262541-53-5		
Mol. weight	20000 Da		
RL-2560	DBCO-mPEG (30kDa)		
alpha-Dibenzoazacyclooctyne-omega-methoxy-poly(ethylene glycol)			
CAS-No.	2262541-53-5		
Mol. weight	30000 Da		

			Product details
RL-4110	DBCO-Suc-SS-COOH		
CAS-No.	2749426-25-1		
Formula	$C_{26}H_{24}N_2O_4S_2$		
Mol. weight	468,59 g/mol		
RL-4310	DBCO-Suc-CBT		
N1-(2-cyanobenzo[d]thiazol-6-yl)-N4-(3-(11,12-didehydro-5,6-dihydro-dibenzo[b,f]azocin-yl)-3-oxopropyl) succinamide			
Formula	$C_{30}H_{23}N_5O_3S$		
Mol. weight	533,61 g/mol		
PEG7070	Tetra(-PEG(11)-DBCO)pentaerythritol		
Formula	$C_{193}H_{288}N_{12}O_{60}$		
Mol. weight	3736,45 g/mol		
PEG7075	Bis-PEG(11)-DBCO		
Formula	$C_{64}H_{82}N_4O_{15}$		
Mol. weight	1147,37 g/mol		
PEG7095	DBCO-PEG(24)-amido-PEG(24)-DSPE		
Formula	$C_{163}H_{299}N_4O_{60}P$		
Mol. weight	3306,13 g/mol		

Products with DACN

		Product details
RL-3600	DACN(Ms)*HCl N-(Mesyl)-4,8-diazacyclononyne hydrochloride CAS-No. 2331322-16-6 Formula $C_8H_{16}N_2O_2S^+HCl$ Mol. weight 202,27*36,46 g/mol	
RL-3610	DACN(Ms,Ns) N-(Mesyl)-N'-(2-nosyl)-4,8-diazacyclononyne CAS-No. 2411082-25-0 Formula $C_{14}H_{17}N_3O_5S_2$ Mol. weight 387,43 g/mol	
RL-2735	DACN(Tos)*HCl N-(p-toluenesulfonyl)-4,8-diazacyclononyne hydrochloride CAS-No. 2331322-18-8 Formula $C_{14}H_{18}N_2O_2S^+HCl$ Mol. weight 278,37*36,46 g/mol	
RL-2720	DACN(Tos,Suc-OH) N-succinoyl-N'-(p-toluenesulfonyl)-4,8-diazacyclononyne CAS-No. 2109751-68-8 Formula $C_{18}H_{22}N_2O_5S$ Mol. weight 378,44 g/mol	
RL-2710	DACN(Tos,Ns) N-(o-nitrobenzenesulfonyl)-N'-(p-toluenesulfonyl)-4,8-diazacyclononyne CAS-No. 1797508-58-7 Formula $C_{20}H_{21}N_3O_6S_2$ Mol. weight 463,53 g/mol	

		Product details
RL-2730 DACN(Tos2) N,N'-bis(p-toluenesulfonyl)-4,8-diazacyclononyne CAS-No. 1797508-57-6 Formula C ₂₁ H ₂₄ N ₂ O ₄ S ₂ Mol. weight 432,56 g/mol		
RL-2737 DACN(Tos2,6-OH) 4,8-Bis(p-toluenesulfonyl)-4,8-diazacyclononyl-6-ol CAS-No. 2109751-74-6 Formula C ₂₁ H ₂₄ N ₂ O ₅ S ₂ Mol. weight 448,55 g/mol		

References:

- A strain-promoted [3 + 2] azide-alkyne cycloaddition for covalent modification of biomolecules in living systems; N. J. Agard, J. A. Prescher, C. R. Bertozzi; **J Am Chem Soc** 2004; **126**: 15046-7. <https://doi.org/10.1021/ja044996f>
- A hydrophilic azacyclooctyne for Cu-free click chemistry; E. M. Sletten, C. R. Bertozzi; **Org Lett** 2008; **10**: 3097-9. <https://doi.org/10.1021/ol801141k>
- Live-cell imaging of cellular proteins by a strain-promoted azide-alkyne cycloaddition; K. E. Beatty, J. D. Fisk, B. P. Smart, Y. Y. Lu, J. Szychowski, M. J. Hangauer, J. M. Baskin, C. R. Bertozzi, D. A. Tirrell; **ChemBiochem** 2010; **11**: 2092-5. <https://doi.org/10.1002/cbic.201000419>
- Bioconjugation with strained alkenes and alkynes; M. F. Debets, S. S. van Berkel, J. Dommerholt, A. T. Dirks, F. P. Rutjes, F. L. van Delft; **Acc Chem Res** 2011; **44**: 805-15. <https://doi.org/10.1021/ar200059z>
- Reactivity of biarylazacyclooctynones in copper-free click chemistry; C. G. Gordon, J. L. Mackey, J. C. Jewett, E. M. Sletten, K. N. Houk, C. R. Bertozzi; **J Am Chem Soc** 2012; **134**: 9199-208. <https://doi.org/10.1021/ja3000936>
- Bioorthogonal labelling of biomolecules: new functional handles and ligation methods; M. F. Debets, J. C. van Hest, F. P. Rutjes; **Org Biomol Chem** 2013; **11**: 6439-55. <https://doi.org/10.1039/c3ob41329b>
- 18F-labeling using click cycloadditions; K. Kettenbach, H. Schieferstein, T. L. Ross; **Biomed Res Int** 2014; **2014**: 361329. <https://doi.org/10.1155/2014/361329>
- Click Chemistry and Radiochemistry: The First 10 Years; J. P. Meyer, P. Adumeau, J. S. Lewis, B. M. Zeglis; **Bioconjug Chem** 2016; **27**: 2791-2807. <https://doi.org/10.1021/acs.bioconjchem.6b00561>
- Strain-Promoted 1,3-Dipolar Cycloaddition of Cycloalkynes and Organic Azides; J. Dommerholt, F. Rutjes, F. L. van Delft; **Top Curr Chem** 2016; **374**: 16. <https://doi.org/10.1007/s41061-016-0016-4>
- A facile preparation of functional cycloalkynes via an azide-to-cycloalkyne switching approach; S. Yoshida, T. Kuribara, H. Ito, T. Meguro, Y. Nishiyama, F. Karaki, Y. Hatakeyama, Y. Koike, I. Kii, T. Hosoya; **Chem Commun (Camb)** 2019; **55**: 3556-3559. <https://doi.org/10.1039/c9cc01113g>
- Selective strain-promoted azide-alkyne cycloadditions through transient protection of bicyclo[6.1.0] nonynes with silver or gold; K. Adachi, T. Meguro, Y. Sakata, K. Igawa, K. Tomooka, T. Hosoya, S. Yoshida; **Chem Commun (Camb)** 2020; **56**: 9823-9826. <https://doi.org/10.1039/d0cc04606j>
- Facile assembly of three cycloalkyne-modules onto a platform compound bearing thiophene S,S-dioxide moiety and two azido groups; T. Meguro, Y. Sakata, T. Morita, T. Hosoya, S. Yoshida; **Chem Commun (Camb)** 2020; **56**: 4720-4723. <https://doi.org/10.1039/d0cc01810d>
- Copper-Free Huisgen Cycloaddition for the 14-3-3-Templated Synthesis of Fusicoccin-Peptide Conjugates; R. Masuda, Y. Kawasaki, K. Igawa, Y. Manabe, H. Fujii, N. Kato, K. Tomooka, J. Ohkanda; **Chem Asian J** 2020; **15**: 742-747. <https://doi.org/10.1002/asia.202000042>

[↑ back to content](#)

The Click Reaction

Amino Acid Derivatives and Related Building Blocks for Click Chemistry

Spermines and Amines for Click Chemistry

Click Reagents for Drug Delivery

Click Chemistry Tools for Proteomics

Carbohydrates for Click Chemistry


Proteolysis Targeting Chimeras (PROTACs®)

Index


1.2.2. 3rd Generation: Inverse Electron-Demand Diels-Alder (IEDDA) Reactions

Click chemistry is frequently the method of choice for site-selective labeling and crosslinking. However, in biological systems, the cytotoxicity of copper used for the classical Cu-promoted 1,3-dipolar cycloaddition may cause major problems. The copper-free strain-promoted alkyne-azide cycloaddition (SPAAC) utilizing cyclooctynes, on the other hand, is limited by its moderate reaction kinetics for the application in live cells, where the concentration of biomolecules is usually low. Another potential drawback of cyclooctynes is the extensive patent coverage of many variants.

Tetrazine ligation presents the option for a copper-free, rapid, and fully bioorthogonal type of Click chemistry. Mechanistically, this reaction proceeds *via* an inverse electron-demand Diels-Alder cycloaddition reaction between a tetrazine and a strained alkene such as *trans*-cyclooctene (TCO), cyclopropane or norbornene, followed by a retro-Diels-Alder reaction under elimination of N₂, the latter rendering the reaction irreversible.




Stability vs. Faster Reaction Kinetics: 6-Me or 6-H Tetrazines



There are two main types of tetrazines that are widely applied: 6-methyl-substituted tetrazines and 6-hydrogen-substituted tetrazines. Methyl-substituted tetrazines exhibit a high stability even when dissolved in aqueous media, while still offering faster reaction kinetics with TCO derivatives than any other bioorthogonal reaction pairs (approx. 1000 M⁻¹s⁻¹). Moreover, they tolerate a wide array of reaction conditions which renders them the prime choice for applications such as protein labeling. Hydrogen-substituted tetrazines, on the other hand, show lower stability and less tolerance to harsh reaction conditions, but offer extremely fast reaction kinetics (up to 30000 M⁻¹s⁻¹) for applications like *in vivo* imaging.

This method excels at very low concentrations (e.g., in biological systems) due to the extremely rapid second order reaction rate constants (between approx. $800 \text{ M}^{-1}\text{s}^{-1}$ and $30000 \text{ M}^{-1}\text{s}^{-1}$). Moreover, the tetrazine-TCO ligation can be performed in aqueous media and has been applied in live cell imaging. These properties make tetrazine Click chemistry the method of choice for labeling or crosslinking biomolecules in living cells.




Choice of Spacer: Alkyl or PEG?

Tetrazines equipped with alkyl spacers are suitable for reactions in organic solvents. For applications in aqueous media, however, PEG-spacers are usually the superior choice. Moreover, tetrazines equipped with PEG-spacers are ideal for the functionalization of proteins since PEGs are known to reduce the aggregation of labeled polypeptides.

In summary, the reaction between a tetrazine (Tz) and a *trans*-cyclooctene (TCO) is the innovative third generation Click reaction that proceeds without the use of copper or other catalysts. It is rapid, fully bioorthogonal, irreversible and excels at very low concentrations.

Products with Tetrazine

		Product details
HAA9470	H-L-Phe(4-MeTz)-OH*TFA	
(S)-2-amino-3-(4-(6-methyl-1,2,4,5-tetrazin-3-yl)phenyl)propanoic acid trifluoroacetic acid salt		
CAS-No.	1698038-85-5 net	 
Formula	$\text{C}_{12}\text{H}_{13}\text{N}_5\text{O}_2 \cdot \text{CF}_3\text{COOH}$	
Mol. weight	259,27*114,02 g/mol	


[↑ back to content](#)

HAA9480 H-L-Phe(4-Azido-PrTz)-OH*TFA

(S)-2-amino-3-(4-(6-(3-azidopropyl)-1,2,4,5-tetrazin-3-yl)phenyl)propanoic acid trifluoroacetic acid salt

Formula $C_{14}H_{16}N_8O_2 \cdot CF_3COOH$

Mol. weight 328,34*114,02 g/mol


**HAA9490 H-L-Phe(3-MeTz)-OH*TFA**

(S)-2-amino-3-(3-(6-methyl-1,2,4,5-tetrazin-3-yl)phenyl)propanoic acid trifluoroacetic acid salt

CAS-No. 2036323-75-6 net

Formula $C_{12}H_{13}N_5O_2 \cdot CF_3COOH$

Mol. weight 259,27*114,02 g/mol


**HAA9500 H-L-Phe(3-BuTz)-OH*TFA**

(S)-2-amino-3-(3-(6-butyl-1,2,4,5-tetrazin-3-yl)phenyl)propanoic acid trifluoroacetic acid salt

CAS-No. 2036323-83-6 net

Formula $C_{15}H_{19}N_5O_2 \cdot CF_3COOH$

Mol. weight 301,35*114,02 g/mol


**HAA9510 H-L-Phe(3-iPrTz)-OH*TFA**

(S)-2-amino-3-(3-(6-isopropyl-1,2,4,5-tetrazin-3-yl)phenyl)propanoic acid trifluoroacetic acid salt

CAS-No. 2421119-12-0 net

Formula $C_{14}H_{17}N_5O_2 \cdot CF_3COOH$

Mol. weight 287,32*114,02 g/mol


**HAA9170 H-L-Lys(MeTz-PhAc)-OH*TFA**

N-(2-(4-(6-methyl-1,2,4,5-tetrazin-3-yl)phenyl)acetyl)-L-lysine TFA salt

CAS-No. 2578384-82-2 (net)

Formula $C_{17}H_{22}N_6O_3 \cdot CF_3COOH$

Mol. weight 358,40*114,02 g/mol


**RL-2140 (Me)Tz-butanoic acid**







4-(6-methyl-1,2,4,5-tetrazin-3-yl)butanoic acid







CAS-No. 1923268-81-8





Formula $C_7H_{10}N_4O_2$

Mol. weight 182,18 g/mol










		Product details
RL-3905	MeTz-PEG(4)-STP	
sodium 2,3,5,6-tetrafluoro-4-((1-(4-(6-methyl-1,2,4,5-tetrazin-3-yl)phenoxy)-3,6,9,12-tetraoxapentadecan-15-oyl)oxy)benzenesulfonate		
Formula	$C_{26}H_{27}F_4N_4NaO_{10}S$	
Mol. weight	686,56 g/mol	
RL-2360	MeTz-Bzl-NH₂*HCl	
Methyltetrazine-benzylamine*HCl		
CAS-No.	1345955-28-3	
Formula	$C_{10}H_{11}N_5*HCl$	
Mol. weight	201,23*36,46 g/mol	
RL-2130	(Me)Tz-benzoic acid	
4-(6-methyl-1,2,4,5-tetrazin-3-yl)benzoic acid		
CAS-No.	1345866-66-1	
Formula	$C_{10}H_8N_4O_2$	
Mol. weight	216,2 g/mol	
RL-2580	Tz-benzoic acid	
4-(1,2,4,5-tetrazin-3-yl)benzoic acid		
CAS-No.	1345866-65-0	
Formula	$C_9H_6N_4O_2$	
Mol. weight	202,17 g/mol	
RL-2300	MeTz-PhAcOH	
Methyltetrazine-phenylacetic acid		
CAS-No.	1380500-88-8	
Formula	$C_{11}H_{10}N_4O_2$	
Mol. weight	230,22 g/mol	
RL-2320	MeTz-PhAc-NHS	
Methyltetrazine-phenylacetyl succinimidyl ester		
CAS-No.	1644644-96-1	
Formula	$C_{15}H_{13}N_5O_4$	
Mol. weight	327,29 g/mol	

			Product details
RL-3915	MeTz-PhAc-Sulfo-NHS		
sodium 1-(2-(4-(6-methyl-1,2,4,5-tetrazin-3-yl)phenyl)acetoxy)-2,5-dioxopyrrolidine-3-sulfonate			
CAS-No.	1821017-46-2		
Formula	C ₁₅ H ₁₂ N ₅ NaO ₇		
Mol. weight	429,34 g/mol		
RL-2230	Bz-(Me)Tz-NHS		
2,5-dioxopyrrolidin-1-yl 5-(4-(6-methyl-1,2,4,5-tetrazin-3-yl)benzylamino)-5-oxopentanoate			
CAS-No.	1454558-58-7		
Formula	C ₂₄ H ₂₅ N ₅ O ₇		
Mol. weight	412,41 g/mol		
RL-2240	Bz-Tz-NHS		
2,5-dioxopyrrolidin-1-yl 5-(4-(1,2,4,5-tetrazin-3-yl)benzylamino)-5-oxopentanoate			
CAS-No.	1244040-64-9		
Formula	C ₁₈ H ₁₈ N ₆ O ₅		
Mol. weight	398,37 g/mol		
RL-2370	MeTz-PEG(4)-NH₂*HCl		
Methyltetrazine-PEG(4)-amine HCL salt			
CAS-No.	1802908-05-9 net		
Formula	C ₁₇ H ₂₅ N ₅ O ₄ *HCl		
Mol. weight	363,41*HCl g/mol		
RL-2310	MeTz-PEG(4)-COOH		
Methyltetrazine-PEG(4)-acid			
CAS-No.	1802907-91-0		
Formula	C ₂₀ H ₂₈ N ₄ O ₇		
Mol. weight	436,56 g/mol		
RL-2330	MeTz-PEG(4)-NHS		
Methyltetrazine-PEG(4)-propanoyl succinimidyl ester			
CAS-No.	1802907-92-1		
Formula	C ₂₄ H ₃₁ N ₅ O ₉		
Mol. weight	533,53 g/mol		




		Product details
RL-2340	MeTz-PEG(4)-mal	
Methyltetrazine-PEG(4)-maleimide		
CAS-No.	1802908-02-6	 
Formula	C ₂₄ H ₃₀ N ₆ O ₇	
Mol. weight	514,53 g/mol	
RL-2250	Bz-Tz-PEG(5)-NHS	
2,5-dioxopyrrolidin-1-yl 1-(4-(1,2,4,5-tetrazin-3-yl)phenyl)-3-oxo-6,9,12,15,18-pentaoxa-2-azahenicosan-21-oate		
CAS-No.	1682653-80-0	 
Formula	C ₂₇ H ₃₆ N ₆ O ₁₀	
Mol. weight	604,61 g/mol	
LS-4280	Biotin-MeTz	
N-(4-(6-methyl-1,2,4,5-tetrazin-3-yl)benzyl)-5-((3aS,4S,6aR)-2-oxohexahydro-1H-thieno[3,4-d]imidazol-4-yl)pentanamide		
CAS-No.	1802883-51-7	 
Formula	C ₂₀ H ₂₅ N ₇ O ₂ S	
Mol. weight	427,53 g/mol	
LS-4290	Biotin-PEG(4)-MeTz	
N-(4-(6-methyl-1,2,4,5-tetrazin-3-yl)benzyl)-1-(5-((3aS,4S,6aR)-2-oxohexahydro-1H-thieno[3,4-d]imidazol-4-yl)pentanamido)-3,6,9,12-tetraoxapentadecan-15-amide		
CAS-No.	1962919-31-8	 
Formula	C ₃₁ H ₄₆ N ₈ O ₇ S	
Mol. weight	674,82 g/mol	

Products with TCO

		Product details
TCO1070	TCO-PEG(3)-NH₂*HCl	
<i>trans</i> -Cyclooctene-PEG(3)-amine		
CAS-No.	2028288-77-7	 
Formula	C ₁₉ H ₃₆ N ₂ O ₅ *HCl	
Mol. weight	372,51*36,46 g/mol	

			Product details
TCO1060	TCO-NH ₂ *HCl		
<i>trans</i> -Cyclooctene-amine hydrochloride			
CAS-No.	1800507-94-1		
Formula	C ₁₂ H ₂₂ N ₂ O ₂ *HCl		
Mol. weight	226,32*36,45 g/mol		
TCO1000	TCO-NHS		
<i>trans</i> -Cyclooctene succinimidyl carbonate			
CAS-No.	1191901-33-3		
Formula	C ₁₃ H ₁₇ NO ₅		
Mol. weight	267,28 g/mol		
TCO1040	TCO-PEG(4)-COOH		
<i>trans</i> -Cyclooctene-PEG(4)-Acid			
CAS-No.	1802913-21-8		
Formula	C ₂₀ H ₃₅ NO ₈		
Mol. weight	417,49 g/mol		
TCO1010	TCO-PEG(4)-NHS		
<i>trans</i> -Cyclooctene-PEG(4)-carboxy succinimidyl ester			
CAS-No.	1621096-79-4		
Formula	C ₂₄ H ₃₈ N ₂ O ₁₀		
Mol. weight	514,57 g/mol		
TCO1050	TCO-PEG(3)-mal		
<i>trans</i> -Cyclooctene-PEG(3)-maleimide			
CAS-No.	1809356-72-6		
Formula	C ₂₆ H ₄₁ N ₃ O ₈		
Mol. weight	523,62 g/mol		
TCO1020	TCO-PEG(12)-NHS		
<i>trans</i> -Cyclooctene-PEG(12)-carboxy succinimidyl ester			
CAS-No.	2185016-39-9		
Formula	C ₄₀ H ₇₀ N ₂ O ₁₈		
Mol. weight	866,99 g/mol		

Products with Norbornene

		Product details
HAA9235 H-L-Lys(Norbornene-methoxycarbonyl)-OH*HCl N-epsilon-(norbornene-methoxycarbonyl)-L-lysine hydrochloride CAS-No. 1378916-76-7 Formula C ₁₅ H ₂₄ N ₂ O ₄ *HCl Mol. weight 296,37*36,46 g/mol		
RL-2080 Norbornene-NHS (Norbornene-2-yl)-N-hydroxysuccinimidylcarbonate CAS-No. 1888335-48-5 Formula C ₁₂ H ₁₃ NO ₅ Mol. weight 251,24 g/mol		
RL-2090 Norbornene-methyl-NHS (Norbornene-2-methyl)-N-hydroxysuccinimidylcarbonate CAS-No. 1986791-87-0 Formula C ₁₃ H ₁₅ NO ₅ Mol. weight 265,26 g/mol		

References:

- *Biomedical applications of tetrazine cycloadditions*; N. K. Devaraj, R. Weissleder; **Acc Chem Res** 2011; **44**: 816-27. <https://doi.org/10.1021/ar200037t>
- *trans-Cyclooctene—a stable, voracious dienophile for bioorthogonal labeling*; R. Selvaraj, J. M. Fox; **Curr Opin Chem Biol** 2013; **17**: 753-60. <https://doi.org/10.1016/j.cbpa.2013.07.031>
- *Inverse electron demand Diels-Alder (IEDDA)-initiated conjugation: a (high) potential click chemistry scheme*; A. C. Knall, C. Slugovc; **Chem Soc Rev** 2013; **42**: 5131-42. <https://doi.org/10.1039/c3cs60049a>
- *The growing impact of bioorthogonal click chemistry on the development of radiopharmaceuticals*; D. Zeng, B. M. Zeglis, J. S. Lewis, C. J. Anderson; **J Nucl Med** 2013; **54**: 829-32. <https://doi.org/10.2967/jnumed.112.115550>
- *Highly accelerated inverse electron-demand cycloaddition of electron-deficient azides with aliphatic cyclooctynes*; J. Dommerholt, O. van Rooijen, A. Borrmann, C. F. Guerra, F. M. Bickelhaupt, F. L. van Delft; **Nat Commun** 2014; **5**: 5378. <https://doi.org/10.1038/ncomms6378>
- *The inverse electron demand Diels-Alder click reaction in radiochemistry*; T. Reiner, B. M. Zeglis; **J Labelled Comp Radiopharm** 2014; **57**: 285-90. <https://doi.org/10.1002/jlcr.3149>
- *3,6-Substituted-1,2,4,5-tetrazines: tuning reaction rates for staged labeling applications*; D. Wang, W. Chen, Y. Zheng, C. Dai, K. Wang, B. Ke, B. Wang; **Org Biomol Chem** 2014; **12**: 3950-5. <https://doi.org/10.1039/c4ob00280f>
- *Click Chemistry and Radiochemistry: The First 10 Years*; J. P. Meyer, P. Adumeau, J. S. Lewis, B. M. Zeglis; **Bioconjug Chem** 2016; **27**: 2791-2807. <https://doi.org/10.1021/acs.bioconjchem.6b00561>
- *Bio-orthogonal Fluorescent Labelling of Biopolymers through Inverse-Electron-Demand Diels-Alder Reactions*; E. Kozma, O. Demeter, P. Kele; **Chembiochem** 2017; **18**: 486-501. <https://doi.org/10.1002/cbic.201600607>
- *Bicyclo[6.1.0]nonyne and tetrazine amino acids for Diels–Alder reactions*; X. Li, Z. Liu, S. Dong; **RSC Advances** 2017; **7**: 44470-44473. <https://doi.org/10.1039/c7ra08136g>

[↑ back to content](#)

- Inverse electron demand Diels-Alder reactions in chemical biology; B. L. Oliveira, Z. Guo, G. J. L. Bernardes; **Chem Soc Rev** 2017; **46**: 4895-4950. <https://doi.org/10.1039/c7cs00184c>
- Advances in Tetrazine Bioorthogonal Chemistry Driven by the Synthesis of Novel Tetrazines and Dienophiles; H. Wu, N. K. Devaraj; **Acc Chem Res** 2018; **51**: 1249-1259. <https://doi.org/10.1021/acs.accounts.8b00062>
- Inverse electron demand Diels-Alder (IEDDA) reactions in peptide chemistry; M. Pagel; **J. Pept. Sci.** 2019; **25**: e3141. <https://doi.org/10.1002/psc.3141>
- Bioorthogonal Fluorescence Turn-On Labeling Based on Bicyclononyne-Tetrazine Cycloaddition Reactions that Form Pyridazine Products; S. J. Siegl, J. Galeta, R. Dzijak, M. Dracinsky, M. Vrabel; **Chempluschem** 2019; **84**: 493-497. <https://doi.org/10.1002/cplu.201900176>
- An Extended Approach for the Development of Fluorogenic trans-Cyclooctene-Tetrazine Cycloadditions; S. J. Siegl, J. Galeta, R. Dzijak, A. Vazquez, M. Del Rio-Villanueva, M. Dracinsky, M. Vrabel; **ChemBiochem** 2019; **20**: 886-890. <https://doi.org/10.1002/cbic.201800711>
- Structural insights into incorporation of norbornene amino acids for click modification of proteins; S. Schneider, M. J. Gattner, M. Vrabel, V. Fluegel, V. Lopez-Carrillo, S. Pril, T. Carell; **ChemBioChem** 2013; **14(16)**: 2114-2118. <https://doi.org/10.1002/cbic.201300435>
- A need for speed: genetic encoding of rapid cycloaddition chemistries for protein labelling in living cells; M. J. Schmidt, D. Summerer; **ChemBioChem** 2012; **13(11)**: 1553-1557. <https://doi.org/10.1002/cbic.201200321>
- A genetically encoded norbornene amino acid for the mild and selective modification of proteins in a copper-free click reaction; E. Kaya, M. Vrabel, C. Deiml, S. Prill, V. S. Fluxa, T. Carell; **Angew. Chem. Int. Ed.** 2012; **51(18)**: 4466-4469. <https://doi.org/10.1002/anie.201109252>
- Amino Acids for Diels-Alder Reactions in Living Cells; T. Plass, S. Milles, C. Koehler, J. Szymanski, R. Mueller, M. Wiessler, C. Schultz, E. A. Lemke; **Angew. Chem. Int. Ed.** 2012; **51(17)**: 4166-4170. <https://doi.org/10.1002/anie.201108231>



We offer custom synthesis of DBCO, DACN, TCO and tetrazine derivatives and related conjugations.


Send us your inquiry to info@iris-biotech.de!

1.2.3. Custom Synthesis of DBCO, Tetrazine and TCO Derivatives

Table 1: A selection of derivatives of tetrazine, TCO and DBCO available by custom synthesis.

DBCO with Cleavable Linkers

e.g., ADC linkers, Dde-based linkers, disulfide-based linkers.



DBCO-PEG-Derivatives

DBCO-PEG-NHS


DBCO-PEG-mal

DBCO-PEG-Bis-Sulfone-Thiol

DBCO-PEG-COOH


DBCO-PEG-NH₂

With mono- or polydisperse PEG.



DBCO-Biotin


With both monodisperse and polydisperse PEG-spacers.



DBCO-Dye

With both monodisperse and polydisperse PEG-spacers.

With the dye of your choice, e.g., ICG, (5)6-carboxyfluorescein.



Tetrazine-PEG-Derivatives


Tetrazine-PEG-NHS

Tetrazine-PEG-mal

Tetrazine-PEG-COOH


Tetrazine-PEG-NH₂

With both monodisperse and polydisperse PEG-spacers.




Tetrazine-Biotin

With both monodisperse and polydisperse PEG-spacers, or alkyl-spacers.




Tetrazine-Dye

With both monodisperse and polydisperse PEG-spacers. With the dye of your choice, e.g., ICG, (5)6-carboxyfluorescein.



TCO-PEG-Derivatives

TCO-PEG-NHS
TCO-PEG-mal
TCO-PEG-COOH
TCO-PEG-NH₂



You want to know more about our custom synthesis capabilities?


Check out our brochure!



2. Amino Acid Derivatives and Related Building Blocks for Click Chemistry

2.1. Recombinant Incorporation of Amino Acids into Proteins


Genetic code expansion is a powerful technology in proteomics, facilitating the site-specific incorporation of noncanonical amino acids (ncAAs) into proteins using the cellular machinery. A wide variety of ncAAs can be incorporated into proteins using this technology that relies on aminoacyl-tRNA synthetase/tRNA pairs.



Certain amino acids such as azidohomoalanine (Aha) can be incorporated into proteins using the cell's native translational apparatus. Aha is a structural analog of methionine (Met), and as such activated by the native methionyl-tRNA synthetase of *Escherichia coli*, replacing methionine in proteins expressed in methionine-depleted bacterial cultures. Aside from being a clickable amino acid, azidohomoalanine is an excellent conformationally sensitive IR probe to study protein folding and protein structure.

In most cases, researchers resort to engineering suitable aminoacyl-tRNA synthetase/tRNA pairs in order to incorporate ncAAs. For example, the usually promiscuous pyrrolysyl-tRNA synthetase (PylRS) machinery can be engineered to accommodate more than 100 ncAAs or α -hydroxy acids into proteins at amber codons, and can be reassigned to other codons such as ochre (UAA) or opal (UGA). Among the most prominent noncanonical amino acids that are routinely incorporated by engineered PylRS/ tRNA^{Pyl} pairs are azido and propargyl analogs of L-lysine, enabling the biochemist to site-specifically introduce an azido or alkyne group into a protein for further Click conjugation.

Recent developments in the field of genetic code expansion include the directed evolution of tRNA synthetases to improve substrate selectivity, as well as the reassignment of further codons to encode ncAAs. Once an azido or alkyne function has been built into the protein sequence, conjugation with a large number of diverse clickable compounds opens up a wide field of possibilities. Alternatively, a protein bearing an azido group can be selectively modified *via* Staudinger ligation (see Fig. 10). Many different applications from therapeutics to diagnostics can be addressed through conjugates with PEG-polymers, dyes, cofactors, antibodies, small molecules, toxins, additional proteins, and peptides.




One recent example of a sophisticated application of genetic code expansion is the site-directed incorporation of two different non-canonical amino acids into human erythropoietin via cell-free protein synthesis (Zemella *et al.*, *Sci. Rep.* 2018). Either *p*-propargyloxyphenylalanine (pPa) or *p*-azido-L-phenyl-alanine (AzF) was incorporated into an erythropoietin amber-mutant (EPO-Amb) via amber suppression in a eukaryotic translationally active lysate. This eukaryotic system also facilitated the glycosylation of EPO which is known to be crucial for its pharmacokinetics. The recombinant EPO variants were subsequently labelled with various fluorophores, as well as functionalized with a PEG of 10 kDa. Similar to glycosylation, the attachment of PEGs has been shown to improve solubility, stability and activity of recombinantly produced EPO (Hoffmann *et al.*, *Mol. Biosyst.* 2016).

The usefulness of the cell-free protein synthesis approach for the incorporation of ncAAs was further demonstrated by the preparation and ligand-free dimerization of functional human epidermal growth factor receptor (EGFR), a complex eukaryotic transmembrane protein (Quast *et al.*, *Sci. Rep.* 2016). EGFR is a receptor tyrosine kinase that dimerizes and autophosphorylates upon binding to its ligand, thereby initiating an intracellular signal transduction cascade. In order to facilitate dimerization in the absence of a ligand, *p*-azido-L-phenylalanine (AzF) was site-selectively incorporated into EGFR, which was verified by Staudinger ligation of a phosphine dye to AzF. Two different EGFR amber mutants that incorporate AzF in the intracellular juxtamembrane domain were synthesized and reacted with a bis-COMBO Click-cross-linking reagent, thereby generating covalently linked receptor dimers.

References:

- Genetic Code Expansion: Inception, Development, Commercialization; M. Manandhar, E. Chun, F. E. Romesberg; *J Am Chem Soc* 2021. <https://doi.org/10.1021/jacs.0c11938>
- Cell-free protein synthesis as a novel tool for directed glycoengineering of active erythropoietin; A. Zemella, L. Thoring, C. Hoffmeister, M. Samalikova, P. Ehren, D. A. Wustenhagen, S. Kubick; *Sci Rep* 2018; **8**: 8514. <https://doi.org/10.1038/s41598-018-26936-x>
- Cell-free synthesis of functional human epidermal growth factor receptor: Investigation of ligand-independent dimerization in Sf21 microsomal membranes using non-canonical amino acids; R. B. Quast, B. Ballion, M. Stech, A. Sonnabend, B. R. Varga, D. A. Wustenhagen, P. Kele, S. M. Schiller, S. Kubick; *Sci Rep* 2016; **6**: 34048. <https://doi.org/10.1038/srep34048>
- Cotranslational incorporation of non-standard amino acids using cell-free protein synthesis; R. B. Quast, D. Mrusek, C. Hoffmeister, A. Sonnabend, S. Kubick; *FEBS Lett* 2015; **589**: 1703-12. <https://doi.org/10.1016/j.febslet.2015.04.041>
- Pyrrolysyl-tRNA synthetase: an ordinary enzyme but an outstanding genetic code expansion tool; W. Wan, J. M. Tharp, W. R. Liu; *Biochim Biophys Acta* 2014; **1844**: 1059-70. <https://doi.org/10.1016/j.bbapap.2014.03.002>
- Genetic encoding and labeling of aliphatic azides and alkynes in recombinant proteins via a pyrrolysyl-tRNA Synthetase/tRNA(CUA) pair and click chemistry; D. P. Nguyen, H. Lusic, H. Neumann, P. B. Kapadnis, A. Deiters, J. W. Chin; *J Am Chem Soc* 2009; **131**: 8720-1. <https://doi.org/10.1021/ja900553w>
- High-throughput screening for methionyl-tRNA synthetases that enable residue-specific incorporation of non-canonical amino acids into recombinant proteins in bacterial cells; T. H. Yoo, D. A. Tirrell; *Angew. Chem. Int. Ed.* 2007; **46**: 5340-3. <https://doi.org/10.1002/anie.200700779>
- Direct charging of tRNACUA with pyrrolysine in vitro and in vivo; S. K. Blight, R. C. Larue, A. Mahapatra, D. G. Longstaff, E. Chang, G. Zhao, P. T. Kang, K. B. Green-Church, M. K. Chan, J. A. Krzycki; *Nature* 2004; **431**: 333. <https://doi.org/10.1038/nature02895>
- Incorporation of azides into recombinant proteins for chemoselective modification by the Staudinger ligation; K. L. Kiick, E. Saxon, D. A. Tirrell, C. R. Bertozzi; *Proc Natl Acad Sci U S A* 2002; **99**: 19-24. <https://doi.org/10.1073/pnas.012583299>
- Global replacement of tryptophan with aminotryptophans generates non-invasive protein-based optical pH sensors; N. Budisa, M. Rubini, J. H. Bae, E. Weyher, W. Wenger, R. Golbik, R. Huber, L. Moroder; *Angew. Chem. Int. Ed.* 2002; **41**: 4066-9. [https://doi.org/10.1002/1521-3773\(20021104\)41:21<4066::AID-ANIE4066>3.0.CO;2-6](https://doi.org/10.1002/1521-3773(20021104)41:21<4066::AID-ANIE4066>3.0.CO;2-6)

			Product details
LS-3405	TCEP*HCl		
3-[bis(2-carboxyethyl)phosphanyl]propanoic acid, hydrochloride			
CAS-No.	51805-45-9		
Formula	C ₉ H ₁₅ O ₆ P*HCl		
Mol. weight	250,19*36,45 g/mol		
			

As part of daily lab procedures, proteins and peptides need to be kept in their reduced state, e.g., for protein analysis, to maintain their activity, to prevent denaturation, to couple them to carriers or payloads, or to inactivate RNAses.

For this purpose, DTT and BME are frequently used, however, they come with certain drawbacks: DTT and BME are easily oxidized by ambient air, they may react with heavy metal ions and interfere with metal ion affinity chromatography (IMAC), they are quenching thiol-reactive reagents like maleimides, and, in the case of BME, it is malodorous.

Discover TCEP as nearly odorless alternative!

2.2. Peptide Synthesis with Azido and Alkyne Amino Acids

Both Boc- and Fmoc-protected derivatives of azido and alkyne amino acids can be readily introduced into peptide sequences by standard SPPS protocols. Such building blocks have found widespread use in techniques such as peptide ligation, bioconjugation, labeling, immobilization and linkerology. Bioconjugation, which is defined as the joining of two biomolecules or the ligation of a synthetic molecule with a biomolecule, stands out in particular among those applications. Targets that are notoriously difficult to access such as glycopeptides and -proteins can be synthesized in a straightforward and chemoselective fashion *via* Click reaction to afford neoglycopeptides and -proteins. Peptides or proteins that aid in the translocation into cells or that facilitate the targeting to certain tissues or organelles may be conjugated to toxins, fluorophores, or oligonucleotides by means of the Click reaction.

Another potential application is the cyclization of peptides *via* Click chemistry. This technique is a well-known approach to stabilize specific conformations in order to optimize peptide binding, and to increase resistance toward proteolytic degradation. If two clickable groups are placed at a suitable distance from each other in a peptide, they can undergo intramolecular cycloaddition with good yields and minimal side reactions. For example, this Click-mediated cyclization may be used to stabilize an α -helical secondary structure when azide and alkyne are located in side-chains at positions *i* and *i*+4, respectively.

Example for a Protocol for Click Reactions in Peptide Synthesis:



Successful protocols have been published applying to 3 μmol peptide in 4 mL $t\text{BuOH}/\text{H}_2\text{O}$ (1:2) with excess of ascorbic acid (40 μmol) and $\text{CuSO}_4 \cdot 5 \text{H}_2\text{O}$ (40 μmol) generating Cu(I) *in situ*. Stirring at room temperature overnight is followed by appropriate chromatographic work up. [Le Chevalier-Isaad *et al.*, *Eur. J. Org. Chem.* 2010]

References:







- *Synthesis and conformational analysis of a cyclic peptide obtained via i to i+4 intramolecular side-chain to side-chain azide-alkyne 1,3-dipolar cycloaddition*; S. Cantel, C. Isaad Ale, M. Scrima, J. J. Levy, R. D. DiMarchi, P. Rovero, J. A. Halperin, A. M. D'Ursi, A. M. Papini, M. Chorev; **J Org Chem** 2008; **73**: 5663-74.
<https://doi.org/10.1021/jo800142s>
- *Side chain-to-side chain cyclization by click reaction*; A. Le Chevalier Isaad, A. M. Papini, M. Chorev, P. Rovero; **J. Pept. Sci.** 2009; **15**: 451-4. <https://doi.org/10.1002/psc.1141>
- *CuI-Catalyzed Azide-Alkyne Intramolecular i-to-(i+4) Side-Chain-to-Side-Chain Cyclization Promotes the Formation of Helix-Like Secondary Structures*; M. Scrima, A. Le Chevalier-Isaad, P. Rovero, A. M. Papini, M. Chorev, A. M. D'Ursi; **Eur. J. Org. Chem.** 2010; **2010**: 446-457. <https://doi.org/10.1002/ejoc.200901157>
- *Improved synthesis and biological evaluation of chelator-modified alpha-MSH analogs prepared by copper-free click chemistry*; N. J. Baumhover, M. E. Martin, S. G. Parameswarappa, K. C. Kloepping, M. S. O'Dorisio, F. C. Pigge, M. K. Schultz; **Bioorg Med Chem Lett** 2011; **21**: 5757-61. <https://doi.org/10.1016/j.bmcl.2011.08.017>
- *„Click“-cyclized (68)Ga-labeled peptides for molecular imaging and therapy: synthesis and preliminary in vitro and in vivo evaluation in a melanoma model system*; M. E. Martin, M. Sue O'Dorisio, W. M. Leverich, K. C. Kloepping, S. A. Walsh, M. K. Schultz; **Recent Results Cancer Res** 2013; **194**: 149-75. https://doi.org/10.1007/978-3-642-27994-2_9







2.3. Azido Amino Acids and Related Derivatives

Azido-Alkyl/Aryl Acids and Alcohols







		Product details
RL-4430 N₃-C₇H₁₄-COOH 8-azido-octanoic acid CAS-No. 217180-76-2 Formula C ₈ H ₁₅ N ₃ O ₂ Mol. weight 185,23 g/mol		
BNN1370 N₃-EDA-Suc-OH Azido-ethylenediamine-succinoyl-OH CAS-No. 2225891-73-4 Formula C ₆ H ₁₀ N ₄ O ₃ Mol. weight 186,17 g/mol		







[↑ back to content](#)

			Product details
RL-4350	N₃-DAPr-Suc-OH		
Azido-propylenediamine-succinoyl-OH			
CAS-No.	929894-58-6		
Formula	C ₇ H ₁₂ N ₄ O ₃		
Mol. weight	200,20 g/mol		
			
RL-4360	N₃-DABu-Suc-OH		
Azido-butylenediamine-succinoyl-OH			
CAS-No.	2226183-50-0		
Formula	C ₈ H ₁₄ N ₄ O ₃		
Mol. weight	214,23 g/mol		
			
HAA6990	N₃-Aca-Aca-OH		
6-(6-azidohexanamido)hexanoic acid			
CAS-No.	866363-71-5		
Formula	C ₁₂ H ₂₂ N ₄ O ₃		
Mol. weight	270,33 g/mol		
			
RL-3650	N₃-Phenylpropionic-OH		
3-(4-azidophenyl)propanoic acid			
CAS-No.	103489-31-2		
Formula	C ₉ H ₉ N ₃ O ₂		
Mol. weight	191,19 g/mol		
			
AAA2190	DAPOA*DCHA		
2-(1,3-diazidopropan-2-yloxy)acetic acid dicyclohexylamine			
CAS-No.	2389064-43-9		
Formula	C ₅ H ₈ N ₆ O ₃ *C ₁₂ H ₂₃ N		
Mol. weight	200,16*181,32 g/mol		
			
HAA2230	N₃-1,4-cis-CHC-OH		
cis-4-Azidocyclohexanecarboxylic acid			
CAS-No.	863222-21-3		
Formula	C ₇ H ₁₁ N ₃ O		
Mol. weight	169,18 g/mol		
			





			Product details
HAA2235	N₃-1,4-<i>trans</i>-CHC-OH		
<i>trans</i> -4-Azidocyclohexanecarboxylic acid			
CAS-No.	1931895-14-5		
Formula	C ₇ H ₁₁ N ₃ O ₂		
Mol. weight	169,18 g/mol		
			
HAA2240	N₃-<i>trans</i>-MCHC-OH		
<i>trans</i> -4-(Azidomethyl)cyclohexanecarboxylic acid			
CAS-No.	170811-10-6		
Formula	C ₈ H ₁₃ N ₃ O ₂		
Mol. weight	183,21 g/mol		
			
HAA3175	N₃-HABA*DCHA (2S)		
(S)-4-azido-2-hydroxybutyric acid dicyclohexylamine			
CAS-No.	959148-55-1		
Formula	C ₄ H ₇ N ₃ O ₃ *C ₁₂ H ₂₃ N		
Mol. weight	145,12*181,32 g/mol		
			
HAA3365	N₃-IsoSer*DCHA (2S)		
(S)-2-Hydroxy-3-azidopropanoic acid dicyclohexylamine			
CAS-No.	1620171-65-4		
Formula	C ₃ H ₅ N ₃ O ₃ *C ₁₂ H ₂₃ N		
Mol. weight	131,09*181,32 g/mol		
			
HAA2245	N₃-PhAc-OH		
(4-Azidophenyl)acetic acid			
CAS-No.	62893-37-2		
Formula	C ₈ H ₇ N ₃ O ₂		
Mol. weight	177,16 g/mol		
			
AAA1970	N₃-Pen-OH		
5-Azido-pentanoic acid			
CAS-No.	79583-98-5		
Formula	C ₅ H ₉ N ₃ O ₂		
Mol. weight	143,14 g/mol		
			

[The Click Reaction](#)
[Amino Acid Derivatives and Related Building Blocks for Click Chemistry](#)
[Spermines and Amines for Click Chemistry](#)
[Click Reagents for Drug Delivery](#)
[Click Chemistry Tools for Proteomics](#)
[Carbohydrates for Click Chemistry](#)
[Proteolysis Targeting Chimeras \(PROTACs*\)](#)
[Index](#)
[↑ back to content](#)


			Product details
AAA1960	N₃-Hx-OH		
6-Azido-hexanoic acid			
CAS-No.	79598-53-1		
Formula	C ₆ H ₁₁ N ₃ O ₂		
Mol. weight	157,17 g/mol		
			
RL-2980	N₃-Aca-OSu		
6-Azidocaproic acid N-hydroxysuccinimidyl ester			
CAS-No.	866363-70-4		
Formula	C ₁₀ H ₁₄ N ₄ O ₄		
Mol. weight	254,24 g/mol		
			
RL-3480	8-Azido-octanoyl-OSu		
8-Azidooctanoic acid N-hydroxysuccinimide ester			
CAS-No.	2576471-56-0		
Formula	C ₁₂ H ₁₈ N ₄ O ₄		
Mol. weight	282,30 g/mol		
			
RL-3200	11-Azidoundecanoic acid		
11-Azido-undecanoic acid			
CAS-No.	118162-45-1		
Formula	C ₁₁ H ₂₁ N ₃ O ₂		
Mol. weight	227,30 g/mol		
			
RL-3170	11-Azido-undecanoyl-OSu		
11-Azidoundecanoic acid N-hydroxysuccinimide ester			
CAS-No.	850080-13-6		
Formula	C ₁₅ H ₂₄ N ₄ O ₄		
Mol. weight	324,38 g/mol		
			
RL-3210	12-Azidododecanoic acid		
12-Azido-dodecanoic acid			
CAS-No.	80667-36-3		
Formula	C ₁₂ H ₂₃ N ₃ O ₂		
Mol. weight	241,33 g/mol		
			







			Product details
RL-3220	12-Azido-dodecanoyl-OSu		
12-Azidododecanoic acid N-hydroxysuccinimide ester			
CAS-No.	2489524-00-5		
Formula	C ₁₆ H ₂₆ N ₄ O ₄		
Mol. weight	338,40 g/mol		
RL-3230	14-Azido-myristic acid		
14-azidotetradecanoic acid			
CAS-No.	176108-61-5		
Formula	C ₁₄ H ₂₇ N ₃ O ₂		
Mol. weight	269,38 g/mol		
RL-3240	16-Azido-palmitic acid		
16-azidohexadecanoic acid			
CAS-No.	112668-54-9		
Formula	C ₁₆ H ₃₁ N ₃ O ₂		
Mol. weight	297,44 g/mol		
RL-3250	18-Azido-stearic acid		
18-azidooctadecanoic acid			
CAS-No.	1529763-58-3		
Formula	C ₁₈ H ₃₅ N ₃ O ₂		
Mol. weight	325,49 g/mol		
RL-2995	4-(Azidomethyl)benzoic acid		
4-Azidomethylbenzoic acid			
CAS-No.	79584-03-5		
Formula	C ₈ H ₇ N ₃ O ₂		
Mol. weight	177,16 g/mol		
RL-2035	ATFB		
4-Azido-2,3,5,6-tetrafluorobenzoic acid			
CAS-No.	122590-77-6		
Formula	C ₇ HF ₄ N ₃ O ₂		
Mol. weight	235,1 g/mol		

[The Click Reaction](#)
[Amino Acid Derivatives and Related Building Blocks for Click Chemistry](#)
[Spermines and Amines for Click Chemistry](#)
[Click Reagents for Drug Delivery](#)
[Click Chemistry Tools for Proteomics](#)
[Carbohydrates for Click Chemistry](#)
[Proteolysis Targeting Chimeras \(PROTACs[®]\)](#)
[Index](#)
[↑ back to content](#)





		Product details
RL-2045	ATFB-NHS	
N-Succinimidyl 4-azido-2,3,5,6-tetrafluorobenzoate		
CAS-No.	126695-58-7	
Formula	C ₁₁ H ₄ F ₄ N ₄ O ₄	
Mol. weight	332,17 g/mol	
PEG5000	N₃-TFBA-O₂Oc	
{2-[2-(4-Azido-2,3,5,6-tetrafluorobenzoyl-amino)ethoxy]ethoxy}acetic acid		
CAS-No.	1993119-45-1	
Formula	C ₁₃ H ₁₂ F ₄ N ₄ O ₅	
Mol. weight	380,25 g/mol	
RL-2990	4-Azidobenzyl alcohol	
(4-azidophenyl)methanol		
CAS-No.	31499-54-4	
Formula	C ₇ H ₇ N ₃ O	
Mol. weight	149,15 g/mol	
RL-4070	2-Azidobenzyl alcohol	
(2-azidophenyl)methanol		
CAS-No.	20615-76-3	
Formula	C ₇ H ₇ N ₃ O	
Mol. weight	149,15 g/mol	

Azido-Alanine and Propionic Acid Derivatives

		Product details
HAA1880	H-L-Aza-OH*HCl hydrate	
(S)-2-Amino-3-azidopropanoic acid hydrochloride hydrate		
CAS-No.	1620171-64-3	
Formula	C ₃ H ₆ N ₄ O ₂ *HCl*nH ₂ O	
Mol. weight	130,11*36,45 g/mol	

		Product details
HAA1885	H-D-Aza-OH*HCl hydrate	
(R)-2-Amino-3-azidopropanoic acid hydrochloride hydrate		
CAS-No.	1379690-01-3	 
Formula	C ₃ H ₆ N ₄ O ₂ *HCl*nH ₂ O	
Mol. weight	130,11*36,45 g/mol	
BAA1820	Boc-L-Aza-OH*CHA	
(S)-2-t-Butyloxycarbonylamino-3-azidopropanoic acid cyclohexylamine		
CAS-No.	2098496-88-7	 
Formula	C ₈ H ₁₆ N ₄ O ₆ *C ₆ H ₁₃ N	
Mol. weight	230,22*99,18 g/mol	
BAA1825	Boc-D-Aza-OH*CHA	
(R)-2-t-Butyloxycarbonylamino-3-azidopropanoic acid cyclohexylamine		
CAS-No.	2098496-96-7	 
Formula	C ₈ H ₁₆ N ₄ O ₆ *C ₆ H ₁₃ N	
Mol. weight	230,22*99,18 g/mol	
FAA1820	Fmoc-L-Aza-OH (solv.)	
(S)-2-(9-Fluorenylmethyloxycarbonylamino)-3-azidopropanoic acid, solvate with DIPE		
CAS-No.	684270-46-0	 
Formula	C ₁₈ H ₁₆ N ₄ O ₄	
Mol. weight	352,34 g/mol	
FAA6870	Fmoc-D-Aza-OH	
(R)-2-(9-Fluorenylmethyloxycarbonylamino)-3-azidopropanoic acid		
CAS-No.	1016163-79-3	 
Formula	C ₁₈ H ₁₆ N ₄ O ₄	
Mol. weight	352,34 g/mol	
FAA9420	Fmoc-L-MeDap(N₃)-OH	
(S)-2-((((9H-fluoren-9-yl)methoxy)carbonyl)(methyl)amino)-3-azidopropanoic acid		
CAS-No.	1263721-08-9	 
Formula	C ₁₉ H ₁₈ N ₄ O ₄	
Mol. weight	366,38 g/mol	






[The Click Reaction](#)
[Amino Acid Derivatives and Related Building Blocks for Click Chemistry](#)
[Spermines and Amines for Click Chemistry](#)
[Click Reagents for Drug Delivery](#)
[Click Chemistry Tools for Proteomics](#)
[Carbohydrates for Click Chemistry](#)
[Proteolysis Targeting Chimeras \(PROTACs[®]\)](#)
[Index](#)
[↑ back to content](#)

		Product details
HAA2130 N₃-L-Dap(Boc)-OH (S)-2-Azido-3-((t-butyloxycarbonyl)amino)propanoic acid CAS-No. 1932432-15-9 Formula C ₈ H ₁₆ N ₄ O ₄ Mol. weight 230,22 g/mol		
HAA2135 N₃-D-Dap(Boc)-OH (R)-2-Azido-3-((t-butyloxycarbonyl)amino)propanoic acid CAS-No. 1630044-08-4 Formula C ₈ H ₁₆ N ₄ O ₄ Mol. weight 230,22 g/mol		
HAA2140 N₃-L-Dap(Fmoc)-OH (S)-2-Azido-3-[(9-fluorenylmethyloxycarbonyl)amino]propanoic acid CAS-No. 880637-82-1 Formula C ₁₈ H ₁₆ N ₄ O ₄ Mol. weight 352,34 g/mol		
HAA2145 N₃-D-Dap(Fmoc)-OH (R)-2-Azido-3-[(9-fluorenylmethyloxycarbonyl)amino]propanoic acid CAS-No. 1807631-13-5 Formula C ₁₈ H ₁₆ N ₄ O ₄ Mol. weight 352,34 g/mol		

References:

- Azidoalanine mutagenicity in *Salmonella*: effect of homologation and alpha-methyl substitution; J. B. Mangold, M. R. Mischke, J. M. LaVelle; **Mutat Res** 1989; **216**: 27-33. [https://doi.org/10.1016/0165-1161\(89\)90020-4](https://doi.org/10.1016/0165-1161(89)90020-4)
- „Click to chelate“: synthesis and installation of metal chelates into biomolecules in a single step; T. L. Mindt, H. Struthers, L. Brans, T. Anguelov, C. Schweinsberg, V. Maes, D. Tourwe, R. Schibli; **J Am Chem Soc** 2006; **128**: 15096-7. <https://doi.org/10.1021/ja066779f>
- Design, synthesis, and biological activity of novel triazole amino acids used to probe binding interactions between ligand and neutral amino acid transport protein SN1; M. Gajewski, B. Seaver, C. S. Esslinger; **Bioorg Med Chem Lett** 2007; **17**: 4163-6. <https://doi.org/10.1016/j.bmcl.2007.05.061>
- Peptide tertiary structure nucleation by side-chain crosslinking with metal complexation and double „click“ cycloaddition; O. Torres, D. Yuksel, M. Bernardina, K. Kumar, D. Bong; **ChemBiochem** 2008; **9**: 1701-5. <https://doi.org/10.1002/cbic.200800040>
- Maintaining biological activity by using triazoles as disulfide bond mimetics; K. Holland-Nell, M. Meldal; **Angew. Chem. Int. Ed.** 2011; **50**: 5204-6. <https://doi.org/10.1002/anie.201005846>

Azido-Homoalanine and Azido-Butanoic Acid Derivatives

		Product details
HAA5730 H-L-Aha-OH*HCl (S)-2-Amino-4-azidobutanoic acid hydrochloride CAS-No. 942518-29-8 Formula $C_4H_8N_4O_2 \cdot HCl$ Mol. weight 144,13*36,45 g/mol		
HAA9280 H-L-Aha-OH 4-Azido-L-homoalanine CAS-No. 120042-14-0 Formula $C_4H_8N_4O_2$ Mol. weight 144,13 g/mol		
HAA1630 H-D-Aha-OH*HCl (R)-2-Amino-4-azidobutanoic acid hydrochloride CAS-No. 1858224-26-6 Formula $C_4H_8N_4O_2 \cdot HCl$ Mol. weight 144,13*36,45 g/mol		
BAA1800 Boc-L-Aha-OH*CHA (S)-2-t-Butyloxycarbonylamino-4-azidobutanoic acid cyclohexylamine CAS-No. 120042-08-2net Formula $C_9H_{16}N_4O_4 \cdot C_6H_{13}N$ Mol. weight 244,25*99,18 g/mol		
BAA1805 Boc-D-Aha-OH*CHA (R)-2-t-Butyloxycarbonylamino-4-azidobutanoic acid cyclohexylamine CAS-No. 1609202-75-6 net Formula $C_9H_{16}N_4O_4 \cdot C_6H_{13}N$ Mol. weight 244,25*99,18 g/mol		

[↑ back to content](#)

The Click Reaction

Amino Acid Derivatives and Related Building Blocks for Click Chemistry

Spermines and Amines for Click Chemistry







Click Reagents for Drug Delivery

Click Chemistry Tools for Proteomics

Carbohydrates for Click Chemistry

Proteolysis Targeting Chimeras (PROTACs[®])




Index

		Product details	
FAA6620	Fmoc-L-Aha-OH		
(S)-2-(9-Fluorenylmethyloxycarbonylamino)-4-azido-butanoic acid			
CAS-No.	942518-20-9		
Formula	C ₁₉ H ₁₈ N ₄ O ₄		
Mol. weight	366,41 g/mol		
FAA6810	Fmoc-D-Aha-OH		
(R)-2-(9-Fluorenylmethyloxycarbonylamino)-4-azido-butanoic acid			
CAS-No.	1263047-53-5		
Formula	C ₁₉ H ₁₈ N ₄ O ₄		
Mol. weight	366,41 g/mol		
ZAA5700	Z-L-Aha-OH*DCHA		
(S)-2-Benzoyloxycarbonylamino-4-azidobutanoic acid dicyclohexylamine			
CAS-No.	1263047-43-3 net		
Formula	C ₁₂ H ₁₄ N ₄ O ₄ *C ₁₂ H ₂₃ N		
Mol. weight	278,26*181,34 g/mol		
HAA2150	N₃-L-Dab(Boc)-OH		
(S)-2-Azido-4-((t-butyloxycarbonyl)amino)butanoic acid			
CAS-No.	1932403-71-8		
Formula	C ₉ H ₁₆ N ₄ O ₄		
Mol. weight	244,25 g/mol		
HAA2155	N₃-D-Dab(Boc)-OH		
(R)-2-Azido-4-((t-butyloxycarbonyl)amino)butanoic acid			
CAS-No.	1922891-74-4		
Formula	C ₉ H ₁₆ N ₄ O ₄		
Mol. weight	244,25 g/mol		
HAA3170	N₃-L-Dab(Fmoc)-OH		
(S)-2-Azido-4-[(9-fluorenylmethyloxycarbonyl)amino]butanoic acid			
CAS-No.	2250436-44-1		
Formula	C ₉ H ₁₆ N ₄ O ₄		
Mol. weight	366,37 g/mol		



References:

- EP2190856 (A2), US2008096819 (A1), WO2009026393 (A2), WO2009026393 (A3), WO2006079364
- Presentation and detection of azide functionality in bacterial cell surface proteins; A. J. Link, M. K. Vink, D. A. Tirrell; **J Am Chem Soc** 2004; **126**: 10598-602. <https://doi.org/10.1021/ja047629c>
- Mild and chemoselective peptide-bond cleavage of peptides and proteins at azido homoalanine; J. W. Back, O. David, G. Kramer, G. Masson, P. T. Kasper, L. J. de Koning, L. de Jong, J. H. van Maarseveen, C. G. de Koster; **Angew. Chem. Int. Ed.** 2005; **44**: 7946-50. <https://doi.org/10.1002/anie.200502431>
- Expanding the diversity of chemical protein modification allows post-translational mimicry; S. I. van Kasteren, H. B. Kramer, H. H. Jensen, S. J. Campbell, J. Kirkpatrick, N. J. Oldham, D. C. Anthony, B. G. Davis; **Nature** 2007; **446**: 1105-9. <https://doi.org/10.1038/nature05757>
- Site-specific modification of *Candida antarctica* lipase B via residue-specific incorporation of a non-canonical amino acid; S. Schoffelen, M. H. Lammermon, M. B. van Eldijk, J. C. van Hest; **Bioconjug Chem** 2008; **19**: 1127-31. <https://doi.org/10.1021/bc800019v>
- Unnatural amino acid incorporation into virus-like particles; E. Strable, D. E. Prasuhn, Jr., A. K. Udit, S. Brown, A. J. Link, J. T. Ngo, G. Lander, J. Quispe, C. S. Potter, B. Carragher, D. A. Tirrell, M. G. Finn; **Bioconjug Chem** 2008; **19**: 866-75. <https://doi.org/10.1021/bc700390r>
- Selective enrichment of azide-containing peptides from complex mixtures; M. A. Nessen, G. Kramer, J. Back, J. M. Baskin, L. E. Smeenk, L. J. de Koning, J. H. van Maarseveen, L. de Jong, C. R. Bertozzi, H. Hiemstra, C. G. de Koster; **J Proteome Res** 2009; **8**: 3702-11. <https://doi.org/10.1021/pr900257z>
- „Clickable“ elastins: elastin-like polypeptides functionalized with azide or alkyne groups; R. L. Teeuwen, S. S. van Berkel, T. H. van Dulmen, S. Schoffelen, S. A. Meeuwissen, H. Zuilhof, F. A. de Wolf, J. C. van Hest; **Chem Commun (Camb)** 2009; 4022-4. <https://doi.org/10.1039/b903903a>




Azido-beta-Homoalanine




		Product details
HAA3970 H-L-Dbu(N₃)-OH*HCl (S)-3-Amino-4-azidobutanoic acid hydrochloride CAS-No. 2389078-78-6 net Formula C ₄ H ₈ N ₄ O ₂ *HCl Mol. weight 144,13*36,45 g/mol		
FAA2035 Fmoc-L-Dbu(N₃)-OH (S)-3-(9-Fluorenylmethyloxycarbonyl)amino-4-azido-butanoic acid CAS-No. 934502-72-4 Formula C ₁₉ H ₁₈ N ₄ O ₄ Mol. weight 366,37 g/mol		
FAA3650 Fmoc-D-Dbu(N₃)-OH (R)-3-(9-Fluorenylmethyloxycarbonyl)amino-4-azido-butanoic acid CAS-No. 1932023-47-6 Formula C ₁₉ H ₁₈ N ₄ O ₄ Mol. weight 366,37 g/mol		

[↑ back to content](#)

		Product details
ZAA1290	Z-L-Dbu(N₃)-OH	
(S)-3-(Benzyloxycarbonyl)amino-4-azido-butanoic acid		
CAS-No.	1932657-23-2	
Formula	C ₁₂ H ₁₄ N ₄ O ₄	
Mol. weight	278,26 g/mol	
		
ZAA1285	Z-D-Dbu(N₃)-OH	
(R)-3-(Benzyloxycarbonyl)amino-4-azido-butanoic acid		
CAS-No.	1931958-82-5	
Formula	C ₁₂ H ₁₄ N ₄ O ₄	
Mol. weight	278,26 g/mol	
		

2-Amino-3-Azido-Butanoic Acid

		Product details
HAA3010	H-Abu(3-N₃)-OH*HCl (2S,3S)	
(2S,3S)-2-amino-3-azidobutanoic acid hydrochloride		
CAS-No.	2737202-68-3	
Formula	C ₄ H ₈ N ₄ O ₂ *HCl	
Mol. weight	144,13*36,45 g/mol	
		
HAA3020	H-Abu(3-N₃)-OH*HCl (2S,3R)	
(2S,3R)-2-amino-3-azidobutanoic acid hydrochloride		
CAS-No.	2737202-63-8	
Formula	C ₄ H ₈ N ₄ O ₂ *HCl	
Mol. weight	144,13*36,45 g/mol	
		
FAA2040	Fmoc-Abu(3-N₃)-OH (2S,3S)	
(2S,3S)-2-(9-Fluorenylmethyloxycarbonyl)amino-3-azido-butanoic acid		
CAS-No.	131669-42-6	
Formula	C ₁₉ H ₁₈ N ₄ O ₄	
Mol. weight	366,37 g/mol	
		

		Product details	
FAA3200	Fmoc-Abu(3-N₃)-OH (2S,3R)		
(2S,3R)-2-(9-Fluorenylmethyloxycarbonyl)amino-3-azido-butanoic acid			
CAS-No.	146306-79-8		
Formula	C ₁₉ H ₁₈ N ₄ O ₄		
Mol. weight	366,37 g/mol		
FAA3540	Fmoc-Abu(3-N₃)-OH (2R,3S)		
(2R,3S)-2-(9-Fluorenylmethyloxycarbonyl)amino-3-azido-butanoic acid			
CAS-No.	1932349-21-7		
Formula	C ₁₉ H ₁₈ N ₄ O ₄		
Mol. weight	366,37 g/mol		
FAA2095	Fmoc-Abu(3-N₃)-OH (2R,3R)		
(2R,3R)-2-(9-Fluorenylmethyloxycarbonyl)amino-3-azido-butanoic acid			
CAS-No.	1229394-75-5		
Formula	C ₁₉ H ₁₈ N ₄ O ₄		
Mol. weight	366,37 g/mol		



**Interested in Click Chemistry
and automated synthesis?**

Watch the recording of our workshop!



The Click Reaction

Amino Acid Derivatives and Related
Building Blocks for Click Chemistry

Spermines and Amines for Click Chemistry

Click Reagents for Drug Delivery

Click Chemistry Tools for Proteomics

Carbohydrates for Click Chemistry

Proteolysis Targeting Chimeras (PROTACs[®])

Index

[↑ back to content](#)

Azido-Masked Amino Function

Azido groups located in amino acid side-chains can be used for various applications. 2-Amino-3-azido-butanoic acid is shown as an example below.




Fig. 11: Possible applications of amino acids bearing azide groups in the side-chain.

A) Azido groups can be used for any type of Click conjugation with any available alkynyl residue forming conjugates with peptides or any other organic molecule.

B) Azido groups may also be used for another prominent type of bioconjugation, namely the Staudinger ligation, which is a further development of the Staudinger reaction. The Staudinger ligation is characterized by high selectivity and a typically rapid and high-yielding turnover. As a biorthogonal reaction, it has been used for the semisynthesis of proteins, for installing posttranslational modifications such as glycosylations, and for DNA labeling.


C) The azido group can be reduced to an amino function and hereby serve as a masked amino group. Prominent methods for the reduction of azido groups include the Staudinger reaction as well as the reduction by DTT. Azido groups are stable towards treatment with piperidine (Fmoc deprotection), $\text{Pd}(0)$ (Alloc removal) and acidic treatment (cleavage of Mtt, Trt or other acid-sensitive groups). However, as it is a pseudohalogenide, care must be taken during coupling steps, as HATU will cause a high degree of racemization. This can be avoided using collidine or other non-nucleophilic bases instead of DIPEA.

Chiral α,β -diamines and diamino acids have increasingly become motifs of interest in organic synthesis owing to their ubiquity in natural products and medicinal agents. For example, these motifs are found in biotin, penicillins, and the antiinfluenza neuraminidase inhibitor Tamiflu. Chiral vicinal diamines and their metal complexes have been employed in stereoselective organic synthesis, in particular as chiral auxiliaries and ligands in catalytic asymmetric synthesis.



References:

- Über neue organische Phosphorverbindungen III. Phosphinmethylderivate und Phosphinimine; H. Staudinger, J. Meyer; **Helv. Chim. Acta** 1919; **2**: 635-646. <https://doi.org/10.1002/hlca.19190020164>
- Reduction of aryl azides by thiols: Implications for the use of photoaffinity reagents; J. V. Staros, H. Bayley, D. N. Standring, J. R. Knowles; **Biochem. Biophys. Res. Commun.** 1978; **80**: 568-572. [https://doi.org/10.1016/0006-291X\(78\)91606-6](https://doi.org/10.1016/0006-291X(78)91606-6)
- Sixty years of Staudinger reaction; Y. G. Gololobov, I. N. Zhmurova, L. F. Kasukhin; **Tetrahedron** 1981; **37**: 437-472. [https://doi.org/10.1016/S0040-4020\(01\)92417-2](https://doi.org/10.1016/S0040-4020(01)92417-2)
- Staudinger ligation: a peptide from a thioester and azide; B. L. Nilsson, L. L. Kiessling, R. T. Raines; **Org Lett** 2000; **2**: 1939-41. <https://doi.org/10.1021/ol0060174>
- A „traceless“ Staudinger ligation for the chemoselective synthesis of amide bonds; E. Saxon, J. I. Armstrong, C. R. Bertozzi; **Org Lett** 2000; **2**: 2141-3. <https://doi.org/10.1021/ol006054v>
- Alpha,beta-diamino acids: biological significance and synthetic approaches; A. Viso, R. Fernandez de la Pradilla, A. Garcia, A. Flores; **Chem Rev** 2005; **105**: 3167-96. <https://doi.org/10.1021/cr0406561>
- Chemoselective peptide cyclization by traceless Staudinger ligation; R. Kleineweischede, C. P. Hackenberger; **Angew. Chem. Int. Ed.** 2008; **47**: 5984-8. <https://doi.org/10.1002/anie.200801514>
- Protein engineering with the traceless Staudinger ligation; A. Tam, R. T. Raines; **Methods Enzymol.** 2009; **462**: 25-44. [https://doi.org/10.1016/S0076-6879\(09\)62002-4](https://doi.org/10.1016/S0076-6879(09)62002-4)
- Bioconjugation via azide-Staudinger ligation: an overview; C. I. Schilling, N. Jung, M. Biskup, U. Schepers, S. Bräse; **Chem Soc Rev** 2011; **40**: 4840-71. <https://doi.org/10.1039/c0cs00123f>
- Update 1 of: alpha,beta-diamino acids: biological significance and synthetic approaches; A. Viso, R. Fernandez de la Pradilla, M. Tortosa, A. Garcia, A. Flores; **Chem Rev** 2011; **111**: PR1-42. <https://doi.org/10.1021/cr100127y>
- The Staudinger Ligation; C. Bednarek, I. Wehl, N. Jung, U. Schepers, S. Bräse; **Chem Rev** 2020; **120**: 4301-4354. <https://doi.org/10.1021/acs.chemrev.9b00665>
- Site-selective traceless Staudinger ligation for glycoprotein synthesis reveals scope and limitations; G. J. Bernardes, L. Linderth, K. J. Doores, O. Boutureira, B. G. Davis; **ChemBiochem** 2011; **12**: 1383-6. <https://doi.org/10.1002/cbic.201100125>




Azido-Citrulline





		Product details
HAA4980 N₃-L-Cit-OH*DCHA (S)-2-Azido-citrulline dicyclohexylamine CAS-No. 1799421-66-1 net Formula C ₆ H ₁₁ N ₅ O ₃ *C ₁₂ H ₂₃ N Mol. weight 201,18*181,32 g/mol		 

Azido-Cysteine


		Product details
HAA2810	N₃-L-Cys(Trt)-OH*CHA	
(R)-2-azido-3-(tritylthio)propanoic acid cyclohexylamine		
CAS-No.	1286670-90-3	
Formula	C ₂₂ H ₁₉ N ₃ O ₂ S*C ₆ H ₁₃ N	
Mol. weight	389,47*99,17 g/mol	
		
FAA9460	Fmoc-L-Cys(Azidoethyl)-OH	
N-(((9H-fluoren-9-yl)methoxy)carbonyl)-S-(2-azidoethyl)-L-cysteine		
CAS-No.	3029657-80-2	
Formula	C ₂₀ H ₂₀ N ₄ O ₄ S	
Mol. weight	412,46 g/mol	
		

Azido-Glycine


		Product details
BAA4895	Boc-Aeg(N₃)-OK	
N-(2-azidoethyl)-N-(tert-butoxycarbonyl)glycine potassium salt		
CAS-No.	2172548-14-8 net	
Formula	C ₉ H ₁₅ KN ₄ O ₄	
Mol. weight	282,34 g/mol	
		
FAA4060	Fmoc-Aeg(N₃)-OH	
N-(9-Fluorenylmethoxycarbonyl)-N-(2-azidoethyl)glycine		
CAS-No.	1935981-35-3	
Formula	C ₁₉ H ₁₈ N ₄ O ₄	
Mol. weight	366,37 g/mol	
		
FAA4055	Fmoc-Abg(N₃)-OH	
N-(9-Fluorenylmethoxycarbonyl)-N-(4-azidobutyl)glycine		
CAS-No.	2250433-81-7	
Formula	C ₂₁ H ₂₂ N ₄ O ₄	
Mol. weight	394,42 g/mol	
		

		Product details
HAA2850	N₃-Gly-Gly-OH*DCHA Azido-glycylglycine dicyclohexylamine CAS-No. 855750-87-7 net Formula C ₄ H ₆ N ₄ O ₃ *C ₁₂ H ₂₃ N Mol. weight 158,12*181,32 g/mol	 
HAA2840	N₃-Gly-Gly-Gly-OH Azido-glycylglycylglycine CAS-No. 1993176-75-2 Formula C ₆ H ₉ N ₅ O ₄ Mol. weight 215,17 g/mol	 
HAA2860	N₃-Gly-Gly-Gly-Gly-Gly-OH CAS-No. 2250433-77-1 Formula C ₁₀ H ₁₅ N ₇ O ₆ Mol. weight 329,27 g/mol	 
HAA9330	N₃-Gly-Aeg(Fmoc)-OH N-(2-((((9H-fluoren-9-yl)methoxy)carbonyl)amino)ethyl)-N-(2-azidoacetyl)glycine CAS-No. 2606227-07-8 Formula C ₂₁ H ₂₁ N ₅ O ₅ Mol. weight 423,43 g/mol	 





Azido-Leucine







		Product details
HAA3350	N₃-L-Leu-OH*BHA (S)-2-azido-4-methylpentanoic acid benzhydrylamine salt CAS-No. 79410-33-6 Formula C ₆ H ₁₁ N ₃ O ₂ *C ₁₃ H ₁₃ N Mol. weight 157,17*183,25 g/mol	 

[The Click Reaction](#)
[Amino Acid Derivatives and Related Building Blocks for Click Chemistry](#)
[Spermines and Amines for Click Chemistry](#)
[Click Reagents for Drug Delivery](#)
[Click Chemistry Tools for Proteomics](#)
[Carbohydrates for Click Chemistry](#)
[Proteolysis Targeting Chimeras \(PROTACs[®]\)](#)
[Index](#)
[↑ back to content](#)

	Product details
HAA2820 N₃-L-Leu-OH*CHA (S)-2-Azido-4-methylpentanoic acid cyclohexylamine CAS-No. 1286670-79-8 Formula C₆H₁₁N₃O₂*C₆H₁₃N Mol. weight 157,17*99,18 g/mol	 

Azido-Lysine

	Product details
HAA9210 H-L-Lys(N₃)-OH N-epsilon-azido-L-lysine CAS-No. 159610-92-1 Formula C₆H₁₂N₄O₂ Mol. weight 172,19 g/mol	 
HAA1625 H-L-Lys(N₃)-OH*HCl N-epsilon-Azido-L-lysine hydrochloride CAS-No. 1454334-76-9 Formula C₆H₁₂N₄O₂*HCl Mol. weight 172,19*36,45 g/mol	 
HAA9340 H-L-Lys(N₃-Gly)-OH*HCl Azidoacetyl-L-Lysine hydrochloride CAS-No. 1198617-82-1 net Formula C₈H₁₅N₅O₃ Mol. weight 229,24 g/mol	 
FAA8855 Fmoc-L-Lys(N₃-Gly)-OH Azidoacetyl-Fmoc-L-Lysine CAS-No. 1198617-89-8 Formula C₂₃H₂₅N₅O₅ Mol. weight 451,48 g/mol	 


		Product details
HAA1890	H-D-Lys(N₃)-OH*HCl	
N-epsilon-Azido-D-lysine hydrochloride		
CAS-No.	2098497-01-7	 
Formula	C ₆ H ₁₂ N ₄ O ₂ *HCl	
Mol. weight	172,19*36,45 g/mol	
HAA2080	H-L-Lys(EO-N₃)-OH*HCl	
(S)-2-amino-6-((2-azidoethoxy)carbonylamino)hexanoic acid hydrochloride		
CAS-No.	1994331-17-7	 
Formula	C ₉ H ₁₇ N ₅ O ₆ *HCl	
Mol. weight	259,26*36,46 g/mol	
HAA9295	H-L-Lys(COCF₂N₃)-OH*HCl	
N6-(2-azido-2,2-difluoroacetyl)-L-lysine		
Formula	C ₈ H ₁₃ F ₂ N ₅ O ₃ *HCl	 
Mol. weight	265,22*36,46 g/mol	
HAA2870	H-(Gly)₃-Lys(N₃)-OH*HCl	
Triglycyl-epsilon-azido-L-lysine hydrochloride		
CAS-No.	2250437-45-5	 
Formula	C ₁₂ H ₂₁ N ₇ O ₅ *HCl	
Mol. weight	343,34*36,45 g/mol	
BAA4900	Boc-L-Lys(N₃)-OK	
Boc-azidolysine potassium salt		
CAS-No.	846549-33-5	 
Formula	C ₁₁ H ₁₉ KN ₄ O ₄	
Mol. weight	310,40 g/mol	
BAA1810	Boc-L-Lys(N₃)-OH*CHA	
N-alpha-t-Butyloxycarbonyl-epsilon-azido-L-lysine cyclohexylamine		
CAS-No.	2098497-30-2	 
Formula	C ₁₁ H ₂₀ N ₄ O ₄ *C ₆ H ₁₃ N	
Mol. weight	272,30*99,18 g/mol	

[The Click Reaction](#)
[Amino Acid Derivatives and Related Building Blocks for Click Chemistry](#)
[Spermines and Amines for Click Chemistry](#)
[Click Reagents for Drug Delivery](#)
[Click Chemistry Tools for Proteomics](#)
[Carbohydrates for Click Chemistry](#)
[Proteolysis Targeting Chimeras \(PROTACs[®]\)](#)
[Index](#)
[↑ back to content](#)

BAA1815 Boc-D-Lys(N₃)-OH*CHA


N-alpha-*t*-Butyloxycarbonyl-epsilon-azido-D-lysine
cyclohexylamine

CAS-No. 1858224-39-1
Formula C₁₁H₂₀N₄O₄*C₆H₁₃N
Mol. weight 272,30*99,18 g/mol


FAA1793 Fmoc-L-Lys(N₃)-OH


N-alpha-(9-Fluorenylmethyloxycarbonyl)-epsilon-azido-L-lysine

CAS-No. 159610-89-6
Formula C₂₁H₂₂N₄O₄
Mol. weight 394,42 g/mol


FAA1835 Fmoc-D-Lys(N₃)-OH


N-alpha-(9-Fluorenylmethyloxycarbonyl)-epsilon-azido-D-lysine

CAS-No. 1198791-53-5
Formula C₂₁H₂₂N₄O₄
Mol. weight 394,42 g/mol


FAA7925 Fmoc-L-Lys(N₃-AEEA)-OH


N2-(((9H-fluoren-9-yl)methoxy)carbo-
nyl)-N6-(2-(2-(2-azidoethoxy)ethoxy)acetyl)-L-lysine

CAS-No. 1236293-83-6
Formula C₂₇H₃₃N₅O₇
Mol. weight 539,59 g/mol


FAA8145 Fmoc-L-Lys(N₃-Aca-DIM)-OH


N-alpha-(9-Fluorenylmethyloxycarbonyl)-N-epsilon-
[6-azido-1-(4,4-dimethyl-2,6-dioxocyclohexylidene)
hexyl]-L-lysine







CAS-No. 2408993-39-3
Formula C₃₅H₄₃N₅O₆
Mol. weight 629,76 g/mol


FAA8595 Fmoc-L-MeLys(N₃)-OH

N-alpha-(9-Fluorenylmethyloxycarbonyl)-N-al-
pha-methyl-epsilon-azido-L-lysine

CAS-No. 1263721-14-7
Formula C₂₂H₂₄N₄O₄
Mol. weight 408,46 g/mol



		Product details
HAA2170	N₃-L-Lys(Boc)-OH (S)-2-Azido-6-[(t-butyloxycarbonyl)amino]hexanoic acid CAS-No. 333366-32-8 Formula C ₁₁ H ₂₀ N ₄ O ₄ Mol. weight 272,3 g/mol	 
HAA2175	N₃-D-Lys(Boc)-OH (R)-2-Azido-6-[(t-butyloxycarbonyl)amino]hexanoic acid CAS-No. 1178899-92-7 Formula C ₁₁ H ₂₀ N ₄ O ₄ Mol. weight 272,3 g/mol	 
HAA2160	N₃-L-Lys(Fmoc)-OH (S)-2-Azido-6-[(9-fluorenylmethyloxycarbonyl)amino]hexanoic acid CAS-No. 473430-12-5 Formula C ₂₁ H ₂₂ N ₄ O ₄ Mol. weight 394,42 g/mol	 
HAA2165	N₃-D-Lys(Fmoc)-OH (R)-2-Azido-6-[(9-fluorenylmethyloxycarbonyl)amino]hexanoic acid CAS-No. 1994300-35-4 Formula C ₂₁ H ₂₂ N ₄ O ₄ Mol. weight 394,42 g/mol	 
HAA2880	N₃-L-Lys(Mtt)-OH (S)-2-Azido-6-[(4-methyltrityl)amino]hexanoic acid CAS-No. 1333231-26-7 Formula C ₂₆ H ₂₈ N ₄ O ₂ Mol. weight 428,53 g/mol	 
FAA8825	Fmoc-L-Lys(COCF₂N₃)-OH N2-(((9H-fluoren-9-yl)methoxy)carbonyl)-N6-(2-azido-2,2-difluoroacetyl)-L-lysine Formula C ₂₃ H ₂₃ F ₂ N ₅ O ₅ Mol. weight 487,46 g/mol	 

[The Click Reaction](#)
[Amino Acid Derivatives and Related Building Blocks for Click Chemistry](#)
[Spermines and Amines for Click Chemistry](#)
[Click Reagents for Drug Delivery](#)
[Click Chemistry Tools for Proteomics](#)
[Carbohydrates for Click Chemistry](#)
[Proteolysis Targeting Chimeras \(PROTACs[®]\)](#)
[Index](#)
[↑ back to content](#)


FAA8830 Fmoc-L-Lys(4-N₃-Z)-OH

N2-(((9H-fluoren-9-yl)methoxy)carbonyl)-N6-(((4-azido-benzyl)oxy)carbonyl)-L-lysine

CAS-No. 1446511-14-3

Formula C₂₉H₂₉N₅O₆

Mol. weight 543,58 g/mol


**HAA9315 H-L-Lys(4-N₃-Z)-OH*HCl**

(2S)-6-(4-Azido-benzylloxycarbonylamino)-2-amino-hexanoic acid hydrochloride

CAS-No. 2084913-49-3

Formula C₁₄H₁₉N₅O₄*HCl

Mol. weight 321,34*36,46 g/mol


**HAA9380 H-L-Lys(2-N₃-Z)-OH**

N6-(((2-azidobenzyl)oxy)carbonyl)-L-lysine

CAS-No. 1131963-69-3

Formula C₁₄H₁₉N₅O₄

Mol. weight 321,34 g/mol


**HAA9370 H-L-Lys(3-N₃-Z)-OH*HCl**

N6-(((3-azidobenzyl)oxy)carbonyl)-L-lysine

CAS-No. 2084913-47-1

Formula C₁₄H₁₉N₅O₄*HCl

Mol. weight 321,34*36,45 g/mol


**FAA8880 Fmoc-L-Lys(2-N₃-Z)-OH**

(2S)-6-(2-Azido-benzylloxycarbonylamino)-2-(9H-fluoren-9-ylmethoxycarbonylamino)-hexanoic acid

CAS-No. 2714331-96-9

Formula C₂₉H₂₉N₅O₆

Mol. weight 543,58 g/mol


**FAA8890 Fmoc-L-Lys(3-N₃-Z)-OH**

(2S)-6-(3-Azido-benzylloxycarbonylamino)-2-(9H-fluoren-9-ylmethoxycarbonylamino)-hexanoic acid






CAS-No. 1836202-27-7

Formula C₂₉H₂₉N₅O₆

Mol. weight 543,58 g/mol



Azido-Ornithine

		Product details
HAA1620 H-L-Orn(N₃)-OH*HCl N-delta-Azido-L-ornithine hydrochloride CAS-No. 156463-09-1n Formula C ₅ H ₁₀ N ₄ O ₂ *HCl Mol. weight 158,16*36,45 g/mol		
HAA1895 H-D-Orn(N₃)-OH*HCl N-delta-Azido-D-ornithine hydrochloride CAS-No. 1858224-08-4 Formula C ₅ H ₁₀ N ₄ O ₂ *HCl Mol. weight 158,16*36,45 g/mol		
BAA1830 Boc-L-Orn(N₃)-OH*CHA N-alpha-t-Butyloxycarbonyl-delta-azido-L-ornithine cyclohexylamine CAS-No. 763139-35-1n Formula C ₁₀ H ₁₈ N ₄ O ₄ *C ₆ H ₁₃ N Mol. weight 258,27*99,18 g/mol		
BAA1835 Boc-D-Orn(N₃)-OH*CHA N-alpha-t-Butyloxycarbonyl-delta-azido-D-ornithine cyclohexylamine CAS-No. 1858224-18-6 Formula C ₁₀ H ₁₈ N ₄ O ₄ *C ₆ H ₁₃ N Mol. weight 258,27*99,18 g/mol		
FAA6880 Fmoc-L-Orn(N₃)-OH N-alpha-(9-Fluorenylmethyloxycarbonyl)-delta-azido-L-ornithine CAS-No. 1097192-04-5 Formula C ₂₀ H ₂₀ N ₄ O ₄ Mol. weight 380,4 g/mol		


[The Click Reaction](#)
[Amino Acid Derivatives and Related Building Blocks for Click Chemistry](#)
[Spermines and Amines for Click Chemistry](#)
[Click Reagents for Drug Delivery](#)
[Click Chemistry Tools for Proteomics](#)
[Carbohydrates for Click Chemistry](#)
[Proteolysis Targeting Chimeras \(PROTACs[®]\)](#)
[Index](#)
[↑ back to content](#)

FAA6890 Fmoc-D-Orn(N₃)-OHN- α -(9-Fluorenylmethyloxycarbonyl)- δ -azido-D-ornithine

CAS-No. 1176270-25-9

Formula C₂₀H₂₀N₄O₄

Mol. weight 380,4 g/mol


**FAA8680 Fmoc-L-MeOrn(N₃)-OH**

(S)-2-((((9H-fluoren-9-yl)methoxy)carbonyl)(methyl)amino)-5-azidopentanoic acid

CAS-No. 2991255-09-3

Formula C₂₁H₂₂N₄O₄


Mol. weight 394,43 g/mol

**HAA2210 N₃-D-Orn(Boc)-OH*CHA**(R)-2-Azido-5-[(*t*-butoxycarbonyl)amino]pentanoic acid cyclohexylamine

CAS-No. 2165877-62-1

Formula C₁₀H₁₈N₄O₄*C₆H₁₃N


Mol. weight 258,27*99,18 g/mol

**HAA9485 N₃-L-Orn(Boc)-ONa**sodium (S)-2-azido-5-[(*tert*-butoxycarbonyl)amino]pentanoate

CAS-No. 1639198-67-6 net

Formula C₁₀H₁₇N₄NaO₄


Mol. weight 280,26 g/mol

**HAA9495 N₃-L-Orn(Boc)-OK**potassium (S)-2-azido-5-[(*tert*-butoxycarbonyl)amino]pentanoate

CAS-No. 1639198-67-6 net

Formula C₁₀H₁₇N₄KO₄

Mol. weight 296,37 g/mol


**HAA2225 N₃-L-Orn(Fmoc)-OH**

(S)-2-Azido-5-[(9-fluorenylmethyloxycarbonyl)amino]pentanoic acid

CAS-No. 1994267-98-9

Formula C₂₀H₂₀N₄O₄





Mol. weight 380,4 g/mol



Reference:

- Application of metal-free triazole formation in the synthesis of cyclic RGD-DTPA conjugates; S. S. van Berkel, A. Dirks, S. A. Meeuwissen, D. L. Pingen, O. C. Boerman, P. Laverman, F. L. van Delft, J. J. Cornelissen, F. P. Rutjes; *Chembiochem* 2008; **9**: 1805-15. <https://doi.org/10.1002/cbic.200800074>

Azido-Proline

		Product details
HAA2125 H-L-Pro(4-N ₃)-OH*HCl (2S,4S) (2S,4S)-4-Azidopyrrolidine-2-carboxylic acid hydrochloride CAS-No. 892128-58-4 Formula C ₅ H ₈ N ₄ O ₂ *HCl Mol. weight 156,14*36,45 g/mol		
HAA3150 H-L-Pro(4-N ₃)-OH*HCl (2S,4R) (2S,4R)-4-Azidopyrrolidine-2-carboxylic acid hydrochloride CAS-No. 1019849-13-8 net Formula C ₅ H ₈ N ₄ O ₂ *HCl Mol. weight 156,14*36,45 g/mol		
HAA3140 H-D-Pro(4-N ₃)-OH*HCl (2R,4S) (2R,4S)-4-Azidopyrrolidine-2-carboxylic acid hydrochloride CAS-No. 2137086-50-9 Formula C ₅ H ₈ N ₄ O ₂ *HCl Mol. weight 156,14*36,45 g/mol		
HAA3190 H-D-Pro(4-N ₃)-OH*HCl (2R,4R) (2R,4R)-4-Azidopyrrolidine-2-carboxylic acid hydrochloride CAS-No. 2737202-69-4 Formula C ₅ H ₈ N ₄ O ₂ *HCl Mol. weight 156,14*36,45 g/mol		


[↑ back to content](#)

BAA1905 Boc-L-Pro(4-N₃)-OH (2S,4S)*cis*-N- α -(*t*-Butyloxycarbonyl)-4-azido-L-proline

CAS-No. 132622-65-2

Formula C₁₀H₁₆N₄O₄


Mol. weight 256,26 g/mol

**BAA1930** Boc-L-Pro(4-N₃)-OH*DCHA (2S,4R)*trans*-N- α -(*t*-Butyloxycarbonyl)-4-azido-L-proline dicyclohexylamine

CAS-No. 132622-68-5 net

Formula C₁₀H₁₆N₄O₄*C₁₂H₂₃N


Mol. weight 256,26*181,32 g/mol

**BAA3110** Boc-D-Pro(4-N₃)-OH*DCHA (2R,4S)*trans*-N- α -(*t*-Butyloxycarbonyl)-4-azido-D-proline dicyclohexylamine

CAS-No. 132622-77-6 net

Formula C₁₀H₁₆N₄O₄*C₁₂H₂₃N


Mol. weight 256,26*181,32 g/mol

**BAA3120** Boc-D-Pro(4-N₃)-OH (2R,4R)*cis*-N- α -(*t*-Butyloxycarbonyl)-4-azido-D-proline

CAS-No. 650601-59-5

Formula C₁₀H₁₆N₄O₄


Mol. weight 256,26 g/mol

**FAA2050** Fmoc-L-Pro(4-N₃)-OH (2S,4S)*cis*-N- α -(9-Fluorenylmethyloxycarbonyl)-4-azido-L-proline

CAS-No. 263847-08-1

Formula C₂₀H₁₈N₄O₄


Mol. weight 378,38 g/mol



**FAA3000** Fmoc-L-Pro(4-N₃)-OH (2S,4R)*trans*-N- α -(9-Fluorenylmethyloxycarbonyl)-4-azido-L-proline

CAS-No. 702679-55-8




Formula C₂₀H₁₈N₄O₄







Mol. weight 378,38 g/mol







		Product details
FAA4630	Fmoc-D-Pro(4-N₃)-OH (2R,4S) <i>trans</i> -N- α -(9-Fluorenylmethyloxycarbonyl)-4-azido-D-proline	 
CAS-No.	2137142-63-1	
Formula	C ₂₀ H ₁₈ N ₄ O ₄	
Mol. weight	378,38 g/mol	
FAA4720	Fmoc-D-Pro(4-N₃)-OH (2R,4R) <i>cis</i> -N- α -(9-Fluorenylmethyloxycarbonyl)-4-azido-D-proline	 
CAS-No.	1378847-51-8	
Formula	C ₂₀ H ₁₈ N ₄ O ₄	
Mol. weight	378,38 g/mol	


Azido-Phenylalanine






		Product details
HAA1850	H-L-Phe(4-N₃)-OH 4-Azido-L-phenylalanine	 
CAS-No.	33173-53-4	
Formula	C ₉ H ₁₀ N ₄ O ₂	
Mol. weight	206,20 g/mol	
HAA2980	H-L-Phe(4-N₃)-OH*HCl 4-Azido-L-phenylalanine hydrochloride	 
CAS-No.	34670-43-4	
Formula	C ₉ H ₁₀ N ₄ O ₂ *HCl	
Mol. weight	206,2*36,45 g/mol	
HAA1855	H-D-Phe(4-N₃)-OH 4-Azido-D-phenylalanine	 
CAS-No.	1241681-80-0	
Formula	C ₉ H ₁₀ N ₄ O ₂	
Mol. weight	206,20 g/mol	

		Product details
HAA1856	H-D-Phe(4-N₃)-OH*HCl 4-Azido-D-phenylalanine hydrochloride CAS-No. 1241681-80-0 Formula C ₉ H ₁₀ N ₄ O ₂ *HCl Mol. weight 206,2*36,45 g/mol	 
HAA4090	H-L-Phe(4-CH₂-N₃)*HCl 4-azidomethyl-L-phenylalanine hydrochloride CAS-No. 1446772-80-0 Formula C ₁₀ H ₁₂ N ₄ O ₂ *HCl Mol. weight 220,23*36,45 g/mol	 
HAA9480	H-L-Phe(4-Azido-PrTz)-OH*TFA (S)-2-amino-3-(4-(6-(3-azidopropyl)-1,2,4,5-tetrazin-3-yl)phenyl)propanoic acid trifluoroacetic acid salt Formula C ₁₄ H ₁₆ N ₈ O ₂ *CF ₃ COOH Mol. weight 328,34*114,02 g/mol	 
BAA1850	Boc-L-Phe(4-N₃)-OH N-alpha-t-Butyloxycarbonyl-4-azido-L-phenylalanine CAS-No. 33173-55-6 Formula C ₁₄ H ₁₈ N ₄ O ₄ Mol. weight 306,32 g/mol	 
BAA1855	Boc-D-Phe(4-N₃)-OH N-alpha-t-Butyloxycarbonyl-4-azido-D-phenylalanine CAS-No. 214630-05-4 Formula C ₁₄ H ₁₈ N ₄ O ₄ Mol. weight 306,32 g/mol	 
BAA4660	Boc-L-Phe(4-CH₂-N₃)-OH N-alpha-t-Butyloxycarbonyl-4-azidomethyl-L-phenylalanine CAS-No. 205127-59-9 Formula C ₁₅ H ₂₀ N ₄ O ₄ Mol. weight 320,35 g/mol	 


			Product details
FAA1905	Fmoc-L-Phe(4-N₃)-OH		
N-alpha-(9-Fluorenylmethyloxycarbonyl)-4-azido-L-phenylalanine			
CAS-No.	163217-43-4		
Formula	C ₂₄ H ₂₀ N ₄ O ₄		
Mol. weight	428,44 g/mol		
			
FAA1910	Fmoc-D-Phe(4-N₃)-OH		
N-alpha-(9-Fluorenylmethyloxycarbonyl)-4-azido-D-phenylalanine			
CAS-No.	1391586-30-3		
Formula	C ₂₄ H ₂₀ N ₄ O ₄		
Mol. weight	428,44 g/mol		
			
FAA7740	Fmoc-L-Phe(4-CH₂-N₃)-OH		
N-alpha-(9-Fluorenylmethyloxycarbonyl)-4-azido-methyl-L-phenylalanine			
CAS-No.	2375587-79-2		
Formula	C ₂₅ H ₂₂ N ₄ O ₄		
Mol. weight	442,47 g/mol		
			
HAA3360	N₃-L-Phe-OH*DCHA		
(S)-2-Azido-3-phenylpropanoic acid dicyclohexylamine			
CAS-No.	79410-36-9		
Formula	C ₉ H ₉ N ₃ O ₂ *C ₁₃ H ₂₃ N		
Mol. weight	191,19*181,32 g/mol		
			

Azido-Tyrosine

			Product details
HAA3940	H-L-Tyr(3-N₃)-OH		
3-Azido-L-tyrosine			
CAS-No.	129960-90-3		
Formula	C ₉ H ₁₀ N ₄ O ₃		
Mol. weight	222,2 g/mol		
			

		Product details
HAA9310	H-L-Tyr(2-azidoethyl)-OH	
O-2-azidoethyl-tyrosine		
CAS-No.	1570523-47-5	
Formula	C ₁₁ H ₁₄ N ₄ O ₃	
Mol. weight	250,26 g/mol	
		
BAA4650	Boc-L-Tyr(2-azidoethyl)-OH	
N-alpha-t-Butyloxycarbonyl-O-(2-azidoethyl)-L-tyrosine		
CAS-No.	1434445-10-9	
Formula	C ₁₆ H ₂₂ N ₄ O ₅	
Mol. weight	350,38 g/mol	
		
BAA2235	Boc-L-Tyr(PEG(3)-N₃)-OH*DCHA	
N-alpha-t-Butyloxycarbonyl-O-(2-(2-(2-azidoethoxy)ethoxy)ethyl)-L-tyrosine dicyclohexylamine		
CAS-No.	1831059-64-3 net	
Formula	C ₂₀ H ₃₀ N ₄ O ₇ *C ₁₂ H ₂₃ N	
Mol. weight	438,47*181,32 g/mol	
		
FAA8535	Fmoc-L-Tyr(2-azidoethyl)-OH	
N-alpha-(9-Fluorenylmethyloxycarbonyl)-O-(2-azidoethyl)-L-tyrosine		
CAS-No.	1454816-10-4	
Formula	C ₂₆ H ₂₄ N ₄ O ₅	
Mol. weight	472,50 g/mol	
		
HAA9215	H-L-Tyr(2-azidoethyl)-OH*HCl	
O-(2-azidoethyl)-L-tyrosine hydrochloride		
CAS-No.	1567845-62-8	
Formula	C ₁₁ H ₁₄ N ₄ O ₃ *HCl	
Mol. weight	250,26*36,46 g/mol	
		


Azido-Valine







	Product details
HAA9285 N₃-L-Val-OH*CHA Azido-L-valine cyclohexylammonium salt CAS-No. 1217462-63-9 Formula C ₅ H ₉ N ₃ O ₂ *C ₆ H ₁₃ N Mol. weight 143,15*99,17 g/mol	 







References:

- Improved solid-phase peptide synthesis method utilizing alpha-azide-protected amino acids; J. T. t. Lundquist, J. C. Pelletier; **Org Lett** 2001; **3**: 781-3. <https://doi.org/10.1021/ol0155485>
- Azide reduction during peptide cleavage from solid support-the choice of thioscavenger?; P. E. Schneggenburger, B. Worbs, U. Diederichsen; **J. Pept. Sci.** 2010; **16**: 10-4. <https://doi.org/10.1002/psc.1202>
- A DOTA-peptide conjugate by copper-free click chemistry; M. E. Martin, S. G. Parameswarappa, M. S. O'Dorisio, F. C. Pigge, M. K. Schultz; **Bioorg Med Chem Lett** 2010; **20**: 4805-7. <https://doi.org/10.1016/j.bmcl.2010.06.111>
- CuI-Catalyzed Azide-Alkyne Intramolecular i-to-(i+4) Side-Chain-to-Side-Chain Cyclization Promotes the Formation of Helix-Like Secondary Structures; M. Scrima, A. Le Chevalier-Isaad, P. Rovero, A. M. Papini, M. Chorev, A. M. D'Ursi; **Eur. J. Org. Chem.** 2010; **2010**: 446-457. <https://doi.org/10.1002/ejoc.200901157>
- Synthesis and conformational analysis of a cyclic peptide obtained via i to i+4 intramolecular side-chain to side-chain azide-alkyne 1,3-dipolar cycloaddition; S. Cantel, C. Isaad Ale, M. Scrima, J. J. Levy, R. D. DiMarchi, P. Rovero, J. A. Halperin, A. M. D'Ursi, A. M. Papini, M. Chorev; **J Org Chem** 2008; **73**: 5663-74. <https://doi.org/10.1021/jo800142s>
- Side chain-to-side chain cyclization by click reaction; A. Le Chevalier Isaad, A. M. Papini, M. Chorev, P. Rovero; **J. Pept. Sci.** 2009; **15**: 451-4. <https://doi.org/10.1002/psc.1141>
- An efficient peptide ligation using azido-protected peptides via the thioester method; H. Katayama, H. Hojo, T. Ohira, Y. Nakahara; **Tetrahedron Lett.** 2008; **49**: 5492-5494. <https://doi.org/10.1016/j.tetlet.2008.07.037>

Azido-Linkers and Azido-ADC Linkers

	Product details
RL-4100 Azido-SS-COOH 3-((2-azidoethyl)disulfanyl)propanoic acid CAS-No. 2228857-32-5 Formula C ₅ H ₉ N ₃ O ₂ S ₂ Mol. weight 207,27 g/mol	 

			Product details
RL-3320	Azido-Pen-SS-COOH		
3-((2-(5-azidopentanamido)ethyl)disulfanyl)propanoic acid			
CAS-No.	2576471-47-9		
Formula	C ₁₀ H ₁₈ N ₄ O ₃ S ₂		
Mol. weight	306,40 g/mol		
HNN1090	N₃-Cystamine*HCl		
Azido-cystamine hydrochloride			
CAS-No.	1807512-40-8 net		
Formula	C ₄ H ₁₀ N ₄ S ₂ *HCl		
Mol. weight	178,28*36,45 g/mol		
RL-4120	Biotin-SS-N₃		
N-(2-((2-azidoethyl)disulfanyl)ethyl)-5-((3aS,4S,6aR)-2-oxohexahydro-1H-thieno[3,4-d]imidazol-4-yl)pentanamide			
CAS-No.	1620523-64-9		
Formula	C ₁₄ H ₂₄ N ₆ O ₂ S ₃		
Mol. weight	404,57 g/mol		
RL-4150	Azido-SS-OpNC		
2-((2-azidoethyl)disulfanyl)ethyl (4-nitrophenyl) carbonate			
CAS-No.	2766027-28-3		
Formula	C ₁₁ H ₁₂ N ₄ O ₅ S ₂		
Mol. weight	344,36 g/mol		
RL-3280	N₃-Pen-Dde		
2-(5-azido-1-hydroxypentylidene)-5,5-dimethylcyclohexane-1,3-dione			
CAS-No.	1867129-38-1		
Formula	C ₁₃ H ₁₉ N ₃ O ₃		
Mol. weight	265,31 g/mol		
RL-3290	N₃-Pen-Dtpp		
5-(5-azido-1-hydroxypentylidene)-1,3-dimethylpyrimidine-2,4,6(1H,3H,5H)-trione			
CAS-No.	1867129-42-7		
Formula	C ₁₁ H ₁₅ N ₅ O ₄		
Mol. weight	281,27 g/mol		

		Product details
ADC1680	6-Azidohexanoyl-Val-Ala-PAB-Cl 6-Azidohexanoyl-valyl-alanyl-(4-aminobenzyl chloride) CAS-No. 2706565-03-7 Formula $C_{21}H_{31}ClN_6O_3$ Mol. weight 450,97 g/mol	 
ADC1710	Azido-PEG(4)-Val-Ala-PAB-Cl Azido-tetraethyleneglycol-propanoyl-valyl-alanyl-(4-aminobenzyl chloride) Formula $C_{26}H_{41}ClN_6O_7$ Mol. weight 585,10 g/mol	 
ADC1290	6-Azidohexanoyl-Val-Ala-PAB 6-azidohexanoyl-valyl-alanyl-(4-aminobenzyl alcohol) CAS-No. 2706564-30-7 Formula $C_{21}H_{32}N_6O_4$ Mol. weight 432,52 g/mol	 
ADC1300	6-Azidohexanoyl-Val-Ala-PAB-PNP 6-azidohexanoyl-valyl-alanyl-(4-aminobenzyl)-(4-nitrophenyl)-carbonate Formula $C_{28}H_{35}N_7O_8$ Mol. weight 597,62 g/mol	 
ADC1120	6-Azidohexanoyl-Val-Cit-PAB 6-azidohexanoyl-valyl-citrullyl-(4-aminobenzyl alcohol) CAS-No. 1613321-02-0 Formula $C_{24}H_{38}N_6O_5$ Mol. weight 518,61 g/mol	 
ADC1130	6-Azidohexanoyl-Val-Cit-PAB-PNP 6-azidohexanoyl-valyl-citrullyl-(4-aminobenzyl)-(4-nitrophenyl)-carbonate CAS-No. 1613321-01-9 Formula $C_{31}H_{41}N_9O_9$ Mol. weight 683,71 g/mol	 

[The Click Reaction](#)
[Amino Acid Derivatives and Related Building Blocks for Click Chemistry](#)
[Spermines and Amines for Click Chemistry](#)
[Click Reagents for Drug Delivery](#)
[Click Chemistry Tools for Proteomics](#)
[Carbohydrates for Click Chemistry](#)
[Proteolysis Targeting Chimeras \(PROTACs[®]\)](#)
[Index](#)
[↑ back to content](#)


ADC1580 Azido-cyclobutane-1,1-dicarboxamide-Ala-PAB

3-azidopropyl-cyclobutane-1,1-dicarboxamide-allyl-(4-aminobenzyl alcohol)

CAS-No. 2576471-45-7

Formula $C_{19}H_{26}N_6O_4$

Mol. weight 402,45 g/mol


**ADC1590 Azido-cyclobutane-1,1-dicarboxamide-Ala-PAB-PNP**

3-azidopropyl-cyclobutane-1,1-dicarboxamide-allyl-(4-aminobenzyl)-(4-nitrophenyl)-carbonate

CAS-No. 2576471-36-6

Formula $C_{26}H_{29}N_7O_8$

Mol. weight 567,55 g/mol


**ADC1480 Azido-cyclobutane-1,1-dicarboxamide-Cit-PAB**

3-azidopropyl-cyclobutane-1,1-dicarboxamide-citryl-(4-aminobenzyl alcohol)

CAS-No. 2576471-33-3

Formula $C_{22}H_{32}N_6O_5$

Mol. weight 488,54 g/mol


**ADC1490 Azido-cyclobutane-1,1-dicarboxamide-Cit-PAB-PNP**

3-azidopropyl-cyclobutane-1,1-dicarboxamide-citryl-(4-aminobenzyl)-(4-nitrophenyl)-carbonate

CAS-No. 2576471-44-6

Formula $C_{29}H_{35}N_7O_9$


Mol. weight 653,64 g/mol

**ADC1330 Azido-PEG(4)-Val-Ala-PAB**

azido-tetraethyleneglycol-propanoyl-valyl-allyl-(4-aminobenzyl alcohol)

Formula $C_{26}H_{42}N_6O_8$


Mol. weight 566,65 g/mol



**ADC1340 Azido-PEG(4)-Val-Ala-PAB-PNP**

azido-tetraethyleneglycol-propanoyl-valyl-allyl-(4-aminobenzyl)-(4-nitrophenyl)-carbonate


Formula $C_{33}H_{45}N_7O_{12}$

Mol. weight 731,75 g/mol




		Product details
ADC1160	Azido-PEG(4)-Val-Cit-PAB	
azido-tetraethyleneglycol-propanoyl-valyl-citrul-lyl-(4-aminobenzyl alcohol)		
CAS-No.	2055024-64-9	 
Formula	C ₂₉ H ₄₈ N ₉ O ₉	
Mol. weight	652,74 g/mol	
ADC1170	Azido-PEG(4)-Val-Cit-PAB-PNP	
azido-tetraethyleneglycol-propanoyl-valyl-citrul-lyl-(4-aminobenzyl)-(4-nitrophenyl)-carbonate		
CAS-No.	1869126-60-2	 
Formula	C ₃₆ H ₅₁ N ₉ O ₁₃	
Mol. weight	817,84 g/mol	






Other Azido Building Blocks

		Product details
HNN1110	N₃-TEMPO	
4-Azido-2,2,6,6-tetramethyl-1-piperidinyloxy		
CAS-No.	63697-61-0	 
Formula	C ₉ H ₁₇ N ₄ O	
Mol. weight	197,26 g/mol	






2.4. Alkyne Amino Acids and Related Derivatives

Propargylating Reagents







		Product details
RL-8670	5HP2O-alkyne	
5-hydroxy-5-methyl-1-(prop-2-yn-1-yl)-1,5-dihydro-2H-pyrrol-2-one		
CAS-No.	2484704-61-0	 
Formula	C ₈ H ₉ NO ₂	
Mol. weight	151,16 g/mol	

			Product details
RL-3945	Mal-Alkyne		
3-(2,5-dioxo-2,5-dihydro-1H-pyrrol-1-yl)-N-(prop-2-yn-1-yl)propanamide			
CAS-No.	548777-19-1		
Formula	C ₁₀ H ₁₀ N ₂ O ₃		
Mol. weight	206,20 g/mol		
PEG2755	Propargyl amine		
CAS-No.	2450-71-7		
Formula	C ₃ H ₅ N		
Mol. weight	55,08 g/mol		
MAA1100	Mal-AMCHC-N-Propargylamide		
<i>trans</i> -4-[(2,5-dioxo-2,5-dihydro-1H-pyrrol-1-yl)methyl]-N-(prop-2-yn-1-yl)cyclohexane-1-carboxamide			
CAS-No.	2027476-42-0		
Formula	C ₁₅ H ₁₈ N ₂ O ₃		
Mol. weight	274,32 g/mol		
LS-4310	DecarboxyBiotin-Alkyne		
(3aS,4S,6aR)-4-(6-heptyn-1-yl)tetrahydro-1H-Thieno[3,4-d]imidazol-2(3H)-one			
CAS-No.	887915-53-9		
Formula	C ₁₂ H ₁₈ N ₂ O ₃ S		
Mol. weight	238,35 g/mol		
RL-3490	Biotin-Propargylamide		
Biotinyl-N-propargylamide			
CAS-No.	773888-45-2		
Formula	C ₁₃ H ₁₉ N ₃ O ₂ S		
Mol. weight	281,37 g/mol		


Alkyne-Alkyl Acids and Alkyne-Aryl Acids

		Product details
RL-8660	Poc-Aoa 2-((((prop-2-yn-1-yloxy)carbonyl)amino)oxy)acetic acid CAS-No. 1877011-25-0 Formula $C_6H_7NO_5$ Mol. weight 173,12 g/mol	 
RL-3330	Alkyne-SS-COOH 3-((2-pent-4-ynamidoethyl)disulfanyl)propanoic acid CAS-No. 2279938-29-1 Formula $C_{10}H_{15}NO_3S_2$ Mol. weight 261,36 g/mol	 
RL-4110	DBCO-Suc-SS-COOH CAS-No. 2749426-25-1 Formula $C_{26}H_{24}N_2O_4S_2$ Mol. weight 468,59 g/mol	 
RL-3410	Photo-Click-Heptanoic acid 2-(3-(but-3-ynyl)-3H-diazirin-3-yl)acetic acid CAS-No. 2049109-24-0 Formula $C_7H_8N_2O_2$ Mol. weight 152,15 g/mol	 
RL-2055	Alkyne-myristic acid 13-Tetradecynoic acid CAS-No. 82909-47-5 Formula $C_{14}H_{24}O_2$ Mol. weight 224,34 g/mol	 




[The Click Reaction](#)
[Amino Acid Derivatives and Related Building Blocks for Click Chemistry](#)
[Spermines and Amines for Click Chemistry](#)
[Click Reagents for Drug Delivery](#)
[Click Chemistry Tools for Proteomics](#)
[Carbohydrates for Click Chemistry](#)
[Proteolysis Targeting Chimeras \(PROTACs[®]\)](#)
[Index](#)
[↑ back to content](#)




			Product details
RL-2060	Alkyne-palmitic acid		
15-Hexadecynoic acid			
CAS-No.	99208-90-9		
Formula	$C_{16}H_{28}O_2$		
Mol. weight	252,39 g/mol		
RL-3720	Photo-Click-Palmitic acid		
11-(3-(but-3-yn-1-yl)-3H-diazirin-3-yl)undecanoic acid			
CAS-No.	2988174-40-7		
Formula	$C_{16}H_{26}N_2O_2$		
Mol. weight	278,40 g/mol		
RL-2065	Alkyne-stearic acid		
17-Octadecynoic acid			
CAS-No.	34450-18-5		
Formula	$C_{18}H_{32}O_2$		
Mol. weight	280,45 g/mol		
RL-3760	Photo-Click-Stearic acid		
13-(3-(but-3-yn-1-yl)-3H-diazirin-3-yl)tridecanoic acid			
CAS-No.	2989205-54-9		
Formula	$C_{18}H_{30}N_2O_2$		
Mol. weight	306,45 g/mol		
PEG1935	Propargyl-NHS		
3-(Prop-2-ynyloxy)propanoic acid succinimidyl ester			
CAS-No.	1174157-65-3		
Formula	$C_{10}H_{11}NO_5$		
Mol. weight	225,2 g/mol		
RL-3460	10-Undecynoyl-OSu		
10-Undecynoic acid N-hydroxysuccinimide ester			
CAS-No.	1006592-57-9		
Formula	$C_{15}H_{21}NO_4$		
Mol. weight	279,34 g/mol		

		Product details
RL-3780	Fmoc-N-propargyl-MPBA	
4-(4-(((9-Fluorenylmethyloxycarbonyl)(propargyl)amino)methyl)-3-methoxyphenoxy)butanoic acid		
CAS-No.	1009362-00-8	
Formula	C ₃₀ H ₂₉ NO ₆	
Mol. weight	499,20 g/mol	



Propargylalanine and Propargyl-Propionic Acid Derivatives

		Product details
FAA4170	Fmoc-L-Dap(Pentynoyl)-OH	
N-alpha-(9-Fluorenylmethyloxycarbonyl)-N-beta-(4-pentynoyl)-L-2,3-diaminopropionic acid		
CAS-No.	2250436-47-4	
Formula	C ₂₃ H ₂₂ N ₂ O ₅	
Mol. weight	406,43 g/mol	
		
FAA2080	Fmoc-alpha-Prg-D-Ala-OH	
(S)-2-[(9-Fluorenylmethyloxycarbonyl)amino]-2-methyl-4-pentynoic acid		
CAS-No.	1198791-58-0	
Formula	C ₂₁ H ₁₉ NO ₄	
Mol. weight	349,38 g/mol	
		
FAA4230	Fmoc-L-Dap(Poc)-OH	
N-alpha-(9-Fluorenylmethyloxycarbonyl)-N-beta-propargyloxycarbonyl-L-2,3-diaminopropionic acid		
CAS-No.	2250437-44-4	
Formula	C ₂₂ H ₂₀ N ₂ O ₆	
Mol. weight	408,41 g/mol	
		




Propargylarginine

FAA7400 Fmoc-L-Arg(Propargyl,Pbf)-OH

N-alpha-(9-Fluorenylmethoxycarbonyl)-N'-(2,2,4,6,7-pentamethyldihydrobenzofuran)-N''-propargyl-5-sulfonyl-L-arginine

Formula $C_{37}H_{42}N_4O_7S$
Mol. weight 686,82 g/mol



Product details




Propargylcysteine

HAA2350 H-L-Cys(Propargyl)-OH*HCl

S-Propargyl-L-cysteine hydrochloride

CAS-No. 3262-64-4
Formula $C_6H_9NO_2S \cdot HCl$
Mol. weight 159,21*36,45 g/mol




Product details



BAA2250 Boc-L-Cys(Propargyl)-OH*DCHA

N-alpha-t-Butyloxycarbonyl-S-propargyl-L-cysteine dicyclohexylamine


CAS-No. 1260119-25-2 net
Formula $C_{11}H_{17}NO_4S \cdot C_{12}H_{23}N$
Mol. weight 259,32*181,32 g/mol



FAA3810 Fmoc-L-Cys(Propargyl)-OH

N-alpha-(9-Fluorenylmethoxycarbonyl)-S-propargyl-L-cysteine






CAS-No. 1354752-76-3
Formula $C_{21}H_{19}NO_4S$
Mol. weight 381,44 g/mol








References:

- Photoinduced addition of glycosyl thiols to alkynyl peptides: use of free-radical thiol-yne coupling for post-translational double-glycosylation of peptides; M. Lo Conte, S. Pacifico, A. Chambery, A. Marra, A. Dondoni; *J Org Chem* 2010; **75**: 4644-7. <https://doi.org/10.1021/jo1008178>
- Neoglycopeptides through direct functionalization of cysteine; C. Vala, F. Chrétien, E. Balentova, S. Lamandé-Langle and Y. Chapleur; *Tetrahedron Letters* 2011; **52**: 17-20. <https://doi.org/10.1016/j.tetlet.2010.10.021>

Propargylglycine

		Product details
BAA3230 Boc-N-(propargyl)-glycine N-alpha-t-Butyloxycarbonyl-N-alpha-propargyl-glycine CAS-No. 158979-29-4 Formula C ₁₀ H ₁₂ NO ₄ Mol. weight 213,23 g/mol		
FAA4950 Fmoc-N-(propargyl)-glycine N-alpha-(9-Fluorenylmethyloxycarbonyl)-N-alpha-propargyl-glycine CAS-No. 1033622-38-6 Formula C ₂₀ H ₁₇ NO ₄ Mol. weight 335,35 g/mol		
HAA7151 H-L-Pra-OH L-Propargylglycine CAS-No. 23235-01-0 Formula C ₅ H ₇ NO ₂ Mol. weight 113,11 g/mol		
HAA6490 H-D-Pra-OH*HCl D-Propargylglycine hydrochloride CAS-No. 23235-03-2 Formula C ₅ H ₈ NO ₂ *HCl Mol. weight 113,11*36,45 g/mol		
HAA9445 H-Sar-Sar-NH-Propargyl*HCl Di-sarcosine-propargylamide hydrochloride Formula C ₉ H ₁₅ N ₃ O ₂ *HCl Mol. weight 197,24*36,45 g/mol		

		Product details
BAA1434	Boc-L-Pra-OH*DCHA	
N-alpha-(t-Butyloxycarbonyl)-L-propargylglycine dicyclohexylamine		
CAS-No.	63039-49-6	
Formula	$C_{10}H_{15}NO_4 \cdot C_{12}H_{23}N$	
Mol. weight	213,23*181,32 g/mol	
		
BAA1377	Boc-D-Pra-OH*DCHA	
N-alpha-t-Butyloxycarbonyl-D-propargylglycine dicyclohexylamine		
CAS-No.	63039-47-4	
Formula	$C_{10}H_{15}NO_4 \cdot C_{12}H_{23}N$	
Mol. weight	213,23*181,32 g/mol	
		
FAA1589	Fmoc-L-Pra-OH	
N-alpha-(9-Fluorenylmethyloxycarbonyl)-L-propargylglycine		
CAS-No.	198561-07-8	
Formula	$C_{20}H_{17}NO_4$	
Mol. weight	335,35 g/mol	
		
FAA1690	Fmoc-D-Pra-OH	
N-alpha-(9-Fluorenylmethyloxycarbonyl)-D-propargylglycine		
CAS-No.	220497-98-3	
Formula	$C_{20}H_{17}NO_4$	
Mol. weight	335,36 g/mol	
		
WAA6025	Fmoc-L-Pra-Wang Resin	
Fmoc-L-Propargylglycine-Wang Resin		
Mesh Size	100-200 mesh	
DVB	1% DVB	
		


Propargyllysine

Product details

HAA2085 H-Lys(Pentynoyl)-OH*HCl

(S)-2-Amino-6-(pent-4-ynamido)hexanoic acid hydrochloride


CAS-No. 1167421-22-8 net
Formula $C_{11}H_{18}N_2O_3 \cdot HCl$
Mol. weight 226,27 *36,5 g/mol



HAA9440 H-Lys(Pentynoyl)-OH

N6-(pent-4-ynoyl)-L-lysine


CAS-No. 1167421-22-8
Formula $C_{11}H_{18}N_2O_3$
Mol. weight 226,28 g/mol



HAA2090 H-Lys(Poc)-OH*HCl

(S)-Amino-6-((prop-2-ynyloxy)carbonylamino)hexanoic acid hydrochloride


CAS-No. 1428330-91-9
Formula $C_{10}H_{16}N_2O_4 \cdot HCl$
Mol. weight 228,25*36,45 g/mol



HAA9390 H-Lys(CO-Ethoxypropargyl)-OH*HCl

(2S)-2-amino-6-([2-(prop-2-yn-1-yloxy)ethoxy]carbonylamino)hexanoic acid


Formula $C_{12}H_{20}N_2O_5 \cdot HCl$
Mol. weight 272,30*36,45 g/mol



BAA1960 Boc-L-Lys(Poc)-OH

(S)-2-(t-Butyloxycarbonylamino)-6-((prop-2-ynyloxy)carbonylamino)hexanoic acid

CAS-No. 1202704-91-3
Formula $C_{15}H_{24}N_2O_6$
Mol. weight 328,36 g/mol



The Click Reaction

Amino Acid Derivatives and Related Building Blocks for Click Chemistry

Spermines and Amines for Click Chemistry

Click Reagents for Drug Delivery

Click Chemistry Tools for Proteomics

Carbohydrates for Click Chemistry

Proteolysis Targeting Chimeras (PROTACs[®])

Index

[↑ back to content](#)


FAA4175 Fmoc-L-Lys(pentynoyl)-OH

N-alpha-(9-Fluorenylmethyloxycarbonyl)-N-epsilon-lon-(4-pentynoyl)-L-lysine

CAS-No. 1159531-18-6

Formula $C_{26}H_{28}N_2O_5$

Mol. weight 448,51 g/mol



FAA8115 Fmoc-L-Lys(Pentynoyl-DIM)-OH

N-alpha-(9-Fluorenylmethyloxycarbonyl)-N-epsilon-lon-[1-(4,4-dimethyl-2,6-dioxocyclohexylidene)pent-4-yn-1-yl]-L-lysine

CAS-No. 2408993-33-7

Formula $C_{34}H_{38}N_2O_6$

Mol. weight 570,69 g/mol



FAA8135 Fmoc-D-Lys(pentynoyl)-OH

N-alpha-(9-Fluorenylmethyloxycarbonyl)-N-epsilon-lon-(4-pentynoyl)-D-lysine

CAS-No. 2576508-18-2

Formula $C_{26}H_{28}N_2O_5$

Mol. weight 448,51 g/mol



FAA3150 Fmoc-L-Lys(Pryoc)-OH

(S)-2-((9-Fluorenylmethyloxy)amino)-6-((prop-2-yn-1-oxycarbonylamino)hexanoic acid

CAS-No. 1584133-25-4

Formula $C_{25}H_{26}N_2O_6$

Mol. weight 450,48 g/mol



FAA9565 Fmoc-D-Lys(Pryoc)-OH

N2-(((9H-fluoren-9-yl)methoxy)carbonyl)-N6-((prop-2-yn-1-yloxy)carbonyl)-D-lysine

CAS-No. 2991236-42-9

Formula $C_{25}H_{26}N_2O_6$

Mol. weight 450,49 g/mol



FAA8995 Fmoc-L-Lys(Hexynoyl)-OH


N2-(((9H-fluoren-9-yl)methoxy)carbonyl)-N6-(hex-5-ynoyl)-L-lysine

CAS-No. 1219440-73-9





Formula $C_{27}H_{30}N_2O_5$



Mol. weight 462,55 g/mol




	Product details
FAA8905 Fmoc-L-Lys(CO-Ethoxypropargyl)-OH (2S)-2-({[(9H-fluoren-9-yl)methoxy]carbonyl}amino)-6-({[2-(prop-2-yn-1-yloxy)ethoxy]carbonyl}amino)hexanoic acid Formula $C_{27}H_{30}N_2O_7$ Mol. weight 494,54 g/mol	 

Propargylphenylalanine


	Product details
HAA4970 H-L-Phe(4-NH-Poc)-OH*HCl 4-(Propargyloxycarbonyl)amino-L-phenylalanine hydrochloride CAS-No. 2576508-14-8 Formula $C_{13}H_{14}N_2O_4 \cdot HCl$ Mol. weight 262,26*36,45 g/mol	 
HAA9220 H-L-Phe(4-Eth)-OH*HCl 4-Ethynyl-L-phenylalanine hydrochloride CAS-No. 188640-63-3 Formula $C_{11}H_{11}NO_2 \cdot HCl$ Mol. weight 189,21*36,46 g/mol	 
BAA3980 Boc-L-Phe(4-NH-Poc)-OH N-alpha-t-Butyloxycarbonyl-4-(propargyloxycarbonyl)amino-L-phenylalanine CAS-No. 2576508-03-5 Formula $C_{18}H_{22}N_2O_6$ Mol. weight 362,38 g/mol	 
BAA4670 Boc-L-Phe(4-Eth)-OH N-alpha-t-Butyloxycarbonyl-4-Ethynyl-L-phenylalanine CAS-No. 169158-05-8 Formula $C_{16}H_{19}NO_4$ Mol. weight 289,33 g/mol	 



		Product details
FAA7720 Fmoc-L-Phe(4-NH-Poc)-OH N-alpha-(9-Fluorenylmethyloxycarbonyl)-4-propargyloxycarbonylamino-L-phenylalanine CAS-No. 2576508-07-9 Formula $C_{28}H_{24}N_2O_6$ Mol. weight 484,5 g/mol		
FAA8530 Fmoc-L-Phe(4-Eth)-OH N-alpha-(9-Fluorenylmethyloxycarbonyl)-4-ethynyl-L-phenylalanine CAS-No. 1228049-41-9 Formula $C_{26}H_{21}NO_4$ Mol. weight 411,46 g/mol		

Propargylproline

		Product details
FAA7130 Fmoc-L-Pro(4-NHPoc)-OH (2S,4S) (2S,4S)-1-(9-Fluorenylmethyloxycarbonyl)-4-(propargyloxycarbonyl)amino-pyrrolidine-2-carboxylic acid CAS-No. 2451202-17-6 Formula $C_{26}H_{22}N_2O_6$ Mol. weight 434,44 g/mol		

Propargylserine



		Product details
HAA2355 H-L-Ser(Propargyl)-OH*HCl O-Propargyl-L-serine hydrochloride CAS-No. 1379150-93-2 Formula $C_6H_9NO_3 \cdot HCl$ Mol. weight 143,14*36,45 g/mol		






		Product details
BAA2260	Boc-L-Ser(Propargyl)-OH*DCHA	
N- α -t-Butyloxycarbonyl-O-propargyl-L-serine dicyclohexylamine		
CAS-No.	145205-94-3 (net)	
Formula	C ₁₁ H ₁₇ NO ₅ *C ₁₂ H ₂₃ N	
Mol. weight	243,26*181,32 g/mol	
		
FAA3820	Fmoc-L-Ser(Propargyl)-OH	
N- α -(9-Fluorenylmethyloxycarbonyl)-O-propargyl-L-serine		
CAS-No.	1354752-75-2	
Formula	C ₂₁ H ₁₉ NO ₅	
Mol. weight	365,38 g/mol	
		







References:

- *Isoxazolyl-serine-based agonists of peroxisome proliferator-activated receptor: design, synthesis, and effects on cardiomyocyte differentiation*; Z. L. Wei, P. A. Petukhov, F. Bizik, J. C. Teixeira, M. Mercola, E. A. Volpe, R. I. Glazer, T. M. Willson, A. P. Kozikowski; **J Am Chem Soc** 2004; **126**: 16714-5. <https://doi.org/10.1021/ja046386l>
- *Lacosamide isothiocyanate-based agents: novel agents to target and identify lacosamide receptors*; K. D. Park, P. Morieux, C. Salome, S. W. Cotten, O. Reamtong, C. Evers, S. J. Gaskell, J. P. Stables, R. Liu, H. Kohn; **J Med Chem** 2009; **52**: 6897-911. <https://doi.org/10.1021/jm9012054>

Propargyltyrosine

		Product details
HAA1970	H-L-Tyr(Propargyl)-OH	
O-Propargyl-L-tyrosine		
CAS-No.	610794-20-2	
Formula	C ₁₂ H ₁₃ NO ₃	
Mol. weight	219,24 g/mol	
		
HAA1971	H-L-Tyr(Propargyl)-OH*HCl	
O-Propargyl-L-tyrosine hydrochloride		
CAS-No.	1919043-11-0	
Formula	C ₁₂ H ₁₃ NO ₃ *HCl	
Mol. weight	219,24*36,45 g/mol	
		

		Product details
BAA2265 Boc-L-Tyr(Propargyl)-OH*DCHA N-alpha- <i>t</i> -Butyloxycarbonyl-O-propargyl-L-tyrosine dicyclohexylamine CAS-No. 340732-79-8 Formula C ₁₇ H ₂₁ NO ₅ Mol. weight 319,35 g/mol		
FAA3830 Fmoc-L-Tyr(Propargyl)-OH N-alpha-(9-Fluorenylmethyloxycarbonyl)-O-propargyl-L-tyrosine CAS-No. 1204595-05-0 Formula C ₂₇ H ₂₃ NO ₅ Mol. weight 441,48 g/mol		
Alkyne-Linker and Alkyne-ADC Linker		
PAA1170 Poc-Cystamine*HCl Prop-2-yn-1-yl (2-((2-aminoethyl)disulfaneyl)ethyl) carbamate CAS-No. 1266354-28-2 net Formula C ₈ H ₁₆ N ₂ O ₂ S ₂ *HCl Mol. weight 234,33*36,45 g/mol		
BNN1320 Boc-Cystamine-Poc Tert-butyl (2-((2-(((prop-2-yn-1-yloxy)carbonyl)amino)ethyl)disulfaneyl)ethyl)carbamate CAS-No. 2171512-56-2 Formula C ₁₃ H ₂₂ N ₂ O ₄ S ₂ Mol. weight 334,45 g/mol		
ADC1310 4-Pentynoyl-Val-Ala-PAB 4-pentynoyl-valyl-alanyl-(4-aminobenzyl alcohol) CAS-No. 1956294-75-9 Formula C ₂₀ H ₂₇ N ₃ O ₄ Mol. weight 373,45 g/mol		

		Product details
ADC1690	4-Pentynoyl-Val-Ala-PAB-Cl 4-Pentynoyl-valyl-alanyl-(4-aminobenzyl chloride) Formula $C_{20}H_{26}ClN_3O_3$ Mol. weight 391,90 g/mol	 
ADC1320	4-Pentynoyl-Val-Ala-PAB-PNP 4-pentynoyl-valyl-alanyl-(4-aminobenzyl)-(4-nitrophenyl)-carbonate CAS-No. 1956294-76-0 Formula $C_{27}H_{30}N_4O_8$ Mol. weight 538,55 g/mol	 
ADC1140	4-Pentynoyl-Val-Cit-PAB 4-pentynoyl-valyl-citrullyl-(4-aminobenzyl alcohol) CAS-No. 2708150-97-2 Formula $C_{23}H_{33}N_5O_5$ Mol. weight 459,54 g/mol	 
ADC1150	4-Pentynoyl-Val-Cit-PAB-PNP 4-pentynoyl-valyl-citrullyl-(4-aminobenzyl)-(4-nitrophenyl)-carbonate Formula $C_{30}H_{36}N_6O_9$ Mol. weight 624,64 g/mol	 
ADC1350	Alkyne-PEG(4)-Val-Ala-PAB propargyl-tetraethyleneglycol-propanoyl-valyl-alanyl-(4-aminobenzyl alcohol) CAS-No. 2348405-90-1 Formula $C_{29}H_{45}N_3O_9$ Mol. weight 579,68 g/mol	 
ADC1720	Alkyne-PEG(4)-Val-Ala-PAB-Cl Propargyl-tetraethyleneglycol-propanoyl-valyl-alanyl-(4-aminobenzyl chloride) Formula $C_{29}H_{44}ClN_3O_8$ Mol. weight 598,13 g/mol	 

[The Click Reaction](#)
[Amino Acid Derivatives and Related Building Blocks for Click Chemistry](#)
[Spermines and Amines for Click Chemistry](#)
[Click Reagents for Drug Delivery](#)
[Click Chemistry Tools for Proteomics](#)
[Carbohydrates for Click Chemistry](#)
[Proteolysis Targeting Chimeras \(PROTACs[®]\)](#)
[Index](#)
[↑ back to content](#)


ADC1360 Alkyne-PEG(4)-Val-Ala-PAB-PNP

propargyl-tetraethyleneglycol-propanoyl-valyl-alanyl-(4-aminobenzyl)-(4-nitrophenyl)-carbonate

CAS-No. 2348405-91-2

Formula $C_{36}H_{48}N_4O_{13}$


Mol. weight 744,79 g/mol

**ADC1180 Alkyne-PEG(5)-Val-Cit-PAB**

propargyl-tetraethyleneglycol-propanoyl-valyl-citru-lyl-(4-aminobenzyl alcohol)

Formula $C_{32}H_{51}N_5O_{10}$


Mol. weight 665,77 g/mol

**ADC1190 Alkyne-PEG(5)-Val-Cit-PAB-PNP**

propargyl-tetraethyleneglycol-propanoyl-valyl-citru-lyl-(4-aminobenzyl)-(4-nitrophenyl)-carbonate

Formula $C_{39}H_{54}N_6O_{14}$

Mol. weight 830,88 g/mol


**ADC1600 Propargyl-cyclobutane-1,1-dicarboxamide-Ala-PAB**

propargyl-cyclobutane-1,1-dicarboxamide-ala-nyl-(4-aminobenzyl alcohol)

CAS-No. 2576471-28-6

Formula $C_{19}H_{23}N_3O_4$

Mol. weight 357,40 g/mol


**ADC1610 Propargyl-cyclobutane-1,1-dicarboxamide-Ala-PAB-PNP**

propargyl-cyclobutane-1,1-dicarboxamide-ala-nyl-(4-aminobenzyl)-(4-nitrophenyl)-carbonate

CAS-No. 2576471-42-4

Formula $C_{26}H_{26}N_4O_8$

Mol. weight 522,51 g/mol


**ADC1500 Propargyl-cyclobutane-1,1-dicarboxamide-Cit-PAB**




propargyl-cyclobutane-1,1-dicarboxamide-citru-lyl-(4-aminobenzyl alcohol)

CAS-No. 2576471-27-5



Formula $C_{22}H_{29}N_5O_5$

Mol. weight 443,50 g/mol



		Product details
ADC1510	Propargyl-cyclobutane-1,1-dicarboxamide-Cit-PAB-PNP propargyl-cyclobutane-1,1-dicarboxamide-citru- lyl-(4-aminobenzyl)-(4-nitrophenyl)-carbonate	 
CAS-No.	2576471-40-2	
Formula	C ₂₉ H ₃₂ N ₆ O ₉	
Mol. weight	608,60 g/mol	
ADC1790	Alkyne-HMPO-OH (2-ethoxy-7-(prop-2-yn-1-yloxy)benzo[d][1,3]dioxol-5-yl) methanol	 
CAS-No.	3028213-63-7	
Formula	C ₁₃ H ₁₄ O ₅	
Mol. weight	250,25 g/mol	
ADC1800	Alkyne-HMPO-PNP (2-ethoxy-7-(prop-2-yn-1-yloxy)benzo[d][1,3]dioxol-5-yl) methyl (4-nitrophenyl) carbonate	 
CAS-No.	3028213-53-5	
Formula	C ₂₀ H ₁₇ NO ₉	
Mol. weight	415,35 g/mol	

Auxiliary Reagents

		Product details
RL-2010	TBTA Tris[(1-benzyl-1H-1,2,3-triazol-4-yl)methyl]amine	 
CAS-No.	510758-28-8	
Formula	C ₃₀ H ₃₀ N ₁₀	
Mol. weight	530,63 g/mol	
RL-2210	THPTA Tris(3-hydroxypropyltriazolylmethyl)amine	 
CAS-No.	760952-88-3	
Formula	C ₁₈ H ₃₀ N ₁₀ O ₃	
Mol. weight	434,51 g/mol	

3. Spermines and Amines for Click Chemistry

Polyamines are aliphatic cations with multiple functions in cell proliferation and differentiation and are essential for normal cell growth and development in eukaryotes. These molecules carry positive charges at their primary and secondary amino groups at physiological pH. Consequently, polyamines bind to various anionic macromolecules including DNA, RNA, acidic phospholipids, and certain proteins. These polycationic alkylamines are involved in various critical cellular functions, such as maintaining chromatin structure, regulating ion-channels, maintaining membrane stability, modulating protein synthesis, and scavenging free radicals. Polyamines also serve as substrates for transglutaminase reactions and for the synthesis of the translational regulator hypusine.

Crucial parts of the biological function of polyamines are the regulation of gene expression by altering DNA structure, the modulation of protein synthesis by binding to RNA, and the modulation of signal transduction pathways. The binding of polyamines to both RNA and DNA leads to conformational changes of those nucleic acids. Polyamines cause the conformational transition of DNA from the B form to the Z form and also cause bending of DNA. Both structural alterations are known to influence transcription. Close to 80% of all polyamines in the cell are associated with RNA, while spermine in particular has been shown to stabilize tRNA structures. Binding of polyamines to RNA causes structural changes that increase the efficiency of protein synthesis.

Polyamines are also known to modulate DNA-protein interactions, e.g., by enhancing the binding of specific gene-regulatory proteins to certain regulatory sequences termed response elements. The poly-amine spermine has been reported to facilitate the binding of estrogen receptor and nuclear factor κ B (NF- κ B) to their respective response elements at 100 to 500 μ M concentrations. Polyamines are also involved in modulating ligand-receptor interactions, for example N-methyl-D-aspartate (NMDA) receptors, which are important for the excitatory synaptic transmission in the brain and spinal cord.






Moreover, polyamines have been implicated as important molecules in virus-host interactions since many viruses utilize and manipulate polyamines for their own replication. Those pathogens depend on the presence of polyamines for numerous aspects of their replication cycles, such as DNA and RNA polymerization, genome packaging, and viral protein translation. Certain viruses even appear to stimulate polyamine synthesis upon infection, a fact that underlines the importance of this class of molecules for the viral life cycle.

The polyamine metabolic pathway and thus polyamine levels are strictly regulated in cells. However, dysregulation of polyamine metabolism is a frequently observed event in cancer. For example, elevated levels of polyamines have been associated with breast, colon, prostate, and skin cancers. Consequently, polyamine synthesis, metabolism, uptake, and function may be promising targets for cancer therapy.

References:

- *Structural analysis of DNA interactions with biogenic polyamines and cobalt(III)hexamine studied by Fourier transform infrared and capillary electrophoresis*; A. A. Ouameur, H. A. Tajmir-Riahi; **J Biol Chem** 2004; **279**: 42041-54. <https://doi.org/10.1074/jbc.M406053200>
- *Polyamines and Their Role in Virus Infection*; B. C. Mounce, M. E. Olsen, M. Vignuzzi, J. H. Connor; **Microbiol. Mol. Biol. Rev.** 2017; **81**: e00029-17. <https://doi.org/10.1128/MMBR.00029-17>

- Targeting polyamine metabolism for cancer therapy and prevention; T. R. Murray-Stewart, P. M. Woster, R. A. Casero, Jr.; *Biochem J* 2016; **473**: 2937-53. <https://doi.org/10.1042/BCJ20160383>
- Polyamine catabolism in carcinogenesis: potential targets for chemotherapy and chemoprevention; V. Battaglia, C. DeStefano Shields, T. Murray-Stewart, R. A. Casero, Jr.; *Amino Acids* 2014; **46**: 511-9. <https://doi.org/10.1007/s00726-013-1529-6>
- Targeting polyamine metabolism and function in cancer and other hyperproliferative diseases; R. A. Casero, Jr., L. J. Marton; *Nat Rev Drug Discov* 2007; **6**: 373-90. <https://doi.org/10.1038/nrd2243>

		Product details
FNN1020	Fmoc-EDA-N₃	
1-[(9-Fluorenylmethoxycarbonyl)amino]-2-azidoethane		
CAS-No.	432507-62-5	
Formula	C ₁₇ H ₁₆ N ₄ O ₂	
Mol. weight	308,33 g/mol	
		
FNN1030	Fmoc-DAP-N₃	
1-[(9-Fluorenylmethoxycarbonyl)amino]-3-azidopropane		
CAS-No.	1021422-85-4	
Formula	C ₁₈ H ₁₈ N ₄ O ₂	
Mol. weight	322,36 g/mol	
		
SNN1170	Spermine(HHNN₃)*3HCl	
N1-(3-Aminopropyl)-N4-(3-azidopropyl)butane-1,4-diamine trihydrochloride		
CAS-No.	1190203-77-0	
Formula	C ₁₀ H ₂₄ N ₆ *3HCl	
Mol. weight	228,34*109,38 g/mol	
		
SNN1230	Spermine(N₃BBB)	
tert-butyl (4-((3-azidopropyl)(tert-butoxycarbonyl)amino)butyl(3-((tert-butoxycarbonyl)amino)propyl) carbamate		
CAS-No.	1190203-80-5	
Formula	C ₂₅ H ₄₈ N ₆ O ₆	
Mol. weight	528,70 g/mol	
		
SNN1210	Spermine(N₃HHN₃)*2TsOH	
N1-(3-Aminopropyl)-N4-(3-azidopropyl)butane-1,4-diamine bistosylate		
CAS-No.	2250433-79-3	
Formula	C ₁₀ H ₂₂ N ₈ *C ₁₄ H ₁₆ O ₆ S ₂	
Mol. weight	254,34*344,40 g/mol	
		

[↑ back to content](#)

The Click Reaction

Amino Acid Derivatives and Related Building Blocks for Click Chemistry

Spermines and Amines for Click Chemistry

Click Reagents for Drug Delivery

Click Chemistry Tools for Proteomics

Carbohydrates for Click Chemistry


Proteolysis Targeting Chimeras (PROTACs[®])

Index

4. Click Reagents for Drug Delivery

4.1. Principles of Polymer Therapeutics

Peptides, proteins, and other biomolecules have a high potential as drugs due to their usually high specificity and efficacy. However, they often show poor pharmacokinetic properties, with their low stability under physiological conditions being a major factor. Since synthetic biomolecules are similar to endogenous molecules found in the human body, they are quickly degraded by enzymes and cleared from the system. Especially peptides and small proteins are susceptible to renal clearance. Additionally, the immunogenic responses and side effects elicited by many drugs, in particular protein drugs, are exacerbated by their hydrophobicity. Conjugation to biocompatible polymers, such as PEG (poly(ethylene glycol)), PGA (poly(glutamic acid)) or POX (poly(2-oxazoline)), increases aqueous solubility of a drug and often drastically enhances its pharmacokinetics at both the whole organism and subcellular level. By using bi- or multifunctional polymers, a linkage between two compounds can be formed or multivalent conjugates generated.



The main advantage of proteins, antibodies, siRNA, mRNA and other biomolecules when applied as drugs is their high specificity in combination with their low side effects, as they usually only interact with their dedicated target. A current focus is the study of targeted drug delivery systems for the controlled delivery and/or release of therapeutic agents. To this end, a biocompatible polymeric carrier is covalently linked to an active agent and a targeting moiety. This recognition part can be a peptide or protein that specifically binds to a certain cellular receptor, or an antibody against a specific antigen on the cell surface. Especially in the context of antibody-drug conjugate, the introduction of a PEG moiety can be very beneficial for pharmacokinetics by attenuating the frequently hydrophobic nature of payload molecules. After internalization of the whole conjugate, the active part (e.g., a nucleic acid or toxin) is released by e.g., variations in pH, temperature, or enzyme concentration. Consequently, the active agent is enabled to exert its function (e.g., the inhibition of a certain enzyme or the initiation of apoptosis) at a high local drug concentration, further increasing drug efficacy. If the active agent and/or the targeting entity are tailored to certain characteristics of an individual or a group of individuals (e.g., specific cancer cell antigens), this approach opens the door to individualized therapies ("personalized medicine").

To enlarge the field of Polymer Therapeutics to new classes of drug molecules, there is a constant search for new types of polymers. One interesting group are homopolymers of amino acids. Typical examples are polylysine and poly- α -glutamic acid (PGA) - polymers that do not exist in nature, but which show good physiological properties due to their similarity to natural proteins. In particular, poly(glutamic acid) has been identified as suitable carrier system.






Fig. 13: Multivalent presentation of a drug on poly- α -(glutamic acid).

PGA shows the ability to conjugate with partners on its N- and C-termini, analogous to the α and ω derivatization of a PEG – poly(ethylene glycol). Additionally, the glutamic acid side-chains may be used for further decoration of the polymer. Therefore, a multivalent presentation of a specific molecule along the polymer chain is possible, which is especially interesting for small molecules (Fig. 13). It is theoretically also possible to PEGylate a small molecule. However, inactivation of the small drug molecule is often the consequence. PGA is hence an ideal carrier for low molecular weight APIs.



4.2. PEGylation Improves Drug Delivery and Pharmacokinetics


Small drug molecules, but also large biomolecules like antibodies suffer from rapid clearance, causing a sharp decrease in plasma concentration of the drug as it is removed from the body. Consequently, drug administration has to be repeated within relatively short time intervals in order to keep its plasma concentration over a certain threshold. Otherwise, immunogenic reactions may be triggered.




PEGylated drugs show decreased rates of renal clearance and reduced immunogenicity. Consequently, plasma half-life of the drug is significantly increased, extending the time intervals between applications of the drug over the course of the treatment. This is due to the following mechanisms:

1. Preventing Degradation and Reducing Immunogenicity:

PEG chains cover the surface of a biopharmaceutical and thus effectively shield it from recognition by the immune system. This PEG layer has characteristics that are rather similar to a solvent, preventing uptake by cells of the retinal endothelial system (macrophage system). Therefore, recognition by the immune system (antibodies, proteases, and other degradation enzymes) is significantly attenuated. The drug stays intact and is not degraded or metabolized during its presence in and journey through the body.



Poly(ethylene glycol) is a polymeric linear structure with repeating polyethylene oxide units.



Mass spectrum of a polyethylene glycol showing the typical signals with a difference of $m/z = 44$

$$\mathcal{D} = \frac{M_w^\circ}{M_n^\circ} \geq 1 \text{ with } M_w^\circ = \frac{\sum N_x M_x^2}{\sum N_x M_x} \text{ and } M_n^\circ = \frac{\sum N_x M_x}{\sum N_x}$$

Whenever there is a distribution of molecular weights, the weight average M_w° is always greater than the number average M_n° and the polydispersity $\mathcal{D}$ is greater than 1.

Fig. 16: Composition of poly(ethylene glycols), typical mass spectrum of a polydisperse PEG, and formula for the polydispersity $\mathcal{D}$.

2. Preventing Excretion:

PEG is very hygroscopic by nature and surrounded by a large solvation sphere of water. Thus, the overall hydrodynamic radius of a biopharmaceutical may be increased by PEGylation by up to an order of magnitude, to a size larger than the diameter of the glomerular capillaries (6 to 12 nm). Consequently, a PEGylated drug can no longer be excreted through the kidneys, and pharmacologic half-life is significantly extended.

Chemical/Physical Properties and Quality Parameters of PEGs, Dispersity

Depending on whether a given PEG consists of a single molecular weight species (a defined number n of repeating units) or of a range of species with an average mass and a distribution of n around a mean value, PEG polymers are referred to as monodisperse or polydisperse, respectively. If the polymer is polydisperse, its mass spectrum will show a range of different molecular weights (Fig. 16).

A measure of the distribution of molecular weights in a polymer is given by the Dispersity $\mathcal{D}$, which is defined as the ratio between the weight average molecular weight M_w and the number average molecular weight M_n . The weight average M_w does not “count” species just by their number but takes into account the total weight of each species and is therefore a much more realistic indicator of the gross mechanical properties of a polymer.

In case of a homogeneous PEG, which consists only of one polymer species with a defined chain length, M_w is equal to M_n , thus the dispersity $\bar{D}$ equals 1 and the compound is referred to as monodisperse. Whenever there is a distribution of molecular weights, the weight average M_w is always greater than the number average M_n , and consequently the dispersity $\bar{D}$ is greater than 1. The dispersity of PEGs typically used in PEGylations ranges between 1.05 and 1.50.

However, whenever a PEGylated drug candidate needs to be approved by EMEA, FDA and other authorities, it is easier and faster if this compound is a defined species with a defined molecular weight. Therefore, the need for large but monodisperse PEGs is increasing.

Summary of Chemical and Physical Properties of PEGs:

- Good solubility in BOTH hydrophilic AND hydrophobic solvents as water, toluene, methylene chloride, and many other organic solvents.
- Insoluble in diethyl ether, hexane, ethylene glycol.
- Insoluble in water at elevated temperature.
- The solubility is influenced by forming derivatives.
- Highly mobile in water with high exclusion volume; large hydrodynamic radius.
- Complex formation with metal cations.
- Can be used to precipitate proteins and nucleic acids.
- Form two-phase system with aqueous solutions of other polymers.
- Non-toxic and FDA approved for use in drug products.

PEGylating Biopharmaceuticals and Small Molecules has the Following Effects:




- Improves solubility of conjugated molecules.
- Renders proteins non-immunogenic and tolerogenic.
- Reduces the rate of renal clearance through the kidney and alters pharmacokinetics.
- Alters electroosmotic flow.
- Increases cell permeability.

References:

- Structure of a PEGylated protein reveals a highly porous double-helical assembly; G. Cattani, L. Voegelé and P. B. Crowley; **Nat Chem** 2015; **7**: 823-8. <https://doi.org/10.1038/nchem.2342>
- Peptide and Protein PEGylation III: Advances in Chemistry and Clinical Applications; F. M. Veronese, J. M. Harris; **Adv. Drug Deliv. Rev.** 2008; **60**: 1-88.
- PEGylation of native disulfide bonds in proteins; S. Brocchini, S. Balan, A. Godwin, J. W. Choi, M. Zloh, S. Shaunak; **Nat Protoc** 2006; **1**: 2241-52. <https://doi.org/10.1038/nprot.2006.346>
- Pegylated Arg-Gly-Asp Peptide: 64Cu Labeling and PET Imaging of Brain Tumor $\alpha v \beta 3$ -Integrin Expression; X. Chen, Y. Hou, M. Tohme, R. Park, V. Khankaldyyan, I. Gonzales-Gomez, J. R. Bading, W. E. Laug, P. S. Conti; **J Nucl Med** 2004; **45**: 1776.
- Chemistry for peptide and protein PEGylation; M. J. Roberts, M. D. Bentley, J. M. Harris; **Adv Drug Deliv Rev** 2002; **54**: 459-76. [https://doi.org/10.1016/S0169-409X\(02\)00022-4](https://doi.org/10.1016/S0169-409X(02)00022-4)
- Peptide and protein PEGylation: a review of problems and solutions; F. M. Veronese; **Biomaterials** 2001; **22**: 405-17. [https://doi.org/10.1016/S0142-9612\(00\)00193-9](https://doi.org/10.1016/S0142-9612(00)00193-9)
- Functionalization of poly(ethylene glycol) and monomethoxy-poly(ethylene glycol); A. F. Bückmann, M. Morr, G. Johansson; **Makromol. Chem.** 1981; **182**: 1379-1384. <https://doi.org/10.1002/macp.1981.021820509>

- New, easily removable poly(ethylene glycol) supports for the liquid-phase method of peptide synthesis; V. N. R. Pillai, M. Mutter, E. Bayer, I. Gatfield; **J. Org. Chem.** 1980; **45**: 5364-5370. <https://doi.org/10.1021/jo01314a032>
- Synthesis and characterization of poly(ethylene glycol) derivatives; J. M. Harris, E. C. Struck, M. G. Case, M. S. Paley, M. Yalpani, J. M. Van Alstine, D. E. Brooks; **J Polymer Sci Polymer Chem Ed** 1984; **22**: 341-352. <https://doi.org/10.1002/pol.1984.170220207>
- Attachment of drugs to polyethylene glycols; S. Zalipsky, C. Gilon, A. Zilkha; **Eur. Polym. J.** 1983; **19**: 1177-1183. [https://doi.org/10.1016/0014-3057\(83\)90016-2](https://doi.org/10.1016/0014-3057(83)90016-2)
- Effect of covalent attachment of polyethylene glycol on immunogenicity and circulating life of bovine liver catalase; A. Abuchowski, J. R. McCoy, N. C. Palczuk, T. van Es, F. F. Davis; **J Biol Chem** 1977; **252**: 3582-6.
- Process design for large-scale purification of formate dehydrogenase from *Candida boidinii* by affinity partitioning; A. Cordes, M.-R. Kula; **J Chromat** 1986; **376**: 375-384. [https://doi.org/10.1016/s0378-4347\(00\)80853-1](https://doi.org/10.1016/s0378-4347(00)80853-1)
- Chemical modification of horseradish peroxidase with ethanal-methoxypolyethylene glycol: Solubility in organic solvents, activity, and properties; P. Wirth, J. Soupe, D. Tritsch, J.-F. Biellmann; **Bioorg. Chem.** 1991; **19**: 133-142. [https://doi.org/10.1016/0045-2068\(91\)90029-o](https://doi.org/10.1016/0045-2068(91)90029-o)
- The Synthesis of Substituted Methoxy-Poly(Ethyleneglycol) Derivatives Suitable for Selective Protein Modification; T. P. Kogan; **Synth. Commun.** 1992; **22**: 2417-2424. <https://doi.org/10.1080/00397919208019100>

4.3. Azido-PEG Derivatives for Click Chemistry


		Product details
RL-3935 Mal-PEG(3)-N₃ N-(2-(2-(2-(2-azidoethoxy)ethoxy)ethoxy)ethyl)-3-(2,5-dioxo-2,5-dihydro-1H-pyrrol-1-yl)propanamide CAS-No. 1858264-36-4 Formula C ₁₅ H ₂₃ N ₅ O ₆ Mol. weight 369,38 g/mol		
RL-4030 PFB-mercaptopropionyl-PEG3-N₃ Perfluorobiphenyl-mercaptopropionyl-PEG(3)-N ₃ Formula C ₂₃ H ₂₁ F ₉ N ₄ O ₄ S Mol. weight 620,49 g/mol		
PEG1088 N₃-PEG(8)-OH alpha-Azido-omega-hydroxy octa(ethylene glycol) CAS-No. 352439-36-2 Formula C ₁₆ H ₃₃ N ₃ O ₈ Mol. weight 395,45 g/mol		

PEG2065 Biotin-TEG-ATFBABiotin-triethylenglycol-(*p*-azido-tetrafluorobenzamide)

CAS-No. 1264662-85-2

Formula $C_{27}H_{37}F_4N_7O_6S$

Mol. weight 663,68 g/mol


**PEG4940 Biotin-PEG(3)-N₃**

11-[D(+)-Biotinylamino]-1-azido-3,6,9-trioxaundecane

CAS-No. 875770-34-6

Formula $C_{18}H_{32}N_6O_5S$

Mol. weight 444,55 g/mol


**PEG6795 Biotin-PEG3-Ph(4-N₃)**

Biotin-PEG3-phenyl azide

CAS-No. 2088238-77-9

Formula $C_{25}H_{37}N_7O_6S$


Mol. weight 563,67 g/mol

**PEG7990 Biotin-PEG(4)-N₃**(3a*S*,4*S*,6a*R*)-4-((28-azido-5,21-dioxo-9,12,15,18-tetraoxa-6,22-diazaoctacosyl)tetrahydro-1*H*-thieno[3,4-*d*]imidazol-2(3*H*)-one

CAS-No. 1006592-62-6

Formula $C_{27}H_{49}N_7O_5S$


Mol. weight 615,79 g/mol

**PEG7960 Biotin-PEG(4)-Dde-N₃**N-(19-azido-15-(4,4-dimethyl-2,6-dioxocyclohexylidene)-3,6,9,12-tetraoxa-16-azanonadecyl)-5-((3a*S*,4*S*,6a*R*)-2-oxohexahydro-1*H*-thieno[3,4-*d*]imidazol-4-yl)pentanamide

CAS-No. 1802907-93-2

Formula $C_{32}H_{53}N_7O_8S$


Mol. weight 695,87 g/mol






**PEG8000 Biotin-PEG(4)-Picolyl-N₃**(3a*S*,4*S*,6a*R*)-4-((1-(6-(azidomethyl)pyridin-3-ylamino)-1,17-dioxo-4,7,10,13-tetraoxa-16-azahenicosan-21-yl)tetrahydro-1*H*-thieno[3,4-*d*]imidazol-2(3*H*)-one

CAS-No. 2222687-71-8

Formula $C_{27}H_{42}N_8O_7S$

Mol. weight 622,74 g/mol



		Product details
PEG7970 Biotin-PEG(4)-Dde-Picolyl-N₃ N-(6-(azidomethyl)pyridin-3-yl)-15-(4,4-dimethyl-2,6-dioxocyclohexylidene)-1-(5-((3aS,4S,6aR)-2-oxohexahydro-1H-thieno[3,4-d]imidazol-4-yl)pentanamido)-3,6,9,12-tetraoxa-16-azanonadecan-19-amide 		
CAS-No. 2055048-42-3 Formula C ₃₈ H ₅₇ N ₉ O ₉ S Mol. weight 815,98 g/mol		
PEG8100 Biotin-PEG(4)-SS-Azide N-(2-((3-((3-azidopropyl)amino)-3-oxopropyl)disulfaneyl)ethyl)-1-(5-((3aS,4S,6aR)-2-oxohexahydro-1H-thieno[3,4-d]imidazol-4-yl)pentanamido)-3,6,9,12-tetraoxapentadecan-15-amide 		
CAS-No. 1260247-52-6 Formula C ₂₉ H ₅₂ N ₈ O ₈ S ₃ Mol. weight 736,96 g/mol		
PEG4330 Biotin-PEG(7)-N₃ alpha-Biotin-omega-azido hepta(ethylene glycol) 		
CAS-No. 1334172-75-6 Formula C ₂₆ H ₄₈ N ₆ O ₉ S Mol. weight 620,76 g/mol		
PEG4340 Biotin-PEG(11)-N₃ [2-(2-aminoethoxy)ethoxy]acetic acid <i>tert</i> -butyl ester*HCl 		
CAS-No. 956494-20-5 Formula C ₃₄ H ₆₄ N ₆ O ₁₃ S Mol. weight 796,97 g/mol		
PEG4350 Biotin-PEG(23)-N₃ alpha-Biotin-omega-azido 23(ethylene glycol) 		
CAS-No. 956494-20-5 Formula C ₅₈ H ₁₁₂ N ₆ O ₂₅ S Mol. weight 1325,6 g/mol		

The Click Reaction

Amino Acid Derivatives and Related Building Blocks for Click Chemistry

Spermines and Amines for Click Chemistry

Click Reagents for Drug Delivery







Click Chemistry Tools for Proteomics







Carbohydrates for Click Chemistry

Proteolysis Targeting Chimeras (PROTACs®)







Index







[↑ back to content](#)

		Product details	
PEG7215	N₃-PEG2-NMe3 mesylate		
2-(2-(2-azidoethoxy)ethoxy)-N,N,N-trimethylethan-1-aminium			
CAS-No.	1447915-25-4 net		
Formula	C ₁₀ H ₂₄ N ₄ O ₅ S		
Mol. weight	312,39 g/mol		
PEG4980	H₂N-PEG(2)-N₃*TosOH		
2-[2-(2-Azidoethoxy)ethoxy]ethanaminium tosylat			
CAS-No.	166388-57-4		
Formula	C ₇ H ₁₄ N ₄ O ₂ *C ₇ H ₈ O ₃ S		
Mol. weight	174,20*172,20 g/mol		
PEG3060	H₂N-PEG(3)-N₃		
1-Amino-11-azido-3,6,9-trioxaundecane			
CAS-No.	134179-38-7		
Formula	C ₈ H ₁₈ N ₄ O ₃		
Mol. weight	218,25 g/mol		
RL-4320	N₃-PEG(3)-NH-Suc		
1-azido-13-oxo-3,6,9-trioxa-12-azahexadecan-16-oic acid			
CAS-No.	1202400-17-6		
Formula	C ₁₂ H ₂₂ N ₄ O ₆		
Mol. weight	318,33 g/mol		
BNN1150	N₃-TOTA		
1-Azido-4,7,10-trioxa-13-tridecanamine			
CAS-No.	1162336-72-2		
Formula	C ₁₀ H ₂₂ N ₄ O ₃		
Mol. weight	246,31 g/mol		
PEG5320	N₃-PEG(4)-NH₂		
14-Azido-3,6,9,12-tetraoxatetradecan-1-amine			
CAS-No.	951671-92-4		
Formula	C ₁₀ H ₂₂ N ₄ O ₄		
Mol. weight	262,31 g/mol		

			Product details
PEG1087	H₂N-PEG(6)-N₃		
alpha-Amino-omega-azido hexa(ethylene glycol)			
CAS-No.	957486-82-7		
Formula	C ₁₄ H ₃₀ N ₄ O ₆		
Mol. weight	350,42 g/mol		
PEG2350	H₂N-PEG(7)-N₃		
alpha-Amino-omega-azido hepta(ethylene glycol)			
CAS-No.	1333154-77-0		
Formula	C ₁₆ H ₃₄ N ₄ O ₇		
Mol. weight	394,46 g/mol		
PEG3050	H₂N-PEG(8)-N₃		
alpha-Amino-omega-azido octa(ethylene glycol)			
CAS-No.	857891-82-8		
Formula	C ₁₈ H ₃₈ N ₄ O ₈		
Mol. weight	438,52 g/mol		
PEG1081	H₂N-PEG(11)-N₃		
alpha-Amino-omega-azido undecae(ethylene glycol)			
CAS-No.	1800414-71-4		
Formula	C ₂₄ H ₅₀ N ₄ O ₁₁		
Mol. weight	570,69 g/mol		
PEG3070	H₂N-PEG(23)-N₃		
alpha-Azido-omega-amino 23(ethylene glycol)			
CAS-No.	2172677-19-7		
Formula	C ₄₈ H ₉₈ N ₄ O ₂₃		
Mol. weight	1099,3 g/mol		
PEG3080	H₂N-PEG(35)-N₃		
alpha-Azido-omega-amino 35(ethylene glycol)			
CAS-No.	749244-38-0		
Formula	C ₇₂ H ₁₄₆ N ₄ O ₃₅		
Mol. weight	1627,94 g/mol		

[The Click Reaction](#)
[Amino Acid Derivatives and Related Building Blocks for Click Chemistry](#)
[Spermines and Amines for Click Chemistry](#)
[Click Reagents for Drug Delivery](#)
[Click Chemistry Tools for Proteomics](#)
[Carbohydrates for Click Chemistry](#)
[Proteolysis Targeting Chimeras \(PROTACs[®]\)](#)
[Index](#)
[↑ back to content](#)

			Product details
PEG3010	H ₂ N-PEG-N ₃ (3 kDa)		
alpha-Amino-omega-azido poly(ethylene glycol)			
Mol. weight	3000 Da		
PEG3030	H ₂ N-PEG-N ₃ (5 kDa)		
alpha-Amino-omega-azido poly(ethylene glycol)			
Mol. weight	5000 Da		
PEG3000	H ₂ N-PEG-N ₃ (10 kDa)		
alpha-Amino-omega-azido poly(ethylene glycol)			
Mol. weight	10000 Da		
PEG3020	H ₂ N-PEG-N ₃ (20 kDa)		
alpha-Amino-omega-azido poly(ethylene glycol)			
Mol. weight	20000 Da		
PEG4960	Boc-NH-PEG(2)-N ₃		
1-(t-Butyloxycarbonyl-amino)-3,6-dioxa-8-octaneazide			
CAS-No.	950683-55-3		
Formula	C ₁₁ H ₂₂ N ₄ O ₄		
Mol. weight	274,32 g/mol		
BAA4920	Boc-Aoa-NH-PEG3-N ₃		
tert-butyl ((14-azido-2-oxo-6,9,12-trioxa-3-azatetra-decyl)oxy)carbamate			
CAS-No.	1803191-52-7		
Formula	C ₁₅ H ₂₉ N ₅ O ₇		
Mol. weight	391,43 g/mol		

			Product details
PEG8160	N₃-PEG(3)-NH-Boc		
t-Butyl N-(2-(2-(2-(2-azidoethoxy)ethoxy)ethoxy)ethyl) carbamate			
CAS-No.	642091-68-7		
Formula	C ₁₃ H ₂₆ N ₄ O ₅		
Mol. weight	318,37 g/mol		
PEG4900	N₃-EEEt-OH		
2-[2-(2-Azidoethoxy)ethoxy]ethanol			
CAS-No.	86520-52-7		
Formula	C ₆ H ₁₃ N ₃ O ₃		
Mol. weight	175,19 g/mol		
PEG3760	N₃-PEG(3)-OH		
alpha-Azido-omega-hydroxy tetra(ethylene glycol)			
CAS-No.	86770-67-4		
Formula	C ₈ H ₁₇ N ₃ O ₄		
Mol. weight	219,24 g/mol		
PEG5300	N₃-PEG(4)-OH		
2-(2-(2-(2-(2-Azidoethoxy)ethoxy)ethoxy)ethoxy) ethanol			
CAS-No.	86770-68-5		
Formula	C ₁₀ H ₂₁ N ₃ O ₅		
Mol. weight	263,29 g/mol		
PEG6720	N₃-PEG(5)-OH		
17-Azido-3,6,9,12,15-pentaoxaheptadecan-1-ol			
CAS-No.	86770-69-6		
Formula	C ₁₂ H ₂₅ N ₃ O ₆		
Mol. weight	307,34 g/mol		
PEG1390	N₃-PEG(12)-OH		
35-Azido-3,6,9,12,15,18,21,24,27,30,33-undecaaxapentatriacontan-1-ol			
CAS-No.	73342-16-2		
Formula	C ₂₄ H ₄₉ N ₃ O ₁₂		
Mol. weight	571,66 g/mol		

[The Click Reaction](#)
[Amino Acid Derivatives and Related Building Blocks for Click Chemistry](#)
[Spermines and Amines for Click Chemistry](#)
[Click Reagents for Drug Delivery](#)
[Click Chemistry Tools for Proteomics](#)
[Carbohydrates for Click Chemistry](#)
[Proteolysis Targeting Chimeras \(PROTACs[®]\)](#)
[Index](#)
[↑ back to content](#)


PEG3770 N₃-PEG(24)-OH

alpha-Azido-omega-hydroxy 24(ethylene glycol)

CAS-No. 73342-16-2

Formula C₄₈H₉₇N₃O₂₄

Mol. weight 1100,29 g/mol


**PEG3780 N₃-PEG(36)-OH**

alpha-Azido-omega-hydroxy 36(ethylene glycol)

CAS-No. 73342-16-2

Formula C₇₂H₁₄₅N₃O₃₆


Mol. weight 1628,92 g/mol

**PEG5350 HO-PEG-N₃ (3 kDa)**

alpha-Hydroxy-omega-azido poly(ethylene glycol)

CAS-No. 73342-16-2


Mol. weight 3000 Da

**PEG5360 HO-PEG-N₃ (5 kDa)**

alpha-Hydroxy-omega-azido poly(ethylene glycol)

CAS-No. 73342-16-2


Mol. weight 5000 Da

**PEG5330 HO-PEG-N₃ (10 kDa)**

alpha-Hydroxy-omega-azido poly(ethylene glycol)

CAS-No. 73342-16-2


Mol. weight 10000 Da







**PEG5340 HO-PEG-N₃ (20 kDa)**

alpha-Hydroxy-omega-azido poly(ethylene glycol)







CAS-No. 73342-16-2







Mol. weight 20000 Da






			Product details
PEG5170	N₃-TOTA-Suc		
1-Azido-4,7,10-trioxa-13-tridecaneamine succinamic acid			
CAS-No.	1993176-74-1		
Formula	C ₁₄ H ₂₆ N ₄ O ₆		
Mol. weight	346,38 g/mol		
PEG5400	N₃-AEEEA*CHA		
11-Azido-3,6,9-trioxaundecanoic acid cyclohexylamine			
CAS-No.	172531-37-2		
Formula	C ₈ H ₁₅ N ₃ O ₅ *C ₆ H ₁₃ N		
Mol. weight	233,22*99,17 g/mol		
PEG7950	N₃-AEEA-OK		
Potassium 8-azido-3,6-dioxaoctanoate			
CAS-No.	882518-90-3		
Formula	C ₆ H ₁₀ KN ₃ O ₄		
Mol. weight	39,10*188,16 g/mol		
PEG5390	N₃-O₂Oc-OtBu		
8-Azido-3,6-dioxaoctanoic acid t-butyl ester			
CAS-No.	251564-45-1		
Formula	C ₁₀ H ₁₉ N ₃ O ₄		
Mol. weight	245,28 g/mol		
PEG2780	N₃-O₂Oc-OH*CHA		
[2-(2-azidoethoxy)ethoxy]acetic acid cyclohexylamine salt			
CAS-No.	2098500-94-6		
Formula	C ₆ H ₁₁ N ₃ O ₄ *C ₆ H ₁₃ N		
Mol. weight	189,17*99,17 g/mol		
PEG2790	N₃-O₂Oc-O₂Oc-OH		
8-(8-Azido-3,6-dioxaoctanoylamido)-3,6-dioxaoctanoic acid			
CAS-No.	1254054-60-8		
Formula	C ₁₂ H ₂₂ N ₄ O ₇		
Mol. weight	334,33 g/mol		

[The Click Reaction](#)
[Amino Acid Derivatives and Related Building Blocks for Click Chemistry](#)
[Spermines and Amines for Click Chemistry](#)
[Click Reagents for Drug Delivery](#)
[Click Chemistry Tools for Proteomics](#)
[Carbohydrates for Click Chemistry](#)
[Proteolysis Targeting Chimeras \(PROTACs*\)](#)
[Index](#)
[↑ back to content](#)



			Product details
RL-4370	N₃-PEG(3)-NH-DIG-OH		
Diglycolic acid PEG3 azide			
CAS-No.	239081-53-9		
Formula	C ₁₂ H ₂₂ N ₄ O ₇		
Mol. weight	334,33 g/mol		
RL-4380	Fmoc-NH-PEG(3)-N₃		
(9H-fluoren-9-yl)methyl (2-(2-(2-azidoethoxy)ethoxy)ethoxy)ethyl)carbamate			
CAS-No.	1172605-58-1		
Formula	C ₂₃ H ₂₈ N ₄ O ₅		
Mol. weight	440,50 g/mol		
PEG2345	N₃-PEG(4)-COOH		
15-Azido-4,7,10,13-tetraoxa-pentadecanoic acid			
CAS-No.	1257063-35-6		
Formula	C ₁₁ H ₂₁ N ₃ O ₆		
Mol. weight	291,3 g/mol		
PEG6915	N₃-PEG(4)-TFP		
2,3,5,6-tetrafluorophenyl 1-azido-3,6,9,12-tetraoxapentadecan-15-oate			
CAS-No.	1807505-33-4		
Formula	C ₁₇ H ₂₁ F ₄ N ₃ O ₆		
Mol. weight	439,36 g/mol		
PEG1400	N₃-PEG(4)-NHS		
15-Azido-4,7,10,13-tetraoxa-pentadecanoic acid succinimidyl ester			
CAS-No.	944251-24-5		
Formula	C ₁₅ H ₂₄ N ₄ O ₈		
Mol. weight	388,37 g/mol		
PEG4170	N₃-PEG(8)-COOH		
alpha-Azido-omega-(propionic acid) octa(ethylene glycol)			
CAS-No.	1214319-92-2		
Formula	C ₁₉ H ₃₇ N ₃ O ₁₀		
Mol. weight	467,51 g/mol		






		Product details
PEG1405 N₃-PEG(8)-NHS 1-Azido-3,6,9,12,15,18,21,24-octaoxaheptacosan-27-oic acid succinimidyl ester CAS-No. 1204834-00-3 Formula C ₂₃ H ₄₀ N ₄ O ₁₂ Mol. weight 564,58 g/mol		
PEG2015 N₃-PEG(9)-COOH O-(2-Azidoethyl)-O-[2-(diglycolyl-amino)ethyl]heptaethylene glycol CAS-No. 846549-37-9 Formula C ₂₂ H ₄₂ N ₄ O ₁₂ Mol. weight 554,59 g/mol		
PEG6895 Aminooxy-PEG(11)-N₃*HCl O-(35-azido-3,6,9,12,15,18,21,24,27,30,33-undecaoxapentatriacontyl)hydroxylamine hydrochloride CAS-No. 2560602-21-1 net Formula C ₂₄ H ₅₀ N ₄ O ₁₂ *HCl Mol. weight 586,68*36,45 g/mol		
PEG4180 N₃-PEG(12)-COOH alpha-Azido-omega-(propionic acid) dodeca(ethylene glycol) CAS-No. 1167575-20-3 Formula C ₂₇ H ₅₃ N ₃ O ₁₄ Mol. weight 643,72 g/mol		
PEG6925 N₃-PEG(12)-TFP 2,3,5,6-tetrafluorophenyl 1-azido-3,6,9,12,15,18,21,24,27,30,33,36-dodecaoxanonatriacantan-39-oate CAS-No. 2735663-71-3 Formula C ₃₃ H ₅₃ F ₄ N ₃ O ₁₄ Mol. weight 791,79 g/mol		
PEG1395 N₃-PEG(12)-NHS 1-Azido-3,6,9,12,15,18,21,24,27,30,33,36-dodecaoxanonatriacantan-39-oic acid succinimidyl ester CAS-No. 1108750-59-9 Formula C ₃₁ H ₅₆ N ₄ O ₁₆ Mol. weight 740,79 g/mol		

[The Click Reaction](#)
[Amino Acid Derivatives and Related Building Blocks for Click Chemistry](#)
[Spermines and Amines for Click Chemistry](#)
[Click Reagents for Drug Delivery](#)
[Click Chemistry Tools for Proteomics](#)
[Carbohydrates for Click Chemistry](#)
[Proteolysis Targeting Chimeras \(PROTACs[®]\)](#)
[Index](#)
[↑ back to content](#)

			Product details
PEG4190	N₃-PEG(24)-COOH		
alpha-Azido-omega-(propionic acid) 24(ethylene glycol)			
CAS-No.	1167575-20-3		
Formula	C ₅₁ H ₁₀₁ N ₃ O ₂₆		
Mol. weight	1172,35 g/mol		
PEG7650	N₃-PEG(24)-TFP		
alpha-Azido-omega-(2,3,5,6-tetrafluorophenyl propionate) 24(ethylene glycol)			
Formula	C ₅₇ H ₁₀₁ F ₄ N ₃ O ₂₆		
Mol. weight	1320,41 g/mol		
PEG7660	N₃-PEG(36)-TFP		
alpha-Azido-omega-(2,3,5,6-tetrafluorophenyl propionate) 36(ethylene glycol)			
CAS-No.	3017961-70-2		
Formula	C ₈₁ H ₁₄₉ F ₄ N ₃ O ₃₈		
Mol. weight	1849,04 g/mol		







Whenever free thiol groups (e.g., from cysteine) are used for conjugation, maleimides are typically the reaction partner of choice. However, maleimides also react with other functional groups, for example -COOH, -OH or -NH₂ which may lead to the formation of unwanted impurities. The iodo group reacts more specifically with thiols, resulting in much cleaner conjugates.







			Product details
PEG7190	Bromoacetamido-PEG(3)-N₃		
Bromoacetamido-tri(ethylene glycol)-azide			
CAS-No.	940005-81-2		
Formula	C ₁₀ H ₁₉ BrN ₄ O ₄		
Mol. weight	339,19 g/mol		
PEG7200	Bromoacetamido-PEG(11)-N₃		
Bromoacetamido-undeca(ethylene glycol)-azide			
CAS-No.	2172677-17-5		
Formula	C ₂₆ H ₅₁ BrN ₄ O ₁₂		
Mol. weight	691,61 g/mol		

			Product details
PEG7210	Bromoacetamido-PEG(23)-N₃		
Bromoacetamido-23(ethylene glycol)-azide			
CAS-No.	2735663-83-7		
Formula	C ₅₀ H ₉₉ BrN ₄ O ₂₄		
Mol. weight	1220,24 g/mol		
PEG3130	I-PEG-N₃ (10 kDa)		
alpha-Iodo-omega-azido poly(ethylene glycol)			
Mol. weight	10000 Da		
PEG3140	I-PEG-N₃ (20 kDa)		
alpha-Iodo-omega-azido poly(ethylene glycol)			
Mol. weight	20000 Da		
PEG3150	I-PEG-N₃ (3 kDa)		
alpha-Iodo-omega-azido poly(ethylene glycol)			
Mol. weight	3000 Da		
PEG3160	I-PEG-N₃ (5 kDa)		
alpha-Iodo-omega-azido poly(ethylene glycol)			
Mol. weight	5000 Da		







Reference:


- *Quantitative reactivity profiling predicts functional cysteines in proteomes*; E. Weerapana, C. Wang, G. M. Simon, F. Richter, S. Khare, M. B. Dillon, D. A. Bachovchin, K. Mowen, D. Baker, B. F. Cravatt; **Nature** 2010; **468**: 790-5. <https://doi.org/10.1038/nature09472>

			Product details
PEG1690	MeO-PEG(4)-N₃		
13-Azido-2,5,8,11-tetraoxa-tridecane			
CAS-No.	606130-90-9		
Formula	C ₉ H ₁₉ N ₃ O ₄		
Mol. weight	233,26 g/mol		
PEG1705	MeO-PEG(8)-N₃		
25-Azido-2,5,8,11,14,17,20,23-octaopentacosane			
CAS-No.	869718-80-9		
Formula	C ₁₇ H ₃₅ N ₃ O ₈		
Mol. weight	409,48 g/mol		
PEG1660	MeO-PEG(12)-N₃		
37-Azido-2,5,8,11,14,17,20,23,26,29,32,35-dodecaoxaheptatriacontane			
CAS-No.	2170098-29-8		
Formula	C ₂₅ H ₅₁ N ₃ O ₁₂		
Mol. weight	585,69 g/mol		
PEG1710	MeO-PEG(24)-N₃		
alpha-Methoxy-omega-azido-24(ethylene glycol)			
CAS-No.	89485-61-0		
Formula	C ₄₉ H ₉₉ N ₃ O ₂₄		
Mol. weight	1114,34 g/mol		
PEG3430	MeO-PEG(36)-N₃		
alpha-Methoxy-omega-azido-36(ethylene glycol)			
CAS-No.	89485-61-0		
Formula	C ₇₃ H ₁₄₇ N ₃ O ₃₆		
Mol. weight	1642,95 g/mol		
PEG1219	MeO-PEG-N₃ (750 Da)		
alpha-Methoxy-omega-azido poly(ethylene glycol)			
Mol. weight	750 Da		

			Product details
PEG1225	MeO-PEG-N₃ (2 kDa)		
alpha-Methoxy-omega-azido poly(ethylene glycol)			
Mol. weight	2000 Da		
PEG2040	MeO-PEG-N₃ (5 kDa)		
alpha-Methoxy-omega-azido poly(ethylene glycol)			
Mol. weight	5000 Da		
PEG2045	MeO-PEG-N₃ (10 kDa)		
alpha-Methoxy-omega-azido poly(ethylene glycol)			
Mol. weight	10000 Da		
PEG2050	MeO-PEG-N₃ (20 kDa)		
alpha-Methoxy-omega-azido poly(ethylene glycol)			
Mol. weight	20000 Da		
PEG4830	Azido-PEG-Si(OMe)₃ (3 kDa)		
alpha-Azido-omega-trimethoxysilyl poly(ethylene glycol)			
Mol. weight	3000 Da		
PEG4835	Azido-PEG-Si(OMe)₃ (5 kDa)		
alpha-Azido-omega-trimethoxysilyl poly(ethylene glycol)			
Mol. weight	5000 Da		





[The Click Reaction](#)
[Amino Acid Derivatives and Related Building Blocks for Click Chemistry](#)
[Spermines and Amines for Click Chemistry](#)
[Click Reagents for Drug Delivery](#)
[Click Chemistry Tools for Proteomics](#)
[Carbohydrates for Click Chemistry](#)
[Proteolysis Targeting Chimeras \(PROTACs[®]\)](#)
[Index](#)
[↑ back to content](#)

			Product details
PEG4840	Azido-PEG-Si(OMe)₃ (10 kDa)		
alpha-Azido-omega-trimethoxysilyl poly(ethylene glycol)			
Mol. weight	10000 Da		
PEG4845	Azido-PEG-Si(OMe)₃ (20 kDa)		
alpha-Azido-omega-trimethoxysilyl poly(ethylene glycol)			
Mol. weight	20000 Da		
RL-4420	N₃-PEG(3)-SQA		
3-((2-(2-(2-azidoethoxy)ethoxy)ethoxy)ethyl)amino)-4-ethoxycyclobut-3-ene-1,2-dione			
CAS-No.	2920089-51-4		
Formula	C ₁₄ H ₂₂ N ₄ O ₆		
Mol. weight	342,35 g/mol		
PEG6655	N₃-PEG-SQA (3 kDa)		
alpha-Azido-omega-squaric acid ethyl ester poly(ethylene glycol)			
Mol. weight	3000 Da		
PEG6660	N₃-PEG-SQA (5 kDa)		
alpha-Azido-omega-squaric acid ethyl ester poly(ethylene glycol)			
Mol. weight	5000 Da		
PEG6645	N₃-PEG-SQA (10 kDa)		
alpha-Azido-omega-squaric acid ethyl ester poly(ethylene glycol)			
Mol. weight	10000 Da		







		Product details
PEG6650	N₃-PEG-SQA (20 kDa)	
alpha-Azido-omega-squaric acid ethyl ester poly(ethylene glycol)		
Mol. weight	20000 Da	
		







[The Click Reaction](#)
[Amino Acid Derivatives and Related Building Blocks for Click Chemistry](#)

4.4. Alkyne-PEG Derivatives for Click Chemistry

		Product details
PEG4950	Biotin-PEG(4)-alkyne	
15-[D(+)-Biotinylamino]-4,7,10,13-tetraoxapentadec-1-yne		
CAS-No.	1262681-31-1	
Formula	C ₂₁ H ₃₅ N ₃ O ₆ S	
Mol. weight	457,58 g/mol	
		
PEG7980	Biotin-PEG(4)-Dde-Alkyne	
N-(15-(4,4-dimethyl-2,6-dioxocyclohexylidene)-3,6,9,12-tetraoxa-16-azanonadec-18-ynyl)-5-((3aS,4S,6aR)-2-oxohexahydro-1H-thieno[3,4-d]imidazol-4-yl)pentanamide		
CAS-No.	1802908-00-4	
Formula	C ₃₂ H ₅₀ N ₄ O ₈ S	
Mol. weight	650,83 g/mol	
		
PEG8010	Biotin-PEG(4)-Alkyne	
(3aS,4S,6aR)-4-(5-(2,1-dioxo-8,11,14,17-tetraoxa-4,20-diazapentacos-1-yn-25-yl)tetrahydro-1H-thieno[3,4-d]imidazol-2(3H)-one		
CAS-No.	1006592-45-5	
Formula	C ₂₄ H ₄₀ N ₄ O ₇ S	
Mol. weight	528,66 g/mol	
		
PEG8110	Biotin-PEG(4)-SS-Alkyne	
N-(2-((3-oxo-3-(prop-2-ynylamino)propyl)disulfany)ethyl)-1-(5-((3aS,4S,6aR)-2-oxohexahydro-1H-thieno[3,4-d]imidazol-4-yl)pentanamido)-3,6,9,12-tetraoxapentadecan-15-amide		
CAS-No.	1260247-54-8	
Formula	C ₂₉ H ₄₉ N ₅ O ₈ S ₃	
Mol. weight	691,92 g/mol	
		

[Spermines and Amines for Click Chemistry](#)
[Click Reagents for Drug Delivery](#)
[Click Chemistry Tools for Proteomics](#)
[Carbohydrates for Click Chemistry](#)
[Proteolysis Targeting Chimeras \(PROTACs[®]\)](#)
[Index](#)
[↑ back to content](#)

			Product details
PEG6815	Alkyne-O₂Oc-OH		
Pentynoyl-8-amino-3,6-dioxaoctanoic acid			
CAS-No.	2134622-80-1		
Formula	C ₁₁ H ₁₇ NO ₅		
Mol. weight	243,26 g/mol		
PEG6805	Biotin-PEG3-NH-Pentynoyl		
alpha-Biotin-omega-pentynoyl 3(ethylene glycol)			
CAS-No.	869082-82-6		
Formula	C ₂₃ H ₃₈ N ₄ O ₆ S		
Mol. weight	498,64 g/mol		
PEG5430	Alkyne-PEG(4)-NH₂		
Alkyne-PEG(4)-amine			
CAS-No.	1013921-36-2		
Formula	C ₁₁ H ₂₁ NO ₄		
Mol. weight	231,29 g/mol		
PEG2960	H₂N-PEG-alkyne (3 kDa)		
alpha-Amino-omega-propargylacetamido poly(ethylene glycol)			
Mol. weight	3000 Da		
PEG2980	H₂N-PEG-alkyne (5 kDa)		
alpha-Amino-omega-propargylacetamido poly(ethylene glycol)			
Mol. weight	5000 Da		
PEG2950	H₂N-PEG-alkyne (10 kDa)		
alpha-Amino-omega-propargylacetamido poly(ethylene glycol)			
Mol. weight	10000 Da		


		Product details	
PEG2970	H₂N-PEG-alkyne (20 kDa)		
alpha-Amino-omega-propargylacetamido poly(ethylene glycol)			
Mol. weight	20000 Da		
PAA1050	Poc-O₂Oc-OH*DCHA		
8-(Popargyloxycarbonyl-amino)-3,6-dioxaoctanoic acid dicyclohexylamine			
CAS-No.	2988660-57-5		
Formula	C ₁₀ H ₁₅ NO ₆ *C ₁₂ H ₂₃ N		
Mol. weight	245,23*181,32 g/mol		
PEG8170	Propargyl-PEG(5)-COOH		
4,7,10,13,16-pentaoxanonadec-18-ynoic acid			
CAS-No.	1245823-51-1		
Formula	C ₁₄ H ₂₄ O ₇		
Mol. weight	304,34 g/mol		
PEG5410	Alkyne-PEG(4)-NHS		
Alkyne-PEG(4)-succinimidyl ester			
CAS-No.	1393330-40-9		
Formula	C ₁₈ H ₂₇ NO ₉		
Mol. weight	401,41 g/mol		
PEG2860	NHS-PEG-alkyne (3 kDa)		
alpha-Succinimidyl ester-omega-propargylacetamido poly(ethylene glycol)			
Mol. weight	3000 Da		
PEG2880	NHS-PEG-alkyne (5 kDa)		
alpha-Succinimidyl ester-omega-propargylacetamido poly(ethylene glycol)			
Mol. weight	5000 Da		

[The Click Reaction](#)
[Amino Acid Derivatives and Related Building Blocks for Click Chemistry](#)
[Spermines and Amines for Click Chemistry](#)
[Click Reagents for Drug Delivery](#)
[Click Chemistry Tools for Proteomics](#)
[Carbohydrates for Click Chemistry](#)
[Proteolysis Targeting Chimeras \(PROTACs*\)](#)
[Index](#)
[↑ back to content](#)

PEG2850 NHS-PEG-alkyne (10 kDa)


alpha-Succinimidyl ester-omega-propargylacetamido poly(ethylene glycol)

Mol. weight 10000 Da

**PEG2870 NHS-PEG-alkyne (20 kDa)**


alpha-Succinimidyl ester-omega-propargylacetamido poly(ethylene glycol)

Mol. weight 20000 Da

**PEG2910 NHS-PEG(NH-Boc)-alkyne (3 kDa)**


alpha-Succinimidyl ester-omega-(N-t-Butyloxycarbonyl-L-propargyl-glycyl) poly(ethylene glycol)

Mol. weight 3000 Da

**PEG2930 NHS-PEG(NH-Boc)-alkyne (5 kDa)**


alpha-Succinimidyl ester-omega-(N-t-Butyloxycarbonyl-L-propargyl-glycyl) poly(ethylene glycol)

Mol. weight 5000 Da

**PEG2900 NHS-PEG(NH-Boc)-alkyne (10 kDa)**


alpha-Succinimidyl ester-omega-(N-t-Butyloxycarbonyl-L-propargyl-glycyl) poly(ethylene glycol)





Mol. weight 10000 Da

**PEG2920 NHS-PEG(NH-Boc)-alkyne (20 kDa)**

alpha-Succinimidyl ester-omega-(N-t-Butyloxycarbonyl-L-propargyl-glycyl) poly(ethylene glycol)







Mol. weight 20000 Da









		Product details
PEG2915	NHS-PEG(NH-Fmoc)-alkyne (3 kDa)	
alpha-Succinimidyl ester-omega-(N-(9-Fluorenylmethyloxycarbonyl)-L-propargyl-glycyl) poly(ethylene glycol)		
Mol. weight	3000 Da	
		
PEG2935	NHS-PEG(NH-Fmoc)-alkyne (5 kDa)	
alpha-Succinimidyl ester-omega-(N-(9-Fluorenylmethyloxycarbonyl)-L-propargyl-glycyl) poly(ethylene glycol)		
Mol. weight	5000 Da	
		
PEG2905	NHS-PEG(NH-Fmoc)-alkyne (10 kDa)	
alpha-Succinimidyl ester-omega-(N-(9-Fluorenylmethyloxycarbonyl)-L-propargyl-glycyl) poly(ethylene glycol)		
Mol. weight	10000 Da	
		
PEG2925	NHS-PEG(NH-Fmoc)-alkyne (20 kDa)	
alpha-Succinimidyl ester-omega-(N-(9-Fluorenylmethyloxycarbonyl)-L-propargyl-glycyl) poly(ethylene glycol)		
Mol. weight	20000 Da	
		

Maleimides are frequently used for conjugation to free thiol groups (e.g., from cysteine). However, maleimides also react with other functional groups such as OH or NH₂, potentially leading to the formation of impurities. The iodo group reacts more specifically with thiols, resulting in much cleaner conjugates.

		Product details
PEG3110	I-PEG-alkyne (3 kDa)	
alpha-Iodo-omega-propargylacetamido poly(ethylene glycol)		
Mol. weight	3000 Da	
		

			Product details
PEG3120	I-PEG-alkyne (5 kDa)		
alpha-Iodo-omega-propargylacetamido poly(ethylene glycol)			
Mol. weight	5000 Da		
PEG3090	I-PEG-alkyne (10 kDa)		
alpha-Iodo-omega-propargylacetamido poly(ethylene glycol)			
Mol. weight	10000 Da		
PEG3100	I-PEG-alkyne (20 kDa)		
alpha-Iodo-omega-propargylacetamido poly(ethylene glycol)			
Mol. weight	20000 Da		
PEG5440	Alkyne-PEG(4)-mal		
Alkyne-PEG(4)-maleimide			
CAS-No.	1609651-90-2		
Formula	C ₁₈ H ₂₆ N ₂ O ₇		
Mol. weight	382,41 g/mol		
PEG2840	MeO-PEG-alkyne (750 Da)		
alpha-Methoxy-omega-propargylacetamido poly(ethylene glycol)			
CAS-No.	1993176-75-2		
Mol. weight	750 Da		
PEG2810	MeO-PEG-alkyne (2 kDa)		
alpha-Methoxy-omega-propargylacetamido poly(ethylene glycol)			
CAS-No.	850180-69-7		
Mol. weight	2000 Da		

		Product details
PEG2830	MeO-PEG-alkyne (5 kDa)	
alpha-Methoxy-omega-propargylacetamido poly(ethylene glycol)		
CAS-No.	850180-69-7	
Mol. weight	5000 Da	
		
PEG2800	MeO-PEG-alkyne (10 kDa)	
alpha-Methoxy-omega-propargylacetamido poly(ethylene glycol)		
CAS-No.	850180-69-7	
Mol. weight	10000 Da	
		
PEG2820	MeO-PEG-alkyne (20 kDa)	
alpha-Methoxy-omega-propargylacetamido poly(ethylene glycol)		
CAS-No.	850180-69-7	
Mol. weight	20000 Da	
		
PEG4810	Alkyne-PEG-Si(OMe)₃ (3 kDa)	
alpha-Propargylacetamido-omega-trimethoxysilyl poly(ethylene glycol)		
Mol. weight	3000 Da	
		
PEG4815	Alkyne-PEG-Si(OMe)₃ (5 kDa)	
alpha-Propargylacetamido-omega-trimethoxysilyl poly(ethylene glycol)		
Mol. weight	5000 Da	
		
PEG4820	Alkyne-PEG-Si(OMe)₃ (10 kDa)	
alpha-Propargylacetamido-omega-trimethoxysilyl poly(ethylene glycol)		
Mol. weight	10000 Da	
		

The Click Reaction

 Amino Acid Derivatives and Related
 Building Blocks for Click Chemistry

Spermines and Amines for Click Chemistry

Click Reagents for Drug Delivery






Click Chemistry Tools for Proteomics

Carbohydrates for Click Chemistry






 Proteolysis Targeting Chimeras (PROTACs[®])

Index







[↑ back to content](#)





			Product details
PEG4825	Alkyne-PEG-Si(OMe)₃ (20 kDa)		
alpha-Propargylacetamido-omega-trimethoxysilyl poly(ethylene glycol)			
Mol. weight	20000 Da		
PEG6570	Alkynyl-PEG-SQA (3 kDa)		
alpha-Pentynyl-omega-squaric acid ethyl ester poly(ethylene glycol)			
Mol. weight	3000 Da		
PEG6575	Alkynyl-PEG-SQA (5 kDa)		
alpha-Pentynyl-omega-squaric acid ethyl ester poly(ethylene glycol)			
Mol. weight	5000 Da		
PEG6560	Alkynyl-PEG-SQA (10 kDa)		
alpha-Pentynyl-omega-squaric acid ethyl ester poly(ethylene glycol)			
Mol. weight	10000 Da		
PEG6565	Alkynyl-PEG-SQA (20 kDa)		
alpha-Pentynyl-omega-squaric acid ethyl ester poly(ethylene glycol)			
Mol. weight	20000 Da		

4.5. Poly(Amino Acids) and Poly(Oxazolines) for Click Chemistry

		Product details
PGA1085 Prg-PGA(20) Propargyl-poly(L-glutamic acid) sodium salt Mol. weight 3000 Da		
PGA1095 Prg-PGA(100) Propargyl-poly(L-glutamic acid) sodium salt Mol. weight 15000 Da		
PGA1125 N₃-PGA(20) Azido-ethyltri(ethylene glycol)-poly(L-glutamic acid) sodium salt Mol. weight 3000 Da		
PGA1130 N₃-PGA(50) Azido-ethyltri(ethylene glycol)-poly(L-glutamic acid) sodium salt Mol. weight 7500 Da		
PGA1135 N₃-PGA(100) Azido-ethyltri(ethylene glycol)-poly(L-glutamic acid) sodium salt Mol. weight 15000 Da		

[↑ back to content](#)

		Product details
PGA1205 <i>n</i> Bu-PGA(200)[Prg(20)] n-Butyl-poly(L-glutamic acid gamma-propargyl amide) sodium salt (10-20mol% propargyl substitution) Mol. weight 30000 Da		
PGA1290 <i>n</i> Bu-PGA(20)[PEG2-N ₃ (10% mod)] n-Butyl-poly(L-glutamic acid gamma-azido-ethyltri(ethylene glycol) amide) sodium salt (10-20mol% azido substitution) Mol. weight 3600 Da		
PGA1295 <i>n</i> Bu-PGA(50)[PEG2-N ₃ (10% mod)] n-Butyl-poly(L-glutamic acid gamma-azido-ethyltri(ethylene glycol) amide) sodium salt (10-20mol% azido substitution) Mol. weight 9100 Da		
PGA1300 <i>n</i> Bu-PGA(100)[PEG2-N ₃ (10% mod)] n-Butyl-poly(L-glutamic acid gamma-azido-ethyltri(ethylene glycol) amide) sodium salt (10-20mol% azido substitution) Mol. weight 18300 Da		
POX1200 Me-PMeOx(50)-N ₃ alpha-Methyl-poly(2-methyl-2-oxazoline)-omega-azide CAS-No. 26375-28-0 Formula CH ₃ (C ₄ H ₇ NO)50N ₃ Mol. weight 4300 Da		
POX1210 Me-PMeOx(100)-N ₃ alpha-Methyl-poly(2-methyl-2-oxazoline)-omega-azide CAS-No. 26375-28-0 Formula CH ₃ (C ₄ H ₇ NO)100N ₃ Mol. weight 8500 Da		

		Product details	
POX2200	Me-PEtOx(50)-N₃		
alpha-Methyl-poly(2-ethyl-2-oxazoline)-omega-azide			
CAS-No.	25805-17-8		
Formula	CH ₃ (C ₅ H ₉ NO)50N ₃		
Mol. weight	5000 Da		
			
POX2210	Me-PEtOx(100)-N₃		
alpha-Methyl-poly(2-ethyl-2-oxazoline)-omega-azide			
CAS-No.	25805-17-8		
Formula	CH ₃ (C ₅ H ₉ NO)100N ₃		
Mol. weight	5000 Da		
			
POX2250	HOOC-PEtOx(50)-N₃		
alpha-Carboxymethyl-poly(2-ethyl-2-oxazoline)-omega-azide			
CAS-No.	25805-17-8		
Formula	HOCOCH ₂ (C ₅ H ₉ NO)50N ₃		
Mol. weight	5000 Da		
			
POX1250	HOOC-PMeOx(50)-N₃		
alpha-Carboxymethyl-poly(2-methyl-2-oxazoline)-omega-azide			
CAS-No.	26375-28-0		
Formula	HOCOCH ₂ (C ₄ H ₇ NO)50N ₃		
Mol. weight	4300 Da		
			


[The Click Reaction](#)
[Amino Acid Derivatives and Related Building Blocks for Click Chemistry](#)
[Spermines and Amines for Click Chemistry](#)
[Click Reagents for Drug Delivery](#)
[Click Chemistry Tools for Proteomics](#)
[Carbohydrates for Click Chemistry](#)
[Proteolysis Targeting Chimeras \(PROTACs[®]\)](#)
[Index](#)
[↑ back to content](#)

5. Click Chemistry Tools for Proteomics


5.1. Indocyanine Green Dyes for Click Chemistry

Indocyanine green (ICG) dye, a material approved by the FDA for various applications, is a powerful tool for imaging in live cells and tissues.




Iris Biotech offers a series of ICG dyes functionalized with various clickable moieties, such as tetrazine, alkyne, azide or DBCO. Moreover, we offer ICG-equipped with different popular functional groups for conjugation, e.g., maleimide, 2-cyanobenzothiazole (CBT), and NHS.



ICG exhibits an absorption maximum in the near infrared region (NIR) at ca. 800 nm with a slight absorption in the visible range, resulting in a low auto-fluorescence. The emission maximum is at 810 nm. This absorption/emission profile allows for tissue-penetrating excitation without causing tissue damage. Consequently, ICG has found use in fields as diverse as angiography, detection of solid tumours and fluorescence image-guided surgery.



ICG Derivatives for Click Chemistry

		Product details
RL-2840	ICG-azide	
Indocyanine green azide		
Formula	$C_{48}H_{56}N_6O_4S$	
Mol. weight	813,06 g/mol	
		
RL-2880	ICG-alkyne	
Indocyanine green alkyne		
CAS-No.	1622335-41-4	
Formula	$C_{48}H_{53}N_3O_4S$	
Mol. weight	768,02 g/mol	
		
RL-2870	ICG-DBCO	
Indocyanine green dibenzotriazacyclooctyne		
Formula	$C_{63}H_{64}N_4O_5S$	
Mol. weight	989,27 g/mol	
		

References:


- Near-infrared fluorescence: application to in vivo molecular imaging; S. A. Hilderbrand, R. Weissleder; *Curr Opin Chem Biol* 2010; **14**: 71-9. <https://doi.org/10.1016/j.cbpa.2009.09.029>
- NIR dyes for bioimaging applications; J. O. Escobedo, O. Rusin, S. Lim, R. M. Strongin; *Curr Opin Chem Biol* 2010; **14**: 64-70. <https://doi.org/10.1016/j.cbpa.2009.10.022>
- In vivo molecular imaging of cancer with a quenching near-infrared fluorescent probe using conjugates of monoclonal antibodies and indocyanine green; M. Ogawa, N. Kosaka, P. L. Choyke, H. Kobayashi; *Cancer Res* 2009; **69**: 1268-72. <https://doi.org/10.1158/0008-5472.CAN-08-3116>

5.2. Clickable Linkers for Selective Protein Labeling




One way to selectively label a protein is to recombinantly express a modified version containing a sequence that selectively reacts with a specific linker type. Two examples of this approach are the His Tag acylation and the HaloTag®.

In the His Tag acylation approach, a N-terminal GlyHis₆ tag attached to a protein of interest selectively reacts with a 4-methoxyphenyl ester, generating an acylated N-terminus. While 4-methoxyphenyl esters are too unreactive to undergo acylation with any other primary amine, a proximal imidazole in the Gly-His₆ sequence acts as a catalyst to facilitate selective acylation of the N-terminal glycine.

For the acylation of lysine, Jensen *et al.* developed the peptide sequence His_n-Lys-His_m (Lys-His tag) that directs the acylation of the designated lysine N-ε amine under mild conditions and with high selectivity over native lysine residues.

		Product details
RL-3010	N₃Ac-OPhOMe	
4-Methoxyphenyl 2-azidoacetate		
CAS-No.	2546513-31-7	
Formula	C ₉ H ₉ N ₃ O ₃	
Mol. weight	207,19 g/mol	
		






The HaloTag® is a protein tag whose sequence can easily be fused to a gene coding for a protein of interest. Functionally, it is a haloalkane dehalogenase that binds and forms covalent bonds to specific halogenated ligands. Those ligands are composed of two parts: a chloroalkane linker that forms the bond with HaloTag® protein, and a functional group or affinity handle. A HaloTag®-containing fusion protein is thus able to selectively label itself with an appropriate haloalkane dehalogenase ligand.

			Product details
RL-3640	Halo-PEG(5)-azide		
1-azido-21-chloro-3,6,9,12,15-pentaoxahenicosane			
CAS-No.	1261238-21-4		
Formula	C ₁₆ H ₃₂ ClN ₃ O ₅		
Mol. weight	381,90 g/mol		
RL-3670	Halo-DBCO		
N-[2-[2-[(6-chlorohexyl)oxy]ethoxy]ethyl]-gamma-oxo-dibenz[b,f]azocine-5(6H)-butanamide			
CAS-No.	1808119-16-5		
Formula	C ₂₉ H ₃₅ ClN ₂ O ₄		
Mol. weight	511,06 g/mol		
RL-3700	Halo-PEG(2)-Azide		
1-Azido-12-chloro-3,6-dioxadodecane			
CAS-No.	2568146-55-2		
Formula	C ₁₀ H ₂₀ ClN ₃ O ₂		
Mol. weight	249,74 g/mol		
RL-3710	Halo-PEG(4)-Azide		
1-Azido-18-chloro-3,6,9,12-tetraoxaoctadecane			
CAS-No.	2989398-83-4		
Formula	C ₁₄ H ₂₈ ClN ₃ O ₄		
Mol. weight	337,85 g/mol		

Reference:

- *Direct pH measurements by using subcellular targeting of 5(and 6-) carboxysemaphthorhodafuor in mammalian cells; H. A. Benink, M. G. McDougall, D. H. Klaubert, G. V. Los; **Biotechniques** 2018; **47(3)**: 769-774.*
<https://doi.org/10.2144/000113220>


Besides HaloTag®, we also offer clickable SNAP-tag® and CLIP-tag™ ligands. The SNAP-tag® is a 20 kDa self-labeling protein tag bearing a cysteine moiety that undergoes an irreversible reaction with synthetic O6-benzylguanine (BG) derivatives, resulting in a covalent thioether bond. The CLIP-tag™ is a modified version of the SNAP-tag®, engineered to react with benzylcytosine (BC) instead of benzylguanine (BG).

			Product details
RL-3930	Alkyne-SNAP		
N-(4-(((2-amino-9H-purin-6-yl)oxy)methyl)benzyl)pent-4-ynamide			
CAS-No.	1104822-07-2		
Formula	C ₁₈ H ₁₈ N ₆ O ₂		
Mol. weight	350,38 g/mol		
RL-3940	Alkyne-PEG(5)-SNAP		
N-(4-(((2-amino-9H-purin-6-yl)oxy)methyl)benzyl)-4,7,10,13,16-pentaoxononadec-18-ynamide			
Formula	C ₂₇ H ₃₆ N ₆ O ₇		
Mol. weight	556,62 g/mol		
RL-3950	Azide-SNAP		
N-(4-(((2-amino-9H-purin-6-yl)oxy)methyl)benzyl)-6-azidohexanamide			
Formula	C ₁₉ H ₂₃ N ₉ O ₂		
Mol. weight	409,45 g/mol		
RL-3960	Azide-PEG(4)-SNAP		
N-(4-(((2-amino-9H-purin-6-yl)oxy)methyl)benzyl)-1-azido-3,5,7,9-tetraoxadodecan-12-amide			
Formula	C ₂₄ H ₃₃ N ₉ O ₆		
Mol. weight	543,59 g/mol		
RL-4010	DBCO-SNAP		
			
Formula	C ₃₄ H ₃₁ N ₇ O ₃		
Mol. weight	585,67 g/mol		




References:

- Site-specific protein labeling with SNAP-Tags; N. B. Cole; **Curr Protoc Protein Sci.** 2013; **73(30)**: 1-30.
<https://doi.org/10.1002/0471140864.ps3001s73>
- Directed evolution of O6-alkylguanine-DNA alkyltransferase for efficient labeling of fusion proteins with small molecules in vivo; A. Juillerat, T. Gronemeyer, A. Keppler, S. Gendreizig, H. Pick, H. Vogel, K. Johnsson; **Chem. Biol.** 2003; **10(4)**: 313-317. [https://doi.org/10.1016/s1074-5521\(03\)00068-1](https://doi.org/10.1016/s1074-5521(03)00068-1)
- Site-specific, Covalent Labeling of Recombinant Antibody Fragments via Fusion to an Engineered Version of 6-O-Alkylguanine DNA Alkyltransferase; F. Kampmeier, M. Ribbert, T. Nachreiner, S. Dembski, F. Beaufils, A. Brecht, S. Barth; **Bioconjugate Chem.** 2009; **20(5)**: 1010-1015. <https://doi.org/10.1021/bc9000257>
- SNAP-Tag Technology: A Useful Tool to Determine Affinity Constants and Other Functional Parameters of Novel Antibody Fragments; J. Niesen, M. Sack, M. Seidel, R. Fendel, S. Barth, R. Fischer, C. Stein; **Bioconjugate Chem.** 2016; **27(8)**: 1931-1941. <https://doi.org/10.1021/acs.bioconjchem.6b00315>
- Snap-, CLIP- and Halo-Tag Labelling of Budding Yeast Cells; F. Stagge, G. Y. Mitronova, V. N. Belov, C. A. Wurm, S. Jakobs; **PLoS ONE** **8(10)**: e78745. <https://doi.org/10.1371/journal.pone.0078745>

In recent years, the so-called cyanobenzothiazole (CBT) click condensation has gained increased interest. It is highly biocompatible and controllable and runs under physiological conditions. In the CBT click reaction, the nitrile function of 2-cyanobenzothiazole (2-CBT) reacts with a β -aminothiol (1,2 aminothiol) and forms a 4,5-dihydro-1,3 thiazole. The β -aminothiol may originate from, e.g., an N-terminal D- or L-cysteine residue, which can be easily introduced recombinantly or during SPPS.




[↑ back to content](#)

		Product details
RL-4280 4-Pentynoyl-Gly-CBT N-(2-((2-cyanobenzo[d]thiazol-6-yl)amino)-2-oxoethyl)pent-4-ynamide Formula $C_{15}H_{12}N_4O_2S$ Mol. weight 312,35 g/mol		
RL-4290 6-Azidopentanoyl-Gly-CBT 5-azido-N-(2-((2-cyanobenzo[d]thiazol-6-yl)amino)-2-oxoethyl)pentanamide Formula $C_{15}H_{15}N_7O_2S$ Mol. weight 357,39 g/mol		
RL-4310 DBCO-Suc-CBT N1-(2-cyanobenzo[d]thiazol-6-yl)-N4-(3-(11,12-didehydro-5,6-dihydro-dibenzo[b,f]azocin-yl)-3-oxopropyl)succinamide Formula $C_{30}H_{23}N_5O_3S$ Mol. weight 533,61 g/mol		




References:

- *N, S-Double labeling of N-terminal cysteines via an alternative conjugation pathway with 2-cyanobenzo-thiazole*; W. Wang, J. Gao; **Org. Chem.** 2020; **85**: 1756-1763. <https://doi.org/10.1021/acs.joc.9b02959>
- *Thiol-cyanobenzothiazole ligation for the efficient preparation of peptide-PNA conjugates*; N. Patil, J. Karas, B. Turner, F. Shabanpoor; **Bioconjugate Chem.** 2019; **30**: 73-799.
<https://doi.org/10.1021/acs.bioconjchem.8b00908>
- *CBT-Cys click reaction for optical bioimaging in vivo*. X. Hu, R. Tang, L. Bai, S. Liu, G. Liang, X. Sun; **View.** 2023; 20220065. <https://doi.org/10.1002/VIW.20220065>

While nitriles in general may react with sulfhydryl groups (e.g., from cysteines) in a rather unspecific way, the advantage of α -cyanopyridines (2-picolinonitrils) is their strong selectivity for 1,2-aminothiols (e.g., provided as N-terminal cysteine or internally as non-canonical amino acid). Thus, this reaction may be used for the synthesis of cyclic peptides or to fuse ligands to peptide chains.



This click-like reaction is biorthogonal, biocompatible, catalyst-free, and selective; it proceeds readily in aqueous solutions at physiological pH and ambient temperature. The formation of the resulting 2-thiazoline (alternative name: 2,5-dihydrothiazole) is driven by the release of ammonia. The reaction product is stable at physiological conditions.

		Product details
RL-8625 N₃-PEG(4)-CINA N-(14-azido-3,6,9,12-tetraoxatetradecyl)-2-cyanoisonicotinamide Formula C ₁₇ H ₂₄ N ₆ O ₅ Mol. weight 392,42 g/mol		
RL-8630 Alkyne-PEG(4)-CINA 2-cyano-N-(3,6,9,12-tetraoxapentadec-14-yn-1-yl)isonicotinamide Formula C ₁₈ H ₂₃ N ₃ O ₅ Mol. weight 361,40 g/mol		
RL-8640 DBCO-PEG(4)-CINA Dibenzoazacyclooctyne-tetra(ethylene glycol)-cyanoisonicotinamide Formula C ₃₆ H ₃₉ N ₅ O ₇ Mol. weight 653,74 g/mol		






Find many more derivatives in our webshop!

References:

- *The Cyanopyridine-Aminothiol Click Reaction: Expanding Horizons in Chemical Biology*; C. Nitsche; **SynLett.** 2024; **35**: A-E. <https://dx.doi.org/10.1055/a-2214-7612>
- *Biocompatible and Selective Generation of Bicyclic Peptides*; S. Ullrich, J. George, A. Coram, R. Morewood, C. Nitsche; **Angew. Chem Int. Ed.** 2022; **61(43)**: e20228400. <https://doi.org/10.1002/anie.202208400>
- *Tobacco Etch Virus protease: A shortcut across biotechnologies*; F. Cesaratto, O. Burrone, G. Petris; **J Biotechnol.** 2016; **231(10)**: 239-249. <https://doi.org/10.1016/j.jbiotec.2016.06.012>


Besides, we offer click-functionalized biotin derivatives. Biotinylation represents a versatile tool for purification, immobilization, and labeling. It is a rather small molecule compared to other labels, thus minimizing its steric impact on the biotinylated carrier itself.

[↑ back to content](#)

		Product details
LS-4660 Biotinyl-DAPe-N₃ N-(5-azidopentyl)-5-((3aS,4S,6aR)-2-oxohexahydro-1H-thieno[3,4-d]imidazol-4-yl)pentanamide CAS-No. 1349190-76-6 Formula C ₁₅ H ₂₆ N ₆ O ₂ S Mol. weight 354,47 g/mol		
LS-4210 Biotin-N₃ N-(3-azidopropyl)-5-((3aS,4S,6aR)-2-oxohexahydro-1H-thieno[3,4-d]imidazol-4-yl)pentanamide CAS-No. 908007-17-0 Formula C ₁₃ H ₂₂ N ₆ O ₂ S Mol. weight 326,42 g/mol		
LS-4300 DecarboxyBiotin-N₃ (3aS,4S,6aR)-4-(5-Azidopentyl)tetrahydro-1H-Thieno[3,4-d]imidazol-2(3H)-one CAS-No. 1260586-88-6 Formula C ₁₀ H ₁₇ N ₅ OS Mol. weight 255,34 g/mol		



Find many more biotin derivatives
in our dedicated brochure!



6. Carbohydrates for Click Chemistry

Glycoconjugates, i.e. glycans linked to proteins or lipids, are an essential part of all living organisms. In higher organisms, but also in lower eukaryotes and some bacteria and archaea, many proteins are post-translationally modified by linking oligosaccharides to amino acid side-chains, forming glycoproteins. Glycosylation is the most complex posttranslational modification and can be observed on membrane proteins, secreted proteins and peptides, or proteins in the cytosol and nucleus.

Glycoconjugates display a multitude of biological effects from protein folding and stabilization, to cell surface interaction through molecular recognition motifs for cell-cell communication, and structural support and protection.

Abnormal glycosylation patterns can be observed in pathological conditions such as neurodegenerative diseases or tumor growth and metastasis. Moreover, glycosylation patterns play a decisive role in the infection pathways of and the immune response against many pathogens, further underlining the importance of this type of modification.

Synthetic glycoconjugates are interesting targets for the investigation of immunogenicity, infection pathways or structure activity relationships, and for the development of novel drugs and vaccines. Carbohydrates functionalized for Click chemistry provide mild and selective access to such glycoconjugates.










Fig. 21: Simplified representation of a eukaryotic cell and its cell surface glycans.





6.1. Galactose Derivatives

			Product details
GBB1445	alpha-GalNAc-N₃		
1-O-(2-Azidoethoxy)-2-acetamido-2deoxy-alpha-D-galactopyranoside			
CAS-No.	195384-42-0		
Formula	C ₁₀ H ₁₈ N ₄ O ₆		
Mol. weight	290,27 g/mol		
			
GBB1370	alpha-GalNAc-TEG-N₃		
1-O-(2-(2-(2-Azidoethoxy)ethoxy)ethoxy)-2-acetamido-2deoxy-alpha-D-galactopyranoside			
CAS-No.	882873-70-3		
Formula	C ₁₄ H ₂₆ N ₄ O ₈		
Mol. weight	378,38 g/mol		
			
GBB1430	beta-Gal-Et-N₃		
1-(2-Azidoethoxy)-beta-D-galactopyranose			
CAS-No.	151651-54-6		
Formula	C ₈ H ₁₅ N ₃ O ₆		
Mol. weight	249,22 g/mol		
			
GBB1380	beta-Gal-TEG-N₃		
1-O-(2-(2-(2-Azidoethoxy)ethoxy)ethoxy)-beta-D-galactopyranoside			
CAS-No.	126765-27-3		
Formula	C ₁₂ H ₂₃ N ₃ O ₈		
Mol. weight	337,33 g/mol		
			
GBB1375	alpha-GalNAc-TEG-Alkyne		
1-O-(2-(2-(2-(Prop-2-ynyloxy)ethoxy)ethoxy)ethoxy)-2-acetamido-2deoxy-alpha-D-galactopyranoside			
Formula	C ₁₇ H ₂₉ NO ₉		
Mol. weight	391,41 g/mol		
			


		Product details
GBB1385	beta-Gal-TEG-Alkyne	 
CAS-No.	1072903-79-7	
Formula	C ₁₅ H ₂₆ O ₉	
Mol. weight	350,36 g/mol	

[The Click Reaction](#)
[Amino Acid Derivatives and Related Building Blocks for Click Chemistry](#)
[Spermines and Amines for Click Chemistry](#)
[Click Reagents for Drug Delivery](#)
[Click Chemistry Tools for Proteomics](#)
[Carbohydrates for Click Chemistry](#)
[Proteolysis Targeting Chimeras \(PROTACs[®]\)](#)
[Index](#)




6.2. Glucose Derivatives

		Product details
GBB1435	beta-Glc-N₃	 
CAS-No.	165331-08-8	
Formula	C ₈ H ₁₅ N ₃ O ₆	
Mol. weight	249,22 g/mol	
GBB1390	beta-Glc-TEG-N₃	 
CAS-No.	156058-83-2	
Formula	C ₁₂ H ₂₃ N ₃ O ₈	
Mol. weight	337,33 g/mol	
GBB1395	beta-Glc-TEG-Alkyne	 
CAS-No.	1072903-76-4	
Formula	C ₁₅ H ₂₆ O ₉	
Mol. weight	350,36 g/mol	
GBB1400	beta-GlcNAc-TEG-N₃	 
1-O-(2-(2-(2-Azidoethoxy)ethoxy)ethoxy)-2-acetamido-2-deoxy-beta-D-glucopyranoside		
CAS-No.	86520-54-9	
Formula	C ₁₄ H ₂₆ N ₄ O ₈	
Mol. weight	378,73 g/mol	




[↑ back to content](#)

		Product details
GBB1405 beta-GlcNAc-TEG-Alkyne		
1-O-(2-(2-(2-(Prop-2-ynyloxy)ethoxy)ethoxy)ethoxy)-2-acetamido-2-deoxy-beta-D-glucopyranoside		
CAS-No.	1884184-44-4	
Formula	C ₁₇ H ₂₉ NO ₉	
Mol. weight	391,41 g/mol	
		

6.3. Mannose Derivatives

		Product details
GBB1440	alpha-Man-N₃	
1-(2-Azidoethoxy)-alpha-D-mannopyranose		
CAS-No.	155196-97-7	
Formula	C ₈ H ₁₅ N ₃ O ₆	
Mol. weight	249,22 g/mol	
		
GBB1420	alpha-Man-TEG-N₃	
1-O-(2-(2-(2-Azidoethoxy)ethoxy)ethoxy)-alpha-D-mannopyranoside		
CAS-No.	246855-76-5	
Formula	C ₁₂ H ₂₃ N ₃ O ₈	
Mol. weight	337,33 g/mol	
		
GBB1425	alpha-Man-TEG-Alkyne	
1-O-(2-(2-(2-(Prop-2-ynyloxy)ethoxy)ethoxy)ethoxy)-alpha-D-mannopyranoside		
CAS-No.	1072903-80-0	
Formula	C ₁₅ H ₂₆ O ₉	
Mol. weight	350,36 g/mol	
		

6.4. Lactose Derivatives

		Product details
GBB1455 beta-Lac-EO-N₃ 2-Azidoethyl 4-O-beta-D-galactopyranosyl-beta-(1->4)-D-glucopyranoside CAS-No. 230286-11-0 Formula C ₁₄ H ₂₅ N ₃ O ₁₁ Mol. weight 411,36 g/mol		
GBB1410 beta-Lac-TEG-N₃ (2-(2-(2-Azidoethoxy)ethoxy)ethyl) 4-O-beta-D-galactopyranosyl-beta-(1->4)-D-glucopyranoside CAS-No. 246855-74-3 Formula C ₁₈ H ₃₃ N ₃ O ₁₃ Mol. weight 499,47 g/mol		
GBB1415 beta-Lac-TEG-Alkyne CAS-No. 1442747-66-1 Formula C ₂₁ H ₃₆ O ₁₄ Mol. weight 512,5 g/mol		

References:

- *Nanovesicles displaying functional linear and branched oligomannose self-assembled from sequence-defined Janus glycodendrimers*; Q. Xiao, M. Delbianco, S. E. Sherman, A. M. Reveron Perez, P. Bharate, A. Pardo-Vargas, C. Rodriguez-Emmenegger, N. Y. Kostina, K. Rahimi, D. Soder, M. Moller, M. L. Klein, P. H. Seeberger, V. Percec; **Proc Natl Acad Sci U S A** 2020; 202003938. <https://doi.org/10.1073/pnas.2003938117>
- *Glycans in drug discovery*; P. Valverde, A. Ardá, N.-C. Reichardt, J. Jiménez-Barbero, A. Gimeno; **MedChemComm** 2019. <https://doi.org/10.1039/c9md00292h>
- *Stepwise orthogonal click chemistry toward fabrication of paclitaxel/galactose functionalized fluorescent nanoparticles for HepG2 cell targeting and delivery*; C. H. Lai, T. C. Chang, Y. J. Chuang, D. L. Tzou, C. C. Lin; **Bio-conjug Chem** 2013; **24**: 1698-709. <https://doi.org/10.1021/bc400219t>
- *Design of multivalent galactosyl carborane as a targeting specific agent for potential application to boron neutron capture therapy*; C. H. Lai, Y. C. Lin, F. I. Chou, C. F. Liang, E. W. Lin, Y. J. Chuang, C. C. Lin; **Chem Commun (Camb)** 2012; **48**: 612-4. <https://doi.org/10.1039/c1cc14447b>
- *Essentials of Glycobiology*; A. Varki, R. Schauer, R. D. Cummings, J. D. Esko, H. H. Freeze, P. Stanley, C. R. Bertozzi, G. W. Hart, M. E. Etzler; 2009.
- *Role of glycosylation in development*; R. S. Haltiwanger, J. B. Lowe; **Annu. Rev. Biochem.** 2004; **73**: 491-537. <https://doi.org/10.1146/annurev.biochem.73.011303.074043>
- *Glycosylation and the immune system*; P. M. Rudd, T. Elliott, P. Cresswell, I. A. Wilson, R. A. Dwek; **Science** 2001; **291**: 2370-6. <https://doi.org/10.1126/science.291.5512.2370>


[↑ back to content](#)

- Glycosylation in cellular mechanisms of health and disease; K. Ohtsubo, J. D. Marth; **Cell** 2006; **126**: 855-67.
<https://doi.org/10.1016/j.cell.2006.08.019>
- The role of selectins in inflammation and disease; K. Ley; **Trends Mol Med** 2003; **9**: 263-8.
- Carbohydrate diversity: synthesis of glycoconjugates and complex carbohydrates; A. Holemann, P. H. Seeberger; **Curr Opin Biotechnol** 2004; **15**: 615-22. <https://doi.org/10.1016/j.copbio.2004.10.001>
- Post-translational modifications in the context of therapeutic proteins; G. Walsh, R. Jefferis; **Nat. biotechnol.** 2006; **24**: 1241-52. <https://doi.org/10.1038/nbt1252>
- Glycomics: a pathway to a class of new and improved therapeutics; Z. Shriver, S. Raguram, R. Sasisekharan; **Nat Rev Drug Discov** 2004; **3**: 863-73. <https://doi.org/10.1038/nrd1521>
- The bittersweet promise of glycobiology; A. Dove; **Nat. biotechnol.** 2001; **19**: 913-7.
<https://doi.org/10.1038/nbt1001-913>
- Carbohydrate and protein immobilization onto solid surfaces by sequential Diels-Alder and azide-alkyne cycloadditions; X. L. Sun, C. L. Stabler, C. S. Cazalis, E. L. Chaikof; **Bioconjug Chem** 2006; **17**: 52-7.
<https://doi.org/10.1021/bc0502311>
- Covalent display of oligosaccharide arrays in microtiter plates; M. C. Bryan, F. Fazio, H. K. Lee, C. Y. Huang, A. Chang, M. D. Best, D. A. Calarese, O. Blixt, J. C. Paulson, D. Burton, I. A. Wilson, C. H. Wong; **J Am Chem Soc** 2004; **126**: 8640-1. <https://doi.org/10.1021/ja048433f>
- Multivalent, bifunctional dendrimers prepared by click chemistry; P. Wu, M. Malkoch, J. N. Hunt, R. Vestberg, E. Kaltgrad, M. G. Finn, V. V. Fokin, K. B. Sharpless, C. J. Hawker; **Chem Commun (Camb)** 2005: 5775-7.
<https://doi.org/10.1039/b512021g>
- A chemoenzymatic approach to glycopeptide antibiotics; H. Lin, C. T. Walsh; **J Am Chem Soc** 2004; **126**: 13998-4003. <https://doi.org/10.1021/ja045147v>



**You cannot find the derivative
you are looking for?**

Get in contact for a custom synthesis!



7. Proteolysis Targeting Chimeras (PROTACs®)

Targeted protein degradation *via* proteolysis-targeting chimeras (PROTACs) is an emerging attempt to cure diseases caused by the irregular expression of certain disease-causing proteins. Such protein degraders act as bifunctional linkers and allow to feed the protein of interest (POI) to the cell's Ubiquitin-Proteasome system, thus, to eliminate the malexpressed proteins. These PROTACs consist of three components: one ligand with high affinity for E3 ubiquitin ligase, another one with high affinity for the POI and an appropriate cross-linker joining both ligands. This linker can also be used to increase the solubility, if needed, e.g., by incorporation of PEGs. The resulting proximity of both, the recruited POI and the E3 ligase, allows the polyubiquitination of the POI by the E3 associated E2 enzyme. This leads to a labeling of the POI for degradation through the proteasome.

"PROTAC® is a registered trademark of Arvinas Operations, Inc., and is used under license."

Mode of action:

1. A cross-linker unites the POI ligand and E3 ligase ligand = PROTAC.
2. The three-component PROTAC recruits the POI and the E2-associated E3 ligase *via* the respective ligands = Ternary complex.
3. Several Ubiquitins are added to lysine residues of the POI = Polyubiquitination.
4. The ubiquitinated POI is degraded by the proteasome.















Fig. 22: Targeted protein degradation via proteolysis-targeting chimeras.

To construct a suitable PROTAC, we provide a variety of E3 ubiquitin ligase ligands in combination with linkers of various length and an elective amino-, carboxyl-, click- or thiol-reactive end ("Partial PROTACs").

Click-Reactive Partial PROTACs

			Product details
PTC1400	Pomalidomide-PEG1-Alkyne		
N-(2-(2,6-Dioxopiperidin-3-yl)-1,3-dioxisoindolin-4-yl)-3-(prop-2-yn-1-yloxy)propanamide			
CAS-No.	2236109-19-4		
Formula	C ₁₉ H ₁₇ N ₃ O ₆		
Mol. weight	383,35 g/mol		
PTC1410	Pomalidomide-PEG2-Alkyne		
N-(2-(2,6-Dioxopiperidin-3-yl)-1,3-dioxisoindolin-4-yl)-3-(2-(prop-2-yn-1-yloxy)ethoxy)propanamide			
CAS-No.	2243000-25-9		
Formula	C ₂₁ H ₂₁ N ₃ O ₇		
Mol. weight	427,41 g/mol		
PTC1420	Pomalidomide-PEG3-Alkyne		
N-(2-(2,6-Dioxopiperidin-3-yl)-1,3-dioxisoindolin-4-yl)-3-(2-(2-(prop-2-yn-1-yloxy)ethoxy)ethoxy)propanamide			
CAS-No.	2236109-20-7		
Formula	C ₂₃ H ₂₅ N ₃ O ₈		
Mol. weight	471,46 g/mol		
PTC1440	Pomalidomide-PEG5-Alkyne		
N-(2-(2,6-Dioxopiperidin-3-yl)-1,3-dioxisoindolin-4-yl)-4,7,10,13,16-pentaoxanonadec-18-ynamide			
Formula	C ₂₇ H ₃₃ N ₃ O ₁₀		
Mol. weight	559,57 g/mol		
PTC1460	(S,R,S)-AHPC-PEG1-Alkyne		
(2S,4R)-1-((S)-3,3-Dimethyl-2-(3-(prop-2-yn-1-yloxy)propanamido)butanoyl)-4-hydroxy-N-(4-(4-methylthiazol-5-yl)benzyl)pyrrolidine-2-carboxamide			
CAS-No.	2242965-71-3		
Formula	C ₂₈ H ₃₆ N ₄ O ₅ S		
Mol. weight	540,67 g/mol		

		Product details
PTC1470 (S,R,S)-AHPC-PEG2-Alkyne (2S,4R)-1-((S)-3,3-Dimethyl-2-(3-(2-(prop-2-yn-1-yloxy)ethoxy)propanamido)butanoyl)-4-hydroxy-N-(4-(4-methylthiazol-5-yl)benzyl)pyrrolidine-2-carboxamide CAS-No. 2242965-72-4 Formula $C_{30}H_{40}N_4O_6S$ Mol. weight 584,73 g/mol		
PTC1480 (S,R,S)-AHPC-PEG3-Alkyne (2S,4R)-1-((S)-2-(tert-Butyl)-4-oxo-7,10,13-trioxa-3-aza-hexadec-15-ynoyl)-4-hydroxy-N-(4-(4-methylthiazol-5-yl)benzyl)pyrrolidine-2-carboxamide CAS-No. 2374122-30-0 Formula $C_{32}H_{44}N_4O_7S$ Mol. weight 628,78 g/mol		
PTC1490 (S,R,S)-AHPC-PEG4-Alkyne (2S,4R)-1-((S)-2-(tert-Butyl)-4-oxo-7,10,13,16-tetraoxa-3-azanonadec-18-ynoyl)-4-hydroxy-N-(4-(4-methylthiazol-5-yl)benzyl)pyrrolidine-2-carboxamide CAS-No. 2267282-21-1 Formula $C_{34}H_{48}N_4O_8S$ Mol. weight 672,83 g/mol		
PTC1500 (S,R,S)-AHPC-PEG5-Alkyne (2S,4R)-1-((S)-2-(tert-Butyl)-4-oxo-7,10,13,16,19-pentaoxa-3-azadocos-21-ynoyl)-4-hydroxy-N-(4-(4-methylthiazol-5-yl)benzyl)pyrrolidine-2-carboxamide CAS-No. 2817805-63-1 Formula $C_{36}H_{52}N_4O_9S$ Mol. weight 716,88 g/mol		
PTC1510 (S,R,S)-AHPC-PEG6-Alkyne (2S,4R)-1-((S)-2-(tert-Butyl)-4-oxo-7,10,13,16,19,22-hexaoxa-3-azapentacos-24-ynoyl)-4-hydroxy-N-(4-(4-methylthiazol-5-yl)benzyl)pyrrolidine-2-carboxamide Formula $C_{38}H_{56}N_4O_{10}S$ Mol. weight 760,94 g/mol		
PTC1520 Pomalidomid- PEG1-N₃ 2-(2-Azidoethoxy)-N-(2-(2,6-dioxopiperidin-3-yl)-1,3-dioxoisindolin-4-yl)acetamide CAS-No. 2133360-04-8 Formula $C_{17}H_{16}N_6O_6$ Mol. weight 400,35 g/mol		

[The Click Reaction](#)
[Amino Acid Derivatives and Related Building Blocks for Click Chemistry](#)
[Spermines and Amines for Click Chemistry](#)
[Click Reagents for Drug Delivery](#)
[Click Chemistry Tools for Proteomics](#)
[Carbohydrates for Click Chemistry](#)
[Proteolysis Targeting Chimeras \(PROTACs*\)](#)
[Index](#)
[↑ back to content](#)


PTC1530 Pomalidomid- PEG2-N₃

2-(2-(2-Azidoethoxy)ethoxy)-N-(2-(2,6-dioxopiperidin-3-yl)-1,3-dioxoisindolin-4-yl)acetamide

CAS-No. 2267306-14-7

Formula C₁₉H₂₀N₆O₇

Mol. weight 444,4 g/mol


**PTC1540 Pomalidomid- PEG3-N₃**

2-(2-(2-(2-azidoethoxy)ethoxy)ethoxy)-N-(2-(2,6-dioxopiperidin-3-yl)-1,3-dioxoisindolin-4-yl)acetamide

CAS-No. 2267306-15-8

Formula C₂₁H₂₄N₆O₈

Mol. weight 488,45 g/mol


**PTC1560 Pomalidomid-C6-PEG3-butyl-N₃**

6-(2-(2-((6-Azidohexyl)oxy)ethoxy)ethoxy)-N-(2-(2,6-dioxopiperidin-3-yl)-1,3-dioxoisindolin-4-yl)hexanamide

CAS-No. 2300178-66-7

Formula C₂₉H₄₀N₆O₈

Mol. weight 600,66 g/mol


**PTC1570 Pomalidomid-C6-PEG1-C3-PEG1-butyl-N₃**

6-((5-((6-Azidohexyl)oxy)pentyl)oxy)-N-(2-(2,6-dioxopiperidin-3-yl)-1,3-dioxoisindolin-4-yl)hexanamide

CAS-No. 2300178-65-6

Formula C₃₀H₄₂N₆O₇


Mol. weight 598,69 g/mol

**PTC1580 Pomalidomid-PEG6-butyl-N₃**

4-azido-N-(2-(2,6-dioxopiperidin-3-yl)-1,3-dioxoisindolin-4-yl)-3,6,9,12,15,18-hexaoxatetracosanamide

Formula C₃₁H₄₄N₆O₁₁

Mol. weight 676,71 g/mol


**PTC1590 (S,R,S)-AHPC-PEG1-N₃**





(2S,4R)-1-((S)-2-(2-(2-Azidoethoxy)acetamido)-3,3-dimethylbutanoyl)-4-hydroxy-N-(4-(4-methylthiazol-5-yl)benzyl)pyrrolidine-2-carboxamide

CAS-No. 2101200-09-1

Formula C₂₆H₃₅N₇O₅S

Mol. weight 557,67 g/mol




		Product details
PTC1600 (S,R,S)-AHPC-PEG2-N₃ (2S,4R)-1-((S)-2-(2-(2-Azidoethoxy)ethoxy)aceta- mido)-3,3-dimethylbutanoyl)-4-hydroxy-N-(4-(4- methylthiazol-5-yl)benzyl)pyrrolidine-2-carboxamide CAS-No. 2010159-45-0 Formula C ₂₈ H ₃₉ N ₇ O ₆ S Mol. weight 601,72 g/mol		
PTC1610 (S,R,S)-AHPC-PEG3-N₃ (2S,4R)-1-((S)-14-azido-2-(tert-butyl)-4-oxo-6,9,12-tri- oxa-3-azatetradecanoyl)-4-hydroxy-N-(4-(4-methylthi- azol-5-yl)benzyl)pyrrolidine-2-carboxamide CAS-No. 1797406-80-4 Formula C ₃₀ H ₄₃ N ₇ O ₇ S Mol. weight 645,77 g/mol		
PTC1640 (S,R,S)-AHPC-PEG6-N₃ (2S,4R)-1-((S)-23-Azido-2-(tert-butyl)-4-oxo- 6,9,12,15,18,21-hexaoxa-3-azatricosanoyl)-4-hydroxy-N- (4-(4-methylthiazol-5-yl)benzyl)pyrrolidine-2-carbox- amide CAS-No. 2086298-71-5 Formula C ₃₆ H ₅₅ N ₇ O ₁₀ S Mol. weight 777,93 g/mol		
PTC1680 (S,R,S)-AHPC-PEG6-butyl-N₃ (2S,4R)-1-((S)-27-Azido-2-(tert-butyl)-4-oxo- 6,9,12,15,18,21-hexaoxa-3-azaheptacosanoyl)-4-hy- droxy-N-(4-(4-methylthiazol-5-yl)benzyl)pyrrolidi- ne-2-carboxamide Formula C ₄₀ H ₆₃ N ₇ O ₁₀ S Mol. weight 834,03 g/mol		

[The Click Reaction](#)
[Amino Acid Derivatives and Related
Building Blocks for Click Chemistry](#)
[Spermines and Amines for Click Chemistry](#)
[Click Reagents for Drug Delivery](#)
[Click Chemistry Tools for Proteomics](#)
[Carbohydrates for Click Chemistry](#)
[Proteolysis Targeting Chimeras \(PROTACs*\)](#)
[Index](#)


**For more PROTAC building blocks,
see our webshop!**


[↑ back to content](#)

In addition to these pre-designed building blocks, we offer custom synthesis of your required ligand-linker combination or “complete PROTAC”. This allows to design a library of slightly different PROTACs in order to find the best combination for your application, as even small changes in ligands and cross-linkers might affect the efficiency of the formation of the ternary complex.



References:

- Bifunctional Molecules beyond PROTACs; S. J. Conway; **J. Med. Chem.** 2020; **63**: 2802-2806.
<https://doi.org/10.1021/acs.jmedchem.0c00293>.
- Targeted protein degradation by PROTACs; T. K. Neklesa, J. D. Winkler, C. M. Crews; **Pharmacol Ther** 2017; **174**: 138-144. <https://doi.org/10.1016/j.pharmthera.2017.02.027>.
- Targeted Protein Degradation: from Chemical Biology to Drug Discovery; P. M. Cromm, C. M. Crews; **Cell Chem Biol** 2017; **24**: 1181-1190. <https://doi.org/10.1016/j.chembiol.2017.05.024>.
- Targeted Protein Degradation by Small Molecules; D. P. Bondeson, C. M. Crews; **Annu Rev Pharmacol Toxicol** 2017; **57**: 107-123. <https://doi.org/10.1146/annurev-pharmtox-010715-103507>.
- Small-Molecule PROTACs: New Approaches to Protein Degradation; M. Toure, C. M. Crews; **Angew Chem Int Ed** 2016; **55**: 1966-73. <https://doi.org/10.1002/anie.201507978>.
- Impact of linker length on the activity of PROTACs; K. Cyrus, M. Wehenkel, E. Y. Choi, H. J. Han, H. Lee, H. Swanson, K. B. Kim; **Mol Biosyst** 2011; **7**: 359-64. <https://doi.org/10.1039/c0mb00074d>.
- Development of potent monoclonal antibody auristatin conjugates for cancer therapy; S. O. Doronina, B. E. Toki, M. Y. Torgov, B. A. Mendelsohn, C. G. Cervený, D. F. Chace, R. L. DeBlanc, R. P. Gearing, T. D. Bovee, C. B. Siegall, J. A. Francisco, A. F. Wahl, D. L. Meyer, P. D. Senter; **Nat Biotechnol** 2003; **21**: 778-84.
<https://doi.org/10.1038/nbt832>

- A versatile acid-labile linker for antibody–drug conjugates; M. C. Finniss, K. S. Chu, C. J. Bowerman, J. C. Luft, Z. A. Haroon, J. M. DeSimone; **Med. Chem. Commun.** 2014; **5**: 1355-1358. <https://doi.org/10.1039/c4md00150h>
- Gemtuzumab ozogamicin, a potent and selective anti-CD33 antibody-calicheamicin conjugate for treatment of acute myeloid leukemia; P. R. Hamann, L. M. Hinman, I. Hollander, C. F. Beyer, D. Lindh, R. Holcomb, W. Hallett, H. R. Tsou, J. Upešlaciš, D. Shochat, A. Mountain, D. A. Flowers, I. Bernstein; **Bioconjug Chem** 2002; **13**: 47-58. <https://doi.org/10.1021/bc010021y>
- Antibody conjugates of 7-ethyl-10-hydroxycamptothecin (SN-38) for targeted cancer chemotherapy; S. J. Moon, S. V. Govindan, T. M. Cardillo, C. A. D'Souza, H. J. Hansen, D. M. Goldenberg; **J. Med. Chem.** 2008; **51**: 6916-26. <https://doi.org/10.1021/jm800719t>
- Novel Silyl Ether-Based Acid-Cleavable Antibody-MMAE Conjugates with Appropriate Stability and Efficacy; Y. Wang, S. Fan, D. Xiao, F. Xie, W. Li, W. Zhong, X. Zhou; **Cancers (Basel)** 2019; **11**: 957. <https://doi.org/10.3390/cancers11070957>



For more infos on PROTACs, see
our brochure on linker technologies!



[↑ back to content](#)

Index

Product code	Product name	Page	Product code	Product name	Page
RL-2130	(Me)Tz-benzoic acid	20	ADC1790	Alkyne-HMPO-OH	82
RL-2140	(Me)Tz-butanoic acid	19	ADC1800	Alkyne-HMPO-PNP	82
PTC1460	(S,R,S)-AHPC-PEG1-Alkyne	133	RL-2055	Alkyne-myristic acid	68
PTC1590	(S,R,S)-AHPC-PEG1-N ₃	135	PEG6815	Alkyne-O ₂ Oc-OH	107
PTC1470	(S,R,S)-AHPC-PEG2-Alkyne	134	RL-2060	Alkyne-palmitic acid	69
PTC1600	(S,R,S)-AHPC-PEG2-N ₃	136	RL-8630	Alkyne-PEG(4)-CINA	124
PTC1480	(S,R,S)-AHPC-PEG3-Alkyne	134	PEG5440	Alkyne-PEG(4)-mal	111
PTC1610	(S,R,S)-AHPC-PEG3-N ₃	136	PEG5430	Alkyne-PEG(4)-NH ₂	107
PTC1490	(S,R,S)-AHPC-PEG4-Alkyne	134	PEG5410	Alkyne-PEG(4)-NHS	108
PTC1500	(S,R,S)-AHPC-PEG5-Alkyne	134	ADC1350	Alkyne-PEG(4)-Val-Ala-PAB	80
PTC1510	(S,R,S)-AHPC-PEG6-Alkyne	134	ADC1720	Alkyne-PEG(4)-Val-Ala-PAB-Cl	80
PTC1680	(S,R,S)-AHPC-PEG6-butyl-N ₃	136	ADC1360	Alkyne-PEG(4)-Val-Ala-PAB-PNP	81
PTC1640	(S,R,S)-AHPC-PEG6-N ₃	136	RL-3940	Alkyne-PEG(5)-SNAP	121
RL-3460	10-Undecynoyl-OSu	69	ADC1180	Alkyne-PEG(5)-Val-Cit-PAB	81
RL-3170	11-Azido-undecanoyl-OSu	35	ADC1190	Alkyne-PEG(5)-Val-Cit-PAB-PNP	81
RL-3200	11-Azidoundecanoic acid	35	PEG4820	Alkyne-PEG-Si(OMe) ₃ (10 kDa)	112
RL-3220	12-Azido-dodecanoyl-OSu	36	PEG4825	Alkyne-PEG-Si(OMe) ₃ (20 kDa)	113
RL-3210	12-Azidododecanoic acid	35	PEG4810	Alkyne-PEG-Si(OMe) ₃ (3 kDa)	112
RL-3230	14-Azido-myristic acid	36	PEG4815	Alkyne-PEG-Si(OMe) ₃ (5 kDa)	112
RL-3240	16-Azido-palmitic acid	36	RL-3930	Alkyne-SNAP	121
RL-3250	18-Azido-stearic acid	36	RL-3330	Alkyne-SS-COOH	68
RL-4070	2-Azidobenzyl alcohol	37	RL-2065	Alkyne-stearic acid	69
RL-2995	4-(Azidomethyl)benzoic acid	36	PEG6560	Alkynyl-PEG-SQA (10 kDa)	113
RL-2990	4-Azidobenzyl alcohol	37	PEG6565	Alkynyl-PEG-SQA (20 kDa)	113
RL-4280	4-Pentynoyl-Gly-CBT	123	PEG6570	Alkynyl-PEG-SQA (3 kDa)	113
ADC1310	4-Pentynoyl-Val-Ala-PAB	79	PEG6575	Alkynyl-PEG-SQA (5 kDa)	113
ADC1690	4-Pentynoyl-Val-Ala-PAB-Cl	80	GBB1445	alpha-GalNAc-N ₃	127
ADC1320	4-Pentynoyl-Val-Ala-PAB-PNP	80	GBB1375	alpha-GalNAc-TEG-Alkyne	127
ADC1140	4-Pentynoyl-Val-Cit-PAB	80	GBB1370	alpha-GalNAc-TEG-N ₃	127
ADC1150	4-Pentynoyl-Val-Cit-PAB-PNP	80	GBB1440	alpha-Man-N ₃	129
RL-8670	5HP20-alkyne	66	GBB1425	alpha-Man-TEG-Alkyne	129
ADC1290	6-Azidohexanoyl-Val-Ala-PAB	64	GBB1420	alpha-Man-TEG-N ₃	129
ADC1680	6-Azidohexanoyl-Val-Ala-PAB-Cl	64	PEG6895	Aminoxy-PEG(11)-N ₃ *HCl	100
ADC1300	6-Azidohexanoyl-Val-Ala-PAB-PNP	64	RL-2035	ATFB	36
ADC1120	6-Azidohexanoyl-Val-Cit-PAB	64	RL-2045	ATFB-NHS	37
ADC1130	6-Azidohexanoyl-Val-Cit-PAB-PNP	64	RL-3960	Azide-PEG(4)-SNAP	121
RL-4290	6-Azidopentanoyl-Gly-CBT	123	RL-3950	Azide-SNAP	121
RL-3480	8-Azido-octanoyl-OSu	35	ADC1580	Azido-cyclobutane-1,1-dicarboxamide-Ala-PAB	65

Product code	Product name	Page	Product code	Product name	Page
ADC1590	Azido-cyclobutane-1,1-dicarboxamide-Ala-PAB-PNP	65	PEG8000	Biotin-PEG(4)-Picolyl-N ₃	91
ADC1480	Azido-cyclobutane-1,1-dicarboxamide-Cit-PAB	65	PEG8100	Biotin-PEG(4)-SS-Azide	92
ADC1490	Azido-cyclobutane-1,1-dicarboxamide-Cit-PAB-PNP	65	PEG8120	Biotin-PEG(4)-SS-DBCO	10
ADC1330	Azido-PEG(4)-Val-Ala-PAB	65	PEG4330	Biotin-PEG(7)-N ₃	92
ADC1710	Azido-PEG(4)-Val-Ala-PAB-Cl	64	PEG6805	Biotin-PEG3-NH-Pentynoyl	107
ADC1340	Azido-PEG(4)-Val-Ala-PAB-PNP	65	PEG6795	Biotin-PEG3-Ph(4-N ₃)	91
ADC1160	Azido-PEG(4)-Val-Cit-PAB	66	RL-3490	Biotin-Propargylamide	67
ADC1170	Azido-PEG(4)-Val-Cit-PAB-PNP	66	RL-4120	Biotin-SS-N ₃	63
PEG4840	Azido-PEG-Si(OMe) ₃ (10 kDa)	105	PEG2065	Biotin-TEG-ATFBA	91
PEG4845	Azido-PEG-Si(OMe) ₃ (20 kDa)	105	LS-4660	Biotinyl-DAPE-N ₃	125
PEG4830	Azido-PEG-Si(OMe) ₃ (3 kDa)	104	PEG7075	Bis-PEG(11)-DBCO	14
PEG4835	Azido-PEG-Si(OMe) ₃ (5 kDa)	104	BAA4895	Boc-Aeg(N ₃)-OK	47
RL-3320	Azido-Pen-SS-COOH	63	BAA4920	Boc-Aoa-NH-PEG3-N ₃	95
RL-4100	Azido-SS-COOH	62	BNN1320	Boc-Cystamine-Poc	79
RL-4150	Azido-SS-OpNC	63	BAA1805	Boc-D-Aha-OH*CHA	40
GBB1430	beta-Gal-Et-N ₃	127	BAA1825	Boc-D-Aza-OH*CHA	38
GBB1385	beta-Gal-TEG-Alkyne	128	BAA1815	Boc-D-Lys(N ₃)-OH*CHA	51
GBB1380	beta-Gal-TEG-N ₃	127	BAA1835	Boc-D-Orn(N ₃)-OH*CHA	54
GBB1435	beta-Glc-N ₃	128	BAA1855	Boc-D-Phe(4-N ₃)-OH	59
GBB1395	beta-Glc-TEG-Alkyne	128	BAA1377	Boc-D-Pra-OH*DCHA	73
GBB1390	beta-Glc-TEG-N ₃	128	BAA3120	Boc-D-Pro(4-N ₃)-OH (2R,4R)	57
GBB1405	beta-GlcNAc-TEG-Alkyne	129	BAA3110	Boc-D-Pro(4-N ₃)-OH*DCHA (2R,4S)	57
GBB1400	beta-GlcNAc-TEG-N ₃	128	BAA1800	Boc-L-Aha-OH*CHA	40
GBB1455	beta-Lac-EO-N ₃	130	BAA1820	Boc-L-Aza-OH*CHA	38
GBB1415	beta-Lac-TEG-Alkyne	130	BAA2250	Boc-L-Cys(Propargyl)-OH*DCHA	71
GBB1410	beta-Lac-TEG-N ₃	130	BAA1810	Boc-L-Lys(N ₃)-OH*CHA	50
LS-4270	Biotin-DBCO	9	BAA4900	Boc-L-Lys(N ₃)-OK	50
LS-4280	Biotin-MeTz	22	BAA1960	Boc-L-Lys(Poc)-OH	74
LS-4210	Biotin-N ₃	125	BAA1830	Boc-L-Orn(N ₃)-OH*CHA	54
PEG4340	Biotin-PEG(11)-N ₃	92	BAA4660	Boc-L-Phe(4-CH2-N ₃)-OH	59
PEG4350	Biotin-PEG(23)-N ₃	92	BAA4670	Boc-L-Phe(4-Eth)-OH	76
PEG4940	Biotin-PEG(3)-N ₃	91	BAA1850	Boc-L-Phe(4-N ₃)-OH	59
PEG4950	Biotin-PEG(4)-alkyne	106	BAA3980	Boc-L-Phe(4-NH-Poc)-OH	76
PEG8010	Biotin-PEG(4)-Alkyne	106	BAA1434	Boc-L-Pra-OH*DCHA	73
RL-2520	Biotin-PEG(4)-DBCO	10	BAA1905	Boc-L-Pro(4-N ₃)-OH (2S,4S)	57
PEG7980	Biotin-PEG(4)-Dde-Alkyne	106	BAA1930	Boc-L-Pro(4-N ₃)-OH*DCHA (2S,4R)	57
PEG8140	Biotin-PEG(4)-Dde-DBCO	10	BAA2260	Boc-L-Ser(Propargyl)-OH*DCHA	78
PEG7960	Biotin-PEG(4)-Dde-N ₃	91	BAA4650	Boc-L-Tyr(2-azidoethyl)-OH	61
PEG7970	Biotin-PEG(4)-Dde-Picolyl-N ₃	92	BAA2235	Boc-L-Tyr(PEG(3)-N ₃)-OH*DCHA	61
LS-4290	Biotin-PEG(4)-MeTz	22	BAA2265	Boc-L-Tyr(Propargyl)-OH*DCHA	79
PEG7990	Biotin-PEG(4)-N ₃	91	BAA3230	Boc-N-(propargyl)-glycine	72

[↑ back to content](#)

The Click Reaction

Amino Acid Derivatives and Related
Building Blocks for Click Chemistry

Spermines and Amines for Click Chemistry

Click Reagents for Drug Delivery

Click Chemistry Tools for Proteomics

Carbohydrates for Click Chemistry

Proteolysis Targeting Chimeras (PROTACs*)

Index

Product code	Product name	Page	Product code	Product name	Page
PEG4960	Boc-NH-PEG(2)-N ₃	95	PEG6760	DBCO-PEG(24)-TFP	12
PEG7200	Bromoacetamido-PEG(11)-N ₃	101	RL-2480	DBCO-PEG(3)-BisSulfonThiol-Linker	9
PEG6800	Bromoacetamido-PEG(12)-DBCO	11	PEG6765	DBCO-PEG(36)-TFP	13
PEG7210	Bromoacetamido-PEG(23)-N ₃	102	RL-8640	DBCO-PEG(4)-CINA	124
PEG6810	Bromoacetamido-PEG(24)-DBCO	12	RL-2500	DBCO-PEG(4)-mal	11
PEG7190	Bromoacetamido-PEG(3)-N ₃	101	RL-2420	DBCO-PEG(4)-NH ₂ *TFA	10
PEG6790	Bromoacetamido-PEG(4)-DBCO	11	RL-2510	DBCO-PEG(4)-OH	10
RL-2230	Bz-(Me)Tz-NHS	21	PEG6740	DBCO-PEG(4)-TFP	11
RL-2240	Bz-Tz-NHS	21	RL-2450	DBCO-PEG(5)-COOH	11
RL-2250	Bz-Tz-PEG(5)-NHS	22	RL-2460	DBCO-PEG(5)-NHS	11
RL-4180	ClICr [®] base compound	6	RL-4010	DBCO-SNAP	121
RL-4190	ClICr [®] -beta-Ala-NH ₂ *TFA	6	RL-4310	DBCO-Suc-CBT	14, 123
RL-4330	ClICr [®] -OSu	6	RL-4110	DBCO-Suc-SS-COOH	14, 68
RL-4200	ClICr [®] -Suc	6	RL-2421	DBCO-Sulfo-PEG(4)-NH ₂	10
RL-3600	DACN(Ms)*HCl	15	LS-4310	DecarboxyBiotin-Alkyne	67
RL-3610	DACN(Ms,Ns)	15	LS-4300	DecarboxyBiotin-N ₃	125
RL-2735	DACN(Tos)*HCl	15	FAA4055	Fmoc-Abg(N ₃)-OH	47
RL-2710	DACN(Tos,Ns)	15	FAA2095	Fmoc-Abu(3-N ₃)-OH (2R,3R)	44
RL-2720	DACN(Tos,Suc-OH)	15	FAA3540	Fmoc-Abu(3-N ₃)-OH (2R,3S)	44
RL-2730	DACN(Tos2)	16	FAA3200	Fmoc-Abu(3-N ₃)-OH (2S,3R)	44
RL-2737	DACN(Tos2,6-OH)	16	FAA2040	Fmoc-Abu(3-N ₃)-OH (2S,3S)	43
AAA2190	DAPOA*DCHA	33	FAA4060	Fmoc-Aeg(N ₃)-OH	47
RL-4020	DBCO-C6-Alkyne	8	FAA2080	Fmoc-alpha-Prg-D-Ala-OH	70
RL-2430	DBCO-COOH	8	FAA6810	Fmoc-D-Aha-OH	41
ADC1620	DBCO-cyclobutane-1,1-dicarboxamide-Ala-PAB	8	FAA6870	Fmoc-D-Aza-OH	38
ADC1630	DBCO-cyclobutane-1,1-dicarboxamide-Ala-PAB-PNP	9	FAA3650	Fmoc-D-Dbu(N ₃)-OH	42
ADC1520	DBCO-cyclobutane-1,1-dicarboxamide-Cit-PAB	9	FAA1835	Fmoc-D-Lys(N ₃)-OH	51
ADC1530	DBCO-cyclobutane-1,1-dicarboxamide-Cit-PAB-PNP	9	FAA8135	Fmoc-D-Lys(pentynoyl)-OH	75
RL-2490	DBCO-mal	8	FAA9565	Fmoc-D-Lys(Pryoc)-OH	75
RL-2540	DBCO-mPEG (10kDa)	13	FAA6890	Fmoc-D-Orn(N ₃)-OH	55
RL-2550	DBCO-mPEG (20kDa)	13	FAA1910	Fmoc-D-Phe(4-N ₃)-OH	60
RL-2560	DBCO-mPEG (30kDa)	13	FAA1690	Fmoc-D-Pra-OH	73
RL-2530	DBCO-mPEG (5kDa)	13	FAA4720	Fmoc-D-Pro(4-N ₃)-OH (2R,4R)	58
RL-2120	DBCO-NH ₂	8	FAA4630	Fmoc-D-Pro(4-N ₃)-OH (2R,4S)	58
RL-2440	DBCO-NHS	8	FNN1030	Fmoc-DAP-N ₃	84
PEG6830	DBCO-PEG(12)-(5)6-carboxyfluorescein	12	FNN1020	Fmoc-EDA-N ₃	84
PEG6770	DBCO-PEG(12)-MAL	12	FAA6620	Fmoc-L-Aha-OH	41
PEG6750	DBCO-PEG(12)-TFP	12	FAA7400	Fmoc-L-Arg(Propargyl,Pbf)-OH	71
PEG7095	DBCO-PEG(24)-amido-PEG(24)-DSPE	14	FAA1820	Fmoc-L-Aza-OH (solv.)	38
PEG6780	DBCO-PEG(24)-MAL	13	FAA9460	Fmoc-L-Cys(Azidoethyl)-OH	47
PEG7460	DBCO-PEG(24)-OMe	12	FAA3810	Fmoc-L-Cys(Propargyl)-OH	71

Product code	Product name	Page	Product code	Product name	Page
FAA4170	Fmoc-L-Dap(Pentynoyl)-OH	70	HAA1895	H-D-Orn(N ₃)-OH*HCl	54
FAA4230	Fmoc-L-Dap(Poc)-OH	70	HAA1855	H-D-Phe(4-N ₃)-OH	58
FAA2035	Fmoc-L-Dbu(N ₃)-OH	42	HAA1856	H-D-Phe(4-N ₃)-OH*HCl	59
FAA8880	Fmoc-L-Lys(2-N ₃ -Z)-OH	53	HAA6490	H-D-Pra-OH*HCl	72
FAA8890	Fmoc-L-Lys(3-N ₃ -Z)-OH	53	HAA3190	H-D-Pro(4-N ₃)-OH*HCl (2R,4R)	56
FAA8830	Fmoc-L-Lys(4-N ₃ -Z)-OH	53	HAA3140	H-D-Pro(4-N ₃)-OH*HCl (2R,4S)	56
FAA8905	Fmoc-L-Lys(CO-Ethoxypropargyl)-OH	76	HAA9280	H-L-Aha-OH	40
FAA8825	Fmoc-L-Lys(COCF ₂ N ₃)-OH	52	HAA5730	H-L-Aha-OH*HCl	40
FAA8995	Fmoc-L-Lys(Hexynoyl)-OH	75	HAA1880	H-L-Aza-OH*HCl hydrate	37
FAA1793	Fmoc-L-Lys(N ₃)-OH	51	HAA2350	H-L-Cys(Propargyl)-OH*HCl	71
FAA8145	Fmoc-L-Lys(N ₃ -Aca-DIM)-OH	51	HAA3970	H-L-Dbu(N ₃)-OH*HCl	42
FAA7925	Fmoc-L-Lys(N ₃ -AEEA)-OH	51	HAA9380	H-L-Lys(2-N ₃ -Z)-OH	53
FAA8855	Fmoc-L-Lys(N ₃ -Gly)-OH	49	HAA9370	H-L-Lys(3-N ₃ -Z)-OH*HCl	53
FAA4175	Fmoc-L-Lys(pentynoyl)-OH	75	HAA9315	H-L-Lys(4-N ₃ -Z)-OH*HCl	53
FAA8115	Fmoc-L-Lys(Pentynoyl-DIM)-OH	75	HAA9390	H-L-Lys(CO-Ethoxypropargyl)-OH*HCl	74
FAA3150	Fmoc-L-Lys(Pryoc)-OH	75	HAA9295	H-L-Lys(COCF ₂ N ₃)-OH*HCl	50
FAA9420	Fmoc-L-MeDap(N ₃)-OH	38	HAA2080	H-L-Lys(EO-N ₃)-OH*HCl	50
FAA8595	Fmoc-L-MeLys(N ₃)-OH	51	HAA9170	H-L-Lys(MeTz-PhAc)-OH*TFA	19
FAA8680	Fmoc-L-MeOrn(N ₃)-OH	55	HAA9210	H-L-Lys(N ₃)-OH	49
FAA6880	Fmoc-L-Orn(N ₃)-OH	54	HAA1625	H-L-Lys(N ₃)-OH*HCl	49
FAA7740	Fmoc-L-Phe(4-CH ₂ -N ₃)-OH	60	HAA9340	H-L-Lys(N ₃ -Gly)-OH*HCl	49
FAA8530	Fmoc-L-Phe(4-Eth)-OH	77	HAA9235	H-L-Lys(Norbornene-methoxycarbonyl)-OH*HCl	24
FAA1905	Fmoc-L-Phe(4-N ₃)-OH	60	HAA9440	H-L-Lys(Pentynoyl)-OH	74
FAA7720	Fmoc-L-Phe(4-NH-Poc)-OH	77	HAA2085	H-L-Lys(Pentynoyl)-OH*HCl	74
FAA1589	Fmoc-L-Pra-OH	73	HAA2090	H-L-Lys(Poc)-OH*HCl	74
WAA6025	Fmoc-L-Pra-Wang Resin	73	HAA1620	H-L-Orn(N ₃)-OH*HCl	54
FAA3000	Fmoc-L-Pro(4-N ₃)-OH (2S,4R)	57	HAA9500	H-L-Phe(3-BuTz)-OH*TFA	19
FAA2050	Fmoc-L-Pro(4-N ₃)-OH (2S,4S)	57	HAA9510	H-L-Phe(3-iPrTz)-OH*TFA	19
FAA7130	Fmoc-L-Pro(4-NHPoc)-OH (2S,4S)	77	HAA9490	H-L-Phe(3-MeTz)-OH*TFA	19
FAA3820	Fmoc-L-Ser(Propargyl)-OH	78	HAA9480	H-L-Phe(4-Azido-PrTz)-OH*TFA	19, 59
FAA8535	Fmoc-L-Tyr(2-azidoethyl)-OH	61	HAA4090	H-L-Phe(4-CH ₂ -N ₃)*HCl	59
FAA3830	Fmoc-L-Tyr(Propargyl)-OH	79	HAA9220	H-L-Phe(4-Eth)-OH*HCl	76
FAA4950	Fmoc-N-(propargyl)-glycine	72	HAA9470	H-L-Phe(4-MeTz)-OH*TFA	18
RL-3780	Fmoc-N-propargyl-MPBA	70	HAA1850	H-L-Phe(4-N ₃)-OH	58
RL-4380	Fmoc-NH-PEG(3)-N ₃	99	HAA2980	H-L-Phe(4-N ₃)-OH*HCl	58
HAA2870	H-(Gly)3-Lys(N ₃)-OH*HCl	50	HAA4970	H-L-Phe(4-NH-Poc)-OH*HCl	76
HAA3020	H-Abu(3-N ₃)-OH*HCl (2S,3R)	43	HAA7151	H-L-Pra-OH	72
HAA3010	H-Abu(3-N ₃)-OH*HCl (2S,3S)	43	HAA3150	H-L-Pro(4-N ₃)-OH*HCl (2S,4R)	56
HAA1630	H-D-Aha-OH*HCl	40	HAA2125	H-L-Pro(4-N ₃)-OH*HCl (2S,4S)	56
HAA1885	H-D-Aza-OH*HCl hydrate	38	HAA2355	H-L-Ser(Propargyl)-OH*HCl	77
HAA1890	H-D-Lys(N ₃)-OH*HCl	50	HAA9310	H-L-Tyr(2-azidoethyl)-OH	61

[↑ back to content](#)

Product code	Product name	Page	Product code	Product name	Page
HAA9215	H-L-Tyr(2-azidoethyl)-OH*HCl	61	RL-2870	ICG-DBCO	7, 118
HAA3940	H-L-Tyr(3-N ₃)-OH	60	RL-3945	Mal-Alkyne	67
HAA1970	H-L-Tyr(Propargyl)-OH	78	MAA1100	Mal-AMCHC-N-Propargylamide	67
HAA1971	H-L-Tyr(Propargyl)-OH*HCl	78	RL-3935	Mal-PEG(3)-N ₃	90
HAA9445	H-Sar-Sar-NH-Propargyl*HCl	72	PEG7465	Me-PEG(12)-DBCO	7
PEG1081	H ₂ N-PEG(11)-N ₃	94	POX2210	Me-PETox(100)-N ₃	116
PEG4980	H ₂ N-PEG(2)-N ₃ *TosOH	93	POX2200	Me-PETox(50)-N ₃	116
PEG3070	H ₂ N-PEG(23)-N ₃	94	POX1210	Me-PMeOx(100)-N ₃	115
PEG3060	H ₂ N-PEG(3)-N ₃	93	POX1200	Me-PMeOx(50)-N ₃	115
PEG3080	H ₂ N-PEG(35)-N ₃	94	PEG1660	MeO-PEG(12)-N ₃	103
PEG1087	H ₂ N-PEG(6)-N ₃	94	PEG1710	MeO-PEG(24)-N ₃	103
PEG2350	H ₂ N-PEG(7)-N ₃	94	PEG3430	MeO-PEG(36)-N ₃	103
PEG3050	H ₂ N-PEG(8)-N ₃	94	PEG1690	MeO-PEG(4)-N ₃	103
PEG2950	H ₂ N-PEG-alkyne (10 kDa)	107	PEG1705	MeO-PEG(8)-N ₃	103
PEG2970	H ₂ N-PEG-alkyne (20 kDa)	108	PEG2800	MeO-PEG-alkyne (10 kDa)	112
PEG2960	H ₂ N-PEG-alkyne (3 kDa)	107	PEG2810	MeO-PEG-alkyne (2 kDa)	111
PEG2980	H ₂ N-PEG-alkyne (5 kDa)	107	PEG2820	MeO-PEG-alkyne (20 kDa)	112
PEG3000	H ₂ N-PEG-N ₃ (10 kDa)	95	PEG2830	MeO-PEG-alkyne (5 kDa)	112
PEG3020	H ₂ N-PEG-N ₃ (20 kDa)	95	PEG2840	MeO-PEG-alkyne (750 Da)	111
PEG3010	H ₂ N-PEG-N ₃ (3 kDa)	95	PEG2045	MeO-PEG-N ₃ (10 kDa)	104
PEG3030	H ₂ N-PEG-N ₃ (5 kDa)	95	PEG1225	MeO-PEG-N ₃ (2 kDa)	104
RL-3670	Halo-DBCO	9, 120	PEG2050	MeO-PEG-N ₃ (20 kDa)	104
RL-3700	Halo-PEG(2)-Azide	120	PEG2040	MeO-PEG-N ₃ (5 kDa)	104
RL-3710	Halo-PEG(4)-Azide	120	PEG1219	MeO-PEG-N ₃ (750 Da)	103
RL-3640	Halo-PEG(5)-azide	120	RL-2360	MeTz-Bzl-NH ₂ *HCl	20
PEG5330	HO-PEG-N ₃ (10 kDa)	97	RL-2310	MeTz-PEG(4)-COOH	21
PEG5340	HO-PEG-N ₃ (20 kDa)	97	RL-2340	MeTz-PEG(4)-mal	22
PEG5350	HO-PEG-N ₃ (3 kDa)	97	RL-2370	MeTz-PEG(4)-NH ₂ *HCl	21
PEG5360	HO-PEG-N ₃ (5 kDa)	97	RL-2330	MeTz-PEG(4)-NHS	21
POX2250	HOOC-PETox(50)-N ₃	116	RL-3905	MeTz-PEG(4)-STP	20
POX1250	HOOC-PMeOx(50)-N ₃	116	RL-2320	MeTz-PhAc-NHS	20
PEG3090	I-PEG-alkyne (10 kDa)	111	RL-3915	MeTz-PhAc-Sulfo-NHS	21
PEG3100	I-PEG-alkyne (20 kDa)	111	RL-2300	MeTz-PhAcOH	20
PEG3110	I-PEG-alkyne (3 kDa)	110	HAA2230	N ₃ -1,4- <i>cis</i> -CHC-OH	33
PEG3120	I-PEG-alkyne (5 kDa)	111	HAA2235	N ₃ -1,4- <i>trans</i> -CHC-OH	34
PEG3130	I-PEG-N ₃ (10 kDa)	102	HAA6990	N ₃ -Aca-Aca-OH	33
PEG3140	I-PEG-N ₃ (20 kDa)	102	RL-2980	N ₃ -Aca-OSu	35
PEG3150	I-PEG-N ₃ (3 kDa)	102	PEG7950	N ₃ -AEEA-OK	98
PEG3160	I-PEG-N ₃ (5 kDa)	102	PEG5400	N ₃ -AEEA*CHA	98
RL-2880	ICG-alkyne	118	RL-4430	N ₃ -C7H14-COOH	32
RL-2840	ICG-azide	118	HNN1090	N ₃ -Cystamine*HCl	63

Product code	Product name	Page	Product code	Product name	Page
HAA2155	N ₃ -D-Dab(Boc)-OH	41	PEG3770	N ₃ -PEG(24)-OH	97
HAA2135	N ₃ -D-Dap(Boc)-OH	39	PEG7650	N ₃ -PEG(24)-TFP	101
HAA2145	N ₃ -D-Dap(Fmoc)-OH	39	PEG8160	N ₃ -PEG(3)-NH-Boc	96
HAA2175	N ₃ -D-Lys(Boc)-OH	52	RL-4370	N ₃ -PEG(3)-NH-DIG-OH	99
HAA2165	N ₃ -D-Lys(Fmoc)-OH	52	RL-4320	N ₃ -PEG(3)-NH-Suc	93
HAA2210	N ₃ -D-Orn(Boc)-OH*CHA	55	PEG3760	N ₃ -PEG(3)-OH	96
RL-4360	N ₃ -DABu-Suc-OH	33	RL-4420	N ₃ -PEG(3)-SQA	105
RL-4350	N ₃ -DAPr-Suc-OH	33	PEG3780	N ₃ -PEG(36)-OH	97
BNN1370	N ₃ -EDA-Suc-OH	32	PEG7660	N ₃ -PEG(36)-TFP	101
PEG4900	N ₃ -EEEEt-OH	96	RL-8625	N ₃ -PEG(4)-CINA	124
HAA9330	N ₃ -Gly-Aeg(Fmoc)-OH	48	PEG2345	N ₃ -PEG(4)-COOH	99
HAA2860	N ₃ -Gly-Gly-Gly-Gly-Gly-OH	48	PEG5320	N ₃ -PEG(4)-NH ₂	93
HAA2840	N ₃ -Gly-Gly-Gly-OH	48	PEG1400	N ₃ -PEG(4)-NHS	99
HAA2850	N ₃ -Gly-Gly-OH*DCHA	48	PEG5300	N ₃ -PEG(4)-OH	96
HAA3175	N ₃ -HABA*DCHA (2S)	34	PEG6915	N ₃ -PEG(4)-TFP	99
AAA1960	N ₃ -Hx-OH	35	PEG6720	N ₃ -PEG(5)-OH	96
HAA3365	N ₃ -IsoSer*DCHA (2S)	34	PEG4170	N ₃ -PEG(8)-COOH	99
HAA4980	N ₃ -L-Cit-OH*DCHA	46	PEG1405	N ₃ -PEG(8)-NHS	100
HAA2810	N ₃ -L-Cys(Trt)-OH*CHA	47	PEG1088	N ₃ -PEG(8)-OH	90
HAA2150	N ₃ -L-Dab(Boc)-OH	41	PEG2015	N ₃ -PEG(9)-COOH	100
HAA3170	N ₃ -L-Dab(Fmoc)-OH	41	PEG6645	N ₃ -PEG-SQA (10 kDa)	105
HAA2130	N ₃ -L-Dap(Boc)-OH	39	PEG6650	N ₃ -PEG-SQA (20 kDa)	106
HAA2140	N ₃ -L-Dap(Fmoc)-OH	39	PEG6655	N ₃ -PEG-SQA (3 kDa)	105
HAA3350	N ₃ -L-Leu-OH*BHA	48	PEG6660	N ₃ -PEG-SQA (5 kDa)	105
HAA2820	N ₃ -L-Leu-OH*CHA	49	PEG7215	N ₃ -PEG2-NMe3 mesylate	93
HAA2170	N ₃ -L-Lys(Boc)-OH	52	RL-3280	N ₃ -Pen-Dde	63
HAA2160	N ₃ -L-Lys(Fmoc)-OH	52	RL-3290	N ₃ -Pen-Dtpp	63
HAA2880	N ₃ -L-Lys(Mtt)-OH	52	AAA1970	N ₃ -Pen-OH	34
HAA9495	N ₃ -L-Orn(Boc)-OK	55	PGA1135	N ₃ -PGA(100)	114
HAA9485	N ₃ -L-Orn(Boc)-ONa	55	PGA1125	N ₃ -PGA(20)	114
HAA2225	N ₃ -L-Orn(Fmoc)-OH	55	PGA1130	N ₃ -PGA(50)	114
HAA3360	N ₃ -L-Phe-OH*DCHA	60	HAA2245	N ₃ -PhAc-OH	34
HAA9285	N ₃ -L-Val-OH*CHA	62	RL-3650	N ₃ -Phenylpropionic-OH	33
PEG2790	N ₃ -O ₂ Oc-O ₂ Oc-OH	98	HNN1110	N ₃ -TEMPO	66
PEG2780	N ₃ -O ₂ Oc-OH*CHA	98	PEG5000	N ₃ -TFBA-O ₂ Oc	37
PEG5390	N ₃ -O ₂ Oc-OtBu	98	BNN1150	N ₃ -TOTA	93
PEG4180	N ₃ -PEG(12)-COOH	100	PEG5170	N ₃ -TOTA-Suc	98
PEG1395	N ₃ -PEG(12)-NHS	100	HAA2240	N ₃ -trans-MCHC-OH	34
PEG1390	N ₃ -PEG(12)-OH	96	RL-3010	N ₃ Ac-OphOMe	119
PEG6925	N ₃ -PEG(12)-TFP	100	PGA1300	<i>n</i> Bu-PGA(100)[PEG2-N ₃ (10% mod)]	115
PEG4190	N ₃ -PEG(24)-COOH	101	PGA1290	<i>n</i> Bu-PGA(20)[PEG2-N ₃ (10% mod)]	115

The Click Reaction

Amino Acid Derivatives and Related
Building Blocks for Click Chemistry

Spermines and Amines for Click Chemistry

Click Reagents for Drug Delivery

Click Chemistry Tools for Proteomics

Carbohydrates for Click Chemistry

Proteolysis Targeting Chimeras (PROTACs*)

Index

[↑ back to content](#)

Product code	Product name	Page	Product code	Product name	Page
PGA1205	<i>n</i> Bu-PGA(200)[Prg(20)]	115	PEG8170	Propargyl-PEG(5)-COOH	108
PGA1295	<i>n</i> Bu-PGA(50)[PEG2-N ₃ (10% mod)]	115	SNN1170	Spermine(HHHN ₃)*3HCl	84
PEG2900	NHS-PEG(NH-Boc)-alkyne (10 kDa)	109	SNN1230	Spermine(N ₃ BBB)	84
PEG2920	NHS-PEG(NH-Boc)-alkyne (20 kDa)	109	SNN1210	Spermine(N ₃ HHN ₃)*2TsOH	84
PEG2910	NHS-PEG(NH-Boc)-alkyne (3 kDa)	109	RL-2010	TBTA	82
PEG2930	NHS-PEG(NH-Boc)-alkyne (5 kDa)	109	LS-3405	TCEP*HCl	31
PEG2905	NHS-PEG(NH-Fmoc)-alkyne (10 kDa)	110	TCO1060	TCO-NH ₂ *HCl	23
PEG2925	NHS-PEG(NH-Fmoc)-alkyne (20 kDa)	110	TCO1000	TCO-NHS	23
PEG2915	NHS-PEG(NH-Fmoc)-alkyne (3 kDa)	110	TCO1020	TCO-PEG(12)-NHS	23
PEG2935	NHS-PEG(NH-Fmoc)-alkyne (5 kDa)	110	TCO1050	TCO-PEG(3)-mal	23
PEG2850	NHS-PEG-alkyne (10 kDa)	109	TCO1070	TCO-PEG(3)-NH ₂ *HCl	22
PEG2870	NHS-PEG-alkyne (20 kDa)	109	TCO1040	TCO-PEG(4)-COOH	23
PEG2860	NHS-PEG-alkyne (3 kDa)	108	TCO1010	TCO-PEG(4)-NHS	23
PEG2880	NHS-PEG-alkyne (5 kDa)	108	PEG7070	Tetra(-PEG(11)-DBCO)pentaerythritol	14
RL-2090	Norbornene-methyl-NHS	24	RL-2210	THPTA	82
RL-2080	Norbornene-NHS	24	RL-2580	Tz-benzoic acid	20
RL-4030	PFB-mercaptopropionyl-PEG3-N ₃	90	ZAA1285	Z-D-Dbu(N ₃)-OH	43
RL-3410	Photo-Click-Heptanoic acid	68	ZAA5700	Z-L-Aha-OH*DCHA	41
RL-3720	Photo-Click-Palmitic acid	69	ZAA1290	Z-L-Dbu(N ₃)-OH	43
RL-3760	Photo-Click-Stearic acid	69			
RL-8660	Poc-Aoa	68			
PAA1170	Poc-Cystamine*HCl	79			
PAA1050	Poc-O ₂ OC-OH*DCHA	108			
PTC1520	Pomalidomid- PEG1-N ₃	134			
PTC1530	Pomalidomid- PEG2-N ₃	135			
PTC1540	Pomalidomid- PEG3-N ₃	135			
PTC1570	Pomalidomid-C6-PEG1-C3-PEG1-butyl-N ₃	135			
PTC1560	Pomalidomid-C6-PEG3-butyl-N ₃	135			
PTC1580	Pomalidomid-PEG6-butyl-N ₃	135			
PTC1400	Pomalidomide-PEG1-Alkyne	133			
PTC1410	Pomalidomide-PEG2-Alkyne	133			
PTC1420	Pomalidomide-PEG3-Alkyne	133			
PTC1440	Pomalidomide-PEG5-Alkyne	133			
PGA1095	Prg-PGA(100)	114			
PGA1085	Prg-PGA(20)	114			
PEG2755	Propargyl amine	67			
ADC1600	Propargyl-cyclobutane-1,1-dicarboxamide-Ala-PAB	81			
ADC1610	Propargyl-cyclobutane-1,1-dicarboxamide-Ala-PAB-PNP	81			
ADC1500	Propargyl-cyclobutane-1,1-dicarboxamide-Cit-PAB	81			
ADC1510	Propargyl-cyclobutane-1,1-dicarboxamide-Cit-PAB-PNP	82			
PEG1935	Propargyl-NHS	69			

Notes

The Click Reaction

Amino Acid Derivatives and Related Building Blocks for Click Chemistry

Spermines and Amines for Click Chemistry

Click Reagents for Drug Delivery

Click Chemistry Tools for Proteomics

Carbohydrates for Click Chemistry

Proteolysis Targeting Chimeras (PROTACs[®])

Index

[↑ back to content](#)

Code of Conduct

As business activity of Iris Biotech GmbH impacts people's lives and health, it must be operated in ethical and correct manner and act with integrity and responsibility. To ensure high ethical standards and fair business practices, Iris Biotech GmbH applies an integrated policy known as its Code of Conduct.

In 2001 Iris Biotech GmbH was founded just at the beginning of the Biotech movement and the first remarkable breakthrough of biotech pharma products. Although the biotech field is rather young compared to other industries we believe on long-term business, a good partnership between our business partners and Iris Biotech GmbH and a good reputation. It is our duty as well as our responsibility to maintain and to extend this over the next generations – based on the principles of an honourable and prudent tradesman which based upon the concept of honourable entrepreneurship.

This Code of Conduct has been developed following the “Voluntary Guidelines for Manufacturers of Fine Chemical Intermediates and Active Ingredients” issued by AIME (Agrochemical & Intermediates Manufacturers in Europe) and the requirements of some of our business associates.

Iris Biotech GmbH commits to hold this Code of Conduct and to include and apply its principles in the management system and the company policies.

Ethics

Iris Biotech GmbH undertakes business in an ethical manner and acts with integrity. All corruption, extortion and embezzlement are prohibited. We do not pay or accept bribes or participate in other illegal inducements in business or government relationships. We conduct our business in compliance with all applicable anti-trust laws. Employees are encouraged to report concerns or illegal activities in the workplace, without threat of reprisal, intimidation or harassment.

Labour

Iris Biotech GmbH is committed to uphold the human rights of workers and to treat them with dignity and respect. Child labour, workplace harassment, discrimination, and harsh and inhumane treatment are prohibited. Iris Biotech GmbH respects the rights of the employees to associate freely, join or not join labour unions, seek representation and join workers' councils. Employees are paid and their working timetable is established according to applicable wage and labour laws. Employees are able to communicate openly with management regarding working conditions without threat of reprisal, intimidation or harassment.

General Policies

Contracts and Secrecy Agreements are binding and the confidential information received is only used for intended purposes. Clear management and organizational structures exist to provide efficient normal working and to address problems quickly. Know-how is protected and intellectual property is respected.

Health and Safety

Iris Biotech GmbH provides a safe and healthy working environment to the employees and protects them from overexposure to chemical and physical hazards. Products are produced, stored and shipped under the guidelines of the relevant chemical and safety legislation. Risks and emergency scenarios are identified and evaluated, and their possible impact is minimized by implementing emergency plans and written procedures. Safety information regarding hazardous materials is available to educate, train and protect workers from hazards. Preventive equipment and facilities maintenance is performed at suitable periods to reduce potential hazards. Employees are regularly trained in health and safety matters and are informed about product properties and risk classification when it is required.

Environment

Iris Biotech GmbH operates in an environmentally responsible and efficient manner, minimizing adverse impacts on the environment. Waste streams are managed to ensure a safe handling, movement, storage, recycling and reuse, before and after being generated. Systems to prevent and mitigate accidental spills and releases to the environment are in place. All required environmental permits and licenses are obtained and their operational and reporting requirements are complied with.

Production and Quality Management

A quality management system following the Good Distribution Practices (GDP rules) of Active Pharmaceutical Ingredients is established covering all the aspects of the worldwide distribution of products. Regular audits are performed to evaluate the efficiency and fulfilling of the quality system. Process controls to provide reproducible product quality are established. There are preventive maintenance procedures to ensure plant reliability and the lowest risk of failure. Staff is trained periodically about GMP and GDP rules. Procedures are established and installations are designed to avoid cross contamination. Batch and analytical records are kept for inspection and audit purposes for suitable periods according guidelines.

Research and Development

Research and development staff education is appropriate to their functional activity and they are trained to develop, optimize and scale-up the processes. Intellectual property is respected and know-how protected. Development of manufacturing processes reflects the principles of the Green Chemistry according to the American Chemical Society Green Chemistry Institute. Animal testing is not used unless alternatives are not scientifically valid or accepted by regulators. If animal testing is carried out, animals are treated so that pain and stress are minimized.

Terms and Conditions of Sales

All orders placed by a buyer are accepted and all contracts are made subject to the terms which shall prevail and be effective notwithstanding any variations or additions contained in any order or other document submitted by the buyer. No modification of these terms shall be binding upon Iris Biotech GmbH unless made in writing by an authorised representative of Iris Biotech GmbH.

Placing of Orders

Every order made by the buyer shall be deemed an offer by the buyer to purchase products from Iris Biotech GmbH and will not be binding on Iris Biotech GmbH until a duly authorised representative of Iris Biotech GmbH has accepted the offer made by the buyer. Iris Biotech GmbH may accept orders from commercial, educational or government organisations, but not from private individuals and Iris Biotech GmbH reserves the right to insist on a written order and/or references from the buyer before proceeding.

There is no minimum order value. At the time of acceptance of an order Iris Biotech GmbH will either arrange prompt despatch from stock or the manufacture/acquisition of material to satisfy the order. In the event of the latter Iris Biotech GmbH will indicate an estimated delivery date. In addition to all its other rights Iris Biotech GmbH reserves the right to refuse the subsequent cancellation of the order if Iris Biotech GmbH expects to deliver the product on or prior to the estimated delivery date. Time shall not be of the essence in respect of delivery of the products. If Iris Biotech GmbH is unable to deliver any products by reason of any circumstances beyond its reasonable control („Force Majeure“) then the period for delivery shall be extended by the time lost due to such Force Majeure. Details of Force Majeure will be forwarded by Iris Biotech GmbH to the buyer as soon as reasonably practicable.

Prices, Quotations and Payments

Prices are subject to change. For the avoidance of doubt, the price advised by Iris Biotech GmbH at the time of the buyer placing the order shall supersede any previous price indications. The buyer must contact the local office of Iris Biotech GmbH before ordering if further information is required. Unless otherwise agreed by the buyer and Iris Biotech GmbH, the price shall be for delivery ex-works. In the event that the buyer requires delivery of the products otherwise than ex-works the buyer should contact the local office of Iris Biotech GmbH in order to detail its requirements. Iris Biotech GmbH shall, at its discretion, arrange the buyer's delivery requirements including, without limitation, transit insurance, the mode of transit (Iris Biotech GmbH reserves the right to vary the mode of transit if any regulations or other relevant considerations so require) and any special packaging requirements (including cylinders). For the avoidance of doubt all costs of delivery and packaging in accordance with the buyer's requests over and above that of delivery in standard packaging ex-works shall be for the buyer's account unless otherwise agreed by both parties. Incoterms 2020 shall apply. Any tax, duty or charge imposed by governmental authority or otherwise and any other applicable taxes, duties or charges shall be for the buyer's account. Iris Biotech GmbH may, on request and where possible, provide quotations for multiple packs or bulk quantities, and non-listed items. Irrespective of the type of request or means of response all quotations must be accepted by the buyer without condition and in writing before an order will be accepted by Iris Biotech GmbH. Unless agreed in writing on different terms, quotations are valid for 30 days from the date thereof. Payment terms are net 30 days from invoice date unless otherwise agreed in writing. Iris Biotech GmbH reserves the right to request advance payment at its discretion. For overseas transactions the buyer shall pay all the banking charges of Iris Biotech GmbH. The buyer shall not

be entitled to withhold or set-off payment for the products for any reason whatsoever. Government/ Corporate Visa and MasterCard (and other such credit cards) may be accepted on approved accounts for payment of the products. Personal credit cards are not acceptable. Failure to comply with the terms of payment of Iris Biotech GmbH shall constitute default without reminder. In these circumstances Iris Biotech GmbH may (without prejudice to any other of its rights under these terms) charge interest to accrue on a daily basis at the rate of 2% per month from the date upon which payment falls due to the actual date of payment (such interest shall be paid monthly). If the buyer shall fail to fulfil the payment terms in respect of any invoice of Iris Biotech GmbH Iris Biotech GmbH may demand payment of all outstanding balances from the buyer whether due or not and/or cancel all outstanding orders and/or decline to make further deliveries or provision of services except upon receipt of cash or satisfactory securities. Until payment by the buyer in full of the price and any other monies due to Iris Biotech GmbH in respect of all other products or services supplied or agreed to be supplied by Iris Biotech GmbH to the buyer (including but without limitation any costs of delivery) the property in the products shall remain vested in Iris Biotech GmbH.

Shipping, Packaging and Returns

The buyer shall inspect goods immediately on receipt and inform Iris Biotech GmbH of any shortage or damage within five days. Quality problems must be notified within ten days of receipt. Goods must not be returned without prior written authorisation of Iris Biotech GmbH. Iris Biotech GmbH shall at its sole discretion replace the defective products (or parts thereof) free of charge or refund the price (or proportionate price) to buyer. Opened or damaged containers cannot be returned by the buyer without the written prior agreement of Iris Biotech GmbH. In the case of agreed damaged containers which cannot be so returned, the buyer assumes responsibility for the safe disposal of such containers in accordance with all applicable laws.

Product Quality, Specifications and Technical Information

Products are analysed in the Quality Control laboratories of Iris Biotech GmbH's production partners by methods and procedures which Iris Biotech GmbH considers appropriate. In the event of any dispute concerning reported discrepancies arising from the buyer's analytical results, determined by the buyer's own analytical procedures, Iris Biotech GmbH reserves the right to rely on the results of own analytical methods of Iris Biotech GmbH. Certificates of Analysis or Certificates of Conformity are available at the discretion of Iris Biotech GmbH for bulk orders but not normally for prepack orders. Iris Biotech GmbH reserves the right to make a charge for such certification. Specifications may change and reasonable variation from any value listed should not form the basis of a dispute. Any supply by Iris Biotech GmbH of bespoke or custom product for a buyer shall be to a specification agreed by both parties in writing. Technical information, provided orally, in writing, or by electronic means by or on behalf of Iris Biotech GmbH, including any descriptions, references, illustrations or diagrams in any catalogue or brochure, is provided for guidance purposes only and is subject to change.

Safety

All chemicals should be handled only by competent, suitably trained persons, familiar with laboratory procedures and potential chemical hazards. The burden of safe use of the products of Iris Biotech GmbH vests in the buyer. The buyer assumes all responsibility for warning his employees, and any persons who might reasonably be expected to come into contact with the products, of all risks to person and property in any way connected with the products and for instructing them in their safe handling and use. The buyer also assumes the responsibility for the safe disposal of all products in accordance with all applicable laws.

Uses, Warranties and Liabilities

All products of Iris Biotech GmbH are intended for laboratory research purposes and unless otherwise stated on product labels, in the catalogue and product information sheet of Iris Biotech GmbH or in other literature furnished to the buyer, are not to be used for any other purposes, including but not limited to use as or as components in drugs for human or animal use, medical devices, cosmetics, food additives, household chemicals, agricultural or horticultural products or pesticides. Iris Biotech GmbH offers no warranty regarding the fitness of any product for a particular purpose and shall not be responsible for any loss or damage whatsoever arising there from. No warranty or representation is given by Iris Biotech GmbH that the products do not infringe any letters patent, trademarks, registered designs or other industrial rights. The buyer further warrants to Iris Biotech GmbH that any use of the products in the United States of America shall not result in the products becoming adulterated or misbranded within the meaning of the Federal Food, Drug and Cosmetic Act (or such equivalent legislation in force in the buyer's jurisdiction) and shall not be materials which may not, under sections 404, 505 or 512 of the Act, be introduced into interstate commerce. The buyer acknowledges that, since the products of Iris Biotech GmbH are intended for research purposes, they may not be on the Toxic Substances Control Act 1976 („TSCA“) inventory. The buyer warrants that it shall ensure that the products are approved for use under the TSCA (or such other equivalent legislation in force in the buyer's jurisdiction), if applicable. The buyer shall be responsible for complying with any legislation or regulations governing the use of the products and their importation into the country of destination (for the avoidance of doubt to include, without limitation, the TSCA and all its amendments, all EINECS, ELINCS and NONS regulations). If any licence or consent of any government or other authority shall be required for the acquisition, carriage or use of the products by the buyer the buyer shall obtain the same at its own expense and if necessary produce evidence of the same to Iris Biotech GmbH on demand. Failure to do so shall not entitle the buyer to withhold or delay payment. Any additional expenses or charges incurred by Iris Biotech GmbH resulting from such failure shall be for the buyer's account. Save for death or personal injury caused by negligence of Iris Biotech GmbH, sole obligation of Iris Biotech GmbH and buyer's exclusive remedy with respect to the products proved to the satisfaction of Iris Biotech GmbH to be defective or products incorrectly supplied shall be to accept the return of said products to Iris Biotech GmbH for refund of the actual purchase price paid by the buyer (or proportionate part thereof), or replacement of the defective product (or part thereof) with alternative product. Iris Biotech GmbH shall have no liability to the buyer under or arising directly or indirectly out of or otherwise in connection with the supply of products by Iris Biotech GmbH to the buyer and/or their re-sale or use by the buyer or for any product, process or services of the buyer which in any way comprises the product in contract tort (including negligence or breach of statutory duty) or otherwise for pure economic loss, loss of profit, business, reputation, depletion of brand, contracts, revenues or anticipated savings or for any special indirect or consequential damage or loss of any nature except as may otherwise be expressly provided for in these terms. All implied warranties, terms and representations in respect of the products (whether implied by statute or otherwise) are excluded to the fullest extent permitted by law. The buyer shall indemnify Iris Biotech GmbH for and against any and all losses, damages and expenses, including legal fees and other costs of defending any action, that Iris Biotech GmbH may sustain or incur as a result of any act or omission by the buyer, its officers, agents or employees, its successors or assignees, its customers or all other third parties, whether direct or indirect, in connection with the use of any product. For the avoidance of doubt and in the event that Iris Biotech GmbH supplies bespoke or custom product to the buyer's design or specification, this indemnity shall extend to include any claim by a third party that the manufacture of the product for the buyer or the use of the product by the buyer infringes the intellectual property rights of any third party.

General

Iris Biotech GmbH shall be entitled to assign or sub-contract all or any of its rights and obligations hereunder. The buyer shall not be entitled to assign, transfer, sub-contract or otherwise delegate any of its rights or obligations hereunder. Any delay or forbearance by Iris Biotech GmbH in exercising any right or remedy under these terms shall not constitute a waiver of such right or remedy. If any provision of these terms is held by any competent authority to be invalid or unenforceable in whole or in part the validity of the other provisions of these terms and the remainder of the provision in question shall not be affected. These terms shall be governed by German Law and the German Courts shall have exclusive jurisdiction for the hearing of any dispute between the parties save in relation to enforcement where the jurisdiction of the German Courts shall be non-exclusive.



Get in Contact



Iris Biotech GmbH
Adalbert-Zoellner-Str. 1
95615 Marktredwitz
Germany

☎ +49 (0) 9231 97121-0
☎ +49 (0) 9231 97121-99
✉ info@iris-biotech.de
🌐 www.iris-biotech.de

Distribution Partners

The list contains the current distributors of Iris Biotech in different regions of the world. The latest list of distribution partners and contact details is available at: www.iris-biotech.de/distribution-partner

China:

Chengdu Yoo Technology Co., Ltd.

Japan:

BizCom Japan, Inc.

Shigematsu & Co., Ltd

Cosmo Bio Co., Ltd.

USA & Canada:

Peptide Solutions LLC

India, Bangladesh, Oman, Sri Lanka, United Arab Emirates:

Sumit Biosciences Pvt Ltd.

Empowering Peptide Innovation