



Iris
Biotech

LINKEROLOGY®



Version: IB6_4

Empowering Peptide Innovation

With this guiding theme in mind, Iris Biotech's mission is to support researchers by supplying

- innovative technologies,
- rare compounds,
- as well as a broad portfolio on standard consumables,

available in flexible quantities from small scale to bulk quantities. To fulfill our dedication "Empowering Peptide Innovation", we are attending various conferences, symposia, and exhibitions each year. This allows us to remain in direct contact with scientists all over the world, both from academia and industry, to exchange knowledge, and to gather new ideas to tackle your current challenges.

Guided by our dedication to provide

- competent service,
- as well as novel substances and
- latest technologies,

Iris Biotech is your trusted partner for the world of peptides, while having strong expertise in associated disciplines. Thus, our portfolio comprises reagents and tools for the synthesis and modification of peptides, e.g., amino acids, resins and solvents but also for related technologies such as drug delivery, linkerology® and life sciences.

Owed to the growing demand for tailor-made compounds, our portfolio is fine-tuned by our custom synthesis service at Iris Biotech Laboratories. Our skilled scientists offer profound expertise in

- *de novo* route development,
- upscaling towards larger scale production,
- as well as synthesis optimization for increased efficiency.

Examples are the synthesis of rare chiral building blocks, unnatural amino acid derivatives, sophisticated orthogonal protecting groups, heterocycles, building blocks for nucleotides, PEGs and PEG-analogs as well as specific linkers for controlled drug delivery and release.



Amino Acids



Building Blocks



Life Sciences



Drug Delivery



Reagents



Resins



Linkerology®



Click Chemistry

Portfolio Overview

Peptide Synthesis and Modification

(Protected) Amino Acids

Standards such as Fmoc-D/L-AAA and Boc-D/L-AAA, Smoc amino acids for peptide synthesis in water, variety of protecting groups (e.g., Pbf, Trt, ^tBu, Bzl, Acn, Mob, SIT, Phacn, Allocam, Mmt), unusual amino acids, fluorinated derivatives, substituted prolines, arginine analogs

Building Blocks

Amino alcohols, amino aldehydes, diamines and hydrazines, (pseudoproline) dipeptides, polyamines and spermines, fatty acid derivatives, peptide nucleic acids (PNAs)

Reagents

Coupling reagents, solvents and scavengers, protecting groups

Resins

Preloaded resins (e.g., based on Trityl, TCP, TentaGel, Methoxybenzhydryl, Merrifield, PAM, Rink, Wang), scavenger resins, hydrazone resins, poly(acrylamide) resins, Cyclover

Linkerology® and Drug Delivery

Linkers for Solid Phase Peptide Synthesis

Cleavable Linkers

Val-Ala-based, Val-Cit-based, disulfide-based, Dde-helping hands, pH-sensitive linkers

Photo-Activatable Linkers

Functionalized Linkers

Clickable linkers, trifunctional linkers, linkers with maleimide function, cross-linkers, selective N-term acylation and biotinylation, 5HP2O

PROTACs

Ligands, linkers & modules

Fullerenes, Poly(2-oxazolines), Dextrans & Plant-Derived Cholesterol

Superparamagnetic Iron Oxide Nanoparticles

Poly-Amino Acids

Poly-Arg, Poly-Glu, Poly-Lys, Poly-Orn, Poly-Sar

PEGylation

Branched PEGylating reagents, (amino-)PEG-acids, PEG-amines & hydrazides & guanidines, reagents for Click-conjugation, Biotin-PEG-reagents, PEG-thiols, PEG-maleimides, other PEGylating reagents

Life Sciences

Biotinylation Reagents

Carbohydrates

Galactose, Glucose, Mannose, Xylose and others

Drug Metabolites

Peptides

Substrates & Inhibitors

E.g., protein kinase inhibitors, substrates for fusion (Halo/Snap/Clip)-tagged proteins

Natural Products

Dyes and Fluorescent Labels

E.g., ICG, AMC, DAPI

Maillard & Amadori Reaction Products

Large portfolio of derivatives useful as standards for food, pharma and cosmetics industry


Vitamins

Custom Synthesis

Your project requires a compound not listed in our portfolio?
Get in contact and inquire about our custom synthesis capabilities.

Our experienced scientists are excited to accept your synthetic challenge!

In such cases, your request undergoes the following stages:



Step-by-Step Analysis

- Customer's demands

Process Evaluation

- Detailed literature review
- Synthetic possibilities

Strategy Development

- Protocol development
- Method development and validation
- Customized synthesis

Quality Consistency

- Identity confirmation
- Purity verification

Our Service Promise

All our services are based on high standards, transparency & documentation, trust, honesty & confidentiality, as well as the required know-how.

High Standards

- Values: sustainability & responsibility
- State-of-the-art equipment & latest technologies
- High quality standards
- Qualified suppliers & regular audits

Transparency & Documentation

- Talk to our specialists – customer care
- Certificates of analysis & origin
- Impurity profiling
- Safety data sheets
- Analytical and process reports

Trust, Honesty & Confidentiality

- Intergenerational business valuing partnerships
- Meeting the customer's expectations
- Integrity towards our customers

Our Know-How

- One-step reactions & complex multi-step synthesis
- Scalability from mg to kg quantities
- Route scouting





Table of Contents

1. The Concept of Antibody-Drug Conjugation (ADC)	1
1.1. Technical and Market Background	1
1.2. Linker Design, Connectivity, Degradability, and Drug-Antibody Ratio (DAR)	5

2. Permanent Linkers	17
2.1. PEG-Based Spacer Molecules	17
2.2. Hydrophobic Spacer Molecules	35
2.3. Permanent Linkers with Maleimide Function	42
2.4. 5HP2O as Maleimide Alternative	49
2.5. Photoactivatable Linkers	51

3. Cleavable Linkers	56
3.1. Valine-Alanine-Based Enzymatically Cleavable Linkers	57
3.2. Valine-Citrulline-Based Enzymatically Cleavable Linkers	65
3.3. β -Glucuronide Enzymatically Cleavable Linkers	71
3.4. PH-Responsive Linkers	72
3.5. Disulfide-Based (Self-Immolative) Linkers	74
3.6. Dde-Based Linkers	84

4. Trifunctional Linkers	101
---------------------------------	------------


5. Cross-Linkers for other Bio Applications	103
5.1. Substrates for Fusion (Halo/Snap/Clip)-Tagged Proteins	103
5.2. Specific His Tag Acylation	109
5.3. Bifunctional Protein Cross-Linkage	112
5.4. Triazolecarbaldehydes for Selective N-Terminal Protein Modification	114
5.5. Proteolysis Targeting Chimeras (PROTACs®)	116
5.6. Site-Selective π -Clamp Mediated Cysteine Arylation	131
6. Preparing Carriers for Conjugation	133
6.1. Antibodies, Antibody Formats and Proteins by (Cell-free) Recombinant Methodologies	134
6.2. Aptamers and other Oligonucleotides	135
6.3. Carbon Compounds	137
6.4. Metals	142
6.5. Metal Oxides	144
6.6. Polymeric Surfaces by Plasma Treatment	145
6.7. Silicates	149
Index	150
Code of Conduct	159
Terms and Conditions of Sales	161

1. The Concept of Antibody-Drug Conjugation (ADC)

1.1. Technical and Market Background

Conjugating highly potent small molecules to vastly target specific biomolecules, like antibodies, has become a modern and sophisticated approach, particularly in the field of cancer therapy. The list of ADCs in clinics continues to grow, bolstered by the success of two pioneers in this field:

Adcetris® (Seattle Genetics) has been approved in 2011 for the treatment of Hodgkin's lymphoma and systemic anaplastic large cell lymphoma (ALCL) and reached \$476.9 million sales per year in 2018. This drug is composed of a monoclonal antibody targeting CD30 conjugated to four molecules of monomethyl auristatin E via a self-immolative linkage (Fig. 1). Reduction of interchain disulfide bonds provides reactive cysteine residues, which are then conjugated with maleimide payload linker systems, yielding the final drug compound.




Kadcyla®, another pioneer in this field, has been approved in 2013 for the treatment of HER-2 positive metastatic breast cancer and reached \$981 million sales per year in 2018. In this case, payloads are conjugated to surface accessible lysines resulting in a heterogeneous modification of the core antibody.

Reference:


→ *Antibody-drug conjugates in tumor therapy*; B. Sammet, C. Steinkuhler, N. Sewald; *Pharm Pat Anal* 2012; 1: 65-73. <https://doi.org/10.4155/ppa.12.4>

ADCs – Mode of Action

The typical mode of action of ADCs is shown in [Fig. 2](#). An ADC circulates in plasma until it reaches the target cell. The antibody portion of an ADC then binds to a cell-surface antigen that is ideally specific to a cancer cell. Upon binding, the ADC-antigen protein complex becomes internalized into the cancer cell. When the complex is degraded, it releases the cytotoxin which then binds to its target to cause cancer cell apoptosis. The linker between antibody and payload is typically either permanent or cleavable by hydrolases, such as the protease cathepsin B, by glucuronidases or through reductive conditions and the presence of glutathione.



This concept is a sophisticated approach combining the high specificity of antibodies with the high potency of (small) drug molecules. The disadvantages of antibodies, like low potency, as well as the drawbacks of small drug molecules, like low specificity accompanied by high toxicity through many side effects, are compensated by the advantages of the other counterpart. A smart synergistic combination of both elements significantly enlarges the narrow therapeutic window of a small drug molecule between minimum (efficacious) and maximum (toxic) dosage ([Fig. 3](#)). ADC drugs expand the therapeutic window, as they can increase efficacy and decrease toxicity in comparison to traditional chemotherapeutic cancer treatments. Targeted delivery to cancer cells increases the amount of dosed drug reaching the tumor, thus lowering the minimum effective dose (MED). The MTD (maximum, highest tolerable dose without serious side effects) is increased, as less drug compound reaches healthy, non-target tissues. Historically, defining MTD was the primary objective of phase 1 oncology trials. More recently, especially for new targeted drugs (including ADCs), emphasis has been placed on determining the recommended phase 2 dose (RP2D), which better captures chronic toxicities emergent after multiple treatment cycles (e.g., edema, effusion, pneumonitis, ocular toxicities) and certain grade 2 side effects (e.g., diarrhea, mucositis, cytopenia, neuropathy, severe fatigue) that may become intolerable over time. However, interpreting small molecule and ADC doses in clinical trials requires specific interpretation and accurate conversion of the doses to a common unit and expansion of the therapeutic window needs to be discussed individually. An appreciation that ADCs do not significantly enhance the MTDs of their payloads may provide insight into several existing observations in this field, like ADCs that feature a common drug linker often encounter similar MTDs because of payload-associated platform toxicities, independent of the target antigen. This highlights that most off-target adverse events are antibody independent.




References:

- L. Anthony, (2019). ADC Landscape Review 2019 [PowerPoint slides]. Retrieved from <http://worldadc-usa.com>
- Design and Synthesis of Tesirine, a Clinical Antibody-Drug Conjugate Pyrrolobenzodiazepine Dimer Payload; A. C. Tiberghien, J. N. Levy, L. A. Masterson, N. V. Patel, L. R. Adams, S. Corbett, D. G. Williams, J. A. Hartley, P. W. Howard; *ACS Med Chem Lett* 2016; **7**: 983-987. <https://doi.org/10.1021/acsmchemlett.6b00062>
- Recent advances of antibody drug conjugates for clinical applications; P. Zhao, Y. Zhang, W. Li, C. Jeanty, G. Xiang, Y. Dong; *Acta Pharmaceutica Sinica B* 2020; **10(9)**: 1589-1600. <https://doi.org/10.1016/j.apsb.2020.04.012>
- The therapeutic window of antibody drug conjugates: A dogma in need of revision; R. Colombo, J. R. Rich; *Cancer Cell* 2022; **40**. <https://doi.org/10.1016/j.ccell.2022.09.016>
- Therapeutic index improvement of antibody-drug conjugates; H.-P. Gerber, S. Gangwar, A. Betts; *Mabs* 2023; **15(1)**: 2230618. <https://doi.org/10.1016/j.ccell.2022.09.016>

While initially only small molecules or short peptides have been used as payloads, the panel of conjugates has opened to chelators for radioactive nuclides and larger biomolecules, such as toxic enzymes. Additional variations have been introduced on the antibody side by utilizing antibody fragment combinations or diabodies (Fig. 4).

Points of conjugation are typically the thiol groups of cysteines, the amino functions of lysines or the N-terminus of a monoclonal antibody. Due to the inherent heterogeneity of conjugation to the multiple amines or cysteines found in mAbs, significant research efforts are directed toward the production of discrete, homogeneous ADC products via site-specific conjugation. This may involve genetic engineering of the mAb to introduce discrete, available cysteines or non-natural amino acids with an orthogonally reactive functional handle such as an aldehyde, ketone, azido, or alkynyl tag. These site-specific approaches increase the homogeneity of ADCs and enable novel bioorthogonal chemistries which utilize, reactive moieties rather than thiols or amines. This broad diversity of applicable linkers can then be utilized leading to improved design in future generations of ADCs.



Tab. 1: Global Approved ADCs (Jan. 2025).

Approval Year	(Trade) Name	Payload	Linker
2011	Adcetris	MMAE	enzyme cleavable (vc)
2013	Kadcyla	DM1	non-cleavable thioether
2017	Besponsa	Calicheamicin	acid cleavable disulfide
2017	Mylotarg	Calicheamicin	acid cleavable disulfide
2018	Lumoxiti	Pseudomonas exotoxin A fragment (PE38)	fused peptide linkage
2019	Enhertu	Deruxtecan	enzyme cleavable (ggfg)
2019	Padcev	MMAE	enzyme cleavable (vc)
2019	Polivy	MMAE	enzyme cleavable (vc)
2020	Trodelvy	SN-38	acid cleavable <i>p</i> -aminobenzyl carbonate
2020	Blenrep	MMAF	non-cleavable (vvi)
2020	Akalux	IR700	non-cleavable
2021	Tivdak	MMAE	enzyme cleavable (vc)
2021	Zynlonta	PBD dimer	enzyme cleavable (va)
2021	Aidixi	MMAE	enzyme cleavable (vc)
2022	ELAHERE	DM4	non-cleavable
expected for 2025	Datopotamab deruxtecan	Deruxtecan	enzyme cleavable (ggfg)
expected for 2025	Patritumab deruxtecan	Deruxtecan	enzyme cleavable (ggfg)
expected for 2025	Telisotuzumab vedotin	MMAE	enzyme cleavable (vc)

[↑ back to content](#)

References:

- Site-specific antibody drug conjugates for cancer therapy; S. Panowski, S. Bhakta, H. Raab, P. Polakis, J. R. Junutula; **MABs** 2014; **6**: 34-45. <https://doi.org/10.4161/mabs.27022>
- Advances in Precision Oncology: Targeted Thorium-227 Conjugates As a New Modality in Targeted Alpha Therapy; U. B. Hagemann, K. Wickstroem, S. Hammer, R. M. Bjerke, S. Zitzmann-Kolbe, O. B. Ryan, J. Karlsson, A. Scholz, H. Hennekes, D. Mumberg, A. S. Cuthbertson; **Cancer Biother Radiopharm** 2020; **35(7)**: 497-510. <https://doi.org/10.1089/cbr.2020.3568>

Background Information

IC₅₀ Inhibitory Concentration Concentration causing 50% of maximal inhibition of the desired activity.	EC₅₀ Effective Concentration Concentration causing 50% of maximal response of the desired effect.	ED₅₀ Effective Dose Dose causing the desired effect in 50% of individuals.
GI₅₀ Growth Inhibition Concentration causing 50% inhibition of cell proliferation/cell growth.	TC₅₀ Toxic Concentration Concentration causing a defined toxic effect in 50% of individuals.	TD₅₀ Toxic Dose Dose causing a defined toxic effect in 50% of individuals.
CC₅₀ Cytotoxic Concentration Concentration killing 50% of cells.	LC₅₀ Lethal Concentration Concentration killing 50% of individuals.	LD₅₀ Lethal Dose Dose killing 50% of individuals.

1.2. Linker Design, Connectivity, Degradability, and Drug-Antibody Ratio (DAR)

Antibody-drug conjugates (ADCs), which combine the specificity, favorable pharmacokinetics, and bio-distribution of a monoclonal antibody (mAb) with the cytotoxic potency of a drug are promising new therapeutics for cancer. Along with the development of mAbs and cytotoxic drugs, the design of the linker is essential, as it impacts the efficacy and tolerability of ADCs. The linker needs to provide sufficient stability during systemic circulation while providing rapid and efficient release of the cytotoxic drug in its active state inside the tumor cells.




Fig. 5: Thiols from cysteines, amines from lysines, azido functions from non-canonical amino acids and specific sequences accessible on the surface of antibodies can be addressed by different chemical or enzymatical conjugation methodologies.




Fig. 6: Drug-antibody ratio (DAR) is an important parameter of an ADC. Low DAR could reduce the antitumor efficacy, while high DAR may affect antibody structure, stability, and antigen binding etc. therefore causing loss of activity. DAR values are also important for the therapeutic index of ADCs. In most ADC drug candidates, the DAR values were maintained at about 2-4. Hence, controlling DAR during ADC preparation is a key procedure. Figure provided by Glycotope.

[↑ back to content](#)

The type of linkage between payload and biomolecule can basically either be permanent or cleavable under certain well-defined circumstances (Fig. 5). As payloads typically are highly cytotoxic, it would be fatal if they were released from their carrier during circulation in plasma. Hence, the linker part should be stable to conditions such as pH, redox potential, presence of proteases in plasma, and all other parameters of plasma. However, after internalization it is favorable that the linker is fragmenting in order to release the drug molecule, ideally in a traceless manner. Conjugations with the antibody can rather easily be achieved using active esters forming amide bonds with lysines, which are usually accessible in a high number on the surface. The resulting conjugate, hence, is rather heterogenic with different numbers of payloads attached at different positions. A more and well-defined drug-antibody ratio (DAR) can be achieved by utilizing the disulfide bridges between heavy and light chains of the antibody.


After reductive cleavage of the disulfide bonds, conjugation chemistry can be performed by different kinds of reactions like conventional maleimides or disulfide bond formation. Heterogeneity can be observed if heavy and light antibody chains do not recombine in the original manner.

A highly accurate and specific DAR with well-defined connectivity can be achieved, if unnatural amino acids, e.g., *p*-azidophenylalanine, can be introduced recombinantly. Click chemistry or other Diels-Aldertype reactions can be used to introduce linkers and payloads. In a similar manner, certain peptide fragments can be added, which serve as substrates for ligases in order to conjugate to appropriate linker-payload conjugates.


References:

- *Linker Technologies for Antibody-Drug Conjugates*; B. Nolting; *Antibody-Drug Conjugates* L. Ducry 2013; **1045**: 71-100. https://doi.org/10.1007/978-1-62703-541-5_5
- *In Vivo Applications of Bioorthogonal Reactions: Chemistry and Targeting Mechanisms*; M. M. A. Mitry, F. Greco, H. M. I. Osborn; *Chemistry* 2023; n/a: e202203942. <https://doi.org/10.1002/chem.202203942>


Cleavage Mechanisms



An ADC travels through plasma until it reaches the target cell. After internalization, the complex degrades and releases the payload even with a stable linker. However, release can be accelerated through implementation of moieties which fragmentize under certain conditions. One of the most commonly used spacers is the bifunctional *p*-aminobenzyl alcohol group, which is linked to the peptide through the amino group forming an amide bond, while amine containing cytotoxic drugs are attached through carbamate functionalities to the benzylic hydroxyl group of the linker. The resulting prodrugs are activated upon protease mediated hydrolysis and cleavage of the amide bond of citrulline to the *p*-aminobenzyl fragment, leading to a 1,6-elimination reaction releasing the unmodified drug, carbon dioxide, and remnants of the linker group (Fig. 7, Fig. 8).



In an extension of the peptide-based linker strategies to provide high ADC stability, β -glucuronic acid-based linkers were developed. Facile release of the active drug is realized through cleavage of the β -glucuronide glycosidic bond by the lysosomal enzyme β -glucuronidase. This enzyme is abundantly present in lysosomes and overexpressed in some tumor types, while its activity outside cells is low. The linker is hydrophilic, stable against circulation, and provides ADCs that are highly active both *in vitro* and *in vivo*.








Besides hydrolases, the presence of oxidoreductases in the lysosome is being utilized for the design of cleavable linkers. Cytochrome P450 oxidoreductase (CPR), nitroquinone oxidoreductase 1 (NQO1), and cellular reductants such as glutathione (GSH) transform reducible fragments like quinone or disulfide to self-immolative intermediates.

Trimethyl Lock

The sterical demand of three closely positioned methyl groups (Fig. 9) favors the cleavage of a carbonyl bond by lacton formation. The acidity of the phenol is sufficient to accelerate lactonization at neutral pH and any residue carrying a hydroxyl or amino function will be unlocked, i.e. tracelessly released. The hydroxy group of phenol can be protected and released by a variety of methodologies. This reaction usually requires no elevated temperature. Hence, it will work nicely at physiological conditions.

			Product details
RL-2960	Acetyl-Trimethyl-Lock		
3-(2-Acetoxy-4,6-dimethylphenyl)-3-methylbutyric acid			
CAS-No.	134098-68-3		
Formula	C ₁₅ H ₂₀ O ₄		
Mol. weight	264,14 g/mol		

		Product details
RL-2950	Fourmethyl-Lock	
3-(2,4-dimethyl-3,6-dioxocyclohexa-1,4-dienyl)-3-methylbutanoic acid		
CAS-No.	133544-77-1	 
Formula	C ₁₃ H ₁₆ O ₄	
Mol. weight	236,26 g/mol	
RL-2940	Fivemethyl-Lock	
3-methyl-3-(2,4,5-trimethyl-3,6-dioxocyclohexa-1,4-dienyl)butanoic acid		
CAS-No.	40662-29-1	 
Formula	C ₁₄ H ₁₈ O ₄	
Mol. weight	250,29 g/mol	
RL-2970	Photo-Trimethyl-Lock	
3-(2-Nitroveratryl-4,6-dimethylphenyl)-3-methylbutyric acid		
CAS-No.	2095134-25-9	 
Formula	C ₂₂ H ₂₇ NO ₇	
Mol. weight	417,45 g/mol	
FAA7190	Fmoc-Spr(oNB)-OH	
N-alpha-(9-Fluorenylmethoxycarbonyl)-beta,beta-dimethyl-(2,4-dimethyl-6-(2-nitrobenzyloxy)phenyl)alanine (rac.)		
CAS-No.	1032400-98-8	 
Formula	C ₃₅ H ₃₄ N ₂ O ₇	
Mol. weight	594,66 g/mol	
FAA7200	Fmoc-Spr(oNv)-OH	
N-alpha-(9-Fluorenylmethoxycarbonyl)-beta,beta-dimethyl-(2-methyl-6-(2-nitroveratryl)phenyl)alanine (rac.)		
CAS-No.	1228829-20-6	 
Formula	C ₃₆ H ₃₆ N ₂ O ₉	
Mol. weight	640,68 g/mol	

References:

- *Trimethyl lock: A trigger for molecular release in chemistry, biology, and pharmacology*; M. N. Levine, R. T. Raines; **Chem. Sci.** 2012; **3**: 2412-2420. <https://doi.org/10.1039/C2SC20536J>
- *Photo-triggered fluorescent labelling of recombinant proteins in live cells*; D. Jung, K. Sato, K. Min, A. Shigenaga, J. Jung, A. Otaka, Y. Kwon; **Chem Commun** 2015; **51**: 9670-3. <https://doi.org/10.1039/C2SC20536J>
- *Detection of DT-diaphorase Enzyme with a ParaCEST MRI Contrast Agent*; I. Daryaei, K. M. Jones, M. D. Pagel; **Chemistry** 2017; **23**: 6514-6517. <https://doi.org/10.1002/chem.201700721>

[↑ back to content](#)

- Syntheses and kinetic studies of cyclisation-based self-immolative spacers; S. Huvelle, A. Alouane, T. Le Saux, L. Jullien, F. Schmidt; **Org Biomol Chem** 2017; **15**: 3435-3443. <https://doi.org/10.1039/c7ob00121e>
- Invention of stimulus-responsive peptide-bond-cleaving residue (Spr) and its application to chemical biology tools; A. Shigenaga, J. Yamamoto, T. Kohiki, T. Inokuma, A. Otake; **J Pept Sci** 2017; **23**: 505-513. <https://doi.org/10.1002/psc.2961>
- Trimethyl Lock: A Multifunctional Molecular Tool for Drug Delivery, Cellular Imaging, and Stimuli-Responsive Materials; O. A. Okoh, P. Klahn; **ChemBioChem** 2018; **19**: 1668-1694. <https://doi.org/10.1002/cbic.201800269>

Disulfide Linkers

Disulfide linkers (Fig. 10) are likely first degraded in the lysosome to generate a cysteine-disulfide catabolite followed by disulfide reduction in the cytosol by cellular reductants such as GSH. The kinetics of reduction can be tailored by neighboring one to four methyl groups next to both sulfurs.




Fig. 10: Mechanism of disulfide bond cleavage in lysosomal compartments.


References:

- Modulating Therapeutic Activity and Toxicity of Pyrrolobenzodiazepine Antibody-Drug Conjugates with Self-Immolative Disulfide Linkers; T. H. Pillow, M. Schutten, S. F. Yu, R. Ohri, J. Sadowsky, K. A. Poon, W. Solis, F. Zhong, G. Del Rosario, M. A. T. Go, J. Lau, S. Yee, J. He, L. Liu, C. Ng, K. Xu, D. D. Leipold, A. V. Kamath, D. Zhang, L. Masterson, S. J. Gregson, P. W. Howard, F. Fang, J. Chen, J. Gunzner-Toste, K. K. Kozak, S. Spencer, P. Polakis, A. G. Polson, J. A. Flygare, J. R. Junutula; **Mol. Cancer Ther.** 2017; **16**: 871-878. <https://doi.org/10.1158/1535-7163.MCT-16-0641>
- Mechanisms of drug release in nanotherapeutic delivery systems; P. T. Wong, S. K. Choi; **Chem Rev** 2015; **115**: 3388-432. <https://doi.org/10.1021/cr5004634>
- Expanded Utility of the beta-Glucuronide Linker: ADCs That Deliver Phenolic Cytotoxic Agents; S. C. Jeffrey, J. De Brabander, J. Miyamoto, P. D. Senter; **ACS Med Chem Lett** 2010; **1**: 277-80. <https://doi.org/10.1021/ml100039h>

Multiple Payloads with one self-immolative Linker

p-Hydroxy- and *p*-amino-benzyl fragments will release payloads by a 1,6-elimination cascade, resulting in chinoide intermediates. Under physiological conditions, they readily add water to reform the aromatic ring structure. In case appropriate carbamate substitutions are also placed on position 2 and 2', a fragmentation will occur in a similar manner as by a 1,4-elimination and release any molecules at these positions (Fig. 11).

One of the major challenges related to anticancer chemotherapy is resistance against anticancer drugs. A strategy to revert the resistance of tumor cells is the combined use of different anticancer drugs.



References:

- A novel connector linkage applicable in prodrug design; P. L. Carl, P. K. Chakravarty, J. A. Katzenellenbogen; **J. Med. Chem.** 1981; **24**: 479-80. <https://doi.org/10.1021/jm00137a001>
- The azaquinone-methide elimination: comparison study of 1,6- and 1,4-eliminations under physiological conditions; R. Erez, D. Shabat; **Org Biomol Chem** 2008; **6**: 2669-72. <https://doi.org/10.1039/b808198k>
- Dendritic chain Dendritic chain reaction: responsive release of hydrogen peroxide upon generation and enzymatic oxidation of methanol; M. Avital- Shmilovici, D. Shabat; **Bioorg Med Chem** 2010; **18**: 3643-7. <https://doi.org/10.1016/j.bmc.2010.02.038>
- ABC transporters as multidrug resistance mechanisms and the development of chemosensitizers for their reversal; C. H. Choi; **Cancer Cell Int** 2005; **5**: 30. <https://doi.org/10.1186/1475-2867-5-30>

It has been reported that payload release can be supported by introducing a *N,N'*-dimethylethane-1,2-diamine bridge between carbamate and payload. After release of carbon dioxide, it will cyclize and form 1,3-dimethylimidazolidin-2-one and liberate the payload from the linker construction (Fig. 12).





Fig. 12: Mechanism of assisted traceless release by *N,N'*-dimethylethane-1,2-diamine.

Besides the benzyl system, other moieties have been used for fragmentation reactions. In Fig. 13 different methods are summarized, which have been studied and published. PG is the protecting group and LG the leaving group belonging to the payload to be released. X needs to be a strong electron-donating group, such as O, X or NH, in order to initiate the elimination cascade. While the 1,6-elimination of a benzyl system tends to be the most common system, *ortho*-benzyl undergoing a 1,4-elimination can alternatively be used, as well as styrene fragments (1,8-elimination). However, neither naphthalene rings nor biphenyl structures (1,10-elimination) work, even with a strongly donating amino group.



References:

- Self-immolative spacers: kinetic aspects, structure-property relationships, and applications; A. Alouane, R. Labruere, T. Le Saux, F. Schmidt, L. Jullien; **Angew Chem Int Ed** 2015; **54**: 7492-509.
<https://doi.org/10.1002/anie.201500088>
- Cleavable linkers in chemical biology; G. Leriche, L. Chisholm, A. Wagner; **Bioorg Med Chem** 2012; **20**: 571-82. <https://doi.org/10.1016/j.bmc.2011.07.048>

Dioxoborolane Cross-Linker:






Fig. 14: Tetramethyldioxoborolane cross-linker fragmentize in the presence of H_2O_2 .

Tetramethyldioxoborolane linkers are used as cross-linkers or masking protecting groups of alcohols or amines. They release their payloads under mild oxidative conditions, as it is in the presence of H_2O_2 . Oxidation initially creates a free phenolate which triggers an initial 1,6-elimination. The second step in this cascade is a 1,4-elimination, followed by a fragmentation under release of the first conjugate. A second 1,4-elimination follows releasing the second conjugated compound. Such moieties are for example being used for the preparation of masked H_2O_2 probes releasing their active fluorophore upon trigger detection. Another application is the generation of liposomes, which disrupt upon H_2O_2 detection and then release the active cargo. Peptides bearing two lysines can be cyclized and masked inactive *via* such a borolane cross-linker. Activity of the bioactive peptide will be restored, if H_2O_2 is present in the tissue environment.

		Product details
RL-4140	TetraMe-Dioxoborolane-(OpNC)2	
(5-methyl-2-((4-(4,4,5,5-tetramethyl-1,3,2-dioxaborolan-2-yl)benzyl)oxy)-1,3-phenylene)bis(methylene)bis(4-nitrophenyl) bis(carbonate)		
CAS-No.	1355342-68-5	
Formula	C ₃₆ H ₃₅ BN ₂ O ₁₃	
Mol. weight	714,49 g/mol	
		

		Product details
RL-4130	TetraMe-Dioxaborolane-OpNC	
4-Nitrophenyl 4-(4,4,5,5-tetramethyl-1,3,2-dioxaborolan-2-yl)benzyl carbonate		
CAS-No.	1254765-89-3	
Formula	C ₂₀ H ₂₂ BN ₂ O ₇	
Mol. weight	399,21 g/mol	
		

References:

- Influence of Linkers within Stimuli-Responsive Prodrugs on Cancer Therapy: A Case of Five Doxorubicin Dimer-Based Nanoparticles; Q. Wang, C. Wang, S. Li, Y. Xiong, H. Wang, Z. Li, J. Wan, X. Yang, Z. Li; **Chem. Mater.** 2022; **34(5)**: 2085-2097. <https://doi.org/10.1021/acs.chemmater.1c03346>
- Reactive Oxygen Species-Responsive Liposomes via Boronate-Caged Phosphatidylethanolamine; J. Lou, M. D. Best; **Bioconjugate Chem.** 2020; **31(9)**: 2220-2230. <https://doi.org/10.1021/acs.bioconjchem.0c00397>
- Magnetic Resonance Imaging of PSMA-Positive Prostate Cancer by a Targeted and Activatable Gd(III) MR Contrast Agent; H. Li, D. Luo, C. Yuan, X. Wang, J. Wang, J. P. Basilion, T. J. Meade; **J. Am. Chem. Soc.** 2021; **143(41)**: 17097-17108. <https://doi.org/10.1021/jacs.1c07377>
- Novel N-Methylated Cyclodepsipeptide Prodrugs for Targeted Cancer Therapy; C. Wu, Z. Cheng, D. Lu, K. Liu, Y. Cheng, P. Wang, Y. Zhou, M. Li, X. Shao, H. Li, W. Su, L. Fang; **J. Med. Chem.** 2021; **64(2)**: 991-1000. <https://doi.org/10.1021/acs.jmedchem.0c01387>
- Structure-Based Identification of Potent Lysine-Specific Demethylase 1 Inhibitor Peptides and Temporary Cyclization to Enhance Proteolytic Stability and Cell Growth-Inhibitory Activity; H. Kitagawa, M. Kikuchi, S. Sato, H. Watanabe, N. Umezawa, M. Kato, Y. Hisamatsu, T. Umehara, T. Higuchi; **J. Med. Chem.** 2021; **64(7)**: 3707-3719. <https://doi.org/10.1021/acs.jmedchem.0c01371>
- A mitochondrial-targetable dual functional near-infrared fluorescent probe to monitor pH and H₂O₂ in living cells and mice; X. Bi, Y. Wang, D. Wang, L. Liu, W. Zhu, J. Zhang, X. Zha; **RSC Adv.** 2020; **10**: 26874-26879. <https://doi.org/10.1039/D0RA03905E>
- GSH Activated Biotin-tagged Near-Infrared Probe for Efficient Cancer Imaging; R. Guo, F. Huang, B. Zhang, Y. Yan, J. Che, Y. Jin, Y. Zhuang, R. Dong, Y. Li, B. Tan, R. Song, Y. Hu, X. Dong, X. Li, N. Lin; **Theranostics** 2019; **9(12)**: 3515-3525. <https://doi.org/10.7150/thno.32742> <https://www.thno.org/v09p3515.htm>
- Photoactivatable Organic Semiconducting Pro-nanoenzymes; J. Li, J. Huang, Y. Lyu, J. Huang, Y. Jiang, C. Xie, K. Pu; **J. Am. Chem. Soc.** 2019; **141(9)**: 4073-4079. <https://doi.org/10.1021/jacs.8b13507>
- Versatile Histochemical Approach to Detection of Hydrogen Peroxide in Cells and Tissues Based on Puromycin Staining; C. Yik-Sham Chung, G. A. Timblin, K. Saijo, C. J. Chang; **J. Am. Chem. Soc.** 2018; **140(19)**: 6109-6121. <https://doi.org/10.1021/jacs.8b02279>
- Facile Fabrication of 10-Hydroxycamptothecin-Backboned Amphiphilic Polyprodrug with Precisely Tailored Drug Loading Content for Controlled Release; X. Zhang, M. Zhang, M. Wang, H. Peng, Q. Hua, L. Ma, B. Wang, H. Wei; **Bioconjugate Chem.** 2018; **29(7)**: 2239-2247. <https://doi.org/10.1021/acs.bioconjchem.8b00238>
- Self immolative dioxetane based chemiluminescent probe for H₂O₂ detection; O. Seven, F. Sozmen, I. S. Turan; **Sens. Actuators B** 2017; **239**: 1318-1324. <https://doi.org/10.1016/j.snb.2016.09.120>
- Reactive Oxygen Species-Responsive Protein Modification and Its Intracellular Delivery for Targeted Cancer Therapy; M. Wang, S. Sun, C. I. Neufeld, B. Perez-Ramirez, Q. Xu; **Angew. Chem. Int. Ed.** 2014; **53(49)**: 13444-13448. <https://doi.org/10.1002/anie.201407234>
- pH and hydrogen peroxide dual responsive supramolecular prodrug system for controlled release of bioactive molecules; Y. Wang, H. Wang, Y. Chen, X. Liu, Q. Jin, J. Ji; **Colloids and Surfaces B: Biointerfaces** 2014; **121**: 189-195. <https://doi.org/10.1016/j.colsurfb.2014.06.024>
- Investigation of self-immolative linkers in the design of hydrogen peroxide activated metalloprotein inhibitors; J. L. Major Jourden, K. B. Daniel, S. M. Cohen; **Chem. Commun.** 2011; **47**: 7968-7970. <https://doi.org/10.1039/C1CC12526E>

[↑ back to content](#)




2. Permanent Linkers







2.1. PEG-Based Spacer Molecules

This class of linkers is considered non-cleavable, meaning linker cleavage and payload release do not depend on the differential properties between plasma and cytoplasmic compartments. Instead, the release of the cytotoxic drug is postulated to occur after internalization of the ADC *via* antigen-mediated endocytosis and delivery to lysosomal compartments, where the antibody is degraded to the level of amino acids through intracellular proteolytic degradation. This process releases a drug derivative, formed by the cytotoxic drug, the linker, and the amino acid residue to which the linker was covalently attached.

The following section displays examples of hetero-bifunctional PEG-based spacer molecules. As payloads are quite often rather hydrophobic, PEG fragments help to solubilize the linker-payload conjugate, which is essential to perform successful conjugation to the antibody. It further helps to increase the solubility in physiological media and to improve the pharmacokinetic properties of the whole ADC construct.

Two of the latest approved ADCs, Trodelvy and Zynlonta, were developed with PEG spacers as part of their linkers to improve solubility and stability *in vivo*.

			Product details
PEG4120	NHS-PEG(2)-NHS		
3,6-Dioxaoctandioic acid bissuccinimidyl ester			
CAS-No.	65869-63-8		
Formula	C ₁₆ H ₂₀ N ₂ O ₁₀		
Mol. weight	400,34 g/mol		
PEG1365	H₂N-PEG(2)-CO-OtBu		
3-(2-(2-Aminoethoxy)ethoxy)propanoic acid t-butyl ester			
CAS-No.	756525-95-8		
Formula	C ₁₁ H ₂₃ NO ₄		
Mol. weight	233,3 g/mol		
PEG6835	Boc-NH-PEG(3)-NH₂*HCl		
<i>tert</i> -butyl (2-(2-(2-(2-aminoethoxy)ethoxy)ethoxy)ethyl)carbamate			
CAS-No.	101187-40-0 net		
Formula	C ₁₃ H ₂₈ N ₂ O ₅ *HCl		
Mol. weight	292,38*36,46 g/mol		

		Product details	
PEG4980	H₂N-PEG(2)-N₃*TosOH 2-[2-(2-Azidoethoxy)ethoxy]ethanaminium tosylat CAS-No. 166388-57-4 Formula C ₇ H ₁₄ N ₄ O ₂ *C ₇ H ₈ O ₃ S Mol. weight 174,20*172,20 g/mol		
PEG4960	Boc-NH-PEG(2)-N₃ 1-(t-Butyloxycarbonyl-amino)-3,6-dioxa-8-octaneazide CAS-No. 950683-55-3 Formula C ₁₁ H ₂₂ N ₄ O ₄ Mol. weight 274,32 g/mol		
PEG1555	mal-PEG(2)-COOH 3-(2-(2-(3-Maleimidopropanamido)ethoxy)ethoxy)propanoic acid CAS-No. 756525-98-1 Formula C ₁₄ H ₂₀ N ₂ O ₇ Mol. weight 328,32 g/mol		
PEG1560	mal-PEG(2)-NHS 3-(2-(2-(3-Maleimidopropanamido)ethoxy)ethoxy)propanoic acid succinimidyl ester CAS-No. 955094-26-5 Formula C ₁₈ H ₂₃ N ₃ O ₉ Mol. weight 425,39 g/mol		
PEG2030	TUDA 3,6,9-Trioxaundecandioic acid CAS-No. 13887-98-4 Formula C ₈ H ₁₄ O ₇ Mol. weight 222,19 g/mol		
BNN1016	Boc-DOOA 1-(t-Butyloxycarbonyl-amino)-3,6-dioxa-8-octaneamine, liq. CAS-No. 153086-78-3 Formula C ₁₁ H ₂₄ N ₂ O ₄ Mol. weight 248,32 g/mol		


FNN1007 Fmoc-DOOA*HCl

1-(9-Fluorenylmethyloxycarbonyl-amino)-3,6-dioxa-8-octaneamine hydrochloride

CAS-No. 868599-73-9

Formula $C_{21}H_{26}N_2O_4 \cdot HCl$

Mol. weight 370,45*36,45 g/mol




RL-4390 Fmoc-NH-PEG(3)-NH₂*HCl

(9H-fluoren-9-yl)methyl (2-(2-(2-aminoethoxy)ethoxy)ethoxy)ethyl carbamate hydrochloride

CAS-No. 906079-91-2

Formula $C_{23}H_{30}N_2O_5 \cdot HCl$

Mol. weight 414,50*36,45 g/mol




ZNN1120 Z-TOTA

1-Benzyloxycarbonyl-4,7,10-trioxa-13-tridecaneamine

CAS-No. 220156-99-0

Formula $C_{18}H_{30}N_2O_5$

Mol. weight 354,45 g/mol




PEG5180 Fmoc-DOOA-DIG-OH

2-(2-(2-((9-Fluorenylmethyloxycarbonyl)amino)ethoxy)ethoxy)ethylamino-diglycolic acid

CAS-No. 669073-64-7

Formula $C_{25}H_{30}N_2O_8$

Mol. weight 486,51 g/mol




RL-4400 Fmoc-NH-PEG(3)-NH-Suc-OH

1-(9H-fluoren-9-yl)-3,17-dioxo-2,7,10,13-tetraoxa-4,16-diazaicosan-20-oic acid

CAS-No. 1653992-32-5

Formula $C_{27}H_{34}N_2O_8$

Mol. weight 514,57 g/mol




RL-4410 Fmoc-NH-PEG(3)-DIG-OH







1-(9H-fluoren-9-yl)-3,17-dioxo-2,7,10,13,19-pentaoxa-4,16-diazahenicosan-21-oic acid

CAS-No. 489427-26-1

Formula $C_{27}H_{34}N_2O_9$

Mol. weight 530,57 g/mol



			Product details
RL-4060	Biotin-DOOA		
Biotinyl-1-amino-3,6-dioxo-8-octanamine			
CAS-No.	138529-46-1		
Formula	C ₁₆ H ₃₀ N ₄ O ₄ S		
Mol. weight	374,50 g/mol		
FAA5730	Fmoc-TTD-DIG-OH		
[N1-(9-Fluorenylmethoxycarbonyl)-1,13-diamino-4,7,10-trioxatridecan-diglycolic acid			
CAS-No.	75-09-2		
Formula	C ₂₉ H ₃₈ N ₂ O ₉		
Mol. weight	558,62 g/mol		
PEG2420	H-O₂Cc-OH		
[2-(2-aminoethoxy)ethoxy]acetic acid			
CAS-No.	134978-97-5		
Formula	C ₆ H ₁₃ NO ₄		
Mol. weight	163,17 g/mol		
PEG7940	H-O₂Cc-OH*HCl		
8-amino-3,6-dioxaoctanoic acid hydrochloride			
CAS-No.	134979-01-4		
Formula	C ₆ H ₁₃ NO ₄ *HCl		
Mol. weight	163,17*36,45 g/mol		
PEG2430	H-O₂Cc-OtBu*HCl		
[2-(2-aminoethoxy)ethoxy]acetic acid <i>tert</i> -butyl ester*HCl			
CAS-No.	2098500-69-5		
Formula	C ₁₀ H ₂₁ NO ₄ *HCl		
Mol. weight	219,28*36,45 g/mol		
AAA1905	Aloc-O₂Cc-OH*DCHA		
8-(Allyloxycarbonyl-amino)-3,6-dioxaoctanoic acid dicyclohexylamine			
CAS-No.	560088-74-6		
Formula	C ₁₀ H ₁₇ NO ₆ *C ₁₂ H ₂₃ N		
Mol. weight	247,11*181,32 g/mol		


PEG8080 Boc-O₂Oc-OH

(2-(2-(*t*-Butyloxycarbonylamino)ethoxy)ethoxy)acetic acid

CAS-No. 108466-89-3

Formula C₁₁H₂₁NO₆

Mol. weight 263,29 g/mol




BAA1466 Boc-O₂Oc-OH*DCHA

8-(*t*-Butyloxycarbonyl-amino)-3,6-dioxaoctanoic acid dicyclohexylammonium salt

CAS-No. 560088-79-1

Formula C₁₁H₂₁NO₆*C₁₂H₂₃N

Mol. weight 263,29*181,32 g/mol




DAA1016 Dde-O₂Oc-OH

8-[(4,4-Dimethyl-2,6-dioxocyclohex-1-ylidene)ethyl-amino]-3,6-dioxaoctanoic acid, {2-[2-(Dde-amino)ethoxy]ethoxy}acetic acid

CAS-No. 1263045-93-7

Formula C₁₆H₂₅NO₆

Mol. weight 327,37 g/mol




FAA1435 Fmoc-O₂Oc-OH

8-(9-Fluorenylmethyloxycarbonyl-amino)-3,6-dioxaoctanoic acid

CAS-No. 166108-71-0

Formula C₂₁H₂₃NO₆

Mol. weight 385,42 g/mol




PEG4870 Mal-O₂Oc-OH

{2-[2-(2,5-Dioxo-2,5-dihydro-1H-pyrrol-1-yl)ethoxy]ethoxy}acetic acid

CAS-No. 173323-23-4

Formula C₁₀H₁₃NO₆

Mol. weight 243,21 g/mol




PEG4650 Mtt-O₂Oc-OH*DEA







N-(4-Methyltrityl)-8-amino-3,6-dioxaoctanoic acid diethylamine







CAS-No. 2098500-66-2







Formula C₂₆H₂₉NO₄*C₄H₁₁N







Mol. weight 419,51*73,14 g/mol















			Product details	
PEG2780	N₃-O₂Oc-OH*CHA			
[2-(2-azidoethoxy)ethoxy]acetic acid cyclohexylamine salt				
CAS-No.	2098500-94-6			
Formula	C ₆ H ₁₁ N ₃ O ₄ *C ₆ H ₁₃ N			
Mol. weight	189,17*99,17 g/mol			
PEG7950	N₃-AEEA-OK			
Potassium 8-azido-3,6-dioxaoctanoate				
CAS-No.	882518-90-3			
Formula	C ₆ H ₁₀ KN ₃ O ₄			
Mol. weight	39,10*188,16 g/mol			
PEG5390	N₃-O₂Oc-OtBu			
8-Azido-3,6-dioxaoctanoic acid t-butyl ester				
CAS-No.	251564-45-1			
Formula	C ₁₀ H ₁₉ N ₃ O ₄			
Mol. weight	245,28 g/mol			
PEG5000	N₃-TFBA-O₂Oc			
{2-[2-(4-Azido-2,3,5,6-tetrafluorobenzoyl-amino)ethoxy]ethoxy}acetic acid				
CAS-No.	1993119-45-1			
Formula	C ₁₃ H ₁₂ F ₄ N ₄ O ₅			
Mol. weight	380,25 g/mol			
RL-4370	N₃-PEG(3)-NH-DIG-OH			
Diglycolic acid PEG3 azide				
CAS-No.	239081-53-9			
Formula	C ₁₂ H ₂₂ N ₄ O ₇			
Mol. weight	334,33 g/mol			
RL-4380	Fmoc-NH-PEG(3)-N₃			
(9H-fluoren-9-yl)methyl (2-(2-(2-(2-azidoethoxy)ethoxy)ethoxy)ethyl)carbamate				
CAS-No.	1172605-58-1			
Formula	C ₂₃ H ₂₈ N ₄ O ₅			
Mol. weight	440,50 g/mol			







			Product details	
PAA1050	Poc-O₂Oc-OH*DCHA			
8-(Popargyloxycarbonyl-amino)-3,6-dioxaoctanoic acid dicyclohexylamine				
CAS-No.	2988660-57-5			
Formula	C ₁₀ H ₁₅ NO ₆ *C ₁₂ H ₂₃ N			
Mol. weight	245,23*181,32 g/mol			
ZAA1186	Z-O₂Oc-OH*DCHA			
8-(Benzoyloxycarbonyl-amino)-3,6-dioxaoctanoic acid dicyclohexylamine				
CAS-No.	560088-84-8			
Formula	C ₁₄ H ₁₉ NO ₆ *C ₁₂ H ₂₃ N			
Mol. weight	297,31*181,32 g/mol			
PEG1810	Fmoc-AEEP			
3-(2-(2-(9-Fluorenylmethyloxycarbonyl)aminoethoxy)ethoxy)propanoic acid				
CAS-No.	872679-70-4			
Formula	C ₂₂ H ₂₅ NO ₆			
Mol. weight	399,44 g/mol			
PEG4970	Fmoc-Ebes			
N-[8-(9-Fluorenylmethyloxycarbonyl)amino-3,6-dioxaoctyl]succinamic acid				
CAS-No.	613245-91-3			
Formula	C ₂₅ H ₃₀ N ₂ O ₇			
Mol. weight	470,51 g/mol			
FAA1568	Fmoc-TTDS-OH			
[N1-(9-Fluorenylmethoxycarbonyl)-1,13-diamino-4,7,10-trioxatridecan-succinamic acid				
CAS-No.	172089-14-4			
Formula	C ₂₉ H ₃₈ N ₂ O ₈			
Mol. weight	542,63 g/mol			
PEG4900	N₃-EEEt-OH			
2-[2-(2-Azidoethoxy)ethoxy]ethanol				
CAS-No.	86520-52-7			
Formula	C ₆ H ₁₃ N ₃ O ₃			
Mol. weight	175,19 g/mol			







		Product details
PEG5400	N₃-AEEEA*CHA	
11-Azido-3,6,9-trioxaundecanoic acid cyclohexylamine		
CAS-No.	172531-37-2	
Formula	C ₈ H ₁₅ N ₃ O ₅ *C ₆ H ₁₃ N	
Mol. weight	233,22*99,17 g/mol	
		
PEG5370	Fmoc-AEEE	
2-(2-(2-(9-Fluorenylmethyloxycarbonyl)aminoethoxy)ethoxy)ethanol		
CAS-No.	560088-66-6	
Formula	C ₂₁ H ₂₅ NO ₅	
Mol. weight	371,43 g/mol	
		
PEG7010	Trt-S-EEE	
S-Trityl-2-(2-(2-mercaptoethoxy)ethoxy)ethanol		
CAS-No.	728033-15-6	
Formula	C ₂₅ H ₂₈ O ₃ S	
Mol. weight	408,55 g/mol	
		
BNN1028	Boc-TOTA	
1-(t-Butyloxycarbonyl-amino)-4,7,10-trioxa-13-tridecanamine, liq.		
CAS-No.	194920-62-2	
Formula	C ₁₅ H ₃₂ N ₂ O ₅	
Mol. weight	320,43 g/mol	
		
FNN1011	Fmoc-TOTA*HCl	
1-(9-Fluorenylmethyloxycarbonyl-amino)-4,7,10-trioxa-13-tridecanamine hydrochloride		
CAS-No.	868599-75-1	
Formula	C ₂₅ H ₃₄ N ₂ O ₅ *HCl	
Mol. weight	442,56*36,45 g/mol	
		
BNN1150	N₃-TOTA	
1-Azido-4,7,10-trioxa-13-tridecanamine		
CAS-No.	1162336-72-2	
Formula	C ₁₀ H ₂₂ N ₄ O ₃	
Mol. weight	246,31 g/mol	
		







			Product details	
PEG5170	N₃-TOTA-Suc			
1-Azido-4,7,10-trioxa-13-tridecaneamine succinamic acid				
CAS-No.	1993176-74-1			
Formula	C ₁₄ H ₂₆ N ₄ O ₆			
Mol. weight	346,38 g/mol			
PEG7860	Boc2-AEEEE			
2-(2-(2-(2-(Di-(<i>t</i> -butylmethyloxycarbonyl))aminoethoxy)ethoxy)ethoxy)ethanol				
CAS-No.	2389064-37-1			
Formula	C ₁₈ H ₃₅ NO ₈			
Mol. weight	393,47 g/mol			
PEG5385	Boc,Z-AEEEE			
2-(2-(2-(2-(Benzyloxycarbonyl- <i>tert</i> -Butylmethyloxycarbonyl)aminoethoxy)ethoxy)ethoxy)ethanol				
CAS-No.	2389064-46-2			
Formula	C ₂₁ H ₃₃ NO ₈			
Mol. weight	427,49 g/mol			
PEG5380	Fmoc-AEEEE			
2-(2-(2-(2-(9-Fluorenylmethyloxycarbonyl)aminoethoxy)ethoxy)ethoxy)ethanol				
CAS-No.	868594-41-6			
Formula	C ₂₃ H ₂₉ NO ₆			
Mol. weight	415,48 g/mol			
PEG6730	Trt-S-EEEE			
S-Trityl-2-(2-(2-(2-(2-mercaptoethoxy)ethoxy)ethoxy)ethoxy)ethanol				
CAS-No.	125607-10-5			
Formula	C ₂₇ H ₃₂ O ₄ S			
Mol. weight	452,61 g/mol			
PEG4130	NHS-PEG(3)-NHS			
3,6,9-Trioxaundecandioic acid bissuccinimidyl ester				
CAS-No.	1314378-16-9			
Formula	C ₁₈ H ₂₄ N ₂ O ₁₁			
Mol. weight	444,39 g/mol			







		Product details
PEG1485	mal-PEG(3)-mal	
Bis-(1,13-(3-maleimidopropionyl)amido)-4,7,10-trioxatridecane		
CAS-No.	756525-89-0	 
Formula	C ₂₆ H ₃₄ N ₄ O ₉	
Mol. weight	522,55 g/mol	
PEG7870	Boc-NH-PEG(3)-NH₂	
1-(t-Butyloxycarbonyl)amino-3,6,9-trioxa-undecan-11-amine		
CAS-No.	101187-40-0	 
Formula	C ₁₃ H ₂₈ N ₂ O ₅	
Mol. weight	292,37 g/mol	
PEG7190	Bromoacetamido-PEG(3)-N₃	
Bromoacetamido-tri(ethylene glycol)-azide		
CAS-No.	940005-81-2	 
Formula	C ₁₀ H ₁₉ BrN ₄ O ₄	
Mol. weight	339,19 g/mol	
PEG4370	Fmoc-NH-PEG(3)-COOH	
12-(9-Fluorenylmethyloxycarbonylamino)-4,7,10-trioxa-dodecanoic acid		
CAS-No.	867062-95-1	 
Formula	C ₂₄ H ₂₉ NO ₇	
Mol. weight	443,49 g/mol	
TCO1050	TCO-PEG(3)-mal	
<i>trans</i> -Cyclooctene-PEG(3)-maleimide		
CAS-No.	1809356-72-6	 
Formula	C ₂₆ H ₄₁ N ₃ O ₈	
Mol. weight	523,62 g/mol	
PEG3760	N₃-PEG(3)-OH	
alpha-Azido-omega-hydroxy tetra(ethylene glycol)		
CAS-No.	86770-67-4	 
Formula	C ₈ H ₁₇ N ₃ O ₄	
Mol. weight	219,24 g/mol	






			Product details
PEG3060	H₂N-PEG(3)-N₃		
1-Amino-11-azido-3,6,9-trioxaundecane			
CAS-No.	134179-38-7		
Formula	C ₈ H ₁₈ N ₄ O ₃		
Mol. weight	218,25 g/mol		
PEG5320	N₃-PEG(4)-NH₂		
14-Azido-3,6,9,12-tetraoxatetradecan-1-amine			
CAS-No.	951671-92-4		
Formula	C ₁₀ H ₂₂ N ₄ O ₄		
Mol. weight	262,31 g/mol		
PEG1221	H-O₂Oc-O₂Oc-OH		
17-Amino-10-oxo-3,6,12,15-tetraoxa-9-azaheptadecan-1-oic acid			
CAS-No.	1143516-05-5		
Formula	C ₁₂ H ₂₄ N ₂ O ₇		
Mol. weight	308,33 g/mol		
BAA1485	Boc-O₂Oc-O₂Oc-OH		
17-(t-Butyloxycarbonyl-amino)-9-aza-3,6,12,15-tetraoxa-10-on-heptadecanoic acid			
CAS-No.	1069067-08-8		
Formula	C ₁₇ H ₃₂ N ₂ O ₉		
Mol. weight	408,45 g/mol		
FAA1787	Fmoc-O₂Oc-O₂Oc-OH		
17-(9-Fluorenylmethyloxycarbonyl-amino)-9-aza-3,6,12,15-tetraoxa-10-on-heptadecanoic acid			
CAS-No.	560088-89-3		
Formula	C ₂₇ H ₃₄ N ₂ O ₉		
Mol. weight	530,58 g/mol		
PEG2790	N₃-O₂Oc-O₂Oc-OH		
8-(8-Azido-3,6-dioxaoctanoylamido)-3,6-dioxaoctanoic acid			
CAS-No.	1254054-60-8		
Formula	C ₁₂ H ₂₂ N ₄ O ₇		
Mol. weight	334,33 g/mol		







		Product details
PEG2770 $\text{H-O}_2\text{Oc-O}_2\text{Oc-O}_2\text{Oc-OH}$ 26-amino-10,19-dioxo-3,6,12,15,21,24-hexaoxa-9,18-diazahexacosan-1-oic acid CAS-No. 2773558-06-6 Formula $\text{C}_{18}\text{H}_{35}\text{N}_3\text{O}_{10}$ Mol. weight 453,48 g/mol		
PEG8060 $\text{H-O}_2\text{Oc-O}_2\text{Oc-O}_2\text{Oc-O}_2\text{Oc-OH}$ 8-amino-3,6-dioxaoctanoic acid tetramer CAS-No. 2773558-66-8 Formula $\text{C}_{24}\text{H}_{46}\text{N}_4\text{O}_{13}$ Mol. weight 598,64 g/mol		
PEG1320 $\text{H}_2\text{N-PEG(4)-OH}$ 2-(2-(2-(2-Aminoethoxy)ethoxy)ethoxy)ethanol CAS-No. 86770-74-3 Formula $\text{C}_8\text{H}_{19}\text{NO}_4$ Mol. weight 193,24 g/mol		
PEG1335 $\text{H}_2\text{N-PEG(4)-NHNH-Boc}$ 15-Amino-4,7,10,13-tetraoxa-pentadecanoyl-N'-(t-butylloxycarbonyl)-hydrazid CAS-No. 1263047-17-1 Formula $\text{C}_{16}\text{H}_{33}\text{N}_3\text{O}_7$ Mol. weight 379,45 g/mol		
PEG1805 $\text{Fmoc-NH-PEG(4)-NHNH-Boc}$ 15-(9-Fluorenylmethyloxycarbonyl)amino-4,7,10,13-tetraoxa-pentadecanoyl-N'-(t-butylloxycarbonyl)hydrazid CAS-No. 1263044-77-4 Formula $\text{C}_{31}\text{H}_{43}\text{N}_3\text{O}_9$ Mol. weight 601,69 g/mol		
PEG7880 $\text{Boc-NH-PEG(4)-NH}_2$ 1-(t-Butyloxycarbonyl)amino-3,6,9,12-tetraoxatetra-decan-14-amine CAS-No. 811442-84-9 Formula $\text{C}_{15}\text{H}_{32}\text{N}_2\text{O}_6$ Mol. weight 336,42 g/mol		







			Product details
PEG1915	Boc-NH-PEG(4)-OH		
2-(2-(2-(2-(t-Butyloxycarbonylamino)ethoxy)ethoxy)ethoxy)ethanol			
CAS-No.	106984-09-2		
Formula	C ₁₃ H ₂₇ NO ₆		
Mol. weight	293,36 g/mol		
PEG1535	HO-PEG(4)-CO-OtBu		
15-Hydroxy-4,7,10,13-tetraoxa-pentadecanoic acid t-butyl ester			
CAS-No.	518044-32-1		
Formula	C ₁₅ H ₃₀ O ₇		
Mol. weight	322,39 g/mol		
PEG7220	HO-PEG(4)-TFP		
Hydroxy-tetra(ethylene glycol)-propionyl 2,3,5,6-tetrafluorophenyl ester			
CAS-No.	2130036-53-0		
Formula	C ₁₇ H ₂₂ F ₄ O ₇		
Mol. weight	414,35 g/mol		
PEG1370	H₂N-PEG(4)-COOH		
15-Amino-4,7,10,13-tetraoxa-pentadecanoic acid			
CAS-No.	663921-15-1		
Formula	C ₁₁ H ₂₃ NO ₆		
Mol. weight	265,3 g/mol		
PEG1375	H₂N-PEG(4)-CO-OtBu		
15-Amino-4,7,10,13-tetraoxa-pentadecanoic acid t-butyl ester			
CAS-No.	581065-95-4		
Formula	C ₁₅ H ₃₁ NO ₆		
Mol. weight	321,41 g/mol		
PEG1920	Boc-NH-PEG(4)-COOH		
15-t-Butyloxycarbonylamino-4,7,10,13-tetraoxa-pentadecanoic acid			
CAS-No.	756525-91-4		
Formula	C ₁₆ H ₃₁ NO ₈		
Mol. weight	365,42 g/mol		


		Product details
PEG2145	Dnp-NH-PEG(4)-COOH	
1-(2,4-Dinitrophenylamino)-3,6,9,12-tetraoxapentadecanoic acid		
CAS-No.	858126-76-8	
Formula	C ₁₇ H ₂₅ N ₃ O ₁₀	
Mol. weight	431,39 g/mol	
		
PEG2150	Dnp-NH-PEG(4)-NHS	
1-(2,4-Dinitrophenylamino)-3,6,9,12-tetraoxapentadecanoic acid succinimidyl ester		
CAS-No.	858126-78-0	
Formula	C ₂₁ H ₂₈ N ₄ O ₁₂	
Mol. weight	528,47 g/mol	
		
PEG4410	Fmoc-NH-PEG(4)-NHS	
15-(9-Fluorenylmethoxycarbonyl)amino-4,7,10,13-tetraoxa-pentadecanoic acid succinimidyl ester		
CAS-No.	1314378-14-7	
Formula	C ₃₀ H ₃₆ N ₂ O ₁₀	
Mol. weight	584,24 g/mol	
		
PEG7810	Fmoc-NH-PEG(4)-TFP	
15-(9-Fluorenylmethoxycarbonyl)amino-4,7,10,13-tetraoxa-pentadecanoic acid (2,3,5,6-tetrafluorophenyl) ester		
CAS-No.	2247993-77-5	
Formula	C ₃₂ H ₃₃ F ₄ NO ₈	
Mol. weight	635,6 g/mol	
		
PEG2161	Mtt-NH-PEG(4)-COOH*TEA	
1-(p-Methyltritylamino)-3,6,9,12-tetraoxapentadecanoic acid triethylammonium salt		
CAS-No.	1310680-33-1 (net)	
Formula	C ₃₁ H ₃₉ NO ₆ *C ₆ H ₁₅ N	
Mol. weight	C ₃₁ H ₃₉ NO ₆ *C ₆ H ₁₅ N g/mol	
		
PEG2230	OPSS-PEG(4)-NHS	
N-[3-(o-Pyridyldisulfido)propanoyl]-15-amino-4,7,10,13-tetraoxa-pentadecanoyl succinimidyl ester		
CAS-No.	1334177-95-5	
Formula	C ₂₃ H ₃₃ N ₃ O ₉ S ₂	
Mol. weight	559,65 g/mol	
		

			Product details
PEG5080	Phth-NO-PEG(4)-NHS		
1-Phthalimidoxy-3,6,9,12-tetraoxapentadecan-15-oic acid succinimide ester			
CAS-No.	1415328-95-8		
Formula	C ₂₃ H ₂₈ N ₂ O ₁₁		
Mol. weight	508,48 g/mol		
PEG1495	Z-NH-PEG(4)-COOH		
15-Benzyloxycarbonylamino-4,7,10,13-tetraoxa-pentadecanoic acid			
CAS-No.	756526-00-8		
Formula	C ₁₉ H ₂₉ NO ₈		
Mol. weight	399,44 g/mol		
PEG8170	Propargyl-PEG(5)-COOH		
4,7,10,13,16-pentaoxonadec-18-ynoic acid			
CAS-No.	1245823-51-1		
Formula	C ₁₄ H ₂₄ O ₇		
Mol. weight	304,34 g/mol		
PEG5410	Alkyne-PEG(4)-NHS		
Alkyne-PEG(4)-succinimide ester			
CAS-No.	1393330-40-9		
Formula	C ₁₈ H ₂₇ NO ₉		
Mol. weight	401,41 g/mol		
PEG5430	Alkyne-PEG(4)-NH₂		
Alkyne-PEG(4)-amine			
CAS-No.	1013921-36-2		
Formula	C ₁₁ H ₂₁ NO ₄		
Mol. weight	231,29 g/mol		
PEG5440	Alkyne-PEG(4)-mal		
Alkyne-PEG(4)-maleimide			
CAS-No.	1609651-90-2		
Formula	C ₁₈ H ₂₆ N ₂ O ₇		
Mol. weight	382,41 g/mol		





			Product details
TCO1040	TCO-PEG(4)-COOH		
<i>trans</i> -Cyclooctene-PEG(4)-Acid			
CAS-No.	1802913-21-8		
Formula	C ₂₀ H ₃₅ NO ₈		
Mol. weight	417,49 g/mol		
TCO1010	TCO-PEG(4)-NHS		
<i>trans</i> -Cyclooctene-PEG(4)-carboxy succinimidyl ester			
CAS-No.	1621096-79-4		
Formula	C ₂₄ H ₃₈ N ₂ O ₁₀		
Mol. weight	514,57 g/mol		
RL-2510	DBCO-PEG(4)-OH		
Dibenzoazacyclooctyne-tetra(ethylene glycol)			
CAS-No.	1416711-60-8		
Formula	C ₂₉ H ₃₆ N ₂ O ₆		
Mol. weight	508,61 g/mol		
RL-2420	DBCO-PEG(4)-NH₂*TFA		
Dibenzoazacyclooctyne-tetra(ethylene glycol)-amine trifluoro acetic acid salt			
CAS-No.	1255942-08-5		
Formula	C ₂₉ H ₃₇ N ₃ O ₆ *C ₂ F ₃ HO ₂		
Mol. weight	523,62*114,02 g/mol		
RL-2421	DBCO-Sulfo-PEG(4)-NH₂		
Dibenzoazacyclooctyne-tetra(ethylene glycol)amine			
CAS-No.	2055198-05-3		
Formula	C ₃₂ H ₄₂ N ₄ O ₁₀ S		
Mol. weight	674,76 g/mol		
PEG1970	HS-PEG(4)-COOH		
15-Mercapto-4,7,10,13-tertaoxa-pentadecanoic acid			
CAS-No.	749247-06-1		
Formula	C ₁₁ H ₂₂ O ₆ S		
Mol. weight	282,35 g/mol		







			Product details
PEG1740	Mmt-S-PEG(4)-COOH		
15-(4-Methoxytrityl)thio-4,7,10,13-tetraoxa-pentadecanoic acid			
CAS-No.	1263047-31-9		
Formula	C ₃₁ H ₃₈ O ₇ S		
Mol. weight	554,69 g/mol		
			
PEG6710	Trt-S-PEG(4)-COOH*H2O		
15-Tritylmercapto-4,7,10,13-tetraoxapentadecanoic acid monohydrate			
CAS-No.	882847-05-4 net		
Formula	C ₃₀ H ₃₆ O ₅ S*H ₂ O		
Mol. weight	524,67*18,01 g/mol		
			
PEG3590	Lipoamide-PEG(4)-OMe		
alpha-Lipoamide-omega-methoxy tetra(ethylene glycol)			
CAS-No.	1334172-66-5		
Formula	C ₁₇ H ₃₃ NO ₅ S ₂		
Mol. weight	395,58 g/mol		
			
PEG2345	N₃-PEG(4)-COOH		
15-Azido-4,7,10,13-tetraoxa-pentadecanoic acid			
CAS-No.	1257063-35-6		
Formula	C ₁₁ H ₂₁ N ₃ O ₆		
Mol. weight	291,3 g/mol		
			
PEG1400	N₃-PEG(4)-NHS		
15-Azido-4,7,10,13-tetraoxa-pentadecanoic acid succinimidyl ester			
CAS-No.	944251-24-5		
Formula	C ₁₅ H ₂₄ N ₄ O ₈		
Mol. weight	388,37 g/mol		
			
PEG5300	N₃-PEG(4)-OH		
2-(2-(2-(2-Azidoethoxy)ethoxy)ethoxy)ethoxy ethanol			
CAS-No.	86770-68-5		
Formula	C ₁₀ H ₂₁ N ₃ O ₅		
Mol. weight	263,29 g/mol		
			







			Product details
RL-2310	MeTz-PEG(4)-COOH		
Methyltetrazine-PEG(4)-acid			
CAS-No.	1802907-91-0		
Formula	C ₂₀ H ₂₈ N ₄ O ₇		
Mol. weight	436,56 g/mol		
RL-2330	MeTz-PEG(4)-NHS		
Methyltetrazine-PEG(4)-propanoyl succinimidyl ester			
CAS-No.	1802907-92-1		
Formula	C ₂₄ H ₃₁ N ₅ O ₉		
Mol. weight	533,53 g/mol		
RL-2340	MeTz-PEG(4)-mal		
Methyltetrazine-PEG(4)-maleimide			
CAS-No.	1802908-02-6		
Formula	C ₂₄ H ₃₀ N ₆ O ₇		
Mol. weight	514,53 g/mol		
RL-4030	PFB-mercaptopropionyl-PEG3-N₃		
Perfluorobiphenyl-mercaptopropionyl-PEG(3)-N ₃			
Formula	C ₂₃ H ₂₁ F ₉ N ₄ O ₄ S		
Mol. weight	620,49 g/mol		
RL-4040	PFB-mercaptopropionyl-AEEA		
Perfluorobiphenyl-mercaptopropionyl-AEEA			
Formula	C ₂₁ H ₁₆ F ₉ NO ₅ S		
Mol. weight	565,41 g/mol		
RL-4050	PFB-mercaptopropionyl-TOTA-Biotin		
Perfluorobiphenyl-mercaptopropionyl-TOTA-Biotin			
Formula	C ₃₅ H ₄₁ F ₉ N ₄ O ₆ S ₂		
Mol. weight	848,84 g/mol		







	Product details
RL-4060 Biotin-DOOA Biotinyl-1-amino-3,6-dioxa-8-octanamine CAS-No. 138529-46-1 Formula $C_{16}H_{30}N_4O_4S$ Mol. weight 374,50 g/mol	 







2.2. Hydrophobic Spacer Molecules







	Product details
RL-3460 10-Undecynoyl-OSu 10-Undecynoic acid N-hydroxysuccinimide ester CAS-No. 1006592-57-9 Formula $C_{15}H_{21}NO_4$ Mol. weight 279,34 g/mol	 
RL-2055 Alkyne-myristic acid 13-Tetradecynoic acid CAS-No. 82909-47-5 Formula $C_{14}H_{24}O_2$ Mol. weight 224,34 g/mol	 
RL-2060 Alkyne-palmitic acid 15-Hexadecynoic acid CAS-No. 99208-90-9 Formula $C_{16}H_{28}O_2$ Mol. weight 252,39 g/mol	 
RL-2065 Alkyne-stearic acid 17-Octadecynoic acid CAS-No. 34450-18-5 Formula $C_{18}H_{32}O_2$ Mol. weight 280,45 g/mol	 




		Product details
BNN1330	DETA(HBH)*2HCl	
<i>tert</i> -butyl bis(2-aminoethyl)carbamate		
CAS-No.	1914917-65-9	
Formula	C ₉ H ₂₁ N ₃ O ₂ *2HCl	
Mol. weight	203,29*72,92 g/mol	
		
BNN1340	DPTA(BHB)*HCl	
di- <i>tert</i> -butyl (azanediylbis(propane-3,1-diyl))dicarbamate		
CAS-No.	82409-03-8	
Formula	C ₁₆ H ₃₃ N ₃ O ₄ *HCl	
Mol. weight	331,46*36,46 g/mol	
		
BNN1350	DETA(BHH*2HCl)	
<i>tert</i> -butyl (2-((2-aminoethyl)amino)ethyl)carbamate dihydrochloride		
CAS-No.	162279-67-6	
Formula	C ₉ H ₂₁ N ₃ O ₂ *2HCl	
Mol. weight	203,29*72,92 g/mol	
		
BNN1380	Boc-EDA-Suc-OH	
Boc,Succinoyl-ethylenediamine		
CAS-No.	891781-87-6	
Formula	C ₁₁ H ₂₀ N ₂ O ₅	
Mol. weight	260,29 g/mol	
		
HAA9300	H-Aca-Aca-OH	
6-(6-Aminohexanamido)hexanoic acid		
CAS-No.	2014-58-6	
Formula	C ₁₂ H ₂₄ N ₂ O ₃	
Mol. weight	244,34 g/mol	
		
BAA4870	Boc-Aca-Aca-OH	
N-Boc-6-(6-Aminohexanamido)hexanoic acid		
CAS-No.	14254-45-6	
Formula	C ₁₇ H ₃₂ N ₂ O ₅	
Mol. weight	344,45 g/mol	
		

			Product details
RL-2980	N₃-Aca-OSu		
6-Azidocaproic acid N-hydroxysuccinimidyl ester			
CAS-No.	866363-70-4		
Formula	C ₁₀ H ₁₄ N ₄ O ₄		
Mol. weight	254,24 g/mol		
			
RL-4350	N₃-DAPr-Suc-OH		
Azido-propylenediamine-succinoyl-OH			
CAS-No.	929894-58-6		
Formula	C ₇ H ₁₂ N ₄ O ₃		
Mol. weight	200,20 g/mol		
			
RL-4360	N₃-DABu-Suc-OH		
Azido-butylenediamine-succinoyl-OH			
CAS-No.	2226183-50-0		
Formula	C ₈ H ₁₄ N ₄ O ₃		
Mol. weight	214,23 g/mol		
			
HAA6990	N₃-Aca-Aca-OH		
6-(6-azidohexanamido)hexanoic acid			
CAS-No.	866363-71-5		
Formula	C ₁₂ H ₂₂ N ₄ O ₃		
Mol. weight	270,33 g/mol		
			
AAA1970	N₃-Pen-OH		
5-Azido-pentanoic acid			
CAS-No.	79583-98-5		
Formula	C ₅ H ₉ N ₃ O ₂		
Mol. weight	143,14 g/mol		
			
AAA1960	N₃-Hx-OH		
6-Azido-hexanoic acid			
CAS-No.	79598-53-1		
Formula	C ₆ H ₁₁ N ₃ O ₂		
Mol. weight	157,17 g/mol		
			

			Product details
BNN1370	N₃-EDA-Suc-OH		
Azido-ethylenediamine-succinoyl-OH			
CAS-No.	2225891-73-4		
Formula	C ₆ H ₁₀ N ₄ O ₃		
Mol. weight	186,17 g/mol		
			
RL-4430	N₃-C7H14-COOH		
8-azidooctanoic acid			
CAS-No.	217180-76-2		
Formula	C ₈ H ₁₅ N ₃ O ₂		
Mol. weight	185,23 g/mol		
			
RL-3480	8-Azido-octanoyl-OSu		
8-Azidooctanoic acid N-hydroxysuccinimide ester			
CAS-No.	2576471-56-0		
Formula	C ₁₂ H ₁₈ N ₄ O ₄		
Mol. weight	282,30 g/mol		
			
RL-3170	11-Azido-undecanoyl-OSu		
11-Azidoundecanoic acid N-hydroxysuccinimide ester			
CAS-No.	850080-13-6		
Formula	C ₁₅ H ₂₄ N ₄ O ₄		
Mol. weight	324,38 g/mol		
			
RL-3200	11-Azidoundecanoic acid		
11-Azido-undecanoic acid			
CAS-No.	118162-45-1		
Formula	C ₁₁ H ₂₁ N ₃ O ₂		
Mol. weight	227,30 g/mol		
			
RL-3210	12-Azidododecanoic acid		
12-Azido-dodecanoic acid			
CAS-No.	80667-36-3		
Formula	C ₁₂ H ₂₃ N ₃ O ₂		
Mol. weight	241,33 g/mol		
			

			Product details
RL-3220	12-Azido-dodecanoyl-OSu		
12-Azidododecanoic acid N-hydroxysuccinimide ester			
CAS-No.	2489524-00-5		
Formula	C ₁₆ H ₂₆ N ₄ O ₄		
Mol. weight	338,40 g/mol		
RL-3230	14-Azido-myristic acid		
14-azidotetradecanoic acid			
CAS-No.	176108-61-5		
Formula	C ₁₄ H ₂₇ N ₃ O ₂		
Mol. weight	269,38 g/mol		
RL-3240	16-Azido-palmitic acid		
16-azidohexadecanoic acid			
CAS-No.	112668-54-9		
Formula	C ₁₆ H ₃₁ N ₃ O ₂		
Mol. weight	297,44 g/mol		
RL-3250	18-Azido-stearic acid		
18-azidooctadecanoic acid			
CAS-No.	1529763-58-3		
Formula	C ₁₈ H ₃₅ N ₃ O ₂		
Mol. weight	325,49 g/mol		
BAA4240	14-(Boc-amino)-myristic acid		
14-((t-Butyloxycarbonyl)amino)tetradecanoic acid			
CAS-No.	2307778-46-5		
Formula	C ₁₉ H ₃₇ NO ₄		
Mol. weight	343,51 g/mol		
FAA8160	14-(Fmoc-amino)-myristic acid		
14-((9-Fluorenylmethyloxycarbonyl)amino)tetradecanoic acid			
CAS-No.	1931109-55-5		
Formula	C ₂₉ H ₃₉ NO ₄		
Mol. weight	465,63 g/mol		


			Product details
BAA3900	16-(Boc-amino)-palmitic acid		
16-((t-Butyloxycarbonyl)amino)hexadecanoic acid			
CAS-No.	135747-73-8		
Formula	C ₂₁ H ₄₁ NO ₄		
Mol. weight	371,55 g/mol		
FAA7460	16-(Fmoc-amino)-palmitic acid		
16-((9-Fluorenylmethyloxycarbonyl)amino)hexadecanoic acid			
CAS-No.	1356220-22-8		
Formula	C ₃₁ H ₄₃ NO ₄		
Mol. weight	493,68 g/mol		
BAA3910	18-(Boc-amino)-stearic acid		
18-((t-Butyloxycarbonyl)amino)octadecanoic acid			
CAS-No.	2389064-45-1		
Formula	C ₂₃ H ₄₅ NO ₄		
Mol. weight	399,61 g/mol		
FAA7450	18-(Fmoc-amino)-stearic acid		
18-((9-Fluorenylmethyloxycarbonyl)amino)octadecanoic acid			
CAS-No.	1199580-37-4		
Formula	C ₃₃ H ₄₇ NO ₄		
Mol. weight	521,73 g/mol		
RL-4020	DBCO-C6-Alkyne		
N-(propargylamidoadipoyl)-dibenzoazacyclooctyne			
Formula	C ₂₄ H ₂₂ N ₂ O ₂		
Mol. weight	370,45 g/mol		
RL-1002	FMPB-Linker		
4-(4'-Formyl-3'-methoxyphenoxy) butanoic acid			
CAS-No.	309964-23-6		
Formula	C ₁₂ H ₁₄ O ₅		
Mol. weight	238,24 g/mol		

			Product details
RL-1008	Br-PAM-Linker		
4-Bromomethylphenyl-acetic acid			
CAS-No.	13737-36-5		
Formula	C ₉ H ₉ BrO ₂		
Mol. weight	229,1 g/mol		
			
RL-1050	H-PAL-Linker		
5-[4-Aminomethyl-3,5-dimethoxyphenoxy]-pentanoic acid			
CAS-No.	125666-66-2		
Formula	C ₁₆ H ₂₁ NO ₅		
Mol. weight	283,31 g/mol		
			
RL-1114	HMPB-Linker		
4-(4-Hydroxymethyl-3-methoxyphenoxy)-butyric acid			
CAS-No.	136849-75-7		
Formula	C ₁₂ H ₁₆ O ₅		
Mol. weight	240,24 g/mol		
			




You could not find the product you are looking for?

Please contact us for a custom synthesis!



2.3. Permanent Linkers with Maleimide Function

Michael addition of a thiol to a maleimide is commonly used for numerous bioconjugations. Several commercial constructs like Brentuximab vedotin, Trastuzumab emtansine, and Cimzia contain a thiol-maleimide adduct. However, this reaction is reversible. During the journey of an appropriate thioether containing drug through physiological media, this bond can break, and fragments are released which might contribute to certain unwanted or even toxic reactions.








However, if the succinimide moiety of a maleimide-thiol conjugate is hydrolyzed, the ring-opened product is fully stabilized towards cleavage (Fig. 15). The rates of ring-opening hydrolysis are greatly accelerated by electron withdrawing N-substituents and good nucleophiles in the proximity of the carbonyl functions. Thus, conjugates made with nucleophilic side-chains and electron-withdrawing maleimides may be purposefully hydrolyzed to their ring-opened counterparts and ensure good *in vivo* stability.







References:







- *Covalent Modification of Biomolecules through Maleimide-Based Labeling Strategies*; K. Renault, J. W. Fredy, P. Y. Renard, C. Sabot; **Bioconjug Chem** 2018; **29**: 2497-2513. <https://doi.org/10.1021/acs.bioconjchem.8b00252>
- *Optimisation of the dibromomaleimide (DBM) platform for native antibody conjugation by accelerated post-conjugation hydrolysis*; M. Morais, J. P. M. Nunes, K. Karu, N. Forte, I. Benni, M. E. B. Smith, S. Caddick, V. Chudasama, J. R. Baker; **Org Biomol Chem** 2017; **15**: 2947-2952. <https://doi.org/10.1039/c7ob00220c>
- *Use of a next generation maleimide in combination with THIOMAB™ antibody technology delivers a highly stable, potent and near homogeneous THIOMAB™ antibody-drug conjugate (TDC)*; J. P. M. Nunes, V. Vassileva, E. Robinson, M. Morais, M. E. B. Smith, R. B. Pedley, S. Caddick, J. R. Baker, V. Chudasama; **RSC Advances** 2017; **7**: 24828-24832. <https://doi.org/10.1039/c7ra04606e>







[↑ back to content](#)







- Long-term stabilization of maleimide-thiol conjugates; S. D. Fontaine, R. Reid, L. Robinson, G. W. Ashley, D. V. Santi; *Bioconjug Chem* 2015; **26**: 145-52. <https://doi.org/10.1021/bc5005262>
- Self-hydrolyzing maleimides improve the stability and pharmacological properties of antibody-drug conjugates; R. P. Lyon, J. R. Setter, T. D. Bovee, S. O. Doronina, J. H. Hunter, M. E. Anderson, C. L. Balasubramanian, S. M. Duniho, C. I. Leiske, F. Li, P. D. Senter; *Nat Biotechnol* 2014; **32**: 1059-62. <https://doi.org/10.1038/nbt.2968>
- Mild method for succinimide hydrolysis on ADCs: impact on ADC potency, stability, exposure, and efficacy; L. N. Tumey, M. Charati, T. He, E. Sousa, D. Ma, X. Han, T. Clark, J. Casavant, F. Loganzo, F. Barletta, J. Lucas, E. I. Graziani; *Bioconjug Chem* 2014; **25**: 1871-80. <https://doi.org/10.1021/bc500357n>







		Product details	
RL-3000	Mal-Et-OH		
N-(2-Hydroxyethyl)maleimide			
CAS-No.	1585-90-6		
Formula	C ₆ H ₇ NO ₃		
Mol. weight	141,12 g/mol		
RL-2780	Mal-NH₂*HCl		
2-Maleimidoethylamine hydrochloride			
CAS-No.	134272-64-3		
Formula	C ₆ H ₈ N ₂ O ₂ *HCl		
Mol. weight	140,14*36,45 g/mol		
RL-2490	DBCO-mal		
Dibenzoazacyclooctyne-maleimide			
CAS-No.	1395786-30-7		
Formula	C ₂₅ H ₂₁ N ₃ O ₄		
Mol. weight	427,45 g/mol		
RL-2500	DBCO-PEG(4)-mal		
Dibenzoazacyclooctyne-tetra(ethylene glycol)-maleimide			
CAS-No.	1480516-75-3		
Formula	C ₃₆ H ₄₂ N ₄ O ₉		
Mol. weight	674,74 g/mol		
RL-2340	MeTz-PEG(4)-mal		
Methyltetrazine-PEG(4)-maleimide			
CAS-No.	1802908-02-6		
Formula	C ₂₄ H ₃₀ N ₆ O ₇		
Mol. weight	514,53 g/mol		



			Product details
PEG4870	Mal-O₂Oc-OH		
{2-[2-(2,5-Dioxo-2,5-dihydro-1H-pyrrol-1-yl)ethoxy]ethoxy}acetic acid			
CAS-No.	173323-23-4		
Formula	C ₁₀ H ₁₃ NO ₆		
Mol. weight	243,21 g/mol		
PEG1555	mal-PEG(2)-COOH		
3-(2-(2-(3-Maleimidopropanamido)ethoxy)ethoxy)propanoic acid			
CAS-No.	756525-98-1		
Formula	C ₁₆ H ₂₀ N ₂ O ₇		
Mol. weight	328,32 g/mol		
PEG1560	mal-PEG(2)-NHS		
3-(2-(2-(3-Maleimidopropanamido)ethoxy)ethoxy)propanoic acid succinimidyl ester			
CAS-No.	955094-26-5		
Formula	C ₁₈ H ₂₃ N ₃ O ₉		
Mol. weight	425,39 g/mol		
MAA1020	Mal-beta-Ala-OSu		
3-(Maleimido)propionic acid N-succinimidyl ester			
CAS-No.	55750-62-4		
Formula	C ₁₁ H ₁₀ N ₂ O ₆		
Mol. weight	266,21 g/mol		
RL-2640	Mal-Bu-NHS		
4-Maleimidobutyric acid-NHS ester			
CAS-No.	80307-12-6		
Formula	C ₁₂ H ₁₂ N ₂ O ₆		
Mol. weight	280,23 g/mol		
RL-2670	Mal-Pen-NHS		
5-Maleimidopentanoic acid-NHS ester			
CAS-No.	103750-03-4		
Formula	C ₁₃ H ₁₄ N ₂ O ₆		
Mol. weight	294,26 g/mol		

			Product details
RL-2660	Mal-Hx-NHS		
6-Maleimidohexanoic acid-NHS ester			
CAS-No.	55750-63-5		
Formula	C ₁₄ H ₁₆ N ₂ O ₆		
Mol. weight	308,29 g/mol		
RL-2690	Mal-PrHx-NHS		
6-(3-Maleimidopropionylamino)-hexanoic acid-NHS ester			
CAS-No.	367927-39-7		
Formula	C ₁₇ H ₂₁ N ₃ O ₇		
Mol. weight	379,36 g/mol		
MAA1040	Mal-L-Dap(Boc)-OH*DCHA		
N-alpha-MaleimidoN-beta-t-butyloxycarbonyl-L-2,3-diaminopropionic acid dicyclohexylamine			
CAS-No.	2004724-16-5		
Formula	C ₁₂ H ₁₆ N ₂ O ₆ *C ₁₂ H ₂₃ N		
Mol. weight	284,27*181,32 g/mol		
MAA1080	Mal-L-Dap(Boc)-OPfp		
N-alpha-MaleimidoN-beta-t-butyloxycarbonyl-L-2,3-diaminopropionic acid pentafluorophenolate			
CAS-No.	1887132-90-2		
Formula	C ₁₈ H ₁₅ F ₅ N ₂ O ₆		
Mol. weight	450,31 g/mol		
MAA1060	Mal-D-Dap(Boc)-OH*DCHA		
N-alpha-MaleimidoN-beta-t-butyloxycarbonyl-D-2,3-diaminopropionic acid dicyclohexylamine			
CAS-No.	2382651-11-6 net		
Formula	C ₁₂ H ₁₆ N ₂ O ₆ *C ₁₂ H ₂₃ N		
Mol. weight	284,27*181,32 g/mol		
BAA6475	Boc-L-Dap(Mal)-OH*DCHA		
(2S)-2-((tert-butoxycarbonyl)amino)-3-(2,5-dioxo-2,5-dihydro-1H-pyrrol-1-yl)propanoic acid dicyclohexylamine			
Formula	C ₁₂ H ₁₆ N ₂ O ₆ *C ₁₂ H ₂₃ N		
Mol. weight	284,27*181,32 g/mol		

		Product details
BAA6480	Boc-L-Dab(Mal)-OH (2S)-2-((tert-butoxycarbonyl)amino)-4-(2,5-dioxo-2,5-dihydro-1H-pyrrol-1-yl)butanoic acid CAS-No. 135631-02-6 Formula $C_{13}H_{18}N_2O_6$ Mol. weight 298,30 g/mol	 
MAA1120	Mal-L-Dap(Boc)-OSu N-alpha-Maleimido-N-beta-Boc-L-2,3-diaminopropionic acid NHS ester CAS-No. 1703778-79-3 Formula $C_{16}H_{19}N_3O_8$ Mol. weight 381,34 g/mol	 
MAA5400	Mal-AMCHC-OH <i>trans</i> -4-(maleimidomethyl)cyclohexane-1-carboxylic acid CAS-No. 69907-67-1 Formula $C_{12}H_{15}NO_4$ Mol. weight 237,25 g/mol	 
MAA1000	Mal-AMCHC-OSu <i>trans</i> -N-Succinimidyl 4-(maleimidomethyl)cyclohexane-1-carboxylate CAS-No. 71875-81-5 Formula $C_{16}H_{18}N_2O_6$ Mol. weight 334,33 g/mol	 
MAA1050	Sulfo-SMCC 4-(N-Maleimidomethyl)cyclohexane-1-carboxylic acid 3-sulfo-N-hydroxysuccinimide ester sodium salt (<i>cis/trans</i> mixture) CAS-No. 92921-24-9 Formula $C_{16}H_{17}N_2NaO_9S$ Mol. weight 436,37 g/mol	 
RL-2650	Mal-cHxHx-NHS 6-[<i>trans</i> -4-(Maleimidomethyl)-cyclohexanoylamino]-hexanoic acid-NHS ester CAS-No. 125559-00-4 Formula $C_{22}H_{29}N_3O_7$ Mol. weight 447,48 g/mol	 

			Product details
RL-2600	3-Mal-Bz-NHS		
3-Maleimidobenzoic acid-NHS ester			
CAS-No.	58626-38-3		
Formula	C ₁₅ H ₁₀ N ₂ O ₆		
Mol. weight	314,25 g/mol		
			
RL-2610	3-Mal-MBz-NHS		
3-(Maleimidomethyl)-benzoic acid-NHS ester			
CAS-No.	91574-36-6		
Formula	C ₁₆ H ₁₂ N ₂ O ₆		
Mol. weight	328,28 g/mol		
			
RL-2620	4-Mal-Bz-NHS		
4-Maleimidobenzoic acid-NHS ester			
CAS-No.	64191-06-6		
Formula	C ₁₅ H ₁₀ N ₂ O ₆		
Mol. weight	314,25 g/mol		
			
RL-2630	4-Mal-MBz-NHS		
4-(Maleimidomethyl)-benzoic acid-NHS ester			
CAS-No.	64987-84-4		
Formula	C ₁₆ H ₁₂ N ₂ O ₆		
Mol. weight	328,28 g/mol		
			
RL-2680	Mal-PhBu-NHS		
4-(4-Maleimidophenyl)-butyric acid-NHS ester			
CAS-No.	79886-55-8		
Formula	C ₁₈ H ₁₆ N ₂ O ₆		
Mol. weight	356,33 g/mol		
			
LS-3350	4-(N-Maleimido)benzophenone		
1-(4-Benzoylphenyl)-1H-pyrrole-2,5-dione			
CAS-No.	92944-71-3		
Formula	C ₁₇ H ₁₁ NO ₃		
Mol. weight	277,28 g/mol		
			


		Product details
RL-3430	Mal-N-Boc-Aeg-NHS N-(t-butoxycarbonyl)-N-(2-(maleinimido)ethyl)glycine N-Hydroxysuccinimidyl ester CAS-No. 2576471-29-7 Formula $C_{17}H_{21}N_3O_8$ Mol. weight 395,37 g/mol	 
RL-3450	Mal-CH₂CH₂-N-(CH₂-COOH)₂ 2,2'-((2-(maleinimido)ethyl)azanediyl)diacetic acid CAS-No. 207612-92-8 Formula $C_{10}H_{12}N_2O_6$ Mol. weight 256,21 g/mol	 
RL-3400	Mal-CH₂CH₂-N(Me)-CH₂-COOH N-(2-(maleinimido)ethyl)-N-methylglycine CAS-No. 2576471-52-6 Formula $C_9H_{12}N_2O_4$ Mol. weight 212,21 g/mol	 
MAA1100	Mal-AMCHC-N-Propargylamide <i>trans</i> -4-[(2,5-dioxo-2,5-dihydro-1H-pyrrol-1-yl)methyl]- N-(prop-2-yn-1-yl)cyclohexane-1-carboxamide CAS-No. 2027476-42-0 Formula $C_{15}H_{18}N_2O_3$ Mol. weight 274,32 g/mol	 
RL-8420	Biotin-NH-NH-Mal 6-(2,5-dioxo-2,5-dihydro-1H-pyrrol-1-yl)-N'-(5- ((3aS,4S,6aR)-2-oxohexahydro-1H-thieno[3,4-d]imidazol-4-yl)pentanoyl)hexanehydrazide CAS-No. 116919-18-7 Formula $C_{20}H_{29}N_5O_5S$ Mol. weight 451,54 g/mol	 
RL-8425	Biocytin-Mal N2-(3-(2,5-dioxo-2,5-dihydro-1H-pyrrol-1-yl)propanoyl)- N6-(5-((3aS,4S,6aR)-2-oxohexahydro-1H-thieno[3,4-d] imidazol-4-yl)pentanoyl)-L-lysine CAS-No. 102849-12-7 Formula $C_{23}H_{33}N_5O_5S$ Mol. weight 523,61 g/mol	 

			Product details
TCOI050	TCO-PEG(3)-mal		
<i>trans</i> -Cyclooctene-PEG(3)-maleimide			
CAS-No.	1809356-72-6		
Formula	C ₂₆ H ₄₁ N ₃ O ₈		
Mol. weight	523,62 g/mol		
PEG1485	mal-PEG(3)-mal		
Bis-(1,13-(3-maleimidopropionyl)amido)-4,7,10-trioxatridecane			
CAS-No.	756525-89-0		
Formula	C ₂₆ H ₃₄ N ₄ O ₉		
Mol. weight	522,55 g/mol		

2.4. 5HP2O as Maleimide Alternative

Despite various possibilities for the modification of and conjugation to side-chain functional groups of histidine, lysine, methionine, tryptophane, and tyrosine, respectively, cysteine remains the most attractive target due to its rare occurrence in natural proteins typically allowing for its site-selective modification upon introduction at a specific position. In terms of cysteine modification, maleimides are the reactive moieties of choice acting *via* a Michael addition reaction. However, maleimides are susceptible to hydrolysis, retro-Michael reactions and/or thiol exchange reactions leading to protein conjugates with overall poor stability.

In this context, starting from furan, the research group around Prof. Annemieke Madder developed 5-hydroxy-1,5-dihydro-2H-pyrrol-2-ones (5HP2Os) – representatives of the 5-hydroxy-pyrrolone family – as maleimide replacement technology avoiding all the above-mentioned maleimide-drawbacks. The substituents R₁ and R₂ – as shown in the scheme below – allow for further derivatization and can be chosen depending on the intended subsequent use.












Interested in more details about 5HP2O
as maleimide replacement technology?

Watch the recording of our online workshop!



		Product details	
RL-8670	5HP2O-alkyne		
5-hydroxy-5-methyl-1-(prop-2-yn-1-yl)-1,5-dihydro-2H-pyrrol-2-one			
CAS-No.	2484704-61-0		
Formula	C ₈ H ₉ NO ₂		
Mol. weight	151,16 g/mol		
			
RL-8675	5HP2O-(CH₂)₄-COOH		
5-(2-hydroxy-2-methyl-5-oxo-2,5-dihydro-1H-pyrrol-1-yl)pentanoic acid			
Formula	C ₁₀ H ₁₅ NO ₄		
Mol. weight	213,23 g/mol		
			
RL-8680	5HP2O-(CH₂)₄-NHS		
2,5-dioxopyrrolidin-1-yl 5-(2-hydroxy-2-methyl-5-oxo-2,5-dihydro-1H-pyrrol-1-yl)pentanoate			
Formula	C ₁₄ H ₁₈ N ₂ O ₆		
Mol. weight	310,31 g/mol		
			
RL-8685	5HP2O-PEG(2)-COOH		
3-(2-(2-(2-hydroxy-2-methyl-5-oxo-2,5-dihydro-1H-pyrrol-1-yl)ethoxy)ethoxy)propanoic acid			
Formula	C ₁₂ H ₁₉ NO ₆		
Mol. weight	273,28 g/mol		
			


		Product details
RL-8690 5HP2O-PEG(2)-NHS 2,5-dioxopyrrolidin-1-yl 3-(2-(2-(2-hydroxy-2-methyl-5-oxo-2,5-dihydro-1H-pyrrol-1-yl)ethoxy)ethoxy)propanoate Formula $C_{16}H_{22}N_2O_8$ Mol. weight 370,36 g/mol		
LS-4670 5HP2O((PEG)2-OH)-(CH2)5-Dansyl 2-(2-(2-hydroxyethoxy)ethoxy)ethyl 3-(1-(5-((5-(dimethylamino)naphthalene-1-sulfonamido)pentyl)-2-hydroxy-5-oxo-2,5-dihydro-1H-pyrrol-2-yl)propanoate Formula $C_{30}H_{43}N_3O_9S$ Mol. weight 621,75 g/mol		
RL-8695 5HP2O((PEG)2-OH)-(CH2)4-NHS 2,5-dioxopyrrolidin-1-yl 5-(2-hydroxy-2-(3-(2-(2-(2-hydroxyethoxy)ethoxy)ethoxy)-3-oxopropyl)-5-oxo-2,5-dihydro-1H-pyrrol-1-yl)pentanoate Formula $C_{22}H_{32}N_2O_{11}$ Mol. weight 500,50 g/mol		

References:

- *Bioconjugation Reagent and Methods*. A. Madder, E. De Geyter, E. Antonatou, S. Smolen, D. Kalaitzakis, D. Vassilikogiannakis; 2020. WO2020/174086A2
- *2-Hydroxypyrrolidin-5-ones for bioconjugation and methods for their production*. A. Madder, E. De Geyter, E. Antonatou, S. Smolen, D. Kalaitzakis, D. Vassilikogiannakis; 2020. WO2020/174086A3
- *5-Hydroxy-pyrrolone based building blocks as maleimide alternatives for protein bioconjugation and single-site multi-functionalization*; E. De Geyter, E. Antonatou, D. Kalaitzakis, S. Smolen, A. Iyer, L. Tack, E. Ongenae, G. Vassilikogiannakis, A. Madder; *Chem. Sci.* 2021; **12**: 5246-5252. <https://doi.org/10.1039/d0sc05881e>







2.5. Photoactivatable Linkers






Irradiation of diazirines with UV light (ca. 350-360 nm) yields a highly reactive carbene species that can undergo insertions into C–C, C–H, O–H, and X–H (X = heteroatom) bonds of neighboring molecules to irreversibly form a covalent bond (Fig. 17). The diazirine moiety is the smallest of all photoreactive groups, so introduction of a diazirine-bearing amino acid into a peptide or protein usually does not impair its biological activity. Further, advantages of diazirine crosslinkers are their stability at room temperature and their relative stability against nucleophiles as well as towards both acidic and basic conditions.



References:

- Protein-polymer conjugation via ligand affinity and photoactivation of glutathione S-transferase; E. W. Lin, N. Boehnke, H. D. Maynard; **Bioconjug Chem** 2014; **25**: 1902-9. <https://doi.org/10.1021/bc500380r>
- Cell-based proteome profiling of potential dasatinib targets by use of affinity-based probes; H. Shi, C. J. Zhang, G. Y. Chen, S. Q. Yao; **J Am Chem Soc** 2012; **134**: 3001-14. <https://doi.org/10.1021/ja208518u>
- Probing protein-protein interactions with a genetically encoded photo-crosslinking amino acid; H. W. Ai, W. Shen, A. Sagi, P. R. Chen, P. G. Schultz; **ChemBioChem** 2011; **12**: 1854-7. <https://doi.org/10.1002/cbic.201100194>
- Proteome profiling reveals potential cellular targets of staurosporine using a clickable cell-permeable probe; H. Shi, X. Cheng, S. K. Sze, S. Q. Yao; **Chem Commun** 2011; **47**: 11306-8. <https://doi.org/10.1039/c1cc14824a>
- Direct interaction between an allosteric agonist pepducin and the chemokine receptor CXCR4; J. M. Janz, Y. Ren, R. Looby, M. A. Kazmi, P. Sachdev, A. Grunbeck, L. Haggis, D. Chinnappen, A. Y. Lin, C. Seibert, T. McMurry, K. E. Carlson, T. W. Muir, S. Hunt, 3rd, T. P. Sakmar; **J Am Chem Soc** 2011; **133**: 15878-81. <https://doi.org/10.1021/ja206661w>
- Aliphatic diazirines as photoaffinity probes for proteins: recent developments; J. Das; **Chem Rev** 2011; **111**: 4405-17. <https://doi.org/10.1021/cr100272z>
- Synthesis and application of photoproline - a photoactivatable derivative of proline; B. VanderMeijden, J. A. Robinson; **Arkivoc** 2011; **2011**: 130-136. <https://doi.org/10.3998/ark.5550190.0012.611>
- Photo-crosslinking of proteins in intact cells reveals a dimeric structure of cyclooxygenase-2 and an inhibitor-sensitive oligomeric structure of microsomal prostaglandin E2 synthase-1; P. O. Hetu, M. Ouellet, J. P. Falgoutyret, C. Ramachandran, J. Robichaud, R. Zamboni, D. Riendeau; **Arch Biochem Biophys** 2008; **477**: 155-62. <https://doi.org/10.1016/j.abb.2008.04.038>
- Covalent capture of phospho-dependent protein oligomerization by site-specific incorporation of a diazirine photo-cross-linker; M. Vila-Perello, M. R. Pratt, F. Tulin, T. W. Muir; **J Am Chem Soc** 2007; **129**: 8068-9. <https://doi.org/10.1021/ja072013j>
- Photo-leucine incorporation reveals the target of a cyclodepsipeptide inhibitor of cotranslational translocation; A. L. MacKinnon, J. L. Garrison, R. S. Hegde, J. Taunton; **J Am Chem Soc** 2007; **129**: 14560-1. <https://doi.org/10.1021/ja076250y>
- Synthesis of photoactive analogues of a cystine knot trypsin inhibitor protein; T. Durek, J. Zhang, C. He, S. B. Kent; **Org Lett** 2007; **9**: 5497-500. <https://doi.org/10.1021/ol702461z>
- Photo-leucine and photo-methionine allow identification of protein-protein interactions in living cells; M. Suchanek, A. Radzikowska, C. Thiele; **Nat Methods** 2005; **2**: 261-7. <https://doi.org/10.1038/nmeth752>

			Product details
RL-2910	Photo-Butylamine		
2-(3-methyl-3H-diazirin-3-yl)ethan-1-amine hydrochloride			
CAS-No.	25055-95-2		
Formula	C ₄ H ₈ N ₃ *HCl		
Mol. weight	99,13*36,45 g/mol		
RL-2890	Photo-Pentanoic acid		
3-(3-methyl-3H-diazirin-3-yl)propanoic acid			
CAS-No.	25055-86-1		
Formula	C ₅ H ₈ N ₂ O ₂		
Mol. weight	128,13 g/mol		
RL-2900	Photo-Hexanoic acid		
4-(3-methyl-3H-diazirin-3-yl)butanoic acid			
CAS-No.	16297-97-5		
Formula	C ₆ H ₁₀ N ₂ O ₂		
Mol. weight	142,16 g/mol		
RL-2920	Photo-Benzoic acid		
4-[3-(Trifluoromethyl)-3H-diazirin-3-yl]benzoic acid			
CAS-No.	85559-46-2		
Formula	C ₉ H ₅ F ₃ N ₂ O ₂		
Mol. weight	230,14 g/mol		
RL-2930	Photo-Benzylamine*HCl		
4-[3-(Trifluoromethyl)-3H-diazirin-3-yl]benzylamine hydrochloride			
CAS-No.	1258874-29-1		
Formula	C ₉ H ₈ N ₃ F ₃ *HCl		
Mol. weight	215,18*36,45 g/mol		
RL-3410	Photo-Click-Heptanoic acid		
2-(3-(but-3-ynyl)-3H-diazirin-3-yl)acetic acid			
CAS-No.	2049109-24-0		
Formula	C ₇ H ₈ N ₂ O ₂		
Mol. weight	152,15 g/mol		

		Product details
RL-3270	Fmoc-AEEP-DIM	
3-(2-(2-(9-Fluorenylmethyl)oxycarbonylaminoethoxy)ethoxy)-1-(4,4-dimethyl-2,6-dioxocyclohexylidene)-propan-1-ol		
CAS-No.	1988771-96-5	
Formula	C ₃₀ H ₃₅ NO ₇	
Mol. weight	521,60 g/mol	
		
RL-3260	Fmoc-Aca-DIM	
6-((9-Fluorenylmethyl)oxycarbonylamino)-1-(4,4-dimethyl-2,6-dioxocyclohexylidene)-hexan-1-ol		
CAS-No.	2379561-08-5	
Formula	C ₂₉ H ₃₃ NO ₅	
Mol. weight	475,58 g/mol	
		
RL-2035	ATFB	
4-Azido-2,3,5,6-tetrafluorobenzoic acid		
CAS-No.	122590-77-6	
Formula	C ₇ HF ₄ N ₃ O ₂	
Mol. weight	235,1 g/mol	
		
RL-2045	ATFB-NHS	
N-Succinimidyl 4-azido-2,3,5,6-tetrafluorobenzoate		
CAS-No.	126695-58-7	
Formula	C ₁₁ H ₄ F ₄ N ₄ O ₄	
Mol. weight	332,17 g/mol	
		
PEG5000	N₃-TFBA-O₂Oc	
{2-[2-(4-Azido-2,3,5,6-tetrafluorobenzoyl-amino)ethoxy]ethoxy}acetic acid		
CAS-No.	1993119-45-1	
Formula	C ₁₃ H ₁₂ F ₄ N ₄ O ₅	
Mol. weight	380,25 g/mol	
		

References:

- Tri- and Tetravalent Photoactivable Cross-Linking Agents; A. Welle, F. Billard, J. Marchand-Brynaert; **Synthesis** 2012; **44**: 2249-2254. <https://doi.org/10.1055/s-0031-1290444>
- Perfluorophenyl azides: new applications in surface functionalization and nanomaterial synthesis; L. H. Liu, M. Yan; **Acc Chem Res** 2010; **43**: 1434-43. <https://doi.org/10.1021/ar100066t>
- Chemistry of Bifunctional Photoprobes; R. S. Pandurangi, P. Lusiak, S. Desai, R. R. Kuntz; **Bioorganic Chemistry** 1998; **26**: 201-212. <https://doi.org/10.1006/bioo.1998.1098>
- Recent Trends in the Evaluation of Photochemical Insertion Characteristics of Heterobifunctional Perfluoroaryl Azide Chelating Agents: Biochemical Implications in Nuclear Medicine; R. S. Pandurangi, S. R. Karra, R. R. Kuntz, W. A. Volkert; **Photochemistry and Photobiology** 1997; **65**: 208-221. <https://doi.org/10.1111/j.1751-1097.1997.tb08547.x>
- Chemistry of Bifunctional Photoprobes. 1. Perfluoroaryl Azido Functionalized Phosphorus Hydrazides as Novel Photoreactive Heterobifunctional Chelating Agents: High Efficiency Nitrene Insertion on Model Solvents and Proteins; R. S. Pandurangi, S. R. Karra, K. V. Katti, R. R. Kuntz, W. A. Volkert; **J Org Chem** 1997; **62**: 2798-2807. <https://doi.org/10.1021/jo961867b>



For more information about photochemistry, discover our dedicated brochure!



3. Cleavable Linkers

Peptidic bonds are expected to have a high serum stability, as lysosomal proteolytic enzymes show reduced activities in blood due to endogenous inhibitors and the unfavorably high pH value of blood compared to lysosomes. This was confirmed by preclinical *in vivo* studies, which revealed half-lives of seven to ten days for peptide linkers. Release of a drug conjugated *via* a peptidyl linker to monoclonal antibodies (mAb) occurs specifically due to the action of lysosomal proteases (e.g., cathepsin and plasmin). These proteases may be present at elevated levels in certain tumor tissues. Therefore, peptide linkers combine greater systemic stability with rapid enzymatic release of the drug in the target cell. Besides Val-Ala, Val-Cit and Phe-Lys, other sequences have been reported as lysosomally cleavable peptides, like Gly-Phe-Leu-Gly and Ala-Leu-Ala-Leu.

References:

- *Star structure of antibody-targeted HPMA copolymer-bound doxorubicin: a novel type of polymeric conjugate for targeted drug delivery with potent antitumor effect*; M. Kovar, J. Strohm, T. Etrych, K. Ulbrich, B. Rihova; **Bioconjug Chem** 2002; **13**: 206-15. <https://doi.org/10.1021/bc010063m>
- *Synthesis of a lipophilic daunorubicin derivative and its incorporation into lipidic carriers developed for LDL receptor-mediated tumor therapy*; A. J. Versluis, E. T. Rump, P. C. Rensen, T. J. Van Berkel, M. K. Bijsterbosch; **Pharm Res** 1998; **15**: 531-7. <https://doi.org/10.1023/a:1011917508056>
- *Influence of a peptide linker on biodistribution and metabolism of antibody-conjugated benzyl-EDTA. Comparison of enzymatic digestion in vitro and in vivo*; M. Studer, L. A. Kroger, S. J. DeNardo, D. L. Kukis and C. F. Meares; **Bioconjug Chem** 1992; **3**: 424-9. <https://doi.org/10.1021/bc00017a012>






Fig. 18: Cyclobutane-1,1-dicarboxamide can replace valine in dipeptide linker systems, resulting in improved ADC selectivity.

Peptide-based ADC linkers, such as Val-Cit or Val-Ala, that are cleaved by lysosomal proteases have shown sufficient stability in serum and effective payload-release in targeted cells. However, the use of peptide-based linkers limits the ability to modulate protease specificity. Furthermore, if the linker can preferentially be hydrolyzed by tumor-specific proteases only, safety margin may improve. In this context, a cyclobutane-1,1-dicarboxamide-containing linker (Fig. 18) replacing valine in other sequences has been invented which is hydrolyzed predominantly by cathepsin B, while the typical valine-citrulline dipeptide linker is rather less. ADCs bearing the nonpeptidic linker are as efficacious and stable *in vivo* as those with the dipeptide linker. Hence, the application of the peptidomimetic linker presents new opportunities for improving the selectivity of ADCs.

Reference:


- *Discovery of Peptidomimetic Antibody-Drug Conjugate Linkers with Enhanced Protease Specificity*; B. Wei, J. Gunzner-Toste, H. Yao, T. Wang, J. Wang, Z. Xu, J. Chen, J. Wai, J. Nonomiya, S. P. Tsai, J. Chuh, K. R. Kozak, Y. Liu, S. F. Yu, J. Lau, G. Li, G. D. Phillips, D. Leipold, A. Kamath, D. Su, K. Xu, C. Eigenbrot, S. Steinbacher, R. Ohri, H. Raab, L. R. Staben, G. Zhao, J. A. Flygare, T. H. Pillow, V. Verma, L. A. Masterson, P. W. Howard, B. Safina; **J. Med. Chem.** 2018; **61**: 989-1000. <https://doi.org/10.1021/acs.jmedchem.7b01430>




3.1. Valine-Alanine-Based Enzymatically Cleavable Linkers

		Product details
ADC1290 6-Azidohexanoyl-Val-Ala-PAB 6-azidohexanoyl-valyl-alanyl-(4-aminobenzyl alcohol) CAS-No. 2706564-30-7 Formula C₂₁H₃₂N₆O₄ Mol. weight 432,52 g/mol		
ADC1300 6-Azidohexanoyl-Val-Ala-PAB-PNP 6-azidohexanoyl-valyl-alanyl-(4-aminobenzyl)-(4-nitrophenyl)-carbonate Formula C₂₈H₃₅N₇O₈ Mol. weight 597,62 g/mol		

Reference:


- *NKT cell-dependent glycolipid-peptide vaccines with potent anti-tumour activity*; R. J. Anderson, B. J. Compton, C. W. Tang, A. Authier-Hall, C. M. Hayman, G. W. Swinerd, R. Kowalczyk, P. Harris, M. A. Brimble, D. S. Larsen, O. Gasser, R. Weinkove, I. F. Hermans, G. F. Painter; **Chem. Sci.** 2015; **6**: 5120-5127. <https://doi.org/10.1039/c4sc03599b>

		Product details
ADC1330 Azido-PEG(4)-Val-Ala-PAB azido-tetraethyleneglycol-propanoyl-valyl-alanyl-(4-aminobenzyl alcohol) Formula C₂₆H₄₂N₆O₈ Mol. weight 566,65 g/mol		


		Product details
ADC1340	Azido-PEG(4)-Val-Ala-PAB-PNP	
azido-tetraethyleneglycol-propanoyl-valyl-alanyl-(4-aminobenzyl)-(4-nitrophenyl)-carbonate		
Formula	$C_{33}H_{45}N_7O_{12}$	
Mol. weight	731,75 g/mol	
ADC1310	4-Pentynoyl-Val-Ala-PAB	
4-pentynoyl-valyl-alanyl-(4-aminobenzyl alcohol)		
CAS-No.	1956294-75-9	
Formula	$C_{20}H_{27}N_3O_4$	
Mol. weight	373,45 g/mol	
ADC1320	4-Pentynoyl-Val-Ala-PAB-PNP	
4-pentynoyl-valyl-alanyl-(4-aminobenzyl)-(4-nitrophenyl)-carbonate		
CAS-No.	1956294-76-0	
Formula	$C_{27}H_{30}N_4O_8$	
Mol. weight	538,55 g/mol	

Reference:

- *Integrin-Targeting Knottin Peptide-Drug Conjugates Are Potent Inhibitors of Tumor Cell Proliferation*; N. Cox, J. R. Kintzing, M. Smith, G. A. Grant, J. R. Cochran; **Angew Chem Int Ed** 2016; **55**: 9894-7.
<https://doi.org/10.1002/anie.201603488>




		Product details
ADC1350	Alkyne-PEG(4)-Val-Ala-PAB	
propargyl-tetraethyleneglycol-propanoyl-valyl-alanyl-(4-aminobenzyl alcohol)		
CAS-No.	2348405-90-1	
Formula	C ₂₉ H ₄₅ N ₃ O ₉	
Mol. weight	579,68 g/mol	
		




		Product details
ADC1360 Alkyne-PEG(4)-Val-Ala-PAB-PNP		
propargyl-tetraethyleneglycol-propanoyl-valyl-alanyl-(4-aminobenzyl)-(4-nitrophenyl)-carbonate		
CAS-No.	2348405-91-2	
Formula	C ₃₆ H ₄₈ N ₄ O ₁₃	
Mol. weight	744,79 g/mol	







References:

- *Exploration of the carmaphycins as payloads in antibody drug conjugate anticancer agents*; J. Almaliti, B. Miller, H. Pietraszkiewicz, E. Glukhov, C. B. Naman, T. Kline, J. Hanson, X. Li, S. Zhou, F. A. Valeriote, W. H. Gerwick; **Eur J Med Chem** 2019; **161**: 416-432. <https://doi.org/10.1016/j.ejmech.2018.10.024>
- *Design and synthesis of novel dual-cyclic RGD peptides for alphavbeta3 integrin targeting*; J. Liu, X. Cheng, X. Tian, D. Guan, J. Ao, Z. Wu, W. Huang, Z. Le; **Bioorg Med Chem Lett** 2019; **29**: 896-900. <https://doi.org/10.1016/j.bmcl.2019.01.043>




		Product details	
ADC1040	Boc-Val-Ala-PAB		
t-Butyloxycarbonyl-valyl-alanyl-4-aminobenzylalcohol			
CAS-No.	1884577-99-4		
Formula	C ₂₀ H ₃₁ N ₃ O ₅		
Mol. weight	393,48 g/mol		
			
ADC1660		Boc-Val-Ala-PAB-Cl	
tert-butyloxycarbonyl-valyl-alanyl-4-aminobenzylchloride			
Formula	C ₂₀ H ₃₀ ClN ₃ O ₄		
Mol. weight	411,93 g/mol		
			
ADC1050		Boc-Val-Ala-PAB-PNP	
t-Butyloxycarbonyl-valyl-alanyl-(4-aminobenzyl)-(4-nitrophenyl)carbonate			
CAS-No.	1884578-00-0		
Formula	C ₂₇ H ₃₄ N ₄ O ₉		
Mol. weight	558,58 g/mol		
			



		Product details
ADC1060	Fmoc-Val-Ala-PAB	
9-Fluorenylmethyloxycarbonyl-valyl-alanyl-4-aminobenzylalcohol		
CAS-No.	1394238-91-5	 
Formula	C ₃₀ H ₃₃ N ₃ O ₅	
Mol. weight	515,61 g/mol	
ADC1670	Fmoc-Val-Ala-PAB-Cl	
9-Fluorenylmethyloxycarbonyl-valyl-alanyl-4-aminobenzylchloride		
CAS-No.	1491136-17-4	 
Formula	C ₃₀ H ₃₂ ClN ₃ O ₄	
Mol. weight	534,05 g/mol	
ADC1070	Fmoc-Val-Ala-PAB-PNP	
9-Fluorenylmethyloxycarbonyl-valyl-alanyl-(4-aminobenzyl)-(4-nitrophenyl)carbonate		
CAS-No.	1394238-92-6	 
Formula	C ₃₇ H ₃₆ N ₄ O ₉	
Mol. weight	680,71 g/mol	
ADC1410	Fmoc-Val-Ala-PAB-NMeCH₂CH₂NMe-Boc	
9-Fluorenylmethyloxycarbonyl-valyl-alanyl-4-aminobenzylloxycarbonyl-((t-butyl methyl(2-methylamino)ethyl)carbamate)		
CAS-No.	1691196-82-3	 
Formula	C ₄₀ H ₅₁ N ₅ O ₈	
Mol. weight	729,86 g/mol	

References:


- *Multivalency Increases the Binding Strength of RGD Peptidomimetic-Paclitaxel Conjugates to Integrin alphaV beta3*; A. Raposo Moreira Dias, A. Pina, A. Dal Corso, D. Arosio, L. Belvisi, L. Pignataro, M. Caruso and C. Gennari; **Chemistry** 2017; **23**: 14410-14415. <https://doi.org/10.1002/chem.201703093>
- *Synthesis and biological evaluation of RGD peptidomimetic-paclitaxel conjugates bearing lysosomally cleavable linkers*; A. Dal Corso, M. Caruso, L. Belvisi, D. Arosio, U. Piarulli, C. Albanese, F. Gasparri, A. Marsiglio, F. Sola, S. Troiani, B. Valsasina, L. Pignataro, D. Donati, C. Gennari; **Chemistry** 2015; **21**: 6921-9. <https://doi.org/10.1002/chem.201500158>
- *Elongated multiple electronic cascade and cyclization spacer systems in activatable anticancer prodrugs for enhanced drug release*; F. M. de Groot, W. J. Loos, R. Koekkoek, L. W. van Berkum, G. F. Busscher, A. E. Seelen, C. Albrecht, P. de Bruijn, H. W. Scheeren; **J Org Chem** 2001; **66**: 8815-30. <https://doi.org/10.1021/jo0158884>

			Product details
ADC1270	MC-Val-Ala-PAB		
6-maleimidohehexanoyl-valyl-alanyl-(4-aminobenzyl alcohol)			
CAS-No.	1870916-87-2		
Formula	C ₂₅ H ₃₄ N ₄ O ₆		
Mol. weight	486,56 g/mol		
ADC1700	MC-Val-Ala-PAB-Cl		
6-Maleimidohehexanoyl-valyl-alanyl-(4-aminobenzyl chloride)			
CAS-No.	2983182-03-0		
Formula	C ₂₅ H ₃₃ ClN ₄ O ₅		
Mol. weight	521,01 g/mol		
ADC1280	MC-Val-Ala-PAB-PNP		
6-maleimidohehexanoyl-valyl-alanyl-(4-aminobenzyl)-(4-nitrophenyl)-carbonate			
CAS-No.	1639939-40-4		
Formula	C ₃₂ H ₃₇ N ₅ O ₁₀		
Mol. weight	651,66 g/mol		

References:




- Improved Methodology for the Synthesis of a Cathepsin B Cleavable Dipeptide Linker, Widely Used in Antibody-Drug Conjugate Research; D. Mondal, J. Ford, K. G. Pinney; **Tetrahedron Lett** 2018; **59**: 3594-3599.
<https://doi.org/10.1016/j.tetlet.2018.08.021>
- Next generation maleimides enable the controlled assembly of antibody-drug conjugates via native disulfide bond bridging; F. F. Schumacher, J. P. Nunes, A. Maruani, V. Chudasama, M. E. Smith, K. A. Chester, J. R. Baker, S. Caddick; **Org Biomol Chem** 2014; **12**: 7261-9. <https://doi.org/10.1039/c4ob01550a>
- Site-Specific Conjugation of Auristatins onto Engineered scFv Using Second Generation Maleimide to Target HER2-positive Breast Cancer in Vitro; N. Aubrey, E. Allard-Vannier, C. Martin, F. Bryden, S. Letast, C. Colas, Z. Lakhri, N. Collinet, I. Dimier-Poisson, I. Chourpa, M. C. Viaud-Massuad, N. Joubert; **Bioconjug Chem** 2018; **29**: 3516-3521. <https://doi.org/10.1021/acs.bioconjchem.8b00668>
- Impact of cathepsin B-sensitive triggers and hydrophilic linkers on in vitro efficacy of novel site-specific antibody-drug conjugates; F. Bryden, C. Martin, S. Letast, E. Lles, I. Vieitez-Villemin, A. Rousseau, C. Colas, M. Brachet-Botineau, E. Allard-Vannier, C. Larbouret, M. C. Viaud-Massuad, N. Joubert; **Org Biomol Chem** 2018; **16**: 1882-1889. <https://doi.org/10.1039/c7ob02780j>




		Product details
ADC1080 Mal-Dap(Boc)-Val-Ala-PAB-PNP		
N-alpha-Maleimido-N-beta- <i>t</i> -butyloxycarbonyl-L-2,3-diaminopropionyl-valyl-alanyl-(4-aminobenzyl)-(4-nitrophenyl)-carbonate		
Formula	C ₃₄ H ₄₀ N ₆ O ₁₂	
Mol. weight	721,71 g/mol	



References:

- *Linker Technologies for Antibody-Drug Conjugates*; B. Nolting; **Antibody-Drug Conjugates** L. Ducry 2013; **1045**: 71-100. https://doi.org/10.1007/978-1-62703-541-5_5
- *Self-hydrolyzing maleimides improve the stability and pharmacological properties of antibody-drug conjugates*; R. P. Lyon, J. R. Setter, T. D. Bovee, S. O. Doronina, J. H. Hunter, M. E. Anderson, C. L. Balasubramanian, S. M. Duniho, C. I. Leiske, F. Li, P. D. Senter; **Nat Biotechnol** 2014; **32**: 1059-62. <https://doi.org/10.1038/nbt.2968>
- *Self-Stabilizing Linker Conjugate*; Lyon R., Doronina S., Bovee T.; Seattle Genetics, Inc.; U.S. Patent No. 9,504,756, 2013
- *In Vivo Antitumor Activity of a Novel Acetazolamide-Cryptophycin Conjugate for the Treatment of Renal Cell Carcinomas*; S. Cazzamalli, E. Figueras, L. Petho, A. Borbely, C. Steinkuhler, D. Neri, N. Sewald; **ACS Omega** 2018; **3**: 14726-14731. <https://doi.org/10.1021/acsomega.8b02350>


		Product details
ADC1390	Mal-beta-Ala-PEG(4)-Val-Ala-PAB	
maleimido-beta-alanyl-tetraethyleneglycol-propionyl-valyl-alanyl-(4-aminobenzyl alcohol)		
CAS-No.	2417003-93-9	
Formula	C ₃₃ H ₄₉ N ₅ O ₁₁	
Mol. weight	691,77 g/mol	
		
ADC1400	Mal-beta-Ala-PEG(4)-Val-Ala-PAB-PNP	
maleimido-beta-alanyl-tetraethyleneglycol-propionyl-valyl-alanyl-(4-aminobenzyl)-(4-nitrophenyl)-carbonate		
CAS-No.	2417003-94-0	
Formula	C ₄₀ H ₅₂ N ₆ O ₁₅	
Mol. weight	856,87 g/mol	
		
ADC1730	Mal-PhAc-Val-Ala-PAB	
(S)-2-(2-(4-(2,5-dioxo-2,5-dihydro-1H-pyrrol-1-yl)phenyl)acetamido)-N-((S)-1-((4-(hydroxymethyl)phenyl)amino)-1-oxopropan-2-yl)-3-methylbutanamide		
Formula	C ₂₇ H ₃₀ N ₄ O ₆	
Mol. weight	506,56 g/mol	
		



ADC1740 Mal-PhAc-Val-Ala-PAB-PNP

4-((S)-2-((S)-2-(2-(4-(2,5-dioxo-2,5-dihydro-1H-pyrrol-1-yl)phenyl)acetamido)-3-methylbutanamido)propanamido)benzyl (4-nitrophenyl) carbonate


Formula $C_{34}H_{33}N_5O_{10}$
Mol. weight 671,66 g/mol



ADC1770 Mal-PhAc-PEG(4)-Val-Ala-PAB

(S)-2-(3-(2-(2-(4-(2,5-dioxo-2,5-dihydro-1H-pyrrol-1-yl)phenyl)acetamido)ethoxy)propanamido)-N-((S)-1-((4-(hydroxymethyl)phenyl)amino)-1-oxopropan-2-yl)-3-methylbutanamide


Formula $C_{32}H_{39}N_5O_8$
Mol. weight 621,69 g/mol



ADC1780 Mal-PhAc-PEG(4)-Val-Ala-PAB-PNP

4-((2S,5S)-15-(4-(2,5-dioxo-2,5-dihydro-1H-pyrrol-1-yl)phenyl)-5-isopropyl-2-methyl-4,7,14-trioxo-10-oxa-3,6,13-triazapentadecanamido)benzyl (4-nitrophenyl) carbonate


Formula $C_{39}H_{42}N_6O_{12}$
Mol. weight 786,80 g/mol



ADC1580 Azido-cyclobutane-1,1-dicarboxamide-Ala-PAB

3-azidopropyl-cyclobutane-1,1-dicarboxamide-allyl-(4-aminobenzyl alcohol)


CAS-No. 2576471-45-7
Formula $C_{19}H_{26}N_6O_4$
Mol. weight 402,45 g/mol



ADC1590 Azido-cyclobutane-1,1-dicarboxamide-Ala-PAB-PNP

3-azidopropyl-cyclobutane-1,1-dicarboxamide-allyl-(4-aminobenzyl)-(4-nitrophenyl)-carbonate


CAS-No. 2576471-36-6
Formula $C_{26}H_{29}N_7O_8$
Mol. weight 567,55 g/mol








ADC1600 Propargyl-cyclobutane-1,1-dicarboxamide-Ala-PAB

propargyl-cyclobutane-1,1-dicarboxamide-allyl-(4-aminobenzyl alcohol)

CAS-No. 2576471-28-6
Formula $C_{19}H_{23}N_3O_4$
Mol. weight 357,40 g/mol








		Product details
ADC1610	Propargyl-cyclobutane-1,1-dicarboxamide-Ala-PAB-PNP	
propargyl-cyclobutane-1,1-dicarboxamide-ala-nyl-(4-aminobenzyl)-(4-nitrophenyl)-carbonate		
CAS-No.	2576471-42-4	
Formula	C ₂₆ H ₂₆ N ₄ O ₈	
Mol. weight	522,51 g/mol	
		
ADC1620	DBCO-cyclobutane-1,1-dicarboxamide-Ala-PAB	
dibenzoazacyclooctyne-cyclobutane-1,1-dicarboxami- de-alanyl-(4-aminobenzyl alcohol)		
CAS-No.	2576471-46-8	
Formula	C ₃₄ H ₃₄ N ₄ O ₅	
Mol. weight	578,66 g/mol	
		
ADC1630	DBCO-cyclobutane-1,1-dicarboxamide-Ala-PAB-PNP	
dibenzoazacyclooctyne-cyclobutane-1,1-dicarboxami- de-alanyl-(4-aminobenzyl)-(4-nitrophenyl)-carbonate		
CAS-No.	2576471-43-5	
Formula	C ₄₁ H ₃₇ N ₅ O ₉	
Mol. weight	743,76 g/mol	
		
ADC1560	Mal-cyclobutane-1,1-dicarboxamide-Ala-PAB	
5-maleimidopentyl-cyclobutane-1,1-dicarboxami- de-alanyl-(4-aminobenzyl alcohol)		
CAS-No.	2576471-41-3	
Formula	C ₂₅ H ₃₂ N ₄ O ₆	
Mol. weight	484,54 g/mol	
		
ADC1570	Mal-cyclobutane-1,1-dicarboxamide-Ala-PAB-PNP	
5-maleimidopentyl-cyclobutane-1,1-dicarboxami- de-alanyl-(4-aminobenzyl)-(4-nitrophenyl)-carbonate		
CAS-No.	2576471-35-5	
Formula	C ₃₂ H ₃₅ N ₅ O ₁₀	
Mol. weight	649,65 g/mol	
		

Reference:

- *Discovery of Peptidomimetic Antibody-Drug Conjugate Linkers with Enhanced Protease Specificity*; B. Wei, J. Gunzner-Toste, H. Yao, T. Wang, J. Wang, Z. Xu, J. Chen, J. Wai, J. Nonomiya, S. P. Tsai, J. Chuh, K. R. Kozak, Y. Liu, S. F. Yu, J. Lau, G. Li, G. D. Phillips, D. Leipold, A. Kamath, D. Su, K. Xu, C. Eigenbrot, S. Steinbacher, R. Ohri, H. Raab, L. R. Staben, G. Zhao, J. A. Flygare, T. H. Pillow, V. Verma, L. A. Masterson, P. W. Howard, B. Safina; **J. Med. Chem.** 2018; **61**: 989-1000. <https://doi.org/10.1021/acs.jmedchem.7b01430>





[↑ back to content](#)

3.2. Valine-Citrulline-Based Enzymatically Cleavable Linkers

		Product details
ADC1090	Mal-Dap(Boc)-Val-Cit-PAB-PNP	
N-alpha-Maleimido-N-beta-t-butyloxycarbonyl-L-2,3-diaminopropionyl-valyl-citrullyl-(4-aminobenzyl)-(4-nitrophenyl)-carbonate		
CAS-No.	2281797-57-5	
Formula	C ₃₇ H ₄₆ N ₈ O ₁₃	
Mol. weight	810,81 g/mol	
ADC1750	Mal-PhAc-Val-Cit-PAB	
(S)-2-((S)-2-(2-(4-(2,5-dioxo-2,5-dihydro-1H-pyrrol-1-yl)phenyl)acetamido)-3-methylbutanamido)-N-(4-(hydroxymethyl)phenyl)-5-ureidopentanamide		
Formula	C ₃₀ H ₃₆ N ₆ O ₇	
Mol. weight	592,65 g/mol	
ADC1760	Mal-PhAc-Val-Cit-PAB-PNP	
4-((S)-2-((S)-2-(2-(4-(2,5-dioxo-2,5-dihydro-1H-pyrrol-1-yl)phenyl)acetamido)-3-methylbutanamido)-5-ureidopentanamido)benzyl (4-nitrophenyl) carbonate		
CAS-No.	3034752-26-3	
Formula	C ₃₇ H ₃₉ N ₇ O ₁₁	
Mol. weight	757,76 g/mol	
ADC1120	6-Azidohexanoyl-Val-Cit-PAB	
6-azidohexanoyl-valyl-citrullyl-(4-aminobenzyl alcohol)		
CAS-No.	1613321-02-0	
Formula	C ₂₄ H ₃₈ N ₈ O ₅	
Mol. weight	518,61 g/mol	
ADC1130	6-Azidohexanoyl-Val-Cit-PAB-PNP	
6-azidohexanoyl-valyl-citrullyl-(4-aminobenzyl)-(4-nitrophenyl)-carbonate		
CAS-No.	1613321-01-9	
Formula	C ₃₁ H ₄₁ N ₉ O ₉	
Mol. weight	683,71 g/mol	



Reference:

- *NKT cell-dependent glycolipid-peptide vaccines with potent anti-tumour activity*; R. J. Anderson, B. J. Compton, C. W. Tang, A. Authier-Hall, C. M. Hayman, G. W. Swinerd, R. Kowalczyk, P. Harris, M. A. Brimble, D. S. Larsen, O. Gasser, R. Weinkove, I. F. Hermans, G. F. Painter; **Chem. Sci.** 2015; **6**: 5120-5127. <https://doi.org/10.1039/c4sc03599b>

		Product details
ADC1160 Azido-PEG(4)-Val-Cit-PAB azido-tetraethyleneglycol-propanoyl-valyl-citrullyl-(4-aminobenzyl alcohol) CAS-No. 2055024-64-9 Formula C ₂₉ H ₄₈ N ₈ O ₉ Mol. weight 652,74 g/mol		
ADC1170 Azido-PEG(4)-Val-Cit-PAB-PNP azido-tetraethyleneglycol-propanoyl-valyl-citrullyl-(4-aminobenzyl)-(4-nitrophenyl)-carbonate CAS-No. 1869126-60-2 Formula C ₃₆ H ₅₁ N ₉ O ₁₃ Mol. weight 817,84 g/mol		
ADC1140 4-Pentynoyl-Val-Cit-PAB 4-pentynoyl-valyl-citrullyl-(4-aminobenzyl alcohol) CAS-No. 2708150-97-2 Formula C ₂₃ H ₃₃ N ₅ O ₅ Mol. weight 459,54 g/mol		
ADC1150 4-Pentynoyl-Val-Cit-PAB-PNP 4-pentynoyl-valyl-citrullyl-(4-aminobenzyl)-(4-nitrophenyl)-carbonate Formula C ₃₀ H ₃₆ N ₆ O ₉ Mol. weight 624,64 g/mol		



Reference:




- *Integrin-Targeting Knottin Peptide-Drug Conjugates Are Potent Inhibitors of Tumor Cell Proliferation*. N. Cox, J. R. Kintzing, M. Smith, G. A. Grant, J. R. Cochran; **Angew. Chem. Int. Ed.** 2016; **55(34)**: 9894-9897. <https://doi.org/10.1002/anie.201603488>

		Product details
ADC1180	Alkyne-PEG(5)-Val-Cit-PAB propargyl-tetraethyleneglycol-propanoyl-valyl-citrullyl-(4-aminobenzyl alcohol)	 
Formula	C ₃₂ H ₅₁ N ₅ O ₁₀	
Mol. weight	665,77 g/mol	
ADC1190	Alkyne-PEG(5)-Val-Cit-PAB-PNP propargyl-tetraethyleneglycol-propanoyl-valyl-citrullyl-(4-aminobenzyl)-(4-nitrophenyl)-carbonate	 
Formula	C ₃₉ H ₅₄ N ₆ O ₁₄	
Mol. weight	830,88 g/mol	

References:


- Exploration of the carmaphycins as payloads in antibody drug conjugate anticancer agents. J. Almaliti, B. Miller, H. Pietraszkiewicz, E. Glukhov, C. B. Naman, T. Kline, J. Hanson, X. Li, S. Zhou, F. A. Valeriote, W. H. Gerwick; *Eur J Med Chem.* 2019; **161**: 416-432. <https://doi.org/10.1016/j.ejmech.2018.10.024>
- Design and synthesis of novel dual-cyclic RGD peptides for αvβ3 integrin targeting. J. Liu, X. Cheng, X. Tian, D. Guan, J. Ao, Z. Wu, W. Huang, Z. Le; *Bioorg Med Chem Lett.* 2019; **29(7)**: 896-900. <https://doi.org/10.1016/j.bmcl.2019.01.043>

		Product details
ADC1020	Boc-Val-Cit-PAB t-Butyloxycarbonyl-valyl-citrullyl-4-aminobenzylalcohol	 
CAS-No.	870487-09-5	
Formula	C ₂₃ H ₃₇ N ₃ O ₆	
Mol. weight	479,59 g/mol	
ADC1010	Boc-Val-Cit-PAB-PNP t-Butyloxycarbonyl-valyl-citrullyl-(4-aminobenzyl)-(4-nitrophenyl)carbonate	 
CAS-No.	870487-10-8	
Formula	C ₃₀ H ₄₀ N ₆ O ₁₀	
Mol. weight	644,67 g/mol	

		Product details
ADC1030	Fmoc-Val-Cit-PAB	
9-Fluorenylmethyloxycarbonyl-valyl-citrullyl-4-aminobenzylalcohol		
CAS-No.	159858-22-7	
Formula	C ₃₃ H ₃₉ N ₃ O ₆	
Mol. weight	601,29 g/mol	
		
ADC1000	Fmoc-Val-Cit-PAB-PNP	
9-Fluorenylmethyloxycarbonyl-valyl-citrullyl-(4-aminobenzyl)-(4-nitrophenyl)carbonate		
CAS-No.	863971-53-3	
Formula	C ₄₀ H ₄₂ N ₅ O ₁₀	
Mol. weight	766,80 g/mol	
		
ADC1240	Fmoc-Val-Cit-PAB-NMeCH₂CH₂NMe-Boc	
9-Fluorenylmethyloxycarbonyl-valyl-citrullyl-4-aminobenzylloxycarbonyl-((t-butyloxymethyl(2-methylamino)ethyl)carbamate)		
CAS-No.	1802297-96-6	
Formula	C ₄₃ H ₅₇ N ₇ O ₉	
Mol. weight	815,95 g/mol	
		

References:

- *Multivalency Increases the Binding Strength of RGD Peptidomimetic-Paclitaxel Conjugates to Integrin αVβ3.* A. R. M. Dias, A. Pina, A. Dal Corso, D. Arosio, L. Belvisi, L. Pignataro, M. Caruso, C. Gennari; **Chem. Eur. J.** 2017; **23(58)**: 14410-14415. <https://doi.org/10.1002/chem.201703093>
- *Synthesis and Biological Evaluation of RGD Peptidomimetic-Paclitaxel Conjugates Bearing Lysosomally Cleavable Linkers.* A. D. Corso, M. Caruso, L. Belvisi, D. Arosio, U. Piarulli, C. Albanese, F. Gasparri, A. Marsiglio, F. Sola, S. Troiani, B. Valsasina, L. Pignataro, D. Donati, C. Gennari; **Chem. Eur. J.** 2015; **21(18)**: 6921-6929. <https://doi.org/10.1002/chem.201500158>
- *Elongated Multiple Electronic Cascade and Cyclization Spacer Systems in Activatable Anticancer Prodrugs for Enhanced Drug Release.* F. M. H. de Groot, W. J. Loos, R. Koekkoek, L. W. A. van Berkom, G. F. Busscher, A. E. Seelen, C. Albrecht, P. Bruijn, H. W. Scheeren; **J. Org. Chem.** 2001; **66(26)**: 8815-8830. <https://doi.org/10.1021/jo0158884>

		Product details
ADC1100	MC-Val-Cit-PAB	
6-maleimidohexanoyl-valyl-citrullyl-(4-aminobenzyl alcohol)		
CAS-No.	159857-80-4	
Formula	C ₂₈ H ₄₀ N ₆ O ₇	
Mol. weight	572,65 g/mol	
		

[↑ back to content](#)

The Concept of Antibody-Drug Conjugation (ADC)

Permanent Linkers


Cleavable Linkers

Trifunctional Linkers

Cross-Linkers for other Bio Applications




Preparing Carriers for Conjugation







Index

		Product details
ADC1110	MC-Val-Cit-PAB-PNP	
6-maleimido-hexanoyl-valyl-citrullyl-(4-aminobenzyl)-(4-nitrophenyl)-carbonate		
CAS-No.	159857-81-5	
Formula	C ₃₅ H ₄₃ N ₇ O ₁₁	
Mol. weight	737,76 g/mol	
		

Reference:

- *Improved Methodology for the Synthesis of a Cathepsin B Cleavable Dipeptide Linker, Widely Used in Antibody-Drug Conjugate Research.* D. Mondal, J. Ford, K. G. Pinney; **Tetrahedron Lett.** 2018; **59(40)**: 3594-3599.
<https://doi.org/10.1016/j.tetlet.2018.08.021>

		Product details
ADC1220	Mal-beta-Ala-PEG(4)-Val-Cit-PAB	
maleimido-beta-alanyl-tetraethyleneglycol-propionyl-valyl-citrullyl-(4-aminobenzyl alcohol)		
CAS-No.	1949793-41-2	
Formula	C ₃₆ H ₅₅ N ₇ O ₁₂	
Mol. weight	777,86 g/mol	
		
ADC1230	Mal-beta-Ala-PEG(4)-Val-Cit-PAB-PNP	
maleimido-beta-alanyl-tetraethyleneglycol-propionyl-valyl-citrullyl-(4-aminobenzyl)-(4-nitrophenyl)-carbonate		
CAS-No.	2003260-12-4	
Formula	C ₄₃ H ₅₈ N ₈ O ₁₆	
Mol. weight	942,96 g/mol	
		
ADC1480	Azido-cyclobutane-1,1-dicarboxamide-Cit-PAB	
3-azidopropyl-cyclobutane-1,1-dicarboxamide-citrullyl-(4-aminobenzyl alcohol)		
CAS-No.	2576471-33-3	
Formula	C ₂₂ H ₃₂ N ₆ O ₅	
Mol. weight	488,54 g/mol	
		

		Product details
ADC1490	Azido-cyclobutane-1,1-dicarboxamide-Cit-PAB-PNP 3-azidopropyl-cyclobutane-1,1-dicarboxamide-citrullyl-(4-aminobenzyl)-(4-nitrophenyl)-carbonate CAS-No. 2576471-44-6 Formula $C_{29}H_{35}N_9O_9$ Mol. weight 653,64 g/mol	 
ADC1500	Propargyl-cyclobutane-1,1-dicarboxamide-Cit-PAB propargyl-cyclobutane-1,1-dicarboxamide-citrullyl-(4-aminobenzyl alcohol) CAS-No. 2576471-27-5 Formula $C_{22}H_{29}N_5O_5$ Mol. weight 443,50 g/mol	 
ADC1510	Propargyl-cyclobutane-1,1-dicarboxamide-Cit-PAB-PNP propargyl-cyclobutane-1,1-dicarboxamide-citrullyl-(4-aminobenzyl)-(4-nitrophenyl)-carbonate CAS-No. 2576471-40-2 Formula $C_{29}H_{32}N_6O_9$ Mol. weight 608,60 g/mol	 
ADC1520	DBCO-cyclobutane-1,1-dicarboxamide-Cit-PAB dibenzoazacyclooctyne-cyclobutane-1,1-dicarboxamide-citrullyl-(4-aminobenzyl alcohol) CAS-No. 2576471-51-5 Formula $C_{37}H_{40}N_6O_6$ Mol. weight 664,75 g/mol	 
ADC1530	DBCO-cyclobutane-1,1-dicarboxamide-Cit-PAB-PNP dibenzoazacyclooctyne-cyclobutane-1,1-dicarboxamide-citrullyl-(4-aminobenzyl)-(4-nitrophenyl)-carbonate CAS-No. 2576471-34-4 Formula $C_{44}H_{43}N_7O_{10}$ Mol. weight 829,85 g/mol	 
ADC1460	Mal-cyclobutane-1,1-dicarboxamide-Cit-PAB 5-maleimidopentyl-cyclobutane-1,1-dicarboxamide-citrullyl-(4-aminobenzyl alcohol) CAS-No. 1799663-03-8 Formula $C_{28}H_{38}N_6O_7$ Mol. weight 570,64 g/mol	 

The Concept of Antibody-Drug Conjugation (ADC)

Permanent Linkers

Cleavable Linkers

Trifunctional Linkers

Cross-Linkers for other Bio Applications

Preparing Carriers for Conjugation

Index

[↑ back to content](#)


ADC1470 Mal-cyclobutane-1,1-dicarboxamide-Cit-PAB-PNP

5-maleimidopentyl-cyclobutane-1,1-dicarboxamide-citrullyl-(4-aminobenzyl)-(4-nitrophenyl)-carbonate

CAS-No. 2204228-34-0

Formula C₃₅H₄₁N₇O₁₁

Mol. weight 735,74 g/mol



Reference:

- *Discovery of Peptidomimetic Antibody–Drug Conjugate Linkers with Enhanced Protease Specificity.* B. Wei, J. Gunzner-Toste, H. Yao, T. Wang, J. Wang, Z. Xu, J. Chen, J. Wai, J. Nonomiya, S. Ping Tsai, J. Chuh, K. R. Kozak, Y. Liu, S. Yu, J. Lau, G. Li, G. D. Phillips, D. Leipold, A. Kamath, D. Su, K. Xu, C. Eigenbrot, S. Steinbacher, R. Ohri, H. Raab, L. R. Staben, G. Zhao, J. A. Flygare, T. H. Pillow, V. Verma, L. A. Masterson, P. W. Howard, B. Safina; **J. Med. Chem.** 2018; **61**(3): 989-1000. <https://doi.org/10.1021/acs.jmedchem.7b01430>



3.3. β-Glucuronide Enzymatically Cleavable Linkers

As an extension of the linkerology® toolbox, the design of linkers with improved stability during systemic circulation is highly desired. As the drug-releasing lysosomal enzyme β-glucuronidase is abundantly present within lysosomes and overexpressed in some tumor types but low outside cells, β-glucuronic acid-based linkers provide the potential for high ADC stability in the systemic circulation and selective intracellular drug release. Especially for ADCs based on highly hydrophobic drugs, the incorporation of the highly hydrophilic β-glucuronides may circumvent the tendency of aggregation. For example, a drug-linker consisting of a β-glucuronide linked to auristatin MMAF was prepared. Rat plasma stability analysis revealed an extrapolated half-life of 81 days, compared with about six days for the corresponding valine-citrulline dipeptide-linked MMAF.



Interested in β-Glucuronide Enzymatically Cleavable Linkers?

Please contact our Custom Synthesis for more details!



References:


- *Expanded Utility of the β -Glucuronide Linker: ADCs That Deliver Phenolic Cytotoxic Agents*; S. C. Jeffrey, J. De Brabander, J. Miyamoto, P. D. Senter; **ACS Med. Chem. Lett.** 2010; **1**: 277-280. <https://doi.org/10.1021/ml100039h>
- *Development and Properties of β -Glucuronide Linkers for Monoclonal Antibody-Drug Conjugates*; S. C. Jeffrey, J. B. Andreyka, S. X. Bernhardt, K. M. Kissler, T. Kline, J. S. Lenox, R. F. Moser, M. T. Nguyen, N. M. Okeley, I. J. Stone, X. Zhang, P. D. Senter; **Bioconjugate Chem.** 2006; **17**: 831-840. <https://doi.org/10.1021/bc0600214>
- *Linker Technologies for Antibody-Drug Conjugates*; B. Nolting; **Antibody-Drug Conjugates L. Ducry** 2013; **1045**: 71-100. https://doi.org/10.1007/978-1-62703-541-5_5

3.4. PH-Responsive Linkers

Various linker types for conjugation are available on the market, responding to and being – ideally tracelessly – cleaved by virtue of a certain trigger, e.g., peptidic linkers are cleaved in the presence of specific enzymes.

Thus, conditions solely present in the affected tissue need to be selected to avoid premature drug release. Compared to healthy cells, abnormal metabolism and proliferation in tumor cells lead to a lowered intracellular pH of around 5.0-6.0, while a lower pH (4.0-5.0) is reported for the lysosomes.

At Iris, we offer a pH-sensitive linker based on an alkyne-functionalized, *para*-nitrophenyl (PNP)-modified 5-(hydroxymethyl)-pyrogallol-orthoester (HMPO).



This construct is stable during circulation in plasma (tested for 24 hours at pH values of 7.4 and 6.6.) while being cleaved tracelessly via 1,6-elimination at pH 5.5. Payloads can be easily coupled: The HMPO-PNP linker ADC1800 reacts with base-activated amides while releasing the nitrophenyl residue.




Fig. 20: Coupling of payload and linker. Reaction conditions: The drug (depicted in red) is dissolved in THF at -78°C. The strong base potassium hexamethyldisilylride (KHMDS) and then [ADC1800 on page 74](#) are added. For the workup, the temperature is allowed to rise to -5 °C over 4 hours. Ammonium chloride is added to quench excess linker, then the product is extracted and purified by chromatography.

While the payload is attached as carbamate, the connection to the carrier protein (usually a tumor-specific antibody) is realized via the propargyl residue through biorthogonal CuAAC click reaction, forming a 1,2,3-triazole. For this, the antibody must carry azido groups, which may be introduced via recombinant biosynthesis, utilizing pyrrolysyl-tRNA synthetase and a suitably modified lysine (e.g., <https://www.iris-biotech.de/HAA2080>), or with an amino-reactive azido linker like, e.g., [PEG1400 on page 33](#) or [RL-2980 on page 37](#).




Fig. 21: Combination of antibody and drug-linker adduct to generate the antibody-drug-conjugate (ADC).

At acidic pH, the linker gets protonated inducing a 1,6-elimination reaction. The oxocarbonyl moiety is released as carbon dioxide, leading to the irreversible and traceless linker cleavage and payload release.






Fig. 22: Payload release at acidic pH in the tumor microenvironment.

While [ADC1800 on page 74](#) as the activated form of the linker reacts directly with the amide-containing drug molecule, we also offer the precursor [ADC1790 on page 74](#) (Alkyne-HMPO-OH) which can be further customized with your linker of choice.


		Product details
ADC1800 Alkyne-HMPO-PNP (2-ethoxy-7-(prop-2-yn-1-yloxy)benzo[d][1,3]dioxol-5-yl) methyl (4-nitrophenyl) carbonate CAS-No. 3028213-53-5 Formula $C_{20}H_{17}NO_9$ Mol. weight 415,35 g/mol		
ADC1790 Alkyne-HMPO-OH (2-ethoxy-7-(prop-2-yn-1-yloxy)benzo[d][1,3]dioxol-5-yl) methanol CAS-No. 3028213-63-7 Formula $C_{13}H_{14}O_5$ Mol. weight 250,25 g/mol		

References:

- A pH-responsive crosslinker platform for antibody-drug conjugate (ADC) targeting delivery; F. Migliorini, E. Cini, E. Dreassi, F. Finetti, G. Ievoli, G. Macri, E. Petricci, E. Rango, L. Trabalzini, M. Taddei; **Chem. Commun.** 2022; **58(75)**: 10532-10535. <https://doi.org/10.1039/D2CC03052G>
- A Self-Immolative Linker for the pH-Responsive Release of Amides; A. Petrini, G. Ievoli, F. Migliorini, M. Taddei, S. Siciliano; **Molecules** 2023; **28(6)**: 2445. <https://doi.org/10.3390/molecules28062445>
- Conjugates of PSMA-binding moieties with cytotoxic agents; WO2024028258
- Compositions and methods for treatment of sexual dysfunction and related diseases, disorders, and conditions; WO2022251699


3.5. Disulfide-Based (Self-Immolative) Linkers

Another chemically labile linkage extensively exploited in the development of antibody-drug conjugates are disulfides. They are stable at physiological pH and are designed to release the drug upon internalization inside cells. The cytosol provides a significantly more reducing environment compared to the extracellular milieu and the presence of cytoplasmic thiol cofactor, such as reduced glutathione (GSH). Additionally, the intracellular enzyme protein disulfide isomerase, or similar enzymes capable of cleaving disulfide bonds, may also contribute to the preferential cleavage of disulfide bonds inside cells. GSH is reported to be present in cells in the concentration range of 0.5-10 mM, compared with a significantly lower concentration of GSH or cysteine in plasma at approximately 5 μ M. This is especially true for tumor cells, where irregular blood flow leads to a hypoxic state, resulting in enhanced activity of reductive enzymes and therefore in even higher glutathione concentrations.






The stability of disulfide bridges can be fine-tuned by adjacent residues (Fig. 23). Methyl groups are bulky enough to have a significant influence on the thermodynamic stability of the disulfide bridge. While one additional methyl group already enhances the stability drastically, two methyl groups make the disulfide bond practically stable towards reductive cleavage. A methylation number of three or four will completely lock the disulfide bridge towards further modifications. As the direct conjugation of cleavable triggers to bioactive agents through disulfide bridges suffers from ineffective cleavage in case of bulky moieties and resulting steric hindrance as well as restricted possibilities for trigger-drug combinations, disulfide based self-immolative linkers (DSILs) provide a robust strategy for selective activation upon disulfide cleavage in the reductive cytoplasmic milieu.

Disulfide-based self-immolative linkers benefit of the reversibility of disulfide-bond formation. Upon oxidation, free thiols form less nucleophilic disulfide bonds, preventing self-immolative fragmentation. However, this process can be reversed in the presence of reducing agents, such as GSH. Those specifications allow for sufficient stability in the extracellular milieu but spontaneous self-immolative reaction within the cytosol upon GSH-mediated disulfide cleavage. Variations in the linker's chemical composition (disulfide ethoxycarbonyl (SSE) vs. disulfide benzyloxycarbonyl (SSB)) result in chemically tunable kinetics of the self-immolative cleavage due to different response rates towards GSH, showing higher rates for SSB-based DSILs compared to SSE-based ones (Fig. 24). Thus, the choice of the linker allows for fine-tuning of the cleavage speed and payload release.



For ease of synthesis, Iris Biotech offers pyridyl disulfides as building blocks for the preparation of disulfide-based self-immolative linkers. Pyridyl disulfides undergo a disulfide exchange reaction with thiol groups to form disulfide bonds over a broad pH range also suitable for physiological pH. During the reaction, a disulfide exchange occurs between the biomolecule's thiol group and the reagent's 2-pyridyl-dithiol group. As a result, pyridine-2-thione is released, which can be followed spectrophotometrically ($\lambda_{\text{max}} = 343 \text{ nm}$) to monitor the progress of the reaction.

			Product details
HNN1090	N₃-Cystamine*HCl		
Azido-cystamine hydrochloride			
CAS-No.	1807512-40-8 net		
Formula	C ₄ H ₁₀ N ₄ S ₂ *HCl		
Mol. weight	178,28*36,45 g/mol		
BNN1170	Boc-Cystamine		
2-(t-Butyloxycarbonylamino)ethyldithio-2'-ethylamine			
CAS-No.	485800-26-8		
Formula	C ₉ H ₂₀ N ₂ O ₂ S ₂		
Mol. weight	252,40 g/mol		
BNN1063	Boc-Cystamine*HCl		
2-(t-Butyloxycarbonylamino)ethyldithio-2'-ethylamine hydrochloride			
CAS-No.	93790-49-9		
Formula	C ₉ H ₂₀ N ₂ O ₂ S ₂ *HCl		
Mol. weight	252,40*36,45 g/mol		

[↑ back to content](#)


BNN1360 Di-Boc-Cystamine

N,N'-Bis-*tert*-butoxycarbonyl-cystamine

CAS-No. 67385-10-8

Formula $C_{14}H_{28}N_2O_4S_2$

Mol. weight 352,51 g/mol




BAA2180 Boc-Cystamine-Suc-OH

4-((2-*t*-Butoxycarbonylaminoethyl)disulfanyl)ethylamino)-4-oxobutanoic acid

CAS-No. 946849-79-2

Formula $C_{13}H_{24}N_2O_5S_2$

Mol. weight 352,47 g/mol




RL-3370 Fmoc-Cystamine*HCl

2-((9-Fluorenylmethoxycarbonylamino)ethyl)disulfanyl-(2-aminoethane) hydrochloride

CAS-No. 2893917-85-4

Formula $C_{19}H_{22}N_2O_2S_2 \cdot HCl$

Mol. weight 374,52*36,45 g/mol




RL-3310 Fmoc-Cystamine-Suc

4-((2-((9-Fluorenylmethoxycarbonylamino)ethyl)disulfanyl)ethylamino)-4-oxobutanoic acid

CAS-No. 946849-80-5

Formula $C_{23}H_{26}N_2O_5S_2$

Mol. weight 474,59 g/mol




RL-4100 Azido-SS-COOH

3-((2-azidoethyl)disulfanyl)propanoic acid

CAS-No. 2228857-32-5

Formula $C_5H_9N_3O_2S_2$

Mol. weight 207,27 g/mol




RL-3320 Azido-Pen-SS-COOH







3-((2-(5-azidopentanamido)ethyl)disulfanyl)propanoic acid

CAS-No. 2576471-47-9

Formula $C_{10}H_{18}N_4O_3S_2$

Mol. weight 306,40 g/mol



			Product details
RL-3330	Alkyne-SS-COOH		
3-((2-pent-4-ynamidoethyl)disulfanyl)propanoic acid			
CAS-No.	2279938-29-1		
Formula	C ₁₀ H ₁₅ NO ₃ S ₂		
Mol. weight	261,36 g/mol		
RL-4090	Mal-SS-COOH		
3-((2-(3-(2,5-dioxo-2,5-dihydro-1H-pyrrol-1-yl)propanamido)ethyl)disulfanyl)propanoic acid			
CAS-No.	2128735-24-8		
Formula	C ₁₂ H ₁₆ N ₂ O ₅ S ₂		
Mol. weight	332,39 g/mol		
RL-4110	DBCO-Suc-SS-COOH		
CAS-No.	2749426-25-1		
Formula	C ₂₄ H ₂₄ N ₂ O ₄ S ₂		
Mol. weight	468,59 g/mol		
RL-2190	Boc-SS-COOH		
3-((2-(tert-butoxycarbonylamino)ethyl)disulfanyl)propanoic acid			
CAS-No.	485800-27-9		
Formula	C ₁₀ H ₁₉ NO ₄ S ₂		
Mol. weight	281,39 g/mol		
RL-2200	Fmoc-SS-COOH		
3-((2-(((9H-fluoren-9-yl)methoxy)carbonylamino)ethyl)disulfanyl)propanoic acid			
CAS-No.	864235-83-6		
Formula	C ₂₆ H ₂₁ NO ₄ S ₂		
Mol. weight	403,52 g/mol		
RL-2810	Boc-AEDI-OH		
2-((2-(tert-Butyloxycarbonylamino)ethyl)disulfanyl)-2-methylpropanoic acid			
CAS-No.	144700-78-7		
Formula	C ₁₁ H ₂₁ NO ₄ S ₂		
Mol. weight	295,42 g/mol		

The Concept of Antibody-Drug Conjugation (ADC)

Permanent Linkers

Cleavable Linkers

Trifunctional Linkers

Cross-Linkers for other Bio Applications

Preparing Carriers for Conjugation

Index

[↑ back to content](#)


RL-2800 Fmoc-AEDI-OH

2-((2-((9-Fluorenylmethyloxycarbonyl)amino)ethyl)disulfanyl)-2-methylpropanoic acid

CAS-No. 1823244-38-7

Formula $C_{21}H_{23}NO_4S_2$

Mol. weight 417,54 g/mol




RL-3510 Boc-NH-SS-OH

2-((2-(*t*-Butyloxycarbonylamino)ethyl)disulfaneyl)ethan-1-ol

CAS-No. 877864-07-8

Formula $C_9H_{19}NO_3S_2$

Mol. weight 253,38 g/mol




RL-3520 Boc-NH-SS-OpNC

2-((2-(*t*-Butyloxycarbonylamino)ethyl)disulfaneyl)ethan-1-yl *p*-nitrophenylcarbonate

CAS-No. 2040301-00-4

Formula $C_{16}H_{22}N_2O_5S_2$

Mol. weight 418,48 g/mol




RL-3530 Fmoc-NH-SS-OH

2-((2-((9-Fluorenylmethyloxycarbonyl)amino)ethyl)disulfaneyl)ethan-1-ol

CAS-No. 2576471-39-9

Formula $C_{19}H_{21}NO_3S_2$

Mol. weight 375,50 g/mol




RL-4150 Azido-SS-OpNC

2-((2-(azidoethyl)disulfanyl)ethyl (4-nitrophenyl) carbonate

CAS-No. 2766027-28-3

Formula $C_{11}H_{12}N_4O_5S_2$

Mol. weight 344,36 g/mol




RL-3540 Fmoc-NH-SS-OpNC







2-((2-((9-Fluorenylmethyloxycarbonyl)amino)ethyl)disulfaneyl)ethan-1-yl *p*-nitrophenylcarbonate

CAS-No. 2576471-53-7

Formula $C_{26}H_{24}N_2O_7S_2$

Mol. weight 540,61 g/mol



			Product details
RL-4160	pNCO-SS-OpNC		
disulfanediy(bis(ethane-2,1-diyl) bis(4-nitrophenyl) bis(carbonate)) CAS-No. 1435972-52-3 Formula $C_{18}H_{16}N_2O_{10}S_2$ Mol. weight 484,45 g/mol			
RL-3560	Boc-NH-SS-Bzl-OH		
4-((2-(t-Butyloxycarbonylamino)ethyl)disulfaneyl) benzylalcohol CAS-No. 2576471-54-8 Formula $C_{16}H_{21}NO_3S_2$ Mol. weight 315,45 g/mol			
RL-3570	Boc-NH-SS-Bzl-OpNC		
4-((2-(t-Butyloxycarbonylamino)ethyl)disulfaneyl) benzyl p-nitrophenylcarbonate CAS-No. 2576471-38-8 Formula $C_{21}H_{24}N_2O_7S_2$ Mol. weight 480,55 g/mol			
RL-3580	Fmoc-NH-SS-Bzl-OH		
4-((2-((9-Fluorenylmethyloxycarbonyl)amino)ethyl)disulfaneyl)benzylalcohol CAS-No. 2064282-26-2 Formula $C_{24}H_{23}NO_3S_2$ Mol. weight 437,57 g/mol			
RL-3300	Biotin-SS-COOH		
3-((2-(5-((3aS,4S,6aR)-2-oxohexahydro-1H-thieno[3,4-d]imidazol-4-yl)pentanamido)ethyl)disulfanyl)propanoic acid CAS-No. 104582-29-8 Formula $C_{15}H_{25}N_3O_4S_3$ Mol. weight 407,57 g/mol			
RL-8415	Biotin-Hx-SS-Py		
5-((3aS,4S,6aR)-2-oxohexahydro-1H-thieno[3,4-d]imidazol-4-yl)-N-(6-(3-(pyridin-2-yl)disulfanyl)propanamido) hexyl)pentanamide CAS-No. 129179-83-5 Formula $C_{24}H_{37}N_5O_3S_3$ Mol. weight 539,77 g/mol			


RL-4120 Biotin-SS-N₃

N-(2-((2-azidoethyl)disulfanyl)ethyl)-5-((3aS,4S,6aR)-2-oxohexahydro-1H-thieno[3,4-d]imidazol-4-yl)pentanamide

CAS-No. 1620523-64-9

Formula C₁₄H₂₄N₆O₂S₃

Mol. weight 404,57 g/mol



PEG8100 Biotin-PEG(4)-SS-Azide

N-(2-((3-((3-azidopropyl)amino)-3-oxopropyl)disulfanyl)ethyl)-1-(5-((3aS,4S,6aR)-2-oxohexahydro-1H-thieno[3,4-d]imidazol-4-yl)pentanamido)-3,6,9,12-tetraoxapentadecan-15-amide

CAS-No. 1260247-52-6

Formula C₂₉H₅₂N₈O₈S₃

Mol. weight 736,96 g/mol



PEG8110 Biotin-PEG(4)-SS-Alkyne

N-(2-((3-oxo-3-(prop-2-ynylamino)propyl)disulfanyl)ethyl)-1-(5-((3aS,4S,6aR)-2-oxohexahydro-1H-thieno[3,4-d]imidazol-4-yl)pentanamido)-3,6,9,12-tetraoxapentadecan-15-amide

CAS-No. 1260247-54-8

Formula C₂₉H₄₉N₅O₈S₃


Mol. weight 691,92 g/mol


PEG8120 Biotin-PEG(4)-SS-DBCO

N-(2-((3-(3-(azadibenzocyclooctyn-1-yl)-3-oxopropylamino)-3-oxopropyl)disulfanyl)ethyl)-1-(5-((3aS,4S,6aR)-2-oxohexahydro-1H-thieno[3,4-d]imidazol-4-yl)pentanamido)-3,6,9,12-tetraoxapentadecan-15-amide

Formula C₄₄H₆₀N₆O₉S₃

Mol. weight 913,18 g/mol



PEG8090 Biotin-PEG(4)-SS-COOH






9,25-dioxo-29-((3aS,4S,6aR)-2-oxohexahydro-1H-thieno[3,4-d]imidazol-4-yl)-12,15,18,21-tetraoxa-4,5-dithia-8,24-diazanonacosan-1-oic acid





CAS-No. 1380166-80-2

Formula C₂₆H₄₆N₄O₉S₃

Mol. weight 654,86 g/mol



			Product details
LS-3570	Biotin-SS-Tyramide		
N-(2-((3-(4-hydroxyphenethylamino)-3-oxopropyl)disulfanyl)ethyl)-5-(2-oxohexahydro-1H-thieno[3,4-d]imidazol-4-yl)pentanamide CAS-No. 678975-20-7 Formula $C_{23}H_{34}N_4O_4S_3$ Mol. weight 526,74 g/mol			
LS-3930	Biotin-PEG(4)-SS-Tyramide		
N-(2-((3-(4-hydroxyphenethylamino)-3-oxopropyl)disulfanyl)ethyl)-1-(5-(2-oxohexahydro-1H-thieno[3,4-d]imidazol-4-yl)pentanamido)-3,6,9,12-tetraoxapentadecan-15-amide Formula $C_{34}H_{55}N_5O_9S_3$ Mol. weight 774,02 g/mol			
LS-3960	Tyramide-SS-amine*HCl		
3-((2-aminoethyl)disulfanyl)-N-(4-hydroxyphenethyl)propanamide hydrochloride CAS-No. 2576471-37-7 net Formula $C_{13}H_{20}N_2O_2S_2 \cdot HCl$ Mol. weight 300,44*36,45 g/mol			
RL-3500	OPSS-OpNC		
2-(2-Pyridithio)ethyl-<i>p</i>-nitrophenylcarbonate CAS-No. 874302-76-8 Formula $C_{14}H_{12}N_2O_5S_2$ Mol. weight 352,38 g/mol			
RL-3890	OPSS-PAB		
2-(pyridin-2-yl)disulfaneyl)ethyl 4-(hydroxymethyl)phenyl)carbamate CAS-No. 2362536-42-1 Formula $C_{15}H_{16}N_2O_3S_2$ Mol. weight 336,42 g/mol			
RL-3920	OPSS-Bzl-PAB		
4-(pyridin-2-yl)disulfaneyl)benzyl 4-(hydroxymethyl)phenyl)carbamate Formula $C_{20}H_{18}N_2O_3S_2$ Mol. weight 398,50 g/mol			

			Product details
RL-4170	2-OPSS-Bzl-OpNC		
O-(4-nitrophenyl) S-(2-(pyridin-2-yl)disulfanyl)benzyl carbonothioate			
CAS-No.	1384425-52-8		
Formula	C ₁₉ H ₁₄ N ₂ O ₄ S ₃		
Mol. weight	430,51 g/mol		
RL-3550	OPSS-Bzl-OpNC		
(4-(pyridin-2-yl)disulfanyl)benzyl p-nitrophenylcarbonate			
CAS-No.	1151989-04-6		
Formula	C ₁₉ H ₁₄ N ₂ O ₅ S ₂		
Mol. weight	414,45 g/mol		
RL-3820	OPSS-PAB-OpNC		
2-(pyridin-2-yl)disulfanyl ethyl (4-(((4-nitrophenoxy) carbonyl)oxy)methyl)phenyl)carbamate			
CAS-No.	2362536-70-5		
Formula	C ₂₂ H ₁₉ N ₃ O ₇ S ₂		
Mol. weight	501,53 g/mol		
RL-3850	OPSS-Bzl-PAB-OpNC		
4-(pyridin-2-yl)disulfanyl benzyl (4-(((4-nitrophenoxy) carbonyl)oxy)methyl)phenyl)carbamate			
Formula	C ₂₇ H ₂₁ N ₃ O ₇ S ₂		
Mol. weight	563,60 g/mol		

References:

- Labeling Carboxyl Groups of Surface-Exposed Proteins Provides an Orthogonal Approach for Cell Surface Isolation; N. E. Ozkan Kucuk, E. Sanal, E. Tan, T. Mitchison, N. Ozlu; **J Proteome Res** 2018; **17**: 1784-1793. <https://doi.org/10.1021/acs.jproteome.7b00825>
- Development of trifunctional probes for glycoproteomic analysis; C. S. Tsai, P. Y. Liu, H. Y. Yen, T. L. Hsu, C. H. Wong; **Chem Commun** 2010; **46**: 5575-7. <https://doi.org/10.1039/c0cc00345j>
- Fluorometric assay for tissue transglutaminase-mediated transamidation activity; C. Gnaccarini, W. Ben-Tahar, W. D. Lubell, J. N. Pelletier, J. W. Keillor; **Bioorg Med Chem** 2009; **17**: 6354-9. <https://doi.org/10.1016/j.bmc.2009.07.031>
- Convenient synthesis of photoaffinity probes and evaluation of their labeling abilities; T. Kan, Y. Kita, Y. Morohashi, Y. Tominari, S. Hosoda, T. Tomita, H. Natsugari, T. Iwatsubo, T. Fukuyama; **Org Lett** 2007; **9**: 2055-8. <https://doi.org/10.1021/ol070376i>
- Facile synthesis toward the construction of an activity probe library for glycosidases; T. H. Shie, Y. L. Chiang, J. J. Lin, Y. K. Li, L. C. Lo; **Carbohydr Res** 2006; **341**: 443-56. <https://doi.org/10.1016/j.carres.2005.12.005>
- A simple photo-affinity labeling protocol; H.-y. Li, Y. Liu, K. Fang, K. Nakanishi; **Chemical Communications** 1999; 365-366. <https://doi.org/10.1039/a809507h>

- *Linker Technologies for Antibody–Drug Conjugates*; B. Nolting; **Antibody–Drug Conjugates L. Ducry** 2013; **1045**: 71–100. https://doi.org/10.1007/978-1-62703-541-5_5
- *Disulfide-Based Self-Immolative Linkers and Functional Bioconjugates for Biological Applications*; Z. Deng, J. Hu, S. Liu; **Macromol Rapid Commun** 2020; **41**: e1900531. <https://doi.org/10.1002/marc.201900531>
- *Reduction-Triggered Transformation of Disulfide-Containing Micelles at Chemically Tunable Rates*; Z. Deng, S. Yuan, R. X. Xu, H. Liang, S. Liu; **Angew. Chem. Int. Ed.** 2018; **57**: 8896–8900. <https://doi.org/10.1002/anie.201802909>
- *Modulated Fragmentation of Proapoptotic Peptide Nanoparticles Regulates Cytotoxicity*; T. Suma, J. Cui, M. Müllner, S. Fu, J. Tran, K. F. Noi, Y. Ju, F. Caruso; **J. Am. Chem. Soc.** 2017; **139**: 4009–4018. <https://doi.org/10.1021/jacs.6b11302>
- *Engineering Intracellular Delivery Nanocarriers and Nanoreactors from Oxidation-Responsive Polymersomes via Synchronized Bilayer Cross-Linking and Permeabilizing Inside Live Cells*; Z. Deng, Y. Qian, Y. Yu, G. Liu, J. Hu, G. Zhang, S. Liu; **J. Am. Chem. Soc.** 2016; **138**: 10452–10466. <https://doi.org/10.1021/jacs.6b04115>

The Concept of Antibody–Drug Conjugation (ADC)

Permanent Linkers

Cleavable Linkers

Trifunctional Linkers

Cross-Linkers for other Bio Applications

Preparing Carriers for Conjugation


Index

3.6. Dde-Based Linkers

The Dde [N-1-(4,4-dimethyl-2,6-dioxocyclohex-1-ylidene)-3-ethyl] protecting group is commonly utilized to protect the side-chain amine groups of lysine, ornithine, 2,4-diaminobutyric acid, and 2,3-diaminopropionic acid. Dde shows orthogonal cleavage conditions to Fmoc (piperidine or DBU) and tBu (TFA) deprotecting protocols and is stable to denaturing washing conditions, while allowing for a mild and selective removal in the presence of other protecting groups using a buffered aqueous solution of hydrazine or hydroxylamine, thus representing a versatile tool for the site-specific modification of peptides. Advantageously, the cleavage can be followed spectrophotometrically since the reaction product of Dde with hydrazine is a chromophoric derivative.

Placing Dde as one terminal group of a linker and a functional group prone for conjugation as the other, or using Dde as the central connective portion of a linker, allows for the creation of new bifunctional linkers that can be selectively and temporarily attached to:

- Appropriately modified biomolecules for binding to streptavidin (with terminal biotin) (*Fig. 25 (A)*), or conjugation to any solid supports, e.g., via Click reaction (*Fig. 25 (B)*).
- Solubilizing tags, e.g., hexa-lysine (“helping-hand linkers”, *Fig. 25 (C)*), oligo-arginine, PEGs (*Fig. 25 (D)*) or other hydrophilic groups improving solubility of hydrophobic peptides or other compounds when being attached to either the N-terminus or any lysine side-chain within a peptide sequence.
- Dyes and any other conjugate for monitoring, diagnostics, targeting or other purposes.



Dde/ivDde linkers are implemented in simple and nearly quantitative steps:

1. Orthogonal deprotection of lysine residues in a peptide or N-terminus or any other amino function of a hydrophobic compound.
2. On-resin incorporation of the linker.
3. Fmoc-SPPS elongation.
4. Cleavage of the peptide from the resin and removal of all side-chain protecting groups.
5. The tagged peptide can be separated from truncated sequences.
6. In-solution cleavage using mild aqueous hydrazine to cleave the Dde linker after purification, streptavidin attachment, NCL-based assembly or another reaction step. The cleavage can be monitored spectroscopically as the resulting pyrazole shows a strong absorption at 290 nm.

Dde/ivDde becomes particularly useful for handling and purification of insoluble and aggregationprone peptides, as any appropriate solubilizing promoting group can be attached to create so-called “helping-hand” linkers that can be removed in a traceless manner (Fig. 26).




Fig. 26: Dde-based “helping-hand” linkers improve the solubility and allow for the purification of hydrophobic peptides. Removal of the handle can be easily monitored spectroscopically at 290 nm.

The Concept of Antibody-Drug Conjugation (ADC)

Permanent Linkers

Cleavable Linkers

Trifunctional Linkers

Cross-Linkers for other Bio Applications

Preparing Carriers for Conjugation

Index


Procedure for Removing Helping-Hands from Peptides (adapted from Jacobsen *et al.*, *J Am Chem Soc* 2016):

10 mL of 2 M hydrazine stock solution (pH 7.5) is being prepared as follows:

1. Weigh 5.7 g Guanidinium chloride and 75 mg DTT into 15 mL Falcon tube.
2. Add 1 mL of 1 M NaH_2PO_4 .
3. Add 2 mL of 10 M hydrazine in water.
4. Add 0.5 mL of 12 M HCl.
5. Dissolve solution by thorough vortexing.
6. Adjust pH to 7.5 by adding concentrated HCl.
7. Fill to a final volume of 10 mL with water.
8. Filter solution using 0.2 μm syringe filter.

Cleavage of the helping-hand can be triggered by equivolume addition of 2 M hydrazine stock solution into the solution of the peptide. After adding the hydrazine solution, subtle adjustment may be necessary to achieve a final solution pH of 7.5. The reaction is normally completed within minutes. Deprotection can be monitored spectrophotometrically at 290 nm.

Despite its widespread use, the biotinylation of proteins for subsequent purification *via* Streptavidin beads bears certain hurdles, e.g., concerning the removal of the proteins from the beads due to the strong binding. One possible improvement is represented by the use of appropriately derivatized Dde-linkers. The connection of such a bifunctional linker with a biotin moiety on the one end, and a clickable group (alkyne, e.g., DBCO) or tyramide on the other, allows for the selective attachment to appropriately modified biomolecules, as well as the mild release of captured proteins from the beads after purification (Fig. 27).



Aside from the commonly used cleavage solution for Dde consisting of 2% hydrazine monohydrate in H₂O, the following procedure may be used in order to ensure full orthogonality between Dde and Fmoc.

Selective Removal of Dde/ivDde using hydroxylamine (adapted from Díaz-Mochón *et al.*, *Org. Lett.* 2004):

1.25 g (1.80 mmol) of NH₂OH·HCl and 0.918 g (1.35 mmol) of Imidazole were suspended in 5 mL NMP, and the mixture sonicated until complete dissolution. This solution can be stirred for at least 2 weeks at -20 °C. Just before reaction, five volumes of this solution were diluted with one volume of alternatively DCM or DMF.

Dde is easy to cleave, but not very robust. Thus, during Fmoc cleavage, Dde might migrate to free lysine ε-amino groups (“scrambling”) or, in rare cases, even to the peptide’s free amino terminus. Especially during the synthesis of longer peptide sequences, a certain extent of Dde is removed during Fmoc cleavage with piperidine.

The sterically more demanding protecting group ivDde is more stable towards piperidine and does usually not migrate to free lysine amino groups. However, for some sequences, total removal of the robust ivDde is hardly possible - especially near the C-terminus or in aggregating sequences.

A newer group for the orthogonal protecting of amino groups reported in literature is MeDmb (methyl dimethylbarbituric acid). At Iris Biotech, we are also offering the sterically more demanding ivDmb. There is no “one-fits-all” when you are to decide on the protecting group selection for orthogonal lysine side-chain protection. When your peptide shows a low tendency for scrambling, Lys(Dde) may be fine. Otherwise, you can switch to more robust ones such as Lys(ivDde) or Lys(ivDmb). If you do not succeed in completely removing the ivDde protecting group, the brand-new ivDmb should be your choice. For more details on the ivDmb protecting group, see our blog post:

<https://www.iris-biotech.de/blog/potm-next-generation-lysine-side-chain-protecting-groups/>

References:

- A novel lysine-protecting procedure for continuous flow solid phase synthesis of branched peptides; B. W. Bycroft, W. C. Chan, S. R. Chhabra, N. D. Hone; **Journal of the Chemical Society, Chemical Communications** 1993: 778-779. <https://doi.org/10.1039/c39930000778>
- An appraisal of new variants of Dde amine protecting group for solid phase peptide synthesis; S. R. Chhabra, B. Hothi, D. J. Evans, P. D. White, B. W. Bycroft, W. C. Chan; **Tetrahedron Letters** 1998; **39**: 1603-1606. [https://doi.org/10.1016/S0040-4039\(97\)10828-0](https://doi.org/10.1016/S0040-4039(97)10828-0)
- Investigation on the stability of the Dde protecting group used in peptide synthesis: migration to an unprotected lysine; K. Augustyns, W. Kraas, G. Jung; **J Pept Res** 1998; **51**: 127-33. <https://doi.org/10.1111/j.1399-3011.1998.tb00630.x>
- Evaluation of ivDde as a quasi-orthogonal protecting group for Fmoc solid-phase peptide synthesis; R. R. Wilhelm, A. Srinivasan, M. A. Schmidt; *Peptides for the New Millennium: Proceedings of the 16th American Peptide Symposium June 26–July 1, 1999, Minneapolis, Minnesota, USA 2000*: 58-59. https://doi.org/10.1007/0-306-46881-6_19
- Synthesis of a chlorothalonil peptide conjugate mimicking protein-bound pesticide residues; H. Hrenn, W. Schwack, W. Seilmeier, H. Wieser; **Tetrahedron Lett** 2003; **44**: 1911-1913. [https://doi.org/10.1016/S0040-4039\(03\)00121-7](https://doi.org/10.1016/S0040-4039(03)00121-7)
- Reaction of 1,3-dimethyl-5-acetyl-barbituric acid (DAB) with primary amines. Access to intermediates for selectively protected spermidines; E. T. da Silva, E. L. S. Lima; **Tetrahedron Lett** 2003; **44**: 3621-3624. [https://doi.org/10.1016/S0040-4039\(03\)00709-3](https://doi.org/10.1016/S0040-4039(03)00709-3)
- Scope and Limitations of Barbituric and Thiobarbituric Amino Acid Derivatives as Protecting Groups for Solid-Phase Peptide Synthesis: Towards a Green Protecting Group; S. Ramkisson, H. H. Al-Rasheed, K. A. Dahlous, B. G. De La Torre, A. El-Faham, F. Albericio; **ChemistrySelect** 2021; **6**: 6626-6630. <https://doi.org/10.1002/slct.202101539>

[↑ back to content](#)

The Concept of Antibody-Drug Conjugation (ADC)

Permanent Linkers

Cleavable Linkers

Trifunctional Linkers

Cross-Linkers for other Bio Applications


Preparing Carriers for Conjugation

Index

LS-4020 Biotin-Dde

2-(1-hydroxy-5-((3aS,4S,6aR)-2-oxohexahydro-1H-thieno[3,4-d]imidazol-4-yl)pentylidene)-5,5-dimethylcyclohexane-1,3-dione


CAS-No. 194038-08-9
Formula $C_{18}H_{26}N_2O_4S$
Mol. weight 366,48 g/mol



LS-4000 Biotin-Dde-Tyramide

2-(1-(4-hydroxyphenethylamino)-5-((3aS,4S,6aR)-2-oxohexahydro-1H-thieno[3,4-d]imidazol-4-yl)pentylidene)-5,5-dimethylcyclohexane-1,3-dione


CAS-No. 2819732-80-2
Formula $C_{26}H_{35}N_2O_4S$
Mol. weight 485,64 g/mol



PEG8130 Biotin-PEG(4)-Dde-Tyramide

N-(15-(4,4-dimethyl-2,6-dioxocyclohexylidene)-18-(4-hydroxyphenyl)-3,6,9,12-tetraoxa-16-aza-octadecyl)-5-((3aS,4S,6aR)-2-oxohexahydro-1H-thieno[3,4-d]imidazol-4-yl)pentanamide


CAS-No. 2814457-48-0
Formula $C_{34}H_{56}N_4O_9S$
Mol. weight 732,93 g/mol



PEG7960 Biotin-PEG(4)-Dde-N₃

N-(19-azido-15-(4,4-dimethyl-2,6-dioxocyclohexylidene)-3,6,9,12-tetraoxa-16-azanonadecyl)-5-((3aS,4S,6aR)-2-oxohexahydro-1H-thieno[3,4-d]imidazol-4-yl)pentanamide


CAS-No. 1802907-93-2
Formula $C_{32}H_{53}N_7O_8S$
Mol. weight 695,87 g/mol








PEG7970 Biotin-PEG(4)-Dde-Picolyl-N₃

N-(6-(azidomethyl)pyridin-3-yl)-15-(4,4-dimethyl-2,6-dioxocyclohexylidene)-1-(5-((3aS,4S,6aR)-2-oxohexahydro-1H-thieno[3,4-d]imidazol-4-yl)pentanamide)-3,6,9,12-tetraoxa-16-azanonadecan-19-amide

CAS-No. 2055048-42-3
Formula $C_{38}H_{57}N_9O_9S$
Mol. weight 815,98 g/mol














		Product details
RL-8700	N₃-L-Glu(Dde)-OtBu	
<i>tert</i> -butyl (S)-2-azido-5-(4,4-dimethyl-2,6-dioxocyclohexylidene)-5-hydroxypentanoate		
Formula	C ₁₇ H ₂₅ N ₃ O ₅	 
Mol. weight	351,40 g/mol	
PEG7980	Biotin-PEG(4)-Dde-Alkyne	
N-(15-(4,4-dimethyl-2,6-dioxocyclohexylidene)-3,6,9,12-tetraoxa-16-azanonadec-18-ynyl)-5-((3aS,4S,6aR)-2-oxohexahydro-1H-thieno[3,4-d]imidazol-4-yl)pentanamide		
CAS-No.	1802908-00-4	 
Formula	C ₃₂ H ₅₀ N ₄ O ₈ S	
Mol. weight	650,83 g/mol	
PEG8140	Biotin-PEG(4)-Dde-DBCO	
N-(15-(4,4-dimethyl-2,6-dioxocyclohexylidene)-19-oxo-19-(azadibenzocyclooctyn-1-yl)-3,6,9,12-tetraoxa-16-azanonadecyl)-5-((3aS,4S,6aR)-2-oxohexahydro-1H-thieno[3,4-d]imidazol-4-yl)pentanamide		
CAS-No.	1807512-43-1	 
Formula	C ₄₇ H ₆₁ N ₅ O ₉ S	
Mol. weight	872,08 g/mol	
RL-3280	N₃-Pen-Dde	
2-(5-azido-1-hydroxypentylidene)-5,5-dimethylcyclohexane-1,3-dione		
CAS-No.	1867129-38-1	 
Formula	C ₁₃ H ₁₉ N ₃ O ₃	
Mol. weight	265,31 g/mol	
RL-3290	N₃-Pen-Dtpp	
5-(5-azido-1-hydroxypentylidene)-1,3-dimethylpyrimidine-2,4,6(1H,3H,5H)-trione		
CAS-No.	1867129-42-7	 
Formula	C ₁₁ H ₁₅ N ₅ O ₄	
Mol. weight	281,27 g/mol	

The Dde derived linker might cleave under mildly acidic and even neutral conditions in the one or the other case. The DTPM derived linker is totally stable under acidic conditions as well as to a wide range of chemical treatments, including particularly harsh sodium methoxide-based deacetylation of chemically introduced glycans.

Reference:

- Combining triazole ligation and enzymatic glycosylation on solid phase simplifies the synthesis of very long glycoprotein analogues; M. Galibert, V. Piller, F. Piller, V. Aucagne, A. F. Delmas; **Chem. Sci.** 2015; **6**: 3617-3623.
<https://doi.org/10.1039/c5sc00773a>

		Product details
FAA8115	Fmoc-L-Lys(Pentynoyl-DIM)-OH	
N-alpha-(9-Fluorenylmethyloxycarbonyl)-N-epsilon-[1-(4,4-dimethyl-2,6-dioxocyclohexylidene)pent-4-yn-1-yl]-L-lysine		
CAS-No.	2408993-33-7	
Formula	C ₃₄ H ₃₈ N ₂ O ₆	
Mol. weight	570,69 g/mol	
		
FAA8145	Fmoc-L-Lys(N₃-Aca-DIM)-OH	
N-alpha-(9-Fluorenylmethyloxycarbonyl)-N-epsilon-[6-azido-1-(4,4-dimethyl-2,6-dioxocyclohexylidene)hexyl]-L-lysine		
CAS-No.	2408993-39-3	
Formula	C ₃₅ H ₄₃ N ₃ O ₆	
Mol. weight	629,76 g/mol	
		
RL-3260	Fmoc-Aca-DIM	
6-((9-Fluorenylmethyl)oxycarbonylamino)-1-(4,4-dimethyl-2,6-dioxocyclohexylidene)-hexan-1-ol		
CAS-No.	2379561-08-5	
Formula	C ₂₉ H ₃₃ NO ₅	
Mol. weight	475,58 g/mol	
		
RL-3270	Fmoc-AEEP-DIM	
3-(2-(2-(9-Fluorenylmethyl)oxycarbonylaminoethoxy)ethoxy)-1-(4,4-dimethyl-2,6-dioxocyclohexylidene)-propan-1-ol		
CAS-No.	1988771-96-5	
Formula	C ₃₀ H ₃₅ NO ₇	
Mol. weight	521,60 g/mol	
		
RL-3470	Fmoc-AEEP-DMB	
(9-Fluorenylmethyloxycarbonyl)amino-PEG(2)-Dtp		
Formula	C ₂₈ H ₃₁ N ₃ O ₈	
Mol. weight	537,57 g/mol	
		

			Product details
DAA1016	Dde-O₂Oc-OH		
8-[(4,4-Dimethyl-2,6-dioxocyclohex-1-ylidene)ethyl-amino]-3,6-dioxaoctanoic acid, {2-[2-(Dde-amino)ethoxy]ethoxy}acetic acid			
CAS-No.	1263045-93-7		
Formula	C ₁₆ H ₂₅ NO ₆		
Mol. weight	327,37 g/mol		
			
PEG8150	Fmoc-PEG(4)-Dde		
1-(9H-Fluorenylmethyloxycarbonylamino)-15-(4,4-dimethyl-2,6-dioxocyclohexylidene)-3,6,9,12-tetraoxapentadecyl-15-ol			
CAS-No.	2093409-87-9		
Formula	C ₃₄ H ₄₃ NO ₉		
Mol. weight	609,71 g/mol		
			
BAA1191	Boc-L-Dab(Dde)-OH		
N-alpha-t-Butyloxycarbonyl-N-gamma-[1-(4,4-dimethyl-2,6-dioxocyclohex-1-ylidene)ethyl]-L-2,4-diaminobutyric acid			
CAS-No.	1263045-50-6		
Formula	C ₁₉ H ₃₀ N ₂ O ₆		
Mol. weight	382,46 g/mol		
			
BAA1171	Boc-D-Dab(Dde)-OH		
N-alpha-t-Butyloxycarbonyl-N-gamma-[1-(4,4-dimethyl-2,6-dioxocyclohex-1-ylidene)ethyl]-D-2,4-diaminobutyric acid			
CAS-No.	1263046-41-8		
Formula	C ₁₉ H ₃₀ N ₂ O ₆		
Mol. weight	382,46 g/mol		
			
BAA1193	Boc-L-Dap(Dde)-OH		
N-alpha-t-Butyloxycarbonyl-N-beta-[1-(4,4-dimethyl-2,6-dioxocyclohex-1-ylidene)ethyl]-L-2,3-diaminopropionic acid			
CAS-No.	1263045-09-5		
Formula	C ₁₈ H ₂₈ N ₂ O ₆		
Mol. weight	368,43 g/mol		
			
BAA1176	Boc-D-Dap(Dde)-OH		
N-alpha-t-Butyloxycarbonyl-N-beta-[1-(4,4-dimethyl-2,6-dioxocyclohex-1-ylidene)ethyl]-D-2,3-diaminopropionic acid			
CAS-No.	1263047-33-1		
Formula	C ₁₈ H ₂₈ N ₂ O ₆		
Mol. weight	368,43 g/mol		
			

The Concept of Antibody-Drug Conjugation (ADC)

Permanent Linkers

Cleavable Linkers

Trifunctional Linkers

Cross-Linkers for other Bio Applications

Preparing Carriers for Conjugation

Index

[↑ back to content](#)


DAA1012 Dde-L-Dap(Fmoc)-OH

N-alpha-(4-4-Dimethyl-2,6-dioxocyclohex-1-ylidene)ethyl-N-beta-(9-fluorenylmethyloxycarbonyl)-L-2,3-diaminopropionic acid

CAS-No. 1263046-98-5

Formula $C_{28}H_{30}N_2O_6$

Mol. weight 490,56 g/mol




DAA1006 Dde-D-Dap(Fmoc)-OH

N-alpha-(4-4-Dimethyl-2,6-dioxocyclohex-1-ylidene)ethyl-N-beta-(9-fluorenylmethyloxycarbonyl)-D-2,3-diaminopropionic acid

CAS-No. 1263046-87-2

Formula $C_{28}H_{30}N_2O_6$

Mol. weight 490,56 g/mol




DAA1010 Dde-L-Dab(Fmoc)-OH

N-alpha-(4-4-Dimethyl-2,6-dioxocyclohex-1-ylidene)ethyl-N-gamma-(9-fluorenylmethyloxycarbonyl)-L-2,4-diaminobutyric acid

CAS-No. 1263045-85-7

Formula $C_{29}H_{32}N_2O_6$

Mol. weight 504,59 g/mol




DAA1004 Dde-D-Dab(Fmoc)-OH

N-alpha-(4-4-Dimethyl-2,6-dioxocyclohex-1-ylidene)ethyl-N-gamma-(9-fluorenylmethyloxycarbonyl)-D-2,4-diaminobutyric acid

CAS-No. 1263046-84-9

Formula $C_{29}H_{32}N_2O_6$

Mol. weight 504,59 g/mol




BAA1197 Boc-L-Orn(Dde)-OH

N-alpha-t-Butyloxycarbonyl-N-delta-[1-(4,4-dimethyl-2,6-dioxocyclohex-1-ylidene)ethyl]-L-ornithine

CAS-No. 1272755-14-2

Formula $C_{20}H_{32}N_2O_6$

Mol. weight 396,49 g/mol




DAA1001 Dde-L-Orn(Aloc)-OH







N-alpha-(4-4-Dimethyl-2,6-dioxocyclohex-1-ylidene)ethyl-N-delta-allyloxycarbonyl-L-ornithine

CAS-No. 1423017-98-4

Formula $C_{19}H_{28}N_2O_6$

Mol. weight 380,44 g/mol



		Product details
DAA1002	Dde-L-Orn(Fmoc)-OH	
N-alpha-(4-4-Dimethyl-2,6-dioxocyclohex-1-ylidene)ethyl-N-delta-(9-fluorenylmethoxycarbonyl)-L-ornithine		
CAS-No.	1423017-87-1	
Formula	C ₃₀ H ₃₄ N ₂ O ₆	
Mol. weight	518,62 g/mol	
		
BAA1286	Boc-L-Lys(Dde)-OH*DCHA	
N-alpha-t-Butyloxycarbonyl-N-epsilon-(4,4-dimethyl-2,6-dioxocyclohex-1-ylidene)ethyl-L-lysine dicyclohexylamine		
CAS-No.	444795-66-8	
Formula	C ₂₁ H ₃₄ N ₂ O ₆ *C ₁₂ H ₂₃ N	
Mol. weight	410,51*181,32 g/mol	
		
BAA1287	Boc-L-Lys(ivDde)-OH	
N-alpha-t-Butyloxycarbonyl-N-epsilon-[1-(4,4-dimethyl-2,6-dioxocyclohex-1-ylidene)-3-methylbutyl]-L-lysine		
CAS-No.	862847-44-7	
Formula	C ₂₄ H ₄₀ N ₂ O ₆	
Mol. weight	452,6 g/mol	
		
DAA1013	Dde-L-Lys(Aloc)-OH*DCHA	
N-alpha-(4-4-Dimethyl-2,6-dioxocyclohex-1-ylidene)ethyl-N-epsilon-allyloxycarbonyl-L-lysine dicyclohexylamine		
CAS-No.	264230-73-1 net	
Formula	C ₂₀ H ₃₀ N ₂ O ₆ *C ₁₂ H ₂₃ N	
Mol. weight	394,47*181,32 g/mol	
		
DAA1014	Dde-L-Lys(Boc)-OH	
N-alpha-(4-4-Dimethyl-2,6-dioxocyclohex-1-ylidene)ethyl-N-epsilon-t-butyloxycarbonyl-L-lysine		
CAS-No.	1189586-14-8	
Formula	C ₂₁ H ₃₄ N ₂ O ₆	
Mol. weight	410,51 g/mol	
		
DAA1015	Dde-L-Lys(Fmoc)-OH	
N-alpha-(4-4-Dimethyl-2,6-dioxocyclohex-1-ylidene)ethyl-N-epsilon-(9-fluorenylmethoxycarbonyl)-L-lysine		
CAS-No.	156648-40-7	
Formula	C ₃₁ H ₃₆ N ₂ O ₆	
Mol. weight	532,64 g/mol	
		

The Concept of Antibody-Drug Conjugation (ADC)

Permanent Linkers

Cleavable Linkers

Trifunctional Linkers

Cross-Linkers for other Bio Applications

Preparing Carriers for Conjugation


Index

[↑ back to content](#)

DAA1017 Dde-D-Lys(Fmoc)-OH

N-alpha-(4,4-Dimethyl-2,6-dioxocyclohex-1-ylidene)ethyl-N-epsilon-(9-fluorenylmethyloxycarbonyl)-D-lysine


CAS-No. 1301706-71-7
Formula $C_{31}H_{36}N_2O_6$
Mol. weight 532,64 g/mol



DAA1018 ivDde-L-Dap(Fmoc)-OH

N-alpha-[(4,4-Dimethyl-2,6-dioxocyclohex-1-ylidene)-3-methylbutyl]-N-beta-(9-fluorenylmethyloxycarbonyl)-L-2,3-diaminopropionic acid


CAS-No. 2389078-71-9
Formula $C_{31}H_{36}N_2O_6$
Mol. weight 532,63 g/mol



DAA1019 ivDde-L-Lys(Fmoc)-OH

N-alpha-[(4,4-Dimethyl-2,6-dioxocyclohex-1-ylidene)-3-methylbutyl]-N-epsilon-(9-fluorenylmethyloxycarbonyl)-L-lysine


CAS-No. 1446752-60-8
Formula $C_{34}H_{42}N_2O_6$
Mol. weight 574,71 g/mol



DAA1030 ivDde-D-Lys(Fmoc)-OH

N-alpha-[(4,4-Dimethyl-2,6-dioxocyclohex-1-ylidene)-3-methylbutyl]-N-epsilon-(9-fluorenylmethyloxycarbonyl)-D-lysine


CAS-No. 2308529-94-2
Formula $C_{34}H_{42}N_2O_6$
Mol. weight 574,71 g/mol



DAA1020 Dde-L-Met-OH

N-alpha-[1-(4,4-dimethyl-2,6-dioxocyclohexylidene)ethyl]-L-methionine


CAS-No. 1435266-87-7
Formula $C_{15}H_{23}NO_4S$
Mol. weight 313,13 g/mol









FAA8690 Fmoc-Aeg(Dde)-OH

(9-Fluorenylmethyloxycarbonyl)-N-(2-((1-(4,4-dimethyl-2,6-dioxocyclohexylidene)ethyl)amino)ethyl)glycine

CAS-No. 2988301-08-0
Formula $C_{29}H_{32}N_2O_6$
Mol. weight 504,58 g/mol




		Product details
FAA8815	Fmoc-L-Abu(3-Dde-amino)-OH (2S,3S) 2-(Fmoc-amino)-3-(Dde-amino)butanoic acid (2S,3S) Formula $C_{29}H_{32}N_2O_6$ Mol. weight 504,58 g/mol	 
FAA1462	Fmoc-L-Dap(Dde)-OH N-alpha-(9-Fluorenylmethyloxycarbonyl)-N-beta-[(4,4-dimethyl-2,6-dioxocyclohex-1-ylidene)ethyl]-2,3-diaminopropionic acid CAS-No. 247127-51-1 Formula $C_{28}H_{30}N_2O_6$ Mol. weight 490,56 g/mol	 
FAA1476	Fmoc-D-Dap(Dde)-OH N-alpha-(9-Fluorenylmethyloxycarbonyl)-N-beta-[(4,4-dimethyl-2,6-dioxocyclohex-1-ylidene)ethyl]-D-2,3-diaminopropionic acid CAS-No. 210830-03-8 Formula $C_{28}H_{30}N_2O_6$ Mol. weight 490,56 g/mol	 
FAA1464	Fmoc-L-Dap(ivDde)-OH N-alpha-(9-Fluorenylmethyloxycarbonyl)-N-beta-[1-(4,4-dimethyl-2,6-dioxocyclohex-1-ylidene)-3-methylbutyl]-L-2,3-diaminopropionic acid CAS-No. 607366-20-1 Formula $C_{31}H_{36}N_2O_6$ Mol. weight 532,64 g/mol	 
FAA1478	Fmoc-D-Dap(ivDde)-OH N-alpha-(9-Fluorenylmethyloxycarbonyl)-N-beta-[1-(4,4-dimethyl-2,6-dioxocyclohex-1-ylidene)-3-methylbutyl]-D-2,3-diaminopropionic acid CAS-No. 1228900-15-9 Formula $C_{31}H_{36}N_2O_6$ Mol. weight 532,64 g/mol	 
FAA1365	Fmoc-L-Dab(Dde)-OH N-alpha-(9-Fluorenylmethyloxycarbonyl)-N-gamma-(4,4-dimethyl-2,6-dioxocyclohex-1-ylidene)ethyl-L-2,4-diaminobutyric acid CAS-No. 235788-61-1 Formula $C_{29}H_{32}N_2O_6$ Mol. weight 504,59 g/mol	 

FAA1318 Fmoc-D-Dab(Dde)-OH

N-alpha-(9-Fluorenylmethyloxycarbonyl)-N-gamma-(4,4-dimethyl-2,6-dioxocyclohex-1-ylidene)ethyl-D-2,4-diaminobutyric acid


CAS-No. 596797-14-7
Formula $C_{29}H_{32}N_2O_6$
Mol. weight 504,59 g/mol



FAA1458 Fmoc-L-Dab(ivDde)-OH

N-alpha-(9-Fluorenylmethyloxycarbonyl)-N-gamma-[1-(4,4-dimethyl-2,6-dioxocyclohex-1-ylidene)-3-methylbutyl]-L-2,4-diaminobutyric acid


CAS-No. 607366-21-2
Formula $C_{32}H_{38}N_2O_6$
Mol. weight 546,67 g/mol



FAA1473 Fmoc-D-Dab(ivDde)-OH

N-alpha-(9-Fluorenylmethyloxycarbonyl)-N-gamma-[1-(4,4-dimethyl-2,6-dioxocyclohex-1-ylidene)-3-methylbutyl]-D-2,4-diaminobutyric acid


CAS-No. 872169-32-3
Formula $C_{32}H_{38}N_2O_6$
Mol. weight 546,67 g/mol



FAA1502 Fmoc-L-Orn(Dde)-OH

N-alpha-(9-Fluorenylmethyloxycarbonyl)-N-delta-[1-(4,4-dimethyl-2,6-dioxocyclohex-1-ylidene)ethyl]-L-ornithine


CAS-No. 269062-80-8
Formula $C_{30}H_{34}N_2O_6$
Mol. weight 518,62 g/mol



FAA2090 Fmoc-D-Orn(Dde)-OH

N-alpha-(9-Fluorenylmethyloxycarbonyl)-N-delta-[1-(4,4-dimethyl-2,6-dioxocyclohex-1-ylidene)ethyl]-D-ornithine


CAS-No. 1419640-31-5
Formula $C_{30}H_{34}N_2O_6$
Mol. weight 518,62 g/mol













FAA1503 Fmoc-L-Orn(ivDde)-OH

N-alpha-(9-Fluorenylmethyloxycarbonyl)-N-delta-[1-(4,4-dimethyl-2,6-dioxocyclohex-1-ylidene)-3-methylbutyl]-L-ornithine solvate

CAS-No. 1198321-33-3
Formula $C_{33}H_{40}N_2O_6$
Mol. weight 560,7 g/mol



		Product details
FAA1493 Fmoc-D-Orn(ivDde)-OH.solv. N-alpha-(9-Fluorenylmethyloxycarbonyl)-N-de-lta-[1-(4,4-dimethyl-2,6-dioxocyclohex-1-ylidene)-3-methylbutyl]-D-ornithine.solv. CAS-No. 1272754-86-5 Formula $C_{33}H_{40}N_2O_6$ Mol. weight 560,7 g/mol		
FAA1390 Fmoc-L-Lys(Dde)-OH N-alpha-(9-Fluorenylmethyloxycarbonyl)-N-ep-silon-(4,4-dimethyl-2,6-dioxocyclohex-1-ylidene)ethyl-L-lysine CAS-No. 150629-67-7 Formula $C_{31}H_{36}N_2O_6$ Mol. weight 532,64 g/mol		
FAA1486 Fmoc-D-Lys(Dde)-OH N-alpha-(9-Fluorenylmethyloxycarbonyl)-N-ep-silon-[(4,4-dimethyl-2,6-dioxocyclohex-1-ylidene)ethyl]-D-lysine CAS-No. 333973-51-6 Formula $C_{31}H_{36}N_2O_6$ Mol. weight 532,64 g/mol		
FAA1488 Fmoc-D-Lys(ivDde)-OH N-alpha-(9-Fluorenylmethyloxycarbonyl)-N-ep-silon-[1-(4,4-dimethyl-2,6-dioxocyclohex-1-ylidene)-3-methylbutyl]-D-lysine CAS-No. 1272755-33-5 Formula $C_{34}H_{42}N_2O_6$ Mol. weight 574,72 g/mol		
FAA1500 Fmoc-L-Lys(ivDde)-OH N-alpha-(9-Fluorenylmethyloxycarbonyl)-N-ep-silon-[1-(4,4-dimethyl-2,6-dioxocyclohex-1-ylidene)-3-methylbutyl]-L-lysine CAS-No. 204777-78-6 Formula $C_{34}H_{42}N_2O_6$ Mol. weight 574,72 g/mol		
FAA8840 Fmoc-L-Lys(MeDmb)-OH (2S)-6-[[1-(1,3-dimethyl-2,4,6-trioxo-1,3-diazinan-5-ylidene)ethyl]amino]-2-(((9H-fluoren-9-yl)methoxy)carbonyl)amino)hexanoic acid Formula $C_{29}H_{32}N_4O_7$ Mol. weight 548,60 g/mol		

		Product details
FAA8845 Fmoc-D-Lys(MeDmb)-OH (2R)-6-[[1-(1,3-dimethyl-2,4,6-trioxo-1,3-diazinan-5-ylidene)ethyl]amino]-2-[[[(9H-fluoren-9-yl)methoxy]carbonyl]amino]hexanoic acid Formula $C_{29}H_{32}N_4O_7$ Mol. weight 548,60 g/mol		
FAA7975 Fmoc-L-Lys(ivDmb)-OH N2-(((9H-fluoren-9-yl)methoxy)carbonyl)-N6-(1-(1,3-dimethyl-2,4,6-trioxotetrahydropyrimidin-5(2H)-ylidene)-3-methylbutyl)-L-lysine Formula $C_{32}H_{38}N_4O_7$ Mol. weight 590,68 g/mol		
FAA1401 Fmoc-L-MeLys(Dde)-OH N-alpha-(9-Fluorenylmethoxycarbonyl)-N-alpha-methyl-N-epsilon-(4,4-dimethyl-2,6-dioxocyclohex-1-ylidene)ethyl-L-lysine CAS-No. 1428229-84-8 Formula $C_{32}H_{38}N_2O_6$ Mol. weight 546,67 g/mol		
FAA7935 Fmoc-L-MeLys(ivDde)-OH N2-(((9H-fluoren-9-yl)methoxy)carbonyl)-N6-(1-(4,4-dimethyl-2,6-dioxocyclohexylidene)-3-methylbutyl)-N2-methyl-L-lysine CAS-No. 1173996-67-2 Formula $C_{35}H_{44}N_2O_6$ Mol. weight 588,75 g/mol		

References:

- Synthesis of hydrophobic insulin-based peptides using a helping hand strategy; M. M. Disotuar, M. E. Petersen, J. M. Nogueira, M. S. Kay, D. H. Chou; **Org Biomol Chem** 2019; **17**: 1703-1708. <https://doi.org/10.1039/c8ob01212a>
- Chemical synthesis of Shiga toxin subunit B using a next-generation traceless "helping hand" solubilizing tag; J. M. Fulcher, M. E. Petersen, R. J. Giesler, Z. S. Cruz, D. M. Eckert, J. N. Francis, E. M. Kawamoto, M. T. Jacobsen, M. S. Kay; **Org Biomol Chem** 2019; **17**: 10237-10244. <https://doi.org/10.1039/c9ob02012h>
- Mapping the Binding Site of BMS-708163 on gamma-Secretase with Cleavable Photoprobes; N. Gertsik, C. W. Am Ende, K. F. Geoghegan, C. Nguyen, P. Mukherjee, S. Mente, U. Seneviratne, D. S. Johnson Y. M. Li; **Cell Chem Biol** 2017; **24**: 3-8. <https://doi.org/10.1016/j.chembiol.2016.12.006>
- A Helping Hand to Overcome Solubility Challenges in Chemical Protein Synthesis; M. T. Jacobsen, M. E. Petersen, X. Ye, M. Galibert, G. H. Lorimer, V. Aucagne, M. S. Kay; **J Am Chem Soc** 2016; **138**: 11775-82. <https://doi.org/10.1021/jacs.6b05719>
- Compounds and methods for purifying peptides produced by solid phase peptide synthesis; Aucagne V., Delmas A.; CNRS; U.S. Patent No. 9,073,969, 2015
- Cleavable trifunctional biotin reagents for protein labelling, capture and release; Y. Yang, S. H. Verhelst; **Chem Commun** 2013; **49**: 5366-8. <https://doi.org/10.1039/c3cc42076k>
- Investigation on the stability of the Dde protecting group used in peptide synthesis: migration to an unprotected lysine; K. Augustyns, W. Kraas, G. Jung; **J Pept Res** 1998; **51**: 127-33. <https://doi.org/10.1111/j.1399-3011.1998.tb00630.x>
- Full Orthogonality between Dde and Fmoc: The Direct Synthesis of PNA-Peptide Conjugates; J. J. Díaz-Mochón, L. Bialy, M. Bradley; **Org. Lett.** 2004; **7**: 1127-1129. <https://doi.org/10.1021/ol049905y>
- Scope and Limitations of Barbituric and Thiobarbituric Amino Acid Derivatives as Protecting Groups for Solid-Phase Peptide Synthesis: Towards a Green Protecting Group; S. Ramkisson, H. H. Al-Rasheed, K. A. Dahlous, B. G. De La Torre, A. El-Faham, F. Albericio; **Chem. Select** 2021; **6(26)**: 6626-6630. <https://doi.org/10.1002/slct.202101539>



Iris
Biotech













Any Questions or Suggestions?

We are there for you – simply choose one of the numerous possibilities to get in touch!

📞 +49 (0) 9231 97121-0
 📠 +49 (0) 9231 97121-99
 ✉ info@iris-biotech.de
 🌐 www.iris-biotech.de

4. Trifunctional Linkers

		Product details
AAA2190 DAPOA*DCHA 2-(1,3-diazidopropan-2-yloxy)acetic acid dicyclohexylamine CAS-No. 2389064-43-9 Formula C ₅ H ₈ N ₆ O ₃ *C ₁₂ H ₂₃ N Mol. weight 200,16*181,32 g/mol		
BNN1330 DETA(HBH)*2HCl <i>tert</i> -butyl bis(2-aminoethyl)carbamate CAS-No. 1914917-65-9 Formula C ₉ H ₂₁ N ₃ O ₂ *2HCl Mol. weight 203,29*72,92 g/mol		
FAA7570 Fmoc2-DAPOA 2-((1,3-bis((9-fluorenylmethyloxycarbonyl)amino)propan-2-yl)oxy)acetic acid CAS-No. 688350-01-8 Formula C ₃₅ H ₃₂ N ₂ O ₇ Mol. weight 592,64 g/mol		
RL-2710 DACN(Tos,Ns) N-(o-nitrobenzenesulfonyl)-N'-(p-toluenesulfonyl)-4,8-diazacyclononyne CAS-No. 1797508-58-7 Formula C ₂₀ H ₂₁ N ₂ O ₅ S ₂ Mol. weight 463,53 g/mol		
RL-2720 DACN(Tos,Suc-OH) N-succinoyl-N'-(p-toluenesulfonyl)-4,8-diazacyclononyne CAS-No. 2109751-68-8 Formula C ₁₈ H ₂₂ N ₂ O ₅ S Mol. weight 378,44 g/mol		


		Product details
RL-2730	DACN(Tos₂) <i>N,N'</i> -bis(<i>p</i> -toluenesulfonyl)-4,8-diazacyclononyne CAS-No. 1797508-57-6 Formula C ₂₁ H ₂₄ N ₂ O ₄ S ₂ Mol. weight 432,56 g/mol	 
RL-3600	DACN(Ms)*HCl <i>N</i> -(Mesyl)-4,8-diazacyclononyne hydrochloride CAS-No. 2331322-16-6 Formula C ₈ H ₁₄ N ₂ O ₂ S*HCl Mol. weight 202,27*36,46 g/mol	 
RL-3610	DACN(Ms,Ns) <i>N</i> -(Mesyl)- <i>N'</i> -(2-nosyl)-4,8-diazacyclononyne CAS-No. 2411082-25-0 Formula C ₁₄ H ₁₇ N ₃ O ₆ S ₂ Mol. weight 387,43 g/mol	 
RL-2735	DACN(Tos)*HCl <i>N</i> -(<i>p</i> -toluenesulfonyl)-4,8-diazacyclononyne hydrochloride CAS-No. 2331322-18-8 Formula C ₁₄ H ₁₈ N ₂ O ₂ S*HCl Mol. weight 278,37*36,46 g/mol	 
HAA9330	N₃-Gly-Aeg(Fmoc)-OH <i>N</i> -(2-((((9H-fluoren-9-yl)methoxy)carbonyl)amino)ethyl)- <i>N</i> -(2-azidoacetyl)glycine CAS-No. 2606227-07-8 Formula C ₂₁ H ₂₁ N ₅ O ₅ Mol. weight 423,43 g/mol	 

5. Cross-Linkers for other Bio Applications


5.1. Substrates for Fusion (Halo/Snap/Clip)-Tagged Proteins

Site-specific protein labeling is a versatile tool for studying protein function and interaction in living cells. Various peptide-based protein tags have been developed. Herein, we are focusing on HaloTag®, SNAP-tag® and CLIP-Tag™, as well as corresponding substrates offered by Iris Biotech. Those labeling systems enable the specific, covalent attachment of in principle any molecule of choice to a protein of interest.

The HaloTag® (Fig. 28) is a 33 kDa self-labeling protein tag derived from the haloalkane dehalogenase DhaA from *Rhodococcus rhodochrous*. Its active site reacts in a nucleophilic attack with chloroalkane linker substrates to form an irreversible bond in the case of the mutated enzyme. The chloroalkane linker can easily be functionalized with a label of choice, e.g., with a fluorophore or biotin. For the wild-type enzyme, this intermediate would be hydrolyzed, leading to the regeneration of the enzyme.



The SNAP-tag® (Fig. 29) is a 20 kDa self-labeling protein tag based on a modified form of the human O6-alkyl-guanine-DNA-alkyltransferase (hAGT), a DNA repair enzyme. A cysteine residue within the SNAP-tag® undergoes an irreversible reaction with synthetic O6-benzylguanine (BG) derivatives, resulting in a covalent thioether bond. The BG moiety can easily be further functionalized with a label of choice, e.g., fluorophore, biotin, generally without affecting the reaction of the substrate with the SNAP-tag®.



The CLIP-tag™ (Fig. 30) (20 kDa) is a modified version of the SNAP-tag, engineered to react with benzylcytosine (BC) instead of benzylguanine (BG). Thus, properties are similar. CLIP-tag™- and SNAP-tag®-fused proteins can be labeled simultaneously in the same cells.




Fig. 30: Schematic illustration of the CLIP-tag™ reaction.







Tab. 2: Table summarizing main properties of HaloTag®, SNAP-tag® and CLIP-tag™.







	HaloTag	SNAP-tag	CLIP-tag
Origin	Haloalkane dehalogenase (<i>Rhodococcus rhodochrous</i>)	Human O6-alkylguanine-DNA-alkyltransferase	Human O6-alkylguanine-DNA-alkyltransferase
Reactivity	Chloroalkane derivatives	O6-benzylguanine derivatives	Benzylcytosine derivatives
Length	297 amino acids	182 amino acids	182 amino acids
Molecular Weight	33.6 kDa	19.4 kDa	19.4 kDa

Iris Biotech offers a biotin- ([RL-3860 on page 108](#)) as well as an ICG-functionalized ([RL-3830 on page 108](#)) SNAP-tag® substrate, as well as the corresponding CLIP-tag™ suitable ([RL-3840 on page 106](#), [RL-3870 on page 106](#)) derivatives. Biotinylated proteins can for example be selectively isolated based on the high affinity towards avidin representing a useful tool for purification, immobilization, and labeling. Indocyanine green (ICG) is a near-infrared fluorescence imaging dye (absorption maximum 800 nm + slight absorption in the visible range; emission maximum 810 nm) approved by the FDA.

As substrates for the HaloTag® various products are offered, e.g., suitable for further functionalization via Click chemistry.

SNAP-tag® is a registered trademark and CLIP-tag™ a trademark of New England Biolabs, Inc. HaloTag® is a registered trademark to Promega Corporation. HaloTag® Technology is proprietary to Promega Corporation. PROTAC® is a registered trademark of Arvinas Operations, Inc., and is used under license.

			Product details
RL-3180	Halo-PEG(2)-Suc		
4-((2-(2-((6-chlorohexyl)oxy)ethoxy)ethyl)amino)-4-oxobutanoic acid			
CAS-No.	1488363-39-8		
Formula	C ₁₄ H ₂₆ ClNO ₅		
Mol. weight	323,81 g/mol		
RL-3640	Halo-PEG(5)-azide		
1-azido-21-chloro-3,6,9,12,15-pentaoxahenicosane			
CAS-No.	1261238-21-4		
Formula	C ₁₆ H ₃₂ ClN ₃ O ₅		
Mol. weight	381,90 g/mol		
RL-3670	Halo-DBCO		
N-[2-[2-[(6-chlorohexyl)oxy]ethoxy]ethyl]-gamma-oxo-dibenz[b,f]azocine-5(6H)-butanamide			
CAS-No.	1808119-16-5		
Formula	C ₂₉ H ₃₅ ClN ₂ O ₄		
Mol. weight	511,06 g/mol		
RL-3680	Halo-PEG(2)-NH₂*HCl		
12-Chloro-3,6-dioxa-dodecan-1-amine hydrochloride			
CAS-No.	1035373-85-3		
Formula	C ₁₀ H ₂₂ ClNO ₂ *HCl		
Mol. weight	223,74*36,46 g/mol		
RL-3690	Halo-PEG(4)-NH₂*HCl		
18-Chloro-3,6,9,12-tetraoxa-octadecan-1-amine hydrochloride			
CAS-No.	1261238-20-3		
Formula	C ₁₄ H ₃₀ ClNO ₄ *HCl		
Mol. weight	311,85*36,46 g/mol		
RL-3700	Halo-PEG(2)-Azide		
1-Azido-12-chloro-3,6-dioxadodecane			
CAS-No.	2568146-55-2		
Formula	C ₁₀ H ₂₀ ClN ₃ O ₂		
Mol. weight	249,74 g/mol		

			Product details
RL-3710	Halo-PEG(4)-Azide		
1-Azido-18-chloro-3,6,9,12-tetraoxaoctadecane			
CAS-No.	2989398-83-4		
Formula	C ₁₄ H ₂₈ ClN ₃ O ₄		
Mol. weight	337,85 g/mol		
RL-3840	ICG-CLIP		
4-(2-(((1E,3E,5E,7E)-7-(3-(6-(((4-aminopyrimidin-2-yl)oxy)methyl)benzyl)amino)-6-oxohexyl)-1,1-dimethyl-1,3-dihydro-2H-benzo[e]indol-2-ylidene)hepta-1,3,5-trien-1-yl)-1,1-dimethyl-1H-benzo[e]indol-3-ium-3-yl)butane-1-sulfonate			
Formula	C ₅₇ H ₆₂ N ₆ O ₅ S		
Mol. weight	943,22 g/mol		
RL-3870	Biotin-Clip		
N-(4-(((4-aminopyrimidin-2-yl)oxy)methyl)benzyl)-6-(5-(((3aS,4S,6aR)-2-oxohexahydro-1H-thieno[3,4-d]imidazol-4-yl)pentanamido)hexanamide)			
CAS-No.	1004524-73-5		
Formula	C ₂₈ H ₃₉ N ₇ O ₄ S		
Mol. weight	569,73 g/mol		
RL-3835	SNAP-acid		
5-(((4-(((2-amino-9H-purin-6-yl)oxy)methyl)benzyl)amino)-5-oxopentanoic acid			
CAS-No.	881663-34-9		
Formula	C ₁₈ H ₂₀ N ₆ O ₄		
Mol. weight	384,40 g/mol		
RL-3930	Alkyne-SNAP		
N-(4-(((2-amino-9H-purin-6-yl)oxy)methyl)benzyl)pent-4-ynamide			
CAS-No.	1104822-07-2		
Formula	C ₁₈ H ₁₈ N ₆ O ₂		
Mol. weight	350,38 g/mol		
RL-3940	Alkyne-PEG(5)-SNAP		
N-(4-(((2-amino-9H-purin-6-yl)oxy)methyl)benzyl)-4,7,10,13,16-pentaoxanonadec-18-ynamide			
Formula	C ₂₇ H ₃₆ N ₆ O ₇		
Mol. weight	556,62 g/mol		

The Concept of Antibody-Drug Conjugation (ADC)

Permanent Linkers

Cleavable Linkers

Trifunctional Linkers

Cross-Linkers for other Bio Applications

Preparing Carriers for Conjugation


Index

[↑ back to content](#)

RL-3950 Azide-SNAP

N-(4-(((2-amino-9H-purin-6-yl)oxy)methyl)benzyl)-6-azidohexanamide


Formula $C_{19}H_{23}N_5O_2$
Mol. weight 409,45 g/mol



RL-3960 Azide-PEG(4)-SNAP

N-(4-(((2-amino-9H-purin-6-yl)oxy)methyl)benzyl)-1-azido-3,5,7,9-tetraoxadodecan-12-amide


Formula $C_{24}H_{33}N_5O_6$
Mol. weight 543,59 g/mol



RL-3970 Mal-SNAP

N-(4-(((2-amino-9H-purin-6-yl)oxy)methyl)benzyl)-3-(2,5-dioxo-2,5-dihydro-1H-pyrrol-1-yl)propanamide


Formula $C_{20}H_{19}N_7O_4$
Mol. weight 421,42 g/mol



RL-3980 Mal-PEG(4)-SNAP

N-(4-(((2-amino-9H-purin-6-yl)oxy)methyl)benzyl)-1-(3-(2,5-dioxo-2,5-dihydro-1H-pyrrol-1-yl)propanamido)-3,6,9,12-tetraoxapentadecan-15-amide


CAS-No. 1151762-31-0
Formula $C_{31}H_{40}N_8O_9$
Mol. weight 668,71 g/mol



RL-3990 OPSS-SNAP

N-(4-(((2-amino-9H-purin-6-yl)oxy)methyl)benzyl)-3-(pyridin-2-yl)disulfaneyl)propanamide


Formula $C_{21}H_{21}N_7O_2S_2$
Mol. weight 467,57 g/mol






RL-4000 OPSS-PEG(4)-SNAP

N-(4-(((2-amino-9H-purin-6-yl)oxy)methyl)benzyl)-1-(3-(pyridin-2-yl)disulfaneyl)propanamido)-3,6,9,12-tetraoxapentadecan-15-amide

CAS-No. 2278175-10-1
Formula $C_{32}H_{42}N_8O_7S_2$
Mol. weight 714,86 g/mol



		Product details
RL-4010	DBCO-SNAP	 
Formula	$C_{34}H_{31}N_7O_3$	
Mol. weight	585,67 g/mol	
RL-3830	ICG-SNAP	 
4-(2-(((1E,3E,5E,7E)-7-(3-(6-(((2-amino-9H-purin-6-yl)oxy)methyl)benzyl)amino)-6-oxohexyl)-1,1-dimethyl-1,3-dihydro-2H-benzo[e]indol-2-ylidene)hepta-1,3,5-trien-1-yl)-1,1-dimethyl-1H-benzo[e]indol-3-ium-3-yl)butane-1-sulfonate		
Formula	$C_{58}H_{62}N_8O_5S$	
Mol. weight	983,25 g/mol	
RL-3860	Biotin-SNAP	 
N-(4-(((2-amino-9H-purin-6-yl)oxy)methyl)benzyl)-6-(5-(((3aS,4S,6aR)-2-oxohexahydro-1H-thieno[3,4-d]imidazol-4-yl)pentanamido)hexanamide		
CAS-No.	471918-16-8	
Formula	$C_{29}H_{39}N_9O_4S$	
Mol. weight	609,75 g/mol	

References:


- Site-specific protein labeling with SNAP-Tags; N. B. Cole; **Curr Protoc Protein Sci.** 2013; **73(30)**: 1-30.
<https://doi.org/10.1002/0471140864.ps3001s73>
- A general method for the covalent labeling of fusion proteins with small molecules in vivo; A. Keppler, S. Gendreizig, T. Gronemeyer, H. Pick, H. Vogel, K. Johnsson; **Nat. Biotechnol.** 2003; **21**: 86-89. <https://doi.org/10.1038/nbt765>
- HaloTag: A Novel Protein Labeling Technology for Cell Imaging and Protein Analysis; g: v: Los; L. P. Encell, M. G. McDougall, D. D. Hartzell, N. Karassina, C. Zimprich, M. G. Wood, R. Learish, R. F. Ohana, M. Urh, D. Simpson, J. Mendez, K. Zimmerman, P. Otto, G. Vidugiris, J. Zhu, A. Darzins, D. H. Klaubert, R. F. Bulleit, K. V. Wood; **ACS Chem. Biol.** 2008; **3(6)**: 373-382. <https://doi.org/10.1021/cb800025k>
- Directed evolution of O6-alkylguanine-DNA alkyltransferase for efficient labeling of fusion proteins with small molecules in vivo; A. Juillerat, T. Gronemeyer, A. Keppler, S. Gendreizig, H. Pick, H. Vogel, K. Johnsson; **Chem. Biol.** 2003; **10(4)**: 313-317. [https://doi.org/10.1016/s1074-5521\(03\)00068-1](https://doi.org/10.1016/s1074-5521(03)00068-1)
- Site-specific, Covalent Labeling of Recombinant Antibody Fragments via Fusion to an Engineered Version of 6-O-Alkylguanine DNA Alkyltransferase; F. Kampmeier, M. Ribbert, T. Nachreiner, S. Dembski, F. Beaufils, A. Brecht, S. Barth; **Bioconjugate Chem.** 2009; **20(5)**: 1010-1015. <https://doi.org/10.1021/bc9000257>
- SNAP-Tag Technology: A Useful Tool to Determine Affinity Constants and Other Functional Parameters of Novel Antibody Fragments; J. Niesen, M. Sack, M. Seidel, R. Fendel, S. Barth, R. Fischer, C. Stein; **Bioconjugate Chem.** 2016; **27(8)**: 1931-1941. <https://doi.org/10.1021/acs.bioconjchem.6b00315>
- The clinical use of indocyanine green as a near-infrared fluorescent contrast agent for image-guided oncologic surgery. Boudewijn E. Schaafsma MD, J. Sven D. Mieog MD, Merlijn Hutteman MSc, Joost R. van der Vorst MD, Peter J.K. Kuppen PhD, Clemens W.G.M. Lwik PhD, John V. Frangioni MD, PhD, Cornelis J.H. van de Velde MD, PhD, Alexander L. Vahrmeijer MD PhD; **Surg. Oncol.** 2011; **104**: 323-332. <https://doi.org/10.1002/jso.21943>
- Degradation kinetics of indocyanine green in aqueous solution. Vishal Saxena, Mostafa Sadoqi, Jun Shao; **J. Pharm. Sci.** 2003; **92**: 2090-2097. <https://doi.org/10.1002/jps.10470>

[↑ back to content](#)

- Stability assessment of indocyanine green within dextran-coated mesocapsules by absorbance spectroscopy. Mohammad Abbas Yaseen; Jie Yu; Michael S. Wong; Bahman Anvari; **Journal of Biomedical Optics** 2007; **12(6)**: 064031. <https://doi.org/10.1117/1.2821423>
- Imaging proteins inside cells with fluorescent tags; G. Crivat, J. W. Taraska; **Trends Biotechnol** 2012; **30(1)**: 8-16. <https://doi.org/10.1016/j.tibtech.2011.08.002>
- Visualizing Biochemical Activities in Living Cells through Chemistry; L. Reymond, G. Lukinavicius, K. Umezawa, D. Maurel, M. A. Brun, A. Masharina, K. Bojkowska, B. Mollwitz, A. Schena, R. Griss, K. Johnsson; **CHIMICA Int. J. Chem.** 2011; **65(11)**: 868-871. <https://doi.org/10.2533/chimia.2011.868>
- Cell Penetration Profiling Using the Chloroalkane Penetration Assay; L. Peraro, K. L. Deprey, M. K. Moser, Z. Zou, H. L. Ball, B. Levine, J. A. Kritzer; **J. Am. Chem. Soc.** 2018; **140(36)**: 11360-11369. <https://doi.org/10.1021/jacs.8b06144>
- HaloTag technology: a versatile platform for biomedical applications; C. G. England, H. Luo, W. Cai; **Bioconjug Chem** 2015; **26**: 975-86. <https://doi.org/10.1021/acs.bioconjchem.5b00191>
- HaloTag: a novel protein labeling technology for cell imaging and protein analysis; G. V. Los, L. P. Encell, M. G. McDougall, D. D. Hartzell, N. Karassina, C. Zimprich, M. G. Wood, R. Learish, R. F. Ohana, M. Urh, D. Simpson, J. Mendez, K. Zimmerman, P. Otto, G. Vidugiris, J. Zhu, A. Darzins, D. H. Klaubert, R. F. Bulleit, K. V. Wood; **ACS Chem Biol** 2008; **3**: 373-382. <https://doi.org/10.1021/cb800025k>
- Self-labelling enzymes as universal tags for fluorescence microscopy, super-resolution microscopy and electron microscopy; V. Liss, B. Barlag, M. Nietschke, M. Hensel; **Scientific Reports** 2016; 5. <https://doi.org/10.1038/srep17740>
- Snap-, CLIP- and Halo-Tag Labelling of Budding Yeast Cells; F. Stagge, G. Y. Mitronova, V. N. Belov, C. A. Wurm, S. Jakobs; **PLoS ONE** 8(10): e78745. <https://doi.org/10.1371/journal.pone.0078745>

5.2. Specific His Tag Acylation

Gly-His tag for selective α-amine acylation:



Lys-His tag for selective ε-amine acylation:




Fig. 31: Schematic illustration of Gly-His tag acylation and Lys-His tag acylation. Functionalized esters can be used for the introduction of the desired moieties, e.g., azides, biotin.

A frequently used tag for protein modification is the so-called His tag. His tags, or polyhistidine tags, comprise a consecutive series of six to ten histidine residues. His tags are widely used for protein purification by immobilized-metal ion affinity chromatography, allowing to extract a protein of interest from thousands of other proteins present in a cell lysate. However, as an inadvertent side-reaction, it is reported that His tagged proteins can undergo N-terminal acylation with gluconolactone (GDL). The discovery of this reaction by Geoghegan and coworkers triggered the inspiration for the development of His tag-based peptide segments for selective acylations. Thus, Jensen *et al.* developed two methods that use poly-histidine sequences to direct the highly selective acylation of proteins, either at the N-terminus or at a specific lysine residue.

The highly selective and efficient N-terminal acylation of proteins is based on a Gly-His₆ segment (Gly-His tag) complemented by the use of 4-methoxy phenyl esters as finely tuned acylating agents, resulting in stable conjugates. Other acylating agents, e.g., lactons, thioesters, N-hydroxysuccinimide ester were also tested but gave only limited N-terminal acylation or low selectivity.

General Conjugation Protocol:

A 35 μ M solution of GH6-protein in 200 mM HEPES buffer at pH 7.5 is incubated with 40 equiv. of azido-acetyl 4-methoxyphenylester for 24 h at 4 °C. The formation of the mono-functionalized product can be observed by ESI-MS and can reach 70% to 90% conversion. A higher conversion rate can be achieved by the addition of two aliquots of 10 equiv. of the acylating agent in the course of the next 48 h.

For the acylation of lysine, Jensen *et al.* developed the peptide sequence His_n-Lys-His_m (Lys-His tag) that directs the acylation of the designated lysine N- ϵ amine under mild conditions and with high selectivity over native lysine residues.

Both methods provide highly selective acylation, yield stable conjugates, and do not require the use of metal ions. In detail, referring to literature, yields for the protein modification with Gly-His are typically 60–80% mono-acylation with 1–5% over-acylation, while for Lys-His it is typically 50–70% monoacylation with 1–8% over-acylation. The Gly-His tag as well as the Lys-His tag maintain the capacity for immobilized metal ion affinity chromatography and have been shown robust for the attachment of azides, fluorophores, and biotin to different proteins, including antibodies.

Mechanistic studies indicate that the very high selectivity of the His tag acylation is based on specific base catalysis, in which a histidine side-chain assists deprotonation during the direct acylation of the glycine α -amine (Fig. 32). The ester preferentially reacts with assistance from histidine side-chain imidazoles since they are not protonated ($pK_a \sim 6.0$) at the pH of the reaction, in contrast to the N-terminal α -amine ($pK_a \sim 7.6$ –8.0) and lysine side-chains ($pK_a \sim 10.5$). The presence of the additional five His residues in the His tag may serve to modulate the basicity of the imidazole nitrogen of the catalytic residue. A recent study has shown that the pK_a values of individual histidine side-chains in a His₆-tag span a range from 4.8–7.5.






Fig. 32: Imidazole rings of neighboring histidines in a His tag catalyze the acylation of a glycyl N-terminus via a base catalyzed mechanism.

		Product details
RL-3010	N ₃ Ac-OPhOMe	
4-Methoxyphenyl 2-azidoacetate		
CAS-No.	2546513-31-7	
Formula	C ₉ H ₉ N ₃ O ₃	
Mol. weight	207,19 g/mol	
		
RL-3100	Biotin-AEEA-OPhOMe	
2-(2-(2-(Biotinamido)ethoxy)ethoxy)acetic acid 4-methoxyphenyl ester		
CAS-No.	2546513-67-9	
Formula	C ₂₃ H ₃₃ N ₃ O ₇ S	
Mol. weight	495,59 g/mol	
		

References:

- Selective Acylation of Proteins at Gly and Lys in His tags; K. J. Jensen, M. B. Thygesen, K. K. Sørensen, S. Wu, T. Treiberg, S. Schoffelen; **ChemBioChem** 2022. <https://doi.org/10.1002/cbic.202200359>
- Site-specific covalent labeling of His tag fused proteins with N-acyl-N-alkyl sulfonamide reagent; V. Thimaradka, J. H. Oh, C. Heroven, A. R. Aricescu, M. Yuzaki, T. Tamura, I. Hamachi; **Bioorg. Med. Chem.** 2021; **15(30)**: 115947. <https://doi.org/10.1016/j.bmc.2020.115947>
- Selective N-terminal acylation of peptides and proteins with a Gly-His tag sequence; M. C. Martos-Maldonado, C. T. Hjuler, K. K. Sorensen, M. B. Thygesen, J. E. Rasmussen, K. Villadsen, S. R. Midtgaard, S. Kol, S. Schoffelen, K. J. Jensen; **Nat Commun** 2018; **9**: 3307. <https://doi.org/10.1038/s41467-018-05695-3>
- Spontaneous alpha-N-6-phosphogluconoylation of a "His tag" in Escherichia coli: the cause of extra mass of 258 or 178 Da in fusion proteins; K. F. Geohegan, H. B. Dixon, P. J. Rosner, L. R. Hoth, A. J. Lanzetti, K. A. Borzilleri, E. S. Marr, L. H. Pezzullo, L. B. Martin, P. K. LeMotte, A. S. McCol, A. V. Kamath, J. G. Stroh; **Anal Biochem.** 1999; **267(1)**: 169-84. <https://doi.org/10.1006/abio.1998.2990>




You are interested in autocatalytic His tags
for the chemical modification of proteins?

Watch the recording of our workshop!




5.3. Bifunctional Protein Cross-Linkage

		Product details
RL-2770	BSSS	
Bis(sulfosuccinimidyl) suberate sodium salt		
CAS-No.	127634-19-9	
Formula	$C_{16}H_{18}N_2Na_2O_{14}S_2$	
Mol. weight	572,43 g/mol	



This molecule ([Fig. 33](#)) carries amino reactive sulfo-NHS esters on both ends and is a water-soluble, homo-bi functional protein cross-linker (spacer length: 11.4 Å). Due to its water solubility, conjugation reactions can conveniently take place at physiological conditions. This 8-atom spacer is non-cleavable and the molecule is not cell membrane permeable. It can be used to prepare antibody-protein conjugates, for crosslinking cell surface proteins, and for covalently binding an antibody to an immobilized Protein A or Protein G resin.



General BSSS Cross-Linking Protocol:

1. Allow vial of BSSS to fully equilibrate to ambient temperature before opening to prevent condensation inside the vial (BSSS is moisture-sensitive).
2. Immediately before use, prepare a 50 mM solution of BSSS, by dissolving 10 mg BSSS in 350 mL of 25 mM sodium phosphate, pH 7.4 (do not use amine containing buffers for the conjugation reaction).
3. Add BSSS solution (20-fold excess cross-linker to protein) to the protein sample so that the final concentration is between 0.5 to 5 mM.
4. Allow the sample to react at room temperature for 45 minutes to 1 hour. Allow slightly longer if reaction must be done on ice (the reaction rate is only slightly slower at low temperatures).
5. Quench any unreacted BSSS with 25 mM to 60 mM Tris and allow to react for 10-15 minutes at room temperature.
6. Desalt sample to remove unreacted BSSS, i.e., by gel filtration, dialysis, etc.

References:

- *Isotope-tagged cross-linking reagents. A new tool in mass spectrometric protein interaction analysis*; D. R. Muller, P. Schindler, H. Towbin, U. Wirth, H. Voshol, S. Hoving, M. O. Steinmetz; *Anal Chem* 2001; **73**: 1927-34. <https://doi.org/10.1021/ac001379a>
- *The Golgi-associated hook3 protein is a member of a novel family of microtubule-binding proteins*; J. H. Walenta, A. J. Didier, X. Liu, H. Kramer; *J Cell Biol* 2001; **152**: 923-34. <https://doi.org/10.1083/jcb.152.5.923>
- *Rapamycin potentiates transforming growth factor beta-induced growth arrest in nontransformed, oncogene-transformed, and human cancer cells*; B. K. Law, A. Chytil, N. Dumont, E. G. Hamilton, M. E. Waltner-Law, M. E. Aakre, C. Covington, H. L. Moses; *Mol Cell Biol* 2002; **22**: 8184-98. <https://doi.org/10.1128/mcb.22.23.8184-8198.2002>
- *Chemical cross-linking and mass spectrometry for mapping three-dimensional structures of proteins and protein complexes*; A. Sinz; *J Mass Spectrom* 2003; **38**: 1225-37. <https://doi.org/10.1002/jms.559>
- *Mapping low-resolution three-dimensional protein structures using chemical cross-linking and Fourier transform ion-cyclotron resonance mass spectrometry*; G. H. Dihazi, A. Sinz; *Rapid Commun Mass Spectrom* 2003; **17**: 2005-14. <https://doi.org/10.1002/rcm.1144>
- *Mapping the topology and determination of a low-resolution three-dimensional structure of the calmodulin-melittin complex by chemical cross-linking and high-resolution FTICRMS: direct demonstration of multiple binding modes*; D. M. Schulz, C. Ihling, G. M. Clore, A. Sinz; *Biochemistry* 2004; **43**: 4703-15. <https://doi.org/10.1021/bi036149f>
- *Selective inactivation of adrenomedullin over calcitonin gene-related peptide receptor function by the deletion of amino acids 14-20 of the mouse calcitonin-like receptor*; D. Koller, L. M. Ittner, R. Muff, K. Husmann, J. A. Fischer, W. Born; *J Biol Chem* 2004; **279**: 20387-91. <https://doi.org/10.1074/jbc.M313058200>
- *Mactinin, a fragment of cytoskeletal alpha-actinin, is a novel inducer of heat shock protein (Hsp)-90 mediated monocyte activation*; S. D. Luikart, A. Panoskaltis-Mortari, T. Hinkel, R. T. Perri, K. Gupta, T. R. Oegema, P. Gupta; *BMC Cell Biol* 2009; **10**: 60. <https://doi.org/10.1186/1471-2121-10-60>
- *The program for processing newly synthesized histones H3.1 and H4*; E. I. Campos, J. Fillingham, G. Li, H. Zheng, P. Voigt, W. H. Kuo, H. Seepany, Z. Gao, L. A. Day, J. F. Greenblatt, D. Reinberg; *Nat Struct Mol Biol* 2010; **17**: 1343-51. <https://doi.org/10.1038/nsmb.1911>
- *Sar1 assembly regulates membrane constriction and ER export*; K. R. Long, Y. Yamamoto, A. L. Baker, S. C. Watkins, C. B. Coyne, J. F. Conway, M. Aridor; *J Cell Biol* 2010; **190**: 115-28. <https://doi.org/10.1083/jcb.201004132>
- *Mps1 directs the assembly of Cdc20 inhibitory complexes during interphase and mitosis to control M phase timing and spindle checkpoint signaling*; J. Maciejowski, K. A. George, M. E. Terret, C. Zhang, K. M. Shokat, P. V. Jallepalli; *J Cell Biol* 2010; **190**: 89-100. <https://doi.org/10.1083/jcb.201001050>
- *Protection against protein aggregation by alpha-crystallin as a mechanism of preconditioning*; J. E. Ferns, C. S. Theisen, E. E. Fibuch, N. W. Seidler; *Neurochem Res* 2012; **37**: 244-52. <https://doi.org/10.1007/s11064-011-0601-4>
- *Identification of IGPR-1 as a novel adhesion molecule involved in angiogenesis*; N. Rahimi, K. Rezazadeh, J. E. Mahoney, E. Hartsough, R. D. Meyer; *Mol Biol Cell* 2012; **23**: 1646-56. <https://doi.org/10.1091/mbc.E11-11-0934>
- *TWEAK-independent Fn14 self-association and NF-kappaB activation is mediated by the C-terminal region of the Fn14 cytoplasmic domain*; S. A. Brown, E. Cheng, M. S. Williams, J. A. Winkles; *PLoS One* 2013; **8**: e65248. <https://doi.org/10.1371/journal.pone.0065248>

- ZRF1 controls the retinoic acid pathway and regulates leukemogenic potential in acute myeloid leukemia; S. Demajo, I. Uribealago, A. Gutierrez, C. Ballare, S. Capdevila, M. Roth, J. Zuber, J. Martin-Caballero, L. Di Croce; **Oncogene** 2014; **33**: 5501-10. <https://doi.org/10.1038/onc.2013.501>
- Ubiquitin Associates with the N-Terminal Domain of Nerve Growth Factor: The Role of Copper(II) Ions; V. Lanza, A. Travaglia, G. Maltieri, R. Fattorusso, G. Grasso, G. Di Natale, V. Zito, G. Arena, D. Milardi, E. Rizzarelli; **Chemistry** 2016; **22**: 17767-17775. <https://doi.org/10.1002/chem.201603650>
- Potentiation of Surface Stability of AMPA Receptors by Sulfhydryl Compounds: A Redox-Independent Effect by Disrupting Palmitoylation; J. Han, H. Zhang, S. Wang, J. Zhou, Y. Luo, L. H. Long, Z. L. Hu, F. Wang, J. G. Chen, P. F. Wu; **Neurochem Res** 2016; **41**: 2890-2903. <https://doi.org/10.1007/s11064-016-2006-x>
- Dramatic Domain Rearrangements of the Cyanobacterial Orange Carotenoid Protein upon Photoactivation; H. Liu, H. Zhang, G. S. Orf, Y. Lu, J. Jiang, J. D. King, N. R. Wolf, M. L. Gross, R. E. Blankenship; **Biochemistry** 2016; **55**: 1003-9. <https://doi.org/10.1021/acs.biochem.6b00013>
- Nuclear Speckle-related Protein 70 Binds to Serine/Arginine-rich Splicing Factors 1 and 2 via an Arginine/Serine-like Region and Counteracts Their Alternative Splicing Activity; C. H. Kim, Y. D. Kim, E. K. Choi, H. R. Kim, B. R. Na, S. H. Im, C. D. Jun; **J Biol Chem** 2016; **291**: 6169-81. <https://doi.org/10.1074/jbc.M115.689414>
- PGL germ granule assembly protein is a base-specific, single-stranded RNase; S. T. Aoki, A. M. Kershner, C. A. Bingman, M. Wickens, J. Kimble; **Proc Natl Acad Sci U S A** 2016; **113**: 1279-84. <https://doi.org/10.1073/pnas.1524400113>

5.4. Triazolecarbaldehydes for Selective N-Terminal Protein Modification

The selective modification of the N-terminus (i.e., an α -amino group) is frequently required when, e.g., payloads are attached to antibodies to generate ADCs, or proteins need to be precisely labeled, e.g., with fluorescent dyes or biotin. Besides, the N-terminus is an attractive site for modification as it represents a unique position in proteins and is typically not involved in folding. Classical reagents for the conjugation of molecules to amines are NHS and TFP esters, iso(thio)cyanates and the reductive amination with aldehydes. However, all of them react rather unspecific.

Herein, we introduce a novel reagent for the specific N-terminal labeling of proteins and peptides: 1-(4-nitrophenyl)-1H-1,2,3-triazole-4-carbaldehyde (1-NP-triazole-4-CHO) ([RL-8650 on page 116](#)).

This nitrophenol-functionalized triazole aldehyde reacts with primary amines under mild conditions and high conversion rates with reported yields of 90% and more. This allows to introduce labels, e.g., for imaging. The coupling reaction may be performed in aqueous solution buffered with 10 mM potassium phosphate or MOPS at pH 7.5, by adding the aldehyde dissolved in DMF and incubating for 4-16 hours. Work-up may be performed by chromatography.




Fig. 34: Activation of a primary amine (e.g., carrying a label = green star) with 1-NP-triazole-4-carbaldehyde. The amine binds to the aldehyde, forming an imine, which undergoes a rearrangement thus integrating the amine nitrogen as part of the triazole. Then, *p*-nitrophenylaniline is eliminated and an activated aldehyde triazole is formed which then can be coupled to the *N*-terminus of a peptide.

With the success of Click chemistry the triazole moiety gained more and more presence in drug design and discovery. Several triazole derivatives with interesting biological activities are already reported in literature making it an attractive structural motive for peptide modification.

The activated triazole aldehyde – bearing the desired label (illustrated as green star in the scheme below) – can be specifically conjugated to the *N*-terminus of a peptide forming a *N*-terminal 4-imidazolinone. The reaction may be performed overnight, in an aqueous buffer at physiological pH, e.g., 10 mM phosphate buffer, at 37 °C.






Fig. 35: Coupling of the aldehyde-functionalized triazole reagent to the *N*-terminus of a peptide forming an imidazolinone.

		Product details
RL-8650	1-NP-triazole-4-CHO	
1-(4-nitrophenyl)triazole-4-carbaldehyde		
CAS-No.	113934-26-2	
Formula	C ₉ H ₆ N ₄ O ₃	
Mol. weight	218,17 g/mol	
		
RL-8655	4-(Diethoxymethyl)-1-NP-triazole	
4-(Diethoxymethyl)-1-(4-nitrophenyl)-1H-1,2,3-triazole		
CAS-No.	1012081-58-1	
Formula	C ₁₃ H ₁₆ N ₄ O ₄	
Mol. weight	292,29 g/mol	
		

References:

- *Tandem synthesis of 1-formyl-1,2,3-triazoles*; J. T. Fletcher, J. A. Christensen, E. M. Villa; **Tetrahedron Lett.** 2017; **58(47)**: 4450-4454. <https://doi.org/10.1016/j.tetlet.2017.10.023>
- *Synthesis and antibiotic activity of a small molecules library of 1,2,3-triazole derivatives*; M. Aufort, J. Herscovici, P. Bouhours, N. Moreau, C. Girard; **Bioorg Med Chem Lett.** 2008; **18(3)**: 1195-1198. <https://doi.org/10.1016/j.bmcl.2007.11.111>
- *Triazolecarbaldehyde reagents for one-step N-terminal protein modification*, A. Onoda, N. Inouse, E. Sumiyoshi, T. Hayashi; **ChemBioChem.** 2019; **21(9)**: 1274-1278. <https://doi.org/10.1002/cbic.201900692>


5.5. Proteolysis Targeting Chimeras (PROTACs®)

Targeted protein degradation ([Fig. 36](#)) via proteolysis-targeting chimeras (PROTACs) is an emerging attempt to cure diseases caused by the irregular expression of certain disease-causing proteins. Such protein degraders act as bifunctional linkers and allow to feed the protein of interest (POI) to the cell's Ubiquitin-Proteasome system, thus, to eliminate the malexpressed proteins. These PROTACs consist of three components: one ligand with high affinity for E3 ubiquitin ligase, another one with high affinity for the protein of interest (POI) and an appropriate cross-linker joining both ligands. This linker can also be used to increase the solubility, if needed, e.g., by incorporation of PEGs. The resulting proximity of both, the recruited POI and the E3 ligase, allows the polyubiquitination of the POI by the E3 associated E2 enzyme. This leads to a labeling of the POI for degradation through the proteasome.

PROTAC® is a registered trademark of Arvinas Operations, Inc., and is used under license.


Mode of action:







1. A cross-linker unites the POI ligand and E3 ligase ligand = PROTAC.
2. The three-component PROTAC recruits the POI and the E2-associated E3 ligase via the respective ligands = Ternary complex.
3. Several Ubiquitins are added to lysine residues of the POI = Polyubiquitination.
4. The ubiquitinated POI is degraded by the proteasome.






To construct a suitable PROTAC, we provide a variety of E3 ubiquitin ligase ligands in combination with linkers of various length and an elective amino-, carboxyl-, click- or thiol-reactive end ("Partial PROTACs").



E3-Ligase Ligands & Negative Controls







		Product details
PTC1000	Pomalidomide	
1,3-dioxo-2-(2,6-dioxopiperidin-3-yl)-4-aminoisoindoline,		
CAS-No.	19171-19-8	
Formula	C ₁₃ H ₁₁ N ₃ O ₄	
Mol. weight	273,24 g/mol	
		





		Product details
PTC1010	N-Methylated pomalidomide	
4-Amino-2-(1-methyl-2,6-dioxopiperidin-3-yl)isoindoline-1,3-dione		
CAS-No.	1352827-50-9	 
Formula	C ₁₄ H ₁₃ N ₃ O ₄	
Mol. weight	287,27 g/mol	
PTC1020	Lenalidomide	
1-Oxo-4-amino-2-(2,6-dioxopiperidin-3-yl)isoindole		
CAS-No.	191732-72-6	 
Formula	C ₁₃ H ₁₃ N ₃ O ₃	
Mol. weight	259,26 g/mol	
PTC1030	(±)-Thalidomide	
(±)-2-(2,6-Dioxo-3-piperidiny)-1H-isoindole-1,3(2H)-dione		
CAS-No.	50-35-1	 
Formula	C ₁₃ H ₁₀ N ₂ O ₄	
Mol. weight	258,23 g/mol	
PTC1040	(S,R,S)-AHPC hydrochloride	
(2S,4R)-1-((S)-2-Amino-3,3-dimethylbutanoyl)-4-hydroxy-N-(4-(4-methylthiazol-5-yl)benzyl)pyrrolidine-2-carboxamide hydrochloride		
CAS-No.	1448189-80-7	 
Formula	C ₂₂ H ₃₀ N ₄ O ₃ S*xHCl	
Mol. weight	430,56 (free base) g/mol	
PTC1050	(S,S,S)-AHPC hydrochloride	
(2S,4S)-1-((S)-2-Amino-3,3-dimethylbutanoyl)-4-hydroxy-N-(4-(4-methylthiazol-5-yl)benzyl)pyrrolidine-2-carboxamide hydrochloride		
CAS-No.	2115897-23-7	 
Formula	C ₂₂ H ₃₀ N ₄ O ₃ S*xHCl	
Mol. weight	430,56 (free base) g/mol	
PTC1060	VH298	
(2S,4R)-1-((S)-2-(1-Cyanocyclopropanecarboxamido)-3,3-dimethylbutanoyl)-4-hydroxy-N-(4-(4-methylthiazol-5-yl)benzyl)pyrrolidine-2-carboxamide		
CAS-No.	2097381-85-4	 
Formula	C ₂₇ H ₃₃ N ₅ O ₄ S	
Mol. weight	523,65 g/mol	

		Product details
PTC1070	Nimbolide	
CAS-No.	25990-37-8	
Formula	C ₂₇ H ₃₀ O ₇	
Mol. weight	466,52 g/mol	
		
PTC1080	Nutlin-3	
(±)-4-[4,5-Bis(4-chlorophenyl)-2-(2-isopropoxy-4-methoxy-phenyl)-4,5-dihydro-imidazole-1-carbonyl]-piperazin-2-one		
CAS-No.	548472-68-0	
Formula	C ₃₀ H ₃₀ Cl ₂ N ₄ O ₄	
Mol. weight	581,49 g/mol	
		
PTC1090	Nutlin-3a	
(-)-4-(4,5-Bis(4-chlorophenyl)-2-(2-isopropoxy-4-methoxyphenyl)-4,5-dihydro-1H-imidazole-1-carbonyl)piperazin-2-one		
CAS-No.	675576-98-4	
Formula	C ₃₀ H ₃₀ Cl ₂ N ₄ O ₄	
Mol. weight	581,49 g/mol	
		


Amino-Reactive Partial PROTACs







		Product details
PTC1100	Pomalidomide-C3-COOH	
4-((2-(2,6-Dioxopiperidin-3-yl)-1,3-dioxoisindolin-4-yl)amino)butanoic acid		
CAS-No.	2225940-47-4	
Formula	C ₁₇ H ₁₇ N ₃ O ₆	
Mol. weight	359,33 g/mol	
		
PTC1110	Pomalidomide-C6-COOH	
7-((2-(2,6-Dioxopiperidin-3-yl)-1,3-dioxoisindolin-4-yl)amino)heptanoic acid		
CAS-No.	2225940-50-9	
Formula	C ₂₀ H ₂₃ N ₃ O ₆	
Mol. weight	401,41 g/mol	
		

		Product details
PTC1120	Pomalidomide-C9-COOH	
10-((2-(2,6-Dioxopiperidin-3-yl)-1,3-dioxoisindolin-4-yl)amino)decanoic acid		
CAS-No.	2243000-24-8	 
Formula	C ₂₃ H ₂₉ N ₃ O ₆	
Mol. weight	443,5 g/mol	
PTC1130	Pomalidomide-PEG2-butyl COOH	
7-(2-(2-((2-(2,6-Dioxopiperidin-3-yl)-1,3-dioxoisindolin-4-yl)amino)-2-oxoethoxy)ethoxy)heptanoic acid		
CAS-No.	2421216-99-9	 
Formula	C ₂₄ H ₂₉ N ₃ O ₉	
Mol. weight	503,5 g/mol	
PTC1140	Pomalidomide-PEG6-butyl COOH	
1-((2-(2,6-Dioxopiperidin-3-yl)-1,3-dioxoisindolin-4-yl)amino)-1-oxo-3,6,9,12,15,18-hexaoxapentacosan-25-oic acid		
Formula	C ₃₂ H ₄₅ N ₃ O ₁₃	 
Mol. weight	679,71 g/mol	
PTC1150	Pomalidomide-PEG1-COOH	
3-(2-((2-(2,6-Dioxopiperidin-3-yl)-1,3-dioxoisindolin-4-yl)amino)ethoxy)propanoic acid		
CAS-No.	2139348-60-8	 
Formula	C ₁₈ H ₁₉ N ₃ O ₇	
Mol. weight	389,36 g/mol	
PTC1160	Pomalidomide-PEG2-COOH	
3-(2-(2-((2-(2,6-Dioxopiperidin-3-yl)-1,3-dioxoisindolin-4-yl)amino)ethoxy)ethoxy)propanoic acid		
CAS-No.	2140807-17-4	 
Formula	C ₂₀ H ₂₃ N ₃ O ₈	
Mol. weight	433,42 g/mol	
PTC1170	Pomalidomide-PEG3-COOH	
3-(2-(2-(2-((2-(2,6-Dioxopiperidin-3-yl)-1,3-dioxoisindolin-4-yl)amino)ethoxy)ethoxy)ethoxy)propanoic acid		
CAS-No.	2138440-82-9	 
Formula	C ₂₂ H ₂₇ N ₃ O ₉	
Mol. weight	477,46 g/mol	

		Product details
PTC1180	Pomalidomide-PEG4-COOH	
1-((2-(2,6-Dioxopiperidin-3-yl)-1,3-dioxisoindolin-4-yl)amino)-3,6,9,12-tetraoxapentadecan-15-oic acid		
CAS-No.	2138440-81-8	
Formula	C ₂₄ H ₃₁ N ₃ O ₁₀	
Mol. weight	521,52 g/mol	
		
PTC1190	Pomalidomide-PEG5-COOH	
1-((2-(2,6-Dioxopiperidin-3-yl)-1,3-dioxisoindolin-4-yl)amino)-3,6,9,12,15-pentaoxaoctadecan-18-oic acid		
CAS-No.	2139348-63-1	
Formula	C ₂₆ H ₃₅ N ₃ O ₁₁	
Mol. weight	565,57 g/mol	
		
PTC1200	Pomalidomide-PEG6-COOH	
1-((2-(2,6-Dioxopiperidin-3-yl)-1,3-dioxisoindolin-4-yl)amino)-3,6,9,12,15,18-hexaoxahenicosan-21-oic acid		
CAS-No.	2225148-49-0	
Formula	C ₂₈ H ₃₉ N ₃ O ₁₂	
Mol. weight	609,62 g/mol	
		
PTC1220	(S,R,S)-AHPC-PEG2-butyl COOH	
(S,R,S)-AHPC-2-2-6-acid,7-(2-(2-(((S)-1-((2S,4R)-4-Hydroxy-2-((4-(4-methylthiazol-5-yl)benzyl)carbamoyl)pyrrolidin-1-yl)-3,3-dimethyl-1-oxobutan-2-yl)amino)-2-oxoethoxy)ethoxy)heptanoic acid		
CAS-No.	2421187-89-3	
Formula	C ₃₃ H ₄₈ N ₄ O ₈ S	
Mol. weight	660,82 g/mol	
		

Carboxy-Reactive Partial PROTACs


		Product details
PTC1230	Pomalidomide-PEG1-NH ₂ hydrochloride	
2-(2-Aminoethoxy)-N-(2-(2,6-dioxopiperidin-3-yl)-1,3-dioxisoindolin-4-yl)acetamide hydrochloride		
CAS-No.	2380273-67-4	
Formula	C ₁₇ H ₁₈ N ₄ O ₆ ·xHCl	
Mol. weight	374,35 (free base) g/mol	
		

		Product details
PTC1240 Pomalidomide-PEG2-NH₂ hydrochloride 2-(2-(2-Aminoethoxy)ethoxy)-N-(2-(2,6-dioxopiperidin-3-yl)-1,3-dioxoisindolin-4-yl)acetamide hydrochloride CAS-No. 2380273-73-2 Formula C ₁₉ H ₂₂ N ₄ O ₇ ·xHCl Mol. weight 418,40 (free base) g/mol		
PTC1250 Pomalidomide-PEG3-NH₂ hydrochloride 2-(2-(2-(2-Aminoethoxy)ethoxy)ethoxy)-N-(2-(2,6-dioxopiperidin-3-yl)-1,3-dioxoisindolin-4-yl)acetamide hydrochloride CAS-No. 2380273-75-4 Formula C ₂₁ H ₂₆ N ₄ O ₈ ·xHCl Mol. weight 462,45 (free base) g/mol		
PTC1260 Pomalidomide-PEG4-NH₂ hydrochloride 14-Amino-N-(2-(2,6-dioxopiperidin-3-yl)-1,3-dioxoisindolin-4-yl)-3,6,9,12-tetraoxatetradecanamide hydrochloride CAS-No. 2331259-45-9 Formula C ₂₃ H ₃₀ N ₄ O ₉ ·xHCl Mol. weight 506,41 (free base) g/mol		
PTC1270 Pomalidomide-PEG5-NH₂ hydrochloride 17-Amino-N-(2-(2,6-dioxopiperidin-3-yl)-1,3-dioxoisindolin-4-yl)-3,6,9,12,15-pentaoxaheptadecanamide hydrochloride CAS-No. 2421217-05-0 Formula C ₂₅ H ₃₄ N ₄ O ₁₀ ·xHCl Mol. weight 550,56 (free base) g/mol		
PTC1280 Pomalidomide-PEG6-NH₂ hydrochloride 20-Amino-N-(2-(2,6-dioxopiperidin-3-yl)-1,3-dioxoisindolin-4-yl)-3,6,9,12,15,18-hexaoxaicosanamide hydrochloride Formula C ₂₇ H ₃₈ N ₄ O ₁₁ ·xHCl Mol. weight 594,61 (free base) g/mol		
PTC1300 Pomalidomide-PEG6-butyl-NH₂ hydrochloride 4-Amino-N-(2-(2,6-dioxopiperidin-3-yl)-1,3-dioxoisindolin-4-yl)acetamide Formula C ₃₁ H ₄₆ N ₄ O ₁₁ ·xHCl Mol. weight 650,72 (free base) g/mol		

PTC1310 (S,R,S)-AHPC-PEG1-NH₂ hydrochloride

(2S,4R)-1-((S)-2-(2-(2-Aminoethoxy)acetamido)-3,3-dimethylbutanoyl)-4-hydroxy-N-(4-(4-methylthiazol-5-yl)benzyl)pyrrolidine-2-carboxamide hydrochloride


CAS-No. 2711076-33-2
Formula C₂₆H₃₇N₅O₅S*xHCl
Mol. weight 531,67 (free base) g/mol



PTC1320 (S,R,S)-AHPC-PEG2-NH₂ hydrochloride

(2S,4R)-1-((S)-2-(2-(2-(2-Aminoethoxy)ethoxy)acetamido)-3,3-dimethylbutanoyl)-4-hydroxy-N-(4-(4-methylthiazol-5-yl)benzyl)pyrrolidine-2-carboxamide hydrochloride


CAS-No. 2097973-72-1
Formula C₂₈H₄₁N₅O₆S*xHCl
Mol. weight 575,72 (free base) g/mol



PTC1330 (S,R,S)-AHPC-PEG3-NH₂ hydrochloride

(2S,4R)-1-((S)-14-Amino-2-(tert-butyl)-4-oxo-6,9,12-trioxa-3-azatetradecanoyl)-4-hydroxy-N-(4-(4-methylthiazol-5-yl)benzyl)pyrrolidine-2-carboxamide hydrochloride


CAS-No. 2097971-11-2
Formula C₃₀H₄₅N₅O₇S*xHCl
Mol. weight 619,77 (free base) g/mol



PTC1340 (S,R,S)-AHPC-PEG4-NH₂ hydrochloride

(2S,4R)-1-((S)-17-Amino-2-(tert-butyl)-4-oxo-6,9,12,15-tetraoxa-3-azaheptadecanoyl)-4-hydroxy-N-(4-(4-methylthiazol-5-yl)benzyl)pyrrolidine-2-carboxamide hydrochloride


CAS-No. 2010159-57-4
Formula C₃₂H₄₉N₅O₈S*xHCl
Mol. weight 663,83 (free base) g/mol







PTC1350 (S,R,S)-AHPC-PEG5-NH₂ hydrochloride


(2S,4R)-1-((S)-20-Amino-2-(tert-butyl)-4-oxo-6,9,12,15,18-pentaoxa-3-azaicosanoyl)-4-hydroxy-N-(4-(4-methylthiazol-5-yl)benzyl)pyrrolidine-2-carboxamide hydrochloride

CAS-No. 2377275-23-3
Formula C₃₄H₅₃N₅O₉S*xHCl
Mol. weight 707,88 (free base) g/mol



		Product details
PTC1360 (S,R,S)-AHPC-PEG6-NH₂ hydrochloride (2S,4R)-1-((S)-23-Amino-2-(<i>tert</i> -butyl)-4-oxo-6,9,12,15,18,21-hexaoxa-3-azatricosanoyl)-4-hydroxy-N-(4-(4-methylthiazol-5-yl)benzyl)pyrrolidine-2-carboxamide hydrochloride 		
CAS-No.	2924759-90-8	
Formula	C ₃₆ H ₅₇ N ₅ O ₁₀ S*xHCl	
Mol. weight	751,93 (free base) g/mol	
PTC1370 (S,R,S)-AHPC-PEG2-butyl-NH₂ hydrochloride (2S,4R)-1-((S)-2-(2-(2-((6-Aminoethyl)oxy)ethoxy)acetamido)-3,3-dimethylbutanoyl)-4-hydroxy-N-(4-(4-methylthiazol-5-yl)benzyl)pyrrolidine-2-carboxamide 		
CAS-No.	2421187-85-9 net	
Formula	C ₃₂ H ₄₉ N ₅ O ₆ S*xHCl	
Mol. weight	631,83 (free base) g/mol	
PTC1380 (S,R,S)-AHPC-C6-PEG3-butyl-NH₂ hydrochloride (2S,4R)-1-((S)-22-Amino-2-(<i>tert</i> -butyl)-4-oxo-10,13,16-trioxa-3-azadocosanoyl)-4-hydroxy-N-(4-(4-methylthiazol-5-yl)benzyl)pyrrolidine-2-carboxamide 		
CAS-No.	2421187-84-8 net	
Formula	C ₃₈ H ₆₁ N ₅ O ₇ S*xHCl	
Mol. weight	731,99 (free base) g/mol	
PTC1390 (S,R,S)-AHPC-PEG6-butyl-NH₂ hydrochloride (2S,4R)-1-((S)-27-Amino-2-(<i>tert</i> -butyl)-4-oxo-6,9,12,15,18,21-hexaoxa-3-azaheptacosanoyl)-4-hydroxy-N-(4-(4-methylthiazol-5-yl)benzyl)pyrrolidine-2-carboxamide 		
Formula	C ₄₀ H ₆₅ N ₅ O ₁₀ S*xHCl	
Mol. weight	808,04 (free base) g/mol	

Click-Reactive Partial PROTACs

		Product details
PTC1400 Pomalidomide-PEG1-Alkyne N-(2-(2,6-Dioxopiperidin-3-yl)-1,3-dioxoisindolin-4-yl)-3-(prop-2-yn-1-yloxy)propanamide 		
CAS-No.	2236109-19-4	
Formula	C ₁₉ H ₁₇ N ₃ O ₆	
Mol. weight	383,35 g/mol	


PTC1410 Pomalidomide-PEG2-Alkyne

N-(2-(2,6-Dioxopiperidin-3-yl)-1,3-dioxoisindolin-4-yl)-3-(2-(prop-2-yn-1-yloxy)ethoxy)propanamide

CAS-No. 2243000-25-9

Formula $C_{21}H_{21}N_3O_7$

Mol. weight 427,41 g/mol




PTC1420 Pomalidomide-PEG3-Alkyne

N-(2-(2,6-Dioxopiperidin-3-yl)-1,3-dioxoisindolin-4-yl)-3-(2-(2-(prop-2-yn-1-yloxy)ethoxy)ethoxy)propanamide

CAS-No. 2236109-20-7

Formula $C_{23}H_{25}N_3O_8$

Mol. weight 471,46 g/mol




PTC1440 Pomalidomide-PEG5-Alkyne

N-(2-(2,6-Dioxopiperidin-3-yl)-1,3-dioxoisindolin-4-yl)-4,7,10,13,16-pentaoxanonadec-18-ynamide

Formula $C_{27}H_{33}N_3O_{10}$

Mol. weight 559,57 g/mol




PTC1460 (S,R,S)-AHPC-PEG1-Alkyne

(2S,4R)-1-((S)-3,3-Dimethyl-2-(3-(prop-2-yn-1-yloxy)propanamido)butanoyl)-4-hydroxy-N-(4-(4-methylthiazol-5-yl)benzyl)pyrrolidine-2-carboxamide

CAS-No. 2242965-71-3

Formula $C_{28}H_{36}N_4O_5S$

Mol. weight 540,67 g/mol




PTC1470 (S,R,S)-AHPC-PEG2-Alkyne

(2S,4R)-1-((S)-3,3-Dimethyl-2-(3-(2-(prop-2-yn-1-yloxy)ethoxy)propanamido)butanoyl)-4-hydroxy-N-(4-(4-methylthiazol-5-yl)benzyl)pyrrolidine-2-carboxamide

CAS-No. 2242965-72-4

Formula $C_{30}H_{40}N_4O_6S$

Mol. weight 584,73 g/mol




PTC1480 (S,R,S)-AHPC-PEG3-Alkyne







(2S,4R)-1-((S)-2-(tert-Butyl)-4-oxo-7,10,13-trioxa-3-aza-hexadec-15-ynoyl)-4-hydroxy-N-(4-(4-methylthiazol-5-yl)benzyl)pyrrolidine-2-carboxamide

CAS-No. 2374122-30-0

Formula $C_{32}H_{44}N_4O_7S$

Mol. weight 628,78 g/mol



		Product details
PTC1490 (S,R,S)-AHPC-PEG4-Alkyne (2S,4R)-1-((S)-2-(<i>tert</i> -Butyl)-4-oxo-7,10,13,16-tetraoxa-3-azanonadec-18-ynoyl)-4-hydroxy-N-(4-(4-methylthiazol-5-yl)benzyl)pyrrolidine-2-carboxamide CAS-No. 2267282-21-1 Formula C ₃₄ H ₄₈ N ₄ O ₈ S Mol. weight 672,83 g/mol		
PTC1500 (S,R,S)-AHPC-PEG5-Alkyne (2S,4R)-1-((S)-2-(<i>tert</i> -Butyl)-4-oxo-7,10,13,16,19-pentaoxa-3-azadocos-21-ynoyl)-4-hydroxy-N-(4-(4-methylthiazol-5-yl)benzyl)pyrrolidine-2-carboxamide CAS-No. 2817805-63-1 Formula C ₃₆ H ₅₂ N ₄ O ₉ S Mol. weight 716,88 g/mol		
PTC1510 (S,R,S)-AHPC-PEG6-Alkyne (2S,4R)-1-((S)-2-(<i>tert</i> -Butyl)-4-oxo-7,10,13,16,19,22-hexaoxa-3-azapentacos-24-ynoyl)-4-hydroxy-N-(4-(4-methylthiazol-5-yl)benzyl)pyrrolidine-2-carboxamide Formula C ₃₈ H ₅₆ N ₄ O ₁₀ S Mol. weight 760,94 g/mol		
PTC1520 Pomalidomid- PEG1-N₃ 2-(2-Azidoethoxy)-N-(2-(2,6-dioxopiperidin-3-yl)-1,3-dioxoisindolin-4-yl)acetamide CAS-No. 2133360-04-8 Formula C ₁₇ H ₁₆ N ₆ O ₆ Mol. weight 400,35 g/mol		
PTC1530 Pomalidomid- PEG2-N₃ 2-(2-(2-Azidoethoxy)ethoxy)-N-(2-(2,6-dioxopiperidin-3-yl)-1,3-dioxoisindolin-4-yl)acetamide CAS-No. 2267306-14-7 Formula C ₁₉ H ₂₀ N ₆ O ₇ Mol. weight 444,4 g/mol		
PTC1540 Pomalidomid- PEG3-N₃ 2-(2-(2-(2-azidoethoxy)ethoxy)ethoxy)-N-(2-(2,6-dioxopiperidin-3-yl)-1,3-dioxoisindolin-4-yl)acetamide CAS-No. 2267306-15-8 Formula C ₂₁ H ₂₄ N ₆ O ₈ Mol. weight 488,45 g/mol		


PTC1560 Pomalidomid-C6-PEG3-butyl-N₃

6-(2-(2-((6-Azidohexyl)oxy)ethoxy)ethoxy)-N-(2-(2,6-dioxopiperidin-3-yl)-1,3-dioxoisindolin-4-yl)hexanamide

CAS-No. 2300178-66-7

Formula C₂₉H₄₀N₆O₈

Mol. weight 600,66 g/mol




PTC1570 Pomalidomid-C6-PEG1-C3-PEG1-butyl-N₃

6-((5-((6-Azidohexyl)oxy)pentyl)oxy)-N-(2-(2,6-dioxopiperidin-3-yl)-1,3-dioxoisindolin-4-yl)hexanamide

CAS-No. 2300178-65-6

Formula C₃₀H₄₂N₆O₇

Mol. weight 598,69 g/mol




PTC1580 Pomalidomid-PEG6-butyl-N₃

4-azido-N-(2-(2,6-dioxopiperidin-3-yl)-1,3-dioxoisindolin-4-yl)-3,6,9,12,15,18-hexaoxatetracosanamide

Formula C₃₁H₄₄N₆O₁₁

Mol. weight 676,71 g/mol




PTC1590 (S,R,S)-AHPC-PEG1-N₃

(2S,4R)-1-((S)-2-(2-(2-Azidoethoxy)acetamido)-3,3-dimethylbutanoyl)-4-hydroxy-N-(4-(4-methylthiazol-5-yl)benzyl)pyrrolidine-2-carboxamide

CAS-No. 2101200-09-1

Formula C₂₆H₃₅N₇O₅S

Mol. weight 557,67 g/mol




PTC1600 (S,R,S)-AHPC-PEG2-N₃

(2S,4R)-1-((S)-2-(2-(2-(2-Azidoethoxy)ethoxy)acetamido)-3,3-dimethylbutanoyl)-4-hydroxy-N-(4-(4-methylthiazol-5-yl)benzyl)pyrrolidine-2-carboxamide

CAS-No. 2010159-45-0

Formula C₂₈H₃₉N₇O₆S

Mol. weight 601,72 g/mol




PTC1610 (S,R,S)-AHPC-PEG3-N₃






(2S,4R)-1-((S)-14-azido-2-(tert-butyl)-4-oxo-6,9,12-trioxo-3-azatetradecanoyl)-4-hydroxy-N-(4-(4-methylthiazol-5-yl)benzyl)pyrrolidine-2-carboxamide

CAS-No. 1797406-80-4

Formula C₃₀H₄₃N₇O₇S

Mol. weight 645,77 g/mol



		Product details
PTC1640 (S,R,S)-AHPC-PEG6-N₃ (2S,4R)-1-((S)-23-Azido-2-(<i>tert</i> -butyl)-4-oxo-6,9,12,15,18,21-hexaoxa-3-azatricosanoyl)-4-hydroxy-N-(4-(4-methylthiazol-5-yl)benzyl)pyrrolidine-2-carboxamide CAS-No. 2086298-71-5 Formula C ₃₆ H ₅₅ N ₇ O ₁₀ S Mol. weight 777,93 g/mol		
PTC1680 (S,R,S)-AHPC-PEG6-butyl-N₃ (2S,4R)-1-((S)-27-Azido-2-(<i>tert</i> -butyl)-4-oxo-6,9,12,15,18,21-hexaoxa-3-azaheptacosanoyl)-4-hydroxy-N-(4-(4-methylthiazol-5-yl)benzyl)pyrrolidine-2-carboxamide Formula C ₄₀ H ₆₃ N ₇ O ₁₀ S Mol. weight 834,03 g/mol		
Thiol-Reactive Partial PROTACs		
		Product details
PTC1690 Pomalidomid-PEG2-butyl-I N-(2-(2,6-Dioxopiperidin-3-yl)-1,3-dioxoisindolin-4-yl)-2-(2-((6-iodohexyl)oxy)ethoxy)acetamide CAS-No. 1835705-72-0 Formula C ₂₃ H ₂₈ IN ₃ O ₇ Mol. weight 585,39 g/mol		
PTC1700 Pomalidomid-C6-PEG3-butyl-I N-(2-(2,6-Dioxopiperidin-3-yl)-1,3-dioxoisindolin-4-yl)-6-(2-(2-((6-iodohexyl)oxy)ethoxy)ethoxy)hexanamide CAS-No. 1835705-70-8 Formula C ₂₉ H ₄₀ IN ₃ O ₈ Mol. weight 685,55 g/mol		
PTC1710 Pomalidomid-C6-PEG1-C3-PEG1-butyl-I N-(2-(2,6-Dioxopiperidin-3-yl)-1,3-dioxoisindolin-4-yl)-6-((5-((6-iodohexyl)oxy)pentyl)oxy)hexanamide CAS-No. 1835705-76-4 Formula C ₃₀ H ₄₂ IN ₃ O ₇ Mol. weight 683,57 g/mol		

[↑ back to content](#)


PTC1720 Pomalidomid-PEG6-butyl-I

N-(2-(2,6-Dioxopiperidin-3-yl)-1,3-dioxoisindolin-4-yl)-24-iodo-3,6,9,12,15,18-hexaoxatetracosanamide

CAS-No. 1835705-74-2

Formula $C_{31}H_{44}IN_3O_{11}$

Mol. weight 761,6 g/mol




PTC1730 (S,R,S)-AHPC-PEG2-butyl-Cl

(2S,4R)-1-((S)-2-(2-((6-Chlorohexyl)oxy)ethoxy)acetamido)-3,3-dimethylbutanoyl)-4-hydroxy-N-(4-(4-methylthiazol-5-yl)benzyl)pyrrolidine-2-carboxamide

CAS-No. 1835705-57-1

Formula $C_{32}H_{47}ClN_4O_6S$

Mol. weight 651,26 g/mol




PTC1750 (S,R,S)-AHPC-C6-PEG3-butyl-Cl

(2S,4R)-1-((S)-2-(tert-Butyl)-22-chloro-4-oxo-10,13,16-trioxa-3-azadocosanoyl)-4-hydroxy-N-(4-(4-methylthiazol-5-yl)benzyl)pyrrolidine-2-carboxamide

CAS-No. 1835705-55-9

Formula $C_{38}H_{59}ClN_4O_7S$


Mol. weight 751,42 g/mol



References:

- *Bifunctional Molecules beyond PROTACs*; S. J. Conway; **J. Med. Chem.** 2020; **63**: 2802-2806.
<https://doi.org/10.1021/acs.jmedchem.0c00293>
- *Targeted protein degradation by PROTACs*; T. K. Neklesa, J. D. Winkler, C. M. Crews; **Pharmacol Ther** 2017; **174**: 138-144. <https://doi.org/10.1016/j.pharmthera.2017.02.027>
- *Targeted Protein Degradation: from Chemical Biology to Drug Discovery*; P. M. Cromm, C. M. Crews; **Cell Chem Biol** 2017; **24**: 1181-1190. <https://doi.org/10.1016/j.chembiol.2017.05.024>
- *Targeted Protein Degradation by Small Molecules*; D. P. Bondeson, C. M. Crews; **Annu Rev Pharmacol Toxicol** 2017; **57**: 107-123. <https://doi.org/10.1146/annurev-pharmtox-010715-103507>
- *Small-Molecule PROTACs: New Approaches to Protein Degradation*; M. Toure, C. M. Crews; **Angew Chem Int Ed** 2016; **55**: 1966-73. <https://doi.org/10.1002/anie.201507978>
- *Impact of linker length on the activity of PROTACs*; K. Cyrus, M. Wehenkel, E. Y. Choi, H. J. Han, H. Lee, H. Swanson, K. B. Kim; **Mol Biosyst** 2011; **7**: 359-64. <https://doi.org/10.1039/c0mb00074d>
- *Development of potent monoclonal antibody auristatin conjugates for cancer therapy*; S. O. Doronina, B. E. Toki, M. Y. Torgov, B. A. Mendelsohn, C. G. Cervený, D. F. Chace, R. L. DeBlanc, R. P. Gearing, T. D. Bovee, C. B. Siegall, J. A. Francisco, A. F. Wahl, D. L. Meyer, P. D. Senter; **Nat Biotechnol** 2003; **21**: 778-84. <https://doi.org/10.1038/nbt832>
- *Gemtuzumab ozogamicin, a potent and selective anti-CD33 antibody-calicheamicin conjugate for treatment of acute myeloid leukemia*; P. R. Hamann, L. M. Hinman, I. Hollander, C. F. Beyer, D. Lindh, R. Holcomb, W. Hallett, H. R. Tsou, J. Upeslakis, D. Shochat, A. Mountain, D. A. Flowers, I. Bernstein; **Bioconjug Chem** 2002; **13**: 47-58.
<https://doi.org/10.1021/bc010021y>

In addition to these pre-designed building blocks, we offer custom synthesis of your required ligand-linker combination (Fig. 37) or “complete PROTAC”. This allows to design a library of slightly different PROTACs in order to find the best combination for your application, as even small changes in ligands and cross-linkers might affect the efficiency of the formation of the ternary complex.




Please contact us for Custom Synthesis of the PROTAC[®] linker fragment of your choice or complete functional PROTAC[®].

5.6. Site-Selective π -Clamp Mediated Cysteine Arylation




For the modification of proteins, cysteine is a suitable choice for bioconjugation because of the high nucleophilicity of its thiol group as well as the low abundance of cysteine residues in the majority of naturally occurring proteins (ca. 1.7%). However, conventional chemical cysteine-based conjugation techniques are not site-specific resulting in a multiple labelling mixture at random positions.

In nature, selective reactions in proteins are triggered by the formation of certain microenvironments by three-dimensional secondary structures of the surrounding polypeptide architecture. Inspired by nature, Pentelute *et al.* envisioned a new strategy for site-selective chemistry, the so called “ π -clamp” arylation. Cysteine embedded in the four-amino-acid sequence (Phe-Cys-Pro-Phe) tunes up the thiol reactivity for site-selective conjugation with perfluoroaromatic reagents. Thus, the π -clamp allows the selective modification of one cysteine site in a protein containing multiple other endogenous cysteine residues.



The reported reaction is site-specific, operational under physiological conditions, enzyme-free, and as efficient as the commonly used azide-alkyne click chemistry. Furthermore, the π -clamp works as part of the N-terminus, the C-terminus, as well as in the middle of a polypeptide chain. Besides, its small size hardly perturbs the target protein's native structure.

In addition to the possibility for site-selective cysteine labeling in any linear peptide, the π -clamp approach should allow for macrocyclization between cysteine thiolates as a last, high-yielding synthetic step either via an exogenously added perfluoroaryl-based linker, or by incorporating non-crosslinked perfluoroaryl-based moieties first, followed by their macrocyclization with dithiol reagents.

		Product details
RL-4030	PFB-mercaptopropionyl-PEG3-N₃ Perfluorobiphenyl-mercaptopropionyl-PEG(3)-N ₃ Formula $C_{23}H_{21}F_9N_4O_4S$ Mol. weight 620,49 g/mol	 
RL-4040	PFB-mercaptopropionyl-AEEA Perfluorobiphenyl-mercaptopropionyl-AEEA Formula $C_{21}H_{16}F_9NO_5S$ Mol. weight 565,41 g/mol	 
RL-4050	PFB-mercaptopropionyl-TOTA-Biotin Perfluorobiphenyl-mercaptopropionyl-TOTA-Biotin Formula $C_{35}H_{41}F_9N_4O_6S_2$ Mol. weight 848,84 g/mol	 







References:

- Site-Specific Small Molecule Labeling of an Internal Loop in JC Polyomavirus Pentamers Using the π -Clamp-Mediated Cysteine Conjugation; J. A. Baccile, P. J. Voorhees, A. J. Chillo, M. Berry, R. Morgenstern, T. J. Schwertfeger, F. M. Rossi, C. D. S. Nelson; **ChemBiochem** 2021; **22(21)**: 3037-3041.
<https://doi.org/10.1002/cbic.202100188>
- Convergent diversity-oriented side-chain macrocyclization scan for unprotected polypeptides; Y. Zou, A. M. Spokoyny, C. Zhang, M. D. Simon, H. Yu, Y.-S. Lin, B. L. Pentelute; **Org. Biomol. Chem.** 2014; **12**: 566-573.
<https://doi.org/10.1039/C3OB42168F>
- Enzymatic „click“ ligation: selective cysteine modification in polypeptides enabled by promiscuous glutathione S-transferase; C. Zhang, A. M. Spokoyny, Y. Zou, M. D. Simon, B. L. Pentelute; **Angew Chem Int Ed Engl.** 2013; **52(52)**: 14001-5. <https://doi.org/10.1002/anie.201306430>
- π -Clamp-mediated cysteine conjugation; C. Zhang, M. Welborn, T. Zhu, N. J. Yang, M. S. Santos, T. Van Voorhis, B. L. Pentelute; **Nat Chem.** 2016; **8**: 120-128. <https://doi.org/10.1038/nchem.2413>

[↑ back to content](#)


6. Preparing Carriers for Conjugation

Tab. 3: Conceptual overview of conjugation technologies.

Carrier	Conjugation Chemistry	Cleavage Mechanism	Fragmentation for Traceless Release	Cargo Functionality
Biopolymers: • Peptides • Proteins • Antibodies • Single Chain • Nanobodies • Camelids • Oligonucleotides • Aptamers	Thioether formation with maleimide Disulfide bond formation Acylation of amines His tag acylation Click conjugation (CuAAC, SPAAC, IEDDA) Enzyme supported conjugation: • HaloTag® • CLIP-Tag™ • SNAP-Tag® • Sequence dependent conjugation (Sortase, Ligase)	Enzymatic hydrolysis: • Val-Ala • Val-Cit • Phe-Lys • Gly-Phe-Leu-Gly • Ala-Leu-Ala-Leu • Cyclobutyl-Ala • Cyclobutyl-Cit • Glucuronic acid 	<i>p</i>-Aminobenzyl <i>p</i>-Hydroxybenzyl <i>p</i>-Mercaptobenzyl  Oxathiolone  Dimethylimidazolidinone 	Primary & secondary amines  Tertiary amines  Alcohols  Phenols  Carboxylic acids 
Carbon: • Nanotubes • Fullerenes	Nitrene addition via photoactivation of perfluoroarylazides			
Metals: • Gold • Silver	Affinity of sulfur to gold and silver			
Metal oxide	Chelate formation	Oxidation: Tetramethyl-dioxaborolane pH: 5-(hydroxymethyl)pyrogallol orthoester (HMPO)		
Plastic polymers: • Teflon • Polyethylene • Polystyrene • Latex	Introduction of reactive functional groups by plasma treatment (e.g., ammonia or acrylic acid)			
Silicates	Affinity between silicon and oxygen			

6.1. Antibodies, Antibody Formats and Proteins by (Cell-free) Recombinant Methodologies

Cell-free protein synthesis, often referred to as *in vitro* translation, is a fast and viable technique which, in comparison to *in vivo* protein synthesis, leads to the production of a target protein in a considerably less laborious way. Cell-free systems are based on lysates derived from *E. coli* or eukaryotic sources and they allow for the synthesis of a broad spectrum of structurally diverse and modified proteins. The path from a DNA template to the protein of interest is reduced to only a few hours of time and additionally, no genetically modified organisms (GMOs) are needed. A variety of proteins, such as certain membrane proteins (e.g., GPCRs), toxins or transcriptional and translational factors, whose synthesis in conventional *in vivo* systems is often associated with difficulties, can be synthesized *in vitro*. Only the target protein is synthesized in cell-free systems, since endogenous mRNA templates are removed during the process of lysate preparation. The open and flexible character of cell-free systems allows the well-tuned adjustment of the reaction environment by supplementing the reaction mixture with co-factors, chaperones, detergents, rare tRNAs and buffers of varying ion composition, depending on the demands of the individual protein. The reaction conditions during protein synthesis significantly impact protein folding, protein activity and functionality.



Cell-free protein synthesis provides a fast access to proteins and antibody formats including introduction of mutations or non-canonical amino acids at defined positions for subsequent conjugation with permanent or self-immolative linkers and payloads. This includes otherwise difficult to synthesize toxic proteins, membrane proteins, labelled proteins, protein-conjugates, and antibody formats.



You need more details about cell-free protein synthesis?

Watch the recording of our workshop!



[↑ back to content](#)

Do you have a certain bioconjugate in mind?

- Membrane protein carrying a conjugation function on a specific position.
- Antibody, single chain or nanobody decorated with azido or alkyne function for initial linker attachment.
- Biomolecule-linker conjugate, ready to load your payload.
- Or the final biomolecule-drug-conjugate.
- Adaptation of your current protein ready for linker attachment.



Consult with us, we will work out the best solution for you!

References:

- Enriched cell-free and cell-based native membrane derived vesicles (nMV) enabling rapid in-vitro electrophysiological analysis of the voltage-gated sodium channel 1.5; Y. Pandey, S. K. Dondapati, S. Kubick; **Biochim Biophys Acta Biomembr** 2023; **1865**: 184144. <https://doi.org/10.1016/j.bbamem.2023.184144>
- Rapid One-Step Capturing of Native, Cell-Free Synthesized and Membrane-Embedded GLP-1R; L. Haueis, M. Stech, E. Schneider, T. Lanz, N. Hebel, A. Zemella, S. Kubick; **Int J Mol Sci** 2023; **24**: 2808. <https://doi.org/10.3390/ijms24032808>
- Evaluation of the Ion Channel Assembly in a Eukaryotic Cell-Free System Focusing on Two-Pore Domain Potassium Channels K2P; J. Ullrich, C. Ohlhoff, S. K. Dondapati, A. Zemella, S. Kubick; **International Journal of Molecular Sciences** 2023; **24**: 6299. <https://doi.org/10.3390/ijms24076299>
- Cell-Free Systems Enable the Production of AB(5) Toxins for Diagnostic Applications; F. Ramm, L. Jack, D. Kaser, J. L. Schlosshauer, A. Zemella, S. Kubick; **Toxins (Basel)** 2022; **14**: 233. <https://doi.org/10.3390/toxins14040233>
- The Potential of Eukaryotic Cell-Free Systems as a Rapid Response to Novel Zoonotic Pathogens: Analysis of SARS-CoV-2 Viral Proteins; F. Ramm, S. K. Dondapati, H. A. Trinh, D. Wenzel, R. M. Walter, A. Zemella, S. Kubick; **Front Bioeng Biotechnol** 2022; **10**: 896751. <https://doi.org/10.3389/fbioe.2022.896751>
- Synthesis of an Anti-CD7 Recombinant Immunotoxin Based on PE24 in CHO and E. coli Cell-Free Systems; S. K. Krebs, M. Stech, F. Jorde, N. Rakotoarinoro, F. Ramm, S. Marinoff, S. Bahrke, A. Danielczyk, D. A. Wustenhagen, S. Kubick; **Int J Mol Sci** 2022; **23**: 13697. <https://doi.org/10.3390/ijms232213697>
- A Cell-free Expression Pipeline for the Generation and Functional Characterization of Nanobodies; L. Haueis, M. Stech, S. Kubick; **Front Bioeng Biotechnol** 2022; **10**: 896763. <https://doi.org/10.3389/fbioe.2022.896763>
- Synthesis of Fluorescently Labeled Antibodies Using Non-Canonical Amino Acids in Eukaryotic Cell-Free Systems; M. Stech, N. Rakotoarinoro, T. Teichmann, A. Zemella, L. Thoring, S. Kubick; **Structural Proteomics: High-Throughput Methods** R. J. Owens 2021: 175-190. https://doi.org/10.1007/978-1-0716-1406-8_9
- Cell-Free Protein Synthesis: A Promising Option for Future Drug Development; S. K. Dondapati, M. Stech, A. Zemella, S. Kubick; **BioDrugs** 2020; **34**: 327-348. <https://doi.org/10.1007/s40259-020-00417-y>
- Accelerating the Production of Druggable Targets: Eukaryotic Cell-Free Systems Come into Focus; L. Thoring, A. Zemella, D. Wustenhagen, S. Kubick; **Methods and Protocols** 2019; **2**: 30. <https://doi.org/10.3390/mps2020030>

6.2. Aptamers and other Oligonucleotides

Terminal azide bearing nucleotides closely resemble natural nucleotides and are the basis of chemoenzymatic oligonucleotide labeling. T7 RNA polymerase, poly(A) polymerase and the terminal deoxynucleotidyl transferase (Tdt) tolerate azide and alkyne decorated nucleotides and incorporate them at a defined position. Any mRNA can be site-specifically labeled without special requirements or altered production protocols. By click labeling at the 3'-end the half-life of the mRNA could be modulated, which might be an option for increased mRNA stability.

Feeding living cells with 5-ethynyl uridine (EU) will lead to EU labelled RNA, which further can be modified *via* Click conjugation, e.g., with fluorescent dyes or linkers carrying small molecules, peptides or proteins. Incorporation into DNA can also occur after some incubation time. This is most likely *via* conversion of the EU ribonucleoside into EdU aided by intracellular ribonucleotide reductases.

With 5-ethynyl-dA-CEP and 5-ethynyl-dU-CEP alkyne moieties can be used under standard conditions for solid phase synthesis of oligonucleotides followed by subsequent Click conjugation.




Fig. 40: Alkyne and azido functions can be implemented at specific positions of oligonucleotide sequences by various building blocks.



With partners we have the possibility to provide modified oligonucleotide sequences ready for conjugation with suitable (self-immolative) linkers and payloads.

[↑ back to content](#)


References:

- Molecular Epitope Determination of Aptamer Complexes of the Multidomain Protein C-Met by Proteolytic Affinity-Mass Spectrometry; L. Lupu, P. Wiegand, N. Huttman, S. Rawer, W. Kleinekofort, I. Shugureva, A. S. Kichkailo, F. N. Tomilin, A. Lazarev, M. V. Berezovski, M. Przybylski; **ChemMedChem** 2020; **15**: 363-369.
<https://doi.org/10.1002/cmdc.201900489>
- Chemoenzymatic Preparation of Functional Click-Labeled Messenger RNA; S. Croce, S. Serdjukow, T. Carell, T. Frischmuth; **Chembiochem : a European journal of chemical biology** 2020; **21**: 1641-1646.
<https://doi.org/10.1002/cbic.201900718>
- Exploring RNA transcription and turnover in vivo by using click chemistry; C. Y. Jao, A. Salic; **Proc Natl Acad Sci U S A** 2008; **105**: 15779-84. <https://doi.org/10.1073/pnas.0808480105>
- Click Reaction on Solid Phase Enables High Fidelity Synthesis of Nucleobase-Modified DNA; F. Tolle, M. Rosenthal, F. Pfeiffer, G. Mayer; **Bioconjug Chem** 2016; **27**: 500-3. <https://doi.org/10.1021/acs.bioconjchem.5b00668>
- Ethynyl side-chain hydration during synthesis and workup of „clickable“ oligonucleotides: bypassing acetyl group formation by triisopropylsilyl protection; S. A. Ingale, H. Mei, P. Leonard, F. Seela; **J Org Chem** 2013; **78**: 11271-82. <https://doi.org/10.1021/jo401780u>
- Synthesis of highly modified DNA by a combination of PCR with alkyne-bearing triphosphates and click chemistry; J. Gierlich, K. Gutsmedl, P. M. Gramlich, A. Schmidt, G. A. Burley, T. Carell; **Chemistry** 2007; **13**: 9486-94.
<https://doi.org/10.1002/chem.200700502>
- Directed DNA metallization; G. A. Burley, J. Gierlich, M. R. Mofid, H. Nir, S. Tal, Y. Eichen, T. Carell; **J Am Chem Soc** 2006; **128**: 1398-9. <https://doi.org/10.1021/ja055517v>
- A versatile modification of on-column oligodeoxynucleotides using a copper-catalyzed oxidative acetylenic coupling reaction; N. Minakawa, Y. Ono, A. Matsuda; **J Am Chem Soc** 2003; **125**: 11545-52.
<https://doi.org/10.1021/ja036055t>
- DNA duplexes stabilized by modified monomer residues: synthesis and stability; D. Graham, J. A. Parkinson, T. Brown; *Journal of the Chemical Society, Perkin Transactions 1* 1998: 1131-1138.
<https://doi.org/10.1039/a707031d>

6.3. Carbon Compounds

Fullerenes and carbon nanotubes are subject of ongoing research as they possess unique geometrical shapes, as well as appealing photochemical, electrochemical, and physical properties. In addition, they act as efficient radical scavenger and antioxidant, as well as nano carrier for gene and drug delivery. Thus, a wide variety of operations can be considered. There are different options for conjugating linkers and payloads to such type of carriers.

- a) Fullerenes carrying phenylbutyric methyl ester (PBM) can be activated by saponification followed by ester or amide formation with appropriate substitutions.
- b) Fullerene derivatives with malonic acid moieties can react readily with nucleophiles, e.g., the amino functions of amino acids, amino-PEGs or other linkers and payloads enabling multiple payload decoration on one single fullerene.
- c) Perhydroxylated fullerenols are a unique class of water-soluble fullerenes. Their alcohol functions can further be derivatized by ester or ether formation, e.g., via Mitsunobo reaction.
- d) Fullerenes and carbon nanotubes carrying no functional group can be prepared for conjugation by photoaffinity labeling reagents such as aryl azides. Upon photolysis (258 nm), N₂ is liberated and a stabilized singlet perfluorophenylnitrene is being formed *in situ*, which reacts with neighboring molecules by insertion and addition reactions in moderate to good yields.




FLL1070 Fullerene C₆₀ (malonic acid)_n

Buckminsterfullerene-n(malonic acid)

Formula C₆₀(C₃H₂O₄)_n

Mol. weight 720,66+(102,05)n g/mol



Product details



The Concept of Antibody-Drug Conjugation (ADC)

Permanent Linkers

Cleavable Linkers

Trifunctional Linkers

Cross-Linkers for other Bio Applications

Preparing Carriers for Conjugation

Index

[↑ back to content](#)


FLL1020 Fullerene C₆₀ (PBM)

Fulleren-phenyl-(4-phenylbutyric acid methyl ester)

CAS-No. 160848-22-6

Formula C₇₂H₁₄O₂

Mol. weight 910,9 g/mol




FLL1010 Fullerene C₆₀ (PBM)2

Fulleren-diphenyl-bis(4-phenylbutyric acid methyl ester)

CAS-No. 1048679-01-1

Formula C₈₄H₂₈O₄

Mol. weight 1104,14 g/mol




FLL1080 Fullerene C₇₀ (malonic acid)n

Buckminsterfullerene-n-(malonic acid)

Formula C₇₀(C₃H₂O₄)_n

Mol. weight 840,77+(102,05)n g/mol




FLL1060 Fullerene C₇₀ (PBM)

Fulleren-phenyl-(4-phenylbutyric acid methyl ester)

CAS-No. 609771-63-3

Formula C₈₂H₁₄O₂

Mol. weight 1031,01 g/mol




FLL1050 Fullerene C₇₀ (PBM)2

Fulleren-diphenyl-bis(4-phenylbutyric acid methyl ester)

Formula C₉₄H₂₈O₄

Mol. weight 1221,25 g/mol




FLL1030 Fullerenol C60

Polyhydroxylated Fullerene

Formula C₆₀(OH)_n

Mol. weight 720,66+(17,01)n g/mol



FLL1090 Fullerenol C70

Polyhydroxylated Fullerene

 Formula $C_{70}(OH)_n$



 Mol. weight $840,77+(17,01)n$ g/mol


Fig. 42: Carbon nanotubes and fullerenes can be surface decorated with carboxylic acid moieties for subsequent derivatization and linker attachment via photolysis of perfluoroarylazides which form in situ stable and reactive nitrenes.

References:

- Fullerene derivatives with amino acids, peptides and proteins: From synthesis to biomedical application; E. I. Pochkaeva, N. E. Podolsky, D. N. Zakusilo, A. V. Petrov, N. A. Charykov, T. D. Vlasov, A. V. Penkova, L. V. Vasina, I. V. Murin, V. V. Sharoyko, K. N. Semenov; **Prog. Solid. State Ch.** 2020; **57**: 100255.
<https://doi.org/10.1016/j.progsolidstchem.2019.100255>
- Fullerene-based delivery systems; H. Kazemzadeh, M. Mozafari; **Drug Discov Today** 2019; **24**: 898-905.
<https://doi.org/10.1016/j.drudis.2019.01.013>
- Fullerenes in biology and medicine; E. Castro, A. H. Garcia, G. Zavala,, L. Echegoyen; **J. Mater. Chem. B.** 2017;
<https://doi.org/10.1039/c7tb00855d>
- Fullerene C60 with cytoprotective and cytotoxic potential: prospects as a novel treatment agent in Dermatology? A. Rondags, W. Yan Yuen, M. F. Jonkman, B. Horváth; **Exp. Dermatol.** 2017; **26** (3): 220-224.
<https://doi.org/10.1111/exd.13172>
- Water-soluble fullerenes for medical applications; I. Raovi; **Mater. Sci. Technol.** 2016; **33**: 777-794.
<https://doi.org/10.1080/02670836.2016.1198114>

- Functionalized Fullerenes in Photodynamic Therapy; Y.-Y. Huang, S. K. Sharma, R. Yin, T. Agrawal, L. Y. Chiang, M. R. Hamblin; **J. Biomed. Nanotechnol.** 2014; **10**: 1918-1936. <https://doi.org/10.1166/jbn.2014.1963>
- Fullerene-biomolecule conjugates and their biomedical applications; X. Yang, A. Ebrahimi, J. Li, Q. Cui; **Int. J. Nanomed.** 2014; **19**: 77-92. <http://dx.doi.org/10.2147/IJN.S52829>
- Fullerenol Nanoparticles: Toxicity and Antioxidant Activity; R. Injac, M. Prijatelj, B. Strukelj; *Oxidative Stress and Nanotechnology: Methods and Protocols*. D. Armstrong, D. J. Bharali 2013: 75-100. https://doi.org/10.1007/978-1-62703-475-3_5
- Medicinal chemistry and pharmacological potential of fullerenes and carbon nanotubes; F. Cataldo, T. Da Ros; Springer Science & Business Media; 2008; 1.
- Medicinal applications of fullerenes; R. Bakry, R. M. Vallant, M. Najam-ul-Haq, M. Rainer, Z. Szabo, C. W. Huck, G. K. Bonn; **Int. J. Nanomedicine** 2007; **2**: 639-649.
- Chapter 7 – Functionalization and application of [60]fullerene; A. Mateo-Alonso, D. Bonifazi, M. Prato; **Carbon Nanotechnology** 2006; 155-189. <https://doi.org/10.1016/B978-044451855-2/50010-3>
- Fullerene derivatives: an attractive tool for biological applications; S. Bosi, T. Da Ros, G. Spalluto, M. Prato; **Eur. J. Med. Chem.** 2003; **38**: 913-923. <https://doi.org/10.1016/j.ejmech.2003.09.005>
- Biological Applications of Fullerenes; A. W. Jensen, S. R. Wilson, D. I. Schuster; **Bioorg. Med. Chem.** 1996; **4**: 767-779. [https://doi.org/10.1016/0968.0896\(96\)00081-8](https://doi.org/10.1016/0968.0896(96)00081-8)



Interested in products and technologies for drug delivery?

Download our brochure
Polymer Therapeutics!




6.4. Metals

Nanotechnology and nanobiotechnology using gold or silver particles are broadly diverse, rapidly expanding areas of study in medical diagnostics and therapeutics, sensorics and chemistry. Gold nanoparticles (AuNPs), particularly, have found a wide range of biomedical and environmental monitoring applications (drug delivery, diagnostics, biosensing, bio-imaging, theranostics, and hazardous chemical sensing) due to their excellent optoelectronic and enhanced physico-chemical properties.


Metal particles, however, are not water soluble without further modification. Due to the soft character of gold and sulfur, thiols readily form strong dative bonds to gold and silver surfaces creating a self-assembled monolayer (SAMs), which modifies surfaces for subsequent coupling of proteins, PEGs and other molecules. The formed nanoparticles show excellent stability and can be stored for years.

The bond between gold or silver surface atoms and monothiols is sensitive to reducing agents such as DTT (Cleland's Reagent), while the disulfide lipoic acid moiety (also known as thiocetic acid) binds far stronger to metal surfaces and is much more resistant towards removal from the metal surface by DTT, TCEP and similar reagents than monothiols.



A variety of protocols exist in the literature for reducing lipoic acid to dihydrolipoic acid (DHLA), which binds instantly to the surface. Typically, tris(2-carboxyethyl)phosphine (TCEP) or sodium borohydride (NaBH_4) are being used as reducing agents. In general, TCEP reduction is carried out in water or aqueous buffer (excluding phosphate buffer, in which TCEP is unstable), in three times or greater molar excess to the lipoic acid derivative, using an incubation temperature of 25 °C to 50 °C, for about 1-2 hours. Each reduction procedure must be optimized for the particular lipoic acid derivative being reduced to the corresponding DHLA derivative.

Lipoamido-PEG-acids and lipoamido-PEG-alcohols can be used as intermediates for further derivatization after attachment to the surface. The density of functional groups on the surface can be tuned by co-coating the bifunctional mercapto- or lipoic-PEG with methoxy-PEG-lipoamides.



References:


- Gold nanoparticle surface engineering strategies and their applications in biomedicine and diagnostics; K. Mahato, S. Nagpal, M. A. Shah, A. Srivastava, P. K. Maurya, S. Roy, A. Jaiswal, R. Singh, P. Chandra; **3 Biotech** 2019; **9**: 57. <https://doi.org/10.1007/s13205-019-1577-z>
- Recent advances in separation and detection methods for thiol compounds in biological samples; T. Toyo'oka; **J Chromatogr B Analyt Technol Biomed Life Sci** 2009; **877**: 3318-30. <https://doi.org/10.1016/j.jchromb.2009.03.034>
- Polyethylene glycol-based bidentate ligands to enhance quantum dot and gold nanoparticle stability in biological media; B. C. Mei, K. Susumu, I. L. Medintz, H. Mattoussi; **Nat Protoc** 2009; **4**: 412-23. <https://doi.org/10.1038/nprot.2008.243>
- Influence of anchoring ligands and particle size on the colloidal stability and in vivo biodistribution of polyethylene glycol-coated gold nanoparticles in tumor-xenografted mice; G. Zhang, Z. Yang, W. Lu, R. Zhang, Q. Huang, M. Tian, L. Li, D. Liang, C. Li; **Biomaterials** 2009; **30**: 1928-36. <https://doi.org/10.1016/j.biomaterials.2008.12.038>
- Bioconjugate Techniques (Third Edition); G. T. Hermanson; 2013: 1146. <https://doi.org/10.1016/C2009-0-64240-9>
- Toward reliable gold nanoparticle patterning on self-assembled DNA nanoscaffold; J. Sharma, R. Chhabra, C. S. Andersen, K. V. Gothelf, H. Yan, Y. Liu; **J Am Chem Soc** 2008; **130**: 7820-1. <https://doi.org/10.1021/ja802853r>
- Modular poly(ethylene glycol) ligands for biocompatible semiconductor and gold nanocrystals with extended pH and ionic stability; B. C. Mei, K. Susumu, I. L. Medintz, J. B. Delehanty, T. J. Mountziaris, H. Mattoussi; **Journal of Materials Chemistry** 2008; **18**: 4949-4958. <https://doi.org/10.1039/b810488c>
- Oriented immobilization of antibodies with GST-fused multiple Fc-specific B-domains on a gold surface; T. H. Ha, S. O. Jung, J. M. Lee, K. Y. Lee, Y. Lee, J. S. Park, B. H. Chung; **Anal Chem** 2007; **79**: 546-56. <https://doi.org/10.1021/ac061639+>
- Design of biotin-functionalized luminescent quantum dots; K. Susumu, H. T. Uyeda, I. L. Medintz, H. Mattoussi; **J Biomed Biotechnol** 2007; **2007**: 90651. <https://doi.org/10.1155/2007/90651>
- Simultaneous determination of alpha-lipoic acid and its reduced form by high-performance liquid chromatography with fluorescence detection; S. Satoh, T. Toyo'oka, T. Fukushima, S. Inagaki; **J Chromatogr B Analyt Technol Biomed Life Sci** 2007; **854**: 109-15. <https://doi.org/10.1016/j.jchromb.2007.04.003>
- Enhanced oligonucleotide-nanoparticle conjugate stability using thioctic acid modified oligonucleotides; J. A. Dougan, C. Karlsson, W. E. Smith, D. Graham; **Nucleic Acids Res** 2007; **35**: 3668-75. <https://doi.org/10.1093/nar/gkm237>
- Biosensing with Luminescent Semiconductor Quantum Dots; K. Sapsford, T. Pons, I. Medintz, H. Mattoussi; **Sensors** 2006; **6**: 925-953.
- Synthesis and reactions of functionalised gold nanoparticles; M. Brust, J. Fink, D. Bethell, D. J. Schiffrin, C. Kiely; **Journal of the Chemical Society, Chemical Communications** 1995: 1655-1656. <https://doi.org/10.1039/c39950001655>
- Self-assembled organic monolayers: model systems for studying adsorption of proteins at surfaces; K. L. Prime, G. M. Whitesides; **Science** 1991; **252**: 1164-7. <https://doi.org/10.1126/science.252.5009.1164>

6.5. Metal Oxides

Nanotechnology and nanobiotechnology using quantum dots, magnetic particles, or metal oxides are broadly diverse, rapidly expanding areas of study in medical diagnostics and therapeutics, sensoric and chemistry. Many metal oxides are not water soluble without further modification.

Compounds containing functional groups which enable noncovalent binding or chelates to oxide surfaces and which carry residues like poly(ethylene glycol) (PEG) or other polymers equip metal oxide nanoparticles with colloidal stability and stealthiness.

Such functional groups are phosphonic acid compounds $R-H_2PO_3$ which display a strong affinity to metal ionic centers and are characterized by their multidentate binding ability. With appropriate residues the self-assembled monolayer (SAM) and multilayers can be designed. Examples have been published of cerium (CeO_2), iron ($\gamma-Fe_2O_3$), aluminum (Al_2O_3), and titanium (TiO_2) oxides of different sizes and morphologies.






References:

- Versatile Coating Platform for Metal Oxide Nanoparticles: Applications to Materials and Biological Science; J. F. Berret, A. Gailliot; **Langmuir** 2022; **38**: 5323-5338. <https://doi.org/10.1021/acs.langmuir.2c00338>
- Efficient modification of metal oxide surfaces with phosphonic acids by spray coating; A. Bulusu, S. A. Paniagua, B. A. MacLeod, A. K. Sigdel, J. J. Berry, D. C. Olson, S. R. Marder, S. Graham; **Langmuir** 2013; **29**: 3935-42. <https://doi.org/10.1021/la303354t>

6.6. Polymeric Surfaces by Plasma Treatment


Plastic polymers, like polyethylene, polystyrene, PTFE, or co-polymers thereof are materials used in many facets of daily life, including tools and devices with biological, human, and pharmacologic applications. A major property of such materials is that they are rather chemically inert to most environmental conditions.

Low pressure plasma treatment offers a new method to decorate such polymers with specific functional groups, such as amines or carboxylic acids, which offer the opportunity to further modify the surface with specific molecules. PEGs, for example, can turn such usually quite hydrophobic surfaces very hydrophilic. Attachments of fluorescent dyes will stain the particles accordingly. Attachment of peptides, proteins, like streptavidin, or antibodies opens the door to numerous biological and pharmacological applications.



One major drawback of working on surfaces is the limitation of analytical methods, as conventional technologies, like mass spectroscopy, chromatographic purification technologies, and other methodologies widely established for small molecules and biologics are not applicable in this case. For surface analytics X-ray photoelectron spectroscopy (XPS) is the method of choice. It is also known as ESCA (Electron Spectroscopy for Chemical Analysis) and is an established method for the analysis of chemical compositions of surfaces. It enables a highly sensitive, quantitative detection of all elements except hydrogen and helium, as well as for the identification of binding and oxidation states on solid surfaces. The method is very sensitive to surfaces, so that even very thin layers can be studied (2 nm to 10 nm depth of information).

The basic principle of XPS is based on the irradiation of a sample surface in vacuum with soft X-rays and analyzing the energy of the emitted photoelectrons. This energy is different from element to element and also depends on the oxidative level of the element. Therefore, analyzing a carbon signal a quantitative determination e.g., between carbon-hydrogen, carbon-nitrogen and carbonyl carbon can be measured. The XPS spectrum results from plotting the number of detected electrons per energy interval against their kinetic energy.



Information gained with X-ray photoelectron spectroscopy:

- Detection limit: 0.01 – 1 at%, sub-monolayer
- Detectable elements: Li – U
- Chemical bonding information
- Quantitative information
- Information depth: 1 nm to 10 nm
- Lateral resolution: ca. 30 μm
- Depth profiling
- Imaging/mapping




Fig. 50: XPS representing signals of the elements sodium, fluorine, oxygen, nitrogen, calcium, carbon, chlorine, and silicon. The peak area allows quantitative determination of each element.




Fig. 51: Expanded XPS spectrum of a carbon signal showing the ratio between C1 (aliphatic carbon), C2 (carbon with single bond to a hetero atom N or O), C3 (carbonyl, amide carbon), C4 (carboxyl, carbonate carbon), and C5 (fluorinated carbon).

Services available from us:

- Beads of polyethylene, polystyrene, teflon or co-polymers thereof should be surface decorated with permanent or cleavable linkers and loaded with small or biomolecules.
- Turning polymer-based surface of devices hydrophilic, biocompatible or non-immunogenic.
- Equipping the surface of membranes, devices or other parts made of a certain material with an alternative or orthogonal property.



Consult with us, we will be happy to carry your application through to market.

References:

- Characterisation of PEGylated PLGA nanoparticles comparing the nanoparticle bulk to the particle surface using UV/vis spectroscopy, SEC, ¹H NMR spectroscopy, and X-ray photoelectron spectroscopy; S. Spek, M. Haeuser, M. M. Schaefer, K. Langer; **Applied Surface Science** 2015; **347**: 378-385.
<https://doi.org/10.1016/j.apsusc.2015.04.071>
- Reaction of human macrophages on protein corona covered TiO₂ nanoparticles; C. F. Borgognoni, M. Mormann, Y. Qu, M. Schaefer, K. Langer, C. Ozturk, S. Wagner, C. Chen, Y. Zhao, H. Fuchs, K. Riehemann; **Nano-medicine** 2015; **11**: 275-82. <https://doi.org/10.1016/j.nano.2014.10.001>
- Engineering Biomaterials Surfaces Using Micropatterning; L. Gagne, G. Laroche; **Advanced Materials Research** 2006; **15-17**: 77-82. <https://doi.org/10.4028/www.scientific.net/AMR.15-17.77>
- Quantification of cation-exchanged zeolites by XPS and EDS: A comparative study; S. Fibikar, M. T. Rinke, A. Schäfer, L. D. Cola; **Microporous and Mesoporous Materials** 2010; **132**: 296-299.
<https://doi.org/10.1016/j.micromeso.2010.02.016>
- Comparison of Atmospheric-Pressure Plasma versus Low-Pressure RF Plasma for Surface Functionalization of PTFE for Biomedical Applications; C. Sarra-Bournet, S. Turgeon, D. Mantovani, G. Laroche; **Plasma Processes and Polymers** 2006; **3**: 506-515. <https://doi.org/10.1002/ppap.200600012>
- Fast element mapping of titanium wear around implants of different surface structures; U. Meyer, M. Buhner, A. Buchter, B. Kruse-Losler, T. Stamm, H. P. Wiesmann; **Clin Oral Implants Res** 2006; **17**: 206-11.
<https://doi.org/10.1111/j.1600-0501.2005.01184.x>



Please find a variety of trimethylsilyl-PEG derivatives in our PEGylation Brochure.



[↑ back to content](#)

The Concept of Antibody-Drug Conjugation (ADC)

Permanent Linkers

Cleavable Linkers

Trifunctional Linkers


Cross-Linkers for other Bio Applications

Preparing Carriers for Conjugation

Index

6.7. Silicates

Glass, quartz and silicates presenting silanol groups at their surface, which can further react with different silane reagents, such as chlorosilanes or alkoxy silanes. It can help to etch the surface prior to coating with acids to remove an outer non-reactive layer and expose silanol groups making them available to form silicon-oxygen-silicon structures and attaching a coating molecule on the surface.



References:

- Low cost and scalable method for modifying surfaces of hollow particles from hydrophilic to hydrophobic; J. Sharma, G. Polizos, D. Hun, K. Nawaz, R. Sahore; *RSC Adv* 2020; **10**: 31065-31069.
<https://doi.org/10.1039/d0ra06114j>
- Silica and Silane based polymer composite coating on glass slide by dip-Coating Method; S. Sriram, R. K. Singh, A. Kumar; *Surfaces and Interfaces* 2020; **19**: 100472.
<https://doi.org/10.1016/j.surfin.2020.100472>

Index

Product code	Product name	Page	Product code	Product name	Page
PTC1040	(S,R,S)-AHPC hydrochloride	118	RL-3240	16-Azido-palmitic acid	39
PTC1750	(S,R,S)-AHPC-C6-PEG3-butyl-Cl	129	BAA3910	18-(Boc-amino)-stearic acid	40
PTC1380	(S,R,S)-AHPC-C6-PEG3-butyl-NH ₂ hydrochloride	124	FAA7450	18-(Fmoc-amino)-stearic acid	40
PTC1460	(S,R,S)-AHPC-PEG1-Alkyne	125	RL-3250	18-Azido-stearic acid	39
PTC1590	(S,R,S)-AHPC-PEG1-N ₃	127	RL-4170	2-OPSS-Bzl-OpNC	83
PTC1310	(S,R,S)-AHPC-PEG1-NH ₂ hydrochloride	123	RL-2600	3-Mal-Bz-NHS	47
PTC1470	(S,R,S)-AHPC-PEG2-Alkyne	125	RL-2610	3-Mal-MBz-NHS	47
PTC1220	(S,R,S)-AHPC-PEG2-butyl COOH	121	RL-8655	4-(Diethoxymethyl)-1-NP-triazole	116
PTC1730	(S,R,S)-AHPC-PEG2-butyl-Cl	129	LS-3350	4-(N-Maleimido)benzophenone	47
PTC1370	(S,R,S)-AHPC-PEG2-butyl-NH ₂ hydrochloride	124	RL-2620	4-Mal-Bz-NHS	47
PTC1600	(S,R,S)-AHPC-PEG2-N ₃	127	RL-2630	4-Mal-MBz-NHS	47
PTC1320	(S,R,S)-AHPC-PEG2-NH ₂ hydrochloride	123	ADC1310	4-Pentynoyl-Val-Ala-PAB	58
PTC1480	(S,R,S)-AHPC-PEG3-Alkyne	125	ADC1320	4-Pentynoyl-Val-Ala-PAB-PNP	58
PTC1610	(S,R,S)-AHPC-PEG3-N ₃	127	ADC1140	4-Pentynoyl-Val-Cit-PAB	66
PTC1330	(S,R,S)-AHPC-PEG3-NH ₂ hydrochloride	123	ADC1150	4-Pentynoyl-Val-Cit-PAB-PNP	66
PTC1490	(S,R,S)-AHPC-PEG4-Alkyne	126	RL-8695	5HP2O((PEG)2-OH)-(CH ₂) ₄ -NHS	51
PTC1340	(S,R,S)-AHPC-PEG4-NH ₂ hydrochloride	123	LS-4670	5HP2O((PEG)2-OH)-(CH ₂) ₅ -Dansyl	51
PTC1500	(S,R,S)-AHPC-PEG5-Alkyne	126	RL-8675	5HP2O-(CH ₂) ₄ -COOH	50
PTC1350	(S,R,S)-AHPC-PEG5-NH ₂ hydrochloride	123	RL-8680	5HP2O-(CH ₂) ₄ -NHS	50
PTC1510	(S,R,S)-AHPC-PEG6-Alkyne	126	RL-8670	5HP2O-alkyne	50
PTC1680	(S,R,S)-AHPC-PEG6-butyl-N ₃	128	RL-8685	5HP2O-PEG(2)-COOH	50
PTC1390	(S,R,S)-AHPC-PEG6-butyl-NH ₂ hydrochloride	124	RL-8690	5HP2O-PEG(2)-NHS	51
PTC1640	(S,R,S)-AHPC-PEG6-N ₃	128	ADC1290	6-Azidoheptanoyl-Val-Ala-PAB	57
PTC1360	(S,R,S)-AHPC-PEG6-NH ₂ hydrochloride	124	ADC1300	6-Azidoheptanoyl-Val-Ala-PAB-PNP	57
PTC1050	(S,S,S)-AHPC hydrochloride	118	ADC1120	6-Azidoheptanoyl-Val-Cit-PAB	65
PTC1030	(±)-Thalidomide	118	ADC1130	6-Azidoheptanoyl-Val-Cit-PAB-PNP	65
RL-8650	1-NP-triazole-4-CHO	116	RL-3480	8-Azido-octanoyl-OSu	38
RL-3460	10-Undecynoyl-OSu	35	RL-2960	Acetyl-Trimethyl-Lock	9
RL-3170	11-Azido-undecanoyl-OSu	38	ADC1790	Alkyne-HMPO-OH	74
RL-3200	11-Azidoundecanoic acid	38	ADC1800	Alkyne-HMPO-PNP	74
RL-3220	12-Azido-dodecanoyl-OSu	39	RL-2055	Alkyne-myristic acid	35
RL-3210	12-Azidododecanoic acid	38	RL-2060	Alkyne-palmitic acid	35
BAA4240	14-(Boc-amino)-myristic acid	39	PEG5440	Alkyne-PEG(4)-mal	31
FAA8160	14-(Fmoc-amino)-myristic acid	39	PEG5430	Alkyne-PEG(4)-NH ₂	31
RL-3230	14-Azido-myristic acid	39	PEG5410	Alkyne-PEG(4)-NHS	31
BAA3900	16-(Boc-amino)-palmitic acid	40	ADC1350	Alkyne-PEG(4)-Val-Ala-PAB	58
FAA7460	16-(Fmoc-amino)-palmitic acid	40	ADC1360	Alkyne-PEG(4)-Val-Ala-PAB-PNP	59

The Concept of Antibody-Drug Conjugation (ADC)

Permanent Linkers

Cleavable Linkers

Trifunctional Linkers

Cross-Linkers for other Bio Applications

Preparing Carriers for Conjugation

Index

[↑ back to content](#)

Product code	Product name	Page	Product code	Product name	Page
RL-3940	Alkyne-PEG(5)-SNAP	106	RL-3300	Biotin-SS-COOH	80
ADC1180	Alkyne-PEG(5)-Val-Cit-PAB	67	RL-4120	Biotin-SS-N ₃	81
ADC1190	Alkyne-PEG(5)-Val-Cit-PAB-PNP	67	LS-3570	Biotin-SS-Tyramide	82
RL-3930	Alkyne-SNAP	106	PEG5385	Boc,Z-AEEEE	25
RL-3330	Alkyne-SS-COOH	78	BAA4870	Boc-Aca-Aca-OH	36
RL-2065	Alkyne-stearic acid	35	RL-2810	Boc-AEDl-OH	78
AAA1905	Aloc-O ₂ Oc-OH*DCHA	20	BNN1170	Boc-Cystamine	76
RL-2035	ATFB	54	BNN1063	Boc-Cystamine*HCl	76
RL-2045	ATFB-NHS	54	BAA2180	Boc-Cystamine-Suc-OH	77
RL-3960	Azide-PEG(4)-SNAP	107	BAA1171	Boc-D-Dab(Dde)-OH	92
RL-3950	Azide-SNAP	107	BAA1176	Boc-D-Dap(Dde)-OH	92
ADC1580	Azido-cyclobutane-1,1-dicarboxamide-Ala-PAB	63	BNN1016	Boc-DOOA	18
ADC1590	Azido-cyclobutane-1,1-dicarboxamide-Ala-PAB-PNP	63	BNN1380	Boc-EDA-Suc-OH	36
ADC1480	Azido-cyclobutane-1,1-dicarboxamide-Cit-PAB	69	BAA1191	Boc-L-Dab(Dde)-OH	92
ADC1490	Azido-cyclobutane-1,1-dicarboxamide-Cit-PAB-PNP	70	BAA6480	Boc-L-Dab(Mal)-OH	46
ADC1330	Azido-PEG(4)-Val-Ala-PAB	57	BAA1193	Boc-L-Dap(Dde)-OH	92
ADC1340	Azido-PEG(4)-Val-Ala-PAB-PNP	58	BAA6475	Boc-L-Dap(Mal)-OH*DCHA	45
ADC1160	Azido-PEG(4)-Val-Cit-PAB	66	BAA1286	Boc-L-Lys(Dde)-OH*DCHA	94
ADC1170	Azido-PEG(4)-Val-Cit-PAB-PNP	66	BAA1287	Boc-L-Lys(ivDde)-OH	94
RL-3320	Azido-Pen-SS-COOH	77	BAA1197	Boc-L-Orn(Dde)-OH	93
RL-4100	Azido-SS-COOH	77	PEG4960	Boc-NH-PEG(2)-N ₃	18
RL-4150	Azido-SS-OpNC	79	PEG7870	Boc-NH-PEG(3)-NH ₂	26
RL-8425	Biocytin-Mal	48	PEG6835	Boc-NH-PEG(3)-NH ₂ *HCl	17
RL-3100	Biotin-AEEA-OPhOMe	111	PEG1920	Boc-NH-PEG(4)-COOH	29
RL-3870	Biotin-Clip	106	PEG7880	Boc-NH-PEG(4)-NH ₂	28
LS-4020	Biotin-Dde	89	PEG1915	Boc-NH-PEG(4)-OH	29
LS-4000	Biotin-Dde-Tyramide	89	RL-3560	Boc-NH-SS-Bzl-OH	80
RL-4060	Biotin-DOOA	20, 35	RL-3570	Boc-NH-SS-Bzl-OpNC	80
RL-8415	Biotin-Hx-SS-Py	80	RL-3510	Boc-NH-SS-OH	79
RL-8420	Biotin-NH-NH-Mal	48	RL-3520	Boc-NH-SS-OpNC	79
PEG7980	Biotin-PEG(4)-Dde-Alkyne	90	BAA1485	Boc-O ₂ Oc-O ₂ Oc-OH	27
PEG8140	Biotin-PEG(4)-Dde-DBCO	90	PEG8080	Boc-O ₂ Oc-OH	21
PEG7960	Biotin-PEG(4)-Dde-N ₃	89	BAA1466	Boc-O ₂ Oc-OH*DCHA	21
PEG7970	Biotin-PEG(4)-Dde-Picolyl-N ₃	89	RL-2190	Boc-SS-COOH	78
PEG8130	Biotin-PEG(4)-Dde-Tyramide	89	BNN1028	Boc-TOTA	24
PEG8110	Biotin-PEG(4)-SS-Alkyne	81	ADC1040	Boc-Val-Ala-PAB	59
PEG8100	Biotin-PEG(4)-SS-Azide	81	ADC1660	Boc-Val-Ala-PAB-Cl	59
PEG8090	Biotin-PEG(4)-SS-COOH	81	ADC1050	Boc-Val-Ala-PAB-PNP	59
PEG8120	Biotin-PEG(4)-SS-DBCO	81	ADC1020	Boc-Val-Cit-PAB	67
LS-3930	Biotin-PEG(4)-SS-Tyramide	82	ADC1010	Boc-Val-Cit-PAB-PNP	67
RL-3860	Biotin-SNAP	108	PEG7860	Boc2-AEEEE	25

Product code	Product name	Page	Product code	Product name	Page
RL-1008	Br-PAM-Linker	41	RL-3260	Fmoc-Aca-DIM	54, 91
PEG7190	Bromoacetamido-PEG(3)-N ₃	26	RL-2800	Fmoc-AEDI-OH	79
RL-2770	BSSS	112	PEG5370	Fmoc-AEEE	24
RL-3600	DACN(Ms)*HCl	102	PEG5380	Fmoc-AEEEE	25
RL-3610	DACN(Ms,Ns)	102	PEG1810	Fmoc-AEEP	23
RL-2735	DACN(Tos)*HCl	102	RL-3270	Fmoc-AEEP-DIM	54, 91
RL-2710	DACN(Tos,Ns)	101	RL-3470	Fmoc-AEEP-DMB	91
RL-2720	DACN(Tos,Suc-OH)	101	FAA8690	Fmoc-Aeg(Dde)-OH	95
RL-2730	DACN(Tos2)	102	RL-3370	Fmoc-Cystamine*HCl	77
AAA2190	DAPOA*DCHA	101	RL-3310	Fmoc-Cystamine-Suc	77
RL-4020	DBCO-C6-Alkyne	40	FAA1318	Fmoc-D-Dab(Dde)-OH	97
ADC1620	DBCO-cyclobutane-1,1-dicarboxamide-Ala-PAB	64	FAA1473	Fmoc-D-Dab(ivDde)-OH	97
ADC1630	DBCO-cyclobutane-1,1-dicarboxamide-Ala-PAB-PNP	64	FAA1476	Fmoc-D-Dap(Dde)-OH	96
ADC1520	DBCO-cyclobutane-1,1-dicarboxamide-Cit-PAB	70	FAA1478	Fmoc-D-Dap(ivDde)-OH	96
ADC1530	DBCO-cyclobutane-1,1-dicarboxamide-Cit-PAB-PNP	70	FAA1486	Fmoc-D-Lys(Dde)-OH	98
RL-2490	DBCO-mal	43	FAA1488	Fmoc-D-Lys(ivDde)-OH	98
RL-2500	DBCO-PEG(4)-mal	43	FAA8845	Fmoc-D-Lys(MeDmb)-OH	99
RL-2420	DBCO-PEG(4)-NH ₂ *TFA	32	FAA2090	Fmoc-D-Orn(Dde)-OH	97
RL-2510	DBCO-PEG(4)-OH	32	FAA1493	Fmoc-D-Orn(ivDde)-OH.solv.	98
RL-4010	DBCO-SNAP	108	FNN1007	Fmoc-DOOA*HCl	19
RL-4110	DBCO-Suc-SS-COOH	78	PEG5180	Fmoc-DOOA-DIG-OH	19
RL-2421	DBCO-Sulfo-PEG(4)-NH ₂	32	PEG4970	Fmoc-Ebes	23
DAA1004	Dde-D-Dab(Fmoc)-OH	93	FAA8815	Fmoc-L-Abu(3-Dde-amino)-OH (2S,3S)	96
DAA1006	Dde-D-Dap(Fmoc)-OH	93	FAA1365	Fmoc-L-Dab(Dde)-OH	96
DAA1017	Dde-D-Lys(Fmoc)-OH	95	FAA1458	Fmoc-L-Dab(ivDde)-OH	97
DAA1010	Dde-L-Dab(Fmoc)-OH	93	FAA1462	Fmoc-L-Dap(Dde)-OH	96
DAA1012	Dde-L-Dap(Fmoc)-OH	93	FAA1464	Fmoc-L-Dap(ivDde)-OH	96
DAA1013	Dde-L-Lys(Aloc)-OH*DCHA	94	FAA1390	Fmoc-L-Lys(Dde)-OH	98
DAA1014	Dde-L-Lys(Boc)-OH	94	FAA1500	Fmoc-L-Lys(ivDde)-OH	98
DAA1015	Dde-L-Lys(Fmoc)-OH	94	FAA7975	Fmoc-L-Lys(ivDmb)-OH	99
DAA1020	Dde-L-Met-OH	95	FAA8840	Fmoc-L-Lys(MeDmb)-OH	98
DAA1001	Dde-L-Orn(Aloc)-OH	93	FAA8145	Fmoc-L-Lys(N ₃ -Aca-DIM)-OH	91
DAA1002	Dde-L-Orn(Fmoc)-OH	94	FAA8115	Fmoc-L-Lys(Pentynoyl-DIM)-OH	91
DAA1016	Dde-O ₂ Oc-OH	21, 92	FAA1401	Fmoc-L-MeLys(Dde)-OH	99
BNN1350	DETA(BHH*2HCl)	36	FAA7935	Fmoc-L-MeLys(ivDde)-OH	99
BNN1330	DETA(HBH)*2HCl	36, 101	FAA1502	Fmoc-L-Orn(Dde)-OH	97
BNN1360	Di-Boc-Cystamine	77	FAA1503	Fmoc-L-Orn(ivDde)-OH	97
PEG2145	Dnp-NH-PEG(4)-COOH	30	PEG4370	Fmoc-NH-PEG(3)-COOH	26
PEG2150	Dnp-NH-PEG(4)-NHS	30	RL-4410	Fmoc-NH-PEG(3)-DIG-OH	19
BNN1340	DPTA(BHB)*HCl	36	RL-4380	Fmoc-NH-PEG(3)-N ₃	22
RL-2940	Fivemethyl-Lock	10	RL-4400	Fmoc-NH-PEG(3)-NH-Suc-OH	19

[↑ back to content](#)

The Concept of Antibody-Drug Conjugation (ADC)

Permanent Linkers

Cleavable Linkers

Trifunctional Linkers

Cross-Linkers for other Bio Applications

Preparing Carriers for Conjugation

Index

Product code	Product name	Page	Product code	Product name	Page
RL-4390	Fmoc-NH-PEG(3)-NH ₂ *HCl	19	RL-1050	H-PAL-Linker	41
PEG1805	Fmoc-NH-PEG(4)-NHNH-Boc	28	PEG1365	H ₂ N-PEG(2)-CO-OtBu	17
PEG4410	Fmoc-NH-PEG(4)-NHS	30	PEG4980	H ₂ N-PEG(2)-N ₃ *TosOH	18
PEG7810	Fmoc-NH-PEG(4)-TFP	30	PEG3060	H ₂ N-PEG(3)-N ₃	27
RL-3580	Fmoc-NH-SS-Bzl-OH	80	PEG1375	H ₂ N-PEG(4)-CO-OtBu	29
RL-3530	Fmoc-NH-SS-OH	79	PEG1370	H ₂ N-PEG(4)-COOH	29
RL-3540	Fmoc-NH-SS-OpNC	79	PEG1335	H ₂ N-PEG(4)-NHNH-Boc	28
FAA1787	Fmoc-O ₂ Oc-O ₂ Oc-OH	27	PEG1320	H ₂ N-PEG(4)-OH	28
FAA1435	Fmoc-O ₂ Oc-OH	21	RL-3670	Halo-DBCO	105
PEG8150	Fmoc-PEG(4)-Dde	92	RL-3700	Halo-PEG(2)-Azide	105
FAA7190	Fmoc-Spr(oNB)-OH	10	RL-3680	Halo-PEG(2)-NH ₂ *HCl	105
FAA7200	Fmoc-Spr(oNV)-OH	10	RL-3180	Halo-PEG(2)-Suc	105
RL-2200	Fmoc-SS-COOH	78	RL-3710	Halo-PEG(4)-Azide	106
FNN1011	Fmoc-TOTA*HCl	24	RL-3690	Halo-PEG(4)-NH ₂ *HCl	105
FAA5730	Fmoc-TTD-DIG-OH	20	RL-3640	Halo-PEG(5)-azide	105
FAA1568	Fmoc-TTDS-OH	23	RL-1114	HMPB-Linker	41
ADC1060	Fmoc-Val-Ala-PAB	60	PEG1535	HO-PEG(4)-CO-OtBu	29
ADC1670	Fmoc-Val-Ala-PAB-Cl	60	PEG7220	HO-PEG(4)-TFP	29
ADC1410	Fmoc-Val-Ala-PAB-NMeCH ₂ CH ₂ NMe-Boc	60	PEG1970	HS-PEG(4)-COOH	32
ADC1070	Fmoc-Val-Ala-PAB-PNP	60	RL-3840	ICG-CLIP	106
ADC1030	Fmoc-Val-Cit-PAB	68	RL-3830	ICG-SNAP	108
ADC1240	Fmoc-Val-Cit-PAB-NMeCH ₂ CH ₂ NMe-Boc	68	DAA1030	ivDde-D-Lys(Fmoc)-OH	95
ADC1000	Fmoc-Val-Cit-PAB-PNP	68	DAA1018	ivDde-L-Dap(Fmoc)-OH	95
FAA7570	Fmoc2-DAPOA	101	DAA1019	ivDde-L-Lys(Fmoc)-OH	95
RL-1002	FMPB-Linker	40	PTC1020	Lenalidomide	118
RL-2950	Fourmethyl-Lock	10	PEG3590	Lipoamide-PEG(4)-OMe	33
FLL1070	Fullerene C ₆₀ (malonic acid)n	138	MAA1100	Mal-AMCHC-N-Propargylamide	48
FLL1020	Fullerene C ₆₀ (PBM)	139	MAA5400	Mal-AMCHC-OH	46
FLL1010	Fullerene C ₆₀ (PBM)2	139	MAA1000	Mal-AMCHC-OSu	46
FLL1080	Fullerene C ₇₀ (malonic acid)n	139	MAA1020	Mal-beta-Ala-OSu	44
FLL1060	Fullerene C ₇₀ (PBM)	139	ADC1390	Mal-beta-Ala-PEG(4)-Val-Ala-PAB	62
FLL1050	Fullerene C ₇₀ (PBM)2	139	ADC1400	Mal-beta-Ala-PEG(4)-Val-Ala-PAB-PNP	62
FLL1030	Fullerenol C60	139	ADC1220	Mal-beta-Ala-PEG(4)-Val-Cit-PAB	69
FLL1090	Fullerenol C70	140	ADC1230	Mal-beta-Ala-PEG(4)-Val-Cit-PAB-PNP	69
HAA9300	H-Aca-Aca-OH	36	RL-2640	Mal-Bu-NHS	44
PEG8060	H-O ₂ Oc-O ₂ Oc-O ₂ Oc-O ₂ Oc-OH	28	RL-3400	Mal-CH ₂ CH ₂ -N(Me)-CH ₂ -COOH	48
PEG2770	H-O ₂ Oc-O ₂ Oc-O ₂ Oc-OH	28	RL-3450	Mal-CH ₂ CH ₂ -N-(CH ₂ -COOH)2	48
PEG1221	H-O ₂ Oc-O ₂ Oc-OH	27	RL-2650	Mal-cHxHx-NHS	46
PEG2420	H-O ₂ Oc-OH	20	ADC1560	Mal-cyclobutane-1,1-dicarboxamide-Ala-PAB	64
PEG7940	H-O ₂ Oc-OH*HCl	20	ADC1570	Mal-cyclobutane-1,1-dicarboxamide-Ala-PAB-PNP	64
PEG2430	H-O ₂ Oc-OtBu*HCl	20	ADC1460	Mal-cyclobutane-1,1-dicarboxamide-Cit-PAB	70

Product code	Product name	Page	Product code	Product name	Page
ADC1470	Mal-cyclobutane-1,1-dicarboxamide-Cit-PAB-PNP	71	RL-2980	N ₃ -Aca-OSu	37
MAA1060	Mal-D-Dap(Boc)-OH*DCHA	45	PEG7950	N ₃ -AEEA-OK	22
ADC1080	Mal-Dap(Boc)-Val-Ala-PAB-PNP	62	PEG5400	N ₃ -AEEEA*CHA	24
ADC1090	Mal-Dap(Boc)-Val-Cit-PAB-PNP	65	RL-4430	N ₃ -C7H14-COOH	38
RL-3000	Mal-Et-OH	43	HNN1090	N ₃ -Cystamine*HCl	76
RL-2660	Mal-Hx-NHS	45	RL-4360	N ₃ -DABu-Suc-OH	37
MAA1040	Mal-L-Dap(Boc)-OH*DCHA	45	RL-4350	N ₃ -DAPr-Suc-OH	37
MAA1080	Mal-L-Dap(Boc)-OPfp	45	BNN1370	N ₃ -EDA-Suc-OH	38
MAA1120	Mal-L-Dap(Boc)-OSu	46	PEG4900	N ₃ -EEET-OH	23
RL-3430	Mal-N-Boc-Aeg-NHS	48	HAA9330	N ₃ -Gly-Aeg(Fmoc)-OH	102
RL-2780	Mal-NH ₂ *HCl	43	AAA1960	N ₃ -Hx-OH	37
PEG4870	Mal-O ₂ Oc-OH	21, 44	RL-8700	N ₃ -L-Glu(Dde)-OtBu	90
PEG1555	mal-PEG(2)-COOH	18, 44	PEG2790	N ₃ -O ₂ Oc-O ₂ Oc-OH	27
PEG1560	mal-PEG(2)-NHS	18, 44	PEG2780	N ₃ -O ₂ Oc-OH*CHA	22
PEG1485	mal-PEG(3)-mal	26, 49	PEG5390	N ₃ -O ₂ Oc-OtBu	22
RL-3980	Mal-PEG(4)-SNAP	107	RL-4370	N ₃ -PEG(3)-NH-DIG-OH	22
RL-2670	Mal-Pen-NHS	44	PEG3760	N ₃ -PEG(3)-OH	26
ADC1770	Mal-PhAc-PEG(4)-Val-Ala-PAB	63	PEG2345	N ₃ -PEG(4)-COOH	33
ADC1780	Mal-PhAc-PEG(4)-Val-Ala-PAB-PNP	63	PEG5320	N ₃ -PEG(4)-NH ₂	27
ADC1730	Mal-PhAc-Val-Ala-PAB	62	PEG1400	N ₃ -PEG(4)-NHS	33
ADC1740	Mal-PhAc-Val-Ala-PAB-PNP	63	PEG5300	N ₃ -PEG(4)-OH	33
ADC1750	Mal-PhAc-Val-Cit-PAB	65	RL-3280	N ₃ -Pen-Dde	90
ADC1760	Mal-PhAc-Val-Cit-PAB-PNP	65	RL-3290	N ₃ -Pen-Dtpp	90
RL-2680	Mal-PhBu-NHS	47	AAA1970	N ₃ -Pen-OH	37
RL-2690	Mal-PrHx-NHS	45	PEG5000	N ₃ -TFBA-O ₂ Oc	22, 54
RL-3970	Mal-SNAP	107	BNN1150	N ₃ -TOTA	24
RL-4090	Mal-SS-COOH	78	PEG5170	N ₃ -TOTA-Suc	25
ADC1270	MC-Val-Ala-PAB	61	RL-3010	N ₃ Ac-OPhOMe	111
ADC1700	MC-Val-Ala-PAB-Cl	61	PEG4120	NHS-PEG(2)-NHS	17
ADC1280	MC-Val-Ala-PAB-PNP	61	PEG4130	NHS-PEG(3)-NHS	25
ADC1100	MC-Val-Cit-PAB	68	PTC1070	Nimbolide	119
ADC1110	MC-Val-Cit-PAB-PNP	69	PTC1080	Nutlin-3	119
RL-2310	MeTz-PEG(4)-COOH	34	PTC1090	Nutlin-3a	119
RL-2340	MeTz-PEG(4)-mal	34, 43	RL-3550	OPSS-Bzl-OpNC	83
RL-2330	MeTz-PEG(4)-NHS	34	RL-3920	OPSS-Bzl-PAB	82
PEG1740	Mmt-S-PEG(4)-COOH	33	RL-3850	OPSS-Bzl-PAB-OpNC	83
PEG2161	Mtt-NH-PEG(4)-COOH*TEA	30	RL-3500	OPSS-OpNC	82
PEG4650	Mtt-O ₂ Oc-OH*DEA	21	RL-3890	OPSS-PAB	82
PTC1010	N-Methylated pomalidomide	118	RL-3820	OPSS-PAB-OpNC	83
HAA6990	N ₃ -Aca-Aca-OH	37	PEG2230	OPSS-PEG(4)-NHS	30
			RL-4000	OPSS-PEG(4)-SNAP	107

The Concept of Antibody-Drug Conjugation (ADC)

Permanent Linkers

Cleavable Linkers

Trifunctional Linkers

Cross-Linkers for other Bio Applications

Preparing Carriers for Conjugation

Index

[↑ back to content](#)

Product code	Product name	Page	Product code	Product name	Page
RL-3990	OPSS-SNAP	107	PTC1190	Pomalidomide-PEG5-COOH	121
RL-4040	PFB-mercaptopropionyl-AEEA	34, 132	PTC1270	Pomalidomide-PEG5-NH ₂ hydrochloride	122
RL-4030	PFB-mercaptopropionyl-PEG3-N ₃	34, 132	PTC1140	Pomalidomide-PEG6-butyl COOH	120
RL-4050	PFB-mercaptopropionyl-TOTA-Biotin	34, 132	PTC1300	Pomalidomide-PEG6-butyl-NH ₂ hydrochloride	122
RL-2920	Photo-Benzoin acid	53	PTC1200	Pomalidomide-PEG6-COOH	121
RL-2930	Photo-Benzylamine*HCl	53	PTC1280	Pomalidomide-PEG6-NH ₂ hydrochloride	122
RL-2910	Photo-Butylamine	53	ADC1600	Propargyl-cyclobutane-1,1-dicarboxamide-Ala-PAB	63
RL-3410	Photo-Click-Heptanoic acid	53	ADC1610	Propargyl-cyclobutane-1,1-dicarboxamide-Ala-PAB-PNP	64
RL-2900	Photo-Hexanoic acid	53	ADC1500	Propargyl-cyclobutane-1,1-dicarboxamide-Cit-PAB	70
RL-2890	Photo-Pentanoic acid	53	ADC1510	Propargyl-cyclobutane-1,1-dicarboxamide-Cit-PAB-PNP	70
RL-2970	Photo-Trimethyl-Lock	10	PEG8170	Propargyl-PEG(5)-COOH	31
PEG5080	Phth-NO-PEG(4)-NHS	31	RL-3835	SNAP-acid	106
RL-4160	pNCO-SS-OpNC	80	MAA1050	Sulfo-SMCC	46
PAA1050	Poc-O ₂ C-OH*DCHA	23	TCO1050	TCO-PEG(3)-mal	26, 49
PTC1520	Pomalidomid- PEG1-N ₃	126	TCO1040	TCO-PEG(4)-COOH	32
PTC1530	Pomalidomid- PEG2-N ₃	126	TCO1010	TCO-PEG(4)-NHS	32
PTC1540	Pomalidomid- PEG3-N ₃	126	RL-4140	TetraMe-Dioxoborolane-(OpNC) ₂	15
PTC1710	Pomalidomid-C6-PEG1-C3-PEG1-butyl-I	128	RL-4130	TetraMe-Dioxoborolane-OpNC	16
PTC1570	Pomalidomid-C6-PEG1-C3-PEG1-butyl-N ₃	127	PEG7010	Trt-S-EEE	24
PTC1700	Pomalidomid-C6-PEG3-butyl-I	128	PEG6730	Trt-S-EEEE	25
PTC1560	Pomalidomid-C6-PEG3-butyl-N ₃	127	PEG6710	Trt-S-PEG(4)-COOH*H ₂ O	33
PTC1690	Pomalidomid-PEG2-butyl-I	128	LS-3960	Tyramide-SS-amine*HCl	82
PTC1720	Pomalidomid-PEG6-butyl-I	129	PTC1060	VH298	118
PTC1580	Pomalidomid-PEG6-butyl-N ₃	127	PEG1495	Z-NH-PEG(4)-COOH	31
PTC1000	Pomalidomide	117	ZAA1186	Z-O ₂ C-OH*DCHA	23
PTC1100	Pomalidomide-C3-COOH	119	ZNN1120	Z-TOTA	19
PTC1110	Pomalidomide-C6-COOH	119			
PTC1120	Pomalidomide-C9-COOH	120			
PTC1400	Pomalidomide-PEG1-Alkyne	124			
PTC1150	Pomalidomide-PEG1-COOH	120			
PTC1230	Pomalidomide-PEG1-NH ₂ hydrochloride	121			
PTC1410	Pomalidomide-PEG2-Alkyne	125			
PTC1130	Pomalidomide-PEG2-butyl COOH	120			
PTC1160	Pomalidomide-PEG2-COOH	120			
PTC1240	Pomalidomide-PEG2-NH ₂ hydrochloride	122			
PTC1420	Pomalidomide-PEG3-Alkyne	125			
PTC1170	Pomalidomide-PEG3-COOH	120			
PTC1250	Pomalidomide-PEG3-NH ₂ hydrochloride	122			
PTC1180	Pomalidomide-PEG4-COOH	121			
PTC1260	Pomalidomide-PEG4-NH ₂ hydrochloride	122			
PTC1440	Pomalidomide-PEG5-Alkyne	125			

Index

Iris Biotech GmbH • Email: info@iris-biotech.de • www.iris-biotech.de • **Empowering Peptide Innovation** 156

Notes

Code of Conduct

As business activity of Iris Biotech GmbH impacts people's lives and health, it must be operated in ethical and correct manner and act with integrity and responsibility. To ensure high ethical standards and fair business practices, Iris Biotech GmbH applies an integrated policy known as its Code of Conduct.

In 2001 Iris Biotech GmbH was founded just at the beginning of the Biotech movement and the first remarkable breakthrough of biotech pharma products. Although the biotech field is rather young compared to other industries we believe on long-term business, a good partnership between our business partners and Iris Biotech GmbH and a good reputation. It is our duty as well as our responsibility to maintain and to extend this over the next generations – based on the principles of an honourable and prudent tradesman which based upon the concept of honourable entrepreneurship.

This Code of Conduct has been developed following the “Voluntary Guidelines for Manufacturers of Fine Chemical Intermediates and Active Ingredients” issued by AIME (Agrochemical & Intermediates Manufacturers in Europe) and the requirements of some of our business associates.

Iris Biotech GmbH commits to hold this Code of Conduct and to include and apply its principles in the management system and the company policies.

Ethics

Iris Biotech GmbH undertakes business in an ethical manner and acts with integrity. All corruption, extortion and embezzlement are prohibited. We do not pay or accept bribes or participate in other illegal inducements in business or government relationships. We conduct our business in compliance with all applicable anti-trust laws. Employees are encouraged to report concerns or illegal activities in the workplace, without threat of reprisal, intimidation or harassment.

Labour

Iris Biotech GmbH is committed to uphold the human rights of workers and to treat them with dignity and respect. Child labour, workplace harassment, discrimination, and harsh and inhumane treatment are prohibited. Iris Biotech GmbH respects the rights of the employees to associate freely, join or not join labour unions, seek representation and join workers' councils. Employees are paid and their working timetable is established according to applicable wage and labour laws. Employees are able to communicate openly with management regarding working conditions without threat of reprisal, intimidation or harassment.

General Policies

Contracts and Secrecy Agreements are binding and the confidential information received is only used for intended purposes. Clear management and organizational structures exist to provide efficient normal working and to address problems quickly. Know-how is protected and intellectual property is respected.

Health and Safety

Iris Biotech GmbH provides a safe and healthy working environment to the employees and protects them from overexposure to chemical and physical hazards. Products are produced, stored and shipped under the guidelines of the relevant chemical and safety legislation. Risks and emergency scenarios are identified and evaluated, and their possible impact is minimized by implementing emergency plans and written procedures. Safety information regarding hazardous materials is available to educate, train and protect workers from hazards. Preventive equipment and facilities maintenance is performed at suitable periods to reduce potential hazards. Employees are regularly trained in health and safety matters and are informed about product properties and risk classification when it is required.

Environment

Iris Biotech GmbH operates in an environmentally responsible and efficient manner, minimizing adverse impacts on the environment. Waste streams are managed to ensure a safe handling, movement, storage, recycling and reuse, before and after being generated. Systems to prevent and mitigate accidental spills and releases to the environment are in place. All required environmental permits and licenses are obtained and their operational and reporting requirements are complied with.

Production and Quality Management

A quality management system following the Good Distribution Practices (GDP rules) of Active Pharmaceutical Ingredients is established covering all the aspects of the worldwide distribution of products. Regular audits are performed to evaluate the efficiency and fulfilling of the quality system. Process controls to provide reproducible product quality are established. There are preventive maintenance procedures to ensure plant reliability and the lowest risk of failure. Staff is trained periodically about GMP and GDP rules. Procedures are established and installations are designed to avoid cross contamination. Batch and analytical records are kept for inspection and audit purposes for suitable periods according guidelines.

Research and Development

Research and development staff education is appropriate to their functional activity and they are trained to develop, optimize and scale-up the processes. Intellectual property is respected and know-how protected. Development of manufacturing processes reflects the principles of the Green Chemistry according to the American Chemical Society Green Chemistry Institute. Animal testing is not used unless alternatives are not scientifically valid or accepted by regulators. If animal testing is carried out, animals are treated so that pain and stress are minimized.

Terms and Conditions of Sales

All orders placed by a buyer are accepted and all contracts are made subject to the terms which shall prevail and be effective notwithstanding any variations or additions contained in any order or other document submitted by the buyer. No modification of these terms shall be binding upon Iris Biotech GmbH unless made in writing by an authorised representative of Iris Biotech GmbH.

Placing of Orders

Every order made by the buyer shall be deemed an offer by the buyer to purchase products from Iris Biotech GmbH and will not be binding on Iris Biotech GmbH until a duly authorised representative of Iris Biotech GmbH has accepted the offer made by the buyer. Iris Biotech GmbH may accept orders from commercial, educational or government organisations, but not from private individuals and Iris Biotech GmbH reserves the right to insist on a written order and/or references from the buyer before proceeding.

There is no minimum order value. At the time of acceptance of an order Iris Biotech GmbH will either arrange prompt despatch from stock or the manufacture/acquisition of material to satisfy the order. In the event of the latter Iris Biotech GmbH will indicate an estimated delivery date. In addition to all its other rights Iris Biotech GmbH reserves the right to refuse the subsequent cancellation of the order if Iris Biotech GmbH expects to deliver the product on or prior to the estimated delivery date. Time shall not be of the essence in respect of delivery of the products. If Iris Biotech GmbH is unable to deliver any products by reason of any circumstances beyond its reasonable control („Force Majeure“) then the period for delivery shall be extended by the time lost due to such Force Majeure. Details of Force Majeure will be forwarded by Iris Biotech GmbH to the buyer as soon as reasonably practicable.

Prices, Quotations and Payments

Prices are subject to change. For the avoidance of doubt, the price advised by Iris Biotech GmbH at the time of the buyer placing the order shall supersede any previous price indications. The buyer must contact the local office of Iris Biotech GmbH before ordering if further information is required. Unless otherwise agreed by the buyer and Iris Biotech GmbH, the price shall be for delivery ex-works. In the event that the buyer requires delivery of the products otherwise than ex-works the buyer should contact the local office of Iris Biotech GmbH in order to detail its requirements. Iris Biotech GmbH shall, at its discretion, arrange the buyer's delivery requirements including, without limitation, transit insurance, the mode of transit (Iris Biotech GmbH reserves the right to vary the mode of transit if any regulations or other relevant considerations so require) and any special packaging requirements (including cylinders). For the avoidance of doubt all costs of delivery and packaging in accordance with the buyer's requests over and above that of delivery in standard packaging ex-works shall be for the buyer's account unless otherwise agreed by both parties. Incoterms 2020 shall apply. Any tax, duty or charge imposed by governmental authority or otherwise and any other applicable taxes, duties or charges shall be for the buyer's account. Iris Biotech GmbH may, on request and where possible, provide quotations for multiple packs or bulk quantities, and non-listed items. Irrespective of the type of request or means of response all quotations must be accepted by the buyer without condition and in writing before an order will be accepted by Iris Biotech GmbH. Unless agreed in writing on different terms, quotations are valid for 30 days from the date thereof. Payment terms are net 30 days from invoice date unless otherwise agreed in writing. Iris Biotech GmbH reserves the right to request advance payment at its discretion. For overseas transactions the buyer shall pay all the banking charges of Iris Biotech GmbH. The buyer shall not

be entitled to withhold or set-off payment for the products for any reason whatsoever. Government/ Corporate Visa and MasterCard (and other such credit cards) may be accepted on approved accounts for payment of the products. Personal credit cards are not acceptable. Failure to comply with the terms of payment of Iris Biotech GmbH shall constitute default without reminder. In these circumstances Iris Biotech GmbH may (without prejudice to any other of its rights under these terms) charge interest to accrue on a daily basis at the rate of 2% per month from the date upon which payment falls due to the actual date of payment (such interest shall be paid monthly). If the buyer shall fail to fulfil the payment terms in respect of any invoice of Iris Biotech GmbH Iris Biotech GmbH may demand payment of all outstanding balances from the buyer whether due or not and/or cancel all outstanding orders and/or decline to make further deliveries or provision of services except upon receipt of cash or satisfactory securities. Until payment by the buyer in full of the price and any other monies due to Iris Biotech GmbH in respect of all other products or services supplied or agreed to be supplied by Iris Biotech GmbH to the buyer (including but without limitation any costs of delivery) the property in the products shall remain vested in Iris Biotech GmbH.

Shipping, Packaging and Returns

The buyer shall inspect goods immediately on receipt and inform Iris Biotech GmbH of any shortage or damage within five days. Quality problems must be notified within ten days of receipt. Goods must not be returned without prior written authorisation of Iris Biotech GmbH. Iris Biotech GmbH shall at its sole discretion replace the defective products (or parts thereof) free of charge or refund the price (or proportionate price) to buyer. Opened or damaged containers cannot be returned by the buyer without the written prior agreement of Iris Biotech GmbH. In the case of agreed damaged containers which cannot be so returned, the buyer assumes responsibility for the safe disposal of such containers in accordance with all applicable laws.

Product Quality, Specifications and Technical Information

Products are analysed in the Quality Control laboratories of Iris Biotech GmbH's production partners by methods and procedures which Iris Biotech GmbH considers appropriate. In the event of any dispute concerning reported discrepancies arising from the buyer's analytical results, determined by the buyer's own analytical procedures, Iris Biotech GmbH reserves the right to rely on the results of own analytical methods of Iris Biotech GmbH. Certificates of Analysis or Certificates of Conformity are available at the discretion of Iris Biotech GmbH for bulk orders but not normally for prepack orders. Iris Biotech GmbH reserves the right to make a charge for such certification. Specifications may change and reasonable variation from any value listed should not form the basis of a dispute. Any supply by Iris Biotech GmbH of bespoke or custom product for a buyer shall be to a specification agreed by both parties in writing. Technical information, provided orally, in writing, or by electronic means by or on behalf of Iris Biotech GmbH, including any descriptions, references, illustrations or diagrams in any catalogue or brochure, is provided for guidance purposes only and is subject to change.

Safety

All chemicals should be handled only by competent, suitably trained persons, familiar with laboratory procedures and potential chemical hazards. The burden of safe use of the products of Iris Biotech GmbH vests in the buyer. The buyer assumes all responsibility for warning his employees, and any persons who might reasonably be expected to come into contact with the products, of all risks to person and property in any way connected with the products and for instructing them in their safe handling and use. The buyer also assumes the responsibility for the safe disposal of all products in accordance with all applicable laws.

Uses, Warranties and Liabilities

All products of Iris Biotech GmbH are intended for laboratory research purposes and unless otherwise stated on product labels, in the catalogue and product information sheet of Iris Biotech GmbH or in other literature furnished to the buyer, are not to be used for any other purposes, including but not limited to use as or as components in drugs for human or animal use, medical devices, cosmetics, food additives, household chemicals, agricultural or horticultural products or pesticides. Iris Biotech GmbH offers no warranty regarding the fitness of any product for a particular purpose and shall not be responsible for any loss or damage whatsoever arising there from. No warranty or representation is given by Iris Biotech GmbH that the products do not infringe any letters patent, trademarks, registered designs or other industrial rights. The buyer further warrants to Iris Biotech GmbH that any use of the products in the United States of America shall not result in the products becoming adulterated or misbranded within the meaning of the Federal Food, Drug and Cosmetic Act (or such equivalent legislation in force in the buyer's jurisdiction) and shall not be materials which may not, under sections 404, 505 or 512 of the Act, be introduced into interstate commerce. The buyer acknowledges that, since the products of Iris Biotech GmbH are intended for research purposes, they may not be on the Toxic Substances Control Act 1976 („TSCA“) inventory. The buyer warrants that it shall ensure that the products are approved for use under the TSCA (or such other equivalent legislation in force in the buyer's jurisdiction), if applicable. The buyer shall be responsible for complying with any legislation or regulations governing the use of the products and their importation into the country of destination (for the avoidance of doubt to include, without limitation, the TSCA and all its amendments, all EINECS, ELINCS and NONS regulations). If any licence or consent of any government or other authority shall be required for the acquisition, carriage or use of the products by the buyer the buyer shall obtain the same at its own expense and if necessary produce evidence of the same to Iris Biotech GmbH on demand. Failure to do so shall not entitle the buyer to withhold or delay payment. Any additional expenses or charges incurred by Iris Biotech GmbH resulting from such failure shall be for the buyer's account. Save for death or personal injury caused by negligence of Iris Biotech GmbH, sole obligation of Iris Biotech GmbH and buyer's exclusive remedy with respect to the products proved to the satisfaction of Iris Biotech GmbH to be defective or products incorrectly supplied shall be to accept the return of said products to Iris Biotech GmbH for refund of the actual purchase price paid by the buyer (or proportionate part thereof), or replacement of the defective product (or part thereof) with alternative product. Iris Biotech GmbH shall have no liability to the buyer under or arising directly or indirectly out of or otherwise in connection with the supply of products by Iris Biotech GmbH to the buyer and/or their re-sale or use by the buyer or for any product, process or services of the buyer which in any way comprises the product in contract tort (including negligence or breach of statutory duty) or otherwise for pure economic loss, loss of profit, business, reputation, depletion of brand, contracts, revenues or anticipated savings or for any special indirect or consequential damage or loss of any nature except as may otherwise be expressly provided for in these terms. All implied warranties, terms and representations in respect of the products (whether implied by statute or otherwise) are excluded to the fullest extent permitted by law. The buyer shall indemnify Iris Biotech GmbH for and against any and all losses, damages and expenses, including legal fees and other costs of defending any action, that Iris Biotech GmbH may sustain or incur as a result of any act or omission by the buyer, its officers, agents or employees, its successors or assignees, its customers or all other third parties, whether direct or indirect, in connection with the use of any product. For the avoidance of doubt and in the event that Iris Biotech GmbH supplies bespoke or custom product to the buyer's design or specification, this indemnity shall extend to include any claim by a third party that the manufacture of the product for the buyer or the use of the product by the buyer infringes the intellectual property rights of any third party.

General

Iris Biotech GmbH shall be entitled to assign or sub-contract all or any of its rights and obligations hereunder. The buyer shall not be entitled to assign, transfer, sub-contract or otherwise delegate any of its rights or obligations hereunder. Any delay or forbearance by Iris Biotech GmbH in exercising any right or remedy under these terms shall not constitute a waiver of such right or remedy. If any provision of these terms is held by any competent authority to be invalid or unenforceable in whole or in part the validity of the other provisions of these terms and the remainder of the provision in question shall not be affected. These terms shall be governed by German Law and the German Courts shall have exclusive jurisdiction for the hearing of any dispute between the parties save in relation to enforcement where the jurisdiction of the German Courts shall be non-exclusive.



Get in Contact



Iris Biotech GmbH
Adalbert-Zoellner-Str. 1
95615 Marktredwitz
Germany

+49 (0) 9231 97121-0
+49 (0) 9231 97121-99
info@iris-biotech.de
www.iris-biotech.de

Distribution Partners

The list contains the current distributors of Iris Biotech in different regions of the world. The latest list of distribution partners and contact details is available at: www.iris-biotech.de/distribution-partner

China:

Chengdu Yoo Technology Co., Ltd.

Japan:

BizCom Japan, Inc.

Shigematsu & Co., Ltd

Cosmo Bio Co., Ltd.

USA & Canada:

Peptide Solutions LLC

India, Bangladesh, Oman, Sri Lanka, United Arab Emirates:

Sumit Biosciences Pvt Ltd.

Empowering Peptide Innovation