



Iris
Biotech

LIGATION TECHNOLOGIES



Version: IB7_4

Empowering Peptide Innovation

With this guiding theme in mind, Iris Biotech's mission is to support researchers by supplying

- innovative technologies,
- rare compounds,
- as well as a broad portfolio on standard consumables,

available in flexible quantities from small scale to bulk quantities. To fulfill our dedication "Empowering Peptide Innovation", we are attending various conferences, symposia, and exhibitions each year. This allows us to remain in direct contact with scientists all over the world, both from academia and industry, to exchange knowledge, and to gather new ideas to tackle your current challenges.

Guided by our dedication to provide

- competent service,
- as well as novel substances and
- latest technologies,

Iris Biotech is your trusted partner for the world of peptides, while having strong expertise in associated disciplines. Thus, our portfolio comprises reagents and tools for the synthesis and modification of peptides, e.g., amino acids, resins and solvents but also for related technologies such as drug delivery, linkerology® and life sciences.

Owed to the growing demand for tailor-made compounds, our portfolio is fine-tuned by our custom synthesis service at Iris Biotech Laboratories. Our skilled scientists offer profound expertise in

- *de novo* route development,
- upscaling towards larger scale production,
- as well as synthesis optimization for increased efficiency.

Examples are the synthesis of rare chiral building blocks, unnatural amino acid derivatives, sophisticated orthogonal protecting groups, heterocycles, building blocks for nucleotides, PEGs and PEG-analogs as well as specific linkers for controlled drug delivery and release.



Amino Acids



Building Blocks



Life Sciences



Drug Delivery



Reagents



Resins



Linkerology®



Click Chemistry

Portfolio Overview

Peptide Synthesis and Modification

(Protected) Amino Acids

Standards such as Fmoc-D/L-AAA and Boc-D/L-AAA, Smoc amino acids for peptide synthesis in water, variety of protecting groups (e.g., Pbf, Trt, ^tBu, Bzl, Acn, Mob, SIT, Phacn, Allocam, Mmt), unusual amino acids, fluorinated derivatives, substituted prolines, arginine analogs

Building Blocks

Amino alcohols, amino aldehydes, diamines and hydrazines, (pseudoproline) dipeptides, polyamines and spermines, fatty acid derivatives, peptide nucleic acids (PNAs)

Reagents

Coupling reagents, solvents and scavengers, protecting groups

Resins

Preloaded resins (e.g., based on Trityl, TCP, TentaGel, Methoxybenzhydryl, Merrifield, PAM, Rink, Wang), scavenger resins, hydrazone resins, poly(acrylamide) resins, Cyclover

Linkerology® and Drug Delivery

Linkers for Solid Phase Peptide Synthesis

Cleavable Linkers

Val-Ala-based, Val-Cit-based, disulfide-based, Dde-helping hands, pH-sensitive linkers

Photo-Activatable Linkers

Functionalized Linkers

Clickable linkers, trifunctional linkers, linkers with maleimide function, cross-linkers, selective N-term acylation and biotinylation, 5HP2O

PROTACs

Ligands, linkers & modules

Fullerenes, Poly(2-oxazolines), Dextrans & Plant-Derived Cholesterol

Superparamagnetic Iron Oxide Nanoparticles

Poly-Amino Acids

Poly-Arg, Poly-Glu, Poly-Lys, Poly-Orn, Poly-Sar

PEGylation

Branched PEGylating reagents, (amino-)PEG-acids, PEG-amines & hydrazides & guanidines, reagents for Click-conjugation, Biotin-PEG-reagents, PEG-thiols, PEG-maleimides, other PEGylating reagents

Life Sciences

Biotinylation Reagents

Carbohydrates

Galactose, Glucose, Mannose, Xylose and others

Drug Metabolites

Peptides

Substrates & Inhibitors

E.g., protein kinase inhibitors, substrates for fusion (Halo/Snap/Clip)-tagged proteins

Natural Products

Dyes and Fluorescent Labels

E.g., ICG, AMC, DAPI

Maillard & Amadori Reaction Products

Large portfolio of derivatives useful as standards for food, pharma and cosmetics industry


Vitamins

Custom Synthesis

Your project requires a compound not listed in our portfolio?
Get in contact and inquire about our custom synthesis capabilities.

Our experienced scientists are excited to accept your synthetic challenge!

In such cases, your request undergoes the following stages:



Step-by-Step Analysis

- Customer's demands

Process Evaluation

- Detailed literature review
- Synthetic possibilities

Strategy Development

- Protocol development
- Method development and validation
- Customized synthesis

Quality Consistency

- Identity confirmation
- Purity verification

Our Service Promise

All our services are based on high standards, transparency & documentation, trust, honesty & confidentiality, as well as the required know-how.

High Standards

- Values: sustainability & responsibility
- State-of-the-art equipment & latest technologies
- High quality standards
- Qualified suppliers & regular audits

Transparency & Documentation

- Talk to our specialists – customer care
- Certificates of analysis & origin
- Impurity profiling
- Safety data sheets
- Analytical and process reports

Trust, Honesty & Confidentiality

- Intergenerational business valuing partnerships
- Meeting the customer's expectations
- Integrity towards our customers

Our Know-How

- One-step reactions & complex multi-step synthesis
- Scalability from mg to kg quantities
- Route scouting



Table of Contents

1. Description of the Technology	1
---	----------

2. Tools for the Formation of Thioesters	7
2.1. Dawson Linkers	7
2.2. Hydrazone Resins	10

3. Building Blocks for Ligation	17
3.1. Cysteine Building Blocks for Ligation	17
3.2. Ligation at the Position of Glycine	35
3.3. Building Blocks and Reagents for Ligation with Seleno Amino Acids	37
3.4. Oxime Ligation using Aminooxy-Amino Acids	41
3.5. Activated Cysteine-Based Protein Ligation (ACPL)	43
3.6. Nitrile-Aminothiol Conjugation (NATC) as Biocompatible Ligation Technique	44

4. Enzyme-mediated Ligation	48
------------------------------------	-----------

5. Examples	50
5.1. Kinetically Controlled Four-Fragment Ligation for the Convergent Chemical Synthesis of Proteins	50
5.2. Three-Fragment Synthesis of a 116mer with Desulfurization	51

Index	52
Code of Conduct	57
Terms and Conditions of Sales	59


1. Description of the Technology

Background

Peptides are typically synthesized via solid phase (peptide) synthesis (SP(P)S). The main advantage in comparison to classical solution phase synthesis is the fast and easy separation of the desired product from excess reagents by filtration. Typically, in SP(P)S, the molecule being synthesized (e.g., a growing peptide chain) is attached to an insoluble solid support, while reagents are added to the suspension. This setup enables the removal of excess reagents and dissolved byproducts via simple disposal of the reaction solution. Filtration followed by washing of the resin with various solvents leaves the purified protected peptide on the carrier. Consequently, an excess of reagents can be employed in order to shorten reaction times and allow for a quantitative turnover of the substrate, which in turn leads to higher yields and less side products.

The Main Advantages of SPPS

1. Fast work-up through easy separation of solid support from dissolved reactants and by products by filtration, and multiple, rapid washing steps.
2. Improved reaction times, turnover and yield by use of excess amounts of reagents.
3. The syntheses can easily be automatized.
4. Toxic or hazardous materials can be handled safely while being attached to the resin.
5. Minimal physical product loss.
6. Pseudo-dilution phenomena on an individual bead can enable cyclization and avoid formation of dimers.



Two reaction steps have to be carried out per amino acid of a peptide chain. I.e., the synthesis of a 50mer requires 100 chemical steps. Assuming, that the yield per chemical step is 98%, which in general means is rather high, there is only 13% yield of the target sequence left (Fig. 1). Nevertheless, within a long sequence, the probability that difficult fragments are appearing is rather high, which will lower the overall yield of the synthesis. In contrast, an extremely high average yield of 99.9% per coupling step will result in 90% global yield. As this is rather unrealistic, new chemical methodologies such as fragment condensation and native chemical ligation must be developed to facilitate access to longer peptides of 100 amino acids and more.

Tab. 1: Final yield in SPPS depending on the yield per step.

# of Steps	Yield per Step/Final Yield			
	98.0%	99.0%	99.5%	99.9%
10	82%	90%	95%	99%
20	67%	82%	90%	98%
30	55%	74%	86%	97%
40	45%	67%	82%	96%
50	36%	61%	78%	95%
60	30%	55%	74%	94%
70	24%	49%	70%	93%
80	20%	45%	67%	92%
90	16%	40%	64%	91%
100	13%	37%	61%	90%
110	11%	33%	58%	90%
120	9%	30%	55%	89%
130	7%	27%	52%	88%
140	6%	24%	50%	87%
150	5%	22%	47%	86%
160	4%	20%	45%	85%
170	3%	18%	43%	84%
180	3%	16%	41%	84%
190	2%	15%	39%	83%
200	2%	13%	37%	82%

Description of the Technology

Tools for the Formation of Thioesters

Building Blocks for Ligation

Enzyme-mediated Ligation


Examples

Index

[↑ back to content](#)

Native Chemical Ligation

Native chemical ligation (NCL) of unprotected peptide segments involves the reaction between a first peptide fragment α -thioester and a second peptide fragment, which carries a cysteine on the N-terminus, to yield a product with a native amide bond at the ligation site. Peptide- α -thioalkyl esters are commonly used because of their ease of preparation. These thioalkyl esters are rather unreactive, so the ligation reaction is catalyzed by *in situ* transthioesterification with thiol additives, which are either a mixture of thiophenol/benzyl mercaptan or other alkanethiols. Despite the use of this thiol additive, ligation reactions typically take 24–48 h.




Requirements

1. Unprotected peptide fragment as C-terminal thioester.
2. Cysteine on N-terminus of a second unprotected peptide fragment.

Limitations

1. Cysteine appears very rarely in natural sequences.
2. Multiple, consecutive ligations are hardly possible.
3. Cysteine is prone to racemization.
4. Slow kinetics, as ligations normally take 24-48 h; hence a good thiol additive, like 4-mercaptophenylacetic acid (MPAA), is necessary.

The main advantage of this ligation technology vs. fragment condensation is the usually good solubility of unprotected peptide fragments compared to fully protected peptides, which have very limited solubility. One main challenge, however, is that the presence of cysteine within a peptide sequence is a prerequisite for NCL, but at the same time cysteine is the least abundant of all proteinogenic amino acids since it occurs with an average frequency of only 1.4% within natural sequences. This becomes even more of an obstacle if multiple consecutive ligations are necessary to construct the target sequence. Hence, sophisticated variations must be implemented in the absence of cysteine.



Tab. 2: General information - Elements in NCL, their electronegativity, atomic radius, and electron configuration.

Element	Electronegativity	Atomic Radius	Electron Configuration
O	3.4	0.60 ⁰ Å, 1.38 ⁻² ₄ Å	[He] 2s ² 2p ⁴
S	2.6	1.04 ⁰ Å, 1.84 ⁻² ₆ Å	[Ne] 3s ² 3p ⁴
Se	2.6	1.16 ⁰ Å, 1.98 ⁻² ₆ Å	[Ar] 3d ¹⁰ 4s ² 4p ⁴

Tab. 3: Substituted alcohols, thiols and selenols and their respective pK_a values.

Structure	Leaving Group	pK_a
	Alkyl Alcohol	16-18
	Alkyl Thiol	10-11
	Alkyl Selenol	7-8
	Phenol	9.89
	Thiophenol	6.6
	Cys S-H	8.3
	Sec Se-H	5.2
	Alkyl Carboxylic Acids	~4

The mechanistic details of NCL have been studied by several groups in order to find appropriate additives. Substituted thiophenols with $pK_a > 6$ were found to best combine the ability to exchange rapidly and completely with thioalkyl esters, and to then act as efficient leaving groups in the reaction of the peptide thioester with the thiol side-chain of an N-terminal cysteine. 4-mercaptophenylacetic acid (MPAA), a non-malodorous, water-soluble thiol, proved to be a highly effective and practical additive. MPAA reacts orders of magnitude faster than other thiols in model studies of NCL and in the synthesis of small proteins, as demonstrated in the synthesis of turkey ovomucoid third domain (OMTKY3). MPAA should find broad use in native chemical ligation and in the total synthesis of proteins.

The need for a cysteine residue on the N-terminus of one of the reacting peptide fragments has motivated the development of β -, γ -, and δ -thiolated variants of other amino acids, as well as thiol-containing auxiliaries that can be used as cysteine

surrogates in ligation chemistry. After ligation reactions at these residues, the thiol auxiliary is desulfurized (usually by means of radical-based protocols) to yield native polypeptide products. However, this transformation is not chemoselective in the presence of other unprotected cysteine residues that might be found elsewhere in the sequence. This limitation of desulfurization chemistry has led chemists to expand the native chemical ligation toolbox to include the 21st amino acid, selenocysteine (Sec), as well as various non-natural selenoamino acids (to date specifically based on proline and phenylalanine). The key advantage of carrying out ligation chemistry with selenoamino acids rather than thioamino acids is that chemoselective deselenization can be performed under mild conditions (typically with a phosphine reductant and a hydrogen atom source) that do not affect unprotected cysteine residues. Payne *et al.* and Dery *et al.* have demonstrated that ligation products can be subjected to oxidative deselenization to afford serine in place of selenocysteine at the ligation junction.




Fig. 4: Differences in properties and capabilities using either sulfur or selenium in native chemical ligation.

References:


- [13] Protein synthesis by chemical ligation of unprotected peptides in aqueous solution; T. W. Muir, P. E. Dawson, S. B. H. Kent; **Solid-Phase Peptide Synthesis** 1997; **289**: Academic Press, 266-298.
- Accelerated Protein Synthesis via One-Pot Ligation-Deselenization Chemistry; N. J. Mitchell, J. Sayers, S. S. Kulkarni, D. Clayton, A. M. Goldys, J. Ripoll-Rozada, P. J. Barbosa Pereira, B. Chan, L. Radom, R. J. Payne; **Chem** 2017; **2**: 703-715. <https://doi.org/10.1016/j.chempr.2017.04.003>
- Chemical synthesis of proteins; B. L. Nilsson, M. B. Soellner, R. T. Raines; **Annu Rev Biophys Biomol Struct** 2005; **34**: 91-118. <https://doi.org/10.1146/annurev.biophys.34.040204.144700>
- Construction of Challenging Proline-Proline Junctions via Diselenide-Selenoester Ligation Chemistry; J. Sayers, P. M. T. Karpati, N. J. Mitchell, A. M. Goldys, S. M. Kwong, N. Firth, B. Chan, R. J. Payne; **J Am Chem Soc** 2018; **140**: 13327-13334. <https://doi.org/10.1021/jacs.8b07877>
- Expressed protein ligation: a general method for protein engineering; T. W. Muir, D. Sondhi, P. A. Cole; **Proc Natl Acad Sci U S A** 1998; **95**: 6705-6710. <https://doi.org/10.1073/pnas.95.12.6705>
- Free-radical-based, specific desulfurization of cysteine: a powerful advance in the synthesis of polypeptides and glycopolypeptides; Q. Wan, S. J. Danishefsky; **Angew Chem Int Ed Engl** 2007; **46**: 9248-9252. <https://doi.org/10.1002/anie.200704195>
- Insights into the mechanism and catalysis of the native chemical ligation reaction; E. C. Johnson, S. B. Kent; **J Am Chem Soc** 2006; **128**: 6640-6646. <https://doi.org/10.1021/ja058344i>
- Native Chemical Ligation of Peptides and Proteins; P. A. Cistrone, M. J. Bird, D. T. Flood, A. P. Silvestri, J. C. J. Hintzen, D. A. Thompson, P. E. Dawson; **Curr Protoc Chem Biol** 2019; **11**: e61. <https://doi.org/10.1002/cpch.61>
- Oxidative Deselenization of Selenocysteine: Applications for Programmed Ligation at Serine; L. R. Malins, N. J. Mitchell, S. McGowan, R. J. Payne; **Angew Chem Int Ed Engl** 2015; **54**: 12716-12721. <https://doi.org/10.1002/anie.201504639>
- Protein synthesis by native chemical ligation: expanded scope by using straightforward methodology; T. M. Hackeng, J. H. Griffin, P. E. Dawson; **Proc Natl Acad Sci U S A** 1999; **96**: 10068-10073. <https://doi.org/10.1073/pnas.96.18.10068>
- Rapid additive-free selenocystine-selenoester peptide ligation; N. J. Mitchell, L. R. Malins, X. Liu, R. E. Thompson, B. Chan, L. Radom, R. J. Payne; **J Am Chem Soc** 2015; **137**: 14011-14014. <https://doi.org/10.1021/jacs.5b07237>
- Selenopeptide chemistry; M. Muttenthaler, P. F. Alewood; **J Pept Sci** 2008; **14**: 1223-1239. <https://doi.org/10.1002/psc.1075>
- Syntheses and biological activities of selenium analogs of alpha-rat atrial natriuretic peptide; T. Koide, H. Itoh, A. Otaka, M. Furuya, Y. Kitajima, N. Fujii; **Chem Pharm Bull (Tokyo)** 1993; **41**: 1596-1600. <https://doi.org/10.1248/cpb.41.1596>
- Synthesis of proteins by native chemical ligation; P. E. Dawson, T. W. Muir, I. Clark-Lewis, S. B. Kent; **Science** 1994; **266**: 776-779. <https://doi.org/10.1126/science.7973629>
- Synthesis of selenocysteine peptides and their oxidation to diselenide-bridged compounds; D. Besse, L. Moroder; **Journal of Peptide Science** 1997; **3**: 442-453. [https://doi.org/10.1002/\(sici\)1099-1387\(199711\)3:6<442::Aid-psc122>3.0.Co;2-2](https://doi.org/10.1002/(sici)1099-1387(199711)3:6<442::Aid-psc122>3.0.Co;2-2)
- Synthetic study on selenocystine-containing peptides; T. Koide, H. Itoh, A. Otaka, H. Yasui, M. Kuroda, N. Esaki, K. Soda, N. Fujii; **Chem Pharm Bull (Tokyo)** 1993; **41**: 502-506. <https://doi.org/10.1248/cpb.41.502>
- Total chemical synthesis of enzymes; S. Kent; **J Pept Sci** 2003; **9**: 574-593. <https://doi.org/10.1002/psc.475>
- Über Peptidsynthesen. 8. Mitteilung Bildung von S-haltigen Peptiden durch intramolekulare Wanderung von Aminoacylresten; T. Wieland, E. Bokelmann, L. Bauer, H. U. Lang, H. Lau; **Justus Liebigs Annalen der Chemie** 1953; **583**: 129-149. <https://doi.org/10.1002/jlac.19535830110>

2. Tools for the Formation of Thioesters

2.1. Dawson Linkers

The straightforward C-terminal modification of peptides assembled on a solid support used to be a significant challenge in peptide and protein chemistry. This was also true for C-terminal thioester peptides, which are important intermediates in the generation of active esters, amides, and hydrazides. Moreover, they are also an essential component of many synthetic strategies for protein synthesis. The most efficient approach for the synthesis of peptidyl thioesters used to be the *in situ* neutralization protocol for Boc solid-phase peptide synthesis using thioester linkage.

Importantly, the linker is a stable amide during chain assembly and the key activation step utilizes the most robust reaction in solid phase peptide synthesis: the acylation of an amine. As a result, the method is compatible with amino acid side-chains and protecting groups commonly used in peptide synthesis.




References:

- Alternative chemistries for the synthesis of thrombospondin-1 type 1 repeats; T. K. Tiefenbrunn, J. Blanco-Canosa, P. E. Dawson; **Biopolymers** 2010; **94**: 405-413. <https://doi.org/10.1002/bip.21486>
- An Efficient Fmoc-SPPS Approach for the Generation of Thioester Peptide Precursors for Use in Native Chemical Ligation; J. B. Blanco-Canosa, P. E. Dawson; **Angew Chem** 2008; **120**: 6957-6961. <https://doi.org/10.1002/ange.200705471>
- Protein synthesis assisted by native chemical ligation at leucine; Z. Harpaz, P. Siman, K. S. Kumar, A. Brik; **ChemBiochem : a European journal of chemical biology** 2010; **11**: 1232-1235. <https://doi.org/10.1002/cbic.201000168>
- A reversible protection strategy to improve Fmoc-SPPS of peptide thioesters by the N-Acylurea approach; S. K. Mahto, C. J. Howard, J. C. Shimko, J. J. Ottesen; **ChemBiochem : a European journal of chemical biology** 2011; **12**: 2488-2494. <https://doi.org/10.1002/cbic.201100472>
- A semisynthesis platform for investigating structure-function relationships in the N-terminal domain of the anthrax Lethal Factor; B. L. Pentelute, A. P. Barker, B. E. Janowiak, S. B. Kent, R. J. Collier; **ACS Chem Biol** 2010; **5**: 359-364. <https://doi.org/10.1021/cb100003r>
- Synthesis of native proteins by chemical ligation; P. E. Dawson, S. B. Kent; **Annu. Rev. Biochem.** 2000; **69**: 923-60. <https://doi.org/10.1146/annurev.biochem.69.1.923>
- Synthesis of Proteins by Native Chemical Ligation; P. E. Dawson, T. W. Muir, I. Clark-Lewis, S. B. Kent; **Science** 1994; **266**(5186): 776-779. <https://doi.org/10.1126/science.7973629>
- Modulation of Reactivity in Native Chemical Ligation through the Use of Thiol Additives; P. E. Dawson, M. J. Chrchill, M. R. Ghadiri, S. B. Kent; **J. Am. Chem. Soc.** 1997; **119**(19): 4325-4329. <https://doi.org/10.1021/ja962656r>
- An Efficient Fmoc-SPPS Approach for the Generation of Thioester Peptide Precursors for Use in Native Chemical Ligation; J. B. Blanco-Canosa, P. E. Dawson; **Angew. Chem., Int. Ed.** 2008; **47**(36): 6851-6855. <https://doi.org/10.1002/anie.200705471>
- A Semisynthesis Platform for Investigating Structure-Function Relationships in the N-Terminal Domain of the Anthrax Lethal Factor; B. L. Pentelute, A. P. Barker, B. E. Janowiak, S. B. Kent, R. J. Collier; **ACS Chem. Biol.** 2010; **5**(4): 359-364. <https://doi.org/10.1021/cb100003r>
- A Reversible Protection Strategy To Improve Fmoc-SPPS of Peptide Thioesters by the N-Acylurea Approach; S. K. Mahto, C. J. Howard, J. C. Shimko, J. J. Ottesen; **ChemBioChem** 2011; **12**(16): 2488-2494. <https://doi.org/10.1002/cbic.201100472>
- Making Ends Meet: Microwave-Accelerated Synthesis of Cyclic and Disulfide Rich Proteins Via In Situ Thioesterification and Native Chemical Ligation; **Int. J. Pept. Res. Ther.** 2013; **19**: 43-54. <https://doi.org/10.1007/s10989-012-9331-y>
- Chemical Protein Synthesis Using a Second-Generation N-Acylurea Linker for the Preparation of Peptide-Thioester Precursors; J. B. Blanco-Canosa, B. Nardone, F. Albericio, P. E. Dawson; **J. Am. Chem. Soc.** 2015; **137**(22): 7197-7209. <https://doi.org/10.1021/jacs.5b03504>
- Native Chemical Ligation of Peptides and Proteins; P. A. Cistrone, M. J. Bird, D. T. Flood, A. P. Silvestri, J. C. J. Hintzen, D. A. Thompson, P. E. Dawson; **Curr. Protoc. Chem. Biol.** 2019; **11**(1): e61. <https://doi.org/10.1002/cpch.61>

FAA3165 Fmoc-Dbz-OH

3-[(9-Fluorenylmethoxycarbonyl)amino]-4-amino-benzoic acid






CAS-No. 1071446-05-3
Formula $C_{22}H_{18}N_2O_4$
Mol. weight 374,39 g/mol



Product details




[↑ back to content](#)


		Product details
FAA3166	Fmoc-MeDbz-OH 3-[(9-Fluorenylmethoxycarbonyl)amino]-4-(methylamino)benzoic acid CAS-No. 1788861-35-7 Formula $C_{23}H_{20}N_2O_4$ Mol. weight 388,44 g/mol	 
FAA3167	Fmoc-Dbz(o-Boc)-OH 3-(Fmoc-amino)-4-(Boc-amino)-benzoic acid CAS-No. 1823479-63-5 Formula $C_{27}H_{26}N_2O_6$ Mol. weight 474,51 g/mol	 
FAA3168	Fmoc-Dbz(o-Alloc)-OH 3-(Fmoc-amino)-4-(Alloc-amino)-benzoic acid CAS-No. 2143465-53-4 Formula $C_{26}H_{22}N_2O_6$ Mol. weight 458,46 g/mol	 
FAA3169	Fmoc-Dbz(o-NO₂)-OH 4-(Fmoc-amino)-3-Nitro-benzoic acid CAS-No. 1342864-48-5 Formula $C_{22}H_{16}N_2O_6$ Mol. weight 404,37 g/mol	 
S-30091	Fmoc-MeDbz-RAM TG S TentaGel S Dawson Fmoc-MeDbz RAM Resin Mesh Size 90 µm Loading 0.15 - 0.25 mmol/g	 

2.2. Hydrazone Resins


Peptide hydrazides can be easily synthesized using hydrazone resins, which are polystyrene resins functionalized with the hydrazone linker (Fig. 6). The hydrazone linker is completely stable in the course of standard Fmoc-SPPS. Moreover, it tolerates treatment with 5% TFA/DCM, thus permitting selective removal of Mtt or similar acid-labile protecting groups. Subsequent application of cleavage cocktails containing neat TFA permits to obtain the desired peptides in high yields and purity.



Synthesized peptide hydrazides can be applied as building blocks for the chemoselective conjugation to an aldehyde or ketone group on a second peptide (or other biomacromolecule) using the hydrazone ligation technique (Fig. 7). The resulting conjugate shows a non-natural hydrazone linkage replacing the natural amide bond. This, however, is in many cases tolerated, and the conjugate displays the same properties as a natural protein fragment.



A peptide hydrazide can be converted into a peptide thioester. Two different protocols (Fig. 8 and Fig. 9) have been published recently. These methods include, but are not limited to, solution phase oxidation of C-terminal hydrazides to the corresponding acylazides. The required oxidants – typically used in large excess – cannot be employed in concert with N-terminal thiazolidines, a common peptide protecting strategy utilized in the synthesis of proteins requiring multiple ligations, aryl amines, or other redox-sensitive residues which can be incorporated *via* SPPS.



A more robust method for the generation of Fmoc-SPPS peptide thioesters compatible with multiple ligations utilizes mild and selective conditions at stoichiometric quantities of reactants. This process is suitable for multiple ligations to be performed sequentially in one pot, hence, has the potential for a general route to synthesize larger and complex proteins via multiple chemical ligation steps.










Fig. 9: General route for synthesizing peptide thioesters utilizing hydrazone resin.

Peptide hydrazides exposed to acetylacetone (acac, 2.5 eq.) in a solution of 6 M guanidinium chloride and 2% thiophenol at acidic conditions will initially form a hydrazone. This cyclizes to pyrazole at low pH and can be replaced with thiols like MPAA resulting in the corresponding peptide thioester.

References:


- Chemical synthesis of proteins using peptide hydrazides as thioester surrogates; J. S. Zheng, S. Tang, Y. K. Qi, Z. P. Wang, L. Liu; **Nat Protoc** 2013; **8**: 2483-2495. <https://doi.org/10.1038/nprot.2013.152>
- Convenient method of peptide hydrazide synthesis using a new hydrazone resin; P. S. Chelushkin, K. V. Polyanichko, M. V. Leko, M. Y. Dorosh, T. Bruckdorfer, S. V. Burov; **Tetrahedron Letters** 2015; **56**: 619-622. <https://doi.org/10.1016/j.tetlet.2014.12.056>
- Leveraging the Knorr Pyrazole Synthesis for the Facile Generation of Thioester Surrogates for use in Native Chemical Ligation; D. T. Flood, J. C. J. Hintzen, M. J. Bird, P. A. Cistrone, J. S. Chen, P. E. Dawson; **Angew Chem Int Ed Engl** 2018; **57**: 11634-11639. <https://doi.org/10.1002/anie.201805191>
- Protein chemical synthesis by ligation of peptide hydrazides; G. M. Fang, Y. M. Li, F. Shen, Y. C. Huang, J. B. Li, Y. Lin, H. K. Cui, L. Liu; **Angew Chem Int Ed Engl** 2011; **50**: 7645-7649. <https://doi.org/10.1002/anie.201100996>
- Selective Coupling at the α -Amino Group of Cysteine Using Transfer Active-ester-condensation Technology to Synthesize a Linear Octadecapeptide; Y. Liao, Y. Kong, N. Hu, Z. Jin, P. Wang; **Chem Lett** 2010; **39**: 196-197. <https://doi.org/10.1246/cl.2010.196>
- A shelf stable Fmoc hydrazine resin for the synthesis of peptide hydrazides; M. J. Bird, P. E. Dawson; **Pept. Sci.** 2022; **144**: e24268. <https://doi.org/10.1002/pep2.24268>
- Selective Peptide Cysteine Manipulation on Demand and Difficult Protein Chemical Synthesis Enabled by Controllable Acidolysis of N,S-Benzylidene Thioacetals; H. Wu, Z. Sun, X. Li; **Angew. Chem. Int. Ed.** 2024; **63**: e202403396. <https://doi.org/10.1002/anie.202403396>

Product details		
PYV1000 Fmoc-NHN=Pyv Resin Fmoc-hydrazono-pyruvyl-aminomethylpolystyrene resin Mesh Size 100-200 mesh Loading > 0.3 mmol/g DVB 1% DVB		
BR-5279 Fmoc-NHNH-2CT Resin Fmoc-hydrazine-2-chlorotrityl resin Mesh Size 100-200 mesh Loading 0,4-1,4 mmol/g DVB 1% DVB		
BR-5280 Fmoc-NHNH-Trt-PEG Resin Fmoc-hydrazine-trityl polyethyleneglycol resin Loading 0,10-0,30 mmol/g DVB 1% DVB		
PEG8260 Fmoc-NHNH-Trt-PEG Resin HL Fmoc-hydrazine-trityl polyethyleneglycol resin		
PEG8265 Fmoc-NHNH-Trt-PEG Resin XV Fmoc-hydrazine-trityl polyethyleneglycol resin		
PYV1100 Fmoc-L-Ala-NHN=Pyv Resin Fmoc-L-alanyl-hydrazono-pyruvyl-aminomethylpolystyrene resin Mesh Size 100-200 mesh Loading > 0.3 mmol/g DVB 1% DVB		

PYV1110 Fmoc-L-Arg(Pbf)-NHN=Pyv Resin

Fmoc-N'-2,2,4,6,7-pentamethyldihydrobenzo-furan-5-sulfonyl-L-arginyl-hydrazono-pyruvyl-amino-methylpolystyrene resin


Mesh Size 100-200 mesh
Loading > 0.3 mmol/g
DVB 1% DVB



PYV1120 Fmoc-L-Asn(Trt)-NHN=Pyv Resin

Fmoc-N-beta-trityl-L-asparaginyl-hydrazono-pyruvyl-aminomethylpolystyrene resin


Mesh Size 100-200 mesh
Loading > 0.3 mmol/g
DVB 1% DVB



PYV1130 Fmoc-L-Asp(OtBu)-NHN=Pyv Resin

Fmoc-L-aspartyl-beta-t-butyl ester-alpha-hydrazono-pyruvyl-aminomethylpolystyrene resin


Mesh Size 100-200 mesh
Loading > 0.3 mmol/g
DVB 1% DVB



PYV1140 Fmoc-L-Cys(Trt)-NHN=Pyv Resin

Fmoc-S-trityl-L-cysteiny-l-hydrazono-pyruvyl-aminomethylpolystyrene resin


Mesh Size 100-200 mesh
Loading > 0.3 mmol/g
DVB 1% DVB



PYV1150 Fmoc-L-Glu(tBu)-NHN=Pyv Resin

Fmoc-L-glutamyl-gamma-t-butyl ester-alpha-hydrazono-pyruvyl-aminomethylpolystyrene resin


Mesh Size 100-200 mesh
Loading > 0.3 mmol/g
DVB 1% DVB















PYV1160 Fmoc-L-Gln(Trt)-NHN=Pyv Resin

Fmoc-N-gamma-trityl-L-glutaminy-l-hydrazono-pyruvyl-aminomethylpolystyrene resin


Mesh Size 100-200 mesh
Loading > 0.3 mmol/g
DVB 1% DVB




			Product details
PYV1170	Fmoc-Gly-NHN=Pyv Resin		
Fmoc-glycyl-hydrazono-pyruvyl-aminomethylpolystyrene resin			
Mesh Size	100-200 mesh		
Loading	> 0.3 mmol/g		
DVB	1% DVB		
PYV1180	Fmoc-L-His(Trt)-NHN=Pyv Resin		
Fmoc-N-trityl-L-histidyl-hydrazono-pyruvyl-aminomethylpolystyrene resin			
Mesh Size	100-200 mesh		
Loading	> 0.3 mmol/g		
DVB	1% DVB		
PYV1190	Fmoc-L-Ile-NHN=Pyv Resin		
Fmoc-L-isoleucyl-hydrazono-pyruvyl-aminomethylpolystyrene resin			
Mesh Size	100-200 mesh		
Loading	> 0.3 mmol/g		
DVB	1% DVB		
PYV1200	Fmoc-L-Leu-NHN=Pyv Resin		
Fmoc-L-leucyl-hydrazono-pyruvyl-aminomethylpolystyrene resin			
Mesh Size	100-200 mesh		
Loading	> 0.3 mmol/g		
DVB	1% DVB		
PYV1210	Fmoc-L-Lys(Boc)-NHN=Pyv Resin		
Fmoc-N-epsilon-t-butyloxycarbonyl-L-lysyl-hydrazono-pyruvyl-aminomethylpolystyrene resin			
Mesh Size	100-200 mesh		
Loading	> 0.3 mmol/g		
DVB	1% DVB		
PYV1220	Fmoc-L-Met-NHN=Pyv Resin		
Fmoc-L-methionyl-hydrazono-pyruvyl-aminomethylpolystyrene resin			
Mesh Size	100-200 mesh		
Loading	> 0.3 mmol/g		
DVB	1% DVB		

		Product details
PYV1230 Fmoc-L-Phe-NHN=Pyv Resin Fmoc-L-phenylalanyl-hydrazono-pyruvyl-aminomethylpolystyrene resin Mesh Size 100-200 mesh Loading > 0.3 mmol/g DVB 1% DVB		
PYV1240 Fmoc-L-Pro-NHN=Pyv Resin Fmoc-L-prolinyl-hydrazono-pyruvyl-aminomethylpolystyrene resin Mesh Size 100-200 mesh Loading > 0.3 mmol/g DVB 1% DVB		
PYV1250 Fmoc-L-Ser(tBu)-NHN=Pyv Resin Fmoc-O-t-butyl-L-seryl-hydrazono-pyruvyl-aminomethylpolystyrene resin Mesh Size 100-200 mesh Loading > 0.3 mmol/g DVB 1% DVB		
PYV1260 Fmoc-L-Thr(tBu)-NHN=Pyv Resin Fmoc-O-t-butyl-L-threonyl-hydrazono-pyruvyl-aminomethylpolystyrene resin Mesh Size 100-200 mesh Loading > 0.3 mmol/g DVB 1% DVB		
PYV1270 Fmoc-L-Trp(Boc)-NHN=Pyv Resin Fmoc-N-t-butyloxycarbonyl-L-tryptophyl-hydrazono-pyruvyl-aminomethylpolystyrene resin Mesh Size 100-200 mesh Loading > 0.3 mmol/g DVB 1% DVB		
PYV1280 Fmoc-L-Tyr(tBu)-NHN=Pyv Resin Fmoc-O-t-butyl-L-tyrosyl-hydrazono-pyruvyl-aminomethylpolystyrene resin Mesh Size 100-200 mesh Loading > 0.3 mmol/g DVB 1% DVB		


		Product details
PYV1290	Fmoc-L-Val-NHN=Pyv Resin	
Fmoc-L-valyl-hydrazono-pyruvyl-aminomethylpolystyrene resin		
Mesh Size	100-200 mesh	
Loading	> 0.3 mmol/g	
DVB	1% DVB	





Interested in our resins?

Download our booklet!



Description of the Technology

Tools for the Formation of Thioesters

Building Blocks for Ligation

Enzyme-mediated Ligation

Examples





Index







[↑ back to content](#)

3. Building Blocks for Ligation

3.1. Cysteine Building Blocks for Ligation

1,3-Thiazolidine-4-carboxylic acid (Thz) is a protected cysteine. It can be converted into free cysteine by reaction with aqueous 0.8 M methoxylamine (MeONH₂).

		Product details
HAA1132 H-L-Thz-OH (R)-Thiazolidine-4-carboxylic acid CAS-No. 34592-47-7 Formula C ₄ H ₇ NO ₂ S Mol. weight 133,16 g/mol		
HAA3840 H-L-MeAla(4-Thz)-OH N-alpha-Methyl-beta-(4-thiazolyl)-L-alanine CAS-No. 2131118-50-6 Formula C ₇ H ₁₀ N ₂ O ₂ S Mol. weight 186,23 g/mol		
BAA1135 Boc-L-Thz-OH (R)-N-t-Butyloxycarbonyl-thiazolidine-4-carboxylic acid CAS-No. 51077-16-8 Formula C ₉ H ₁₅ NO ₄ S Mol. weight 233,29 g/mol		
BAA1186 Boc-D-Thz-OH (S)-N-(t-Butyloxycarbonyl)-thiazolidine-4-carboxylic acid CAS-No. 63091-82-7 Formula C ₉ H ₁₅ NO ₄ S Mol. weight 233,29 g/mol		

		Product details	
FAA1427	Fmoc-L-Thz-OH		
(R)-N-(9-Fluorenylmethyloxycarbonyl)-thiazolidine-4-carboxylic acid			
CAS-No.	133054-21-4		
Formula	C ₁₉ H ₁₇ NO ₄ S		
Mol. weight	355,42 g/mol		
FAA1495	Fmoc-D-Thz-OH		
(S)-N-α-(9-Fluorenylmethyloxycarbonyl)-thiazolidine-4-carboxylic acid			
CAS-No.	198545-89-0		
Formula	C ₁₉ H ₁₇ NO ₄ S		
Mol. weight	355,42 g/mol		
FAA1437	Fmoc-L-Thz(Me₂)-OH		
(R)-N-(9-Fluorenylmethyloxycarbonyl)-2,2-dimethyl-thiazolidine-4-carboxylic acid			
CAS-No.	873842-06-9		
Formula	C ₂₁ H ₂₁ NO ₄ S		
Mol. weight	383,46 g/mol		
FAA9200	Fmoc-L-Thz(Dmp)-OH		
(4R)-3-(((9H-fluoren-9-yl)methoxy)carbonyl)-2-(2,4-dimethoxyphenyl)thiazolidine-4-carboxylic acid			
CAS-No.	2648642-22-0		
Formula	C ₂₇ H ₂₅ NO ₆ S		
Mol. weight	491,56 g/mol		
FAA9320	Fmoc-L-Lys(Boc-Thz)-OH		
N2-(((9H-fluoren-9-yl)methoxy)carbonyl)-N6-((R)-3-(tert-butoxycarbonyl)thiazolidine-4-carboxyl)-L-lysine			
Formula	C ₃₀ H ₃₇ N ₃ O ₇ S		
Mol. weight	583,70 g/mol		
FAA9330	Fmoc-L-Dab(Boc-Thz)-OH		
(S)-2-(((9H-fluoren-9-yl)methoxy)carbonyl)amino)-4-((R)-3-(tert-butoxycarbonyl)thiazolidine-4-carboxamido)butanoic acid			
CAS-No.	2968514-54-5		
Formula	C ₂₈ H ₃₃ N ₃ O ₇ S		
Mol. weight	555,65 g/mol		

Description of the Technology

Tools for the Formation of Thioesters







Building Blocks for Ligation







Enzyme-mediated Ligation

Examples

Index

[↑ back to content](#)

			Product details
FAA9340	Fmoc-L-Lys(4-Thz, Boc)-OH		
(2S)-2-((((9H-fluoren-9-yl)methoxy)carbonyl)amino)-4-(3-(tert-butoxycarbonyl)thiazolidin-5-yl)butanoic acid CAS-No. 1240666-28-7 Formula $C_{27}H_{32}N_2O_6S$ Mol. weight 512,62 g/mol			
HAA1017	H-D-Cys-OH*HCl*H₂O		
D-Cysteine hydrochloride monohydrate CAS-No. 32443-99-5 Formula $C_3H_7NO_2S \cdot HCl \cdot H_2O$ Mol. weight 121,2*36,45*18,01 g/mol			
HAA3530	H-L-Cys-NH₂*HCl		
L-Cysteine amide hydrochloride CAS-No. 16359-98-1 Formula $C_3H_8N_2OS \cdot HCl$ Mol. weight 120,17*36,45 g/mol			
HAA1574	H-L-Cys(Bzl)-OH		
S-Benzyl-L-cysteine CAS-No. 3054-01-1 Formula $C_{10}H_{13}NO_2S$ Mol. weight 211,29 g/mol			
HAA6110	H-D-Cys(Bzl)-OH		
S-Benzyl-D-cysteine CAS-No. 23032-53-3 Formula $C_{10}H_{13}NO_2S$ Mol. weight 211,29 g/mol			
HAA6070	H-L-Cys(Acm)-OH*HCl		
S-(Acetyl-aminomethyl)-L-cysteine hydrochloride CAS-No. 28798-28-9 Formula $C_6H_{12}N_2O_3S \cdot HCl$ Mol. weight 192,24*36,45 g/mol			

		Product details
HAA6080	H-L-Cys(Bzl)-OMe*HCl S-Benzyl-L-cysteine methyl ester hydrochloride CAS-No. 16741-80-3 Formula $C_{11}H_{15}NO_2S \cdot HCl$ Mol. weight 225,31*36,45 g/mol	 
HAA6090	H-L-Cys(MBzl)-OH S-(4-Methylbenzyl)-L-cysteine CAS-No. 42294-52-0 Formula $C_{11}H_{15}NO_2S$ Mol. weight 225,3 g/mol	 
HAA9270	H-L-Cys(MDNPE)-OH 1-[4',5'-(methylenedioxy)-2'-nitrophenyl]ethyl]-L-cysteine CAS-No. 1551078-43-3 Formula $C_{12}H_{14}N_2O_6S$ Mol. weight 314,31 g/mol	 
HAA3500	H-D-Cys(Mmt)-OH S-p-methoxytrityl-D-cysteine CAS-No. 926935-33-3 Formula $C_{23}H_{23}NO_3S$ Mol. weight 393,5 g/mol	 
HAA6100	H-L-Cys(Mob)-OH S-(4-Methoxybenzyl)-L-cysteine CAS-No. 2544-31-2 Formula $C_{11}H_{15}NO_3S$ Mol. weight 241,3 g/mol	 
HAA3510	H-L-Cys(Npys)-OH*HCl S-(3-nitro-2-pyridylthio)-L-cysteine hydrochloride CAS-No. 108807-66-5 Formula $C_8H_9N_3O_4S_2 \cdot HCl$ Mol. weight 275,30*36,45 g/mol	 

Description of the Technology

Tools for the Formation of Thioesters







Building Blocks for Ligation







Enzyme-mediated Ligation

Examples

Index

[↑ back to content](#)

			Product details
HAA6160	H-L-Cys(Trt)-OH		
S-Trityl-L-cysteine			
CAS-No.	2799-07-7		
Formula	$C_{22}H_{21}NO_2S$		
Mol. weight	363,48 g/mol		
			
HAA3520	H-D-Cys(Trt)-OMe*HCl		
S-trityl-D-cysteine methyl ester hydrochloride			
CAS-No.	1020369-32-7		
Formula	$C_{23}H_{23}NO_2S^*HCl$		
Mol. weight	377,50*36,45 g/mol		
			
HAA1995	H-L-Cys(Trt)-OtBu*HCl		
S-Trityl-L-cysteine t-butyl ester hydrochloride			
CAS-No.	158009-03-1		
Formula	$C_{26}H_{29}NO_2S^*HCl$		
Mol. weight	419,58*36,45 g/mol		
			
HAA1560	H-L-Cys(Trt)-NH2		
S-Trityl-L-cysteine amide			
CAS-No.	166737-85-5		
Formula	$C_{22}H_{22}N_2OS$		
Mol. weight	362,49 g/mol		
			
HAA6120	H-D-Cys(Trt)-OH		
S-Trityl-D-cysteine			
CAS-No.	25840-82-8		
Formula	$C_{22}H_{21}NO_2S$		
Mol. weight	363,48 g/mol		
			
HAA2100	H-D-Cys(Trt)-OtBu*HCl		
S-Trityl-D-cysteine t-butyl ester hydrochloride			
CAS-No.	439089-10-8		
Formula	$C_{26}H_{29}NO_2S^*HCl$		
Mol. weight	419,58*36,45 g/mol		
			

		Product details
HAA6140	H-L-Cys(StBu)-OH	
S-Thio-t-butyl-L-cysteine		
CAS-No.	30044-51-0	
Formula	C ₇ H ₁₅ NO ₂ S ₂	
Mol. weight	209,32 g/mol	
		
HAA6150	H-L-Cys(tBu)-OH*HCl	
S-t-Butyl-L-cysteine hydrochloride		
CAS-No.	2481-09-6	
Formula	C ₇ H ₁₅ NO ₂ S*HCl	
Mol. weight	177,26*36,45 g/mol	
		
HAA9320	H-L-Cys(oNv)-OH	
S-(4,5-dimethoxy-2-nitrobenzyl)-L-cysteine		
CAS-No.	214633-68-8	
Formula	C ₁₂ H ₁₆ N ₂ O ₆ S	
Mol. weight	316,33 g/mol	
		
AAA1300	Ac-L-Cys-OH	
N-alpha-Acetyl-L-cysteine		
CAS-No.	616-91-1	
Formula	C ₅ H ₉ NO ₃ S	
Mol. weight	163,19 g/mol	
		
AAA2015	Aloc-L-Cys(Trt)-OH	
N-alpha-Allyloxycarbonyl-S-trityl-L-cysteine		
CAS-No.	96865-72-4	
Formula	C ₂₆ H ₂₅ NO ₄ S	
Mol. weight	447,55 g/mol	
		
HAA2810	N₃-L-Cys(Trt)-OH*CHA	
(R)-2-azido-3-(tritylthio)propanoic acid cyclohexylamine		
CAS-No.	1286670-90-3	
Formula	C ₂₂ H ₁₉ N ₃ O ₂ S*C ₆ H ₁₃ N	
Mol. weight	389,47*99,17 g/mol	
		

Description of the Technology

Tools for the Formation of Thioesters







Building Blocks for Ligation







Enzyme-mediated Ligation

Examples

Index

[↑ back to content](#)

		Product details
BAA1083	Boc-L-Cys-OH	
N-alpha-t-Butyloxycarbonyl-L-cysteine		
CAS-No.	20887-95-0	 
Formula	C ₈ H ₁₅ NO ₄ S	
Mol. weight	221,27 g/mol	
BAA1170	Boc-D-Cys-OH	
N-alpha-t-Butyloxycarbonyl-D-cysteine		
CAS-No.	149270-12-2	 
Formula	C ₈ H ₁₅ NO ₄ S	
Mol. weight	221,27 g/mol	
BAA1078	Boc-L-Cys(Acm)-OH	
N-alpha-t-Butyloxycarbonyl-S-(acetyl-amino-methyl)-L-cysteine		
CAS-No.	19746-37-3	 
Formula	C ₁₁ H ₂₀ N ₂ O ₅ S	
Mol. weight	292,36 g/mol	
BAA1510	Boc-L-Cys(Acm,O)-OH	
N-alpha-t-Butyloxycarbonyl-S-(acetyl-amino-methyl)-S-oxo-L-cysteine		
CAS-No.	75893-04-8	 
Formula	C ₁₁ H ₂₀ N ₂ O ₆ S	
Mol. weight	308,35 g/mol	
BAA1079	Boc-L-Cys(Bzl)-OH	
N-alpha-t-Butyloxycarbonyl-S-benzyl-L-cysteine		
CAS-No.	5068-28-0	 
Formula	C ₁₅ H ₂₁ NO ₄ S	
Mol. weight	311,38 g/mol	
BAA5410	Boc-D-Cys(Bzl)-OH	
N-alpha-t-Butyloxycarbonyl-S-benzyl-D-cysteine		
CAS-No.	102830-49-9	 
Formula	C ₁₅ H ₂₁ NO ₄ S	
Mol. weight	311,38 g/mol	

		Product details
BAA5510	Boc-L-Cys(Fm)-OH	
N-alpha- <i>t</i> -Butyloxycarbonyl-S-(9-fluorenyl-methyl)-L-cysteine		
CAS-No.	84888-35-7	
Formula	C ₂₂ H ₂₅ NO ₄ S	
Mol. weight	399,51 g/mol	
		
BAA1080	Boc-L-Cys(MBzl)-OH	
N-alpha- <i>t</i> -Butyloxycarbonyl-S-(4-methyl-benzyl)-L-cysteine		
CAS-No.	61925-77-7	
Formula	C ₁₆ H ₂₃ NO ₄ S	
Mol. weight	325,43 g/mol	
		
BAA5420	Boc-D-Cys(MBzl)-OH	
N-alpha- <i>t</i> -Butyloxycarbonyl-S-(4-methyl-benzyl)-D-cysteine		
CAS-No.	61925-78-8	
Formula	C ₁₆ H ₂₃ NO ₄ S	
Mol. weight	325,43 g/mol	
		
BAA1081	Boc-L-Cys(Mob)-OH	
N-alpha- <i>t</i> -Butyloxycarbonyl-S-(4-methoxy-benzyl)-L-cysteine		
CAS-No.	18942-46-6	
Formula	C ₁₆ H ₂₃ NO ₅ S	
Mol. weight	341,43 g/mol	
		
BAA1860	Boc-L-Cys(Npys)-OH	
N-alpha- <i>t</i> -Butyloxycarbonyl-S-(3-nitro-2-pyridylthio)-L-cysteine		
CAS-No.	76880-29-0	
Formula	C ₁₃ H ₁₇ N ₃ O ₆ S ₂	
Mol. weight	375,42 g/mol	
		
BAA6390	Boc-L-Cys(Phacm)-OH	
N-alpha- <i>t</i> -Butyloxycarbonyl-S-(Phenylacetyl-amino-methyl)-L-cysteine		
CAS-No.	57084-73-8	
Formula	C ₁₇ H ₂₄ N ₂ O ₅ S	
Mol. weight	368,45 g/mol	
		

Description of the Technology

Tools for the Formation of Thioesters







Building Blocks for Ligation







Enzyme-mediated Ligation

Examples

Index

[↑ back to content](#)

		Product details
BAA3140	Boc-L-Cys(Ph)-OH	
N-alpha-t-Butyloxycarbonyl-S-phenyl-L-cysteine		
CAS-No.	163705-28-0	
Formula	C ₁₆ H ₁₉ NO ₄ S	
Mol. weight	297,37 g/mol	
BAA1082		
Boc-L-Cys(tBu)-OH		
N-alpha-t-Butyloxycarbonyl-S-t-butyl-L-cysteine		
CAS-No.	56976-06-8	
Formula	C ₁₇ H ₂₃ NO ₄ S	
Mol. weight	277,37 g/mol	
BAA6415		
Boc-L-Cys(StBu)-OH		
N-(tert-butoxycarbonyl)-S-(tert-butylthio)-L-cysteine		
CAS-No.	30044-61-2	
Formula	C ₁₇ H ₂₃ NO ₄ S ₂	
Mol. weight	309,44 g/mol	
BAA1084		
Boc-L-Cys(Trt)-OH		
N-alpha-t-Butyloxycarbonyl-S-trityl-L-cysteine		
CAS-No.	21947-98-8	
Formula	C ₂₇ H ₂₉ NO ₄ S	
Mol. weight	463,59 g/mol	
BAA5000		
Boc-D-Cys(Trt)-OH		
N-alpha-t-Butyloxycarbonyl-S-trityl-D-cysteine		
CAS-No.	87494-13-1	
Formula	C ₂₇ H ₂₉ NO ₄ S	
Mol. weight	463,59 g/mol	
IAD1040		
Boc-L-Ser[Fmoc-L-Cys(Trt)]-OH		
O-(N-(((9H-fluoren-9-yl)methoxy)carbonyl)-S-trityl-L-cysteinyl)-N-(tert-butoxycarbonyl)-L-serine		
CAS-No.	944283-13-0	
Formula	C ₄₅ H ₄₄ N ₂ O ₈ S	
Mol. weight	772,9	
		

		Product details
IAD2040	Boc-L-Thr[Fmoc-L-Cys(Trt)]-OH O-(N-(((9H-fluoren-9-yl)methoxy)carbonyl)-S-trityl-L-cysteiny)-N-(tert-butoxycarbonyl)-L-threonine CAS-No. 944283-30-1 Formula $C_{46}H_{46}N_2O_8S$ Mol. weight 786,93 g/mol	 
FAA1362	Fmoc-L-Cys-OH*H₂O N-alpha-(9-Fluorenylmethoxycarbonyl)-L-cysteine monohydrat CAS-No. 135248-89-4 Formula $C_{18}H_{17}NO_4S \cdot H_2O$ Mol. weight 343,40*18,01 g/mol	 
FAA1470	Fmoc-D-Cys-OH*H₂O N-alpha-(9-Fluorenylmethoxycarbonyl)-D-cysteine monohydrat CAS-No. 157355-80-1 Formula $C_{18}H_{17}NO_4S \cdot H_2O$ Mol. weight 343,4*18,01 g/mol	 
FAA1980	Fmoc-L-Cys-NH₂ N-alpha-(9-Fluorenylmethoxycarbonyl)-L-cysteine amide CAS-No. 623177-62-8 Formula $C_{18}H_{18}N_2O_3S$ Mol. weight 342,41 g/mol	 
FAA5150	Fmoc-L-Cys(Aapam)-OH N-alpha-(9-Fluorenylmethoxycarbonyl)-S-((4-(allyloxycarbonylamino)phenylacetylaminomethyl)-L-cysteine CAS-No. 1946783-89-6 Formula $C_{31}H_{31}N_3O_7S$ Mol. weight 589,66 g/mol	 
FAA1506	Fmoc-L-Cys(Acm)-OH N-alpha-(9-Fluorenylmethoxycarbonyl)-S-(acetylaminomethyl)-L-cysteine CAS-No. 86060-81-3 Formula $C_{21}H_{22}N_2O_5S$ Mol. weight 414,48 g/mol	 

Description of the Technology

Tools for the Formation of Thioesters







Building Blocks for Ligation







Enzyme-mediated Ligation

Examples

Index

[↑ back to content](#)

			Product details
FAA6230	Fmoc-D-Cys(Acm)-OH		
N-alpha-(9-Fluorenylmethyloxycarbonyl)-S-(acetyl-aminomethyl)-D-cysteine			
CAS-No.	168300-88-7		
Formula	C ₂₁ H ₂₂ N ₂ O ₅ S		
Mol. weight	414,48 g/mol		
			
FAA6720	Fmoc-L-MeCys(Acm)-OH		
N-alpha-(9-Fluorenylmethyloxycarbonyl)-N-alpha-methyl-S-(acetyl-aminomethyl)-L-cysteine			
CAS-No.	481642-19-7		
Formula	C ₂₂ H ₂₄ N ₂ O ₅ S		
Mol. weight	428,5 /mol		
			
FAA4751	Fmoc-L-Cys(Ac-OtBu)-OH*DCHA		
N-alpha-(9-Fluorenylmethyloxycarbonyl)-S-(t-butoxycarbonylmethyl)-L-cysteine dicyclohexylamine			
CAS-No.	269730-62-3 net		
Formula	C ₂₄ H ₂₇ NO ₆ S*C ₁₂ H ₂₃ N		
Mol. weight	457,54*181,32 g/mol		
			
FAA7610	Fmoc-L-Cys(Allocam)-OH		
N-alpha-(9-Fluorenylmethyloxycarbonyl)-S-((allyloxycarbonylamino)methyl)-L-cysteine			
CAS-No.	232953-09-2		
Formula	C ₂₃ H ₂₄ N ₂ O ₆ S		
Mol. weight	456,51 g/mol		
			
FAA6270	Fmoc-L-Cys(Bzl)-OH		
N-alpha-(9-Fluorenylmethyloxycarbonyl)-S-benzyl-L-cysteine			
CAS-No.	53298-33-2		
Formula	C ₂₅ H ₂₃ NO ₄ S		
Mol. weight	433,52 g/mol		
			
FAA6940	Fmoc-L-Cys(Ddm)-OH		
N-alpha-(9-Fluorenylmethyloxycarbonyl)-S-di(4-methoxyphenyl)methyl-L-cysteine			
CAS-No.	1403825-56-8		
Formula	C ₃₃ H ₃₁ NO ₆ S		
Mol. weight	569,67 g/mol		
			

		Product details
FAA3180	Fmoc-L-Cys(S-DMP)-OH	
N-alpha-(9-Fluorenylmethyloxycarbonyl)-S-(2,6-dimethoxythiophenol)-L-cysteine		
CAS-No.	1403834-73-0	
Formula	C ₂₆ H ₂₅ NO ₆ S ₂	
Mol. weight	511,61 g/mol	
		
FAA3190	Fmoc-L-Cys(Dpm)-OH	
N-alpha-(9-Fluorenylmethyloxycarbonyl)-S-diphenylmethyl-L-cysteine		
CAS-No.	247595-29-5	
Formula	C ₃₁ H ₂₇ NO ₄ S	
Mol. weight	509,62 g/mol	
		
FAA5650	Fmoc-D-Cys(Dpm)-OH	
N-alpha-(9-Fluorenylmethyloxycarbonyl)-S-diphenylmethyl-D-cysteine		
CAS-No.	2389078-16-2	
Formula	C ₃₁ H ₂₇ NO ₄ S	
Mol. weight	509,62 g/mol	
		
FAA7945	Fmoc-L-Cys(MDNPE)-OH	
N-(((9H-fluoren-9-yl)methoxy)carbonyl)-S-(1-(6-nitrobenzo[d][1,3]dioxol-5-yl)ethyl)-L-cysteine		
Formula	C ₂₇ H ₂₄ N ₂ O ₈ S	
Mol. weight	536,56 g/mol	
		
FAA1714	Fmoc-L-Cys(MBzl)-OH	
N-alpha-(9-Fluorenylmethyloxycarbonyl)-S-(4-methylbenzyl)-L-cysteine		
CAS-No.	136050-67-4	
Formula	C ₂₆ H ₂₅ NO ₄ S	
Mol. weight	447,53 g/mol	
		
FAA1030	Fmoc-L-Cys(Mmt)-OH	
N-alpha-(9-Fluorenylmethyloxycarbonyl)-S-p-methoxytrityl-L-cysteine		
CAS-No.	177582-21-7	
Formula	C ₃₈ H ₃₃ NO ₅ S	
Mol. weight	615,74 g/mol	
		

Description of the Technology

Tools for the Formation of Thioesters







Building Blocks for Ligation







Enzyme-mediated Ligation

Examples

Index

[↑ back to content](#)

			Product details
FAA1614	Fmoc-D-Cys(Mmt)-OH		
N-alpha-(9-Fluorenylmethyloxycarbonyl)-S-p-methoxytrityl-D-cysteine			
CAS-No.	1198791-73-9		
Formula	C ₃₈ H ₃₃ NO ₃ S		
Mol. weight	615,74 g/mol		
			
FAA4845	Fmoc-alpha-Me-L-Cys(Mmt)-OH		
N-alpha-(9-Fluorenylmethyloxycarbonyl)-alpha-methyl-S-(4-methoxytrityl)-L-cysteine			
CAS-No.	1198791-74-0		
Formula	C ₃₉ H ₃₅ NO ₃ S		
Mol. weight	629,76 g/mol		
			
FAA1715	Fmoc-L-Cys(Mob)-OH		
N-alpha-(9-Fluorenylmethyloxycarbonyl)-S-(4-methoxybenzyl)-L-cysteine			
CAS-No.	141892-41-3		
Formula	C ₂₆ H ₂₅ NO ₃ S		
Mol. weight	463,55 g/mol		
			
FAA4155	Fmoc-L-Cys(Msbh)-OH		
N-alpha-(9-Fluorenylmethyloxycarbonyl)-S-(4,4'-dimethylsulfinylbenzhydryl)-L-cysteine			
CAS-No.	1584646-97-8		
Formula	C ₃₃ H ₃₁ NO ₆ S ₃		
Mol. weight	633,80 g/mol		
			
FAA8150	Fmoc-D-Cys(Msbh)-OH		
N-alpha-(9-Fluorenylmethyloxycarbonyl)-S-(4,4'-dimethylsulfinylbenzhydryl)-D-cysteine			
Formula	C ₃₃ H ₃₁ NO ₆ S ₃		
Mol. weight	633,80 g/mol		
			
FAA8420	Fmoc-L-Cys(NDBF)-OH		
N-alpha-(9-Fluorenylmethyloxycarbonyl)-S-(1-(3-nitro-dibenzofuran-2-yl)-ethyl)-L-cysteine			
CAS-No.	1895883-28-9		
Formula	C ₃₂ H ₂₆ N ₂ O ₇ S		
Mol. weight	582,62 g/mol		
			

		Product details
FAA3970	Fmoc-L-Cys(oNv)-OH	
N-alpha-(9-Fluorenylmethyloxycarbonyl)-S-(2-nitroveratryl)-L-cysteine		
CAS-No.	214633-71-3	
Formula	C ₂₇ H ₂₆ N ₂ O ₆ S	
Mol. weight	538,57 g/mol	
		
FAA6910	Fmoc-L-Cys(Phacm)-OH	
N-alpha-(9-Fluorenylmethyloxycarbonyl)-S-Phenylacetylaminomethyl)-L-cysteine		
CAS-No.	159680-21-4	
Formula	C ₂₇ H ₂₆ N ₂ O ₅ S	
Mol. weight	490,57 g/mol	
		
FAA3710	Fmoc-D-Cys(Phacm)-OH	
N-alpha-(9-Fluorenylmethyloxycarbonyl)-S-(phenylacetylaminomethyl)-D-cysteine		
CAS-No.	1565818-55-4	
Formula	C ₂₇ H ₂₆ N ₂ O ₅ S	
Mol. weight	490,57 g/mol	
		
FAA8495	Fmoc-L-Cys(SIT)-OH	
N-alpha-(9-Fluorenylmethyloxycarbonyl)-S-(sec-isoamyl mercaptan)-L-cysteine		
CAS-No.	2545642-31-5	
Formula	C ₂₃ H ₂₇ NO ₄ S ₂	
Mol. weight	445,59 g/mol	
		
FAA1575	Fmoc-L-Cys(StBu)-OH	
N-alpha-(9-Fluorenylmethyloxycarbonyl)-S-(t-butylthio)-L-cysteine		
CAS-No.	73724-43-3	
Formula	C ₂₂ H ₂₅ NO ₄ S ₂	
Mol. weight	431,57 g/mol	
		
FAA1965	Fmoc-D-Cys(StBu)-OH	
N-alpha-(9-Fluorenylmethyloxycarbonyl)-S-(t-butylthio)-D-cysteine		
CAS-No.	501326-55-2	
Formula	C ₂₂ H ₂₅ NO ₄ S ₂	
Mol. weight	431,57 g/mol	
		

Description of the Technology

Tools for the Formation of Thioesters







Building Blocks for Ligation







Enzyme-mediated Ligation

Examples

Index

[↑ back to content](#)

		Product details
FAA3340	Fmoc-L-MeCys(S-tBu)-OH N-alpha-(9-Fluorenylmethyloxycarbonyl)-N-alpha-methyl-S-(t-butylthio)-L-cysteine CAS-No. 1013096-03-1 Formula $C_{23}H_{27}NO_4S_2$ Mol. weight 445,59 g/mol	 
FAA1716	Fmoc-L-Cys(tBu)-OH N-alpha-(9-Fluorenylmethyloxycarbonyl)-S-t-butyl-L-cysteine CAS-No. 67436-13-9 Formula $C_{22}H_{25}NO_4S$ Mol. weight 399,51 g/mol	 
FAA4160	Fmoc-L-Cys(Thp)-OH N-alpha-(9-Fluorenylmethyloxycarbonyl)-S-tetrahydropyranyl-L-cysteine CAS-No. 1673576-83-4 Formula $C_{23}H_{25}NO_5S$ Mol. weight 427,15 g/mol	 
FAA1040	Fmoc-L-Cys(Trt)-OH N-alpha-(9-Fluorenylmethyloxycarbonyl)-S-trityl-L-cysteine CAS-No. 103213-32-7 Formula $C_{37}H_{31}NO_4S$ Mol. weight 585,71 g/mol	 
FAA3570	Fmoc-L-MeCys(Trt)-OH N-alpha-(9-Fluorenylmethyloxycarbonyl)-N-alpha-methyl-S-trityl-L-cysteine CAS-No. 944797-51-7 Formula $C_{38}H_{33}NO_4S$ Mol. weight 599,74 g/mol	 
FAA1035	Fmoc-D-Cys(Trt)-OH N-alpha-(9-Fluorenylmethyloxycarbonyl)-S-trityl-D-cysteine CAS-No. 167015-11-4 Formula $C_{37}H_{31}NO_4S$ Mol. weight 585,71 g/mol	 

		Product details
FAA4840	Fmoc-alpha-Me-L-Cys(Trt)-OH N-alpha-(9-Fluorenylmethyloxycarbonyl)-alpha-methyl-S-trityl-L-cysteine CAS-No. 725728-43-8 Formula $C_{38}H_{33}NO_4S$ Mol. weight 599,74 g/mol	 
FAA4850	Fmoc-alpha-Me-D-Cys(Trt)-OH N-alpha-(9-Fluorenylmethyloxycarbonyl)-alpha-methyl-S-trityl-D-cysteine CAS-No. 725728-37-0 Formula $C_{38}H_{33}NO_4S$ Mol. weight 599,74 g/mol	 
FAA5670	Fmoc-L-Cys(Trt)-OMe N-alpha-(9-Fluorenylmethyloxycarbonyl)-S-trityl-L-cysteine methyl ester CAS-No. 245088-56-6 Formula $C_{38}H_{33}NO_4S$ Mol. weight 599,74 g/mol	 
ZAA1161	Z-L-Cys(Bzl)-OH N-alpha-Benzylloxycarbonyl-S-benzyl-L-cysteine CAS-No. 3257-18-9 Formula $C_{18}H_{19}NO_4S$ Mol. weight 345,42 g/mol	 
ZAA1310	Z-L-Cys(Trt)-OH N-alpha-Benzylloxycarbonyl-S-trityl-L-cysteine CAS-No. 26311-04-6 Formula $C_{30}H_{27}NO_4S$ Mol. weight 497,60 g/mol	 
PSI1450	Fmoc-L-Ala-L-Cys[PSI(Dmp,H)pro]-OH (S)-3-(N-(9-Fluorenylmethyloxycarbonyl)-L-alanyl)-2-(2,4-dimethoxyphenyl)thiazolidine-4-carboxylic acid CAS-No. 2022956-37-0 Formula $C_{30}H_{30}N_2O_7S$ Mol. weight 562,63 g/mol	 

Description of the Technology

Tools for the Formation of Thioesters







Building Blocks for Ligation




Enzyme-mediated Ligation

Examples

Index

[↑ back to content](#)

		Product details
PSI1470 Fmoc-L-Asp(tBu)-L-Cys[PSI(Dmp,H)pro]-OH (S)-3-(N-(9-Fluorenylmethyloxycarbonyl)-beta-t-butyl-L-aspartyl)-2-(2,4-dimethoxyphenyl)thiazolidine-4-carboxylic acid CAS-No. 1359754-16-7 Formula $C_{35}H_{39}N_2O_9S$ Mol. weight 662,75 g/mol		
PSI1580 Fmoc-L-Cys(Trt)-L-Cys(Psi(Dmp,H)pro)-OH (R)-3-(N-(9-Fluorenylmethyloxycarbonyl)-S-trityl-L-cysteiny)-2-(2,4-dimethoxyphenyl)thiazolidine-4-carboxylic acid CAS-No. 2022956-75-6 Formula $C_{49}H_{44}N_2O_7S_2$ Mol. weight 837,01 g/mol		
PSI1490 Fmoc-L-Glu(tBu)-L-Cys[PSI(Dmp,H)pro]-OH (S)-3-(N-(9-Fluorenylmethyloxycarbonyl)-gamma-t-butyl-L-glutamyl)-2-(2,4-dimethoxyphenyl)thiazolidine-4-carboxylic acid CAS-No. 2565804-44-4 Formula $C_{36}H_{40}N_2O_9S$ Mol. weight 676,78 g/mol		
PSI1440 Fmoc-Gly-L-Cys[PSI(Dmp,H)pro]-OH (S)-3-(N-(9-Fluorenylmethyloxycarbonyl)-glycyl)-2-(2,4-dimethoxyphenyl)thiazolidine-4-carboxylic acid CAS-No. 1926163-05-4 Formula $C_{29}H_{28}N_2O_7S$ Mol. weight 548,61 g/mol		
PSI1510 Fmoc-L-Leu-L-Cys[PSI(Dmp,H)pro]-OH (S)-3-(N-(9-Fluorenylmethyloxycarbonyl)-L-leucyl)-2-(2,4-dimethoxyphenyl)thiazolidine-4-carboxylic acid CAS-No. 1926163-06-5 Formula $C_{33}H_{36}N_2O_7S$ Mol. weight 604,71 g/mol		
PSI1520 Fmoc-L-Lys(Boc)-L-Cys[PSI(Dmp,H)pro]-OH (S)-3-(N-(9-Fluorenylmethyloxycarbonyl)-N-epsilon-tert-butyloxycarbonyl-L-lysyl)-2-(2,4-dimethoxyphenyl)thiazolidine-4-carboxylic acid CAS-No. 1926163-07-6 Formula $C_{38}H_{45}N_3O_9S$ Mol. weight 719,84 g/mol		

		Product details
PSI1560	Fmoc-L-Tyr(tBu)-L-Cys[PSI(Dmp,H)pro]-OH	 
(S)-3-(N-(9-Fluorenylmethyloxycarbonyl)-O-t-butyl-L-thyrosyl)-2-(2,4-dimethoxyphenyl)thiazolidine-4-carboxylic acid		
Formula	C ₄₀ H ₄₂ N ₂ O ₈ S	
Mol. weight	710,84 g/mol	
PSI1570	Fmoc-L-Val-L-Cys[PSI(Dmp,H)pro]-OH	 
(S)-3-(N-(9-Fluorenylmethyloxycarbonyl)-L-valyl)-2-(2,4-dimethoxyphenyl)thiazolidine-4-carboxylic acid		
CAS-No.	1926163-08-7	
Formula	C ₃₂ H ₃₄ N ₂ O ₇ S	
Mol. weight	590,69 g/mol	
SAA1110	Smoc-L-Cys(Trt)-OH	 
N-(((2,7-disulfo-9H-fluoren-9-yl)methoxy)carbonyl)-S-trityl-L-cysteine potassium salt		
CAS-No.	2442552-68-1	
Formula	C ₃₇ H ₂₉ K ₂ NO ₁₀ S ₃	
Mol. weight	822,01 g/mol	

Description of the Technology

Tools for the Formation of Thioesters

Building Blocks for Ligation

Enzyme-mediated Ligation


Examples

Index



Interested in Smoc amino acids?

Download our flyer!



[↑ back to content](#)


3.2. Ligation at the Position of Glycine


On average, cysteine occurs at a very low frequency in natural peptide sequences, a fact that frequently causes challenges when planning NCL strategies. On the other hand, glycine is a much more common amino acid, and being the simplest amino acid in nature, lends itself ideally for the construction of an NCL auxiliary.

The building block *t*Bu-SS-Photo(Fmoc)-Gly-OH is based on a glycine scaffold that is *N*-alkylated with a photolabile thiol-bearing auxiliary. This auxiliary mimics the action of an *N*-terminal cysteine's sulfhydryl group in NCL, and can be tracelessly removed following ligation, leaving only a glycine behind. Being able to perform NCL at the position of glycine residues enables the synthetic chemist with a much higher flexibility for the design of peptides and proteins. Following SPPS, this building block is attached to the *N*-terminus of a peptide sequence in lieu of a glycine residue. The auxiliary's Fmoc-protected amino functionality can subsequently be deprotected and functionalized, e.g., with a PEG, which is useful for increasing the solubility of peptide fragments, and for facilitating their purification by precipitation with EtOH/Et₂O. This is particularly valuable if the peptide's amino acid side-chains are supposed to be further derivatized post-SPPS.

Reference:

- A PEGylated photocleavable auxiliary mediates the sequential enzymatic glycosylation and native chemical ligation of peptides; C. Bello, S. Wang, L. Meng, K. W. Moremen, C. F. Becker; *Angew Chem Int Ed Engl* 2015; **54**: 7711-7715. <https://doi.org/10.1002/anie.201501517>


		Product details
PAA2000 tBu-SS-Photo(Fmoc)-Gly-OH		
Photocleavable-NCL-auxiliary-Gly-OH		
CAS-No.	1994388-93-0	
Formula	C ₃₆ H ₄₄ N ₄ O ₉ S ₂	
Mol. weight	740,89 g/mol	
		



3.3. Building Blocks and Reagents for Ligation with Seleno Amino Acids

General Procedure for the Synthesis of Selenoesters

The peptide residue after Fmoc-SPPS is dissolved in anhydrous DMF and cooled to 0 °C. Diphenyl diselenide (DPDS) (30 eq. in DMF) and Bu₃P (30 eq.) are subsequently added. After 3 h at 0 °C, the solvent is removed *in vacuo*.



One-Pot Additive-Free Diselenide-Selenoester Ligation-Deselenization Reactions

Conditions of Additive-free Ligation:

2.5 mM final concentration of diselenide dimer in 6 M Gdn*HCl, and 0.1 M Na₂HPO₄ (pH 7.2; reduced to 6.2–6.5 upon addition to peptide fragments).







Conditions of One-Pot Deselenization:

Hexane extraction (×5) followed by the addition of 0.25 M TCEP, 25 mM DTT in 6 M Gdn*HCl, and 0.1 M Na₂HPO₄ (pH 5–6).

Reference:

→ Accelerated Protein Synthesis via One-Pot Ligation-Deselenization Chemistry; N. J. Mitchell, J. Sayers, S. S. Kulkarni, D. Clayton, A. M. Goldys, J. Ripoll-Rozada, P. J. Barbosa Pereira, B. Chan, L. Radom, R. J. Payne; *Chem* 2017; **2**: 703-715. <https://doi.org/10.1016/j.chempr.2017.04.003>

Our portfolio contains a variety of selenocysteine (Sec) derivatives as well as selenazolidine carboxylic acids (Sez derivatives). Sez can be deprotected and converted to Sec by treatment with O-methylhydroxylamine (MeONH₂) at pH 4 or by using Cu(II) salts.

		Product details
HAA9255	H-L-Sec(oNv)-OH*TFA Dimethoxynitrobenzyl selenocysteine TFA salt CAS-No. 1644398-13-9 Formula $C_{12}H_{16}N_2O_6Se \cdot CF_3COOH$ Mol. weight 363,24*114,02 g/mol	 
HAA9360	H-L-Sec(MDNPE)-OH Se-(Methyl-o-nitropiperonyl)-selenocysteine CAS-No. 2235373-47-2 Formula $C_{12}H_{14}N_2O_6Se$ Mol. weight 361,21 g/mol	 
HAA9230	H-L-Sec(MDNPE)*TFA (2R)-2-amino-3-((1-(6-nitrobenzo[d][1,3]dioxol-5-yl)ethyl)selanyl)propanoic acid trifluoroacetate CAS-No. 2235373-48-3 Formula $C_{12}H_{14}N_2O_6Se \cdot CF_3CO_2H$ Mol. weight 361,22*114,02 g/mol	 
HAA9350	(H-L-Sec-OH)2 L-Selenocysteine, (H-Sec)2, (H-L-Sec)2 CAS-No. 29621-88-3 Formula $C_6H_{12}N_2O_4Se_2$ Mol. weight 334,11 g/mol	 
HAA9465	H-L-Sec(oNB)-OH*HCl (R)-2-amino-3-((2-nitrobenzyl)selanyl)propanoic acid CAS-No. 324582-23-2 net Formula $C_{10}H_{12}N_2O_4Se \cdot HCl$ Mol. weight 303,18*36,46 g/mol	 
HAA9475	H-L-Sec(NPE)-OH*HCl (2R)-2-amino-3-((1-(2-nitrophenyl)ethyl)selanyl)propanoic acid Formula $C_{11}H_{14}N_2O_4Se \cdot HCl$ Mol. weight 317,02*36,46 g/mol	 

Description of the Technology

Tools for the Formation of Thioesters







Building Blocks for Ligation







Enzyme-mediated Ligation

Examples

Index

[↑ back to content](#)

		Product details
BAA3680 (Boc-L-Sec)₂ N-α-t-Butyloxycarbonyl-L-selenocystine CAS-No. 877754-71-7 Formula C ₁₆ H ₂₈ N ₂ O ₈ Se ₂ Mol. weight 534,35 g/mol		
BAA3760 Boc-L-Sec(Mob)-OH N-α-t-Butyloxycarbonyl-Se-(4-methoxybenzyl)-L-selenocysteine CAS-No. 959415-39-5 Formula C ₁₆ H ₂₃ NO ₅ Se Mol. weight 388,32 g/mol		
BAA4830 Boc-L-Sec(Mob)-OPfp N-α-tert-Butoxycarbonyl-4-methoxybenzyl-L-selenocysteine pentafluorophenyl ester CAS-No. 1257525-48-6 Formula C ₂₂ H ₂₃ F ₅ NO ₅ Se Mol. weight 554,38 g/mol		
FAA8705 Fmoc-L-Sec(Mob)-OH N-α-(9-Fluorenylmethyloxycarbonyl)-Se-(4-methoxybenzyl)-L-selenocysteine CAS-No. 150308-80-8 Formula C ₂₆ H ₂₅ NO ₅ Se Mol. weight 510,46 g/mol		
FAA8760 Fmoc-L-Sec(Mob)-OPfp N-α-(9-Fluorenylmethyloxycarbonyl)-L-4-methoxybenzyl selenocysteine pentafluorophenyl ester CAS-No. 939431-43-3 Formula C ₃₂ H ₂₄ F ₅ NO ₅ Se Mol. weight 676,51 g/mol		
FAA8710 Fmoc-D-Sec(Mob)-OH N-α-(9-Fluorenylmethyloxycarbonyl)-Se-(4-methoxybenzyl)-D-selenocysteine CAS-No. 2987041-11-0 Formula C ₂₆ H ₂₅ NO ₅ Se Mol. weight 510,46 g/mol		

		Product details
FAA8465	Fmoc-L-Sec(Xan)-OH	
Fmoc-Se-xanthyl-L-selenocysteine		
CAS-No.	1639843-35-8	
Formula	C ₃₁ H ₂₅ NO ₅ Se	
Mol. weight	570,49 g/mol	
		
FAA8600	Fmoc-D-Sec(Xan)-OH	
Fmoc-Se-xanthyl-D-selenocysteine		
CAS-No.	2988660-46-2	
Formula	C ₃₁ H ₂₅ NO ₅ Se	
Mol. weight	570,49 g/mol	
		
BAA4880	Boc-L-Sez-OH	
Boc selenazolidine carboxylic acid		
CAS-No.	1841180-44-6	
Formula	C ₉ H ₁₅ NO ₄ Se	
Mol. weight	280,19 g/mol	
		
FAA8860	Fmoc-L-Sez-OH	
Fmoc selenazolidine carboxylic acid		
CAS-No.	1985651-74-8	
Formula	C ₁₉ H ₁₇ NO ₄ Se	
Mol. weight	402,31 g/mol	
		
LS-3405	TCEP*HCl	
3-[bis(2-carboxyethyl)phosphanyl]propanoic acid, hydrochloride		
CAS-No.	51805-45-9	
Formula	C ₉ H ₁₅ O ₆ P*HCl	
Mol. weight	250,19*36,45 g/mol	
		
RL-1020	DTT (racemic)	
DL-Dithiothreitol		
CAS-No.	3483-12-3	
Formula	C ₄ H ₁₀ O ₂ S ₂	
Mol. weight	154,25 g/mol	
		

Description of the Technology

Tools for the Formation of Thioesters

Building Blocks for Ligation

Enzyme-mediated Ligation

Examples

Index

[↑ back to content](#)

3.4. Oxime Ligation using Aminoxy-Amino Acids

Replacing the primary amino group of an amino acid by an aminoxy moiety leads to an increase in nucleophilicity. Thus, after completion of the peptide synthesis and deprotection of the aminoxy-function, chemoselective reactions with carbonyl compounds under formation of a kinetically stable oxime bond can be performed. Compared to imines, oximes display much higher stability toward hydrolysis. This increased stability is explained by the α -effect provided by the heteroatom adjacent to the sp^2 nitrogen.

This simple reaction can be utilized for the synthesis of large proteins by fragment assembly = oxime ligation. Successive deprotection of aldehyde and/or aminoxy moieties allows for the controlled and sequential ligation of peptide fragments.

The oxime linkage is reported as good peptidomimetic instead of the amide bond. The low usage of this ligation type is most probably due to the possible hydrolysis of the oxime linkage as well as the formation of a mixture of products due to *E/Z*-isomerization.










Fig. 12: Illustration of oxime ligation.

References:

- Site-specific cross-linking of proteins to DNA via a new biorthogonal approach employing oxime ligation; S. S. Pujari, Y. Zhang, S. Ji, M. D. Distefano, N. Y. Tretyakova; **Chem. Commun.** 2018; **54**: 6296-6299. <https://doi.org/10.1039/C8CC01300D>
- 44. Amino-oxy-derivatives. Part I. Some α -amino-oxy-acids and α -amino-oxy-hydrazides. D. Mchale, J. Green, P. Mamalis; **J. Chem. Soc.** 1960; 225-229. <https://doi.org/10.1039/JR9600000225>
- SAR by Oxime-Containing Peptide Libraries: Application to Tsg101 Ligand Optimization; F. Liu, A. G. Stephen, A. A. Waheed, M. J. Aman, E. O. Freed, R. J. Fisher, T. R. Burke; **ChemBioChem** 2008; **9(12)**: 2000-2004. <https://doi.org/10.1002/cbic.200800281>
- A Versatile Set of Aminoxy Amino Acids for the Synthesis of Neoglycopeptides; M. R. Carrasco, R. T. Brown; **J. Org. Chem.** 2003; **68**: 8853-8858. <https://doi.org/10.1021/jo034984x>

		Product details
FAA1461	Fmoc-L-Dap(Boc-Aoa)-OH N-alpha-(9-Fluorenylmethyloxycarbonyl)-N-beta-(t-butyloxycarbonyl-aminoxy-acetyl)-L-2,3-diaminopropionic acid CAS-No. 600153-12-6 Formula $C_{25}H_{29}N_3O_8$ Mol. weight 499,53 g/mol	 
FAA8445	Fmoc-Hcan(Boc)-OH (S) (S)-2-(Fmoc-amino)-5-(Boc-aminoxy)pentanoic acid CAS-No. 204844-15-5 Formula $C_{25}H_{30}N_2O_7$ Mol. weight 470,51 g/mol	 
FAA1995	Fmoc-L-His(3-Me)-OH N-alpha-(9-Fluorenylmethyloxycarbonyl)-N ₃ -methyl-L-histidine CAS-No. 252049-16-4 Formula $C_{22}H_{21}N_3O_4$ Mol. weight 391,42 g/mol	 
FAA4370	Fmoc-L-Lys(Boc-Aoa)-OH N-alpha-(9-Fluorenylmethyloxycarbonyl)-N-epsilon-(t-butyloxycarbonyl)aminoxyacetyl-L-lysine CAS-No. 757960-24-0 Formula $C_{28}H_{35}N_3O_8$ Mol. weight 541,59 g/mol	 
FAA4700	Fmoc-L-Lys(Mtt-Aoa)-OH N-alpha-(9-Fluorenylmethyloxycarbonyl)-N-epsilon-(4-methyltrityl)aminoxyacetyl-L-lysine CAS-No. 2250436-45-2 Formula $C_{43}H_{43}N_3O_6$ Mol. weight 697,82 g/mol	 
FAA8450	Fmoc-AAHA(Boc)-OH (S) (S)-2-(Fmoc-amino)-6-(Boc-aminoxy)hexanoic acid CAS-No. 357278-11-6 Formula $C_{26}H_{32}N_2O_7$ Mol. weight 484,54 g/mol	 

Description of the Technology

Tools for the Formation of Thioesters



Building Blocks for Ligation

Enzyme-mediated Ligation

Examples


Index


[↑ back to content](#)

			Product details
FAA8455	Fmoc-L-<i>trans</i>-Hyp(NHBoc)-OH		
	Fmoc-4-(Boc-aminooxy)-proline (2S,4R)		
CAS-No.	1015426-45-5		
Formula	C ₂₅ H ₂₈ N ₂ O ₇		
Mol. weight	468,50 g/mol		
FAA8460	Fmoc-L-<i>cis</i>-Hyp(NHBoc)-OH		
	Fmoc-4-(Boc-aminooxy)-proline (2S,4S)		
CAS-No.	1015426-31-9		
Formula	C ₂₅ H ₂₈ N ₂ O ₇		
Mol. weight	468,50 g/mol		

3.5. Activated Cysteine-Based Protein Ligation (ACPL)

This method is based on the direct activation of cysteine by the small molecule cyanating reagent 2-nitro-5-thiocyanatobenzoic acid (NTCB). The transferred cyanide then forms a thiocyanate which undergoes a reversible intramolecular addition with the cysteine N-amide to generate a five-membered 1-acyl-2-iminothiazolidine intermediate which can be reacted with nucleophiles. When the nucleophile is hydrazine, the afforded peptide hydrazide can then undergo ligation, either *via* reaction with an aldehyde/ketone (hydrazone ligation) or *via* transformation to a peptide thioester and subsequent chemical ligation, both as described above. Without addition of a nucleophilic amine, the formed five-membered ring undergoes hydrolysis.



		Product details
RL-4080	NTCB	
2-nitro-5-thiocyanatobenzoic acid		
CAS-No.	30211-77-9	
Formula	C ₈ H ₄ N ₂ O ₄ S	
Mol. weight	224,19 g/mol	
		

Reference:

→ Expressed Protein Ligation without Intein; Y. Qiao, G. Yu, K. C. Kratch, X. Aria Wang, W. Wei Wang, S. Z. Leeuwon, S. Xu, J. S. Morse, W. Ray Liu; *J. Am. Chem. Soc.* 2002; **124**: 7057-7054. <https://doi.org/10.1021/jacs.0c00252>

3.6. Nitrile-Aminothiol Conjugation (NATC) as Biocompatible Ligation Technique

While nitriles in general may react with sulfhydryl groups (e.g., from cysteines) in a rather unspecific way, the advantage of α -cyanopyridines (2-picolinonitrils) is their strong selectivity for 1,2-aminothiols (e.g., provided as N-terminal cysteine, internally as non-canonical amino acid, or by conversion of thiazolidines). The so-called nitrile-aminothiol conjugation (NATC) proceeds rapidly under physiological pH conditions without the need for any added catalyst and leads to the formation of a thiazolidine ring.




















Fig. 14: The reaction between a 2-cyanopyridine and a 1,2-aminothiol forms a 2-thiazoline.

		Product details
FAA9370	Fmoc-L-3-(2-cyano-4-pyridyl)-alanine	
(S)-2-((((9H-fluoren-9-yl)methoxy)carbonyl)amino)-3-(2-cyanopyridin-4-yl)propanoic acid		
CAS-No.	2245755-79-5	
Formula	C ₂₄ H ₁₉ N ₃ O ₄	
Mol. weight	413,43 g/mol	
		

[↑ back to content](#)

		Product details
FAA9375	Fmoc-L-3-(2-cyano-3-pyridyl)-alanine (S)-2-((((9H-fluoren-9-yl)methoxy)carbonyl)amino)-3-(2-cyanopyridin-3-yl)propanoic acid CAS-No. 2245755-77-3 Formula $C_{24}H_{19}N_3O_4$ Mol. weight 413,43 g/mol	 
FAA9380	Fmoc-L-Dap(2-CINA)-OH (S)-2-((((9H-fluoren-9-yl)methoxy)carbonyl)amino)-3-(2-cyanoisonicotinamido)propanoic acid Formula $C_{25}H_{20}N_4O_5$ Mol. weight 456,46 g/mol	 
FAA9385	Fmoc-L-Dap(6-CNA)-OH (S)-2-((((9H-fluoren-9-yl)methoxy)carbonyl)amino)-3-(6-cyanonicotinamido)propanoic acid Formula $C_{25}H_{20}N_4O_5$ Mol. weight 456,46 g/mol	 
FAA9390	Fmoc-L-Dab(2,6-DCP)-OH (S)-2-((((9H-fluoren-9-yl)methoxy)carbonyl)amino)-4-((2,6-dicyanopyridin-4-yl)amino)butanoic acid CAS-No. 2968514-51-2 Formula $C_{26}H_{21}N_5O_4$ Mol. weight 467,49 g/mol	 
FAA9395	Fmoc-L-Cys(2,6-DCP)-OH N-((((9H-fluoren-9-yl)methoxy)carbonyl)-S-(2,6-dicyanopyridin-4-yl)-L-cysteine CAS-No. 2968514-50-1 Formula $C_{25}H_{18}N_4O_4S$ Mol. weight 470,50 g/mol	 
RL-8625	N₃-PEG(4)-CINA N-(14-azido-3,6,9,12-tetraoxatetradecyl)-2-cyanoisonicotinamide Formula $C_{17}H_{24}N_6O_5$ Mol. weight 392,42 g/mol	 

			Product details
RL-8630	Alkyne-PEG(4)-CINA		
2-cyano-N-(3,6,9,12-tetraoxapentadec-14-yn-1-yl)isonicotinamide			
Formula	$C_{18}H_{23}N_3O_5$		
Mol. weight	361,40 g/mol		
RL-8635	Biotin-PEG(4)-CINA		
2-cyano-N-(16-oxo-20-((3aS,4S,6aR)-2-oxohexahydro-1H-thieno[3,4-d]imidazol-4-yl)-3,6,9,12-tetraoxa-15-azaicosyl)isonicotinamide			
Formula	$C_{27}H_{40}N_6O_7S$		
Mol. weight	592,71 g/mol		
RL-8640	DBCO-PEG(4)-CINA		
Dibenzoazacyclooctyne-tetra(ethylene glycol)-cyanoisonicotinamide			
Formula	$C_{36}H_{39}N_3O_7$		
Mol. weight	653,74 g/mol		
FAA9315	Fmoc-L-Lys(Boc-Cys(Trt))-OH		
N2-(((9H-fluoren-9-yl)methoxy)carbonyl)-N6-(N-(tert-butoxycarbonyl)-S-trityl-L-cysteinyl)-L-lysine			
CAS-No.	587854-43-1		
Formula	$C_{48}H_{51}N_3O_7S$		
Mol. weight	814,01 g/mol		
FAA9320	Fmoc-L-Lys(Boc-Thz)-OH		
N2-(((9H-fluoren-9-yl)methoxy)carbonyl)-N6-((R)-3-(tert-butoxycarbonyl)thiazolidine-4-carbonyl)-L-lysine			
Formula	$C_{30}H_{37}N_3O_7S$		
Mol. weight	583,70 g/mol		
FAA9325	Fmoc-L-Dab(Boc-Cys(Trt))-OH		
(S)-2-(((9H-fluoren-9-yl)methoxy)carbonyl)amino-4-((R)-2-((tert-butoxycarbonyl)amino)-3-(tritylthio)propanamido)butanoic acid			
CAS-No.	2968514-52-3		
Formula	$C_{46}H_{47}N_3O_7S$		
Mol. weight	785,96 g/mol		


		Product details
FAA9330	Fmoc-L-Dab(Boc-Thz)-OH	
(S)-2-((((9H-fluoren-9-yl)methoxy)carbonyl)amino)-4-((R)-3-(tert-butoxycarbonyl)thiazolidine-4-carboxamido)butanoic acid		
CAS-No.	2968514-54-5	
Formula	C ₂₈ H ₃₃ N ₃ O ₇ S	
Mol. weight	555,65 g/mol	
FAA9335	Fmoc-L-Lys(5-STrt, Boc)-OH	
(2S)-2-((((9H-fluoren-9-yl)methoxy)carbonyl)amino)-6-((tert-butoxycarbonyl)amino)-5-(tritylthio)hexanoic acid		
CAS-No.	1240666-29-8	
Formula	C ₄₅ H ₄₆ N ₂ O ₆ S	
Mol. weight	742,93 g/mol	
FAA9340	Fmoc-L-Lys(4-Thz, Boc)-OH	
(2S)-2-((((9H-fluoren-9-yl)methoxy)carbonyl)amino)-4-(3-(tert-butoxycarbonyl)thiazolidin-5-yl)butanoic acid		
CAS-No.	1240666-28-7	
Formula	C ₂₇ H ₃₂ N ₂ O ₆ S	
Mol. weight	512,62 g/mol	

References:



- *Advancing Nitrile-Aminothiol Strategy for Dual and Sequential Bioconjugation*; V. J. Thombare, Y. Wu, K. Pamulapati, M. Han, J. Tailhades, M. J. Cryle, K. D. Roberts, T. Velkov, J. Li, N. A. Patil; **Chem. Eur. J.** 2024; **30**(46): e202401674. <https://doi.org/10.1002/chem.202401674>
- *The Cyanopyridine-Aminothiol Click Reaction: Expanding Horizons in Chemical Biology*; C. Nitsche; **SynLett.** 2024; **35**: A-E. <https://dx.doi.org/10.1055/a-2214-7612>




4. Enzyme-mediated Ligation

The transpeptidase activity of sortase can be used to produce fusion proteins *in vitro*. The enzyme recognition motif (LPXTG) is added to the C-terminus of a protein of interest, while an oligo(glycine) motif is added to the N-terminus of the second protein. Upon addition of sortase A, the two peptides are covalently linked through a native peptide bond while losing one glycine residue.



Here we present several building blocks, which can be used for convenient N-terminal oligo(glycine) design of ligation fragments.

			Product details
HAA2860	N₃-Gly-Gly-Gly-Gly-Gly-OH		
CAS-No.	2250433-77-1		
Formula	C ₁₀ H ₁₅ N ₇ O ₆		
Mol. weight	329,27 g/mol		
HAA2840	N₃-Gly-Gly-Gly-OH		
Azido-glycylglycylglycine			
CAS-No.	1993176-75-2		
Formula	C ₆ H ₉ N ₅ O ₄		
Mol. weight	215,17 g/mol		

		Product details
HAA2850 N₃-Gly-Gly-OH*DCHA Azido-glycylglycine dicyclohexylamine CAS-No. 855750-87-7 net Formula C ₄ H ₆ N ₄ O ₃ *C ₁₂ H ₂₃ N Mol. weight 158,12*181,32 g/mol		
FDP1030 Fmoc-Gly-Gly-OH N-alpha-(9-Fluorenylmethyloxycarbonyl)-glycylglycine CAS-No. 35665-38-4 Formula C ₁₉ H ₁₈ N ₂ O ₅ Mol. weight 354,36 g/mol		
HAA2870 H-(Gly)₃-Lys(N₃)-OH*HCl Triglycyl-epsilon-azido-L-lysine hydrochloride CAS-No. 2250437-45-5 Formula C ₁₂ H ₂₁ N ₇ O ₅ *HCl Mol. weight 343,34*36,45 g/mol		

References:

- Assembly of Oligoglycine Layers on Mica Surface; S. V. Tsygankova, A. A. Chinarev, A. B. Tuzikov, I. S. Zaitsev, N. Severin, A. A. Kalachev, J. P. Rabe, N. V. Bovin; **Journal of Biomaterials and Nanobiotechnology** 2011; **02**: 91-97.
<https://doi.org/10.4236/jbnt.2011.21012>
- Biantennary oligoglycines and glyco-oligoglycines self-associating in aqueous medium; S. V. Tsygankova, A. A. Chinarev, A. B. Tuzikov, N. Severin, A. A. Kalachev, J. P. Rabe, A. S. Gambaryan, N. V. Bovin; **Beilstein J Org Chem** 2014; **10**: 1372-1382. <https://doi.org/10.3762/bjoc.10.140>
- DNA-bending finger: artificial design of 6-zinc finger peptides with polyglycine linker and induction of DNA bending; M. Imanishi, Y. Hori, M. Nagaoka, Y. Sugiura; **Biochemistry** 2000; **39**: 4383-4390.
<https://doi.org/10.1021/bi992989b>
- Effects of the length of a glycine linker connecting the N- and C-termini of a circularly permuted dihydrofolate reductase; M. Iwakura, T. Nakamura; **Protein Eng** 1998; **11**: 707-713. <https://doi.org/10.1093/protein/11.8.707>
- Linkers in the structural biology of protein-protein interactions; V. P. Reddy Chichili, V. Kumar, J. Sivaraman; **Protein Sci** 2013; **22**: 153-167. <https://doi.org/10.1002/pro.2206>
- Optimizing the stability of single-chain proteins by linker length and composition mutagenesis; C. R. Robinson, R. T. Sauer; **Proc Natl Acad Sci U S A** 1998; **95**: 5929-5934. <https://doi.org/10.1073/pnas.95.11.5929>
- Sortase-mediated protein ligation: a new method for protein engineering; H. Mao, S. A. Hart, A. Schink, B. A. Pollok; **J Am Chem Soc** 2004; **126**: 2670-2671. <https://doi.org/10.1021/ja039915e>
- Sortase-tag expressed protein ligation: combining protein purification and site-specific bioconjugation into a single step; R. Warden-Rothman, I. Caturegli, V. Popik, A. Tsourkas; **Anal Chem** 2013; **85**: 11090-11097.
<https://doi.org/10.1021/ac402871k>


5. Examples

5.1. Kinetically Controlled Four-Fragment Ligation for the Convergent Chemical Synthesis of Proteins

The convergent synthesis of proteins by the chemical ligation of multiple peptide segments is achieved by a kinetically controlled four-fragment ligation. This was made possible by the controlled extension of the key Cys-peptide intermediate by ligation at either the N- or C-termini. For more information, please refer to the following publication:

Reference:

→ *Kinetically controlled ligation for the convergent chemical synthesis of proteins*; D. Bang, B. L. Pentelute, S. B. Kent; *Angew Chem Int Ed Engl* 2006; **45**: 3985-3988.



5.2. Three-Fragment Synthesis of a 116mer with Desulfurization

The discovery of native chemical ligation, the regio- and chemoselective coupling of two unprotected peptide segments, facilitated the synthesis of polypeptides with more than 200 amino acids. However, this approach initially relied on the presence of at least one cysteine residue in the sequence at a convenient position. Hence, postligation-desulfurization protocols were developed that allowed for ligation at amino acid residues other than cysteine. Please refer to the following review for a comprehensive overview of “the development and recent progress on the chemical synthesis of peptides and proteins encompassing postligation-desulfurization at alanine, valine, lysine, threonine, leucine, proline, arginine, aspartic acid, glutamate, phenylalanine, glutamine, and tryptophan”.

Reference:

→ *Postligation-Desulfurization: A General Approach for Chemical Protein Synthesis*; J. Ma, J. Zeng, Q. Wan; **Protein Ligation and Total Synthesis II** 2015; (Ed.: L. Liu): Springer International Publishing, Cham, 57-101.




Fig. 17: Example of a synthesis of a peptide via ligation of three segments and thiazolidine as temporarily protected cysteine.

Index

Product code	Product name	Page	Product code	Product name	Page
BAA3680	(Boc-L-Sec)2	39	FAA8150	Fmoc-D-Cys(Msbh)-OH	29
HAA9350	(H-L-Sec-OH)2	38	FAA3710	Fmoc-D-Cys(Phacm)-OH	30
AAA1300	Ac-L-Cys-OH	22	FAA1965	Fmoc-D-Cys(StBu)-OH	31
RL-8630	Alkyne-PEG(4)-CINA	46	FAA1035	Fmoc-D-Cys(Trt)-OH	32
AAA2015	Aloc-L-Cys(Trt)-OH	22	FAA1470	Fmoc-D-Cys-OH*H ₂ O	26
RL-8635	Biotin-PEG(4)-CINA	46	FAA8710	Fmoc-D-Sec(Mob)-OH	39
BAA5410	Boc-D-Cys(Bzl)-OH	23	FAA8600	Fmoc-D-Sec(Xan)-OH	40
BAA5420	Boc-D-Cys(MBzl)-OH	24	FAA1495	Fmoc-D-Thz-OH	18
BAA5000	Boc-D-Cys(Trt)-OH	25	FAA3168	Fmoc-Dbz(o-Alloc)-OH	9
BAA1170	Boc-D-Cys-OH	23	FAA3167	Fmoc-Dbz(o-Boc)-OH	9
BAA1186	Boc-D-Thz-OH	17	FAA3169	Fmoc-Dbz(o-NO ₂)-OH	9
BAA1078	Boc-L-Cys(Acm)-OH	23	FAA3165	Fmoc-Dbz-OH	8
BAA1510	Boc-L-Cys(Acm,O)-OH	23	FDP1030	Fmoc-Gly-Gly-OH	49
BAA1079	Boc-L-Cys(Bzl)-OH	23	PSI1440	Fmoc-Gly-L-Cys[PSI(Dmp,H)pro]-OH	33
BAA5510	Boc-L-Cys(Fm)-OH	24	PYV1170	Fmoc-Gly-NHN=Pyv Resin	14
BAA1080	Boc-L-Cys(MBzl)-OH	24	FAA8445	Fmoc-Hcan(Boc)-OH (S)	42
BAA1081	Boc-L-Cys(Mob)-OH	24	FAA9375	Fmoc-L-3-(2-cyano-3-pyridyl)-alanine	45
BAA1860	Boc-L-Cys(Npys)-OH	24	FAA9370	Fmoc-L-3-(2-cyano-4-pyridyl)-alanine	44
BAA3140	Boc-L-Cys(Ph)-OH	25	PSI1450	Fmoc-L-Ala-L-Cys[PSI(Dmp,H)pro]-OH	33
BAA6415	Boc-L-Cys(StBu)-OH	25	PYV1100	Fmoc-L-Ala-NHN=Pyv Resin	12
BAA1082	Boc-L-Cys(tBu)-OH	25	PYV1110	Fmoc-L-Arg(Pbf)-NHN=Pyv Resin	13
BAA1084	Boc-L-Cys(Trt)-OH	25	PYV1120	Fmoc-L-Asn(Trt)-NHN=Pyv Resin	13
BAA1083	Boc-L-Cys-OH	23	PYV1130	Fmoc-L-Asp(OrtBu)-NHN=Pyv Resin	13
BAA3760	Boc-L-Sec(Mob)-OH	39	PSI1470	Fmoc-L-Asp(tBu)-L-Cys[PSI(Dmp,H)pro]-OH	33
BAA4830	Boc-L-Sec(Mob)-OPfp	39	FAA8460	Fmoc-L-cis-Hyp(NHBoc)-OH	43
IAD1040	Boc-L-Ser[Fmoc-L-Cys(Trt)]-OH	25	FAA9395	Fmoc-L-Cys(2,6-DCP)-OH	45
BAA4880	Boc-L-Sez-OH	40	FAA5150	Fmoc-L-Cys(Aapam)-OH	26
IAD2040	Boc-L-Thr[Fmoc-L-Cys(Trt)]-OH	26	FAA4751	Fmoc-L-Cys(Ac-OrtBu)-OH*DCHA	27
BAA1135	Boc-L-Thz-OH	17	FAA1506	Fmoc-L-Cys(Acm)-OH	26
RL-8640	DBCO-PEG(4)-CINA	46	FAA7610	Fmoc-L-Cys(Allocam)-OH	27
RL-1020	DTT (racemic)	40	FAA6270	Fmoc-L-Cys(Bzl)-OH	27
FAA8450	Fmoc-AAHA(Boc)-OH (S)	42	FAA6940	Fmoc-L-Cys(Ddm)-OH	27
FAA4850	Fmoc-alpha-Me-D-Cys(Trt)-OH	32	FAA3190	Fmoc-L-Cys(Dpm)-OH	28
FAA4845	Fmoc-alpha-Me-L-Cys(Mmt)-OH	29	FAA1714	Fmoc-L-Cys(MBzl)-OH	28
FAA4840	Fmoc-alpha-Me-L-Cys(Trt)-OH	32	FAA7945	Fmoc-L-Cys(MDNPE)-OH	28
FAA6230	Fmoc-D-Cys(Acm)-OH	27	FAA1030	Fmoc-L-Cys(Mmt)-OH	28
FAA5650	Fmoc-D-Cys(Dpm)-OH	28	FAA1715	Fmoc-L-Cys(Mob)-OH	29
FAA1614	Fmoc-D-Cys(Mmt)-OH	29	FAA4155	Fmoc-L-Cys(Msbh)-OH	29

[↑ back to content](#)

Description of the Technology

Tools for the Formation of Thioesters

Building Blocks for Ligation

Enzyme-mediated Ligation

Examples

Index

Product code	Product name	Page	Product code	Product name	Page
FAA8420	Fmoc-L-Cys(NDBF)-OH	29	PYV1240	Fmoc-L-Pro-NHN=Pyv Resin	15
FAA3970	Fmoc-L-Cys(oNv)-OH	30	FAA8705	Fmoc-L-Sec(Mob)-OH	39
FAA6910	Fmoc-L-Cys(Phacm)-OH	30	FAA8760	Fmoc-L-Sec(Mob)-OPfp	39
FAA3180	Fmoc-L-Cys(S-DMP)-OH	28	FAA8465	Fmoc-L-Sec(Xan)-OH	40
FAA8495	Fmoc-L-Cys(SIT)-OH	30	PYV1250	Fmoc-L-Ser(tBu)-NHN=Pyv Resin	15
FAA1575	Fmoc-L-Cys(SiBu)-OH	30	FAA8860	Fmoc-L-Sez-OH	40
FAA1716	Fmoc-L-Cys(tBu)-OH	31	PYV1260	Fmoc-L-Thr(tBu)-NHN=Pyv Resin	15
FAA4160	Fmoc-L-Cys(Thp)-OH	31	FAA9200	Fmoc-L-Thz(Dmp)-OH	18
PSI1580	Fmoc-L-Cys(Trt)-L-Cys(Psi(Dmp,H)pro)-OH	33	FAA1437	Fmoc-L-Thz(Me2)-OH	18
PYV1140	Fmoc-L-Cys(Trt)-NHN=Pyv Resin	13	FAA1427	Fmoc-L-Thz-OH	18
FAA1040	Fmoc-L-Cys(Trt)-OH	31	FAA8455	Fmoc-L-trans-Hyp(NHBoc)-OH	43
FAA5670	Fmoc-L-Cys(Trt)-OMe	32	PYV1270	Fmoc-L-Trp(Boc)-NHN=Pyv Resin	15
FAA1980	Fmoc-L-Cys-NH ₂	26	PSI1560	Fmoc-L-Tyr(tBu)-L-Cys[PSI(Dmp,H)pro]-OH	34
FAA1362	Fmoc-L-Cys-OH*H ₂ O	26	PYV1280	Fmoc-L-Tyr(tBu)-NHN=Pyv Resin	15
FAA9390	Fmoc-L-Dab(2,6-DCP)-OH	45	PSI1570	Fmoc-L-Val-L-Cys[PSI(Dmp,H)pro]-OH	34
FAA9325	Fmoc-L-Dab(Boc-Cys(Trt))-OH	46	PYV1290	Fmoc-L-Val-NHN=Pyv Resin	16
FAA9330	Fmoc-L-Dab(Boc-Thz)-OH	18, 47	FAA3166	Fmoc-MeDbz-OH	9
FAA9380	Fmoc-L-Dap(2-CINA)-OH	45	S-30091	Fmoc-MeDbz-RAM TG S	9
FAA9385	Fmoc-L-Dap(6-CNA)-OH	45	PYV1000	Fmoc-NHN=Pyv Resin	12
FAA1461	Fmoc-L-Dap(Boc-Aoa)-OH	42	BR-5279	Fmoc-NHNH-2CT Resin	12
PYV1160	Fmoc-L-Gln(Trt)-NHN=Pyv Resin	13	BR-5280	Fmoc-NHNH-Trt-PEG Resin	12
PSI1490	Fmoc-L-Glu(tBu)-L-Cys[PSI(Dmp,H)pro]-OH	33	PEG8260	Fmoc-NHNH-Trt-PEG Resin HL	12
PYV1150	Fmoc-L-Glu(tBu)-NHN=Pyv Resin	13	PEG8265	Fmoc-NHNH-Trt-PEG Resin XV	12
FAA1995	Fmoc-L-His(3-Me)-OH	42	HAA2870	H-(Gly)3-Lys(N ₃)-OH*HCl	49
PYV1180	Fmoc-L-His(Trt)-NHN=Pyv Resin	14	HAA6110	H-D-Cys(Bzl)-OH	19
PYV1190	Fmoc-L-Ile-NHN=Pyv Resin	14	HAA3500	H-D-Cys(Mmt)-OH	20
PSI1510	Fmoc-L-Leu-L-Cys[PSI(Dmp,H)pro]-OH	33	HAA6120	H-D-Cys(Trt)-OH	21
PYV1200	Fmoc-L-Leu-NHN=Pyv Resin	14	HAA3520	H-D-Cys(Trt)-OMe*HCl	21
FAA9340	Fmoc-L-Lys(4-Thz, Boc)-OH	19, 47	HAA2100	H-D-Cys(Trt)-OtBu*HCl	21
FAA9335	Fmoc-L-Lys(5-STrt, Boc)-OH	47	HAA1017	H-D-Cys-OH*HCl*H ₂ O	19
PSI1520	Fmoc-L-Lys(Boc)-L-Cys[PSI(Dmp,H)pro]-OH	34	HAA6070	H-L-Cys(Acm)-OH*HCl	19
PYV1210	Fmoc-L-Lys(Boc)-NHN=Pyv Resin	14	HAA1574	H-L-Cys(Bzl)-OH	19
FAA4370	Fmoc-L-Lys(Boc-Aoa)-OH	42	HAA6080	H-L-Cys(Bzl)-OMe*HCl	20
FAA9315	Fmoc-L-Lys(Boc-Cys(Trt))-OH	46	HAA6090	H-L-Cys(MBzl)-OH	20
FAA9320	Fmoc-L-Lys(Boc-Thz)-OH	18, 46	HAA9270	H-L-Cys(MDNPE)-OH	20
FAA4700	Fmoc-L-Lys(Mtt-Aoa)-OH	42	HAA6100	H-L-Cys(Mob)-OH	20
FAA6720	Fmoc-L-MeCys(Acm)-OH	27	HAA3510	H-L-Cys(Npys)-OH*HCl	20
FAA3340	Fmoc-L-MeCys(S-tBu)-OH	31	HAA9320	H-L-Cys(oNv)-OH	22
FAA3570	Fmoc-L-MeCys(Trt)-OH	31	HAA6140	H-L-Cys(SiBu)-OH	22
PYV1220	Fmoc-L-Met-NHN=Pyv Resin	14	HAA6150	H-L-Cys(tBu)-OH*HCl	22
PYV1230	Fmoc-L-Phe-NHN=Pyv Resin	15	HAA1560	H-L-Cys(Trt)-NH ₂	21

Product code	Product name	Page
HAA6160	H-L-Cys(Trt)-OH	21
HAA1995	H-L-Cys(Trt)-OtBu*HCl	21
HAA3530	H-L-Cys-NH ₂ *HCl	19
HAA3840	H-L-MeAla(4-Thz)-OH	17
HAA9230	H-L-Sec(MDNPE)*TFA	38
HAA9360	H-L-Sec(MDNPE)-OH	38
HAA9475	H-L-Sec(NPE)-OH*HCl	38
HAA9465	H-L-Sec(oNB)-OH*HCl	38
HAA9255	H-L-Sec(oNv)-OH*TFA	38
HAA1132	H-L-Thz-OH	17
HAA2860	N ₃ -Gly-Gly-Gly-Gly-Gly-OH	48
HAA2840	N ₃ -Gly-Gly-Gly-OH	48
HAA2850	N ₃ -Gly-Gly-OH*DCHA	49
HAA2810	N ₃ -L-Cys(Trt)-OH*CHA	22
RL-8625	N ₃ -PEG(4)-CINA	45
RL-4080	NTCB	44
SAA1110	Smoc-L-Cys(Trt)-OH	34
PAA2000	tBu-SS-Photo(Fmoc)-Gly-OH	35
LS-3405	TCEP*HCl	40
ZAA1161	Z-L-Cys(Bzl)-OH	32
ZAA1310	Z-L-Cys(Trt)-OH	32

Description of the Technology

Tools for the Formation of Thioesters

Building Blocks for Ligation

Enzyme-mediated Ligation

Examples

Index

Notes

Code of Conduct

As business activity of Iris Biotech GmbH impacts people's lives and health, it must be operated in ethical and correct manner and act with integrity and responsibility. To ensure high ethical standards and fair business practices, Iris Biotech GmbH applies an integrated policy known as its Code of Conduct.

In 2001 Iris Biotech GmbH was founded just at the beginning of the Biotech movement and the first remarkable breakthrough of biotech pharma products. Although the biotech field is rather young compared to other industries we believe on long-term business, a good partnership between our business partners and Iris Biotech GmbH and a good reputation. It is our duty as well as our responsibility to maintain and to extend this over the next generations – based on the principles of an honourable and prudent tradesman which based upon the concept of honourable entrepreneurship.

This Code of Conduct has been developed following the “Voluntary Guidelines for Manufacturers of Fine Chemical Intermediates and Active Ingredients” issued by AIME (Agrochemical & Intermediates Manufacturers in Europe) and the requirements of some of our business associates.

Iris Biotech GmbH commits to hold this Code of Conduct and to include and apply its principles in the management system and the company policies.

Ethics

Iris Biotech GmbH undertakes business in an ethical manner and acts with integrity. All corruption, extortion and embezzlement are prohibited. We do not pay or accept bribes or participate in other illegal inducements in business or government relationships. We conduct our business in compliance with all applicable anti-trust laws. Employees are encouraged to report concerns or illegal activities in the workplace, without threat of reprisal, intimidation or harassment.

Labour

Iris Biotech GmbH is committed to uphold the human rights of workers and to treat them with dignity and respect. Child labour, workplace harassment, discrimination, and harsh and inhumane treatment are prohibited. Iris Biotech GmbH respects the rights of the employees to associate freely, join or not join labour unions, seek representation and join workers' councils. Employees are paid and their working timetable is established according to applicable wage and labour laws. Employees are able to communicate openly with management regarding working conditions without threat of reprisal, intimidation or harassment.

General Policies

Contracts and Secrecy Agreements are binding and the confidential information received is only used for intended purposes. Clear management and organizational structures exist to provide efficient normal working and to address problems quickly. Know-how is protected and intellectual property is respected.

Health and Safety

Iris Biotech GmbH provides a safe and healthy working environment to the employees and protects them from overexposure to chemical and physical hazards. Products are produced, stored and shipped under the guidelines of the relevant chemical and safety legislation. Risks and emergency scenarios are identified and evaluated, and their possible impact is minimized by implementing emergency plans and written procedures. Safety information regarding hazardous materials is available to educate, train and protect workers from hazards. Preventive equipment and facilities maintenance is performed at suitable periods to reduce potential hazards. Employees are regularly trained in health and safety matters and are informed about product properties and risk classification when it is required.

Environment

Iris Biotech GmbH operates in an environmentally responsible and efficient manner, minimizing adverse impacts on the environment. Waste streams are managed to ensure a safe handling, movement, storage, recycling and reuse, before and after being generated. Systems to prevent and mitigate accidental spills and releases to the environment are in place. All required environmental permits and licenses are obtained and their operational and reporting requirements are complied with.

Production and Quality Management

A quality management system following the Good Distribution Practices (GDP rules) of Active Pharmaceutical Ingredients is established covering all the aspects of the worldwide distribution of products. Regular audits are performed to evaluate the efficiency and fulfilling of the quality system. Process controls to provide reproducible product quality are established. There are preventive maintenance procedures to ensure plant reliability and the lowest risk of failure. Staff is trained periodically about GMP and GDP rules. Procedures are established and installations are designed to avoid cross contamination. Batch and analytical records are kept for inspection and audit purposes for suitable periods according guidelines.

Research and Development

Research and development staff education is appropriate to their functional activity and they are trained to develop, optimize and scale-up the processes. Intellectual property is respected and know-how protected. Development of manufacturing processes reflects the principles of the Green Chemistry according to the American Chemical Society Green Chemistry Institute. Animal testing is not used unless alternatives are not scientifically valid or accepted by regulators. If animal testing is carried out, animals are treated so that pain and stress are minimized.

Terms and Conditions of Sales

All orders placed by a buyer are accepted and all contracts are made subject to the terms which shall prevail and be effective notwithstanding any variations or additions contained in any order or other document submitted by the buyer. No modification of these terms shall be binding upon Iris Biotech GmbH unless made in writing by an authorised representative of Iris Biotech GmbH.

Placing of Orders

Every order made by the buyer shall be deemed an offer by the buyer to purchase products from Iris Biotech GmbH and will not be binding on Iris Biotech GmbH until a duly authorised representative of Iris Biotech GmbH has accepted the offer made by the buyer. Iris Biotech GmbH may accept orders from commercial, educational or government organisations, but not from private individuals and Iris Biotech GmbH reserves the right to insist on a written order and/or references from the buyer before proceeding.

There is no minimum order value. At the time of acceptance of an order Iris Biotech GmbH will either arrange prompt despatch from stock or the manufacture/acquisition of material to satisfy the order. In the event of the latter Iris Biotech GmbH will indicate an estimated delivery date. In addition to all its other rights Iris Biotech GmbH reserves the right to refuse the subsequent cancellation of the order if Iris Biotech GmbH expects to deliver the product on or prior to the estimated delivery date. Time shall not be of the essence in respect of delivery of the products. If Iris Biotech GmbH is unable to deliver any products by reason of any circumstances beyond its reasonable control („Force Majeure“) then the period for delivery shall be extended by the time lost due to such Force Majeure. Details of Force Majeure will be forwarded by Iris Biotech GmbH to the buyer as soon as reasonably practicable.

Prices, Quotations and Payments

Prices are subject to change. For the avoidance of doubt, the price advised by Iris Biotech GmbH at the time of the buyer placing the order shall supersede any previous price indications. The buyer must contact the local office of Iris Biotech GmbH before ordering if further information is required. Unless otherwise agreed by the buyer and Iris Biotech GmbH, the price shall be for delivery ex-works. In the event that the buyer requires delivery of the products otherwise than ex-works the buyer should contact the local office of Iris Biotech GmbH in order to detail its requirements. Iris Biotech GmbH shall, at its discretion, arrange the buyer's delivery requirements including, without limitation, transit insurance, the mode of transit (Iris Biotech GmbH reserves the right to vary the mode of transit if any regulations or other relevant considerations so require) and any special packaging requirements (including cylinders). For the avoidance of doubt all costs of delivery and packaging in accordance with the buyer's requests over and above that of delivery in standard packaging ex-works shall be for the buyer's account unless otherwise agreed by both parties. Incoterms 2020 shall apply. Any tax, duty or charge imposed by governmental authority or otherwise and any other applicable taxes, duties or charges shall be for the buyer's account. Iris Biotech GmbH may, on request and where possible, provide quotations for multiple packs or bulk quantities, and non-listed items. Irrespective of the type of request or means of response all quotations must be accepted by the buyer without condition and in writing before an order will be accepted by Iris Biotech GmbH. Unless agreed in writing on different terms, quotations are valid for 30 days from the date thereof. Payment terms are net 30 days from invoice date unless otherwise agreed in writing. Iris Biotech GmbH reserves the right to request advance payment at its discretion. For overseas transactions the buyer shall pay all the banking charges of Iris Biotech GmbH. The buyer shall not

be entitled to withhold or set-off payment for the products for any reason whatsoever. Government/ Corporate Visa and MasterCard (and other such credit cards) may be accepted on approved accounts for payment of the products. Personal credit cards are not acceptable. Failure to comply with the terms of payment of Iris Biotech GmbH shall constitute default without reminder. In these circumstances Iris Biotech GmbH may (without prejudice to any other of its rights under these terms) charge interest to accrue on a daily basis at the rate of 2% per month from the date upon which payment falls due to the actual date of payment (such interest shall be paid monthly). If the buyer shall fail to fulfil the payment terms in respect of any invoice of Iris Biotech GmbH Iris Biotech GmbH may demand payment of all outstanding balances from the buyer whether due or not and/or cancel all outstanding orders and/or decline to make further deliveries or provision of services except upon receipt of cash or satisfactory securities. Until payment by the buyer in full of the price and any other monies due to Iris Biotech GmbH in respect of all other products or services supplied or agreed to be supplied by Iris Biotech GmbH to the buyer (including but without limitation any costs of delivery) the property in the products shall remain vested in Iris Biotech GmbH.

Shipping, Packaging and Returns

The buyer shall inspect goods immediately on receipt and inform Iris Biotech GmbH of any shortage or damage within five days. Quality problems must be notified within ten days of receipt. Goods must not be returned without prior written authorisation of Iris Biotech GmbH. Iris Biotech GmbH shall at its sole discretion replace the defective products (or parts thereof) free of charge or refund the price (or proportionate price) to buyer. Opened or damaged containers cannot be returned by the buyer without the written prior agreement of Iris Biotech GmbH. In the case of agreed damaged containers which cannot be so returned, the buyer assumes responsibility for the safe disposal of such containers in accordance with all applicable laws.

Product Quality, Specifications and Technical Information

Products are analysed in the Quality Control laboratories of Iris Biotech GmbH's production partners by methods and procedures which Iris Biotech GmbH considers appropriate. In the event of any dispute concerning reported discrepancies arising from the buyer's analytical results, determined by the buyer's own analytical procedures, Iris Biotech GmbH reserves the right to rely on the results of own analytical methods of Iris Biotech GmbH. Certificates of Analysis or Certificates of Conformity are available at the discretion of Iris Biotech GmbH for bulk orders but not normally for prepack orders. Iris Biotech GmbH reserves the right to make a charge for such certification. Specifications may change and reasonable variation from any value listed should not form the basis of a dispute. Any supply by Iris Biotech GmbH of bespoke or custom product for a buyer shall be to a specification agreed by both parties in writing. Technical information, provided orally, in writing, or by electronic means by or on behalf of Iris Biotech GmbH, including any descriptions, references, illustrations or diagrams in any catalogue or brochure, is provided for guidance purposes only and is subject to change.

Safety

All chemicals should be handled only by competent, suitably trained persons, familiar with laboratory procedures and potential chemical hazards. The burden of safe use of the products of Iris Biotech GmbH vests in the buyer. The buyer assumes all responsibility for warning his employees, and any persons who might reasonably be expected to come into contact with the products, of all risks to person and property in any way connected with the products and for instructing them in their safe handling and use. The buyer also assumes the responsibility for the safe disposal of all products in accordance with all applicable laws.

Uses, Warranties and Liabilities

All products of Iris Biotech GmbH are intended for laboratory research purposes and unless otherwise stated on product labels, in the catalogue and product information sheet of Iris Biotech GmbH or in other literature furnished to the buyer, are not to be used for any other purposes, including but not limited to use as or as components in drugs for human or animal use, medical devices, cosmetics, food additives, household chemicals, agricultural or horticultural products or pesticides. Iris Biotech GmbH offers no warranty regarding the fitness of any product for a particular purpose and shall not be responsible for any loss or damage whatsoever arising there from. No warranty or representation is given by Iris Biotech GmbH that the products do not infringe any letters patent, trademarks, registered designs or other industrial rights. The buyer further warrants to Iris Biotech GmbH that any use of the products in the United States of America shall not result in the products becoming adulterated or misbranded within the meaning of the Federal Food, Drug and Cosmetic Act (or such equivalent legislation in force in the buyer's jurisdiction) and shall not be materials which may not, under sections 404, 505 or 512 of the Act, be introduced into interstate commerce. The buyer acknowledges that, since the products of Iris Biotech GmbH are intended for research purposes, they may not be on the Toxic Substances Control Act 1976 („TSCA“) inventory. The buyer warrants that it shall ensure that the products are approved for use under the TSCA (or such other equivalent legislation in force in the buyer's jurisdiction), if applicable. The buyer shall be responsible for complying with any legislation or regulations governing the use of the products and their importation into the country of destination (for the avoidance of doubt to include, without limitation, the TSCA and all its amendments, all EINECS, ELINCS and NONS regulations). If any licence or consent of any government or other authority shall be required for the acquisition, carriage or use of the products by the buyer the buyer shall obtain the same at its own expense and if necessary produce evidence of the same to Iris Biotech GmbH on demand. Failure to do so shall not entitle the buyer to withhold or delay payment. Any additional expenses or charges incurred by Iris Biotech GmbH resulting from such failure shall be for the buyer's account. Save for death or personal injury caused by negligence of Iris Biotech GmbH, sole obligation of Iris Biotech GmbH and buyer's exclusive remedy with respect to the products proved to the satisfaction of Iris Biotech GmbH to be defective or products incorrectly supplied shall be to accept the return of said products to Iris Biotech GmbH for refund of the actual purchase price paid by the buyer (or proportionate part thereof), or replacement of the defective product (or part thereof) with alternative product. Iris Biotech GmbH shall have no liability to the buyer under or arising directly or indirectly out of or otherwise in connection with the supply of products by Iris Biotech GmbH to the buyer and/or their re-sale or use by the buyer or for any product, process or services of the buyer which in any way comprises the product in contract tort (including negligence or breach of statutory duty) or otherwise for pure economic loss, loss of profit, business, reputation, depletion of brand, contracts, revenues or anticipated savings or for any special indirect or consequential damage or loss of any nature except as may otherwise be expressly provided for in these terms. All implied warranties, terms and representations in respect of the products (whether implied by statute or otherwise) are excluded to the fullest extent permitted by law. The buyer shall indemnify Iris Biotech GmbH for and against any and all losses, damages and expenses, including legal fees and other costs of defending any action, that Iris Biotech GmbH may sustain or incur as a result of any act or omission by the buyer, its officers, agents or employees, its successors or assignees, its customers or all other third parties, whether direct or indirect, in connection with the use of any product. For the avoidance of doubt and in the event that Iris Biotech GmbH supplies bespoke or custom product to the buyer's design or specification, this indemnity shall extend to include any claim by a third party that the manufacture of the product for the buyer or the use of the product by the buyer infringes the intellectual property rights of any third party.

General

Iris Biotech GmbH shall be entitled to assign or sub-contract all or any of its rights and obligations hereunder. The buyer shall not be entitled to assign, transfer, sub-contract or otherwise delegate any of its rights or obligations hereunder. Any delay or forbearance by Iris Biotech GmbH in exercising any right or remedy under these terms shall not constitute a waiver of such right or remedy. If any provision of these terms is held by any competent authority to be invalid or unenforceable in whole or in part the validity of the other provisions of these terms and the remainder of the provision in question shall not be affected. These terms shall be governed by German Law and the German Courts shall have exclusive jurisdiction for the hearing of any dispute between the parties save in relation to enforcement where the jurisdiction of the German Courts shall be non-exclusive.



Get in Contact



Iris Biotech GmbH
Adalbert-Zoellner-Str. 1
95615 Marktredwitz
Germany

+49 (0) 9231 97121-0
+49 (0) 9231 97121-99
info@iris-biotech.de
www.iris-biotech.de

Distribution Partners

The list contains the current distributors of Iris Biotech in different regions of the world. The latest list of distribution partners and contact details is available at: www.iris-biotech.de/distribution-partner

China:

Chengdu Yoo Technology Co., Ltd.

Japan:

BizCom Japan, Inc.

Shigematsu & Co., Ltd

Cosmo Bio Co., Ltd.

USA & Canada:

Peptide Solutions LLC

India, Bangladesh, Oman, Sri Lanka, United Arab Emirates:

Sumit Biosciences Pvt Ltd.

Empowering Peptide Innovation