


SMOC-AMINO ACIDS

*Innovative Peptide Synthesis using
Water Instead of Organic Solvents*



The Clean Peptide Technology

Replace organic solvents with water during peptide synthesis by using 2,7-disulfo-9-fluorenylmethoxycarbonyl (Smoc) amino acids.

pages 1 – 4

Compatible with
water-swellaable resins.

pages 1 – 4

Pronounced fluorescence
allows real-time monitoring.

page 3

Purification by ion exchange
chromatography.

pages 1, 2



Version: IF13_3

Smoc-Amino Acids

Innovative Peptide Synthesis using Water Instead of Organic Solvents

For the production of synthetic peptides, every year tens of thousands of tons of organic solvents are required in the chemical, cosmetic and pharmaceutical industries. According to the European Chemicals Directive REACH, these solvents are classified as substances of very high concern and their use is associated with significant risks for health and the environment. Following modern regulatory rules, “undesirable” solvents such as DCM, DMF, NMP and THF, which are frequently recommended and employed during SPPS, should be replaced. At the same time, the enormous consumption of solvents and reagents leads to high production costs.

Thus, despite being considered as advanced and efficient technique for peptide production, solid-phase peptide synthesis (SPPS) is associated with severe drawbacks. Therefore, ongoing attempts aim at developing alternative approaches using solvents, which are permitting to reduce the risks for environment and human health.

The Clean Peptide Technology Using Smoc-Amino Acids – eco-friendly peptide manufacturing

Benefits:


- 2,7-disulfo-9-fluorenylmethoxycarbonyl (Smoc) amino acids are water-compatible
→ replacement of organic solvents during peptide synthesis
- Certain side-chains can remain unprotected
→ more atom-efficient
- Smoc-amino acids show a pronounced fluorescence
→ real-time monitoring during peptide synthesis ($\lambda_{\text{ex}} = 280 \text{ nm}$; $\lambda_{\text{em}} = 340 \text{ nm}$)
- As a solid support, a water-swelling resin can be used



You need more details? Watch the recording of our online workshop about the Smoc technology and its use for peptide synthesis in water.



Methodology:



Cleavage of the Smoc-group can be achieved using aqueous piperidine, piperazine, sodium hydroxide, ethanolamine, and ammonia to liberate the N-terminus within 5-15 minutes at ambient temperature together with the formation of the respective disulfonated dibenzofulvene and the products of water or base addition.

Recommended Coupling Protocols:

	Method 1) <i>in situ</i> Oxyma coupling	Method 2) NHS with pre-activation
Preparation & Coupling	3.0 eq. Smoc amino acid, 5.5 eq. EDC, 3.5 eq. Oxyma, 3.0 eq. NaHCO ₃ in 30% isopropanol or MeCN in water 45 min or 2x25 min for double coupling	NHS ester formation: 3.0 eq. Smoc amino acid, 5.5 eq. EDC, 3.5 eq. NHS in water (pH range 5-6) 25 min For the coupling, add the prepared NHS mixture to the amine, adjust the pH to 8.0 and readjust over time for 15 min.
Wash	2x with water	2x with water
Deprotection	1M NaOH (5 min + 10 min) If ester side-chains are present in the se- quence use another base, e.g., piperazine or 4-methylpiperidine.	1M NaOH (5 min + 10 min) If ester side-chains are present in the se- quence use another base, e.g., piperazine or 4-methylpiperidine.
Wash	2x with water	2x with water


Note:

- Rink amide resin needs to be dry before cleavage from the solid support, otherwise side reactions occur.
- Recommended resins: TentaGel® or other resins suitable for polar solvents, e.g., Li resin



Comparison Boc – Fmoc – Smoc







<i>Boc SPPS</i>	<i>Fmoc SPPS</i>	<i>Smoc SPPS</i>
Boc-Arg(Tos)-OH	Fmoc-Arg(Pbf)-OH	Smoc-Arg-OH
Boc-Asn(Xan)-OH	Fmoc-Asn(Trt)-OH	Smoc-Asn-OH
Boc-Asp(OBzl)-OH	Fmoc-Asp(OtBu)-OH	Smoc-Asp(OtBu)-OH
Boc-Cys(Acm)-OH	Fmoc-Cys(Trt)-OH	Smoc-Cys(Trt)-OH
Boc-Gln(Xan)-OH	Fmoc-Gln(Trt)-OH	Smoc-Gln-OH
Boc-Glu(OBzl)-OH	Fmoc-Glu(OtBu)-OH	Smoc-Glu(OtBu)-OH
Boc-His(Dnp)-OH	Fmoc-His(Trt)-OH	Smoc-His-OH*
Boc-Lys(Cbz)-OH	Fmoc-Lys(Boc)-OH	Smoc-Lys(Boc)-OH
Boc-Ser(Bzl)-OH	Fmoc-Ser(tBu)-OH	Smoc-Ser(tBu)-OH
Boc-Thr(Bzl)-OH	Fmoc-Thr(tBu)-OH	Smoc-Thr(tBu)-OH
Boc-Trp(For)-OH	Fmoc-Trp(Boc)-OH	Smoc-Trp-OH*
Boc-Tyr(Bzl)-OH	Fmoc-Tyr(tBu)-OH	Smoc-Tyr-OH*
* also available with side-chain protection		

Exemplary Fluorescence Monitoring:



In comparison to Fmoc SPPS, no increased racemization levels were observed during Smoc SPPS. Regarding the prevention of aspartimide formation during Smoc SPPS in water, reduced temperature is recommended for the deprotection of sequences prone to aspartimide formation.

		Product details
SAA1010	Smoc-L-Ala-OH	
(((2,7-disulfo-9H-fluoren-9-yl)methoxy)carbonyl)-L-alanine potassium salt		
CAS-No.	2442552-59-0	
Formula	C ₁₈ H ₁₅ K ₂ NO ₁₀ S ₂	
Mol. weight	547,63 g/mol	
		
SAA1230	Smoc-beta-Ala-OH	
(((2,7-disulfo-9H-fluoren-9-yl)methoxy)carbonyl)-beta-alanine potassium salt		
CAS-No.	2337407-19-7	
Formula	C ₁₈ H ₁₅ K ₂ NO ₁₀ S ₂	
Mol. weight	547,63 g/mol	
		

		Product details
SAA1050	Smoc-L-Arg-OH (((2,7-disulfo-9H-fluoren-9-yl)methoxy)carbonyl)-L-arginine potassium salt CAS-No. 2337407-38-0 Formula $C_{21}H_{22}K_2N_4O_{10}S_2$ Mol. weight 632,74 g/mol	 
SAA1080	Smoc-L-Asn-OH (((2,7-disulfo-9H-fluoren-9-yl)methoxy)carbonyl)-L-asparagine potassium salt CAS-No. 2337407-22-2 Formula $C_{19}H_{16}K_2N_2O_{11}S_2$ Mol. weight 590,66 g/mol	 
SAA1130	Smoc-L-Asp(OtBu)-OH (((2,7-disulfo-9H-fluoren-9-yl)methoxy)carbonyl)-L-aspartic-acid-beta-t-butyl-ester potassium salt CAS-No. 2337407-41-5 Formula $C_{23}H_{23}K_2NO_{12}S_2$ Mol. weight 647,75 g/mol	 
SAA1110	Smoc-L-Cys(Trt)-OH N-(((2,7-disulfo-9H-fluoren-9-yl)methoxy)carbonyl)-S-trityl-L-cysteine potassium salt CAS-No. 2442552-68-1 Formula $C_{37}H_{29}K_2NO_{10}S_3$ Mol. weight 822,01 g/mol	 
SAA1120	Smoc-L-Glu(OtBu)-OH (((2,7-disulfo-9H-fluoren-9-yl)methoxy)carbonyl)-L-glutamic-acid-gamma-t-butyl-ester potassium salt CAS-No. 2442552-71-6 Formula $C_{24}H_{25}K_2NO_{12}S_2$ Mol. weight 661,78 g/mol	 
SAA1070	Smoc-L-Gln-OH (((2,7-disulfo-9H-fluoren-9-yl)methoxy)carbonyl)-L-glutamine potassium salt CAS-No. 2337407-39-1 Formula $C_{20}H_{18}K_2N_2O_{11}S_2$ Mol. weight 604,68 g/mol	 


SAA1000 Smoc-Gly-OH

(((2,7-disulfo-9H-fluoren-9-yl)methoxy)carbonyl)glycine potassium salt

CAS-No. 2337407-26-6

 Formula $C_{17}H_{13}K_2NO_{10}S_2$

Mol. weight 533,60 g/mol



SAA1140 Smoc-L-His-OH

(((2,7-disulfo-9H-fluoren-9-yl)methoxy)carbonyl)-L-histidine potassium salt

CAS-No. 2442552-74-9 (net)

 Formula $C_{21}H_{17}K_2N_3O_{10}S_2$

Mol. weight 613,69 g/mol



SAA1220 Smoc-L-His(Trt)-OH

N-alpha-(((2,7-disulfo-9H-fluoren-9-yl)methoxy)carbonyl)-N-tau-trityl-L-histidine potassium salt

CAS-No. 2442552-76-1 (net)

 Formula $C_{40}H_{31}K_2N_3O_{10}S_2$

Mol. weight 856,02 g/mol



SAA1030 Smoc-L-Ile-OH

(((2,7-disulfo-9H-fluoren-9-yl)methoxy)carbonyl)-L-isoleucine potassium salt

CAS-No. 2337407-24-4

 Formula $C_{21}H_{21}K_2NO_{10}S_2$

Mol. weight 589,71 g/mol



SAA1040 Smoc-L-Leu-OH

(((2,7-disulfo-9H-fluoren-9-yl)methoxy)carbonyl)-L-leucine potassium salt

CAS-No. 2337407-36-8

 Formula $C_{21}H_{21}K_2NO_{10}S_2$

Mol. weight 589,71 g/mol



SAA1190 Smoc-L-Lys(Boc)-OH







N6-(tert-butoxycarbonyl)-N2-(((2,7-disulfo-9H-fluoren-9-yl)methoxy)carbonyl)-L-lysine potassium salt






CAS-No. 2442552-82-9

 Formula $C_{26}H_{30}K_2N_2O_{12}S_2$

Mol. weight 704,84 g/mol



		Product details
SAA1100	Smoc-L-Met-OH (((2,7-disulfo-9H-fluoren-9-yl)methoxy)carbonyl)-L-methionine potassium salt CAS-No. 2442552-84-1 Formula $C_{20}H_{19}K_2NO_{10}S_3$ Mol. weight 607,75 g/mol	 
SAA1060	Smoc-L-Phe-OH (((2,7-disulfo-9H-fluoren-9-yl)methoxy)carbonyl)-L-phenylalanine potassium salt CAS-No. 2442552-86-3 Formula $C_{24}H_{19}K_2NO_{10}S_2$ Mol. weight 623,73 g/mol	 
SAA1150	Smoc-L-Pro-OH (((2,7-disulfo-9H-fluoren-9-yl)methoxy)carbonyl)-L-proline potassium salt CAS-No. 2337407-20-0 Formula $C_{20}H_{17}K_2NO_{10}S_2$ Mol. weight 573,67 g/mol	 
SAA1170	Smoc-L-Ser(tBu)-OH O-(tert-butyl)-N-(((2,7-disulfo-9H-fluoren-9-yl)methoxy)carbonyl)-L-serine potassium salt CAS-No. 2337407-37-9 (net) Formula $C_{22}H_{23}K_2NO_{11}S_2$ Mol. weight 619,74 g/mol	 
SAA1160	Smoc-L-Thr(tBu)-OH O-(tert-butyl)-N-(((2,7-disulfo-9H-fluoren-9-yl)methoxy)carbonyl)-L-threonine potassium salt CAS-No. 2442552-94-3 Formula $C_{23}H_{25}K_2NO_{11}S_2$ Mol. weight 633,77 g/mol	 
SAA1180	Smoc-L-Trp-OH (((2,7-disulfo-9H-fluoren-9-yl)methoxy)carbonyl)-L-tryptophan potassium salt CAS-No. 2442552-96-5 Formula $C_{26}H_{20}K_2N_2O_{10}S_2$ Mol. weight 662,77 g/mol	 

		Product details
SAA1210	Smoc-L-Trp(Boc)-OH	
1-(<i>tert</i> -butoxycarbonyl)-N-α-(((2,7-disulfo-9H-fluoren-9-yl)methoxy)carbonyl)-L-tryptophan potassium salt		
CAS-No.	2442552-98-7	 
Formula	C ₃₁ H ₂₈ K ₂ N ₂ O ₁₂ S ₂	
Mol. weight	762,88 g/mol	
SAA1090	Smoc-L-Tyr-OH	
(((2,7-disulfo-9H-fluoren-9-yl)methoxy)carbonyl)-L-tyrosine potassium salt		
CAS-No.	2337407-33-5	 
Formula	C ₂₆ H ₁₉ K ₂ NO ₁₁ S ₂	
Mol. weight	639,73 g/mol	
SAA1200	Smoc-L-Tyr(OtBu)-OH	
O-(<i>tert</i> -butyl)-N-(((2,7-disulfo-9H-fluoren-9-yl)methoxy)carbonyl)-L-tyrosine potassium salt		
CAS-No.	2442553-00-4	 
Formula	C ₂₈ H ₂₇ K ₂ NO ₁₁ S ₂	
Mol. weight	695,84 g/mol	
SAA1020	Smoc-L-Val-OH	
(((2,7-disulfo-9H-fluoren-9-yl)methoxy)carbonyl)-L-valine potassium salt		
CAS-No.	2442553-02-6	 
Formula	C ₂₀ H ₁₉ K ₂ NO ₁₀ S ₂	
Mol. weight	575,69 g/mol	
SAA1240	Smoc-Cl	
9-(((chlorocarbonyl)oxy)methyl)-9H-fluorene-2,7-disulfonic acid		
CAS-No.	1899144-09-2	 
Formula	C ₁₅ H ₁₁ ClO ₈ S ₂	
Mol. weight	418,83 g/mol	

References:

- Method for peptide synthesis and apparatus for carrying out a method for solid phase synthesis of peptides; S. Knauer, T. M. L. Roese, O. Avrutina, H. Kolmar, C. Uth; 2016, WO 2016 050764.
- Sustainable Peptide Synthesis Enabled by a Transient Protecting Group; S. Knauer, N. Koch, C. Uth, R. Meusinger, O. Avrutina, H. Kolmar; **Angew. Chem. Int. Ed.** 2020; **59(31)**: 12984-12990. <https://doi.org/10.1002/anie.202003676>
- Improved method for preparing peptides; S. Knauer; WO 2019 101939.
- Method for preparing peptides; S. Knauer; WO 2019 101940.
- Novel amino-Li resin for water-based solid-phase peptide synthesis; C. Uth, S. Englert, O. Avrutina, H. Kolmar, S. Knauer; **J. Pept. Sci.** 2023; **29(12)**: e3527. <https://doi.org/10.1002/psc.3527>



Looking for the right resin?

Download our booklet!



Empowering Peptide Innovation