

CHO | 360-HCP ELISA Kit

BioGenes Berlin

**Superior
Detection for
CHO HCPs**



Ready-to-use **CHO** | 360-HCP ELISA Kits reliable HCP monitoring during drug development and manufacturing



Features

- Different antibody sets by four antigen preparations with distinct HCP pattern
- High specificity and sensitivity - qualified according to ICH guidelines
- Superior performance compared to other commercially available products

Diverse Antigen Preparations and Differentiated Immunization

Our enhanced generic 360-HCP ELISA approach provides four kits, generated by parallel immunization of goats and rabbits using total or fractionated HCP obtained from mock fermentation of CHO-K1 and CHO-S.

**Identify the best performing kit for your process
with the Starter Set!**

Type
A

Type
B

Type
C

Type
D