

ZetaSep Protein A Resolute ZetaSep Affinity FPLC column with 1 mL bed volume	
Product Code:	ZS-1425-M001.0
Quality Parameter	Specification
Packing evaluation asymmetry: (no double peak, no shoulder)	$0.7 \leq A_s \leq 1.9$
Asymmetry LOT homogeneity:	$CV \leq 8 \%$
Notice:	The use of this product is strictly limited to trained personnel for professional manufacturing, laboratory, or research purposes. Final Fitness-For-Use must be determined by and is the sole responsibility of the end-user.