
Nitrocellulose Blotting Membrane

Type: Pre-cut sheets

Product code : 107A-109A,107B-109B

Features:

- Available in pre-cut sheets and rolls with two different pore sizes
- Uniform pore structure
- Compatible with colorimetric and radiolabeling detection methods
- Molecular fixation by baking, drying, alkali treatment and UV crosslinking is possible
- High strength and durability, preventing distortion or contamination in multiple reprobings
- Capable of binding wide range of fragment sizes
- Inherently charged

Applications:

- Northern and Southern blotting
- Electrotransfer
- Capillary blot

Specifications:

Pore sizes (μ)	Pre-cut Sheets Dimensions (cm x cm)		
	15 x 15	20 x 20	30 x 30
0.22	SF107A	SF108A	SF109A
0.45	SF107B	SF108B	SF109B

Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt. Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

Revision No.: 01 / 2020