



### Prima-Trio – LA1006 : Technical Specification

<b>Specifications</b>	
Manufacturer	<b>Himedia</b>
Product Name	<b>Prima - Trio</b>
Code No	<b>LA1006</b>
Sample Capacity	3x(32x0.2ml)
Peltier Quantity	6
Block Material	Aluminum Alloy
Temperature Range(°C)	0-100°C
Heating Rate (°C/sec)	4.5
Cooling Rate (°C/sec)	4
Control Accuracy(°C)	≤ ±0.3
Uniformity (°C)	≤ ±0.4
Control Mode	Block,Tube
Hot-lid Temperature (°C)	30-110°C
Gradient Range (°C)	1-30°C
Minimum Gradient Temp. (°C)	30°C
Maximum Gradient Temp. (°C)	100°C
Display Screen	8 inch, 800 X 600 Pixels, TFT, TOUCH
Programs Saving	10000 files +(USB Flash)
Power-off Recovery	Yes
Temperature Time Increase/Decrease Progressively	Yes
Pause Function	Yes
Tm Calculator	Yes
Edit File While Running	Yes
Support for USB Mouse	Yes
Support for Touch Screen	Yes
Connection	USB,Ethernet
Network Control	Yes
Auto-open Lid	No

<b>Specifications</b>	
PC Control Software	Yes
Input Voltage	AC85-264V, 47-63Hz
Input Power	600W
Footprint (DxWxH) cm	39 X 27X 25.5
Weight (Kg)	9
Certification	CE Certified
Warranty Period	2 years



HiMediaLaboratories™

**HiMedia Laboratories Pvt. Limited**

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est.  
 LBS Marg, Mumbai - 400 086, India Tel: 00-91-22-6147 1919, 2500 3747, 2500 0970  
 Fax: 6147 1920, 2500 5764, 2500 2468, 2500 2286 www.himedialabs.com  
 Email : atc@himedialabs.com

